

## HEARTSTART MRx MONITOR/DEFIBRILLATOR



### PRODUCT DESCRIPTION

No compromises. That's how we created the HeartStart MRx, a combination transport monitor, cardiograph, manual defibrillator, and AED. HeartStart MRx unites our industry-leading patient monitoring capabilities with our patented biphasic defibrillation therapy in one comprehensive, easy-to-use, rugged, and lightweight device, optimized to meet a variety of user needs.

With carefully arranged controls, connections, and menus, HeartStart MRx efficiently provides comprehensive diagnostic measurements without compromising critical, lifesaving therapy. Monitoring and resuscitation functions are clearly delineated, with monitoring located on one side of the device and defibrillation, pacing, and cardioversion on the other. Patient status and care data --- including numeric vital signs and up to four waves, as well as text prompts, alarms, an event timer, event markers, and battery status indicators --- are logically organized on the device's large, backlit, color display.

It is our mission to produce the finest quality medical equipment, enabling you to deliver effective care and ultimately save lives. We strive to meet the evolving needs of the healthcare community with outstanding technology and innovation, while always focusing on what's most important - positive patient outcomes. So, whether it's the long-life batteries, automated operational checks, or a data management solution, HeartStart MRx has the capabilities you need and the performance you demand - no compromises.

### Features

#### Standard Features

- ▶ ST/AR Basic algorithm for arrhythmia detection
- ▶ ECG monitoring through monitoring electrodes and defibrillation pads
- ▶ Synchronized cardioversion
- ▶ Adjustable ECG size and autogain
- ▶ Manual and AED operation
- ▶ SMART Biphasic waveform for defibrillation therapy
- ▶ Large 4-wave color display
- ▶ Strip chart printer
- ▶ Individual, adjustable volume of QRS beeper, voice prompts, and alerts
- ▶ Event summary
- ▶ Configuration mode
- ▶ Service mode
- ▶ Operational checks
- ▶ Automated self-tests with "ready-for-use" indicator
- ▶ Lithium ion battery with fuel gauge
- ▶ ECG out capability

#### Optional Features

- ▶ SpO<sub>2</sub> with Fourier Artifact Suppression Technology (FAST)
- ▶ Noninvasive Blood Pressure
- ▶ Microstream™ EtCO<sub>2</sub>
- ▶ Noninvasive Pacing
- ▶ 12-Lead ECG with Philips interpretation algorithm
- ▶ 75 mm Printer

#### Standard Accessories

- ▶ AC power module
- ▶ Lithium ion battery with fuel gauge
- ▶ Hands-free multifunction electrode cable
- ▶ 3-Lead ECG cable
- ▶ Disposable monitoring electrodes

- ▶ Strip chart recorder paper
- ▶ Defibrillator test load
- ▶ Documentation CD containing Instructions for Use (User's guide), User training workbook, and Application notes
- ▶ Quick reference cards

#### Package

**M3535A** HeartStart MRx Monitor/Defibrillator

#### Ordering Option Information

|            |                                                                      |
|------------|----------------------------------------------------------------------|
| <b>A01</b> | SpO <sub>2</sub>                                                     |
| <b>A02</b> | SpO <sub>2</sub> and NBP                                             |
| <b>A03</b> | SpO <sub>2</sub> , NBP, and EtCO <sub>2</sub>                        |
| <b>A04</b> | EtCO <sub>2</sub>                                                    |
| <b>B01</b> | External Pacing                                                      |
| <b>B02</b> | 12-Lead ECG Acquisition                                              |
| <b>B04</b> | Wider Printer                                                        |
| <b>C01</b> | Standard External Paddles                                            |
| <b>C02</b> | Water Resistant External Paddles                                     |
| <b>C03</b> | Data Card                                                            |
| <b>C04</b> | Carrying Case                                                        |
| <b>C05</b> | Additional Battery                                                   |
| <b>C07</b> | Barrel style Pad Cable (Replacement for Standard Pad Cable)          |
| <b>C15</b> | 5-Lead ECG Cable                                                     |
| <b>LP1</b> | Instructions for Use                                                 |
| <b>LP2</b> | User Training Video (English only)                                   |
| <b>LP3</b> | User Training DVD (English only)                                     |
| <b>SM1</b> | Service Manual (English only)                                        |
| <b>SM3</b> | Service Training Video (English only)                                |
| <b>W01</b> | One-Year On-Site Warranty                                            |
| <b>W24</b> | Five-Year Bench Warranty with Loaner (U.S. & Canada only)            |
| <b>WA1</b> | Five-Year Biomed Warranty (U.S., Canada & Australia only)            |
| <b>W08</b> | Three-Year Bench Warranty with Loaner (Australia & New Zealand only) |

## Upgrades/Supplies/Accessories

### Upgrades

|               |                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------|
| <b>M3530A</b> | SpO <sub>2</sub>                                                                                           |
| <b>M3531A</b> | NBP                                                                                                        |
| <b>M3532A</b> | EtCO <sub>2</sub>                                                                                          |
| <b>M3533A</b> | Pacing                                                                                                     |
| <b>M3534A</b> | 12-Lead ECG<br>Option B02 - Acquisition<br>Option B03 - Transmission<br>Option B04 - Wider Printer         |
| <b>M4760A</b> | Handle and Cap Plate (for Pads)                                                                            |
| <b>M5527A</b> | External Paddles with Paddle Tray<br>Option C01 - Standard Paddles<br>Option C02 - Water Resistant Paddles |

### External Paddles

|               |                                  |
|---------------|----------------------------------|
| <b>M3542A</b> | Standard External Paddles        |
| <b>M3543A</b> | Water Resistant External Paddles |
| <b>M4745A</b> | Sterilizable External Paddles    |

### Internal Paddles

|                     |                                |
|---------------------|--------------------------------|
| <b>M1741A</b>       | 7.5 cm Switchless              |
| <b>M1742A</b>       | 6.0 cm Switchless              |
| <b>M1743A</b>       | 4.5 cm Switchless              |
| <b>M1744A</b>       | 2.8 cm Switchless              |
| <b>M4741A</b>       | 7.5 cm Switched                |
| <b>M4742A</b>       | 6.0 cm Switched                |
| <b>M4743A</b>       | 4.5 cm Switched                |
| <b>M4744A</b>       | 2.8 cm Switched                |
| <b>989803127121</b> | Large Disposable Switched      |
| <b>989803127131</b> | Medium Disposable Switched     |
| <b>989803127141</b> | Small Disposable Switched      |
| <b>989803127151</b> | Large Disposable Switchless    |
| <b>989803127161</b> | Medium Disposable Switchless   |
| <b>989803127171</b> | Small Disposable Switchless    |
| <b>M4740A</b>       | Internal Paddles Adapter Cable |

## Multifunction Electrode Pads

### White barrel connector

|               |                       |
|---------------|-----------------------|
| <b>M3501A</b> | Defib Adult, AAMI     |
| <b>M3502A</b> | Defib Adult, IEC      |
| <b>M3503A</b> | Defib Pediatric, IEC  |
| <b>M3504A</b> | Defib Pediatric, AAMI |

### Grey plug connector

|               |                                         |
|---------------|-----------------------------------------|
| <b>M3713A</b> | Adult Plus                              |
| <b>M3716A</b> | Adult Radiolucent                       |
| <b>M3717A</b> | Pediatric Plus                          |
| <b>M3718A</b> | Adult Radiotransparent/Reduced Skin     |
| <b>M3719A</b> | Pediatric Radiotransparent/Reduced Skin |
| <b>DP2</b>    | Self-adhesive adult electrode pads      |
| <b>DP6</b>    | Self-adhesive adult electrode pads      |

### Pads Cable

|                 |                                           |
|-----------------|-------------------------------------------|
| <b>M3507A</b>   | Defib Hands-free Pads Cable, barrel style |
| <b>M3508A</b>   | Defib Hands-free, plug style              |
| <b>05-10200</b> | Pads Adapter (use with M3507A)            |

## ECG Monitoring Electrodes

|               |                                                   |
|---------------|---------------------------------------------------|
| <b>M2202A</b> | High-Tack Foam; 5 electrodes/pack (60 packs/case) |
|---------------|---------------------------------------------------|

## ECG Cables

### 12-Lead Cable Set

|               |                                                                                          |
|---------------|------------------------------------------------------------------------------------------|
| <b>M3525A</b> | 10-lead ECG Trunk Cable, 12-pin Connector (for 3-Lead and 12-Lead use)                   |
| <b>M3526A</b> | 3-lead ECG Set and Plug with Snap (AAMI)                                                 |
| <b>M3527A</b> | Add 7-lead ECG Set for 12-Lead use (AAMI)                                                |
| <b>M3528A</b> | 3-lead ECG Set and Plug with Snap (IEC)                                                  |
| <b>M3529A</b> | Add 7-lead ECG Set for 12-Lead use (IEC)                                                 |
| <b>M5530A</b> | Combiner Plug for 3-wire Lead Set for use with M3526A/M3528A                             |
| <b>M1949A</b> | 10-lead ECG Patient Trunk Cable, 12-pin ECG Input Connector (for 5-Lead and 12-Lead use) |
| <b>M1968A</b> | 10-electrode Cable Set, Extremities, Grabber (use with M1976A) (AAMI)                    |
| <b>M1976A</b> | 10-electrode Cable Set, Chest, Grabber (use with M1968A) (AAMI)                          |
| <b>M1971A</b> | 10-electrode Cable Set, Extremities, Grabber (use with M1978A) (IEC)                     |
| <b>M1978A</b> | 10-electrode Cable Set, Chest, Grabber (use with M1971A) (IEC)                           |

### ECG Cables

#### 3-Lead Cable Set

|               |                                       |
|---------------|---------------------------------------|
| <b>M1500A</b> | 3-lead ECG Trunk Cable (AAMI)         |
| <b>M1605A</b> | 3-lead ECG Lead Set with Snaps (AAMI) |
| <b>M1510A</b> | 3-lead ECG Trunk Cable (IEC)          |
| <b>M1615A</b> | 3-lead ECG Lead Set with Snaps (IEC)  |

#### 5-Lead Cable Set

|               |                                       |
|---------------|---------------------------------------|
| <b>M1520A</b> | 5-lead ECG Trunk Cable (AAMI)         |
| <b>M1625A</b> | 5-lead ECG Lead Set with Snaps (AAMI) |
| <b>M1530A</b> | 5-lead ECG Trunk Cable (IEC)          |
| <b>M1635A</b> | 5-lead ECG Lead Set with Snaps (IEC)  |

### Sync Cables

|               |                           |
|---------------|---------------------------|
| <b>M1783A</b> | Sync Cable 8 ft. (2.4 m)  |
| <b>M5526A</b> | Sync Cable 25 ft. (7.6 m) |

### SpO<sub>2</sub> Sensors/Cables

|               |                                                                        |
|---------------|------------------------------------------------------------------------|
| <b>M1191A</b> | Reusable SpO <sub>2</sub> Sensor - Adult Finger 6.5 ft. (2 m)          |
| <b>M1192A</b> | Reusable SpO <sub>2</sub> Sensor - Pediatric/Small Adult Finger        |
| <b>M1193A</b> | Reusable SpO <sub>2</sub> Sensor - Adult Thumb                         |
| <b>M1194A</b> | Reusable SpO <sub>2</sub> Sensor - Adult Ear Clip                      |
| <b>M1195A</b> | Reusable SpO <sub>2</sub> Sensor - Infant                              |
| <b>M1941A</b> | SpO <sub>2</sub> Extension Cable, 2 m (6.5 ft.)                        |
| <b>M1943A</b> | Nellcor SpO <sub>2</sub> Sensor Adapter Cable, 1 m (use with M1903/4B) |

### NBP

#### Interconnect Cables

|               |                              |
|---------------|------------------------------|
| <b>M1598B</b> | Adult Pressure 5 ft. (1.5 m) |
| <b>M1599B</b> | Adult Pressure 10 ft. (3 m)  |

#### Reusable Blood Pressure Cuffs

|               |                                                                               |
|---------------|-------------------------------------------------------------------------------|
| <b>40400A</b> | Reusable NBP Cuff Kit, 3 sizes (pediatric, adult, large adult)                |
| <b>40400B</b> | Reusable NBP Cuff Kit, 5 sizes (infant, pediatric, adult, large adult, thigh) |
| <b>40401A</b> | Traditional Reusable NBP Cuff - Infant                                        |
| <b>40401B</b> | Traditional Reusable NBP Cuff - Pediatric                                     |

### NBP

|               |                                               |
|---------------|-----------------------------------------------|
| <b>40401C</b> | Traditional Reusable NBP Cuff - Adult         |
| <b>40401D</b> | Traditional Reusable NBP Cuff - Large Adult   |
| <b>40401E</b> | Traditional Reusable NBP Cuff - Thigh         |
| <b>M4552A</b> | Antimicrobial Reusable NBP Cuff - Infant      |
| <b>M4553A</b> | Antimicrobial Reusable NBP Cuff - Pediatric   |
| <b>M4554A</b> | Antimicrobial Reusable NBP Cuff - Small Adult |
| <b>M4555A</b> | Antimicrobial Reusable NBP Cuff - Adult       |
| <b>M4557A</b> | Antimicrobial Reusable NBP Cuff - Large Adult |
| <b>M4559A</b> | Antimicrobial Reusable NBP Cuff - Thigh       |
| <b>M1572A</b> | Multi-Patient Comfort Cuffs - Pediatric       |
| <b>M1573A</b> | Multi-Patient Comfort Cuffs - Small Adult     |
| <b>M1574A</b> | Multi-Patient Comfort Cuffs - Adult           |
| <b>M1575A</b> | Multi-Patient Comfort Cuffs - Large Adult     |

### Disposable Blood Pressure Cuffs

|               |                                                          |
|---------------|----------------------------------------------------------|
| <b>M4572A</b> | Soft Single-Patient Disposable Cuff - Infant             |
| <b>M4573A</b> | Soft Single-Patient Disposable Cuff - Pediatric          |
| <b>M4574A</b> | Soft Single-Patient Disposable Cuff - Small Adult        |
| <b>M4575A</b> | Soft Single-Patient Disposable Cuff - Adult              |
| <b>M4576A</b> | Soft Single-Patient Disposable Cuff - Adult X-Long       |
| <b>M4577A</b> | Soft Single-Patient Disposable Cuff - Large Adult        |
| <b>M4578A</b> | Soft Single-Patient Disposable Cuff - Large Adult X-Long |
| <b>M4579A</b> | Soft Single-Patient Disposable Cuff - Thigh              |

### EtCO<sub>2</sub>

#### Intubated Circuits

|               |                                                                  |
|---------------|------------------------------------------------------------------|
| <b>M1920A</b> | Filter Line Set - Adult/Pediatric (25 sets/case)                 |
| <b>M1921A</b> | Filter H Set - Humidified Adult/Pediatric (25 sets/case)         |
| <b>M1923A</b> | Filter H Set - Humidified Infant/Neonatal (yellow, 25 sets/case) |

### Non-Intubated Dual Purpose Circuits (CO<sub>2</sub> + O<sub>2</sub>)

|               |                                |
|---------------|--------------------------------|
| <b>M2520A</b> | Smart CapnoLine - Pediatric    |
| <b>M2521A</b> | Smart CapnoLine - Intermediate |
| <b>M2522A</b> | Smart CapnoLine - Adult        |

### EtCO<sub>2</sub>

#### Non-Intubated Single Purpose Circuits (CO<sub>2</sub>)

|               |                                |
|---------------|--------------------------------|
| <b>M2524A</b> | Smart CapnoLine - Pediatric    |
| <b>M2525A</b> | Smart CapnoLine - Intermediate |
| <b>M2526A</b> | Smart CapnoLine - Adult        |

### Power

|                     |                                                                             |
|---------------------|-----------------------------------------------------------------------------|
| <b>M3538A</b>       | Lithium Ion Battery with fuel gauge                                         |
| <b>M3539A</b>       | AC Power Module                                                             |
| <b>989803135301</b> | 2-Bay Battery Support System for Lithium Ion Batteries                      |
| <b>989803135331</b> | 4-Bay Battery Support System for Lithium Ion Batteries                      |
| <b>989803135341</b> | 4-Bay Battery Support System for Sealed Lead Acid and Lithium Ion Batteries |

### Paper

|                     |                                              |
|---------------------|----------------------------------------------|
| <b>40457C</b>       | 50 mm Chemical Thermal, Gray Grid (10 rolls) |
| <b>40457D</b>       | 50 mm Chemical Thermal, Gray Grid (80 rolls) |
| <b>989803138171</b> | 75 mm Chemical Thermal, Gray Grid (10 rolls) |
| <b>989803138181</b> | 75 mm Chemical Thermal, Gray Grid (80 rolls) |

### Color Handle

|               |                            |
|---------------|----------------------------|
| <b>M5521A</b> | Green and color label set  |
| <b>M5522A</b> | Blue and color label set   |
| <b>M5523A</b> | Yellow and color label set |
| <b>M5524A</b> | Rose and color label set   |
| <b>M5525A</b> | Gray and color label set   |

### Miscellaneous

|               |                                                                 |
|---------------|-----------------------------------------------------------------|
| <b>M1781A</b> | Test Load for use with M3507A Pad Cable                         |
| <b>M3725A</b> | Test Load for use with M3508A Pad Cable                         |
| <b>M3541A</b> | Carrying Case (includes 3 accessory pouches and shoulder strap) |
| <b>M3544A</b> | Data Card Tray                                                  |
| <b>M3545A</b> | Data Card and Tray                                              |
| <b>M3537A</b> | Bedrail Hook mount                                              |

## Specifications

### General

**Dimensions with pads:** 12.4 in. (W) x 8.3 in. (D) x 11.7 in. (H)  
(31.5 cm x 21.0 cm x 29.5 cm)

**Dimensions with paddles:** 13.4 in. (W) x 8.3 in. (D) x 13.6 in. (H)  
(34.0 cm x 21.0 cm x 34.5 cm)

**Weight:** 13.2 lbs. (5.99 kg) including pads, pads cable, full roll of paper, and battery. Incremental weight of external standard paddles and paddle tray is 2.5 lbs. (1.1 kg). Additional battery weighs 1.6 lbs. (0.73 kg)

### Defibrillator

**Waveform:** Biphasic Truncated Exponential. Waveform parameters adjusted as a function of patient impedance

**Shock Delivery:** Via multifunction defib electrode pads, or paddles

### Delivered Energy Accuracy:

| Selected Energy | Nominal Delivered Energy vs. Patient Impedance |     |     |     |     |     |     | Accuracy |
|-----------------|------------------------------------------------|-----|-----|-----|-----|-----|-----|----------|
|                 | Load Impedance                                 |     |     |     |     |     |     |          |
|                 | 25                                             | 50  | 75  | 100 | 125 | 150 | 175 |          |
| 1J              | 1.2                                            | 1.3 | 1.2 | 1.1 | 1.0 | 0.9 | 0.8 | ±2J      |
| 2J              | 1.8                                            | 2.0 | 2.0 | 1.9 | 1.7 | 1.6 | 1.5 | ±2J      |
| 3J              | 2.8                                            | 3.0 | 3.0 | 3.1 | 3.0 | 2.9 | 2.7 | ±2J      |
| 4J              | 3.7                                            | 4.0 | 4.0 | 4.1 | 4.2 | 4.2 | 4.0 | ±2J      |
| 5J              | 4.6                                            | 5.0 | 5.1 | 5.1 | 5.2 | 5.2 | 5.0 | ±2J      |
| 6J              | 5.5                                            | 6.0 | 6.1 | 6.2 | 6.3 | 6.3 | 6.1 | ±2J      |
| 7J              | 6.4                                            | 7.0 | 7.1 | 7.2 | 7.3 | 7.3 | 7.1 | ±2J      |
| 8J              | 7.4                                            | 8.0 | 8.1 | 8.2 | 8.4 | 8.3 | 8.1 | ±2J      |
| 9J              | 8.3                                            | 9.0 | 9.1 | 9.3 | 9.4 | 9.4 | 9.1 | ±2J      |
| 10J             | 9.2                                            | 10  | 10  | 10  | 10  | 10  | 10  | ±2J      |
| 15J             | 14                                             | 15  | 15  | 15  | 16  | 16  | 15  | ±15%     |
| 20J             | 18                                             | 20  | 20  | 21  | 21  | 21  | 20  | ±15%     |
| 30J             | 28                                             | 30  | 30  | 31  | 31  | 31  | 30  | ±15%     |
| 50J             | 46                                             | 50  | 51  | 51  | 52  | 52  | 50  | ±15%     |
| 70J             | 64                                             | 70  | 71  | 72  | 73  | 73  | 71  | ±15%     |
| 100J            | 92                                             | 100 | 101 | 103 | 104 | 104 | 101 | ±15%     |
| 120J            | 110                                            | 120 | 121 | 123 | 125 | 125 | 121 | ±15%     |
| 150J            | 138                                            | 150 | 152 | 154 | 157 | 156 | 151 | ±15%     |
| 170J            | 156                                            | 170 | 172 | 175 | 177 | 177 | 172 | ±15%     |
| 200J            | 184                                            | 200 | 202 | 206 | 209 | 209 | 202 | ±15%     |

**Charge Time:** Less than 5 seconds to 200 Joules with a new, fully charged Lithium Ion battery pack at 25°C. Less than 15 seconds when operating without a battery, using the M3539A AC power module alone at 90% rated mains voltage.

### Patient Impedance Range:

Minimum: 15 Ohm (internal defibrillation); 25 Ohm (external defibrillation)

Maximum: 180 Ohm

Note: Actual functional range may exceed the above values.

### Manual Mode

**Manual Output Energy (Selected):** 1-10, 15, 20, 30, 50, 70, 100, 120, 150, 170, 200J; maximum energy limited to 50J for internal defibrillation

**Controls:** On/Off Therapy Knob, Charge, Shock, Sync, Print, Mark Event, ECG Lead Select, Alarm Pause, Event Review

**Energy Selection Control:** Front panel therapy knob

**Charge Control:** Front panel button; button on external paddles

**Shock Control:** Front panel button; buttons on external or switched internal paddles

**Synchronized Control:** Front panel SYNC button

**Indicators:** Text prompts, audio alerts, QRS beeper, battery status, Ready For Use, external power/battery charging, Sync mode

**Armed Indicators:** Charging tone, charged tone, flashing shock button, and energy level indicated on display

### AED Mode

**AED Energy Profile:** 150 Joules nominal into a 50 Ohm test load

**Text and Voice Prompts:** Extensive text/audible messages guide user through configured protocol

**AED Controls:** On/Off, Shock

**Indicators:** Monitor display messages and prompts, voice prompts, battery status, Ready For Use, external power/battery charging

**Armed Indicators:** Flashing shock button, charged tone, and voice prompts

**ECG Analysis:** Evaluates patient ECG and signal quality to determine if a shock is appropriate and evaluates connection impedance for proper defibrillation pad contact

**Shockable Rhythms:** Ventricular fibrillation and certain ventricular tachycardias, including ventricular flutter and polymorphic ventricular tachycardia

**Shock Advisory Algorithm Sensitivity and Specificity:** Meets AAMI DF-

39

### ECG and Arrhythmia Monitoring

**Inputs:** Up to four (4) ECG waves may be viewed on display and up to 2 waves printed simultaneously. Lead I, II, or III is obtained through the 3-lead ECG cable. With a 5-lead cable, leads aVR, aVL, aVF, and V can also be obtained. Pads ECG is obtained through 2 multifunction defibrillation electrode pads.

**Lead Fault:** Leads Off message and dashed line appear on the display if an electrode or lead wire becomes disconnected. Lead Off indicator in wave sector.

**Pad Fault:** Dashed line appears on the display if a pad becomes disconnected.

**Heart Rate Display:** Digital readout on display from 15 to 300 bpm, with an accuracy of  $\pm 10\%$

**Heart Rate/Arrhythmia Alarms:** HR, Asystole, VFIB/VTACH, VTACH, Extreme Tachy, Extreme Brady, PVC rate, Pacer not capture, Pacer not pacing

### Hands Free Defibrillation Cable Length:

M3507A = 7 ft. (2.2 m)

M3508A = 7 ft. (2.2 m)

**ECG Cable Length:** 9 ft. (2.7 m)

**Common Mode Rejection:** Greater than 90 dB measured per AAMI standard for cardiac monitors (EC 13)

**ECG Size:** 2.5, 5, 10, 20, 40 mm/mV, autogain

### Frequency Response:

AC Line Filter: 60 Hz or 50 Hz

3-lead, 5-lead, and Pads:

Pads ECG for Display: Monitor (0.15-40 Hz) or EMS (1-30 Hz)

Pads ECG for Printer: Monitor (0.15-40 Hz) or EMS (1-30 Hz)

Leads ECG for Display: Monitor (0.15-40 Hz) or EMS (1-30 Hz)

Leads ECG for Printer: Diagnostic (0.05-150 Hz) or Monitor (0.15-40 Hz) or EMS (1-30 Hz)

12-lead:

ECG for Display: (0.05 - 150 Hz), (0.05 - 40 Hz), (0.15 - 40 Hz)

ECG for Report: (0.05 - 150 Hz), (0.05 - 40 Hz), (0.15 - 40 Hz), (0.05 - 150 Hz)

### Patient Isolation (defibrillation proof):

ECG: Type CF

SpO<sub>2</sub>: Type CF

EtCO<sub>2</sub>: Type CF

NBP: Type CF

External Defib: Type BF

Internal Defib: Type CF

### Display

**Size:** 8.4 in. diagonal (128 mm x 171 mm)

**Type:** TFT Color LCD

**Resolution:** 640 x 480 pixels (VGA)

**Sweep Speed:** 25 mm/s nominal (stationary trace; sweeping erase bar) for ECG and SpO<sub>2</sub>; 6.25 mm/s for CO<sub>2</sub>

**Wave Viewing Time:** 5 seconds (ECG)

### Battery

**Type:** 6.0 Ah, 14.8 V, rechargeable, Lithium Ion

**Dimensions:** 6.5 in. (H) x 3.8 in. (W) x 1.6 in. (D) (165 mm x 95 mm x 42 mm)

**Weight:** 1.6 lb. (0.73 kg)

**Charge Time:** Approximately 3 hours to 100%. Approximately 2 hours to 80%, indicated by battery gauge indicator. Charging the battery at temperatures above 45°C may degrade battery life.

**Capacity:** At least 5 hours of monitoring with ECG, SpO<sub>2</sub>, and CO<sub>2</sub> monitored continuously and NBP measured every 15 minutes, or at least 50 200-joule shocks (on a new, fully charged battery at room temperature, 25° C). At least 3.5 hours while pacing at 180 ppm at 160 mA and monitoring as described above

**Battery Indicators:** Battery gauge on battery, capacity indicator on display; flashing RFU indicator, chirp, and LOW BATTERY message appears on display for low battery condition\*

\*Low battery condition triggered with at least 10 minutes of monitoring time and 6 maximum energy discharges remain (with a new battery at room temperature, 25° C)

**Storage:** Storing the battery for extended periods at temperatures above 40°C will reduce battery capacity and degrade battery life.

### Thermal Array Printer

**Continuous ECG Strip:** The Print key starts and stops the strip. The printer can be configured to run real time or with a 10-second delay. The strip prints the primary ECG lead with event annotations and measurements.

**Auto Printing:** The printer can be configured to automatically print on Marked Events, Charge, Shock, and Alarm. When an alarm condition occurs, the unit prints the primary ECG wave and the alarming wave, if configured.

**Reports:** The following can be printed: Event Summary, 12-Lead, Operational Check, Configuration, Status Log, and Device Information

**Speed:** 25 or 50 mm/s with an accuracy of  $\pm 5\%$

**Amplitude Accuracy:**  $\pm 5\%$  or  $\pm 40$  uV, whichever is greater

**Paper Size:** 50 mm (W) by 30 m (100 ft.) (L)  
75 mm (W) by 30 m (100 ft.) (L)

### Noninvasive Pacing

**Waveform:** Monophasic Truncated Exponential

**Current Pulse Amplitude:** 10 mA to 175 mA (5 mA increments); accuracy 10% or 5 mA, whichever is greater

**Pulse Width:** 40 ms with  $\pm 10\%$  accuracy

**Rate:** 30 ppm to 180 ppm (10 ppm increments); accuracy  $\pm 1.5\%$

**Modes:** Demand or Fixed Rate

**Refractory Period:** 340 msec (30 to 80 ppm); 240 msec (90 to 180 ppm)

### SpO<sub>2</sub> Pulse Oximetry

**Range:** SpO<sub>2</sub>: 0-100%; pulse rate: 30 to 300 bpm

**Resolution:** 1%

**Display Update Period:** 1 sec. typical numeric update rate

**Accuracy with:**

M1191A sensor - 1 standard deviation 70% to 100%,  $\pm 2.5\%$

M1192A sensor - 1 standard deviation 70% to 100%,  $\pm 2.5\%$

M1193A sensor - 1 standard deviation 70% to 100%,  $\pm 2.5\%$

M1194A sensor - 1 standard deviation 70% to 100%,  $\pm 4.0\%$

M1195A sensor - 1 standard deviation 70% to 100%,  $\pm 4.0\%$

NELLCOR sensors - 1 standard deviation 80% to 100%,  $\pm 3.0\%$

**Pulse Rate Accuracy:** 2% or 1 bpm (whichever is greater)

**Alarm Range:**

Low Limit: 50 to 99% (Adult/Pediatric)

High Limit: 51 to 100% (Adult/Pediatric)

**Alarm Delay:** 10 seconds

### Noninvasive Blood Pressure

**Pressure Range:**

Systolic: 40-260 mmHg

Diastolic: 20-200 mmHg

**Initial Pressure:** 160 mmHg (Adult); 120 mmHg (Pediatric)

**Maximum Pressure:** 280 mmHg

**Overpressure Safety Limits:** Maximum of 300 mmHg

**Cuff Inflation Time:** 50 second maximum (pediatric or adult)

**Accuracy:**  $\pm 3$  mmHg

**Alarm Range:**

Systolic high limit: 35 - 270 (Adult), 35 - 180 (Pediatric)

Systolic low limit: 30 - 265 (Adult), 30 - 175 (Pediatric)

Diastolic high limit: 15 - 245 (Adult), 15 - 150 (Pediatric)

Diastolic low limit: 10 - 240 (Adult), 10 - 145 (Pediatric)

Mean high limit: 25 - 255 (Adult), 25 - 160 (Pediatric)

Mean low limit: 20 - 250 (Adult), 20 - 155 (Pediatric)

**Calibration Schedule:** yearly or every 10,000 cycles

**Auto Mode Repetition Time:** 1, 2.5, 5, 10, 15, 30, 60, or 120 minutes

**Measurement Time:** Auto/manual mode: 30 seconds (average) @ HR > 60 bpm, 170 seconds (maximum)

**Interconnect Tube Length:**

M1598B Connect tubing 5 ft. (1.5 m)

M1599B Connect tubing 10 ft. (3 m)

### End-Tidal Carbon Dioxide

**Range:** 0 to 99 mmHg

**Resolution:** 1mmHg (0.1kPa)

**Accuracy:** For values between 0 and 38 mmHg:  $\pm 2$  mmHg. For values between 39 and 99 mmHg:  $\pm 5\%$  of reading + 0.08% for every 1 mmHg (above 40 mmHg). Values read at sea level after  $\geq 20$  min. warm up

**Alarm Range:**

Low Limit: 10 to 94 mmHg (Adult/Pediatric)

High Limit: 20 to 95 mmHg (Adult/Pediatric)

**Calibration Schedule:** yearly or every 4,000 hours

**Sample Size:** 50 ml per min.

### Airway Respiration Rate

**Range:** 0 to 150 rpm

**Resolution:** 1 rpm

**Accuracy:**

0 to 40 rpm  $\pm 1$  rpm

41 to 70 rpm  $\pm 2$  rpm

71 to 100rpm  $\pm 3\%$

101 to 150rpm  $\pm 5\%$

**Alarm Range:**

Low Limit: 0 to 99 rpm (Adult/Pediatric)

High Limit: 10 to 100 rpm (Adult/Pediatric)

**Apnea Alarm Time:** 10-40 seconds, in increments of 5

### Calibration Gas for CO<sub>2</sub> Measurement System

**Ingredients:** 5% Carbon Dioxide, 21% Oxygen, 74% Nitrogen

**Cylinder Size:** BD

**Method of Preparation:** Gravimetric

**Blend Tolerance:** 0.03%

**Accuracy:** 0.03% absolute

**Moisture:** 10 PPM Maximum

**Expiration Period:** 2 years

**Pressure:** 144 PSIG, Volume: 10L

### 12-Lead ECG

**Inputs:** With a 10-Lead cable, leads I, II, III, aVR, aVL, aVF, V/C1-V/C6 can be obtained. All 12-Lead ECG waves can be viewed on the display simultaneously. All 12 leads can be printed on the strip chart printer in 3x4 format.

**ECG Bandwidth Filters:**

0.15 - 40 Hz

0.05 - 40 Hz

0.05 - 150 Hz

**12-Lead ECG Transmission:** Data card to a PC

**Receiver Type:** HeartStart Event Review Pro, TraceMaster ECG Management System Version A.04.02, and TraceMaster Vue

### Patient Data Storage

**Internal Event Summary:** Up to 8 hours of continuous ECG waveforms and events and up to 62 megabytes of patient data

**Data Card Event Summary:** Up to 8 hours of continuous ECG waveforms and events and up to 62 megabytes of patient data

## Environmental

**Temperature:** 0° C to 45°C operating, -20° to 70°C storage

**Humidity:** Up to 95% relative humidity

**Altitude:**

Operating: 0 to 15,000 ft. (0 to 4,500 m)

Storage: 0 to 15,000 ft. (0 to 4,500 m)

**Shock:**

Operating: Half-sine waveform, duration < 3 ms, acceleration > 145 g, 1 time on all six faces

Storage: Trapezoidal waveform, acceleration  $\geq 30$  g, velocity change=742 cm/s  $\pm 10\%$  on all six faces

**Vibration:**

Operating: Random vibration, 0.30 grms, 5-500 Hz for  $\geq 10$  minutes/axis

PSD=0.0002 g<sup>2</sup>/Hz from 5 to 350 Hz, -6 dB/octave slope from 350 to 500 Hz

Storage: Random vibration, 2.41 grms, 5-500 Hz for  $\geq 10$  minutes/axis

PSD=0.02 g<sup>2</sup>/Hz from 5 to 100 Hz, -6 dB/octave slope from 100 to 137 Hz, 0 dB/octave slope from 137 to 350 Hz, -6 dB/octave slope from 350 to 500 Hz

**Solids Resistance:** IP2x

**Water Resistance:** Meets EN 60601-2-4. Water testing performed with cables connected to the device

**EMC:** Complies with the requirements of standard EN 60601-1-2:2001

**Safety:** Complies with the requirements of applicable safety standards

**Other Considerations:** Device not suitable for use in the presence of concentrated oxygen or a flammable anesthetic mixture with air, oxygen, or nitrous oxide

**Mode of Operation:** Continuous

**AC Power Module:**

Input: 100-240 VAC, 50-60 Hz, 1-0.46 A (Class 1)

Output: 18V, 5A, 90W

**Battery:** 14.8 V Rechargeable, Lithium Ion

**Philips Medical Systems is part  
of Royal Philips Electronics**

**UNITED STATES:**

Philips Medical Systems  
Cardiac and Monitoring Systems  
3000 Minuteman Road  
Andover, MA 01810  
Tel: (800) 934-7372

**CANADA:**

Philips Medical Systems  
2660 Matheson Blvd. E.  
Mississauga, Ontario L4W 5M2  
Tel: (800) 291-6743

**EUROPE, MIDDLE EAST AND AFRICA:**

Philips Medizinsysteme Böblingen GmbH  
Cardiac and Monitoring Systems  
Hewlett-Packard Str. 2  
71034 Böblingen  
Germany  
Fax: (+49) 7031 463 1552

**LATIN AMERICA HEADQUARTERS:**

Philips Medical Systems  
5200 Blue Lagoon Drive  
9th Floor  
Miami, FL 33126  
Tel: (305) 267-4220

**ASIA PACIFIC HEADQUARTERS:**

Philips Medical Systems  
26/F Hopewell Centre  
17 Kennedy Road  
Wanchai  
Hong Kong  
Tel: (+852) 2821 5888  
Fax: (+852) 2527 7136

[www.medical.philips.com](http://www.medical.philips.com)

© 2004 Koninklijke Philips Electronics N.V. All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder. Koninklijke Philips Electronics N.V. reserves the right to make changes in specifications or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. Published December 2004  
4522 981 98961