

# MEKANT

*Ventilator System*

## User Manual

---



MEK ICS Co., Ltd.  
5F, A-BLDG, WooLimLion, #144-3, Sangdaewon-dong,  
Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea 462-725  
TEL : +82-70-7119-2500  
FAX : +82-31-735-2761  
URL : <http://www.mek-ics.com>

CE . 47 .



## Definitions

This operator's manual uses three special indicators to convey information of a specific nature. They include as follow:



### Warning

Indicate a condition that can endanger the patient or the ventilator operator.



### Caution

Indicate a condition that can damage the equipment.



### NOTE :

Indicates points of particular emphasis that make operation of the ventilator more efficient or convenient.



<b>Section 1 : Introduction .....</b>	<b>1</b>
1-1. Warnings and Cautions .....	2
1-2. General description .....	5
1-2-1. Front picture and the explanation .....	5
1-2-2. Back side picture and explanation .....	7
1-2-3. Left side picture and explanation .....	9
1-2-4. Right side picture and explanation .....	10
1-3. Specifications .....	11
1-3-1. Features .....	11
1-3-2. Electrical .....	11
1-3-3. Functional Specifications .....	12
1-3-4. Pneumatic .....	16
1-3-5. Physical data .....	16
1-3-6. Environmental data .....	16
1-4. Standards and safety .....	17
1-5. Symbols .....	18
1-6. Maintenance .....	21
 <b>Section 2 : Installation and set up .....</b>	 <b>23</b>
2-1. Installation .....	23
2-2. External electrical supply .....	23
2-2-1. AC description .....	23
2-2-2. External battery .....	25
2-3. Internal battery .....	26
2-3-1. Internal battery Indicator .....	27
2-3-2. How to change internal battery .....	28
2-4. Oxygen and Air supply .....	29
2-5. Filters .....	30
2-6. Transportation .....	31
2-7. Fuse replacement instruction .....	32



<b>Section 3 : Operating instructions</b>	<b>35</b>
3-1. Front panel	35
3-1-1. Rotational Control Knob	35
3-1-2. Keypad	36
3-1-3. LED Indicators	37
3-2. Operating	37
3-2-1. Operating procedure and how to use	37
3-3. Display	40
3-3-1. Status information	41
3-3-2. Waveform and alarm status	42
3-3-3. Menu and setting information	43
3-3-4. TOUCH information	44
3-3-4-1. VENT TOUCH screen	45
3-3-4-2. Waveform screen	46
3-3-4-3. Menu screen	47
3-3-4-4. Setting screen	48
3-4. Start <b>MV2000</b>	49
3-4-1. VENT MODE	50
3-4-1-1. VACV	51
3-4-1-2. VSIMV	52
3-4-1-3. PACV	56
3-4-1-4. PSIMV	58
3-4-1-5. SPONT	61
3-4-1-6. tBiLEVEL	63
3-4-1-7. NIV	66
3-4-1-8. PRVC	68
3-4-1-9. AUTO	71
3-4-2. Monitoring	72
3-4-3. Graphics	73
3-4-4. Tools	76
3-4-5. Events	77
3-4-6. SYSTEM	78
3-4-6-1. SETUP	78
3-4-6-1-1. BWF	79
3-4-6-1-2. Nebu. Time	79
3-4-6-1-3. O <sub>2</sub> SENSOR	79
3-4-6-1-4. Trend Init.	79
3-4-6-1-5. BTPS	80
3-4-6-1-6. Sound Vol.	80
3-4-6-1-7. Altitude	80
3-4-6-1-8. Air Temp.	80



3-4-6-2. PATIENT .....	81
3-4-6-3. Date .....	82
3-4-6-4. Calibration .....	83
3-4-6-4-1. SELF CALIBRATION .....	84
3-4-6-4-2. O <sub>2</sub> CALIBRATION .....	84
3-4-6-4-3. INS-EXP CALIBRATION .....	84
3-4-6-4-4. EtCO <sub>2</sub> CALIBRATION .....	85
3-4-7. ALARM SET .....	86

## Section 4 : Alarms ..... 87

4-1. General Information .....	87
4-2. Alarm Message .....	88

<b>Appendix</b> .....	90
A. Operating theory .....	90
B. Communications .....	92
C. Parts accessories .....	95

## Service Calls ..... 97



## Section 1

<b>Section 1 : Introduction</b>	1
1-1. Warnings and Cautions	2
1-2. General description	5
1-2-1. Front picture and the explanation	5
1-2-2. Back side picture and explanation	7
1-2-3. Left side picture and explanation	9
1-2-4. Right side picture and explanation	10
1-3. Specifications	11
1-3-1. Features	11
1-3-2. Electrical	11
1-3-3. Functional Specifications	12
1-3-4. Pneumatic	16
1-3-5. Physical data	16
1-3-6. Environmental data	16
1-4. Standards and safety	17
1-5. Symbols	18
1-6. Maintenance	21

## Section 1 : Introduction

This operator's manual explains the way of installation and usage for MEKANT. Notice for use explains general information to be careful while using the product and other notices are written with each function. Users are recommended to be aware of notice for use and all menus before using the product.

This operator's manual is allowed to use by doctors or authorized persons who are permitted to use, or designed to use by qualified and trained persons for using this products.

Service, test, correction and main continuous working will be written in service manual of MEKANT Ventilator. However, service manual is allowed to be possessed by authorized persons by MEKICS.



## 1-1. Warnings and Cautions



### Warning

- Please make sure to use MEKANT by authorized persons and read the manual before use. Check the setting figure of various sensors and devices before use. Users shall check out whether the environment is OK before use.
- Please be aware of the manual provided by the manufacturer when it is used or kept.
- Do not use it near inflammable anesthetic to avoid any possible explosion.
- Do not disassemble or open the main body to avoid electrical impulse. Disassembly of the equipment shall be allowed by service staffs approved by MEKICS Co., Ltd.
- Users shall be careful with circuit joint to be connected with patients.
- Always shall be careful with leakage of electrical current and consult with hospital technical engineer if any doubt.
- In case of abnormal working or any damage in the equipment, it shall not be used for patients.
- Do not touch patient, table or the equipment if a defibrillator is used.
- Disposable stuffs and disposable circuit for patients shall not be reused.
- Reusable stuffs used for patients shall be sterilized or washed for reuse.
- Check out the level before setting alarm. Alarm limitation shall be done for the current selected group. Alarm setting can be changed by level change.
- Do not correct calibration while used for patient. Please correct it after disconnection from patient if any doubt.
- Insulate power source in case that the equipment is wet and call a hospital technical engineer.
- Please connect earth cable to earth connection port to protect patients and the product.
- Do not disassemble the product in any case.
- Do not install the various kinds of sensor to patients too strongly and not to be entangled.
- Do not reuse it if you find this symbol (ⓧ) on the accessory, which means disposable stuff.
- MEKANT shall be tested regularly if it is working in normal condition.
- Accessories are aseptic products, so the package shall not be torn before use and products shall be free from any defects or cracks.



**Caution**

- MEKANT as respiration assistant equipment for adult helps for normal breathing by setting breathing way, breathing volume, breathing times, oxygen density, breathing rate etc in order to solve hypoxia to be appeared from emphysema of the lungs patients or respiratory disease patients. It is respiratory assistant equipment with various alarm function and safety device for safe respiration management, but above all, should be used together with other clinical information for patients' safety and diagnosis.
- Avoid the place where there is moisture.
- Do not handle the equipment with wet hands.
- Avoid the place under direct rays.
- Use the equipment within the limits of temperature 10 ~ 40 centigrade and humidity 30 ~ 75%
- Do not place the equipment near electronic heater.
- Do not place the equipment under high humidity and bad ventilation.
- Do not place the equipment where it is possible to take vibration or sudden shock.
- Do not place the equipment where to be damaged by explosive gas or chemicals.
- Protect the equipment from dust or metal scrap.
- Pull out power code with plug grasped, instead of cable.
- Check out the main equipment and sensor at the same time after cleaning the equipment. Do not use if there is any damage in the equipment.
- Make the cable and sensor clean after use. Clean them by towel with alcohol soaked and do not use benzene or gasoline.
- Do not use the equipment in dangerous place. Always be careful with the possibility of bad influence by environment.
- It is possible to be flew more current than allowed when the assistant device is used. Consult with hospital engineer or authorized local dealer.
- Check out the basic auxiliary device regularly. Do not use a damaged or electrical unstable device. Consult with hospital engineer if any doubt.
- Do not use the goods beyond standard certified one. This product shall be used under standard certification. Manufacturer does not take responsibility for the problem occurred while using the product beyond standard certification.
- The product can be disassembled by educated and authorized staffs. If a user violates this item, after sales service is not provided.



- Do not install the various kinds of sensor to patients too strongly and not to be entangled.
- Do not place the product where it is possible to fall to a patient.
- Do not move the equipment when power cord or high pressure hose is connected. It will be harmed for patient or the equipment has an error by itself
- If a measured figure is turned out not to be correct, consult with hospital engineer or authorized local dealer.
- Connect the equipment isolated from other devices electronically for patients' electrical safety.
- Accessories of MEKANT shall be used with the ones provided by MEKICS or MEKICS's designated dealer or listed in the manual. MEKICS or MEKICS's designated dealer provide standard certified products only.
- MEKANT shall be tested regularly if it is working well.
- Turn off main switch located backside of the product after use.
- Place the product in stable position.
- Use the sensor and cable provided and designated by manufacturer, and consult with hospital engineer or authorized local dealer if you want to use other products.
- Replace the disposable sensor after use.
- Set up the equipment when it is needed.
- Minimize patient's move when it is used.
- Make high pressure hose and cable clean after use. Clean them by towel with alcohol soaked and do not use benzene or gasoline.
- The Oxygen Air pressure which is impressed pneumatic block should be between 30PSI and 40PSI. Otherwise, it would affect the equipment.
- Move the equipment with fastened to avoid any possible damage to other products or people by the monitor or cable.
- Basically, MEKANT is designed not to have functional disorder by surrounded electronic wave. But, it is possible to be occurred any disorder due to high frequency of surrounding. It may be possible to calculate wrong figure, lose present figure or show noisy wave. Please ask MEKICS technical service center or local distributor if you need help.
- Please turn off the power to recharge the battery when it is not used.
- MEKANT ventilator system is a product conforming to **Class II a**.



## 1-2. General Description

MEKANT as respiration assistant equipment for adult helps for normal breathing by setting breathing way, breathing volume, breathing times, oxygen density, breathing rate etc in order to solve hypoxia to be appeared from emphysema of the lungs patients or respiratory disease patients.

It is a respiration assistant equipment with various alarm function and safety device for safe respiration management, but above all, should be used together with other clinical information for patients' safety and diagnosis

### 1-2-1. Front picture and the explanation




Figure 1-1. Front picture



Part	Description
Ⓐ CIRCUIT ARM SOCKET	Install Circuit Arm for patient.
Ⓑ LCD MONITOR	Show the condition of various ventilators and patients.
Ⓒ KEY FRONT	Key button and encoder.
Ⓓ INSPIRATION PORT	Joint hole to connect to patient.
Ⓔ EXPIRATION PORT	Joint hole to connect from patient.
Ⓕ EXPIRATION LEVER	Lock or unlock lever for Expiration Assembly.
Ⓖ CARRIAGE GRIP	Using to move the ventilator.
Ⓗ NEBULIZER PORT	Connection to spray medicines

Table 1-1. Front picture description

**Warning**

Be careful because heated exhalation assembly due to moisture during exhalation.



## 1-2-2. Back side picture and the explanation




Figure 1-2. Back side picture



Part	Description
㉑ <b>SPEAKER</b>	For alarm and sound.
㉒ <b>VGA</b>	VGA cable connection hole to show patient condition on the other Monitor.
㉓ <b>COMM</b>	Cable connection hole to connect PC or central connection.
㉔ <b>TO PNEUMATIC</b>	Connection hole to connect to the main body.
㉕ <b>AC INLET</b>	Connection hole for AC power input.
㉖ <b>EXTERNAL BATTERY</b>	Connection hole for external battery.
㉗ <b>POWER SWITCH</b>	Power on/off switch for common use.
㉘ <b>O<sub>2</sub> INLET</b>	Connection hole to connect to main body using high pressure hose to approve O <sub>2</sub> gas.
㉙ <b>Air INLET</b>	Connection hole to connect to the main body, using high pressure hose to approve air gas. Device to filter the moisture in the Air.
㉚ <b>TANK RELEAF</b>	Outlet to exhaust the over pressed air when the pressure of air tank is over than 40 PSI.
㉛ <b>BATTERY DOOR</b>	Door to replace the internal battery.
㉜ <b>TO MONITOR</b>	Connection for monitor.
㉝ <b>COMM</b>	Connection hole for pneumatic communication.

Table 1-2. Back side picture description

**Warning**

You should be careful with battery's polarity when it is exchanged.  
Products can be damaged when polarity is connected wrongly.

It should be careful of battery's terminal to be short on the metal part



1-2-3. Left side picture and the explanation




Figure 1-3. Left side picture

Part	Description
a LOCKER	Locker to disassemble main body with cart. Pull to unlock and push to lock.

Table 1-3. Left side picture description



1-2-4. Right side picture and the explanation




Figure 1-4. Right side picture

Part	Description
a O <sub>2</sub> CHANGE DOOR	Door to be used for exchanging O <sub>2</sub> Sensor.
b INSPIRATION RELIEF	Safety valve to exhaust by software, hardware and mechanically when the equipment is out of order.

Table 1-4. Right side picture description



Warning

Do not stop up safety valve to exhaust by software, hardware and mechanically when the equipment is out of order.



## 1-3. Specifications

### 1-3-1. Features

User friendly POP-up Menu with rotary and push/pull Encoder.

100% O<sub>2</sub>, Inspiration Pause, Expiration Pause, Key-Lock ,  
Manual and/or Auto Nebulizer

Dual Alarm ( Speaker / Buzzer ) system

External battery ( for Mobile application )

### 1-3-2. Electrical Specifications

AC input voltage	100 - 120 VAC, 10A 50/60 Hz
	200 - 240 VAC, 5A 50/60 Hz

AC input FUSE Rating	250V/6.3A
----------------------	-----------

Power consumption	84 W max
-------------------	----------

External DC input voltage	12V @ 4A Pb battery
---------------------------	---------------------

Internal battery	Battery ; 12V / 4A Pb battery Operating time: 3 hours minimum . Recharge time; 4 Hours max. Please change it after 3 years use. It can be used for 30 minutes normally. On the other hand, it will be more than 30 minutes if power consumption is small.
------------------	---

Communication	Central interface ; RS-232 serial interface Baud rate 115,200 bps, Upgrade (Main board, Pneumatic board)
---------------	--

Display	12.1" LCD monitor, 800*600,
---------	--------------------------------



---

## 1-3-3. Functional Specifications

### **MODE**

V-ACV	Volume Assist Controlled Ventilation Controlling patient's mechanical artificial respiration based on air volume and adjusting inspiration and exhalation time by force.
V-SIMV	Volume based Synchronized Intermittent Mandatory Ventilation Controlling artificial respiration based on air volume in accordance with patient's voluntary respiration.
SPONT	Spontaneous Ventilation Give pressure support centering on voluntary respiration.
P-ACV	Pressure based Assistance Controlled Ventilation Controlling patient's mechanical artificial respiration based on pressure.
P-SIMV	Pressure based Synchronized Intermittent Mandatory Ventilation Controlling artificial respiration based on pressure in accordance with patient's voluntary respiration.
tBiLEVEL	Same as SPONT and exist Upper PEEP (Refer to SPONT mode)
NIV	Non Invasive Ventilation During artificial respiration by mask type, not tubing type, P-ACV working in allowance of air leakage. (Refer to P-ACV mode)
PRVC	Pressure Regulated Volume Controlled Ventilation Breaths are delivered mandatorily with a constant flow to assure present volumes.
AUTO	Start quickly according to body weight. (Working in P-ACV mode)



**SETTING PARAMETER**

	<b><u>Setting range</u></b>	<b><u>Condition</u></b>	<b><u>Default</u></b>
BW	3 kg ~ 150 kg		10 ml/kg
Pressure	5 ~ 80 cmH <sub>2</sub> O	±(1.1 + 4% of the actual reading) cmH <sub>2</sub> O	5 cmH <sub>2</sub> O
Plateau Time (pause)	0 ~ 2 sec		0.1 sec
PEEP	0 ~ 45 cmH <sub>2</sub> O	±(1.1 + 4% of the actual reading) cmH <sub>2</sub> O	5 cmH <sub>2</sub> O
High PEEP	5 ~ 45 cmH <sub>2</sub> O	±(1.1 + 4% of the actual reading) cmH <sub>2</sub> O	+15 cmH <sub>2</sub> O
Low PEEP	0 ~ 45 cmH <sub>2</sub> O	±(1.1 + 4% of the actual reading) cmH <sub>2</sub> O	+5 cmH <sub>2</sub> O
High PEEP time	0.2 ~ 24 sec	±10%	1 sec
Low PEEP time	Infinite		
	<b><u>Setting range</u></b>	<b><u>Condition</u></b>	<b><u>Default</u></b>
Inspiration time	0.2 ~ 9.9 sec	±10%	1 sec
Fend	25 ~ 100 %		100 %
TRISE	0.1 ~ 0.5	±0.05	0.1 sec
EnSENSE %	10 ~ 60 %		30 %
ExSENSE %	10 ~ 60 %		30 %
Pressure Support	0 ~ 45 cmH <sub>2</sub> O	±(1.1 + 4% of the actual reading) cmH <sub>2</sub> O	0 cmH <sub>2</sub> O
Trigger sensitivity	0.5 ~ 20 cmH <sub>2</sub> O	±(1.1 + 4% of the actual reading) cmH <sub>2</sub> O	5 cmH <sub>2</sub> O
Pressure			
Trigger sensitivity	0.5 ~ 20 L/min	±15% (V <sub>MIN</sub> > 10 lpm)	5 lpm
Flow		±0.4 lpm (V <sub>IMIN</sub> ≤ 10 lpm)	
Volume	5 ml ~ 2500 ml	±15% (V <sub>TIDAL</sub> > 100 ml) ±20 ml (V <sub>TIDAL</sub> ≤ 100 ml)	10 ml/kg
Rate	2 ~ 120 bpm	±2 bpm	20 bpm
O <sub>2</sub> %	21 ~ 100 %	±2%	30 %



	<u>Setting range</u>	<u>Default</u>
Nebulizer	10 ~ 180 min Limited 200 ml Flow trigger base flow cannot use	off
Inspiration pause	Resistance measurement	off
Expiration pause	Compliance measurement	off
Manual Inspiration	Setting tidal volume	off
100% O <sub>2</sub>	3 minutes	off

### AUTO MODE(Emergency Start Ventilation)

	<u>Setting range</u>
I : E	1:2
O <sub>2</sub>	60 %
Mode	PACV
Inspiration Pressure	25 cmH <sub>2</sub> O
Inspiration time	1 sec
Expiration time	2 sec
Respiration rate	20 bpm
PEEP	3 cmH <sub>2</sub> O
Trigger type	Pressure
Trigger sensitivity	3 cmH <sub>2</sub> O
Rise %	50 %
Alarm limit	Auto set
E-sense	30 %



**Alarms**

	<b><u>Setting range</u></b>
PRESSURE HIGH	Δ cmH <sub>2</sub> O ~ 12.0 cmH <sub>2</sub> O
PRESSURE LOW	0 cmH <sub>2</sub> O ~ Δ.0 cmH <sub>2</sub> O
VTIDAL HIGH	Δ ml ~ 2Δ.0 ml ,OFF
VTIDAL LOW	0 ml ~ 2Δ.0 ml
PEEP HIGH	6 cmH <sub>2</sub> O ~ 4Δ cmH <sub>2</sub> O
PEEP LOW	0 cmH <sub>2</sub> O ~ 7 cmH <sub>2</sub> O
O <sub>2</sub> HIGH	1% ~ 100%, OFF
O <sub>2</sub> LOW	0% ~ 100%
R-R HIGH	3 ~ 1Δ.0 BPM
R-R LOW	2 ~ 149 BPM
VMIN HIGH	0.1 ~ Δ.0 LPM
VMIN LOW	0 ~ 49.9 LPM
AIR FAIL	-
O <sub>2</sub> FAIL	-
AIR LEAK	0 ml ~ Δ.0 ml
AIR OBST	-
APNEA	OFF, 10 ~ 60 sec
VENT INOP	-
O <sub>2</sub> SENSOR ERR	-
CIRCUIT OPEN	-
BATTERY MODE	-
LOW BATTERY	DC power loss
COM ERROR	



### 1-3-3. Pneumatic Specifications

O <sub>2</sub> and air supply Pressure range	35 ~ 90 psi
	Flow 150 L/min (STPD, dry required)
	100 % O <sub>2</sub> Suction button– auto calibration
Oxygen sensor life	10,000 hours or 1 Years of use, normal
Safety pressure Pneumatic	30 psi +/- 10 %
Internal pressure Limit	
Inspiration Pressure Limit	85 cmH <sub>2</sub> O +/- 10% ;
	Mechanical release valve
Double safety mechanism for over pressure	

### 1-3-4. Physical Data

Dimensions of device	38 W x 43 D x 69 H cm
	(Width x depth x height of device)
Weight of device	35 Kg (with Battery)
Dimensions of Cart	49 W x 56 D x 77 H cm
	(Width x depth x height of device)
Weight of Cart	28 kg
Connectors	
Inspiratory limb connector	ISO 22mm conical male
Expiratory limb connector	ISO 22mm conical male
Air and Oxygen inlets	DISS male / female



### 1-3-6. Environmental Data

Operating environment	Temperature : 10 to 40 °C (50 ~ 104°F)
Atmosphere	700 to 1060 mbar (10.2 ~ 15.4 psi)
Altitude	-443 ~ 3280 m (-1450 ~ 10,000 ft)
Storage environment	
Temperature	-10 to 60 °C
Relative humidity	10 to 85 %
Atmosphere	500 to 1060 mbar (7.3 ~ 15.4 psi)
Altitude	up to 6560 m (20,000 ft)
Oxygen/Air inlet supply pressure	35 ~ 90 psi

### 1-4. Standard and Safety requirements














MEKANT is in conformity with the following tables. Obtainment of license and required testing has been completed.

EN60601-1:1990	General requirements for Safety
EN60601-1-1:2001	Electromechanical Safety
EN60601-1-2:2001	Electromagnetic Compatibility Requirement and tests
IEC 60601-2-12:2001	Particular requirements for the safety of lung ventilators for medical use










## 1-4. Symbols and Labels

Descriptions for marking on MEKANT ventilator and its operational manuals.

No	Name	Symbol	Explanation
1	POWER ON/OFF		Provision of power by switching on and off the system. IEC 601-1
2	Equipotential-ity		It is indicated terminal which do not need ground electric potential to make equipotential by connecting with each part of system or equipment. IEC 601-1, IEC 617-5.21
3	Attention, Consult accompanying documents		Indicates the necessity to refer to the operational manual prior to equipment operation. IEC 601-1
4	stand-by		Indicates the stand-by mode of the system. IEC 617-5.0.9
5	Start (of action)		Indicate the activation of the equipment. IEC 617-5.1.4
6	Direct current		Indicates appropriateness for the equipment against direct current and designate the connecting part. IEC 601-1, IEC 617-5.31
7	Battery check		Battery capacity indicator
			<div>   </div> <div> FULL Charging 80% Charging </div>
			<div>   </div> <div> 50% Charging 20% Charging </div>
			<div>   </div> <div> Low Charging </div>



No	Name	Symbol	Explanation
1	AC Power plug		Indicates the AC plug without ground earth. IEC ۴۱۷-۵۵۳۴Pr
9	Battery		It should be indicated rating sticker that equipment is suitable for only direct current (D.C). Also, indicate it is related connecting part. IEC ۶۰۱-۱, IEC ۴۱۷-۵۰۳۱
1۰	EXT. BATT	EXT. BATT	Indicates external battery for vehicle use.
11	Manual control		Indicates manual operation of the apparatus. ISO ۷۰۰۰-۰۰۰۹۶
1۲	CLEAR		Clears the menu tree from the display.
1۳	Dangerous voltage		Warns the high voltage. IEC ۶۰۱-۱, IEC ۴۱۷-۵۰۳۶, ISO ۳۸۶۴
1۴	۱۰۰% O <sub>۲</sub> MODE	O <sub>۲</sub> ۱۰۰%	Indicate ۱۰۰% O <sub>۲</sub> MODE. This mode can be ON / OFF by O <sub>۲</sub> ۱۰۰% KEY
1۵	NEBULIZER MODE		Indicate NEBULIZER MODE. This mode can be ON / OFF by NEBULIZER KEY, and the operating time is set from 10 to 180 minutes.
1۶	MANUAL INSPIRATION KEY	MANUAL 	Supply the respiration to the patient when pressing the MANUAL INSPIRATION KEY.






No	Name	Symbol	Explanation
17	INSPIRATION HOLD		Increase INSPIRATION TIME to 1 second.
18	EXPIRATION HOLD		Increase EXPIRATION TIME to 2 second.
19	ALARM SOUND OFF Or MUTE		Indicate ALARM SOUND OFF or MUTE. This function can be ON / OFF by pressing 2 MIN KEY.

Table 1-5. Description of Symbols



## ١-٤. Maintenance

The following table shows the periodic check point and maintenance method for MEKANT ventilator system and the circuit system.

Period	Subject	Check point
Daily	Bacteria Filters - Inspiratory Filter - Expiratory Filter	Check if all filters are in good position.
	Patient Circuit System	Check the connection, leakage, Blockage and moisture.
	Collector Vial, Gas Supply Water Traps	Check if there is moisture and remove it if necessary.
Minimum Every half month or when changed the patient	Patient Circuit System, all Bacteria Filters and Humidification Chamber	Exchange everything, and cleaning and sterilizing used one.
Every Half year	All system	Conduct EST
Every year or After ١٠٠ times autoclaving	All Bacteria Filters	Exchange to new filter and scrap used one
Between ١ year and ٢ year	Oxygen Sensor	Exchange to new one
When necessary	Oxygen Sensor	Conduct Oxygen sensor calibration by pressing ١٠٠% O <sub>2</sub> /CAL.
Every ١٠,٠٠٠ hours	Various components	Exchange appropriate components using ١٠,٠٠٠ hours PM kit.

Table 1-6. Check and Maintenance Period



**Caution**

User should check the installation status and problem of the filters and patient circuit system with naked eyes every day. In case of any problem, replace or exchange the parts in problem right away.

When patient is changed or one patient uses it for a long period, exchange the patient circuit system including bacteria filter in every 15 days, and clean and sterilize the used one.

All the components of patient circuit system excluding bacteria filter can be sterilized by any way among ETO(Ethylene Oxide) gas sterilization, Steam Autoclave and Chemical sterilization. But, only Steam Autoclave is available for bacteria filter.

Sterilization shall be done for about 20 minutes at 132 degree centigrade. Exchange to new bacteria filter after one year use or 100 times of Autoclave use.

Oxygen sensor is valid for 2 years, but it can be shortened in case of using at the condition of high FiO<sub>2</sub> or high temperature.

---



<b>Section 2 : Installation and set up</b>	22
2-1. Installation	23
2-2. External electrical supply	23
2-2-1. AC description	23
2-2-2. External battery	25
2-3. Internal battery	26
2-3-1. Internal battery Indicator	27
2-3-2. How to change internal battery	28
2-4. Oxygen and Air supply	29
2-5. Filters	30
2-6. Transportation	31
2-7. Fuse replacement instruction	32

## Section 2 : INSTALLATION AND SET UP

### 2-1. Installation

When installing and moving the equipment, it is required to stay shock free.

### 2-2. External electrical supply

MEKANT is capable to receive both AC and DC power supply.

#### 2-2-1. AC description

External AC connector is located on the rear side of the apparatus.

Refer to Figure 2-1.




Figure 2-1. External AC and DC Power Supply Sockets



**Warning**

Prior to applying AC or DC to the ventilator, the line frequency has to be adjusted to the local standard otherwise it might cause damages to the apparatus.

Power Cord has to be plugged in to the concentric plug on the ground.

Replace covered wire in case of wear and tear.

Please apply the ground cable in order to protect the patient and the apparatus.

Maintain the same level of input & output voltage.

---

**Caution**


Cable Clips are used in order to affix the AC code and also to prevent disconnection.

---

When connected to DC or AC, LED display will be turned on to indicate the connection.



### 2-2-2. External battery



Pin	Name	Direction	Explanation
①	Not connected		Not used
②	GND	→	Pins for the ground.
③	+VCC	←	Pins to connect external battery
④	Not connected		Not used

Table 2-2. Pin Setting of MV2000 Ventilator External Battery



#### Warning

Watch for polarity when connecting the battery.

Note that when external battery is used, internal battery will not be charged.

It should be used for rating voltage, 12V --- 7A and Electric current should be more than 1A.



#### Caution

Cover the Port Pin by putting the cap when the external battery is not in use.

Do not make pin to be shorted because internal Battery is damaged when EXT. BATT is opened, or each pin is shorted.

Please put a cap next to terminal when EXT. BATT is not used.



#### NOTE :

Fuse is built in the machine.

ATO® Fuse (Fast-Acting Type), Catalog Number 257 015, Ampere Rating (A) 5A, Body color code : blue

Vehicle outlet power supply is used.



## 2-3. Internal battery

Internal battery for MV 2000 is used for the provision of power when moving the patients. Recharge of the battery is done while AC power connected as the apparatus simultaneously utilizes the power from the same source of AC. In case of power failure, the operation will be done by internal battery.



### Warning

When replacing the battery, watch for polarity. Red for +, Black for -. Inappropriate connection might cause damages to the apparatus.

Short-circuit might cause explosion of the battery or harm human body due to high voltage.

It should be used for rating voltage, 12V --- 7A and Electric current should be more than 1A.

We recommend that you should contact the distributor or CS team if you find any problems in the internal battery.



### Caution

When difficulties of power supply arise or when transferring the patients, the battery needs to be charged. Replace the battery if it is in deficiency.

In case of the lowest battery voltage, there may be an error in the measured value.

In case of low battery, use the machine after connection to AC common electric power.



### NOTE :

Fuse is built in the machine.

ATO® Fuse (Fast-Acting Type), Catalog Number 257 015, Ampere Rating (A) 5A, Body color code : blue

Internal Battery may be affected by the number of times for discharge, surrounding temperature, charging voltage. Generally it can be used more than 2 years.



### ㄣ-ㄣ-1. Internal battery indicator

It displays battery capability. This is displayed on the upper right hand corner of the LCD display. The battery capacity will be displayed in 4 main statuses.







Battery Indicator	Level Status	Explanation
	Full Charging	
	80% Charging	
	60% Charging	
	40% Charging	
 	Low Charging	In case of discharging status, the system will be automatically turn off to protect the battery 4 minutes later.

Table ㄣ-3-1. Internal battery Indicator



2-3-2. How to change internal battery



No.	Figure	Explanation
①		Unscrew 4 screws on the lower rear cover that is for Battery
②		You can see the internal battery. After changing battery, assemble the battery cover.

Table 2-3-2. How to change internal battery



Warning

When replacing the battery, watch for polarity.

Red for +. Black for -.

Inappropriate connection might cause damages to the apparatus.

Short-circuit might cause explosion of the battery or harm human body due to high voltage.



Caution

All of bolt is locked when turned clockwise. All of bolt is loosed when turned counter clockwise.



## 2-4. Oxygen and Air supply

MEKANT is the system to provide oxygen and air to patients. Oxygen and Air are provided through the adaptor end at the rear by pressure (25 ~ 90 PSI) and get mixed to have appropriate concentration.


No.	Figure	Explanation
a		<p>1. O<sub>2</sub> consists of green adaptor and green high pressure hose.</p> <p>2. Air input consists of a yellow color.</p>

Table 2-4. Connection way of high pressure hose



### Warning

Check the oxygen and air connector to be free from oil to prevent from explosion.

When you supply or stop oxygen and air, be careful when it's assembled and disassembled because high pressure hose can harm to human body due to inside pressure.



### Caution

Check and see if the Oxygen and Air pressure is set at exact pressure range prior to ventilator application.

When high pressure hose is used or moved, pay attention to wiring not to be tripped over persons.



## ۲-۴. Filters

Use ۲ filters between MEKANT Product and patient for the safety of patients. Before supplying for patient INS. Filter is used, and exhalation filter is used at the place for exhaled air from patient's body in order not to supply polluted air to patient.



### Warning

Contaminated gas might be in exposure to patients if the air was provided without applying filter.

Do not reuse the accessories when ☒ symbol is found since it represents its category as disposable.

Replace the disposable sensor with the new one after single use.

Reusable materials which are used for patients shall be sterilized and washed for reuse.



### Caution

Filters must be MEKICS genuine or designated products. If you wish to use other makes, please consult biomedical engineers at the hospitals, or local authorized dealer.

Please use quality certified product, otherwise manufacturer will not be responsible for any problems caused by this.

All accessories are aseptic; therefore the wrapping has to be intact before opening and use.



2-6. Transportation


No.	Figure	Explanation
①		①. Move the product by handle.

Table 2-6. Handle of the Cart



Warning

When moving the apparatus, please make sure that the LCD Monitor is fixed properly since this might cause damage to the Device.

Please move the product when components and power cords are in order.



Caution

Before moving the system, unlock the locker on the caster.



2-1. Fuse replacement instruction




Figure 2-5. AC INLET FUSE

No.	Name	Figure	Explanation
a	INLET FUSE	A diagram showing the steps to change the fuse. It includes a sequence of four images: 1. Pulling the upper part of the inlet toward you. 2. Removing the existing fuse. 3. Inserting the new fuse. 4. Closing the inlet by doing the steps backward.	<p>1. Pull the upper part of inlet toward you.</p> <p>2. Fuses are located where two upper arrows are and those are and those lower 2 arrows indicate the direction when inserted.</p> <p>3. Remove the existing fuse and insert the new fuse as shown in the picture after checking the rating.</p> <p>4. Close the inlet by doing the steps backward.</p>

Table 2-3. How to change fuse



<b>Section 3 : Operating instructions</b>	<b>35</b>
3-1. Front panel	35
3-1-1. Rotational Control Knob	35
3-1-2. Keypad	36
3-1-3. LED Indicators	37
3-2. Operating	37
3-2-1. Operating procedure and how to use	37
3-3. Display	40
3-3-1. Status information	41
3-3-2. Waveform and alarm status	42
3-3-3. Menu and setting information	43
3-3-4. TOUCH information	44
3-3-4-1. VENT TOUCH screen	45
3-3-4-2. Waveform screen	46
3-3-4-3. Menu screen	47
3-3-4-4. Setting screen	48
3-4. Start MV2000	49
3-4-1. VENT MODE	50
3-4-1-1. VACV	51
3-4-1-2. VSIMV	53
3-4-1-3. PACV	56
3-4-1-4. PSIMV	58
3-4-1-5. SPONT	61
3-4-1-6. tBILEVEL	63
3-4-1-7. NIV	66
3-4-1-8. PRVC	68
3-4-1-9. AUTO	71
3-4-2. Monitoring	72
3-4-3. Graphics	73
3-4-4. Tools	76
3-4-5. Events	77



2-4-6. SYSTEM .....	۷۸
2-4-6-1. SETUP .....	۷۸
2-4-6-1-1. BWF .....	۷۹
2-4-6-1-2. Nebu. Time .....	۷۹
2-4-6-1-3. O <sub>2</sub> SENSOR .....	۷۹
2-4-6-1-4. Trend Init. ....	۷۹
2-4-6-1-5. BTPS .....	۸۰
2-4-6-1-6. Sound Vol. ....	۸۰
2-4-6-1-7. Altitude .....	۸۰
2-4-6-1-8. Air Temp. ....	۸۰
2-4-6-2. PATIENT .....	۸۱
2-4-6-3. Date .....	۸۲
2-4-6-4. Calibration .....	۸۳
2-4-6-4-1. SELF CALIBRATION .....	۸۴
2-4-6-4-2. O <sub>2</sub> CALIBRATION .....	۸۴
2-4-6-4-3. INS-EXP CALIBRATION .....	۸۴
2-4-6-4-4. EtCO <sub>2</sub> CALIBRATION .....	۸۵
2-4-7. ALARM SET .....	۸۶



Section 3 : Operating instructions

3-1. Front panel

3-1-1. Rotational control knob

Rotating and Pressing the Encoder allows a user to navigate through and make changes to the display element.




Figure 3-1. Encoder



### ۳-۱-۲ Keypad




Figure ۳-1-2. Keypad

Part Number	Picture	Description
(a)		<b>Start Switch, Standby</b> ۱) If you press START switch for 0.5 sec, power would be on. In case of STOP OFF you want, switch should be pressed for longer than 10 seconds in Vent operating. ۲) Switch should be pressed for longer than 5 seconds in Vent Ready status. ۳) In case of STANDBY you want, switch should be pressed in Vent operating.
(b)		<b>MANUAL INSPIRATION:</b> This function which is for a respiration according to Ventilator's default is operated once by pressing the key.
(c)		<b>100% O<sub>2</sub> :</b> If you press this button, 100% O <sub>2</sub> is provided for 3 minutes.
(d)		<b>CLEAR :</b> It is used for clearance of menu when menu is controlled by Knob







Part Number	Picture	Description
e		<b>ACCEPT RUN</b> : It is used when Ventilator is working as its final value. If ACCEPT key is not pressed, setting value would not be accepted on the change of MODE and VENT SETTING. This is for the safe.
f		<b>ALARM 2 MIN</b> : This makes ALARM remove for 2 minutes.
g		<b>ALARM RESET</b> : This activates when alarm is reset. And previous alarm lamp is removed and the status of NORMAL will be.
j		<b>Encoder</b> : Rotary switch for selecting function on the pop-up menu.

Table 3-1-2. Keypad

## 3-1-3. LED Indicators



Part Number	Picture	Description
h		<b>AC POWER ON LED</b> : LED is on when machine is connected to AC, otherwise, LED is OFF.
i		<b>BATTERY ON LED</b> : Battery is charged as the machine is working when it is connected AC. If it is not, battery ON LED is lighted because the machine is operating by internal battery. If both AC POWER ON LED and BATTERY ON LED are OFF, battery should be charged.

Table 3-1-3. LED Indicators



3-2. Operating





Caution

Buzzer will be sounded when main power is on normally. Otherwise, check something wrong on the main power. (power, battery, fuse)

ALARM LAMP (Warning or INOP) is on and buzzer is sounded if ventilator has defect on working .

Connect the ventilator to the power is earthed.

3-2-1. Operating procedure and how to use

No.	Picture	Description
1		1) Connect high pressure hose and accessories.
2		1) Connect earthed cable to equipotential port and earthed AC power. 2) Connect power cord to an electric outlet and AC inlet. 3) Switch on POWER on the back side.





<p>3</p>		<p>1) Switch on START on the Front Panel. System Fail sound will be disappeared in case power is on normally.</p> <p>2) Wait for the calibration of each sensor. It is presented normal INITIAL, and then normal ventilator READY is appeared.</p>
<p>4</p>		<p>1) Set the present time and date.. (In an emergency, can be omitted.)</p> <p>2) Enter patient's bed number etc. (In an emergency, can be omitted.)</p> <p>3) In Standby Screen, Select patient's ages and input patient's height. Press START button.</p> <p>4) In top right corner, press Modes button and Choose mode suitable for the patient. Input other value as prescription. Adjust value by pressing and rotating.</p> <p>5) Select ACCEPT button.</p> <p>6) Adjust alarms properly.</p> <p>7) Check the wave form and figure on the monitor normally.</p>

Figure 3-2-1. Operating procedure and how to use



3-3. Display





Figure 3-3. LCD Display

Explanation		Reference
Measurement information	Measurement of ventilator or patient factor.	3-3-1
Waveform and alarm status	Graph and alarm statue of ventilator or patient factor.	3-3-2
Menu and Setting information	Menu button and Setting value according to respiratory mode.	3-3-3

Table 3-3. Description of display explanation



### 3-3-1. Status and measurement information

It displays measurement information of ventilator or patient factor. It shows different tables of measurement information to each Mode. By default,  $V_{E\text{ TIDAL}}$ ,  $V_{E\text{ MIN}}$ ,  $R_{\text{RESP}}$ ,  $P_{\text{PEAK}}$  and  $P_{\text{MEAN}}$  be displayed. Press  buttons in bottom left of the screen. It can be change different values.

Measurement items are displayed by the page as follows.

Page number	Displayed items
1	$P_{\text{PEAK}}$ , $P_{\text{PAUSE}}$ , $P_{\text{MEAN}}$
2	PEEP, AutoPEEP, $P_{\text{A.1}}$
3	$V_{\text{I TIDAL}}$ , $V_{\text{E TIDAL}}$ , $V_{\text{E MIN}}$
4	PEEPH, PEEPL, $F_{\text{PEAK}}$
5	$R_{\text{RESP}}$ , $R_{\text{RSPONT}}$ , $V_{\text{EMIN SPONT}}$
6	$T_{\text{I}}$ , $T_{\text{E}}$ , I : E
7	EXP.FLOW, $RSB_{\text{I}}$ , $O_2$
8	CL, RA
9	HR, $SpO_2$
10	$EtCO_2$ , Ins, $R_{\text{RESP}}$


Table 3-3-1. Displayed items by the page

Displayed status can be confirmed in Monitoring Menu overall. If you want more information, refer to [ 3-4-2. Monitoring].



### ㉓-㉓-㉓. Waveform and alarm status

It displays WAVE, ALARM STATUS and MESSAGE to ventilator or patient factor. Refer to Section ㉓ for alarm information.





3-3-3. Menu and Setting information

It displays respiratory mode of VENT Mode. It shows different table of setting table to each Mode.




Figure 3-3-3. Menu and Setting information Display explanation




Part Number	Display	Explanation
①	Menu	Refer to [3-3-4-3. Menu Screen]
②	Setting information	Refer to [3-3-4-4. SETTING TOUCH]
③		It displays AC Power. (It maybe changed Battery Status)
④		TOUCH ON/OFF
⑤		It displays ratio of I:E

Table 3-3-3. Setting information Display explanation



### 3-3-3. TOUCH SCREEN

MV2000(Touch) can be operated by screen touch. It displays the Touch Area on screen. It is possible to touch the area within outlined area.




Figure 3-3-3. TOUCH Screen information

Part Number	Display	Explanation	Ref.
(a)	VENT TOUCH	Select to change the vent mode.	3-3-4-1
(b)	Waveform	If you TOUCH waveform screen, the waveform can be converted to other waveforms.	3-3-4-2
(c)	Status information	Can be change different values. Refer to [3-3-1. Status information]	-
(d)	Menu	Select to change the menu setting with touch.	3-3-4-3
(e)	Setting screen	During vent mode operating, you can change the vent setting value.	3-3-4-4

Table 3-3-4. TOUCH Screen information explanation



### 3-3-4-1. VENT TOUCH SCREEN

You can go directly to the breathing mode settings in [VENT] TOUCH area




Figure 3-3-4-1(1). VENT MODE SETTING SCREEN

If you choose VENT MODE, the following screen will appear. After entering the setting, touch [ACCEPT] button, chosen mode will start. Refer to [3-4. MV2000 Start]




Figure 3-3-4-1(2). VENT MODE setting screen (V-ACV Mode)



3-3-2. Waveform screen

The border area is the waveform screen




Figure 3-3-2. Waveform screen information

If you TOUCH waveform screen, you can be change the measured graph.

Refer to [3-3-3. Graphics]



### 3-3-4. MENU SCREEN

As follow; TOUCH Menu screen or use encoder, for each menu screen appears.



Figure 3-3-4. MENU SCREEN information

For a description of each menu are as follows.

Part Number	Explanation	Ref.
Monitoring	Overall, the measured data shows for each mode.	3-4-2
Graphics	Can set Screen Layout and Trend Timing.	3-4-3
Tools	Use Insp. hold and Exp. hold.	3-4-4
Events	Can check Event occurred.	3-4-5
System	Setup preferences, patient information, date modified, calibration settings are configured.	3-4-6
Alarms	Set the alarms.	3-4-7

Table 3-3-4. MENU SCREEN explanation



### 3-3-3. SETTING screen

In breathing mode of operation, you can select the area of the SETTING screen (figure 3-3-4-3) and change the vent variable. If you re-click that variable, change the breathing variable that is applied immediately.

Monitoring	Graphics	Tools	Events	System	Alarms
V TIDAL mL 450 150 300	RRESP BPM 67 2 18	O2 % 100 21 21	PEEP cmH2O 45 0 3	T PAUSE sec 2.0 0.0 0.1	F PEAK LPM 180 5 40

Figure 3-3-3. SETTING TOUCH screen information (V-ACV mode)

Also, depending on the operating mode automatically changes SETTING screens.

VENT MODE	SETTING screen list
V-ACV	V Tidal, R RESP, T PAUSE, F PEAK, PEEP, O 2
V-SIMV	V Tidal, R RESP, T PAUSE, F PEAK, PEEP, O 2
P-ACV	P INSP, R RESP, T INSP, T RISE, PEEP, O 2
P-SIMV	P INSP, R RESP, T INSP, T RISE, PEEP, O 2
SPONT	P SUPP, T RISE, APNEA, T RIG, PEEP, O 2
tBiLEVEL	P SUPP, R RESP, TH PEEP, PEEP H, PEEP L, O 2
NIV	P INSP, R RESP, T INSP, T RIG, PEEP, O 2
PRVC	V Tidal, R RESP, T INSP, T RIG, PEEP, O 2
AUTO	P INSP, R RESP, T INSP, T RIG, PEEP, O 2

Table 3-3-3. SETTING screen explanation



3-4. Start MV 10 . .

To start MV2000, select the patient's condition and enter the patient's height in the Standby screen, and press START button to begin. Press Modes button in the top right corner to set the VENT mode. Also TOUCH ① area, move immediately to selected mode setting screen.




Figure 3-4. Standby screen

Explanation	
	Add a new patient or load Recently the patient setting.
	If the patient's height decisions, according to age and gender will automatically calculate IBW.(Ideal Body Weight).
	Move immediately to selected mode setting screen.
	If you enable Alarm AUTO, Start button is automatically set by the Alarm.
	Start MV2000.
	Power off MV2000.

Table 3-4. Standby screen explanation



### 3-4-1. VENT MODE

MV2000 has the ventilation mode which are SIMV(VSIMV,PSIMV), ACV(VACV, PACV, other company calls CMV to Mandatory Mode), Spontaneous Mode, tBi-Level, NIV, PRVC and AUTO Mode.




Figure 3-4-1. VENT MODE menu

VENT MODE	Operating MODE	Explanation	Ref.
V-ACV	V-ACV	Volume based Assistance Controlled Ventilation It controls Patient Ventilation by force from air volume. And it controls expiration time and inspiration time by force.	3-4-1-1
V-SIMV	V-SIMV	Volume based Synchronized Intermittent Mandatory Ventilation. It controls Patient Ventilation to air volume from Synchronized Intermittent Mandatory Ventilation	3-4-1-2
P-ACV	P-ACV	Pressure based Assistance Controlled Ventilation It controls Patient Ventilation by force from air pressure.	3-4-1-3
P-SIMV	P-SIMV	Pressure based Synchronized Intermittent Mandatory It controls Patient Ventilation to air pressure from Synchronized Intermittent Mandatory Ventilation.	3-4-1-4
SPONT	SPONT	Spontaneous Ventilation. It supports spontaneous breathing.	3-4-1-5
tBiLEVEL	tBiLEVEL	It is same to SPONT. And it has Upper PEEP. (Refer to SPONT)	3-4-1-6
NIV	NIV	Non Invasive Ventilation It is mask type, not tubing type. It work P-ACV allow to leakage. (Refer to P-ACV)	3-4-1-7
PRVC	PRVC	Pressure Regulated Volume Controlled Ventilation Breath are delivered mandatorily with a constant flow to assure present volumes.	3-4-1-8
AUTO	P-ACV	Quick Start according to the weight of patient. (P-ACV Controlled)	3-4-1-9

Figure 3-4-1. VENT MODE menu explanation



### ۳-۴-۱-۱. V<sub>ACV</sub>

V<sub>ACV</sub> (Volume based Assistance Controlled Ventilation) :

It controls Patient Ventilation by force from tidal volume. And it controls expiration time and inspiration time by force. Vi<sub>TIDAL</sub> means inspiration volume to supply patient.




Figure ۳-۴-۱-۱. V-ACV SETTING screen



Item	Setting Range	Unit	Step	Explanation
V <sub>TIDAL</sub>	50 ~ 2500	mL	50	It decides volume of 1 breath. <b>NOTE :</b> It set the volume of patient weight. It is 10mL per 1kg. (You can set the volume per 1kg in [SYSTEM - SETUP - BWF] Menu.)
R <sub>RESP</sub>	2 ~ 120	BPM	1	It set the respiratory rate per 1 minute.
T <sub>PAUSE</sub>	0.0 ~ 2.0	sec	0.1	Inspiration Pause Time. It set the inspiration pause time.
P <sub>EEP</sub>	0 ~ 45	cmH <sub>2</sub> O	1	It set the PEEP pressure.( Base pressure of patient)
O <sub>2</sub>	21 ~ 100	%	1	It set the percent of O <sub>2</sub> to patient.
F <sub>PEAK</sub>	50 ~ 110	LPM	1	It set the flow-peak (maximum velocity) to patient.
F <sub>END</sub>	25 ~ 100	%	50	It set percent of Flow (supply velocity) / Fpeak (supply maximum velocity) to End-tidal.
E <sub>nSENS</sub>	10 ~ 60	%	10	It set the sensitivity of one breath completion of patient during to setting of percent of patient expiration volume.
TRIGGER TYPE	PRESSURE /FLOW	-	-	It sets the type of trigger. It is Pressure or Flow type.
F(P) <sub>TRIG</sub>	0.5 ~ 2.0	cmH <sub>2</sub> O /LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
S <sub>IGH</sub> MODE	OFF, 30, 60, 90, 120	RPB	30	It sets breathing rate timing in sigh mode. Don't support over 1500mL volume.
ACCEPT	-	-	-	Execute V-ACV Mode.
CANCEL	-	-	-	Go back to previous menu

Table 3-4-1-1. V-ACV SETTING screen explanation

**A/C (Volume Control) - V-ACV**

If MEKANT perceives spontaneous breathing it works the Assisted Mandatory Breath. If not it works Controlled Mandatory Breath.

When MEKANT works in Volume Control Ventilation it decides inspiratory time for setting of tidal volume and peak flow. Also inspiratory time is influenced by flow pattern and plateau.

In the case of I:E ratio decides it for factor of inspiratory time and respiratory rate.

When you change the factor of I:E ratio and inspiratory time, you can know cycle time, inspiratory time, expiratory time and I:E ratio through breath timing bar.



### 3-4-1-2. V<sub>SIMV</sub>

V<sub>SIMV</sub> (Volume based Synchronized Intermittent Mandatory Ventilation) :

It controls Patient Ventilation to tidal volume from Synchronized Intermittent Mandatory Ventilation. V<sub>TIDAL</sub> means inspiration volume to supply patient.



#### Warning

MEKANT synchronized spontaneous breathing or in case of without spontaneous breathing signal MEKANT ventilate by pressure.




Figure 3-4-1-2. V-SIMV SETTING screen



Item	Setting Range	Unit	Step	Explanation
V <sub>TIDAL</sub>	50 ~ 2500	mL	5	It decides volume of 1 breath. <b>NOTE :</b> It set the volume of patient weight. It is 10mL per 1kg. (You can set the volume per 1kg in [SYSTEM - SETUP - BWF] Menu.)
R <sub>RESP</sub>	2 ~ 120	BPM	1	It sets the respiratory rate per 1 minute.
T <sub>PAUSE</sub>	0.0 ~ 2.0	Sec	0.1	Inspiration Pause Time It set the inspiration time.
P <sub>EEP</sub>	0 ~ 45	cmH <sub>2</sub> O	1	It set the PEEP pressure. ( Base pressure of patient)
O <sub>r</sub>	21 ~ 100	%	1	It sets the percent of O <sub>r</sub> to patient.
F <sub>PEAK</sub>	50 ~ 180	LPM	1	It sets the flow-peak (maximum velocity) to patient.
F <sub>END</sub>	25 ~ 100	%	5	It sets percent of Flow (supply velocity) / Fpeak (supply maximum velocity) to End-tidal.
E <sub>nSENS</sub>	10 ~ 60	%	10	It sets the sensitivity of one breath completion of patient during to setting of percent of patient expiration volume.
E <sub>xSENS</sub>	10 ~ 60	%	10	It sets the sensitivity of one breath completion of patient during to setting of percent of patient inspiration volume.
TRIGGER TYPE	PRESSURE /FLOW	-	-	It sets the type of trigger. It is Pressure or Flow type.
F(P) <sub>TRIG</sub>	0.5 ~ 2.0	cmH <sub>2</sub> O /LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
P <sub>SUPP</sub>	0 ~ 45	cmH <sub>2</sub> O	1	It sets the pressure to patient during spontaneous breathing.
ACCEPT	-	-	-	Execute V <sub>SIMV</sub> Mode.
CANCEL	-	-	-	Go back to previous menu


Table 3-4-1-2. V-SIMV SETTING screen explanation



### **SIMV (Volume Control or Pressure Control) only Mandatory Breath**

SIMV is a mixed ventilator mode that allows both mandatory and spontaneous breaths. SIMV mode differs to other company. SIMV of MEKANT guarantees one of mandatory breath per SIMV cycle. This mandatory breath works by patient initiated mandatory (PIM or Assisted Mandatory) breath or ventilator initiated mandatory (VIM or Controlled Mandatory).

As the below figure shows, each SIMV breath cycle ( $T_b$ ) has two parts. One is the mandatory interval ( $T_m$ ) and the other is Spontaneous Interval ( $T_s$ ). PIM is delivered, the  $T_m$  interval ends and the ventilator switches to the spontaneous interval ( $T_s$ ). If a PIM is not delivered, the ventilator delivers a VIM at the end of the mandatory interval, and then switches to the spontaneous interval.



In the SIMV mode of MEKANT, SIMV Breath Cycle is decided by respiratory rate of setting. If respiratory rate is 12 times, SIMV Breath Cycle is 5 seconds. Mandatory Interval is defined as whichever is less  $0.7 \times \text{SIMV Breath Cycle}$ , or 1.0 seconds. Mandatory interval can be shorter according to PIM breath.



### 3-4-1-3. P<sub>ACV</sub>

P<sub>ACV</sub>(Pressure based Assistance Controlled Ventilation) :

It controls Patient Ventilation by force from air pressure.




Figure 3-4-1-3. P-ACV SETTING screen

Item	Setting Range	Unit	Step	Explanation
P <sub>INSP</sub>	5 ~ 15	cmH <sub>2</sub> O	1	It sets the pressure to patient during inspiration.
T <sub>INSP</sub>	0.2 ~ 9.9	sec	0.1	It sets the inspiration time.
R <sub>RESP</sub>	2 ~ 120	BPM	1	It sets the respiratory rate per 1 minute.
P <sub>EEP</sub>	0 ~ 45	cmH <sub>2</sub> O	1	It sets the PEEP pressure. (Base pressure of patient)
O <sub>2</sub>	21 ~ 100	%	1	It sets the percent of O <sub>2</sub> to patient.
T <sub>RISE</sub>	0.1 ~ 0.5	sec	0.1	It sets the flow speed of gas to patient.
E <sub>nSENS</sub>	10 ~ 60	%	10	It sets the sensitivity of one breath completion of patient during to setting of percent of patient expiration volume.
TRIGGER TYPE	PRESSURE /FLOW	-	-	It sets the type of trigger. It is Pressure or Flow type.
F(P) <sub>TRIG</sub>	0.5 ~ 20	cmH <sub>2</sub> O /LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
ACCEPT	-	-	-	Execute PACV Mode
CANCEL	-	-	-	Go back to previous menu


Table 3-10. P<sub>ACV</sub> Menu Explanation



**A/C (Pressure Control) PACV**


When MEKANT provide pressure control ventilation of mandatory breath, it draws inspiratory pressure instead of tidal volume and it sets I:E ratio or inspiratory time instead of peak flow same as other company.

But as the below figure shows, it is differ to slope of increase – baseline (PEEP) to Inspiratory pressure – according to setting of flow acceleration percent (FAP). According to an increment of FAP setting value, slope of pressure is rapid. If not slope of pressure is gentle. If you show slope of flow, you know that FAP is the change of flow volume to patient. FAP sets PCV and PSV each other.



MEKANT keeps on pressure regularly during inspiration same as other company. So it has Ramp shape of flow pattern and cannot use plateau function.

MEKANT is more activity and elaborateness for control of airway pressure because it has exhalation valve and active exhalation valve. So it controls pressure to exhalation valve in patient's talk, cough. And it allows spontaneous breath. It is special function only MEKANT. It explain figure of next page.





### 3-4-1-4. P<sub>SIMV</sub>

**P<sub>SIMV</sub>**(Pressure based Synchronized Intermittent Mandatory Ventilation) :

This mode is that MEKANT control artificial respiration by spontaneous breath and pressure.



#### Warning

MEKANT synchronized spontaneous breathing or in case of without spontaneous breathing signal MEKANT ventilate by pressure.




Figure 3-4-1-4. P-SIMV SETTING screen



Item	Setting Range	Unit	Step	Explanation
P <sub>INSP</sub>	5 ~ 15	cmH <sub>2</sub> O	1	Setting the patient pressure to once breath.
T <sub>INSP</sub>	0.2 ~ 9.9	sec	0.1	Setting the inspiration Time.
R <sub>RESP</sub>	2 ~ 120	BPM	1	Setting number of breath per minute.
P <sub>EEP</sub>	0 ~ 45	cmH <sub>2</sub> O	1	Setting PEEP pressure. (Base Pressure)
O <sub>2</sub>	21 ~ 100	%	1	Setting oxygen density.
TRISE	0.1 ~ 0.5	sec	0.1	It sets the flow speed of gas to patient
E <sub>nSENS</sub>	10 ~ 60	%	10	Setting sensitivity of expiration by ratio of expiration capacity.
E <sub>xSENS</sub>	10 ~ 60	%	10	Setting sensitivity of expiration by ratio of inspiration capacity.
TRIGGER type	PRESSURE /FLOW	cmH <sub>2</sub> O /LPM	-	It set the type of trigger. It is Pressure or Flow type.
TRIGGER	0.5 ~ 20	cmH <sub>2</sub> O /LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
PPS	0 ~ 45	cmH <sub>2</sub> O	1	Setting pressure of patient when it is working spontaneous breathing.
ACCEPT	-	-	-	Execute P-SIMV Mode.
CANCEL	-	-	-	Go back to previous menu


Table 11-11. P<sub>SIMV</sub> Menu Explanation



### **SIMV (Volume Control or Pressure Control) only Mandatory Breath**

SIMV mode is able to work Mandatory Breath and Spontaneous Breath. MEKANT's SIMV mode is difference with other equipment. MEKANT's SIMV mode guarantees at least one Mandatory Breath per each cycle. This Mandatory Breath is one of the alternative Patient Initiated Mandatory (PIM or Assisted Mandatory) Breath or Ventilator Initiated Mandatory (VIM or Controlled Mandatory) Breath.

As bellow figure each SIMV Breath Cycle has two parts one is Mandatory Interval ( $T_m$ ) and other one is Spontaneous Interval ( $T_s$ ). After receive PIM Breath signal the Mandatory Interval was terminated and then to start Spontaneous Interval ( $T_s$ ). If PIM Breath does not happen until Mandatory Interval terminate It will finish Mandatory Interval and then change to Spontaneous Interval with VIM Breath.



MEKANT SIMV mode's SIMV Breath Cycle depend on number of breath. For example if set Respiratory Rate is 12 SIMV Breath Cycle is 5 sec. Mandatory Interval is taken smaller value maximum 1.0 sec or  $0.7 \times \text{SIMV Breath Cycle}$ . Mandatory Interval will be shorter by PIM Breath.



### ۳-۴-۱-۵. SPONT

**SPONT**(Spontaneous Ventilation) :  
It supports spontaneous breathing.



#### Warning

When MEKANT is working SPONT mode to apnea status during setting time, it change mode(V-ACV) of Apnea Backup setting automatically. It does change back SPONT mode, when 3 times spontaneous breaths within 10 seconds in apnea backup mode.




Figure ۳-۴-۱-۵. SPONT SETTING screen



Item	Setting Range	Unit	STEP	Explanation
P <sub>EEP</sub>	0 ~ 40	cmH <sub>2</sub> O	1	It sets the PEEP pressure. (Base pressure of patient)
O <sub>r</sub>	21 ~ 100	%	1	It sets the percent of O <sub>r</sub> to patient.
E <sub>xSENS</sub>	10 ~ 60	%	10	It sets the sensitivity of one breath completion of patient during to setting of percent of patient inspiration volume.
TRIGGER TYPE	PRESSURE /FLOW	-	-	It sets the type of trigger. It is Pressure or Flow type.
F(P) <sub>TRIG</sub>	0.5 ~ 20	cmH <sub>2</sub> O /LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
P <sub>SUPP</sub>	0 ~ 40	cmH <sub>2</sub> O	1	It sets the pressure to patient during spontaneous breathing.
APNEA	OFF ~ 60	sec	10	If patient has not breath of setting time it works Apnea mode. It does not detect breath of patient during "OFF"
ACCEPT	-	-	-	Execute SPONT Mode.
CANCEL	-	-	-	Go back to previous menu

Table 9-9. SPONT Menu Explanation

**SPONT mode only**

If spontaneous breath reaches to sensitivity, it makes pressure between setting of sensitivity and 1 cmH<sub>2</sub>O through flow supply according to patient flow demand. If Pressure is bigger than 1 cmH<sub>2</sub>O, it starts expiration to opening Exhalation valve. If patient has apnea breath, MEKANT occurs Apnea alarm and start Apnea ventilation (VACV Mode).

It works setting of Respiratory Rate, Tidal Volume, Inspiratory Pressure, Peak Flow, Inspiratory Time, FiO<sub>2</sub> and so on.

It works same as pressure support ventilation. Also if patient has apnea breath, MEKANT occurs Apnea alarm and start Apnea ventilation. But it change previous spontaneous mode if patient has 2 times spontaneous breaths within 10 seconds.



### 3-4-1-6. tBiLEVEL

**tBiLEVEL:** It works same as SPONT mode and it has Upper PEEP.




Figure 3-4-1-6. tBiLEVEL SETTING screen

Item	Setting Range	Unit	STEP	Explanation
T <sub>H</sub> PEEP	0.2 ~ 24	sec	0.1	Setting High PEEP Time. (Upper pressure holding time).
R <sub>RESP</sub>	2 ~ 30	BPM	1	Setting patient's number of breath.
PEEP <sub>HIGH</sub>	0 ~ 40	cmH <sub>2</sub> O	1	Setting High PEEP. (Upper pressure).
PEEP <sub>LOW</sub>	0 ~ 40	cmH <sub>2</sub> O	1	Setting Low PEEP. (Base Pressure of patient).
O <sub>2</sub>	21 ~ 100	%	1	Setting oxygen density of air
E <sub>xSENS</sub>	10 ~ 60	%	1.0	Setting sensitivity of expiration by ratio of inspiration capacity.
TRIGGER TYPE	PRESSURE / FLOW	-	-	It sets the type of trigger. It is Pressure or Flow type.
F(P) <sub>TRIG</sub>	0.5 ~ 2.0	cmH <sub>2</sub> O / LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
P <sub>SUPP</sub>	0 ~ 40	cmH <sub>2</sub> O	1	Setting pressure of patient when it is working spontaneous breathing.
ACCEPT	-	-	-	Execute tBiLEVEL Mode.
CANCEL	-	-	-	Go back to previous menu

Table 3-4-1-6. tBiLEVEL SETTING Explanation




**tBILEVEL Mandatory Breath**

MEKANT's tBILEVEL mode is complex mode it may happened Mandatory Breath and Spontaneous Breath. Because it works as SIMV mode but it difference because Spontaneous range not fixed. So, anytime it is able to Spontaneous Breath.


tBILEVEL has two level of PEEP to supply patient. The upper PEEP is PEEP<sub>HIGH</sub> and lower PEEP is PEEP<sub>LOW</sub>. Mandatory Breath's Inspiration is from PEEP<sub>LOW</sub> to PEEP<sub>HIGH</sub> and Expiration is from PEEP<sub>HIGH</sub> to PEEP<sub>LOW</sub>. So tBILEVEL mode's Mandatory Breath is Pressure Control Ventilation. Each change time of before and after small zone is changing zone, this zone purpose is to detect patient breath if it detect patient breath it would adjust little early or late to start change inspiration and expiration for synchronizing with patient breathing.

Setting parameters are Respiratory Rate, PEEP<sub>HIGH</sub>, PEEP<sub>LOW</sub> and so on, In case of Pressure Support Ventilation there is more add Support Pressure, FAP, Essens and so on.



MEKANT's tBILEVEL mode supplies 5 cmH<sub>2</sub>O Pressure in Spontaneous Breath. If user want to more high pressure support it will use Pressure Support Ventilation setting. If Support Pressure + PEEP<sub>L</sub> value is higher than PEEP<sub>H</sub> It will be working pressure support ventilation difference than Spontaneous breath works Pressure Support Ventilation. Refer to the figure below.





tBILEVEL mode is able to work Spontaneous Breath any time it reduce fighting between ventilator and patient. It can reduce sedation so patient gets over more sort

All of Spontaneous Breath supplies Pressure Support Ventilation so patient feel more comfortable. If it using opposite I:E ratio it means very sort Expiratory Time it is Airway Pressure Release Ventilation(APRV) status. This APRV is one of best way to apply ARDS (Adult Respiratory Distress Syndrome).





### 3-4-1-7. NIV

NIV : Mask type ventilation, it works similar P-ACV mode with air leakage.




Figure 3-4-1-7. NIV SETTING screen

Item	Setting Range	Unit	STEP	Explanation
P <sub>INS</sub>	0 ~ 10	cmH <sub>2</sub> O	1	Setting the patient pressure to once breath.
T <sub>INS</sub>	0.2 ~ 9.9	sec	0.1	Setting the inspiration Time
R <sub>RESP</sub>	2 ~ 106	BPM	1	Setting number of breath per minute.
P <sub>EEP</sub>	0 ~ 40	cmH <sub>2</sub> O	1	Setting PEEP pressure. (Base Pressure)
O <sub>2</sub>	21 ~ 100	%	1	Setting oxygen density.
TRIGGER TYPE	PRESSURE / FLOW	-	-	It set the type of trigger. It is Pressure or Flow type.
F(P) <sub>TRIG</sub>	0.5 ~ 2.0	cmH <sub>2</sub> O / LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
ACCEPT	-	-	-	Execute NIV Mode.
CANCEL	-	-	-	Go back to previous menu

Figure 3-4-1-7. NIV SETTING screen explanation




**■ NIV** (None Invasive Ventilation)

Mask type ventilation, it works similar P-ACV mode with air leakage. If air leak alarm occurs frequently, adjust inspiration and expiration leak volume.

It can set [Alarm] => [AIR LEAK] menu.

NIV mode need to prepare unexpected leak air, at the same time, must supply a stable pressure. During exhausted time, supply airway need to be feedback quickly.

By maintain flow trigger performance, it can be make stable PSV.





### ㄣ-ㄝ-ㄧ-ㄨ. PRVC

**PRVC** (Pressure Regulated Volume Controlled Ventilation) :

Breaths are delivered mandatorily to assure present volumes, with a constant inspiratory pressure continuously adapting to the patient's condition. The flow pattern is decelerating.




Figure ㄣ-ㄝ-ㄧ-ㄨ. PRVC SETTING screen




Item	Setting Range	Unit	STEP	Explanation
V <sub>TIDAL</sub>	50 ~ 2500	ml	5	It decides volume of 1 breath. <b>NOTE :</b> It set the volume of patient weight. It is 10mL per 1kg. (You can set the volume per 1kg in [SYSTEM - SETUP - BWF] Menu.)
T <sub>INSP</sub>	0.2 ~ 9.9	sec	0.1	Setting the inspiration Time
R <sub>RESP</sub>	2 ~ 120	BPM	1	It sets the respiratory rate per 1 minute.
P <sub>EEP</sub>	0 ~ 45	cmH <sub>2</sub> O	1	It sets the PEEP pressure.( Base pressure of patient)
O <sub>2</sub>	21 ~ 100	%	1	It sets the percent of O <sub>2</sub> to patient.
TRISE	0.1 ~ 0.5	sec	0.1	It sets the flow speed of gas to patient
E <sub>nSENS</sub>	10 ~ 60	%	10	It sets the sensitivity of one breath completion of patient during to setting of percent of patient expiration volume.
TRIGGER TYPE	PRESSURE /FLOW	-	-	It sets the type of trigger. It is Pressure or Flow type.
F(P) <sub>TRIG</sub>	0.5 ~ 20	cmH <sub>2</sub> O /LPM	0.5	Setting to detect spontaneous breathing signal, There are pressure and Flow type. User can select one.
P <sub>RESSURE LIMIT</sub>	2 ~ 10	cmH <sub>2</sub> O	20	It sets the patient pressure limit to set over pressure.
ACCEPT	-	-	-	Execute PRVC Mode
CANCEL	-	-	-	Go back to previous menu

Table 3-4-1-8. PRVC SETTING screen explanation




**PRVC (Pressure Regulated Volume Controlled Ventilation)**

Despite of the advantage of the pressure mode, a long time or a sudden change in status of the patient, or pure breathing changes for a long time, the PRVC mode resolve this instability mode.



Initial Operation of PRVC ( 1<sup>st</sup> Breathing )



Pressure Regulation Operation at Compliance change of Patient



## 3-4-1-9. AUTO

**AUTO :**

Emergency Quick Start function it requirement patient body weight. It works PACV Mode. This menu purpose is that in case of emergency patient user have to start quickly so it has minimum stability. All parameters are setup default and user can execute with only body weight. When "Auto Mode", on the top left corner, [AUTO] is flashing.




Figure 3-4-1-9. AUTO SETTING screen

Item	Setting Range	Unit	STEP	Explanation
ACCEPT	-	-	-	Execute AUTO Mode
CANCEL	-	-	-	Go back to previous menu

Table 3-4-1-9. AUTO SETTING screen explanation



3-4-2. Monitoring

Overall, the measured data shows for each mode.

Table 1 : Display the contents about ventilator.

Table 2 : Display HR, SpO<sub>2</sub>, EtCO<sub>2</sub>, Ins, R<sub>RESP</sub>.




Figure 3-4-2-(1). Monitoring screen(Table 1)




Figure 3-4-2-(2). Monitoring screen(Table 2)



3-4-3. Graphics

By using Graphics menu, you can change each Screen Layout and Trend Timing. Refer to [Figure 3-4-3(1)].





Figure 3-4-3(1). Graphics screen

You can change the update period of Trend screen.  
When set to [1 RESP], Trend data will be updated each respiratory.



When Layout select, the screen arrangement will be change as follow.

Item	Explanation
	<p><b>Layout 1</b></p> <p>Pressure, Flow, Volume graphs are displayed.</p>
	<p><b>Layout 2</b></p> <p>Two graphs among Pressure, Flow, Volume are displayed.</p> <p>Two graphs among P-V, V-F, P-F are displayed.</p>
	<p><b>Layout 3</b></p> <p>Two graphs among Pressure, Flow, Volume are displayed.</p> <p>One graph among P-V, V-F, P-F are displayed.</p> <p>Two trend graphs are displayed.</p>
	<p><b>Layout 4</b></p> <p>Two graphs among Pressure, Flow, Volume are displayed.</p> <p>Four trend graphs are displayed.</p>



The method of selecting a graph, just touch the screen.  
The TOUCH areas are [Figure 3-4-3(2)] as follow;




Figure 3-4-3(2). Graphics TOUCH areas

Item	Explanation
Ⓐ	Waveform graph display. When waveform TOUCH, select menu is displayed.
Ⓑ	Loop graph display. When Loop graph TOUCH, select menu is displayed.
Ⓒ	Trend graph display. When Trend graph TOUCH, select menu is displayed.

Table 3-4. Graphics menu explanation



3-4-4. Tools

This menu is to measure [Inspiration hold time] and [Expiration hold time].




Figure 3-4-4. Tools screen

Item	Explanation
Insp. hold	It measures for 1 second the compliance and resistance.
Exp. hold	After EXPIRATION 1 second holding, it measures the AUTO PEEP.

Table 3-4-4. Tools Explanation



3-4-4. Events

When an event occurs, Events screen can be found at.




Figure 3-4-4. Events screen

Item	Explanation
Setting	When alarm setting change, events will be displayed.
Alarm	When alarms appear, events will be displayed.
ALL	All events will be displayed.

Table 3-4-4. Events explanation



### 3-4-6. SYSTEM

SYSTEM MENU is made up SETUP menu, PATIENT menu, DATE menu and CALIBRATION menu.




Figure 3-4-6. SETUP MENU screen configuration

#### 3-4-6-1. SETUP

It is setting for user. There are BWF, Nebulizer time, O<sub>2</sub> sensor ON/OFF, Initial Trend, BTPS ON/OFF, Sound Volume, Altitude, Air temperature.

Item	Explanation	Ref.
BWF	Setting body weight of patient in mL/kg.	3-4-6-1-1
Nebu. Time	Setting operating time of nebulizer( 1 ~ 140 minutes)	3-4-6-1-2
O <sub>2</sub> Sensor	Set to use the O <sub>2</sub> Sensor.	3-4-6-1-3
Trend Init	Initialize Trend graph and Trend data.	3-4-6-1-4
BTPS	BTPS ON/OFF switch.	3-4-6-1-5
Sound Vol.	Adjust sound volume.	3-4-6-1-6
Altitude	Set the altitude.	3-4-6-1-7
Air Temp	Set the air temperature.	3-4-6-1-8

Table 3-4-6-1. SETUP MENU screen explanation



### 3-4-6-1-1. BWF

Set body weight of patient in mL/kg.

Item	Setting Range	Unit	STEP	Explanation
BWF	0 ~ 10	mL/kg	1	Volume / Body weight

Table 3-4-6-1-1. B/W FACTOR explanation

### 3-4-6-1-2. Nebu. Time

Set nebulizer operating time in minute.

Item	Setting Range	Unit	STEP	Explanation
Nebulizer Time	10 ~ 140	min	10	Setting operating time of nebulizer.

Table 3-4-6-1-2. Nebu. Time explanation

### 3-4-6-1-3. O<sub>2</sub> SENSOR

Set to use the O<sub>2</sub> Sensor ON/OFF.

- ON : It is setting O<sub>2</sub> sensor is normal.
- OFF : Set O<sub>2</sub> sensor is abnormal or O<sub>2</sub> Gas is empty.

Item	Setting Range	Unit	STEP	Explanation
O <sub>2</sub> SENSOR	ON / OFF	-	-	Set to use the O <sub>2</sub> Sensor ON/OFF.

Table 3-4-6-1-3. O<sub>2</sub> SENSOR explanation

### 3-4-6-1-4. Trend Init.

Initialize Trend graph and Trend data.



### 3-4-6-1-4. BTPS

BTPS is "Body Temperature and Pressure Saturated with Water Vapor".

In respiratory physiology lung volumes and flows are standardized to barometric pressure at sea level.

Item	Setting Range	Unit	STEP	Explanation
BTPS	ON / OFF	-	-	BTPS ON/OFF switch.

Table 3-4-6-1-4. BTPS explanation

### 3-4-6-1-5. Sound Vol.

Adjust sound and beep volume.

Item	Setting Range	Unit	STEP	Explanation
Sound Vol.	10 ~ 100	%	10	Set alarm sound volume.

Table 3-4-6-1-5. Sound Vol. explanation

### 3-4-6-1-6. Altitude

Set altitude height in meter.

Item	Setting Range	Unit	STEP	Explanation
ALTITUDE	0 ~ 1000	METER	20	Set height of altitude.

Table 3-4-6-1-6. Altitude explanation

### 3-4-6-1-7. Air Temp.

Set air temperature in centigrade.

Item	Setting Range	Unit	STEP	Explanation
Air Temp	0 ~ 40	℃	1	Set air temperature of centigrade.

Table 3-4-6-1-7. Air Temperature explanation



3-4-6-2. PATIENT

Input the patient's information.




Figure 3-4-6-2. PATIENT screen

Item	Mode	Explanation
BED NUM.	-	Set BED Number
ID	-	Set the patient's ID number.
IP	-	External communication with a LAN connection, Set the IP address number.
PORT	-	External communication with a LAN connection, Set the IP PORT number.
PassWord	-	Enter the password(77) for IP and PORT.

Table 3-4-6-2. Patient screen explanation



3-4-6-3. Date

Set date and time.

There are main and pneumatic SW version information, Ventilator operating time.




Figure 3-4-6-3. **Date** screen



### 3-4-6-4. Calibration

You can perform the Calibration of Ventilator, O<sub>2</sub> sensor, Ins-Exp, EtCO<sub>2</sub>.




Figure 3-4-6-4. Calibration screen

Item	Explanation	Ref.
SELF	It is calibration mode for flow sensors.	3-4-6-4-1
O <sub>2</sub>	Calibration O <sub>2</sub> Sensor	3-4-6-4-2
INS-EXP	It is calibration mode for striking balance between inspiration flow sensor and expiration flow sensor.	3-4-6-4-3
EtCO <sub>2</sub>	Calibration EtCO <sub>2</sub> Sensor	3-4-6-4-4

Table 3-4-6-4. Calibration explanation



**3-4-6-1. SELF CALIBRATION**

It is calibration mode for flow sensors. It has 2 flow sensor calibrations. One is offset calibration of Inspiration flow sensor and the other is offset calibration of expiration flow sensor.

CAUTION : Before calibration, MV 00 must be separated from patient.

Item	Range	Explanation
SELF CALIBRATION	START / OFF	Inspiration and expiration flow sensor offset calibration.

Table 3-4-6-1. SELF CALIBRATION explanation

**3-4-6-2. O<sub>2</sub> CALIBRATION**

Calibrate O<sub>2</sub> Sensor. Push the encoder button and O<sub>2</sub> calibration start. After few seconds O<sub>2</sub> calibration is stopped. Measured O<sub>2</sub> sensor calibration data is saved in memory.

Item	Range	Explanation
O <sub>2</sub> CALIBRATION	START / OFF	Calibration O <sub>2</sub> Sensor

Table 3-4-6-2. O<sub>2</sub> CALIBRATION explanation**3-4-6-3. INS-EXP CALIBRATION**

It is calibration mode for striking balance between inspiration flow sensor and expiration flow sensor.

Item	Range	Explanation
INS-EXP CALIBRATION	START / OFF	It is calibration mode for striking balance between inspiration flow sensor and expiration flow sensor. Turn the encoders to the left and right, Ins Volume and Exp Volume is adjusted to match.

Table 3-4-6-3. INS-EXP CALIBRATION explanation

Turn the encoders to the left and right, you should be adjusted to match the measured data and setting data.

For more accurate INS-EXP CALIBRATION, operated the ventilator for 2 hours and then INS-EXP CALIBRATION can get better results.

The important things to note that [EXP. GAS] setting is set to DRY GAS state, and then start INS-EXP CALIBRATION. If [EXP. GAS] be set HUMID GAS, less volume is measured.

BF means Base Flow. BF is the difference of inspiratory flow and expiratory flow.



**3-4-6-4. EtCO<sub>r</sub> CALIBRATION**

If you have EtCO<sub>2</sub> module, Calibration is performed.

Item	Range	Explanation
EtCO <sub>r</sub> CALIBRATION	START / OFF	CALIBRATION EtCO <sub>r</sub>

Table 3-4-6-4. EtCO<sub>2</sub> CALIBRATION explanation



### 3-4-1. ALARM SET

To adjust alarm levels for each measured value, Turn the encoders to the left and right, and then set the value.




Figure 3-4-1. Alarm setting screen

Item	Setting Range	Unit	STEP	Explanation
V <sub>E</sub> TIDAL HIGH	20 ~ 200, OFF	mL	10	Upper limit of patient's capacity.
V <sub>E</sub> TIDAL LOW	OFF, 0 ~ 200	mL	10	Lower limit of patient's capacity.
V <sub>E</sub> MIN HIGH	0.1 ~ 50.0	LPM	0.1	Upper limit of patient's breath of capacity.
V <sub>E</sub> MIN LOW	0.0 ~ 49.9	LPM	0.1	Lower limit of patient's breath of capacity.
R <sub>RESP</sub> HIGH	2 ~ 100	BPM	1	Upper limit of number of breath per minute.
R <sub>RESP</sub> LOW	2 ~ 100	BPM	1	Lower limit of number of breath per minute.
P <sub>PEAK</sub> HIGH	0 ~ 120, OFF	cmH <sub>2</sub> O	1	Upper limit of patient's pressure.
P <sub>PEAK</sub> LOW	OFF, 0 ~ 120	cmH <sub>2</sub> O	1	Lower limit of patient's pressure.
O <sub>2</sub> HIGH	0 ~ 100, OFF	%	1	Upper limit of oxygen's density.
O <sub>2</sub> LOW	OFF, 0 ~ 100	%	1	Lower limit of oxygen's density.
AIR LEAK	0 ~ 500	mL	10	Difference breath of capacity between Ins. Volume and Exp. Volume.

Table 3-4-1. ALARM screen explanation



## Section 4

<b>Section ㉞ : Alarms</b>	㉞㉞
㉞-㉞. General Information	㉞㉞
㉞-㉞. Alarm Message	㉞㉞

## Section ㉞ : Alarms

## ㉞-㉞. General Information

All alarm indicates condition of each situation on LED display window. LCD monitor shows their warning message.




Figure ㉞-㉞. Alarm screen



## ㉞-㉟. Alarm Message

If alarm message is indicated, check the current mode and solve the problem according to it.

Item	Alarm range
PRESSURE HIGH	5 cmH <sub>2</sub> O ~ 120 cmH <sub>2</sub> O
PRESSURE LOW	0 cmH <sub>2</sub> O ~ 50 cmH <sub>2</sub> O
VTIDAL HIGH	5 ml ~ 2500 ml ,OFF
VTIDAL LOW	0 ml ~ 2500 ml
PEEP HIGH	5 cmH <sub>2</sub> O ~ 45 cmH <sub>2</sub> O
PEEP LOW	0 cmH <sub>2</sub> O ~ 7 cmH <sub>2</sub> O
O <sub>2</sub> HIGH	1% ~ 100%, OFF
O <sub>2</sub> LOW	0% ~ 100%
R-R HIGH	3 ~ 150 BPM
R-R LOW	2 ~ 149 BPM
VMIN HIGH	0.1 ~ 50 LPM
VMIN LOW	0 ~ 49.9 LPM
AIR FAIL	-
O <sub>2</sub> FAIL	-
AIR LEAK	0 ml ~ 500 ml
AIR OBST	-
APNEA	OFF, 10 ~ 60 sec
VENT INOP	-
O <sub>2</sub> SENSOR ERR	-
CIRCUIT OPEN	-
BATTERY MODE	-
LOW BATTERY	DC power loss
COM ERROR	

- ㉞) AIR FAIL / O<sub>2</sub> FAIL  
Supplied pressure of Air and O<sub>2</sub> is low. Check the Air and O<sub>2</sub> line.
- ㉟) AIR OBST(Air Obstructed tube)  
It is problem between breathing circuit and patient. Check the breathing circuit.
- ㊱) VENT INOP(Vent in- operation)  
The ventilator is abnormal status. Check the all ventilator system.



- ㄹ) O<sub>2</sub> SENSOR ERR  
It is problem O<sub>2</sub> sensor. Check O<sub>2</sub> sensor.
- ㄷ) CIRCUIT OPEN  
It is problem between breathing circuit and patient. Check the breathing circuit.
- ㄹ) BATTERY MODE  
There is not connected to AC power. DC power is working.
- ㄺ) COM ERROR  
Between DSC and Pneumatic module, communication error is occurred. Check the all ventilator system.



## Appendix A

<b>Appendix</b> .....	90
A. Operating theory .....	90
B. Communications .....	92
C. Parts accessories .....	95

## A. Operating theory

Pneumatic System of MEKANT is composed of 2 pcs of gas circuit. One is for Oxygen, the other is for Air. Proportional solenoid Valve (PSOLs) is the core of Pneumatic System. The PSOLs is controlled accurately by Flow sensor and Microprocessor so that provide accuracy respiration to patient.

MEKANT's Pneumatic System is divided into inspiratory module, expiratory module and patient circuit system. Inspiratory module is subdivided into Gas supply, Flow control, Safety Valve and Inspiration monitoring system.




Figure A-1. PNEUMATIC Block Diagram

### Blending Part

Gas supply system provides the equipment high pressure air (35~90 psi) and Oxygen from external source of supply, and sense that supplying pressure is not adequate. The high supplied pressure air and Oxygen is decompressed into 10 psi through each regulator. Adjustment pressure of Regulator should be measured by pressure valve (SV1, SV2).



**Inspiration Part**

Flow control system is the most important part in Pneumatic System. This is composed of Flow sensor (DP<sub>I</sub>) and PSOLs (PV<sub>I</sub>). This system adjusts Tidal volume, Peak Flow, Flow pattern, Respiratory Rate, Oxygen Percentage. Control of PSOLs is done on measured value by Microprocessor. Flow sensor is Differential Pressure type.

**Safety System**

This is divided in Safety Valve section and Inspiration monitoring section. Safety valve section release excess pressure in the equipment and make a patient breathe with room air as a passageway.

Safety valve open and release pressure in case interior pressure reach to 115 cmH<sub>2</sub>O. And it is a role of passageway to connect to room air in case equipment is defect or no power, patient's expiratory limb is closed.

Inspiration monitoring section monitor gas pressure and oxygen concentration provided the patient. Monitor is measured by Absolute Pressure Transducer (PA) which is attached inspiratory Pressure Transducer (P<sub>I</sub>). Oxygen is done by Oxygen sensor.

**Expiration Part**

This section is for Exhalation. PEEP/CPAP is controlled by Exhalation Valve during exhalation. Exhalation Valve of MEKANT is Active type and is similar with PSOL. This active exhalation valve makes tBiLEVEL mode possible and prevent pressure increasing suddenly.

Flow sensor (DP<sub>E</sub>) measure Flow and Volume of patient's expiration gas, and Expiratory pressure Transducer(P<sub>E</sub>) measure pressure.

**Patient Circuit System**


Patient circuit system transits gas which is adjusted in Flow control system to patient, and makes patient's exhalation gas to exhalation system.

There is each Bacteria Filter on the connecting part of patient circuit system and equipment. So they make patient circuit system separate from the ventilator. It makes gas to be provided the patient heat and humidify in case humidifier is used.



## B. Communications

### VGA




Pin	Name	Direction	Description
①	RED	→	Red Video (V <sub>ohm</sub> , .V <sub>p-p</sub> )
②	GREEN	→	Green Video(V <sub>ohm</sub> , .V <sub>p-p</sub> )
③	BLUE	→	Blue Video(V <sub>ohm</sub> , .V <sub>p-p</sub> )
④	ID <sub>Y</sub>	←	Not connected
⑤	GND	—	Ground
⑥	RGND	—	Red Ground
⑦	GGND	—	Green Ground
⑧	BGND	—	Blue Ground
⑨	KEY		Not connected
⑩	SGND	—	Sync Ground
⑪	ID <sub>•</sub>	←	Not connected
⑫	ID <sub>^</sub> or SDA	←	Not connected
⑬	HSYNC	→	Horizontal Sync
	CSYNC	→	Composite Sync
⑭	VSYNC	←	Vertical Sync
⑮	ID <sub>Y</sub> or SCL		Not connected

Table B-1. VGA serial port pin out



## RS-232(serial) port

A 9-pin male connector configured as data terminal equipment (DTE). Allowable current is 0.5A at 10V DC (maximum).



Pin	Name	Direction	Description
①	CD(Carrier Detect)		Not connected
②	RxD(Receive Data)	←	Receive data (RxD)
③	TxD(Transmit Data)	→	Transmit data (TxD)
④	DTR(Data Terminal Ready)		Not connected
⑤	SG(Signal Ground)	—	Ground(GND)
⑥	DSR(Data Set Ready)		Not connected
⑦	RTS(Request to Send)		Not connected
⑧	CTS(Clear to Send)		Not connected
⑨	RI(Ring Indicate)		Not connected


Table B-2. RS-232 serial port pin out




Communication capabilities (cont)

NOTE :

The connection of accessories or equipment to the ventilator’s RS-232 and remote alarm ports requires electrical separation of conductive earth (ground) connections between different items of equipment of a system. The system cable assembly shielding can interconnect the ventilator and other equipment, increasing the risk of excessive enclosure leakage current from external equipment. For safe connection, the shield conductor must have separation of 4mm between the equipment and ventilator conductors, as show below:



LCD MONITOR POWER SUPPLY Connector



Pin	Name	Direction	Description
①	NC		Not connected
②	GND		GND
③	+VCC(+12V/5A)		12V/ Max 5A, OUTPUT LCD MONITOR POWER SUPPLY
④	NC		Not connected










Table B-3. LCD Monitor Power Connector port pin out












## Appendix C

## C. Parts accessories

## ■ STANDARD

Item	Photo	Part Number	Q'ty	Description
Circuit Arm		14.4.35 2.	1	
High Pressure Hose		1414... 2.	1	Air Gas / Yellow
		1414... 3.	1	O <sub>2</sub> Gas / Green
Display Panel Cable		1424.16 6.	1	
O <sub>2</sub> Sensor		14.5.1 ..	1	Inside the Factory
PEEP Silicone Plate		14.6.16 3.	2	
Test Lung		14.6.16 2.	1	
Operation Manual		1417.3 8.	1	
Power Cord		14.6.05 3.	1	3m, 220VAC



Ground Cable		۱۴۱۳.۰۲ ۰۰	۱	
Fuse		۱۴۱۳.۰۱ ۷۰	۲	INLET ۵.۰T ۶.۳A ۲۵.۷۲۰mm
Wrench Driver		۱۴.۸.۰۷ ۲۰	2	
Battery		۱۴.۸.۰۷ ۲۰	1	
Circuit Arm Fix Knob		۱۴.۴.۱۰ ۰۰	1	
53010(O-Ring)		۱۴۲۱.۰۴ ۰۰	12	O-Ring for PEEP Module
53030(O-Ring)		۱۴۲۱.۰۴ ۱۰	3	O-Ring for PEEP Module
Monitor Fix Bolt		۱۴.۸.۰۳ ۳۰	3	M5x12mm Wrench Bolt(For Monitor)
Body Cover		۱۵.۰.۰۱ ۶۰	1	



## Contact Us

### Service calls

Followings are telephone numbers and addresses for contacting various service, product supplies and sales personnel.

Reaching open a service call with MEKICS Co., Ltd., contact the numbers listed below.

MEKICS Co., Ltd.

△F, A-BLDG, WooLimLion, #14-2, Sangdaewon-dong,  
Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea 462-725

Tel : + 82-70-7119-2500 Fax : + 82-31-735-2761

URL : [www.mek-ics.com](http://www.mek-ics.com)

e-mail : [service@mek-ics.com](mailto:service@mek-ics.com)

※ In the event of a malfunction or failure, contact Service Dept. of MEKICS Co., Ltd. along with the model name, serial number, date of purchase and explanation of failure.

### ***Authorized Representative***

Pedro Pestana, Lda

Nif : 507662407

Fax : +351214414138

Telef. : +351914881444

Adress: Rua Bernardo Marques, n.º 3, 2.º Dt.º

2770-199 Paco de Arcos