

Service Manual

DIAP Communications Protocol

This publication contains data proprietary to Mindray DS USA, Inc. It is provided solely for use by qualified personnel to interface external devices with certain Mindray devices. Mindray DS USA, Inc. assumes no responsibility for any adverse interaction resulting from the use of this protocol, or for any regulatory issues that may be raised by connecting devices. The details of this protocol are subject to change at any time, without prior notice.

© Mindray DS USA, Inc., 2008. All rights reserved. Contents of this publication may not be reproduced in any form without permission of Mindray DS USA, Inc.

Table of Contents

| | | |
|----------|---|-----|
| I | DIAP | 1 |
| 1.1 | Introduction..... | 1 |
| 1.2 | DIAP | 1 |
| 1.2.1 | Connection Type..... | 1 |
| 1.2.2 | Physical Protocol | 1 |
| 1.2.3 | Session Protocol..... | 1 |
| 1.2.4 | Command/Reply Legend | 2 |
| 1.3 | Commands and Replies..... | 6 |
| 1.4 | Error Handling | 19 |
| 1.4.1 | Field Level Errors | 19 |
| 1.4.2 | Command Level Errors..... | 19 |
| A | Appendix A..... | A-1 |
| A.1 | CRC Algorithm..... | A-1 |
| A.2 | CRC Test Program..... | A-3 |
| B | Appendix B | B-1 |
| B.1 | DIAP Modem Support (Optional) | B-1 |
| B.1.1 | Overview..... | B-1 |
| B.1.2 | Physical Implementation..... | B-1 |
| B.1.3 | DIAP Slave Requirements | B-1 |
| B.1.4 | DIAP Slave Supported Modems | B-1 |
| B.1.5 | DIAP Slave Modem Initialization Strings | B-1 |
| B.1.6 | DIAP Master Requirements..... | B-2 |
| B.1.7 | DIAP Master Supported Modems..... | B-2 |

This page intentionally left blank

1 *DIAP*

1.1 Introduction

This document defines the DIAP for new patient monitor external data communication design.

1.2 DIAP

The DIAP was designed to be simple, yet extensible for new features without compromising programs written by third parties.

1.2.1 Connection Type

The DIAP is a point-to-point, hierarchical (one or more into one host) protocol. All communications are initiated by the host. No common access media types, including serial multi-point connections, are supported.


It is a request/reply, non-continuous protocol. Only numeric data (no waveforms) are supported.

1.2.2 Physical Protocol

The physical interface for the DIAP is asynchronous serial, with a baud rate of 9600 and 19200, 8 bits, 1 stop bit, no parity.

1.2.3 Session Protocol

The messages must be/are encapsulated as follows:



- All Commands are case-insensitive ASCII.
- All Replies are mixed case as indicated in next section.
- Timeout: No reply for 10 seconds.
- Maximum message size for Commands is 255 total bytes.
- Maximum message size for Replies is 511 bytes.

Where:

| Item | Description | Type | Length |
|--------|--|--|--------|
| HEADER | DIAPxxx where: xxx is a channel number from 000 to 999 - it is always 000 for a patient monitor - it is always 001 to 999 for a central monitor | literal ASCII character + ASCII decimal numeric | 4+3 |
| < , > | are delimiter characters | literal ASCII string | 1, 1 |
| CRC | is a cyclical redundancy check as described in the section beginning on page 2-1. | ASCII hex numeric | 4 |
| LF | is an ASCII line feed character. | literal ASCII character | 1 |

1.2.4 Command/Reply Legend

This is the *legend* to the Command/Reply list. It indicates how the arguments for the Commands and Replies are defined. See “Commands and Replies”(section 1.3).

- Items designated as ‘ASCII decimal numeric’ are padded with leading (left) zeros.
- Items designated as ‘ASCII hex numeric’ are padded with leading (left) zeros.
- Items designated as ‘literal ASCII string’ are padded with trailing (right) spaces.
- Items designated as ‘ASCII string’ are padded with trailing (right) spaces.
- Items designated as ‘ASCII string’ will be truncated by the monitor/central as necessary.
- Items designated as ‘ASCII string’ may contain A-Z, a-z, 0-9, comma, period, dash, space, question mark.

| Item | Description | Type | Length | Units |
|-------------|--|-----------------------|--------|-------|
| XXX | generic ASCII number | ASCII decimal numeric | 3 | - |
| S | numeric sign | ASCII ‘+’ or ‘-’ | 1 | - |
| productName | Accutorr Passport Passport ST VISA Passport 2 Spectrum Spectrum OR Trio Passport V Valiant* VSeries * Valiant = Accutorr V | literal ASCII string | 32 | - |
| mVersion | monitor software version | ASCII string | 3 | - |

| | | | | |
|-----------------------|---|-----------------------|---------------|------|
| mRevision | monitor software revision | ASCII decimal numeric | 3 | - |
| pVersion | protocol format version=000 | ASCII decimal numeric | 3 | - |
| pRevision | protocol software revision** ** Valiant and Passport V = 000 VSeries = 000 | ASCII decimal numeric | 3 | - |
| patRoom1 | patient room, line 1 | ASCII string | 4 | - |
| patRoom2 | patient room, line 2 | ASCII string | 4 | - |
| patName1 | patient name, line 1 | ASCII string | 32 | - |
| patName2 | patient name, line 2 | ASCII string | 32 | - |
| unit Name | unit name | ASCII string | 32 | - |
| co2ins | co2 inspired value | ASCII decimal numeric | 3 | torr |
| co2et | co2 end-tidal value | ASCII decimal numeric | 3 | torr |
| o2ins | o2 inspired value (10x) | ASCII decimal numeric | 4 | % |
| o2et | o2 end-tidal value (10x) | ASCII decimal numeric | 4 | % |
| agtins | agent inspired value (100x) | ASCII decimal numeric | 5 | % |
| aget | agent end-tidal value (100x) | ASCII decimal numeric | 5 | % |
| n2oins | n2o inspired value (10x) | ASCII decimal numeric | 4 | % |
| n2oet | n2o end-tidal value (10x) | ASCII decimal numeric | 4 | % |
| sys | systolic value | ASCII decimal numeric | 3 + 1 sign | mmHg |
| dia | diastolic value | ASCII decimal numeric | 3 + 1 sign | mmHg |
| map | map value | ASCII decimal numeric | 3 + 1 sign | mmHg |
| status | status bits | ASCII hex numeric | 4 | - |
| heartRate Source | NONE ECG SpO ₂ NIBP P1 P2 P3 P4 | literal ASCII string | 16 | - |
| Respiration Source | NONE ECG CO ₂ PAW TV | literal ASCII string | 16 | - |
| beatType | NONE | literal ASCII string | 16 | - |

| | | | | |
|-------------|---|--|-----------------------------|---------------------------|
| rhythm | NONE | literal ASCII string | 16 | - |
| temp | temperature (10x) (100 X for Accutorr Plus) | ASCII decimal numeric | 4 + 1 sign | °C |
| sec | seconds | ASCII decimal numeric | 2 | - |
| min | minutes | ASCII decimal numeric | 2 | - |
| hour | hours | ASCII decimal numeric | 2 | - |
| day | day of month | ASCII decimal numeric | 2 | - |
| month | month of year | ASCII decimal numeric | 2 | - |
| year | year | ASCII decimal numeric | 4 | - |
| weekday | Mon Tue Wed Thu Fri Sat Sun | literal ASCII string | 3 | - |
| timeInSecs | elapsed time (since midnight, 01- Jan-1990) | ASCII decimal numeric | 10 | - |
| Lo | low alarm value or "OFF" | ASCII decimal numeric or literal ASCII string | 3 + 1 sign or 1 space | - |
| Hi | high alarm value or "OFF" | ASCII decimal numeric or literal ASCII string | 3 + 1 sign or 1 space | - |
| pawpValue | pulmonary artery wedge pressure | ASCII decimal numeric | 3 + 1 sign | mmHg |
| coValue | cardiac output value (100x) | ASCII decimal numeric | 4 | liters/min |
| ciValue | cardiac index value (10x) | ASCII decimal numeric | 3 | liters/min/m ² |
| agentName | ISO ENF HAL DES SEV | literal ASCII string | 3 | - |
| stValue | ECG s-t value (10x) | ASCII decimal numeric | 2 + 1 sign | mm |
| timeInMins | time in minutes | ASCII decimal numeric | 3 | minutes |
| apneaAlarm | alarm delay for Apnea or "OFF" | ASCII decimal numeric or literal ASCII string | 3 | seconds |
| histString1 | Item#,Date,Time,Call (001,XX/XX/XX,24:00,Text) | ASCII String | 4+9+6+1 3 | - |

| | | | | |
|---------------|--|--|-----------------------------|------------------------|
| histString100 | Item#,Date,time,Call (100,XX/XX/XX,24:00,Text) | ASCII String | 4+9+6+1 3 | - |
| height | Height in centimeters (10x) | ASCII decimal numeric | 4 | cm |
| weight | Weight in kilograms (100x) | ASCII decimal numeric | 5 | kg |
| bsa | Body Surface Area using Dubois Algorithm BSA = height(cm)0.725 x weight(kg)0.425 x 0.007184 (1000x) | ASCII decimal numeric | 5 | Dubois Algorithm |
| cpp | Cerebral Perfusion Pressure calculated as the difference between the Art Mean and the ICP Mean pressure | ASCII decimal numeric | 3 + 1 sign | mmHg |
| mode | Monitor mode of "Standby", "Demo" or "Normal" is returned | ASCII string | 8 | - |
| vtValue | insp/exp Tidal Volume | ASCII decimal numeric | 4 | ml |
| mvValue | insp/exp Minute Volume (10x) | ASCII decimal numeric | 3 | ml |
| ieRatio | Ratio of insp time to exp time | ASCII string | 5 | - |
| pValue | Spirometry pressure | ASCII decimal numeric | 3 + 1 sign or 1 space | cmH ₂ O |
| complValue | Dynamic Airway compliance value | ASCII decimal numeric | 3 | ml/ cmH ₂ O |
| macValue | Minimum Alveolar Concentration (10x) | ASCII decimal numeric | 2 | - |
| vtLo | low alarm value or "OFF" of Tidal Volume | ASCII decimal numeric or literal ACSII string | 4 | ml |
| vtHi | low alarm value or "OFF" of Tidal Volume | ASCII decimal numeric or literal ACSII string | 4 | ml |
| mvLo | low alarm value or "OFF" of Flow Minute Volume (10x) | ASCII decimal numeric or literal ACSII string | 3 + 1 sign or 1 space | ml |
| mvHi | high alarm value or "OFF" of Flow Minute Volume (10x) | ASCII decimal numeric or literal ACSII string | 3 + 1 sign or 1 space | ml |
| spSensorSize | Spirometry Sensor Size Adult Ped ----- | Literal ACSII string | 5 | - |
| rawValue | Dynamic Airway Resistance | ASCII decimal numeric | 2 | cmH ₂ O/l/s |

| | | | | |
|-------|-----------------------------------|-----------------------|---|-------------------|
| bisef | BIS Spectral Edge Frequency (10x) | ASCII decimal numeric | 3 | Hz |
| sbis | Standard Deviation of BIS (10x) | ASCII decimal numeric | 3 | N/A |
| semg | Standard Deviation of EMG (10x) | ASCII decimal numeric | 3 | dB |
| svi | Stroke Volume Index (10x) | ASCII decimal numeric | 4 | ml/m ² |
| edvi | End Diastolic Volume Index (10x) | ASCII decimal numeric | 4 | ml/m ² |
| esvi | End Systolic Volume Index (10x) | ASCII decimal numeric | 4 | ml/m ² |
| o2eii | Oxygen Extraction Index (10x) | ASCII decimal numeric | 3 | % |

1.3 Commands and Replies

The following is a list of alternating **Commands** and **Replies**. They are encapsulated in the format above. There are no space characters within the **Commands** or **Replies** - they are indicated in the table for clarity only.

- Multiple **Commands** may be sent in one message; they must be separated by semi-colons (;); the **Replies** will be returned in a single message separated by semi-colons. Having replies returned as a single message insures that the information in the replies are all synchronized with respect to one another (atomic transaction).

Example:

Command:

```
DIAP000<t1;t2;deltaT>
```

Reply:

```
DIAP000<t1=42;t2=39;deltaT=3>
```

The value for deltaT will reflect what the values were for t1 and t2 at the time the message was received.

Command:

```
DIAP000<t1>
```

Reply:

```
DIAP000<t1=42>
```

Command:

```
DIAP000<t2>
```

Reply:

```
DIAP000<t2=39>
```

Command:

```
DIAP000<deltaT> (meanwhile, t1 changed from 42 to 40)
```

Reply:

```
DIAP000<deltaT=1>
```

The value of deltaT reflects the most current temperature values and not necessarily those from previous commands.

- If a particular imbedded **Command** is invalid, the particular **Reply** will be “UNSUPPORTED FEATURE”.

Example

Command:

```
DIAP000<ecgHR;p3>
```

Reply:

```
DIAP000<ecgHR=081;p3=UNSUPPORTED FEATURE>
```

- Command keywords ending with a colon (:) will control the DIAP slave.

Example

Command:

```
DIAP001<patient:A10 ,ICU3,John Doe ,>
```

Reply:

```
DIAP001<patient=A10 ,ICU3,John Doe ,>
```

-Indicating that the controlling Command has been accepted by the DIAP Slave.

NOTE: There are 24 trailing spaces after the name “John Doe”, for a total of 32 characters.

- The commands that are supported by a product are indicated by a checkmark (✓) in the product's column(s). All other commands will receive the **UNSUPPORTED FEATURE** reply.
- Protocol changes that affect Command or Reply or Message formats will cause **pVersion** to be changed; changes that affect whether a feature is supported on a given product will cause **pRevision** to change for that product.

This page intentionally left blank.

| Type | Command | Reply | Passport Rev. 009 | VISA/ Ambulatory | VISA/ Passport | Accutorr Plus REV. G | Passport 2 | Spectrum | Spectrum OR | Trio | Passport V Ver 1.03 | Passport V Ver 1.04 | Accutorr V | VSeries | Notes |
|-------|--------------|---------------------|-------------------|------------------|----------------|----------------------|------------|----------|-------------|------|---------------------|---------------------|------------|---------|-------|
| Alarm | alarmHR | alarmHR=lo, hi | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | |
| Alarm | alarmSpo2 | alarmSpo2=lo, hi | √ | | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | |
| Alarm | alarmP1Sys | alarmP1Sys=lo, hi | √ | | √ | | √ | √ | √ | √ | √ | √ | | √ | |
| Alarm | alarmP1Dia | alarmP1Dia=lo, hi | √ | | √ | | √ | √ | √ | √ | √ | √ | | √ | |
| Alarm | alarmP1Map | alarmP1Map=lo, hi | | | | | √ | √ | √ | √ | √ | √ | | √ | |
| Alarm | alarmP2Sys | alarmP2Sys=lo, hi | | | √ | | √ | √ | √ | | √ | √ | | √ | |
| Alarm | alarmP2Dia | alarmP2Dia=lo, hi | | | √ | | √ | √ | √ | | √ | √ | | √ | |
| Alarm | alarmP2Map | alarmP2Map=lo, hi | √ | | √ | | √ | √ | √ | | √ | √ | | √ | |
| Alarm | alarmP3Sys | alarmP3Sys=lo, hi | | | | | | √ | √ | | | | | √ | |
| Alarm | alarmP3Dia | alarmP3Dia=lo, hi | | | | | | √ | √ | | | | | √ | |
| Alarm | alarmP3Map | alarmP3Map=lo, hi | | | | | | √ | √ | | | | | √ | |
| Alarm | alarmP4Sys | alarmP4Sys=lo, hi | | | | | | √ | √ | | | | | √ | |
| Alarm | alarmP4Dia | alarmP4Dia=lo, hi | | | | | | √ | √ | | | | | √ | |
| Alarm | alarmP4Map | alarmP4Map=lo, hi | | | | | | √ | √ | | | | | √ | |
| Alarm | alarmP5Sys | alarmP5Sys=lo, hi | | | | | | | | | | | | √ | |
| Alarm | alarmP5Dia | alarmP5Dia=lo, hi | | | | | | | | | | | | √ | |
| Alarm | alarmP5Map | alarmP5Map=lo, hi | | | | | | | | | | | | √ | |
| Alarm | alarmP6Sys | alarmP6Sys=lo, hi | | | | | | | | | | | | √ | |
| Alarm | alarmP6Dia | alarmP6Dia=lo, hi | | | | | | | | | | | | √ | |
| Alarm | alarmP6Map | alarmP6Map=lo, hi | | | | | | | | | | | | √ | |
| Alarm | alarmP7Sys | alarmP7Sys=lo, hi | | | | | | | | | | | | √ | |
| Alarm | alarmP7Dia | alarmP7Dia=lo, hi | | | | | | | | | | | | √ | |
| Alarm | alarmP7Map | alarmP7Map=lo, hi | | | | | | | | | | | | √ | |
| Alarm | alarmP8Sys | alarmP8Sys=lo, hi | | | | | | | | | | | | √ | |
| Alarm | alarmP8Dia | alarmP8Dia=lo, hi | | | | | | | | | | | | √ | |
| Alarm | alarmP8Map | alarmP8Map=lo, hi | | | | | | | | | | | | √ | |
| Alarm | alarmICPxSys | alarmICPxSys=lo, hi | | | | | | | | | | | | | |
| Alarm | alarmICPxDia | alarmICPxDia=lo, hi | | | | | | | | | | | | | |
| Alarm | alarmICPxMap | alarmICPxMap=lo, hi | | | | | | | | | | | | | |
| Alarm | alarmNibpSys | alarmNibpSys=lo, hi | √ | | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | |
| Alarm | alarmNibpDia | alarmNibpDia=lo, hi | √ | | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | |

| Type | Command | Reply | Passport Rev. 009 | VISA/ Ambulatory | VISA/ Passport | Accutorr Plus REV. G | Passport 2 | Spectrum | Spectrum OR | Trio | Passport V Ver 1.03 | Passport V Ver 1.04 | Accutorr V | VSeries | Notes |
|-------|--------------|-----------------------|-------------------|------------------|----------------|----------------------|------------|----------|-------------|------|---------------------|---------------------|------------|---------|-------|
| Alarm | alarmNibpMap | alarmNibpMap=lo, hi | | | | √ | √ | √ | √ | √ | √ | √ | √ | √ | |
| Alarm | alarmT1 | alarmT1=lo, hi | | | | | √ | √ | √ | √ | √ | √ | | √ | |
| Alarm | alarmT2 | alarmT2=lo, hi | | | | | | √ | √ | √ | | | | √ | |
| Alarm | alarmT3 | alarmT3=lo, hi | | | | | | | | | | | | √ | |
| Alarm | alarmICT | alarmICT=lo, hi | | | | | | | | | | | | | |
| Alarm | alarmDeltaT | alarmDeltaT=lo, hi | | | | | | √ | √ | | | | | √ | |
| Alarm | alarmResp | alarmResp=lo, hi | √ | | √ | | √ | √ | √ | √ | √ | √ | | √ | |
| Alarm | alarmCo2 | alarmCo2=lo, hi | √ | | √ | | | | | √ | | | | | |
| Alarm | alarmCo2Et | alarmCo2Et=lo, hi | | | | | √3 | √3 | √3 | √ | | √4 | | √ | |
| Alarm | alarmCo2Ins | alarmCo2Ins=lo, hi | | | | | √3 | √3 | √3 | √ | | √4 | | √ | |
| Alarm | alarm02 | alarm02=lo, hi | | | | | | | | √ | | | | | |
| Alarm | alarm02Et | alarm02Et=lo, hi | | | | | √ | √ | √ | √ | | √ | | √ | |
| Alarm | alarm02Ins | alarm02Ins=lo, hi | | | | | √ | √ | √ | √ | | √ | | √ | |
| Alarm | alarmAgent | alarmAgent=lo, hi | | | | | | | | √ | | | | | |
| Alarm | alarmHallIns | alarmHallIns=lo, hi | | | | | √ | √ | √ | √ | | √4 | | √ | |
| Alarm | alarmHalEt | alarmHalEt=lo, hi | | | | | √ | √ | √ | √ | | √4 | | √ | |
| Alarm | alarmIsoIns | alarmIsoIns=lo, hi | | | | | √ | √ | √ | √ | | √4 | | √ | |
| Alarm | alarmIsoEt | alarmIsoEt=lo, hi | | | | | √ | √ | √ | √ | | √4 | | √ | |
| Alarm | alarmEnflIns | alarmEnflIns=lo, hi | | | | | √ | √ | √ | √ | | √4 | | √ | |
| Alarm | alarmEnflEt | alarmEnflEt=lo, hi | | | | | √ | √ | √ | √ | | √4 | | √ | |
| Alarm | alarmSevoIns | alarmSevoIns=lo, hi | | | | | √ | √ | √ | √ | | √4 | | √ | |
| Alarm | alarmSevoEt | alarmSevoEt=lo, hi | | | | | √ | √ | √ | √ | | √4 | | √ | |
| Alarm | alarmDesIns | alarmDesIns=lo, hi | | | | | √ | √ | √ | √ | | √4 | | √ | |
| Alarm | alarmDesEt | alarmDesEt=lo, hi | | | | | √ | √ | √ | √ | | √4 | | √ | |
| Alarm | alarmN2o | alarmN2o=lo, hi | | | | | √ | √ | √ | √ | | √ | | √ | |
| Alarm | alarmApnea | alarmApnea=apneaAlarm | √ | | √ | | √ | √ | √ | √ | √ | √ | | √ | |
| Alarm | alarmTblood | alarmTblood=lo, hi | | | | | | √ | √ | | | | | √ | |
| Alarm | Baud Rate | 9600 | √ | √ | √ | √ | √ | √ | √ | √ | | √ | √ | | |
| Alarm | Baud Rate | 19200 | | | | | √ | √ | √ | √ | | √ | | | |
| Alarm | alarmvtinsp | alarmvtinsp=vtLo,vtHi | | | | | | | | √ | | | | | |
| Alarm | alarmvtexp | alarmvtexp=vtLo,vtHi | | | | | | | | √ | | | | | |

| Type | Command | Reply | Passport Rev. 009 | VISA/ Ambulatory | VISA/ Passport | Accutorr Plus REV. G | Passport 2 | Spectrum | Spectrum OR | Trio | Passport V Ver 1.03 | Passport V Ver 1.04 | Accutorr V | VSeries | Notes |
|-------|-------------|-----------------------|-------------------|------------------|----------------|----------------------|------------|----------|-------------|------|---------------------|---------------------|------------|---------|-------|
| Alarm | alarmmvinsp | alarmmvinsp=mvLo,mvHi | | | | | | | √ | | | | | | |
| Alarm | alarmmvexp | alarmmvexp=mvLo,mvHi | | | | | | | √ | | | | | | |
| Alarm | alarmplat | alarmplat=lo,hi | | | | | | | √ | | | | | | |
| Alarm | alarmpeep | alarmpeep=lo,hi | | | | | | | √ | | | | | | |
| Alarm | alarmpip | alarmpip=lo,hi | | | | | | | √ | | | | | | |
| Alarm | alarmpmean | alarmpmean=lo,hi | | | | | | | √ | | | | | | |
| Alarm | alambis | alambis=lo,hi | | | | | | | √ | | | | | √ | |
| Alarm | alarmartsys | alarmartsys=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmartdia | alarmartdia=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmartmap | alarmartmap=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmpasys | alarmpasys=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmpadia | alarmpadia=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmpamap | alarmpamap=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmcvpsys | alarmcvpsys=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmcvpdia | alarmcvpdia=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmcvpmap | alarmcvpmap=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmuasys | alarmuasys=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmuadia | alarmuadia=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmuamap | alarmuamap=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmlvsys | alarmlvsys=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmlvdia | alarmlvdia=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmlvmap | alarmlvmap=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmlasys | alarmlasys=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmladia | alarmladia=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmlamap | alarmlamap=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmrasys | alarmrasys=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmradia | alarmradia=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmramap | alarmramap=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmicpsys | alarmicpsys=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmicpdia | alarmicpdia=lo,hi | | | | | | √ | √ | | | √ | | √ | |
| Alarm | alarmicpmap | alarmicpmap=lo,hi | | | | | | √ | √ | | | √ | | √ | |

| Type | Command | Reply | Passport Rev. 009 | VISA/ Ambulatory | VISA/ Passport | Accutorr Plus REV. G | Passport 2 | Spectrum | Spectrum OR | Trio | Passport V Ver 1.03 | Passport V Ver 1.04 | Accutorr V | VSeries | Notes |
|---------------|---|---|-------------------|------------------|----------------|----------------------|---------------|----------|-------------|------|---------------------|---------------------|------------|---------|---|
| Alarm | alarmicpxsys | alarmicpxsys=lo,hi | | | | | | | | | | | | | |
| Alarm | alarmicpxdia | alarmicpxdia=lo,hi | | | | | | | | | | | | | |
| Alarm | alarmicpxmap | alarmicpxmap=lo,hi | | | | | | | | | | | | | |
| Alarm | alarmsvo2 | alarmsvo2=lo,hi | | | | | | | | | | | | √ | |
| Alarm | alarmcco | alarmcco=lo,hi | | | | | | | | | | | | √ | |
| Alarm | alarmcci | alarmcci=lo,hi | | | | | | | | | | | | √ | |
| Alarm | alarmrso2c1 | alarmrso2c1=lo,hi | | | | | | | | | | | | √ | |
| Alarm | alarmrso2c2 | alarmrso2c2=lo,hi | | | | | | | | | | | | √ | |
| Alarm | alarmrso2c3 | alarmrso2c3=lo,hi | | | | | | | | | | | | √ | |
| Alarm | alarmrso2c4 | alarmrso2c4=lo,hi | | | | | | | | | | | | √ | |
| Control | patient:patRoom1, patRoom2, part Name1, partName2 | patient=patRoom1,patRoom2, part Name1,partName2 | | √ | √ | | Not Supported | | | | | | | | available only when patient is discharged |
| Control | unit:unitName | unit=unitName | | √ | √ | | | | | | | | | | |
| Informational | htcm | htcm=XXXX | | | | | | √ | √ | | √ | √ | | √ | |
| Informational | wtkg | wtkg=XXXXX | | | | | | √ | √ | | | | | √ | |
| Informational | bsa | bsa=XXXXX Dubois | | | | | | √ | √ | | | | | √ | |
| Informational | connect | connect=productName, pVersion,pRevision | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ |
| Informational | disconnect | disconnect= | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ |
| Informational | product | product=productName, mVersion,mRevision | | √ | √ | | | | | | | | | | |
| Informational | patient | patient=patRoom1,patRoom2, patName1,patName2 | | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | |
| Informational | unit | unit=unitName | | √ | √ | | | | | | | | | | |
| Informational | Mode | Mode | | | | | | | | | | | | | Currently not implemented |
| Informational | spsensorsize | spsensorsize=spSensorSize | | | | | | | √ | | | | | | |
| Physiological | ecgResp | ecgResp=XXX | | | | | | | | √ | | | | √ | |
| Physiological | ecgSt1 | ecgSt1=stValue | | | | | √1 | √1 | √1 | | √ | √ | | | |
| Physiological | ecgSt2 | ecgSt2=stValue | | | | | √1 | √1 | √1 | | √ | √ | | | |
| Physiological | ecgST3 | ecgST3=stValue | | | | | √1 | √1 | √1 | | √ | √ | | | |
| Physiological | ecgHR | ecgHR=XXX | | | | | | | | √ | | | | √ | |

| Type | Command | Reply | Passport Rev. 009 | VISA/ Ambulatory | VISA/ Passport | Accutorr Plus REV. G | Passport 2 | Spectrum | Spectrum OR | Trio | Passport V Ver 1.03 | Passport V Ver 1.04 | Accutorr V | VSeries | Notes |
|---------------|---------------------|--|-------------------|------------------|----------------|----------------------|------------|----------|-------------|------|---------------------|---------------------|------------|---------|--------------------------|
| Physiological | ecgBeatType | ecgBeatType=beatTime | | | | | | | | √ | | | | | |
| Physiological | ecgRhythm | ecgRhythm=rhythm | | | | | | | | √ | | | | | |
| Physiological | co | co=coValue,ciValue, timeInSeconds | | | | | | √ | √ | | | | | √ | |
| Physiological | p1 | p1=sys,dia,map | √ | | √ | | √ | √ | √ | √ | √ | √ | | √ | |
| Physiological | p1HR | p1HR=XXX | | | | | | | | | | | | √ | |
| Physiological | p2 | p2=sys,dia,map | √ | | √ | | √ | √ | √ | | √ | √ | | √ | |
| Physiological | p2HR | p2HR=XXX | | | | | | | | | | | | √ | |
| Physiological | p3 | p3=sys,dia,map | | | | | | √ | √ | | | | | √ | |
| Physiological | p3HR | p3HR=XXX | | | | | | | | | | | | √ | |
| Physiological | p4 | p4=sys,dia,map | | | | | | √ | √ | | | | | √ | |
| Physiological | p4HR | p4HR=XXX | | | | | | | | | | | | √ | |
| Physiological | p5 | p5=sys,dia,map | | | | | | | | | | | | √ | |
| Physiological | p5HR | p5HR=XXX | | | | | | | | | | | | √ | |
| Physiological | p6 | p6=sys,dia,map | | | | | | | | | | | | √ | |
| Physiological | p6HR | p6HR=XXX | | | | | | | | | | | | √ | |
| Physiological | p7 | p7=sys,dia,map | | | | | | | | | | | | √ | |
| Physiological | p7HR | p7HR=XXX | | | | | | | | | | | | √ | |
| Physiological | p8 | p8=sys,dia,map | | | | | | | | | | | | √ | |
| Physiological | p8HR | p8HR=XXX | | | | | | | | | | | | √ | |
| Physiological | icpx | icpx=sys,dia,map | | | | | | | | | | | | | |
| Physiological | icpxHR | icpxHR=XXX | | | | | | | | | | | | | |
| Physiological | historyLine | historyLine=histStringX | | √ | √ | | | | | | | | | | most current line number |
| Physiological | historyFile | historyFile=histString1, histString2, ...histString100 | | √ | √ | | | | | | | | | | All available lines |
| Physiological | pawp | Pawp=pawpValue, timeInSeconds | | | | | √ | √ | √ | √ | | | | √ | |
| Physiological | t1 | t1=temp | √ | | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | |
| Physiological | t2 | t2=temp | | | | | | √ | √ | √ | | | | √ | |
| Physiological | t3 | t3=temp | | | | | | | | | | | | √ | |
| Physiological | ict | ict=temp | | | | | | | | | | | | | |
| Physiological | deltaT ² | deltaT=temp | | | | | | √ | √ | | | | | √ | |

| Type | Command | Reply | Passport Rev. 009 | VISA/ Ambulatory | VISA/ Passport | Accutorr Plus REV. G | Passport 2 | Spectrum | Spectrum OR | Trio | Passport V Ver 1.03 | Passport V Ver 1.04 | Accutorr V | VSeries | Notes |
|---------------|-----------------|------------------------------|-------------------|------------------|----------------|----------------------|------------|----------|-------------|------|---------------------|---------------------|------------|---------|----------------|
| Physiological | tBlood | tBlood=temp | | | | | | √ | √ | | | | | √ | |
| Physiological | hr | hr=XXX | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | |
| Physiological | hrSource | hrSource=heartRateSource | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | |
| Physiological | hrBeatType | hrBeatType=beatType | | | | | | | | √ | | | | | |
| Physiological | hrRythm | hrRythm=rhythm | | | | | | | | √ | | | | | |
| Physiological | nibp | nibp=sys, dia, map | √ | | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | |
| Physiological | nibpStatus | nibpStatus=status | √ | | | | | | | √ | | | | | |
| Physiological | nibpHr | nibpHr=XXX | √ | | | √ | | | | √ | | | √ | | |
| Physiological | nibpElapsedTime | nibpElapsedTime=timeInMins | √ | | √ | | √ | √ | √ | √ | √ | √ | | √ | |
| Physiological | nibpStartTime | nibpStartTime=timeInSecs | | | | | | | | | | | | | |
| Physiological | nibpStopTime | nibpStopTime=timeInSecs | | | | √ | | | | | | | √ | | |
| Physiological | resp | resp=XXX | √ | | √ | | √ | √ | √ | √ | √ | √ | | √ | |
| Physiological | respSource | respSource=respirationSource | √ | | √ | | √ | √ | √ | √ | √ | √ | | √ | |
| Physiological | co2 | co2=co2Ins,co2Et | √ | | √ | | √ | √ | √ | √ | | √4 | | √ | co2Ins="- - -" |
| Physiological | co2Resp | co2Resp=XXX | | | | | | | | √ | | | | √ | |
| Physiological | spo2 | spo2=XXX | √ | | √ | √ | √ | √ | √ | √ | √ | √ | √ | √ | |
| Physiological | spo2HR | spo2HR=XXX | √ | | | √ | √ | √ | √ | | √ | √ | √ | √ | |
| Physiological | o2 | o2=o2Ins, o2Et | | | | | √ | √ | √ | √ | √ | √ | | √ | |
| Physiological | agent | agent=agtIns, agtEt | | | | | √ | √ | √ | √ | √ | √ | | √ | |
| Physiological | agentType | agentType=agentName | | | | | √ | √ | √ | √ | √ | √ | | √ | |
| Physiological | n2o | n2o=n2oIns, n2oEt | | | | | √ | √ | √ | √ | √ | √ | | √ | |
| Physiological | art | art=sys,dia,map | | | | | | √ | √ | | | | √ | √ | |
| Physiological | arthr | arthr=XXX | | | | | | √ | √ | | √ | √ | | √ | |
| Physiological | pa | pa=sys,dia,map | | | | | | √ | √ | | | | √ | √ | |
| Physiological | pahr | pahr=XXX | | | | | | √ | √ | | √ | √ | | √ | |
| Physiological | cvp | cvp=sys,dia,map | | | | | | √ | √ | | | | √ | √ | |
| Physiological | cvphr | cvphr=XXX | | | | | | √ | √ | | √ | √ | | √ | |
| Physiological | ua | ua=sys,dia,map | | | | | | √ | √ | | | | √ | √ | |
| Physiological | uahr | uahr=XXX | | | | | | √ | √ | | √ | √ | | √ | |
| Physiological | lv | lv=sys,dia,map | | | | | | √ | √ | | | | √ | √ | |
| Physiological | lvhr | lvhr=XXX | | | | | | √ | √ | | √ | √ | | √ | |

| Type | Command | Reply | Passport Rev. 009 | VISA/ Ambulatory | VISA/ Passport | Accutorr Plus REV. G | Passport 2 | Spectrum | Spectrum OR | Trio | Passport V Ver 1.03 | Passport V Ver 1.04 | Accutorr V | VSeries | Notes |
|---------------|---------|------------------|-------------------|------------------|----------------|----------------------|------------|----------|-------------|------|---------------------|---------------------|------------|---------|------------------|
| Physiological | la | la=sys,dia,map | | | | | | √ | √ | | | √ | | | |
| Physiological | lahr | lahr=XXX | | | | | | √ | √ | | | | | | |
| Physiological | ra | ra=sys,dia,map | | | | | | √ | √ | | | √ | | | |
| Physiological | rahr | rahr=XXX | | | | | | √ | √ | | | | | | |
| Physiological | icp | icp=sys,dia,map | | | | | | √ | √ | | | √ | | | |
| Physiological | icphr | icphr=XXX | | | | | | √ | √ | | | | | | |
| Physiological | cpp | cpp=XXX | | | | | | √ | √ | | | √ | | | |
| Physiological | ci | ci=ciValue | | | | | | √ | √ | | | | | | |
| Physiological | vtinsp | vtinsp=vtValue | | | | | | | | | | | | | |
| Physiological | vtexp | vtexp=vtValue | | | | | | | | | | | | | |
| Physiological | mvinsp | mvinsp=mvValue | | | | | | | | | | | | | |
| Physiological | mvexp | mvexp=mvValue | | | | | | | | | | | | | |
| Physiological | ieratio | ieratio=ieRatio | | | | | | | | | | | | | |
| Physiological | pplat | pplat=pValue | | | | | | | | | | | | | |
| Physiological | peep | peep=pValue | | | | | | | | | | | | | |
| Physiological | compl | compl=complValue | | | | | | | | | | | | | |
| Physiological | pip | pip=pValue | | | | | | | | | | | | | |
| Physiological | pmean | pmean=pValue | | | | | | | | | | | | | |
| Physiological | mac | mac=macValue | | | | | | | | | | | | | |
| Physiological | bis | bis=XX | | | | | | | | | | | | | |
| Physiological | sqi | sqi=XX | | | | | | | | | | | | | BIS Left Channel |
| Physiological | emg | emg=XX | | | | | | | | | | | | | BIS Left Channel |
| Physiological | bislr | bislr=XX | | | | | | | | | | | | | |
| Physiological | bislbc | bislbc=XX | | | | | | | | | | | | | |
| Physiological | bislsef | bislsef=bisef | | | | | | | | | | | | | |
| Physiological | bisltp | bisltp=XX | | | | | | | | | | | | | |
| Physiological | sbisl | sbisl=sbis | | | | | | | | | | | | | |
| Physiological | semgl | semgl=semg | | | | | | | | | | | | | |
| Physiological | bisr | bisr=XX | | | | | | | | | | | | | |
| Physiological | sqir | sqir=XX | | | | | | | | | | | | | |

| Type | Command | Reply | Passport Rev. 009 | VISA/ Ambulatory | VISA/ Passport | Accutorr Plus REV. G | Passport 2 | Spectrum | Spectrum OR | Trio | Passport V Ver 1.03 | Passport V Ver 1.04 | Accutorr V | VSeries | Notes |
|---------------|----------|-----------------|-------------------|------------------|----------------|----------------------|------------|----------|-------------|------|---------------------|---------------------|------------|---------|-------|
| Physiological | emgr | emgr=XX | | | | | | | | | | | | √ | |
| Physiological | bislr | bislr=XX | | | | | | | | | | | | √ | |
| Physiological | bisrbc | bisrbc=XX | | | | | | | | | | | | √ | |
| Physiological | bisrsef | bisrsef=bisef | | | | | | | | | | | | √ | |
| Physiological | bisrtp | bisrtp=XX | | | | | | | | | | | | √ | |
| Physiological | sbisr | sbisr=sbis | | | | | | | | | | | | √ | |
| Physiological | semgr | semgr=semg | | | | | | | | | | | | √ | |
| Physiological | bisasym | bisasym=XX | | | | | | | | | | | | √ | |
| Physiological | cco | cco=coValue | | | | | | | | | | | | √ | |
| Physiological | cci | cci=ciValue | | | | | | | | | | | | √ | |
| Physiological | svo2 | svo2=XX | | | | | | | | | | | | √ | |
| Physiological | sv | sv=XXXX | | | | | | | | | | | | √ | |
| Physiological | svi | svi=svi | | | | | | | | | | | | √ | |
| Physiological | edv | edv=XXX | | | | | | | | | | | | √ | |
| Physiological | edvi | edvi=XXXX | | | | | | | | | | | | √ | |
| Physiological | ccostat | ccostat=coValue | | | | | | | | | | | | √ | |
| Physiological | ccostat | ccostat=ciValue | | | | | | | | | | | | √ | |
| Physiological | svr | svr=XXXX | | | | | | | | | | | | √ | |
| Physiological | svri | svri=XXXX | | | | | | | | | | | | √ | |
| Physiological | ccosnr | ccosnr=XX | | | | | | | | | | | | √ | |
| Physiological | svv | svv=XX | | | | | | | | | | | | √ | |
| Physiological | rvef | rvef=XX | | | | | | | | | | | | √ | |
| Physiological | rvefstat | rvefstat=XX | | | | | | | | | | | | √ | |
| Physiological | edvstat | edvstat=XXX | | | | | | | | | | | | √ | |
| Physiological | edvistat | edvistat=XXXX | | | | | | | | | | | | √ | |
| Physiological | esv | esv=XXX | | | | | | | | | | | | √ | |
| Physiological | esvi | esvi=XXXX | | | | | | | | | | | | √ | |
| Physiological | svstat | svstat=XXXX | | | | | | | | | | | | √ | |
| Physiological | svistat | svistat=svi | | | | | | | | | | | | √ | |
| Physiological | rso2c1 | rso2c1=XX | | | | | | | | | | | | √ | |
| Physiological | rso2bc1 | rso2bc1=XX | | | | | | | | | | | | √ | |

| Type | Command | Reply | Passport Rev. 009 | VISA/ Ambulatory | VISA/ Passport | Accutorr Plus REV. G | Passport 2 | Spectrum | Spectrum OR | Trio | Passport V Ver 1.03 | Passport V Ver 1.04 | Accutorr V | VSeries | Notes |
|---------------|----------------|-----------------------|-------------------|------------------|----------------|----------------------|------------|----------|-------------|------|---------------------|---------------------|------------|---------|-------|
| Physiological | rso2auc1 | rso2auc1=XX | | | | | | | | | | | | √ | |
| Physiological | rso2c2 | rso2c2=XX | | | | | | | | | | | | √ | |
| Physiological | rso2bc2 | rso2bc2=XX | | | | | | | | | | | | √ | |
| Physiological | rso2auc2 | rso2auc2=XX | | | | | | | | | | | | √ | |
| Physiological | rso2c3 | rso2c3=XX | | | | | | | | | | | | √ | |
| Physiological | rso2bc3 | rso2bc3=XX | | | | | | | | | | | | √ | |
| Physiological | rso2auc3 | rso2auc3=XX | | | | | | | | | | | | √ | |
| Physiological | rso2c4 | rso2c4=XX | | | | | | | | | | | | √ | |
| Physiological | rso2bc4 | rso2bc4=XX | | | | | | | | | | | | √ | |
| Physiological | rso2auc4 | rso2auc4=XX | | | | | | | | | | | | √ | |
| Physiological | pvcpermin | pvcpermin=XX | | | | | | | | | | | | √ | |
| Physiological | satseconds | satseconds=XX | | | | | | | | | | | | √ | |
| Physiological | perfusionindex | perfusionindex=XX | | | | | | | | | | | | √ | |
| Physiological | pausespermin | pausespermin=XX | | | | | | | | | | | | √ | |
| Physiological | app | app=XX | | | | | | | | | | | | √ | |
| Physiological | cepp | cepp=XX | | | | | | | | | | | | √ | |
| Physiological | stLeadI | stLeadI=stValue | | | | | | | | | | | | √ | |
| Physiological | stLeadII | stLeadII= stValue | | | | | | | | | | | | √ | |
| Physiological | stLeadIII | stLeadIII= stValue | | | | | | | | | | | | √ | |
| Physiological | stLeadAVR | stLeadAVR= stValue | | | | | | | | | | | | √ | |
| Physiological | stLeadAVF | stLeadAVF= stValue | | | | | | | | | | | | √ | |
| Physiological | stLeadAVL | stLeadAVL= stValue | | | | | | | | | | | | √ | |
| Physiological | stLeadV1 | stLeadV1= stValue | | | | | | | | | | | | √ | |
| Physiological | stLeadV2 | stLeadV2= stValue | | | | | | | | | | | | √ | |
| Physiological | stLeadV3 | stLeadV3= stValue | | | | | | | | | | | | √ | |
| Physiological | stLeadV4 | stLeadV4= stValue | | | | | | | | | | | | √ | |
| Physiological | stLeadV5 | stLeadV5= stValue | | | | | | | | | | | | √ | |
| Physiological | stLeadV6 | stLeadV6= stValue | | | | | | | | | | | | √ | |
| Physiological | stRefLeadI | stRefLeadI= stValue | | | | | | | | | | | | √ | |
| Physiological | stRefLeadII | stRefLeadII= stValue | | | | | | | | | | | | √ | |
| Physiological | stRefLeadIII | stRefLeadIII= stValue | | | | | | | | | | | | √ | |

| Type | Command | Reply | Passport Rev. 009 | VISA/ Ambulatory | VISA/ Passport | Accutorr Plus REV. G | Passport 2 | Spectrum | Spectrum OR | Trio | Passport V Ver 1.03 | Passport V Ver 1.04 | Accutorr V | VSeries | Notes |
|---------------|--------------|-----------------------|-------------------|------------------|----------------|----------------------|------------|----------|-------------|------|---------------------|---------------------|------------|---------|-------|
| Physiological | stRefLeadAVR | stRefLeadAVR= stValue | | | | | | | | | | | | √ | |
| Physiological | stRefLeadAVF | stRefLeadAVF= stValue | | | | | | | | | | | | √ | |
| Physiological | stRefLeadAVL | stRefLeadAVL= stValue | | | | | | | | | | | | √ | |
| Physiological | stRefLeadV1 | stRefLeadV1= stValue | | | | | | | | | | | | √ | |
| Physiological | stRefLeadV2 | stRefLeadV2= stValue | | | | | | | | | | | | √ | |
| Physiological | stRefLeadV3 | stRefLeadV3= stValue | | | | | | | | | | | | √ | |
| Physiological | stRefLeadV4 | stRefLeadV4= stValue | | | | | | | | | | | | √ | |
| Physiological | stRefLeadV5 | stRefLeadV5= stValue | | | | | | | | | | | | √ | |
| Physiological | stRefLeadV6 | stRefLeadV6= stValue | | | | | | | | | | | | √ | |
| Physiological | do2 | do2=XXXX | | | | | | | | | | | | | |
| Physiological | vo2 | vo2=XX | | | | | | | | | | | | | |
| Physiological | scvo2 | scvo2=XX | | | | | | | | | | | | √ | |
| Physiological | o2ei | o2ei=o2ei | | | | | | | | | | | | | |
| Physiological | raw | raw=rawValue | | | | | | | | | | | | √ | |

- 1 *In 12-Lead mode this is the absolute value otherwise this the delta value.*
- 2 *deltaT is defined as the absolute value of t1 - t2.*
- 3 *Not supported when using the Gas Module II.*
- 4 *Depends upon type of CO2 gas module installed*

1.4 Error Handling

1.4.1 Field Level Errors

Invalid data fields in replies are represented by two, three, or more dashes (— or —, etc.), based on the expected value length.

1.4.2 Command Level Errors

If an error occurs, an error reply is returned to the host in place of the normal reply.

These **Command Level Errors** are indicated below (literal ASCII string, prefaced with a %):

| Error | Description |
|---------------------------------|--|
| UNABLE TO CONNECT | The monitor cannot respond. |
| NO OPEN DELIMITER | The monitor could not find the character '<' (open delimiter) in the message received, making the message not interpretable (syntax error). |
| NO CLOSE DELIMITER | The monitor could not find the character '>' (close delimiter) in the message received, making the message not interpretable (syntax error). |
| NO LINE FEED | The monitor could not find the linefeed character in the message received, making the message not interpretable (syntax error). |
| INVALID CRC CHAR | The monitor detected a non-hexadecimal CRC character (syntax error). |
| UNDEFINED ERROR | The monitor detected an undefinable syntax error. |
| RESPONSE TOO LARGE | In a multiple command message, the response is greater than the maximum size response buffer. (Separate the multiple command message into two or more messages). |
| REQUEST CORRUPT | The monitor detected an invalid CRC. |
| CHANNEL OUT OF RANGE | The central monitor could not reply because the channel number requested is invalid. |
| BAD MSG HEADER | The message header did not contain "DIAP" as the first four characters |
| COMMAND TOO LARGE | The message size exceeded the maximum of 255 |
| IBP LABEL NOT APPLIED | The monitor cannot return a value because the requested IBP label has not been applied to IBP1, IBP2, IBP3 or IBP4. |
| IBP LABEL NOT CURRENT HR SOURCE | The monitor cannot return a value because the request IBP label has not been set as the current heart rate source. |

This page intentionally left blank.

A Appendix A

A.1 CRC Algorithm

The algorithm used to create the CRC follows:

```
const unsigned short crcTable[256] = {  
0x0000, 0x0108, 0x0210, 0x0318, 0x0420, 0x0528, 0x0630, 0x0738,  
0x0840, 0x0948, 0x0A50, 0x0B58, 0x0C60, 0x0D68, 0x0E70, 0x0F78,  
0x1081, 0x1189, 0x1291, 0x1399, 0x14A1, 0x15A9, 0x16B1, 0x17B9,  
0x18C1, 0x19C9, 0x1AD1, 0x1BD9, 0x1CE1, 0x1DE9, 0x1EF1, 0x1FF9,  
0x2102, 0x220A, 0x2312, 0x221A, 0x2522, 0x242A, 0x2732, 0x263A,  
0x2942, 0x284A, 0x2B52, 0x2A5A, 0x2D62, 0x2C6A, 0x2F72, 0x2E7A,  
0x3183, 0x308B, 0x3393, 0x329B, 0x35A3, 0x34AB, 0x37B3, 0x36BB,  
0x39C3, 0x38CB, 0x3BD3, 0x3ADB, 0x3DE3, 0x3CEB, 0x3FF3, 0x3EFB,  
0x4204, 0x430C, 0x4414, 0x411C, 0x4624, 0x472C, 0x4434, 0x453C,  
0x4A44, 0x4B4C, 0x4854, 0x495C, 0x4E64, 0x4F6C, 0x4C74, 0x4D7C,  
0x5285, 0x538D, 0x5095, 0x519D, 0x56A5, 0x57AD, 0x54B5, 0x55BD,  
0x5AC5, 0x5BCD, 0x58D5, 0x59DD, 0x5EE5, 0x5FED, 0x5CF5, 0x5DFD,  
0x6306, 0x620E, 0x6116, 0x601E, 0x6726, 0x662E, 0x6536, 0x643E,  
0x6B46, 0x6A4E, 0x6956, 0x685E, 0x6F66, 0x6E6E, 0x6D76, 0x6C7E,
```

| | | | | | | | |
|---------|---------|---------|---------|---------|---------|---------|---------|
| 0x7387, | 0x728F, | 0x7197, | 0x709F, | 0x77A7, | 0x76AF, | 0x75B7, | 0x74BF, |
| 0x7BC7, | 0x7ACF, | 0x79D7, | 0x78DF, | 0x7FE7, | 0x7EEF, | 0x7DF7, | 0x7CFF, |
| 0x8408, | 0x8500, | 0x8618, | 0x8710, | 0x8028, | 0x8120, | 0x8238, | 0x8330, |
| 0x8C48, | 0x8D40, | 0x8E58, | 0x8F50, | 0x8868, | 0x8960, | 0x8A78, | 0x8B70, |
| 0x9489, | 0x9581, | 0x9699, | 0x9791, | 0x90A9, | 0x91A1, | 0x92B9, | 0x93B1, |
| 0x9CC9, | 0x9DC1, | 0x9ED9, | 0x9FD1, | 0x98E9, | 0x99E1, | 0x9AF9, | 0x9BF1, |
| 0xA50A, | 0xA402, | 0xA71A, | 0xA612, | 0xA12A, | 0xA022, | 0xA33A, | 0xA232, |
| 0xAD4A, | 0xAC42, | 0xAF5A, | 0xAE52, | 0xA96A, | 0xA862, | 0xAB7A, | 0xAA72, |
| 0xB58B, | 0xB483, | 0xB79B, | 0xB693, | 0xB1AB, | 0xB0A3, | 0xB3BB, | 0xB2B3, |
| 0xBDCB, | 0xBCC3, | 0xBFDB, | 0xBED3, | 0xB9EB, | 0xB8E3, | 0xBBFB, | 0BAF3, |
| 0xC60C, | 0xC704, | 0xC41C, | 0xC514, | 0xC22C, | 0xC324, | 0xC03C, | 0xC134, |
| 0xCE4C, | 0xCF44, | 0xCC5C, | 0xCD54, | 0xCA6C, | 0xCB64, | 0xC87C, | 0xC974, |
| 0xD68D, | 0xD785, | 0xD49D, | 0xD595, | 0xD2AD, | 0xD3A5, | 0xD0BD, | 0xD1B5, |
| 0xDECD, | 0xDFC5, | 0xDCDD, | 0xDDD5, | 0xDAED, | 0xDBE5, | 0xD8FD, | 0xD9F5, |
| 0xE70E, | 0xE606, | 0xE51E, | 0xE416, | 0xE32E, | 0xE226, | 0xE13E, | 0xE036, |
| 0xEF4E, | 0xEE46, | 0xED5E, | 0xEC56, | 0xEB6E, | 0xEA66, | 0xE97E, | 0xE876, |
| 0xF78F, | 0xF687, | 0xF59F, | 0xF497, | 0xF3AF, | 0xF2A7, | 0xF1BF, | 0xF0B7, |
| 0xFFCF, | 0xFEC7, | 0xFDDF, | 0xFCD7, | 0xFBEB, | 0xFAE7, | 0xF9FF, | 0xF8F7, |

```
};  
  
unsigned short Crc (unsigned short lastCrc, unsigned long sizeofBuffer,  
unsigned char *buffer) {  
    unsigned long i;  
    unsigned char index;  
    for (i = 0L; i < sizeofBuffer; i++) {  
        index = buffer[i] ^ ((unsigned char)(lastCrc & 0x00FF));  
        lastCrc = ((lastCrc >> 8) & 0x00FF) ^ (crcTable[index]);  
    }  
    return (lastCrc);  
}
```

A.2 CRC Test Program

The following program can be used to test the above:

```
#include <stdio.h>

int main(void) {
    unsigned char buffer [] = {0, 0, 0, 0, 0, 0x77, 0xCF};
    printf("CRC(Expected 0x73F3) = 0x%04X\n", Crc(0x0, 7, buffer));
    return(0);
}
```

NOTE: For the example above, the value for the last CRC parameter should be initialized to zero (0).

This page intentionally left blank.

B *Appendix B*

B.1 DIAP Modem Support (Optional)

B.1.1 Overview

Optional support for a modem attached to a DIAP Slave is described here.

B.1.2 Physical Implementation

A standard external modem is connected to the DIAP Slave, such as a Passport, Point of View, VISA, etc. The same modem type is connected to a DIAP Master, such as a personal computer, workstation, minicomputer, etc. A public telephone switching system connects the modems.

B.1.3 DIAP Slave Requirements

The Mindray equipment must have the following additions made to implement DIAP Modem Support:

- Configuration Dialog for Modem Connected (yes/no choice, default no)
- Configuration Dialog for Modem Speaker On (yes/no choice, default no)
- User Dialog for Reset Modem (button)

Upon initialization, the DIAP Slave UART is set per section 1.2. The modem initialization string is then sent (section **Error! Reference source not found.**). The modem initialization string is re-sent anytime the Reset Modem button is depressed.

B.1.4 DIAP Slave Supported Modems

The DIAP Slave shall support the following V.32-compliant or higher modems:

- Hayes Accura 144

B.1.5 DIAP Slave Modem Initialization Strings

The DIAP Slave communicates with the modem as follows:

| Slave Control | ASCII Strings Sent | Modem Returns | Description |
|--|--------------------------------------|---------------|---|
| configuration set to Modem Connected (yes) | (>1 sec pause) +++ (>1 sec pause) | (nothing) | "Force Command Mode" (Hayes standard for modem attention) |
| Modem Speaker On (no) | ATZ (carriage return) | OK | "Reset Modem" |
| | AT&F (carriage return) | OK | "Recall Factory Configuration" |
| and/or | AT&K4 (carriage return) | OK | "Enable XON/XOFF Local Flow Control" |
| Reset Modem active | ATM0 (carriage return) | OK | "Speaker Off" |
| | ATS0 = 1 (carriage return) | OK | "Answer in 1 Ring" |

NOTE: The above assumes that the modem is new or reset to default factory configuration.

NOTE: If the configuration Modem Speaker On = 'yes', then do not send the ATM0 command above.

B.1.6 DIAP Master Requirements

The Master computer must initialize the local serial port (setting communications parameters per section 2.2) and initialize the locally- attached modem (including modem reset, setting up off-hook and telephone dialing to the DIAP Slave modem). It also configures the modem in a manner similar to the DIAP Slave modem (XON/XOFF flow control).

B.1.7 DIAP Master Supported Modems

The DIAP Master shall support any V.32-complaint or higher modem.

*Mindray DS USA, Inc. • 800 MacArthur Boulevard • Mahwah, NJ 07430 • USA •
Dom. Customer Service: 1.800.288.2121 • Intl. Customer Service: +1.201.995.8000 •
Dom. Fax: 1.800.926.4275 • Intl. Fax: +1.201.995.8680 • www.mindray.com*

*Mindray Medical Netherlands B.V. • Drs. W. van Royenstraat 8 • P.O. Box 26 • 3870 CA Hoevelaken • The
Netherlands • Tel: +31 33 25 44 911 • Fax: +31 33 25 37 621*

*Mindray (UK) Limited • 3 Percy Road • St. John's Park • Huntingdon • Cambridgeshire PE29 6SZ • United Kingdom •
Tel: 01480 416840 • Fax: 01480 436588*

*Mindray Medical France SARL • Europarc Créteil • 123, Chemin des Bassins •
94035 Créteil Cedex • France • Tel: (0)1.45.13.91.50 • Fax: (0)1.45.13.91.51*

*Mindray Medical Germany GmbH • Zwischen den Bächen 4 • 64625 Bensheim • Deutschland •
Tel: +49.6251.17524-0 • Fax: +49.6251.17524-20*

*Mindray Medical International Ltd. • 2813 Office Tower, Convention Plaza • No 1 Harbour Road • Wanchai • Hong
Kong • Tel: +852 2793 5596 • Fax: +852 2344 8824*

*Medstar Importação e Exportação Ltda • Av. Vereador José Diniz, 3300 • São Paulo, SP • CEP 04804-000 • Brazil •
Tel: 55 11 2872-3385 • Fax: 55 11 2872-3385*

