
Cosmed does not assume the liability for interpretation mistakes of this documentation or for casual or consequential damages in connection with the provision, representation or use of this documentation.

No parts of this manual may be reproduced or transmitted in any form without the express permission of COSMED Srl.

COSMED Software can be installed only in one device.

Excel is a registered trademark of Microsoft Corporation.

DBIII is a registered trademark of Bordland International Inc.

Lotus 123 is a registered trademark of Lotus Development Corporation .

**microQuark User manual, V Edition
September 2003**

Copyright © 2001 COSMED

Copyright © 2003 COSMED

COSMED Srl - Italy

<http://www.cosmed.it>

Part N. C02268-02-91

Table of contents

| | |
|--|-----------|
| Getting started | 7 |
| Important notices | 8 |
| Intended use..... | 8 |
| Warnings..... | 8 |
| Contraindication | 11 |
| Contraindications for the Spirometer tests..... | 11 |
| Absolute contraindications | 11 |
| Relative contraindications | 11 |
| Contraindications for Bronchial provocation tests..... | 11 |
| Absolute contraindications | 12 |
| Relative contraindications | 12 |
| Environmental condition of use | 13 |
| Safety and conformity | 14 |
| Safety..... | 14 |
| EMC | 14 |
| Quality Assurance..... | 14 |
| Medical Device Directive (CE mark) | 14 |
| Keynotes | 15 |
| Typographic keynotes..... | 15 |
| Graphic keynotes | 15 |
| Systems Overview | 16 |
| Before starting | 17 |
| Checking the packing contents | 17 |
| microQuark standard packaging | 17 |
| Warranty registration | 18 |
| Register the product via software | 18 |
| How to contact COSMED | 18 |
| Complain, feedback and suggestions..... | 18 |
| PC configuration required | 19 |
| Technical features | 20 |
| Measurements | 21 |
| Measured parameters | 22 |

| | |
|--|-----------|
| FVC - Forced Vital Capacity..... | 22 |
| VC/IVC - Slow Vital Capacity and Ventilatory pattern | 23 |
| MVV - Maximum Voluntary Ventilation..... | 23 |
| Bronchoprovocation Response..... | 23 |
| <hr/> | |
| Installation | 25 |
| Prepare the microQuark | 26 |
| Software installation | 27 |
| Installing the software | 27 |
| Run the software..... | 27 |
| PC port configuration | 27 |
| Software main features | 28 |
| Display | 28 |
| Tool bar | 28 |
| Show/hide the toolbar..... | 28 |
| Dialog windows..... | 28 |
| Use of the keyboard | 28 |
| Use of the mouse..... | 28 |
| Scroll bars..... | 29 |
| On-line help..... | 29 |
| Software version..... | 29 |
| <hr/> | |
| Calibration | 31 |
| The calibration program | 32 |
| Running the Calibration program | 32 |
| Log file | 32 |
| Turbine calibration | 33 |
| Checking the system signals | 35 |
| The control panel..... | 35 |
| Using the control panel | 35 |
| <hr/> | |
| Database Management | 37 |
| Settings | 38 |
| Graphs | 38 |
| Serial port..... | 38 |
| Units of measurements | 39 |
| Using extra fields | 39 |

| | |
|----------------------------------|-----------|
| Customise the fields..... | 39 |
| Patient's database | 40 |
| Patient Card | 40 |
| Visit Card | 41 |
| Test Card | 41 |
| Import/export a Tests card | 41 |
| Diagnosis Database..... | 41 |
| Archive maintenance | 43 |
| Reorganise the archive..... | 43 |
| Delete the archive | 43 |
| Backup and restore..... | 43 |
| Backup..... | 43 |
| Restore..... | 44 |

Spirometry **45**

| | |
|---|-----------|
| Setting spirometry options | 46 |
| Spirometry | 46 |
| Automatic Interpretation..... | 46 |
| Quality control..... | 47 |
| Parameters manager | 47 |
| Predicted values manager | 48 |
| Predicteds set | 48 |
| Set the current predicted | 49 |
| Formula definition | 49 |
| Page set-up..... | 50 |
| Spirometry tests | 52 |
| Recommendations for spirometry tests..... | 52 |
| Forced Vital Capacity (pre) | 53 |
| Perform a FVC (pre) test..... | 53 |
| Test encouragement | 54 |
| Perform the FVC test with the encouragement | 54 |
| Slow Vital Capacity..... | 55 |
| Perform a SVC test | 55 |
| Maximum Voluntary Ventilation | 57 |
| Perform a MVV test..... | 57 |
| Bronchial Provocation Test | 58 |
| Bronchodilator test..... | 58 |

| | |
|---|-----------|
| Methacholine and Histamine Bronchial provocation Tests | 58 |
| Perform the test | 59 |
| Bronchial Provocation protocols Database..... | 60 |
| Enter a new Bronchial provocation protocol in the archive..... | 61 |
| Viewing results | 62 |
| Tests of the current patient | 62 |
| Delete a test..... | 62 |
| Printing results..... | 63 |
| Printing Reports..... | 63 |
| Printing the active window | 64 |
| Printing a series of reports..... | 64 |
| Electronic reports (*.pdf)..... | 64 |
| Export data | 65 |
| Export a test | 65 |

| | |
|--|-----------|
| System maintenance | 67 |
| System maintenance..... | 68 |
| Cleaning and disinfection | 68 |
| Cleaning the turbine flowmeter | 69 |
| Precautions during the cleaning of the turbine | 69 |
| Suggested disinfection solutions | 70 |
| Inspections..... | 70 |

| | |
|--|-----------|
| Appendix | 71 |
| Service - Warranty..... | 72 |
| Warranty and limitation of liability | 72 |
| Return goods policy for warranty or non warranty repair..... | 73 |
| Repair Service Policy | 73 |
| Privacy Information | 75 |
| Personal data treatment and purposes..... | 75 |
| How your personal data are treated | 75 |
| The consent is optional, but | 75 |
| Holder of the treatment..... | 75 |
| Customer rights | 76 |
| Converting factors configuration | 77 |
| ATS 94 recommendations..... | 78 |

| | |
|---------------------------------------|-----------|
| ATS recommendations | 78 |
| Predicted values..... | 79 |
| ERS93 | 79 |
| Reference Adult:..... | 79 |
| Reference Paediatric:..... | 79 |
| KNUDSON 83 | 79 |
| Reference Adult/ Paediatric:..... | 79 |
| ITS (White race) | 79 |
| Reference Adult/ Paediatric:..... | 79 |
| ITS (Black race)..... | 79 |
| Reference Adult/ Paediatric:..... | 79 |
| LAM | 80 |
| Reference Adult/ Paediatric:..... | 80 |
| Multicéntrico de Barcelona | 80 |
| Reference Adult/ Paediatric:..... | 80 |
| NHANES III | 80 |
| Reference Adult/ Paediatric:..... | 80 |
| Automatic diagnosis (algorithm) | 80 |
| Quality Control Messages | 81 |
| References..... | 82 |



Getting started

Important notices

Intended use

microQuark is an electrical medical device designed to perform pulmonary function tests. It is to be used by physicians or by trained personnel on a physician responsibility.

Caution: Federal law restricts this device to sale by or on the order of a physician.

This equipment has been conceived with the aim of providing an auxiliary instrument allowing:

- the formulation of lung pathology diagnosis;
- important studies concerning human physiology;
- the collection of important information in sport medicine.

No responsibility attaches COSMED Srl for any accident happened after a wrong use of the device, such as:

- use by non qualified people;
- non respect of the device intended use;
- non respect of the hereunder reported precautions and instructions.

Warnings

The device, the program algorithms and the presentation of measured data have been developed according to the specifications of ATS (American Thoracic Society) and ERS (European Respiratory Society). Other international references have been followed when these were not available. All bibliography references are reported in Appendix.

The present handbook has been developed with respect of the European Medical Device Directive requirements which sort microQuark within Class II a.

It is recommended to read carefully the following precautions before putting the device into operation.

The precautions reported below are of fundamental importance to assure the safety of all COSMED equipment users.

1. This user manual is to be considered as a part of the medical device and should always be kept on hand.
2. Safety, measure accuracy and precision can be assured only:

-
- using the accessories described in the manual or given with the device. Actually non recommended accessories can affect safety unfavourable. Before using non recommended accessories it is necessary to get in touch with the manufacturer;
 - ordinary equipment maintenance, inspections, disinfection and cleaning are performed in the way and with the frequency described;
 - any modification or fixing is carried out by qualified personnel;
 - the environmental conditions and the electrical plants where the device operates are in compliance with the specifications of the manual and the present regulations concerning electrical plants. In particular grounding reliability and leakage current suppression can only be assured when the device three – wire receptacle is connected to a yellow - green return connected to earth ground. Attempting to defeat the proper connection of the ground wire is dangerous for users and equipment.
3. Before powering the system, check the power cables and the plugs. Damaged electrical parts must be replaced immediately by authorised personnel.
 4. Cleaning residue, particulates, and other contaminates (including pieces of torn or broken components) in the breathing circuit pose a safety risk to the patient during testing procedures. Aspiration of contaminates can potentially be life-threatening. You must follow all the cleaning procedures in System Maintenance, and you must thoroughly inspect the components after cleaning and before each patient test.
 5. This device is not suitable for use in presence of flammable anaesthetics. It is not an AP nor an APG device (according to the EN 60 601-1 definitions).
 6. Keep the device away from heat and flame source, flammable or inflammable liquids or gases and explosive atmospheres.
 7. In accordance with their intended use microQuark is not to be handled together with other medical devices unless it is clearly declared by the manufacturer itself.

-
8. It is recommended to use a computer with electromagnetic compatibility CE marking and with low radiation emission displays.
 9. It is necessary to make the PC, connected to the microQuark, compliant with EN 60601-1 by means of an isolation transformer.
 10. Graphical symbols used in accordance to present specifications are described here below:



Equipment type B (EN60601-1)



Danger: high temperature



OFF



ON



Protective earth ground



Alternating current

Contraindication

The physical strain to execute the respiratory manoeuvre is contraindicated in case of some symptoms or pathology. The following list is not complete and must be considered as a piece of mere information.

Contraindications for the Spirometer tests

Absolute contraindications

For FVC, VC and MVV tests:

- Post-operating state from thoracic surgery

For FVC tests:

- Severe instability of the airways (such as a destructive bronchial emphysema)
- Bronchial non-specific marked hypersensitivity
- Serious problems for the gas exchange (total or partial respiratory insufficiency)

Relative contraindications

For FVC tests:

- spontaneous post-pneumothorax state
- arterial-venous aneurysm
- strong arterial hypertension
- pregnancy with complications at the 3rd month.

For MVV test:

- hyperventilation syndrome

Contraindications for Bronchial provocation tests

The bronchial provocation tests must be executed according to the doctor's discretion. There are not data that reveal specific contraindication for the bronchial provocation test through inhalation.

The modern standard processes have been revealing secure in several clinical studies. However it is recommendable to respect the following contraindications:

Absolute contraindications

- Serious bronchial obstruction (FEV1 in adults)
- Recent myocardium infarct
- Recent vascular-cerebral accident
- Known arterial aneurysm
- Incapacity for understanding the provocation test procedures and its implications.

Relative contraindications

- Bronchial obstruction caused by the respiratory manoeuvre.
- Moderate or serious bronchial obstruction. For ex. Predicted value FEV1 less than 1.51 in men and predicted value FEV1 in women less than 1.21.
- Recent infection in the superior air tracts
- During the asthmatic re-acuteing
- Hypertension
- Pregnancy
- A pharmacology treatment epilepsy

Environmental condition of use

COSMED units have been conceived for operating in medically utilised rooms without potential explosion hazards.

The units should not be installed in vicinity of x-ray equipment, motors or transformers with high installed power rating since electric or magnetic interferences may falsify the result of measurements or make them impossible. Due to this the vicinity of power lines is to be avoided as well.

Cosmed equipment are not AP not APG devices (according to EN 60601-1): they are not suitable for use in presence of flammable anaesthetic mixtures with air, oxygen or nitrogen protoxide.

If not otherwise stated in the shipping documents, Cosmed equipment have been conceived for operating under normal environmental temperatures and conditions [IEC 601-1(1988)/EN 60 601-1 (1990)].

- Temperature range 10°C (50°F) and 40°C (104°F).
- Relative humidity range 20% to 80%
- Atmospheric Pressure range 700 to 1060 mBar
- Avoid to use it in presence of noxious fumes or dusty environment and near heat sources.
- Do not place near heat sources.
- Cardiopulmonary resuscitation emergency equipment accessible.
- Adequate floor space to assure access to the patient during exercise testing.
- Adequate ventilation in the room.

Safety and conformity

Safety

IEC 601-1 (1988) /EN 60 601-1 (1990);

Find reported below the complete classification of the device:

- Class I type B device
- Protection against water penetration: IP00, ordinary equipment unprotected against water penetration
- Non sterile device
- Device not suitable in the presence of flammable anaesthetics;
- Continuous functioning equipment;

EMC

The system meets the EMC Directive 89/336

EN 60601-1-2

EN 55011 Class B (emission), IEC 1000-4-2, IEC 1000-4-3, IEC 1000-4-4

Quality Assurance

UNI ISO 9001 (Registration n° 387 Cermet)

Medical Device Directive (CE mark)

MDD 93/42/EEC (Notified Body 0476).

Class IIa

Keynotes

Here are the keynotes used to make the manual easier to read.


Typographic keynotes

These are the typographic keynotes used in the manual.

| Style | Description |
|-----------------|---|
| Bold | indicates a control or a key to be pressed. |
| <i>“Italic”</i> | indicates a messages shown by the firmware. |

Graphic keynotes

These are the graphic keynotes used in the manual.

| Illustration | Description |
|---|--|
|  | shows the button to click in the software to activate the related feature. |

Systems Overview

microQuark is an instrument designed for lung function screening; the core of the system is the “intelligent” flowmeter that, connected through the serial port (RS232), turns any Personal Computer (laptop or desktop) in a complete spirometric lab.

The system is composed by the turbine flowmeter, the measurement and data elaboration device (lightweight and ergonomic), the communication cable and by the Software pack.

Before starting

Before operating the microQuark system we strongly recommend to check the equipment and register you as a customer.

Checking the packing contents

Make sure that the package contains the items listed below. In case of missing or damaged parts, please contact Cosmed technical assistance.

microQuark standard packaging

| Code | Qty | Description |
|---------------|-----|------------------------------|
| C00960-01-04 | 1 | microQuark unit |
| C02235-01-05 | 1 | Turbine |
| A 662 100 001 | 2 | Nose clips |
| C01739-02-35 | 1 | PC Software |
| C00137-01-20 | 20 | Paediatric paper mouthpieces |
| C00136-01-20 | 20 | Adult paper mouthpieces |
| C00063-01-20 | 1 | Conic mouthpiece |
| C00214-01-20 | 1 | Paediatric adapter |
| A 362 300 004 | 1 | Serial cable RS232 USB Power |
| C00067-02-94 | 1 | Registration card |
| C02268-02-91 | 1 | User manual |
| C01999-02-DC | 1 | Conformity declaration |

Warranty registration

Before using the system, please take a moment to fill in the registration form and the warranty and return them to COSMED, by doing this you are eligible to the customers assistance service.

For further information, please refer to the enclosed registration and warranty form. If the form is not enclosed in the packaging, please contact directly COSMED.

Register the product via software

Together with the PC software, a registration software is supplied. With this software it is possible to fill in an electronic form with the customer information.

1. To run the software, double click on the icon **Registration** or select **Registration...** from ? menu.
2. Type the requested information and click **Send...** to send the form via e-mail to COSMED.

How to contact COSMED

For any information you may need, please contact the manufacturer directly at the following address:

COSMED S.r.l.

Via dei Piani di Monte Savello, 37

P.O. Box n. 3

00040 - Pavona di Albano

Rome - ITALY

Voice: +39 (06) 931.5492

Fax: +39 (06) 931.4580

email: customersupport@cosmed.it

Internet: <http://www.cosmed.it>

Complain, feedback and suggestions

If you have any complain, feedback information or suggestion, please inform us at complain@cosmed.it.

PC configuration required

- Pentium 133 MHz.
- Windows 95, 98, XP.
- 16 Mb RAM .
- 3.5 drive.
- VGA, SVGA monitor.
- USB port.
- Serial Port RS 232 available.
- Any Mouse and Printer compatible with the MS Windows™ operative system.
- PC conform to European Directive 89/336 EMC

Technical features

| | |
|---------------------|--------------------------------|
| Flowmeter | Bidirectional digital turbine |
| Flow Range: | 0.03 - 20 l/s |
| Volume Range: | 12 l |
| Accuracy: | ± 3% or 50 ml |
| Resistance @12 l/s: | < 0.7 cmH ₂ O/l/sec |
| Mouthpieces: | Ø 31 and Ø22 mm |
| Serial port: | RS232C |
| Dimensions: | 150 x 45 x 53 mm |
| Weight: | 77g |

Measurements

Measured parameters

FVC - Forced Vital Capacity

| Symbol | UM | Parameter |
|-------------|-------|--|
| FVC | l | Forced Expiratory Vital Capacity |
| FEV1 | l | Forced Expiratory Volume in 1 sec |
| FEV1/FVC% % | | FEV1 as a percentage of FVC |
| PEF | l/sec | Peak Expiratory Flow |
| FEV0.5 | l | Forced Expiratory Volume in 0.5 sec |
| FEV6 | l | Forced Expiratory Volume in 6 sec |
| FEV1/FEV6 % | | FEV1 as a percentage of FEV6 |
| FEV6/FVC% % | | FEV6 as a percentage of FVC |
| Best FVC | l | Best Forced Expiratory Vital Capacity |
| Best FEV1 | l | Best Forced Expiratory Volume in 1 sec |
| Best PEF | l/sec | Best Peak Expiratory Flow |
| Vmax25% | l/sec | Expiratory Flow @25% of the FVC |
| Vmax50% | l/sec | Expiratory Flow @50% of the FVC |
| Vmax75% | l/sec | Expiratory Flow @75% of the FVC |
| FEF25-75% | l/sec | Mid-exp flow between 25-75%FVC |
| FET100% | sec | Forced expiratory time |
| FEV2 | l | Forced Expiratory Volume in 2 sec |
| FEV3 | l | Forced Expiratory Volume in 3 sec |
| FEV2/FVC% % | | FEV2 as a percentage of FVC |
| FEV3/FVC% % | | FEV3 as a percentage of FVC |
| FEV1/VC% % | | Tiffenau index |
| FEF50-75% | l/sec | Mid-exp flow between 50-75%FVC |
| FEF75-85% | l/sec | Mid-exp flow between 75-85%FVC |
| FEF0.2-1.2% | l/sec | Mid-exp flow between 0.2 l - 1.2 l |
| FiVC | L | Inspiratory Forced Vital Capacity |
| FiF25-75% | l/sec | Forced mid-inspiratory flow |
| FiV1 | l/sec | Forced Inspiratory Volume in 1 sec |
| PIF | l/sec | Peak Inspiratory Flow |
| VEXT | ml | Extrapolated Volume (back extrapolation) |
| PEFT | msec | Time to PEF (10% - 90%) |

VC/IVC - Slow Vital Capacity and Ventilatory pattern

| Symbol | UM | Parameter |
|---------------|-----------|-----------------------------------|
| EVC | l | Expiratory Vital Capacity |
| IVC | l | Inspiratory Vital Capacity |
| ERV | l | Expiratory Reserve Volume |
| IRV | l | Inspiratory Reserve Volume |
| IC | l | Inspiratory Capacity |
| VE | l/min | Expiratory Minute Ventilation |
| Vt | l | Tidal Volume |
| Rf | 1/min | Respiratory Frequency |
| Ti | sec | Duration of Inspiration |
| Te | sec | Duration of Expiration |
| Ttot | sec | Duration of Total breathing cycle |
| Ti/Ttot | — | Ti/Ttot ratio |
| Vt/ti | l/sec | Vt/ti ratio |

MVV - Maximum Voluntary Ventilation

| Symbol | UM | Parameter |
|---------------|-----------|-------------------------------|
| MVV | l/min | Maximum Voluntary Ventilation |
| MVt | l | Tidal Volume (during MVV) |
| MRf | 1/min | Maximum Respiratory frequency |
| MVVt | sec | MVV duration time |

Bronchoprovocation Response

| Symbol | UM | Parameter |
|---------------|-----------|--|
| FallFEV1 | % | Fall in FEV1 from baseline or post diluent |
| FallVmax50% | % | Fall in Vmax50% from ref. |
| P10 | — | Provocative dose causing FEV1 to fall 10% |
| P15 | — | Provocative dose causing FEV1 to fall 15% |
| P20 | — | Provocative dose causing FEV1 to fall 20% |



Installation


Prepare the microQuark

Power supply is provided by PC through the keyboard/mouse PS2 plug. The unit is turned on/off automatically by the PC.

1. Insert the turbine and be sure to push the turbine up to touch the end of the holder.




2. Connect the microQuark to the PC, through the RS232 (data) and USB (power supply) port, as shown in the following picture.



Software installation

Installing the software

1. Select **Run...** from Windows **Start** menu.
2. In the Command line box, type **a:\install** (assuming the disk is in drive A:).
3. Click on **OK** (or press **ENTER** key).
4. The program will load up a dialog box and ask for a directory where to be installed.
5. The installation software will present a dialog box in which you have to select the device. Select **microQuark**.



6. When the installation is over, the program will advise you with a message indicating that the installation has been successfully completed, click on **End**.

Notice: The software is copy-protected. Install the software from the original disk.

Run the software

1. In the Windows **Start** menu, open the Program Group in which the software was installed.
2. Click the **Quark PFT** icon.

PC port configuration

The first time the software is used, it is necessary to configure the communication port with the PC (USB, COM1, COM2,...).

For further details, see the chapter *Database management*.

Software main features

Display

The program may contain several windows. The active window is highlighted with a different colour of the caption. Some functions of the program are "active window" sensitive (Print, right key of the mouse).

Tool bar

Many of the functions that may be selected from the menu can be activated more rapidly by clicking with the mouse on the corresponding icon in the tool bar.

Positioning the mouse cursor on one of the buttons of the toolbar (if the option Hints is enabled), the description of the corresponding function is shown in a label.

Show/hide the toolbar

Select **Toolbar** from **Options** menu in order to show or hide the toolbar.

Dialog windows

The typical operating environment of Microsoft Windows is the Dialog box. This window is provided with a series of fields in which input the information.

Use of the keyboard

- To move the cursor among fields, press the **Tab** key until you reach the desired field.
- Press the **Enter** key to confirm the information input on the dialog box or press the **Esc** key to cancel changes.

Use of the mouse

- To move the cursor among fields, move the mouse on the desired field and left-click.
- Click on the **OK** button with the Left button of the mouse to confirm the information input on the dialog box or click on **Cancel** button to cancel changes.

Scroll bars

Some windows are provided with scroll bars that help to see data exceeding the window space available.

- To move the scroll bar row by row click the scroll arrows at the end of the scroll bars
- To move the scroll bar page by page click on the grey area at both sides of the scroll fields

On-line help

COSMED Help is a complete on-line reference tool that you can use at any time. Help is especially useful when you need information quickly or when the user manual is not available. Help contains a description of each command and dialog box, and explains many procedures for accomplishing common tasks.

To get the Help on line, press the **F1** key.

Software version

To know the software version and the serial number of the software, select **About Quark PFT...** from **Help** menu.



Calibration

The calibration program

Running the Calibration program



Start the program and choose **Calibration** from the **Test** Menu. The software runs the Calibration software and the main menu changes accordingly.

Log file

The program creates and updates as default the calibration log file, containing the conditions and the results of all the calibrations performed by the user.

To access the file select **File/Report File...** from the calibration program.

Turbine calibration

microQuark is calibrated by COSMED. ATS recommends a daily calibration of the turbine. However if it is correctly maintained, turbine retains its precision for longer periods. We advice to calibrate the turbine daily to detect malfunctioning.

Note: if you are using a slow PC, we recommend to set an higher refresh time.

In order to calibrate the turbine:




1. Connect the turbine flowmeter to the calibration syringe.
2. Select **Calibration** from **Test** menu.
3. Select **Reference values** from the **Calibration** menu and enter the syringe volume, if different to the displayed one.
4. Select **Turbine** from the **Calibration** menu.
5. When the **Calibration Turbine** dialog box appears with the syringe piston initially pushed all the way in, move the piston in and out for 5 inspiratory strokes and 5 expiratory strokes in order to get the first values appearing on the screen. Then move the syringe piston for other 10 strokes (IN and EX).
6. At each of the 10 steps the software displays the results of the manoeuvre and the percentage error in the reading.

| Results | | | | | |
|---------|-------|------|------|-------|------|
| Exp. | % | Gain | Ins. | % | Gain |
| 3002 | +0.07 | 1046 | 2993 | -0.23 | 1027 |
| 2985 | -0.50 | 1048 | 2995 | -0.17 | 1028 |
| 2972 | -0.93 | 1052 | 3010 | +0.33 | 1026 |
| 2993 | -0.23 | 1052 | 2993 | -0.23 | 1027 |
| 3019 | +0.63 | 1051 | 2984 | -0.53 | 1027 |
| 2992 | -0.27 | 1051 | 2990 | -0.33 | 1028 |
| 2994 | -0.20 | 1051 | 3004 | +0.13 | 1027 |
| 3009 | +0.30 | 1050 | 3008 | +0.27 | 1027 |
| 2988 | -0.40 | 1051 | 2996 | -0.13 | 1027 |
| 2974 | -0.87 | 1051 | 3006 | +0.20 | 1027 |

Move the calibration syringe...

Cancel Help

-
7. At the end of this operation, the software displays the new calibration factors. Press **OK** to store the new value.



Turbine calibration results [X]

| | | | |
|------------|-----------------------------------|---------------------------------------|--|
| Type (mm): | <input type="text" value="26"/> | <input type="button" value="OK"/> | <input type="button" value="Default"/> |
| Gain Exp.: | <input type="text" value="1000"/> | <input type="button" value="Cancel"/> | <input type="button" value="Help"/> |
| Gain Ins.: | <input type="text" value="1000"/> | | |

The 3 litres calibration syringe can be purchased directly from COSMED (P/N: C00600-01-11).


***Note:** If a bacterial filter is used for the tests, do use it also during the turbine calibration.*

Checking the system signals


The control panel

The **Control Panel**, which can be activated from the **Calibration/Control panel...** menu item, is a useful tool to check the main hardware functions of microQuark.

By using the controls on Control Panel you are able to read the signals acquired by the system both as voltages and processed data.



Using the control panel






Database Management

Settings

The software allows to configure some options selecting **Configure** from the **Option** menu.



Graphs

All the graphs visualised and/or printed can be customised in colours and appearance.

1. Select the desired colours of the curves (5 curves max can be overlapped on the same graph).
2. Enable or disable the **Grid** option.
3. Enable or disable the **Show Info Test** option.

Serial port

You must select the serial port RS 232 that will be used to connect the microQuark with the PC.

To select the serial port, click on the proper **COM** button (the selected port must be different from the mouse one).

Units of measurements

It is possible to configure the units of measurements, weight and height, for printing and viewing.

To select the units of measurement click on **cm/Kg** or **in/lb** according to the desired format.

Using extra fields

The Patient's database is organised in 3 different cards (Patient card, Visit Card and Test card.) where it is possible to store the information about patients and visits .

Besides the standard information, it is possible to customise some fields (user free fields), entering and labelling measurements coming from other devices.

The customisable free fields are:

- 3 fields in the Patient Card (Patient's information)
- 3 fields in the Visit Card (information about the visits)
- 3 fields (2 numeric) in the Test card information about Test)

Customise the fields

In the group **User free fields** type the desired text in the 9 fields available.

Patient's database

The Patients database consists of a Patient Card, a Visit Card and a Test Card in which are listed all tests performed by the patient.

Select **Archive Navigator** from the **File** menu or press the button by side.

The Archive Navigator

Patient Card

Name 1001 Last 1001
ID 1001 Male Birth 01/01/1975

Visit Card

Date: 09/10/1995
Height (cm) 180. Weight (Kg) 70. Age 20

Smoke
 No Yes Ex Years 0 Cig/Day 0

Description
Diagnosis
Notes

Tests to Display

09/10/1995 - FVC 6.11
FEV1 4.41
PEF 6.85

Drug DRUG
Diag Borderline Obstruction

Quality Control
Blow out longer (FET100% < 6 sec), chec

OK Help



Note: after having deleted a record (patient, visit or test), it is recommended to reorganize the archive in order to free disk space.

Patient Card

It collects all the information of a patient (first name, last name, date of birth) which remain the same for each visit. For each patient there is only one Patient Card, which is created the first time the Patient performs a test.

To move within the database use the following buttons:



Move to the first patient in the archive



Move to the previous patient in the archive



Move to the next patient in the archive



Move to the last patient in the archive



Find a patient in the archive



Enter a new patient in the archive



Delete current patient from the archive



Edit the current patient card

Visit Card

It collects all information relative to the visit (diagnosis, visit description...) and to the patient information subject to change between one visit and another (height, weight, smoke). Each patient can be related to several Visit Cards provided they have been created in different days. Before carrying out any spirometric test it is necessary to create a new Visit Card or to open the today's Visit Card.

To move within the database use the following buttons:



Move to the first visit in the archive



Move to the previous visit in the archive



Move to the next visit in the archive



Move to the last visit in the archive



Find a visit in the archive



Enter a new visit card in the archive



Delete current visit card from the archive



Edit the current visit card

Test Card

It contains all the information about the test.

To move within the database use the following buttons:



Delete current test from the archive.



Edit the current test

Import/export a Tests card

This function allows to import /export a test card with the respective visit and patient card.

1. Select the patient.



2. Choose the tests to be exported and press the key by side. All data will be imported/exported in the XPO file format (Cosmed proprietary).

Diagnosis Database

The program allows to manage a diagnosis database, whose records are composed by a diagnosis ID code and a string of text.

The report of the visits can be done either by typing the desired text in the field “Diagnosis” of the Visit Card or, more quickly, retrieving from the diagnosis database the desired one.

If you want to insert, modify or delete a diagnosis from the database select **Database Diagnosis...** from the **File** menu.

Archive maintenance

The software allows to manage files selecting **Archive** from the **File** menu.

It is advisable to perform the archive reorganisation every month, in order to free space on the hard disk and/or to correct possible errors present within the database.

It is possible also that your have no more hard disk space. So, you have to delete all the data. In this case, it is useful to perform the initialising.

Reorganise the archive

1. Select **Reorganize archive** from the **File** menu.
2. Wait for the end of the operation before performing any other function.

Delete the archive


1. Select **Initialize Archive** from the **File** menu.
2. Wait for the end of the operation before performing any other function.

Backup and restore

It is strongly recommended to backup files, a warning message will be displayed monthly. This function allows the user to restore the data if the PC or the HD will not work anymore.

Backup

1. Select **Backup archive** from the **File** menu.




2. Selecting the destination path with the **Browse** key or press **New** to create a new directory. Press **OK** to confirm.

-
3. In the dialog box it will appear an estimate of the number of floppy you need in order to back up the archives. Press **OK**.

Restore

1. Select **Restore archive** from the **File** menu.




2. On the **Restore** dialog box specify the drive source and press **OK**, a dialog box will appear indicating all data of the backup processed.

Spirometry

Setting spirometry options

The software allows to configure some options selecting **Configure** from the **Option** menu.

Spirometry



Automatic Interpretation

microQuark has the function of interpreting each test performed by a patient visualising an automatic diagnosis. The algorithm has been calculated basing on “Lung Function Testing: selection of reference values and interpretative strategies, A.R.R.D. 144/1991:1202-1218”. The automatic diagnosis is calculated at the end of the FVC if:

- the automatic diagnosis option is enabled.
- the patient’s anthropometric data allow the calculation of the LLN (Lower Limit of Normal range).
- at least one FVC test has been performed.

To enable/disable the automatic diagnosis:


1. Click on **Enable Automatic Interpretation**.
2. Select the LLN (Lower Limit of Normal Range) criteria among the ATS ($LLN = Pred - 0.674 * SD$), ERS ($LLN = Pred - 1.647 * SD$) or 80%Pred ($LLN = Pred * 0.8$) specifications.

Quality control

microQuark allows a quality test control. The calculation has been carried out referring to “Spirometry in the Lung Health Study: Methods and Quality Control, A.R.R.D. 1991; 143:1215-1223”. The messages concerning the quality control are shown at the end of the test.

To enable/disable the quality control, click on **Enable Quality Control** checkbox.

Parameters manager



The program allows to calculate a huge number of parameters; it is advisable, in order to simplify the analysis of the results, to view, to print and to sort the desired parameters only. Select the menu item **Options/Parameters...**

View

Move the parameters to view into the *Selected parameters* list.

Print

Move the parameters to print into the *Selected parameters* list.

Sort

Drag the parameter up or down with the mouse.


Customise

Add, modify and delete custom parameters.



If it is necessary to restore the default parameters press the button in the left corner of the window to initialise the parameters database.

Predicted values manager



The program contains a preset of predicted equations, but the user is allowed to customise its own predicted sets. Select **Predicteds...** from **Options menu**.

The window is divided into two forms: **Predicteds set** and **Formula definition**.

Predicteds set

This form allows the user to manage the set of predicted. The following information define a set:

Name: identifies the set and cannot be duplicated;

Description: free field;

Age: the adult predicted start since this age.

To enter a new set of predicted click on the **New** button. The field **Name** must be filled and must be unique. To stop without saving click on the **Cancel** button. To save the set, click on the **Save** button.

To delete a set of predicted click on the **Delete** button. If a set is deleted, also the associated formulae are deleted.

It is possible to generate a new set of predicted with the same attributes and the same formulae of the selected one. To do this click on the **Copy...** button and specify a new Name.

To import a set of predicted click on the **Import...** button and select a file of Predicteds files type.

To export a set of predicted click on the **Export...** button.

In the list **Set current predicted**s choose the current predicted for printing and viewing.

Set the current predicted

microQuark allows to calculate the predicted values according to the following configurable sets:

| Adults | Paediatrics |
|---------------|--------------------|
| ERS 93 | Zapletal |
| Knudson83 | Knudson83 |
| ITS white | ITS white |
| ITS black | ITS black |
| LAM | LAM |
| MC Barcelona | MC Barcelona |
| Nhanes III | Nhanes III |

Select the desired choice in the group **Predicted**.

Formula definition

The screenshot shows the 'Predicteds' dialog box with the 'Formula definition' tab active. The 'Predicteds set' dropdown is set to '232'. The 'Description' field is empty. The 'Use the predicted's formulae:' radio button is selected. The '...of the customized formulae:' section has 'Male' selected. There are two rows of input fields: 'Young' and 'Adult', each with a 'Formula' field and a 'Standard Deviation' field. At the bottom are buttons for 'Copy', 'Paste', 'Parameter...', 'Save', 'Delete', 'OK', 'Cancel', 'Apply', and 'Help'.

This form allows the user to manage the formulae associated to a set of predicted's.

Select the set of predicted's from the list **Predicted**s set.

To insert a new parameter click on the **New...** button.

The parameter formulae can be:

- calculated according to the predicted's in the list **Use the predicted's formulae**;


- customised by the user with the option **...or the customised formulae**.

The **Delete** button deletes the selected parameter.

The **Copy** button stores the selected parameter in memory.

The **Paste** button inserts a new parameter from the one copied. If the name is not unique, the user is asked whether to specify a new name or to replace the existing parameter.

Page set-up



Select **Page Setup...** from the **File** menu.

Header All the printouts carried out by the program are preceded by 3 rows of customisable header (usually they contain the name and the address of the Hospital using the spirometer).


Data Patient and visit information are printed below the header. These data are reported on 3 columns and 5 rows. the user may configure the disposition, change and eventually cancel the fields, as he prefers.

Margins Configures the print margins from the borders of the paper. The unit of measure is decided in **Units of measurements**.

Footer Configures information at the bottom of the page.

Printed file name

Defines the automatic name to be assigned to the pdf file, if the report will be printed in this format.



In the example it has been set to create a filename composed by <date of the test> followed by <last name> and <first name>.

Spirometry tests



Note: Read carefully the contraindications in Chapter 1.

Once completed the phases of the introduction of the patient's data and the visit data, it is possible to carry out the spirometric tests.

microQuark allows to perform the following tests:

| Key | Test |
|------------|---|
| FVC pre | Forced Vital Capacity |
| FVC post | Forced Vital Capacity after bronchial stimulation |
| SVC | Slow Vital Capacity |
| MVV | Maximum Voluntary Ventilation |

Before performing any test make sure that:

1. microQuark is properly connected to your PC and the selected serial port (COM1, COM2) corresponds to the one effectively use.
2. The name shown on the status bar corresponds to the patient who is to carrying out the tests.
3. The today's visit card exists.

Recommendations for spirometry tests

- The patient should wear the nose clips
- The turbine has been recently calibrated (ATS recommends a daily calibration)
- The paper mouthpiece or the antibacterial filter is properly connected to the flowmeter through the corresponding adapter

For hygienic reasons, we strongly recommend the use of a bacterial filter.

If a kid must perform the test it is recommended to enable the encouragement function which shows exactly the manoeuvre of the FVC test.

Forced Vital Capacity (pre)

FVC is a reference test to verify obstructive (airflow limitations) and restrictive disorders (lung volume limitations). To achieve good test results it is fundamental a good manoeuvre (quality control messages, real time plots ...)

The main parameters measured during FVC tests are:

| | |
|-----------|--------------------------------------|
| FVC | Forced Vital Capacity |
| FEV1 | Forced Expiratory Volume in 1 second |
| FEV1/FVC% | FEV1 as a percentage of FVC |
| PEF | Peak Expiratory Flow |
| FEF25-75% | Forced mid-Expiratory Flow |

The two representative plots are the Flow/Volume and Volume/Time loops.

By comparing FVC, FEV1 and FEV1/FVC% values the software allows an automatic interpretation concerning the levels of obstructive and/or restrictive disorders.

Perform a FVC (pre) test



1. Select **Forced Vital Capacity pre** from the **Test** menu and wait for the green led is prompted on the right side of the screen.
2. Explain the manoeuvre to the patient and press the **F2** key.
3. Wait some seconds and perform the test.
4. After having performed the test, press **F3** or wait for the automatic end (5 seconds without flow), so that the software displays the F/V and V/t graphs, the main parameters, and the predicted values.
5. In order to visualise the F/V and V/t graph and the main parameters press the following buttons:



view Flow Volume graph



view Volume Time graph



view data of the test



6. Press **Alt+F3** to stop the acquisition discarding the results.
7. Repeat the test until it is correctly performed (ATS recommends 3 times).




- Press **Exit** to visualise the test list carried out during current session together with the results of the main parameters.
- Select the test that you want to save (the Best Test according to the ATS criteria is highlighted as default) and press **OK**.

Test encouragement

During FVC manoeuvre you might experience some lack of collaboration with kids or with other patients. In this case you may find a good help in using the encouragement software tool.

Perform the FVC test with the encouragement

- Select **Encouragement** from **View** menu.
- Perform the test as explained in the previous paragraph.



Slow Vital Capacity

Important test for assessing COPD (chronic obstructive pulmonary disease) patients affected by this disease might present a the Slow Vital Capacity could be higher than the Forced one (FVC).

The main parameters measured during SVC tests are:

- EVC Expiratory Slow Vital Capacity
- IVC Inspiratory Slow Vital Capacity
- ERV Expiratory Reserve Volume
- IRV Inspiratory Reserve Volume

If the inspiratory/expiratory maximal manoeuvre is preceded by a some breaths at tidal volume the software allows to measure the Respiratory Pattern, represented by the following parameters:

- VE Ventilation per minute
- Vt Tidal volume
- Rf Respiratory frequency
- Ttot Breath time
- Ti/Ttot Inspiratory time/Ttot
- Vt/Ti Vt/Ti

Perform a SVC test



1. Select **Slow Vital Capacity** from the **Test** menu and wait for the green led is prompted on the right side of the screen.



2. Press **F2** and instruct the Patient to breath normally until the message “carry out...” is prompted; then ask to perform a Slow Vital Capacity (deep inhalation, maximal slow expiration and deep inhalation again).



3. Press **F3** or wait for automatic interruption (5 seconds without flow) in order to visualise the V/t graph together with the main parameters compared to the predicted values

4. To visualise the V/t graph and the main parameters press the follow buttons:



view Volume Time graph



view data of the test



5. Press **Alt+F3** to stop the acquisition discarding the results.



6. Repeat the test until it is correctly performed (ATS recommends 3 times).
7. Press **Exit** to visualise the test list carried out during current session together with the results of the main parameters.
8. Select the test that you want to save (the Best Test according to the ATS criteria is highlighted by default) and press **OK**.

The reference for the ERV calculation is displayed on the V/T graph.

Maximum Voluntary Ventilation

Test for assessing the maximum ventilatory capacity. In the past, it was commonly performed during routine PF tests, however its clinical use declined over the years. Today MVV test is most commonly performed as part of the exercise tolerance tests, where it is used as an index of maximum ventilatory capacity. Test consists in breathing in and out deeply and rapidly for 12, 15 seconds. The expired volume during this short period is then extrapolated

The most important measured parameter is the following:

MVV Maximum Voluntary Ventilation

Perform a MVV test



1. Select **Maximum Voluntary Ventilation** from the **test** menu and wait for the green led is prompted on the right side of the screen.



2. Press **F2** and make the Patient breath as much deeply and rapidly as possible for at least 12 seconds.



3. Press **F3** or wait for automatic interruption (5 seconds without flow) in order to visualise the V/t graph together with the main parameters compared to the predicted values

4. To visualise the V/t graph and the main parameters press the follow buttons:



view Volume Time graph



view data of the test



5. Press **Alt+F3** to stop the acquisition discarding the results.

6. Repeat the test until it is correctly performed (ATS recommends 3 times).



7. Press **Exit** to visualise the test list carried out during current session together with the results of the main parameters.

8. Select the test that you want to save (the Best Test according to the ATS criteria is highlighted as default) and press **OK**.

Bronchial Provocation Test

Bronchodilator test



Note: Read carefully the contraindications in Chapter 1.

Bronchodilators are administered routinely in the PFT laboratory to determine whether airflow obstruction is reversible. Bronchodilators increase airway calibre by relaxing airway smooth muscle.

The test consists of comparing results between the reference FVC (FVC PRE) and the FVC POST performed after the administration of the drug. Increasing value of 13-15% of FEV₁, respect to the basal value (FVC Pre) is considered as a reversible condition.

Main parameters are the following:

DFEV₁%pre Change of FEV₁ as a percentage of test PRE

DFVC%pre Change of FVC as a percentage of test PRE

DPEF%pre Change of PEF as a percentage of test PRE

Some authors states that the above mentioned parameters are too dependent from the FVC Pre, hence latest reference (ERS93, [A comparison of six different ways of expressing the bronchodilating response in asthma and COPD; reproducibility and dependence of pre bronchodilator FEV₁: E. Dompeling, C.P. van Schayck et Al; ERJ 1992, 5, 975-981]) recommend the following parameters:

DFEV₁%pred Change of FVC as percentage of predicted value

DFEV₁%poss Change of FEV₁ as percentage of possible value

Methacholine and Histamine Bronchial provocation Tests

The most common indication for performing methacholine and histamine bronchial challenges is to diagnose hyperresponsive airways. Some patients demonstrate normal baseline pulmonary function despite complaints of “tightness” wheezing, cough, and a little or not response to bronchodilator. Other patients demonstrate spirometric improvement after use of bronchodilator have diurnal variation in peak flows. In this groups aerosolised bronchial challenges are used to confirm a diagnosis of Asthma.

We can summarise the use of the test as follows:

1. Diagnose asthma
2. Confirm a diagnosis of asthma

-
3. Document the severity of hyperresponsiveness
 4. Follow changes in hyperresponsiveness

When patients with hyperresponsive airways inhale certain pharmacologic agents (i.e. Methacholine or histamine) the airways respond by constricting.

Test consists of executing repeated FVC following the pharmacologic agents inhalation according to an established protocol. The fall of the FEV1 parameter is used to calculate the bronchial hyperresponsiveness. The most important parameter is the PD20 that is amount of drug (mg/ml) that causes a reduction of 20% of the FEV1 respect the basal value (without drug).

Main parameters are:

- P10 Dose that causes a 10% fall of FEV1.
- P15 Dose that causes a 15% fall of FEV1.
- P20 Dose that causes a 20% fall of FEV1.

The representative plot is the *Dose/response curve*, showing the percentage variation of FEV1 versus the Drug dose in logarithmic scale.

The program assumes as the **baseline test** the best **FVC pre** carried out during the today's visit. You can change the reference pre test editing the **Post** test.

The name of the drug, its quantity and its unit of measurement, can be typed immediately before any **FVC post** manoeuvre (manual protocol) or can be stored in a database of bronchoprovocation (**File/Bronchial Provocation protocols Database...**).

Perform the test

(During 1st step only) select **Protocol...** from the **Test** menu and choose the name of the bronchoprovocation protocol that you are going to use (**manual protocol** if you want to type the information about the agent before any manoeuvre)



1. Select **FVC post** from the **Test** menu.
2. Select an existing protocol or click on "manual protocol", and wait the green leds turned on.
3. Press **F2**, or the button by side, to start the test.
4. Press **F3**, or the button by side, to achieve the test.

-
- In order to visualise the V/t graph and the main parameters press the follow buttons:



view Flow Volume graph



view data of the test



view bronchial provocation response



- Press **Alt+F3** to stop the acquisition discarding the results.
- Repeat the test until it is correctly performed (ATS recommends 3 times).
- Press **Exit** to visualise the test list carried out during current session together with the results of the main parameters.
- Select the test that you want to save (the Best Test according to the ATS criteria is highlighted as default) and press **OK**.



Bronchial Provocation protocols Database


The response to a bronchoprovocator is usually assessed in terms of change in the FEV₁, vital capacity or airways resistance on the basis of serial measurements (FVC manoeuvres) in which the results of the initial test constitute the reference values. The international literature proposes several standardised protocols in order to address the methodological issues of the various available techniques.

The possibility to store a bronchoprovocation protocol in a database is useful to simplify and automate the sequence of operations that the Physician need to execute during the bronchoprovocation tests.

The typical sequence of activities to carry out a bronchoprovocation test are:

- Typing and storing a bronchoprovocation protocol in the database (usually only once).
- Selection of protocol among the list of the ones already present in the database before carrying out the FVC post tests (the selection of “manual protocol” allows to execute the test fully manually).
- Performing the Post tests.

Enter a new Bronchial provocation protocol in the archive

1. Select **Bronchoprov. protocols database** from the **File** menu.
2. Type the Protocol name, the Bronchoprovocator name and the unit of measurement in the proper input fields.
3. If the bronchoprovocator has a cumulative effect select the cumulative check button.
4. Enter the quantities for each step and press the button .

Viewing results

All the visualisation functions refer to the test carried out by the Current Patient, whose name is indicated on the left-side of the status bar.

To view tests results:



1. Select the **Patients** from the **File** menu
2. Select the patient corresponding to the test you want to view.
3. Select in the list box of the tests up to 5 tests of the kind (FVC, VC/IVC, or MVV) and press **OK**.

To switch between graph and or data use the following buttons on the toolbar:



view Flow Volume graph (F5)



view Volume Time graph (F6)



view data of the test (F7)



view bronchial provocation response.

If you need more than one visualisation meantime use the **New Window** function from the **Window** menu.

If you need to display a list of visits:

- Select **Visits list...** from the **File** menu.
- Type the name of the Company and/or the time interval desired or simply confirm for the complete list.

Tests of the current patient

If a **current patient** has been selected you can quickly view his tests selecting **Test current patient...** from the **View** menu.

Delete a test



1. Select **Patients** from the **File** menu or press the button by side.
2. Select the test that you want to eliminate from the list of the tests referred to the Current Patient and press **Delete**.

Printing results

You can print out in three different ways:


- printing the Report
- printing the Active Window
- printing a series of reports

Printing Reports



To print a report of the current visit, select **Print report...** from **File** menu. The software will choose automatically the best performed test.

The standard Report is composed by 1, 2 or 3 pages depending if you wish to printout the FVC data and the graphs together on the first page or if you wish to printout the bronchoprovocation response.



Print Report - Choose what to print

FVC graph Multi-breath N₂ Wash-out
 Response Single-breath O₂
 Single-breath CO diffusing capacity
 Single-breath CO ... (no breath hold)
 CO diffusing capacity Steady State
 Respiratory drive

FVC POST tests:

| Test # | Drug | Dose | FEV1 |
|--------|--------------|------|------|
| 7 | Methacholine | 0.08 | 5.58 |
| 8 | Methacholine | 0.23 | 5.13 |
| 9 | Methacholine | 0.54 | 5.05 |
| 10 | Methacholine | 1.16 | 4.56 |
| 11 | Methacholine | 2.41 | 4.28 |
| 12 | Methacholine | 4.91 | 3.80 |

ATS
 One page (no ATS) Preview

Buttons: OK, Cancel, Help

- Selecting the option **One page (no ATS)** the report will contain, on one page, the F/V and V/t graphs of the best test, overlapped on the **FVC Post**, the patient data, the notes, the diagnosis and the test results.
- Otherwise the report will contain two pages, the first with the patient data, the graphs and the diagnosis, and the second one with the measured parameters, according to the ATS recommendations.
- The 3rd page will contain the bronchoprovocation response.

Select the desired options:

- FVC graph** Prints the F/V and V/t curves for the best FVC test.
- One page (no ATS)** Prints data and graphs on the first page.
- Response** Prints the bronchoprovocator response.
- FVC post** Prints data and graphs for the Post FVC test (the test can be selected among the test performed in the current visit).
- Preview** Views a report preview on the screen.

Printing the active window



This printout function is only enabled when the active window (title bar highlighted) is one of the following objects:

- Any kind of Graph.
- Numeric data
- List of visit

To print the active window select **Print Active window** from **File** menu.

Printing a series of reports

Sometimes it is useful to printout automatically a series of reports (all tests carried out with the employees, all tests carried out in the today's session).

To print out proceed as follows:

1. Select **Visit List** from the **File** menu
2. Set the criteria of the visits to be added in the list (from, to,...)
3. Select **Print Report** from the **File** menu.

Electronic reports (*.pdf)

If an Adobe PDF writer "Printer Driver" is installed and set as the default printer, it is possible to store the printout report automatically in any location of the HD or eventually LAN paths according to a customizable filename format.

It is possible to define the created filename format selecting **File/Page Set up...** (see Page set-up).

Export data

With this function you can export the test data in 4 different formats:


- *.txt (ASCII)
- *.xls (Microsoft Excel)
- *.wk1 (Lotus 123)
- *.xpo (Cosmed)

Export a test

1. Select **Export tests** from the **File** menu.
2. Select the test to export from the list box and press **OK**.
3. Type the name and the format of the file in the dialog **Save as**. If the ASCII format is selected, the Text button in the dialog box Save as allows you to configure the separators for character based files.

With the *.xpo Cosmed file format it is possible to import data from another Quark archive. Press **OK** to confirm.

4. Select the folder for the export and type the file name. Press **OK** to confirm. A status bar will show the file creation.



System maintenance

System maintenance

All service operations which are not specified in this user manual should be performed by qualified personnel in accordance with the service handbook (to be required to the manufacturer).

All materials used in the construction of the microQuark are non toxic and pose no safety risk to the patient or operator.

Prior to the device cleaning, disinfection and inspection it is necessary to switch off the device and to disconnect adapters from the supply mains.

In order to guarantee the highest accuracy of measurements we recommend you to disinfect the turbine periodically.

Cleaning and disinfection

Cleaning and disinfecting instructions are of fundamental importance to control infections and assure patient safety. In fact aspiration of residue, particles and contaminated agents are life – threatening.

In this handbook we strongly recommend you to follow the rules worked out by ATS and ERS (see: "Lung Volume Equipment and Infection Control" – ERS/ATS WORKSHOP REPORT SERIES, European Respiratory Journal 1997; 10: 1928 – 1932), which are summarised and adapted for COSMED products as follows:

- Accessible internal as well as external surfaces of equipment exposed to expired gas should be washed and disinfected prior to testing of subsequent patients.
- Disinfecting should ideally be performed by heat sterilisation, but gas or liquid sterilisation can be used if the equipment is well cleaned first (no droplets of saliva/sputum remain).
- Disposable gloves should be worn when handling mouthpieces, when cleaning equipment exposed to saliva or sputum and especially when drawing blood.
- Laboratory staff should wash hands prior to testing of each patient.
- Adopt particular precautions when testing patients with recognised high – risk communicable diseases (e.g. tuberculosis, multidrug – resistant staphylococcus). In these


cases, the clinical need for such testing should justify the risks.

During the disinfection:

- do not use alcohol or other liquids containing glutaraldehyde on the exterior surfaces of the equipment. Actually they can damage polycarbonates plastics and may produce unhealthy substances.
- do not use abrasive powders or glass cleaners containing alcohol or ammonia on the plexiglas components of the equipment
- do not steam autoclave any parts of the equipment unless it is clearly specified.
- do not immerse the optoelectronic reader.



Warning: Do not use alcoholic solutions for the turbine, otherwise there can be damages to the plastic material.



Cleaning the turbine flowmeter

It is necessary to disinfect periodically the turbine for sanitary measures or/and for the correct device function.

As disinfecting solution it is suggested Sodium hypochlorite 5% (bleach).

The disinfecting procedure is easy and may be effected every time the user needs, keeping attention to some precautions:

1. Take out the turbine.
2. Dip it in a disinfectant solution (non alcoholic based) for about 1 hour.
3. Rinse the turbine in a vessel, filled of clean water, shaking gently to remove the disinfectant (do not clean the turbine by putting it under running water!).
4. Let it dry to air.
5. After cleaning the turbine, check if the turbine propeller rotates freely even with a low speed air flow.
6. Connect the turbine to the reader.

Precautions during the cleaning of the turbine

1. Do not expose the turbine to high heat and do not put it under running water.
2. Do not ever dip the optoelectronic reader in any kind of solution, the liquid infiltration would damage the internal circuit.
3. Do not use alcoholic solutions to clean the turbine.

Suggested disinfection solutions

| | |
|----------------|---------------------|
| Helipur H Plus | Braun Melsungen AG |
| Gigasept FF | Schulke & Mayr GmbH |
| Dismozon pur | Bode Chemie GmbH |
| TETA-S | Fresenius AG |
| CIDEX | Johnson & Johnson |

Inspections

The equipment requires easy inspections to be carried out in order to assure a proper electrical and mechanical safety level in the years.

These inspections are highly recommended after a rough use of the equipment or after a period of storage in unfavourable environmental conditions.

Referring to the electrical safety, it is important to check the conditions of insulation materials of cables, plugs and any other visible part by means of simple inspection, when the equipment is switched off and adapters (or electrical feeders) are disconnected from the supply mains.

Extract the turbine from the unit and verify, by inspection, that the turbine axis fits correctly its seats and the blade is strongly fastened on the axis itself (it can be useful to shake slightly the turbine in order to note any anomalous movement).

Check if there are any torn or broken components in the breathing circuits: remember that they can create safety risk to patients during tests.



Appendix

Service - Warranty

Warranty and limitation of liability

COSMED provides a one (1) year limited warranty from the date of the original sale of COSMED products. All COSMED products are guaranteed to be free from defect upon shipment. COSMED's liability for products covered by this warranty is limited exclusively to replacement, repair, or issuance of a credit for the cost of a defective product, at the sole discretion of COSMED. COSMED shall not be liable under the foregoing warranty unless (i) COSMED is promptly notified in writing by Buyer upon discovery of defect; (ii) the defective product is returned to COSMED, transportation charges prepaid by Buyer, (iii) the defective product is received by COSMED no later than four weeks after the last day of the one (1) year limited warranty period; and (iv) COSMED's examination of the defective product establishes, to COSMED's exclusive satisfaction, that such defect was not caused by misuse, neglect, improper installation, unauthorised repair or alteration, or accident. If the product is manufactured by a third-party, COSMED shall make available for the Buyer's benefit only those warranties which COSMED has received from the third-party manufacturer(s). COSMED hereby specifically disclaims any and all warranties and/or liabilities arising from defect(s) and/or damage(s) to and/or caused by products manufactured by third-party manufacturers. Buyer must obtain written authorisation from COSMED prior to the repair or alteration of COSMED product(s). Failure of Buyer to obtain such written authorisation shall void this warranty.

COSMED hereby specifically disclaims any and all other warranties of any kind, whether express or implied, in fact or by law, including, but without limitation, any and all warranties of merchantability and/or fitness for a particular purpose.

COSMED shall not be liable for special, indirect and/or consequential damages, nor for damages of any kind arising from the use of any COSMED's products, whether said products are used alone or in combination with other products or substances.

Determination of the suitability of any of COSMED's product(s) furnished hereunder for the use contemplated by Buyer is the sole risk and responsibility of Buyer, and COSMED has no responsibility in connection therewith. Buyer assumes all risks

and liabilities for loss, damage or injury to persons or property of Buyer or others arising out of the use or possession of COSMED's products.

The limited warranty as herein above set forth shall not be enlarged, diminished, modified or affected by, and no obligation or liability shall arise or grow out of, the renderings of technical advice or service by COSMED, its agents or employees in connection with Buyer's order or use of the product(s) furnished hereunder.

Return goods policy for warranty or non warranty repair

Goods shipped to COSMED for repair are subject to the following conditions:

1. Goods may only be returned after your receipt of a **Service Return Number (SRN)** from COSMED S.r.l.
2. Place your SRN report and Packing List outside the package.
3. Goods returned must be shipped with freight and insurance charges prepaid. **Collect shipments will not be accepted.**
4. The following list of goods are not eligible for return unless proven defective.
 - Special order items
 - Expendable products
 - Goods held over 30 days from COSMED's invoice date.
 - Used goods not in original shipping containers.
 - Goods which have been altered or abused in any way.
5. The following parts are not covered by warranty:
 - consumables
 - fragile glass or plastic parts
 - rechargeable batteries
 - damages due to use of the device not conforming to the indication reported in this manual

Repair Service Policy

Goods returned to seller for Non-Warranty repair will be subject to conditions 1, 2, 3, 4.

The returned goods need to re-enter COSMED together with the customs documents (Pro-forma Invoice and Customs Paper) as requested by the Italian law.

- The shipment has to be qualified as a Temporary Export.

-
- All the goods returned to COSMED without the customs papers will not be accepted.

For European Community members:

Pro-Forma invoice complete with:

- Number
- Description of the goods
- Quantity
- Serial Number
- Value in €
- Number of parcel
- Gross weight
- Net weight
- Reason for resent (i.e. Resent for repair)

In case you should send the system for repair please contact the nearest service centre or contact COSMED at the following address:

COSMED S.r.l.

Via dei Piani di Monte Savello 37

P.O. Box 3

00040 Pavona di Albano - Rome, Italy

tel. +39 (06) 9315492

fax +39 (06) 9314580

E-mail: customersupport@cosmed.it

For USA customers only please contact:

COSMED USA Inc

1808 North Halsted Street

Chicago, IL 60614 USA

Phone: +1 (312) 642-7222

Fax: +1 (312) 642 7212

email: usa.sales@cosmed.it

To ensure that you receive efficient technical assistance, please specify as precisely as possible the nature of the problem as it is specified on the assistance information form.

We advise you to save the original packaging. You may need it in case to ship the unit to a technical assistance centre.

Privacy Information

Dear Customer,

we inform you that your personal data are gathered and will be used by Cosmed Srl in conformity with the requirements of the Italian privacy law (Decreto Legislativo 196/2003). We believe it is important for you to know how we treat your personal data.

Personal data treatment and purposes

We request and process your personal data:

- a. to place an order, register a product, request a service, answer a survey, enter a contest, correspond with us (all of the above, in the following: “service”) and, if necessary, to supply the Competent Authorities with the required information;
- b. in order to define your commercial profile;
- c. in order to use your commercial profile for own marketing and advertising purposes;
- d. for accounting purposes, including e-mailing of commercial invoices;
- e. for providing your information to selected business partners (also abroad), in order to supply the service;

How your personal data are treated

Your personal data will be stored in electronic format, and protected at the best from destruction, loss (even accidental), not authorized accesses, not allowed treatment or use not in conformity with the purposes above listed.

The consent is optional, but...

If you deny the consent, we regret we cannot supply the service.

Holder of the treatment

The holder of the treatment is Cosmed Srl, Via dei Piani di Monte Savello 37, Pavona di Albano Laziale (RM). The responsible of the personal data treatment is indicated in the documentation stored by Cosmed Srl itself.

Customer rights

In accordance with art.7 of the Law, you can:

- a. obtain confirmation of the existence of your personal data and their communication in intelligible form;
- b. obtain:
 - updating, correction or integration of your data;
 - deletion or transformation in anonymous form of your personal data;
- c. deny your consent to the treatment of your personal data;

These rights can be exercised directly requesting in writing to the holder of the treatment.

Converting factors configuration



You can edit the parameters shown in Control Panel by selecting **Control Panel** from the **Calibration** menu in the calibration program, then pressing the button by side.

You might configure the following options:

Name: identify the parameter

Unit of meas.: unit of measurement

Base line and Gain: factors used to convert the acquired raw data (mV) into the final format according to $Y=(mV-BL)*Gain$. The value entered for gain must be multiplied by 1000 (for Gain=1, enter 1000).

Precision: the number of decimals shown as **0**

ATS 94 recommendations

Reference: “Standardization of Spirometry: 1994 Update”
“American J. Respiratory Critical Care Medicine”, Vol. 152,
1107-1136; 1995.

ATS recommendations

Volume range: 8l (BTPS)
Flow range: ± 14 l/sec
Volume accuracy: $\pm 3\%$ or < 50 ml
Flow accuracy: $\pm 5\%$ or < 200 ml/sec
Flowmeter resistance: < 1.5 cmH₂O da 0 a 14 l/sec

Reproducibility: the 2 largest of 3 acceptable FEV1 and FVC values should be within 5% or 150 ml.

The end of test: no change in volume for 1 second with at least 6 seconds of collected volume.

Accumulation time: the maximum time allowed for volume accumulation during the VC manoeuvre should be at least 30 seconds and at least 15 seconds during the FVC.

The spirometer should be store at least 8 FVC manoeuvres.

FEV1 should be calculated by using the “back extrapolation” method to detect the start of the test, extrapolated volume must not be higher then 5% FVC or 150ml.

The graphic resolution of the printed report must be as in the following:

Volume: 10 mm/l
Flow: 5 mm/l/sec
Time: 20 mm/sec
F/V ratio: 2:1

The total number of error (FVC e FEV1 $> \pm 3.5\%$, FEF25-75% $> 5.5\%$) during the measurement of the 24 standard waveforms must be lower than 4.

Predicted values

ERS93

Reference Adult:

Standardized Lung Function Testing: Official Statement of the European Respiratory Society, The European Respiratory Journal Volume 6, Supplement 16, March 1993.

Reference Paediatric:

Compilation of reference values for lung function measurements in children: Ph. H. Quanjer, J. Stocks, G.Polgar, M. Wise, J. Karlberg, G. Borsboom; ERJ 1989, 2, Supp.4,184s-261s.

KNUDSON 83

Reference Adult/ Paediatric:

Changes in the Normal Maximal Expiratory Flow-Volume Curve with Growth and Anging: J. Knudson, D. Lebowitz, J. Holdberg, B. Burrows; ARRD 1983; 127:725-734

Note: SD@FEV1/FVC e FEV1/VC da ERS93

ITS (White race)

Reference Adult/ Paediatric:

Intermountain Thoracic Society: Clinical Pulmonary Function Testing, second edition (1984) pp 101, 144

Note: SD@FEV1/FVC e FEV1/VC da ERS93

ITS (Black race)

Reference Adult/ Paediatric:

Intermountain Thoracic Society: Clinical Pulmonary Function Testing, second edition (1984) pp 101, 144

Note: SD@FEV1/FVC e FEV1/VC da ERS93

LAM

Reference Adult/ Paediatric:

A survey of ventilatory capacity in Chinese subjects in Hong Kong: Lam Kwok-Kwong, Pang Shing et Al. Annals of Human Biology, 1982, vol. 9, No. 5, 459-472.

Note: $SD@FEV1/FVC$ e $FEV1/VC$ da ERS93

Multicéntrico de Barcelona

Reference Adult/ Paediatric:

Spirometric reference values from a Mediterranean population: J. Roca, J. Sanchis, A. Agusti-Vidal, F. Segarra, D. Navajas. R. Rodriguez-Roisin, P. Casan, S. Sans. Bull. Eur. Physiopathol. Respir. 1986, 22, 217-224.

NHANES III

Reference Adult/ Paediatric:

Spirometric reference values from a sample of the general US population: John L. Hankinson, John. R. Odencrantz and Kathleen B. Fedan. Am J Respir Critr Care Med 1999, 159, 1798-187.

Automatic diagnosis (algorithm)

Reference: “Lung Function Testing: selection of reference values and interpretative strategies”, A.R.R.D., 144/ 1991:1202-1218.

$LLN = Pred - 0.674 * SD$ (ATS, 50° percentile)

$LLN = Pred - 1.647 * SD$ (ERS, 95° percentile)

$LLN = Pred * 0.8$ (80%Pred)

| Message interpretation | Criterion |
|--|------------------------------------|
| Normal Spirometry | FVC and $FEV1/FVC > LLN$ |
| Obstructive abnormality (it may be physiological) | $\% Pred FEV1 \geq 100$ |
| Obstructive abn.: mild | $\% Pred FEV1 < 100$ and ≥ 70 |
| Obstructive abn.: moderate | $\% Pred FEV1 < 70$ and ≥ 60 |
| Obstructive abn.: mod. severe | $\% Pred FEV1 < 60$ and ≥ 50 |

| | |
|-------------------------------|----------------------------|
| Obstructive abn.: severe | % Pred FEV1 < 50 and >= 34 |
| Obstructive abn.: very severe | % Pred FEV1 < 34 |
| Restrictive abn.: mild | FVC<LLN and %Pred FVC>=70 |
| Restrictive abn.: moderate | % Pred FVC < 70 and >= 60 |
| Restrictive abn.: mod. severe | % Pred FVC < 60 and >= 50 |
| Restrictive abn.: severe | % Pred FVC < 50 and >= 34 |
| Restrictive abn.: very severe | % Pred FVC < 34 |

Quality Control Messages

Reference: Spirometry in the Lung Health Study: Methods and Quality Control, ARRD 1991; 143:1215-1223.

| Message | Criterion |
|----------------------|--------------------------------------|
| Start faster | VEXT >5% of the FVC and >150ml |
| Blast out harder | PEFT >120 msec |
| Avoid coughing | 50% drop in the flow in first second |
| Blow out longer | FET100% <6 sec. |
| Blow out more air | flow >0.2l/s within 20 ml of FVC |
| Blow out harder | dPEF<10% |
| Take a deeper breath | dFVC<200ml and 5% best FVC |
| Blow out faster | dFEV1<200ml and 5% FEV1 |
| That was a good test | No errors |
| FVC reproducible | diff. 2 max FVC within 0.2 l |
| FEV1 reproducible | diff. 2 max FEV1 within 0.2 l |
| PEF reproducible | diff. 2 max PEF within 10 % |
| MVV time too short | MVV time less than 12 sec |

References

ATS '94: "Standardization of Spirometry: 1994 Update", American J. Respiratory Critical Care Medicine, Vol. 152, 1107-1136; 1995

ERS '93: "Standardized Lung Function Testing: Official Statement of the European Respiratory Society", The European Respiratory Journal Volume 6, Supplement 16, March "

"Standardization of Spirometry: 1987 Update", American Review of Respiratory Disease, Vol. 136, 1285-1289; 1987

Lung function", J.E. Cotes, Blackwell scientific publications

"Guidelines for Clinical Exercises Testing Laboratories", I.L. Pina, G.J. Balady, P. Hanson, A.J. Labovitz, D.W. Madonna, J. Myers. American Heart Association. 1995; 91, 912.

"Office spirometry", R.E. Hyatt - P.L. Enright, Lea & Febiger