



---

## Ultra-Low Freezer Operation Manual

### i.Series®



#### **i.Series**

---

iUF118, iUF126 (Version A)

---

iUF116, iUF124 (Version A)

---

HELMER SCIENTIFIC  
14400 Bergen Boulevard  
Noblesville, IN 46060 USA



PH +1.317.773.9073  
FAX +1.317.773.9082  
USA and Canada 800.743.5637

**360172-A/B**

**ISO 13485:2003 CERTIFIED**

# Document History

Revision	Date	CO	Supersession	Revision Description
A	10 APR 2014*	9339	n/a	Initial release.
B	12 DEC 2014	10128	B supercedes A	Captured product changes prior to product launch

\* Date submitted for Change Order review. Actual release date may vary.

# Contents

---

<b>Section I: General Information</b>	<b>4</b>
<b>1 About this Manual</b>	<b>4</b>
1.1 Intended Audience	4
1.2 Model References	4
1.3 Copyright and Trademark	4
<b>2 Safety</b>	<b>4</b>
2.1 Safety Definitions	4
2.2 Product Labels	5
2.3 Avoiding Injury	5
<b>3 General Recommendations</b>	<b>6</b>
3.1 Intended Use	6
3.2 General Use	6
3.3 Initial Loading	6
<b>4 Specifications</b>	<b>6</b>
<b>5 Compliance</b>	<b>8</b>
5.1 Regulatory Compliance	8
5.2 Electromagnetic Compliance	8
<b>Section II: Installation</b>	<b>9</b>
<b>6 Location Requirements</b>	<b>9</b>
6.1 Placement	9
<b>7 Install Components</b>	<b>9</b>
7.1 Install Rear Stand-Offs	9
7.2 Install AC Power Cord Retainer	10
7.3 Install Shelves	11
7.4 Chart Recorder (Optional)	12
7.4.1 Install and Change Chart Paper	12
<b>Section III: Operation</b>	<b>13</b>
<b>8 Initial Start Up</b>	<b>13</b>
<b>9 Normal Operation</b>	<b>14</b>
9.1 Change Temperature Setpoint	14
9.2 Active Alarms	14
9.3 Mute Active Alarms	15
9.4 Set Alarm Parameters	15
9.5 Access Control	16
9.5.1 Setup	16
9.5.2 Open Freezer with Access Control	17
9.6 i.C <sup>3</sup> ® Icon Reference Guide	17

**Section IV: Maintenance . . . . . 18**

**10 Maintenance Schedule . . . . . 18**

**Section V: Components . . . . . 19**

**11 Exterior . . . . . 19**

    11.1 Front . . . . . 19

    11.2 Rear Panel . . . . . 20

**12 Interior . . . . . 21**

    12.1 Chamber . . . . . 21

## Section I: General Information

### 1 About this Manual

#### 1.1 Intended Audience

This manual is intended for use by end users of the freezer and authorized service technicians.

#### 1.2 Model References

Generic references are used throughout this manual to group models that contain similar features. For example, “iUF models” refers to all models of that size (iUF116, iUF118, iUF124, iUF126). This manual covers all ultra-low freezers, which may be identified singly, by their size, or by their respective “Series.”

#### 1.3 Copyright and Trademark

Helmer®, i.Series®, i.C<sup>3</sup>®, and Rel.i™ are registered trademarks or trademarks of Helmer, Inc. in the United States of America. Copyright © 2014 Helmer, Inc. All other trademarks and registered trademarks are the property of their respective owners.

Helmer, Inc., doing business as (DBA) Helmer Scientific and Helmer.

### 2 Safety

The operator or technician performing maintenance or service on Helmer Scientific products must (a) inspect the product for abnormal wear and damage, (b) choose a repair procedure which will not endanger his/her safety, the safety of others, the product, or the safe operation of the product, and (c) fully inspect and test the product to ensure the maintenance or service has been performed properly.

#### 2.1 Safety Definitions

The following general safety alerts appear with all safety statements within this manual. Read and abide by the safety statement that accompanies the safety alert symbol.



#### **WARNING**

The safety statement that follows this safety alert symbol indicates a hazardous situation which, if not avoided, could result in serious injury.



#### **CAUTION**

The safety statement that follows this safety alert symbol indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.



#### **NOTICE**

The safety statement that follows this safety alert symbol indicates a situation which, if not avoided, could result in damage to the product or stored inventory.

## 2.2 Product Labels

The following general safety and information alerts appear on the product to identify potential hazards to the operator or service technician.



Caution: Safety hazard to operator or service technician



Caution: Unlock all casters



Caution: Electrocution/shock hazard



Earth / ground terminal



Caution: Electrostatic discharge (ESD) hazard



Protective earth / ground terminal

## 2.3 Avoiding Injury

- ▶ Review safety instructions before installing, using, or maintaining the equipment.
- ▶ Before moving unit, remove contents from the chamber.
- ▶ Before moving unit, ensure door is closed and latched, and casters are unlocked and free of debris.
- ▶ Before moving unit, disconnect the AC power cord and secure the cord.
- ▶ When moving unit, use assistance from a second person.
- ▶ Never physically restrict any moving component.
- ▶ Avoid removing electrical service panels and access panels unless so instructed.
- ▶ Use appropriate gloves when handling cold internal components and stored inventory.
- ▶ Keep hands away from pinch points when closing the door.
- ▶ Avoid sharp edges when working inside the electrical compartment and refrigeration compartment.
- ▶ Ensure biological materials are stored at recommended temperatures determined by standards, literature, or good laboratory practices.
- ▶ Proceed with caution when adding and removing samples from the freezer.
- ▶ Total freezer weight (including contents) is not to exceed 1400 lbs (635 kg).
- ▶ Individual shelf load is not to exceed 160 lbs (73 kg).
- ▶ Use supplied power cord only.
- ▶ Using the equipment in a manner not specified by Helmer may impair the protection provided by the equipment.
- ▶ Decontaminate parts prior to sending for service or repair. Contact Helmer or your distributor for decontamination instructions and a Return Authorization Number.
- ▶ Ensure biological materials are stored safely, in accordance with all applicable organizational, regulatory, and legal requirements.
- ▶ The freezer is not considered to be a storage cabinet for flammable or hazardous materials.

## 3 General Recommendations

### 3.1 Intended Use

Helmer ultra-low freezers are intended to provide a controlled temperature environment at ultra-low temperatures required for the storage of biological materials, pharmaceuticals, and reagents used in a research or clinical laboratory.

The devices referenced in this manual are intended to be operated by personnel who have procedures in place for meeting FDA, AABB, or any other applicable regulations for the processing and storage of biological materials, pharmaceuticals, and reagents.

### 3.2 General Use

Allow freezer to come to room temperature before switching power on.

**NOTE** During initial startup, high temperature alarm may activate while freezer reaches operating temperature. You may want to temporarily change the volume on the audible alarm during this period.

**NOTE** This unit is not a “rapid-freezing” device. Freezing large quantities of liquid, or high-water content items, will temporarily increase the chamber temperature and will cause the compressors to operate for a prolonged period of time.

### 3.3 Initial Loading

Allow chamber temperature to stabilize at the setpoint before storing pre-frozen product.

## 4 Specifications

	iUF116	iUF118	iUF124	iUF126
<b>Internal Dimensions (W x H x D)</b>				
<b>Standard/English (in)</b>	23.1 x 49.5 x 23.3	23.1 x 54.1 x 23.3	34.4 x 49.5 x 23.3	34.4 x 54.1 x 23.3
<b>Metric (mm)</b>	587 x 1257 x 592	587 x 1374 x 592	874 x 1257 x 592	874 x 1374 x 592
<b>External Dimensions (W x H x D) <sup>(1)</sup></b>				
<b>Standard/English (in)</b>	33.8 x 78.2 x 34.8	28.9 x 78.2 x 34.8	45.1 x 78.2 x 34.8	40.2 x 78.2 x 34.8
<b>Metric (mm)</b>	859 x 1986 x 884	734 x 1986 x 884	1146 x 1986 x 884	1021 x 1986 x 884
<b>Overall Dimensions (W x H x D) <sup>(2)</sup></b>				
<b>Standard/English (in)</b>	37.1 x 78.2 x 37.5	32.5 x 78.2 x 37.5	48.4 x 78.2 x 37.5	43.8 x 78.2 x 37.5
<b>Metric (mm)</b>	942 x 1986 x 953	826 x 1986 x 953	1229 x 1986 x 953	1113 x 1986 x 953
<b>Physical</b>				
<b>Weight</b>	607 lbs (275 kg)	622 lbs (282 kg)	704 lbs (319 kg)	725 lbs (328 kg)
<b>Interior Volume</b>	16 ft <sup>3</sup> / 453 L	18 ft <sup>3</sup> / 510 L	24 ft <sup>3</sup> / 680 L	26 ft <sup>3</sup> / 736 L
<b>Refrigeration System</b>				
<b>High Stage Refrigerant</b>	R-404A, CFC/HCFC-free			
<b>Low Stage Refrigerant</b>	R-508B, CFC/HCFC-free and R601 natural refrigerant			
<b>High Stage Compressor</b>	1.5 HP, air-cooled			
<b>Low Stage Compressor</b>	1.5 HP, air-cooled			
<b>High Stage Initial Charge</b>	48 oz. (1361 g)			
<b>Low Stage Initial Charge (R508B)</b>	13.5 oz. (383 g)		16 oz. (454 g)	

	iUF116	iUF118	iUF124	iUF126
Low Stage Initial Charge (R601)	0.56 oz (26 g)		0.56 oz (26 g)	
Operational				
Default Setpoint	-80 °C (-112 °F)			
Temperature Control Range	-50 °C to -86 °C (-58 °F to -123 °F)			

(1) Includes casters.

(2) Includes casters, handle, i.C<sup>3</sup> bezel, and door hinges.

	iUF116	iUF118	iUF124	iUF126
Cabinet				
Insulation	Ecomate® insulating foam (non-ODP, non-GWP, and VOC-exempt blowing agent foam)			
Vacuum-Insulated Panels <sup>(3)</sup>	-		-	
Wall Thickness	5.0" (127 mm)	2.7" (69 mm)	5.0" (127 mm)	2.7" (69 mm)
Door Thickness	2.7" (69 mm)			
Internal Compartments	4	5	4	5
External Material	Galvannealed steel with bacteria-resistant powder-coated finish			
Internal Material	Galvannealed steel with bacteria-resistant powder-coated finish			
Shelves	Stainless steel			
Maximum Shelf Load	160 lbs (73 kg)			
External Port	2, standard (top-left corner, rear of cabinet; bottom-left corner, rear of cabinet)			
Temperature Chart Recorder	Optional, 4" (102 mm), 7-day inkless, pressure-sensitive chart paper, backup battery			
Vacuum Break Port	Standard (heated)			
Electrical				
Input Voltage	208/230 V, 60 Hz			
Voltage Tolerance	±10%			
Circuit Breakers	12 A (quantity 2)			
Current Draw <sup>(4)</sup>	11.0 A			
Energy Consumption	18.5 kWh/day	18.5 kWh/day	19 kWh/day	19 kWh/day
Power Source	15 A dedicated circuit			
Boost Cut-In Voltage <sup>(5)</sup>	195 V			
Monitoring System Battery	12 V, 7 Ah rechargeable sealed lead acid battery			
Chart Recorder Battery	Powered from the monitoring system battery			
Temperature Monitor Accuracy				
Chart Recorder	±0.5 °C (0.9 °F)			
i.C <sup>3</sup> Monitor	±0.5 °C (0.9 °F) at setpoint			
Control and Monitoring				
Interface	i.C <sup>3</sup> ® combined monitoring and control interface, 7" color LCD touchscreen			
Alarms	<b>User-configurable:</b> High chamber temperature, low chamber temperature, power failure, door open (time), high ambient temperature, low ambient temperature  <b>Non-configurable:</b> Compressor temperature, condenser temperature, clean filter, CO <sub>2</sub> / LN <sub>2</sub> active, low battery, no battery, refrigeration system, communication failure, sensor failure			
Remote Alarm Interface	Dry contacts (configurable as normally-open or normally-closed)			
Remote Alarm Capacity	0.5 A at 30 V (RMS); 1.0 A at 60 V (DC)			



	iUF116	iUF118	iUF124	iUF126
<b>Environmental</b>				
<b>Operating Standards</b>	<ul style="list-style-type: none"> <li>▶ Indoor use only</li> <li>▶ Altitude (maximum): 2000 m</li> <li>▶ Ambient temperature range: 15 °C to 32 °C</li> <li>▶ Relative humidity (maximum for ambient temperature): up to 80% for temperatures up to 25 °C, decreasing linearly to 53% at 32 °C</li> </ul>			

- (3) Vacuum-insulated panels are included in cabinet walls on indicated models. All models feature vacuum-insulated panels in the exterior door.
- (4) Current draw is measured in full-load Amperes.
- (5) Power conditioning (voltage boost) is an optional feature. If facility voltage is consistently at or below 195V, voltage boost must be installed in order to protect the compressors.



**CAUTION** Maximum shelf load is not to exceed 160 lbs (73 kg).



**CAUTION**

- ▶ The interface on the remote alarm monitoring system is intended for connection to the end user's central alarm system(s) that uses normally-open or normally-closed dry contacts.
- ▶ If an external power supply exceeding 30 V (RMS) or 60 V (DC) is connected to the remote alarm monitoring system's circuit, the remote alarm will not function properly; may be damaged; or may result in injury to the user.

**NOTE** In the event of a power failure, the power failure alarm condition is transmitted through the remote alarm contacts.

## 5 Compliance

### 5.1 Regulatory Compliance

This product is certified to applicable UL and CSA standards by a NRTL.

Sound level is less than 70 dB(A).

### 5.2 Electromagnetic Compliance

This device is suitable for use in a specific electromagnetic environment. The end user of this device is responsible for ensuring the device is used in compliance with the following European Union directives and standards regarding EMC (electromagnetic compliance):

EMC Directive (2004/108/EC) Standards:

- ▶ EN 55011:2009
- ▶ EN 61000-3-2:2006
- ▶ EN 61000-3-3:2008
- ▶ EN 61000-6-1:2007

## Section II: Installation

### 6 Location Requirements

- ▶ Has a dedicated 15 A grounded circuit with dedicated single point receptical meeting the electrical requirements listed on the product specification label.
- ▶ Is clear of direct sunlight, high temperature sources, and heating and air conditioning vents.
- ▶ Minimum 8" (203 mm) above, and minimum 4" (102 mm) behind.
- ▶ Meets limits specified for ambient temperature and relative humidity.

#### 6.1 Placement



##### **WARNING**

- ▶ To prevent tipping, ensure door is closed and latched, and casters are unlocked and free of debris before moving freezer.
- ▶ The freezer is extremely heavy. Helmer recommends that two people work together to move the freezer.

- 1 Ensure all casters are unlocked and door is closed and latched.
- 2 Roll freezer into place and lock casters.
- 3 Adjust leveling feet as necessary to ensure freezer is level.

### 7 Install Components

#### 7.1 Install Rear Stand-Offs



##### **WARNING**

The rear stand-offs include a hole to accept a threaded fastener for anchoring the freezer to a wall. The rear stand-offs do not provide a secure means to anchor the freezer to the wall that can be considered resistant to seismic events.

##### **NOTE**

- ▶ Installation of the rear stand-offs is optional.
- ▶ Anchoring the freezer to the wall is optional.
- ▶ Hardware to anchor the freezer to the wall is not provided with the freezer. The end user is responsible for determining the best method to anchor the freezer to the wall.

##### **Required tools:**

- ▶ 9/16" box wrench



##### *Rear stand-offs*

- 1 Align the holes in the stand-offs with the corresponding threaded holes on the back of the freezer.
- 2 Insert the 3/8" hex head cap screws through the holes in the stand-offs.
- 3 Hand-thread the cap screws into the threaded holes.
- 4 Using a 9/16" open-ended wrench, tighten the cap screws.

## 7.2

## Install AC Power Cord Retainer



**WARNING** Use supplied power cord only.



**NOTICE** Do not position the freezer where it will prevent access to the power cord disconnect, at the power outlet on the wall.

**Required tools:**

- #2 Phillips screwdriver

**Install the retainer:**

- 1 Insert the power cord into the receptacle on rear of the cabinet.
- 2 Install the power cord retainer.
  - a Slide the retainer upward, engaging the groove in the power plug with the slot in the retainer.
  - b Align the holes in the retainer with the corresponding holes on the cabinet.
  - c Insert the screws with lock washers through the retainer and into the holes in the cabinet.
  - d Using a #2 Phillips screwdriver, tighten the screws.



*Power cord retainer.*

## 7.3

**Install Shelves****Install the shelf clips:**

- 1 Open the chamber door and all inner doors.
- 2 Install shelf clips on the shelf standards at the marked locations.

*Shelf clip.***NOTE**

Shelf clips must be installed so the horizontal section is oriented upward.


**Install the shelves:**

- 1 Starting with the bottom shelf, insert the shelf into the chamber at an angle.
- 2 Rotate the shelf so it sets flat on the shelf clips.

*Installed shelves.*

- 3 Working from the bottom to the top, install the remaining shelves as described in steps 1 and 2.
- 4 Close the inner doors and the chamber door.

## 7.4 Chart Recorder (Optional)




*Chart recorder with paper installed.*

Label	Description	Function
A	Left and Right Arrow buttons	Adjust settings and stylus position
B	LED	Indicates status of chart recorder in operating mode, or selected temperature range in paper change mode
C	Chart change button	Adjust position of stylus when changing chart paper, or run a test pattern
D	Stylus	Mark temperature line on paper
E	Reset button (located behind chart paper)	Restart chart recorder

### 7.4.1 Install and Change Chart Paper

- 1 Press and hold **C** button. When stylus begins to move left, release button. The LED flashes to indicate current temperature range.
- 2 When stylus stops moving, remove chart knob then move knob up and away from chart paper.
- 3 Place new chart paper on chart recorder.
- 4 Gently lift stylus and rotate paper so current time line corresponds to time line groove.



- 5 Hold chart paper and reinstall chart knob.

**NOTE** For accurate temperature reading, ensure that current time is aligned with time line groove when chart knob is tightened.

- 6 Press and hold **C** button. When the stylus begins to move right, release the button.
- 7 Confirm the stylus is marking the temperature correctly.

## Section III: Operation

### 8 Initial Start Up


- 1 Plug the power cord into a 15 A grounded circuit with a dedicated single point receptical.
- 2 Switch battery ON/OFF switch **ON**.
- 3 Switch AC ON/OFF switch **ON**.

---

**NOTE** The i.C<sup>3</sup> monitoring and control system will take approximately three minutes to boot up.

---

- 4 Touch the **Home** button.




---

**NOTE** Active alarms are displayed on the Home screen. If an alarm condition other than High Temperature occurs, refer to the service manual for troubleshooting.

---




- 5 If an alarm activates, temporarily mute the alarm by touching the **Mute** button.

## 9 Normal Operation

The i.C<sup>3</sup> Home screen displays temperature and alarm information, and provides icons for reaching other functions of the i.C<sup>3</sup>.



Home screen.



Home screen temperature graph screensaver (touch to return to Home screen).

### 9.1 Change Temperature Setpoint



- ▶ Enter the Settings password.
- ▶ Touch + or – on the spin box to change the value.

**NOTE**

- ▶ Default Settings password is 1234
- ▶ Default setpoint is -80.0 °C

### 9.2 Active Alarms



Home screen with no alarms.



Home screen with active alarm.

Alarm	Description
High Temperature	Chamber temperature reading is above high temperature alarm setpoint
Low Temperature	Chamber temperature reading is below low temperature alarm setpoint
Door Open	Door is open beyond user-specified duration
High Ambient	Ambient temperature reading is above high ambient alarm setpoint
Low Ambient	Ambient temperature reading is below low ambient alarm setpoint
Clean Filter	Condenser filter is obstructed
CO <sub>2</sub> / LN <sub>2</sub> Active	CO <sub>2</sub> / LN <sub>2</sub> backup refrigeration system is active

Low Battery	Backup battery voltage is low
No Battery	Backup battery voltage is too low or battery is disconnected
Power Failure	Power to unit has been disrupted
Sensor Failure	Sensor not functioning properly
Refrigeration System	<ul style="list-style-type: none"> <li>▶ Refrigerant pressure is too high</li> <li>▶ High stage compressor temperature is above the upper limit</li> <li>▶ Low stage compressor temperature is above the upper limit</li> <li>▶ High stage compressor has failed</li> <li>▶ Low stage compressor has failed</li> </ul>
Communication Failure Messages 1, 2, 3	<ol style="list-style-type: none"> <li>1 Communication lost between i.C<sup>3</sup> display board and control board</li> <li>2 Communication lost between i.C<sup>3</sup> display board and internal system memory</li> <li>3 Corrupt database</li> </ol>
Emergency Mode	Chamber temperature sensor has failed (or is failing intermittently) and refrigeration system is operating at 100% duty cycle

### 9.3 Mute Active Alarms


Audible alarms may be muted by touching the **Mute** button to set delay.



Unmuted    Muted

### 9.4 Set Alarm Parameters

The following alarm settings may be changed by the operator. The setpoint for temperature alarms may be changed (where applicable), as well as the time delay between when the alarm condition commences and when the visual and audible alarms are initiated.



Alarm	Description	Default Setpoint	Default Time Delay
High Temperature	Chamber temperature reading is above high temperature alarm setpoint	-70.0 °C	0 minutes
Low Temperature	Chamber temperature reading is below low temperature alarm setpoint	-90.0 °C	0 minutes
Power Failure	Power to unit has been disrupted	n/a	1 minute
Sensor Failure	i.C <sup>3</sup> control and monitoring system has lost communication with a temperature sensor	n/a	0 minutes
Door Open (Time)	Door is open beyond user-specified duration	n/a	1 minute



Alarm	Description	Default Setpoint	Default Time Delay
High Ambient	Ambient temperature reading is above high ambient alarm setpoint	30.0 °C	15 minutes
Low Ambient	Ambient temperature reading is below low ambient alarm setpoint	15.0 °C	15 minutes


## 9.5 Access Control

Allows user-specific secure access to the freezer.

### 9.5.1 Setup

> > > Access Setup


Configure and manage user-specific accounts to allow controlled access to the freezer.



- ▶ Enter the supervisor PIN to set up Access Control.
- ▶ Initial factory supervisor PIN = 5625

#### NOTE

- ▶ The supervisor PIN can not be deleted, and should be changed to prevent unauthorized user ID setup. The supervisor PIN does not allow access to the unit. At least one user ID must be set up to gain access to the unit.
- ▶ The keys provided with the freezer may be used to lock or unlock the exterior door.
- ▶ When setting up Access Control user IDs, ensure the key lock is in the locked position to prevent unauthorized access to the freezer.



#### Add a user profile:

- 1 Touch the Add User button.
  - ▶ The alphanumeric keyboard is displayed.
- 2 Enter the user ID for the new user profile.
- 3 Touch to store the user ID.

- ▶ The alphanumeric keyboard is displayed.
- 4 Enter the four-digit PIN for the new user profile.
- 5 Touch to store the PIN.
  - ▶ The User ID and PIN for the new user profile is displayed in the table.

**NOTE** For complete information on setting up, editing, and deleting Access Control user profiles, refer to the i.C<sup>3</sup> User Guide for i.Series Ultra-Low Freezers.

## 9.5.2 Open Freezer with Access Control



- ▶ Enter a valid PIN using the keypad.

## 9.6 i.C<sup>3</sup>® Icon Reference Guide

	Home		Mute		Brightness Setting
	Event Log		Download		Scroll Arrows
	Settings		Upload		Access Control
	i.C <sup>3</sup> Applications		Temperature Graph		Access Log
	Back Arrow		Information Log		Contact Helmer
	Alarm Conditions		Icon Transfer		Battery Power

## Section IV: Maintenance

### 10 Maintenance Schedule

Maintenance tasks should be completed according to the following schedule. Refer to the service manual and the i.C<sup>3</sup> User Guide for more detail on the various tasks.

**NOTE** These are recommended minimum requirements. Regulations for your organization or physical conditions at your organization may require maintenance items to be performed more frequently (as dictated by facility standard operating procedures), or only by designated service personnel.

Task	Frequency			
	3 months	1 year	2 years	As Needed
Verify the monitor/chamber temperature sensor accuracy. Calibrate the sensor if necessary.				
Verify the ambient temperature sensor accuracy. Calibrate the sensor if necessary.				
Test the High and Low chamber and Ambient Temperature alarms.				
Test the Power Failure alarm (as required by your organization's protocols).				
Test the Door Open alarm.				
Inspect electrical components and wiring terminals in the electrical box for discoloration. Contact Helmer Technical Service if any discoloration is found.				
Inspect and clean the condenser filter.				
Replace the i.C <sup>3</sup> backup battery				
Defrost and clean the chamber, exterior door gasket, and inner doors.				



#### NOTICE

- ▶ Inspect and clean the condenser filter as directed in the maintenance schedule, or when prompted by the i.C<sup>3</sup> control and monitoring system.
- ▶ The Clean Filter alarm monitors the condition of the air filter as a safety measure. The alarm is designed to warn if the filter media becomes clogged such that freezer operation and product integrity will be affected.
- ▶ The Clean Filter alarm could indicate a failure of the condenser fan.


#### NOTE

- ▶ During a power failure, the backup battery provides power to the monitoring system, power failure alarm, and chart recorder (if equipped). If the backup battery is not functioning, the power failure alarm will not be activated.
- ▶ If the backup battery does not provide power to the monitoring system during the power failure alarm test, replace the battery.
- ▶ During a power failure, the Access Control lock will continue to secure the door. To access the freezer during a power failure, the override key must be used.

Section V: Components

11 Exterior

11.1 Front



Left: Front features. Right: Side panel detail.

Label	Description
A	i.C³ user interface
B	USB port
C	Door handle with key lock, electronic lock, and padlock hasp
D	Condenser grill and filter media
E	Side panel (refer to side panel detail)
F	Temperature Chart recorder (optional)
G	Caster
H	AC ON/OFF power switch
I	Circuit breakers
J	Monitoring system backup battery ON/OFF switch
Not shown	Leveling feet

11.2 Rear Panel




Rear features.

Label	Description
A	Condenser discharge grill
B	RJ-45 Ethernet port
C	USB port
D	RS-232 serial port
E	Remote alarm interface
F	LN <sub>2</sub> / CO <sub>2</sub> backup system interface
G	Power connector
Not shown	Product specification label

12 Interior

12.1 Chamber



Interior features.



**CAUTION** Keep hands away from pinch points when closing the door.



- NOTICE**
- ▶ When removing or replacing storage racks, do not allow the storage rack to set on the top edge of a partially-open interior door.
  - ▶ To avoid damage to the interior door hinges, do not apply upward or downward force to the interior doors.

Label	Description
A	Inner door
B	Inner door hinge
C	Inner door gasket
D	Chamber temperature sensor cover
E	Chamber temperature sensor wire cover
Not shown	Chamber / chart recorder temperature sensor (behind cover)
F	Shelf standard
G	Shelf clip
H	Inner door mullion gasket
I	Shelf
J	Inner door retaining clip
K	Inner door catch

END OF MANUAL

HELMER SCIENTIFIC  
14400 Bergen Boulevard  
Noblesville, IN 46060 USA

PH +1.317.773.9073  
FAX +1.317.773.9082  
[www.helmerinc.com](http://www.helmerinc.com)

