

TECHNICAL MANUAL

MASTER PCA




TABLE OF MODIFICATIONS

Information in this Technical manual only refers to devices belonging to the Master PCA.

- File reference :NT 0850
- Revision date :May 1999
- Applicability from serial N°16230001

| Date | Revision number | Pages | Modifications |
|-------------|------------------------|--------------|----------------------|
| May 1999 | rev AO | All | Creation |

Please note: No responsibility whatsoever can be taken by Fresenius Vial S.A for any fundamental change to product specifications (specifications, performance ratings, etc.) made by non-company technician. Small repairs may be carried out by the appropriate maintenance team, at the discretion of the end user and subject to his or her responsibility. We nevertheless recommend that service personnel first attend a training course organised by Fresenius Vial S.A.

It is possible that this document contains errors or typing mistakes. Changes may be made at regular intervals, for inclusion in subsequent editions.. Thank you for your understanding.

COPYRIGHT © 1998 by Fresenius Vial S.A. This technical manual may not be reproduced in whole or in part without the written consent of Fresenius Vial S.A.

Fresenius Vial S.A. - head office : Le Grand Chemin - 38590 Brézins (FRANCE) - With directory and supervision board, capital : 90128000 F - SIREN Grenoble B 408 720 282.

TABLE OF CONTENT

| | |
|--|-----------|
| 1. OVERVIEW..... | 3 |
| 1.1. General description | 3 |
| 1.2. Block diagram | 3 |
| 1.3. Precautions before use | 3 |
| 1.4. Internal safety features..... | 4 |
| 1.5. Technical characteristics | 4 |
| 1.5.1. Device characteristics | 4 |
| 1.5.2. Biological specifications | 7 |
| 1.5.3. Mechanical characteristics | 7 |
| 1.5.4. Dimensions / Weight | 7 |
| 1.5.5. Electrical characteristics..... | 7 |
| 1.5.6. Electronic characteristics..... | 7 |
| 1.5.7. Master PCA Operator's Guide..... | 7 |
| 1.5.8. Components used for manufacturing | 7 |
| 1.5.9. Compliance | 7 |
| 1.5.10. Registrations cards..... | 7 |
| 2. ELECTRONIC BOARDS | 8 |
| 2.1. CPU BOARD | 8 |
| 2.1.1. Functional description | 8 |
| 2.1.2. Regulation - Power supply..... | 8 |
| 2.1.3. Watch dog | 9 |
| 2.1.4. Communication modules..... | 10 |
| 2.1.5. Error message..... | 12 |
| 2.1.6. BUS I2C | 14 |
| 2.1.7. Master/ Pilot link | 15 |
| 2.1.8. Command inputs and visualisation..... | 15 |
| 2.1.9. Description of connectors..... | 19 |
| 2.1.10. Electronic layout | 21 |
| 2.1.11. Installation layout..... | 21 |
| 2.2. ALARM LEDS BOARD..... | 22 |
| 2.2.1. Description of the connector J1 | 22 |
| 2.2.2. Electrical layout - | 22 |
| 2.2.3. Installation layout..... | 22 |
| 2.3. INFUSION LEDS BOARD | 22 |
| 2.3.1. Description of the connector J1 | 22 |
| 2.3.2. Electronic layout | 22 |
| 2.3.3. Installation layout..... | 22 |
| 2.4. FLEXIBLE LINK BOARD | 23 |
| 2.4.1. Description of the connectors..... | 23 |
| 2.4.2. Electronic layout | 24 |
| 2.4.3. Implantation layout | 24 |
| 3. CONFIGURATIONS, CALIBRATIONS AND CHECK | 24 |
| 3.1. Configurations | 24 |
| 3.1.1. Configuration EPROM..... | 25 |
| 3.1.2. Flash EPROM configuration | 25 |
| 3.2. Calibrations | 25 |
| 3.3. Check | 25 |
| 3.3.1. Electrical safety tests..... | 25 |
| 3.3.2. Integrated tests..... | 25 |
| 3.3.3. Test mode | 25 |


| | |
|--|-----------|
| 4. REPLACING SUB-ASSEMBLIES..... | 29 |
| 4.1. Replacing the electronic circuit boards..... | 29 |
| 4.2. Dismounting the support | 29 |
| 4.3. Replacing the flexible circuit | 29 |
| 4.4. Replacing the handle | 29 |
| 5. MAINTENANCE | 30 |
| 5.1. Cleaning and disinfecting | 30 |
| 5.2. Storage | 30 |
| 5.3. Servicing | 30 |
| 5.4. Regular inspections - Preventive maintenance | 30 |
| 5.4.1. Before using checking | 30 |
| 5.4.2. Preventive maintenance | 30 |
| 5.4.3. Internal historical saving battery | 31 |
| 6. ANNEX 1 : ILLUSTRATED PARTS LIST | 33 |
| 6.1. Traceability | 33 |
| 6.1.1. Introduction..... | 33 |
| 6.1.2. traceability table..... | 34 |
| 6.2. General view | 36 |
| 6.3. Mechanical parts list..... | 37 |
| 6.3.1. Aluminium support..... | 37 |
| 6.3.2. Front panel | 37 |
| 6.3.3. cover..... | 38 |
| 6.3.4. Locker..... | 38 |
| 6.4. Electronic parts list..... | 39 |
| 6.4.1. Programmation pass | 39 |
| 6.4.2. Control and command | 39 |
| 6.4.3. Electronic..... | 39 |
| 6.4.4. Patient hand switch | 39 |
| 6.4.5. Connectors | 39 |
| 7. ANNEX 2 : INSTALLATION AND ELECTRONIC LAYOUTS | 40 |
| 7.1. CPU BOARD..... | 40 |
| 7.2. ALARM LEDS BOARD | 40 |
| 7.3. INFUSION LEDS BOARD | 40 |
| 7.4. FLEXIBLE LINK BOARD | 40 |
| 8. ADDENDA AND INFORMATIONS..... | 41 |
| 9. USEFUL ADDRESSES..... | 43 |

1. OVERVIEW


1.1. General description

Information in this document only refers to devices belonging to the Master PCA. No responsibility whatsoever can be taken by Vial Medical for any fundamental change to product specifications (specifications, performance ratings, etc.) made by non-company personnel.

Routine repairs may be carried out by the appropriate maintenance team, at the discretion of the end user and subject to his or her responsibility. We nevertheless recommend that Maintenance personnel first attend a training course organised by Vial Medical.

The symbol  visible on the condensed instruction guide of the device, recommends the Operator Guide should be completely read, in accordance with the EN 60 601-1 Standard.

1.2. Block diagram



1.3. Precautions before use

Please consult Operator's guide

1.4. Internal safety features

The Master PCA is equipped with a continuous functions inspection system activated as soon as the Master is switched ON. Any failure or anomaly in the procedure is immediately detected. Nevertheless, qualified staff of your organisation or our After Sales Service should always be notified in case of abnormal functioning (see useful addresses, chapter 9)

1.5. Technical characteristics

The Master PCA is controlled by a keyboard and a rotary knob. The parameters input by the user, the infusion parameters, and parameters describing the instrument status may be visualised by means of a graphic LCD screen.

1.5.1. Device characteristics

1.5.1.1. Accuracy on average delivered flow rate.

- Accuracy of the $\pm 1 \%$
- Accuracy on the internal diameter of the syringes $\pm 2 \%$

1.5.1.2. Average accuracy on the bolus

This table is valid for a Pilot Anæsthésia or a Pilot C.

| Syringe | BD PLASTIPAK 20 cc | BD PLASTIPAK 50 CC |
|---------------|--------------------|--------------------|
| Volume | 0,2 ml à 0,8 ml | 0,2 ml à 2 ml |
| average error | < - 0,2 % | < - 0,15 % |

Used tubing: SE 1400 S

These measurements are performed in compliance with PrEN60 601-2-24 comply.

1.5.1.3. Bolus volumes and flow rates

| Syringe | BD PLASTIPAK 20 cc | BD PLASTIPAK 50 CC |
|--------------------|--------------------|--------------------|
| Volume | 0,2 ml to 10 ml | 0,2 ml to 30 ml |
| Injection duration | Flash to 15 min | Flash to 15 min |
| Minimum flow rate | 0,8 ml/h | 0,8 ml/h |

Used tubing: SE 1400 S


The maximum flow rate is limited by the performances of Pilot C and Anæsthésia. See connected Pilot operator's guide to know the flow rate limitation.

1.5.1.4. Exactness and accuracy of the displayed values of massic flow rate

According to concentration and massic flow rate, the Master PCA rounds the result of volumic flow rate at the first lower value programmable on the Pilot.

Conversion diagram :

The following diagram allows to determine the error versus the programmed value



1.5.1.5. Exactness and accuracy of the displayed mass values of bolus

As for the massic flow rate, the mass of the bolus is under rounded.

- Error on bolus mass..... < - 0,5 %

1.5.1.6. Accuracy on time

- Error on the locking time. < 1 s

1.5.1.7. Syringes list

| | 50cc / 60 cc | 20 cc |
|----------------------|--------------|-------|
| B-D Perfusion | • | |
| B-D Plastipak | • | • |
| Braun Omnifix | • | • |
| Braun Perfusor | • | • |
| Didactic line France | • | |
| Didactic Perfusion | • | |
| Dispomed | • | |
| Dipomed type P | • | |
| Fresenius Injectomat | • | |

| | 50cc / 60 cc | 20 cc |
|---------------------|--------------|-------|
| Fresenius P Spritze | • | |
| Sherwood Monoject | • | • |
| Ivac | • | |
| Map Gliss LL | • | |
| Map Pic indolor | • | |
| Terumo | • | • |
| Tutoject type T | • | |
| Zeneca PFS | • | |
| | | |

1.5.1.8. Pressure limit

The selection values of the pressure limit threshold are changeable by configuration from the key board. (See operator's guide of the connected Pilot for the operating procedure and the display accuracy). The threshold values or the pressure limits applied by default are the ones programmed on the Pilot.

1.5.1.9. Pressure management

- Occlusion pre-alarm and alarms.
- Disconnection / pressure drop alarm.

1.5.1.10. Occlusion alarm delay according to the infusion flow rate

Correspondences tables between occlusion alarm delay and programmed flow rate for three selectable pressure limits in three pre-adjusted threshold mode. These tables shows obviously that it is very interesting to choose, as soon as possible, the lowest alarm threshold to get the short test alarm delay.

The variable threshold mode allows to select continuously the best adapted value.

| Syringe | Flow rate | Low limit 300 mmHg | Medium limit 500 mmHg | High limit 900 mmHg |
|--------------------|-----------|-----------------------|--------------------------|------------------------|
| Bd Plastipak 50 cc | 1 ml/h | 16'20 | 30'25 | 49'50 |
| | 5 ml/h | 3'40 | 5'19 | 11'30 |
| Bd Plastipak 20 cc | 1 ml/h | 6'20 | 12'35 | 17'45 |
| | 5 ml/h | 45'' | 1'45 | 3'05 |

1 bar = 750 mmHg = 1000 hPa

Used tubing: SE 1400 S

1.5.1.11. Bolus volume at occlusion ending.

When occlusion is detected, the anti-bolus system is activated. By measuring the pressing dynamic strength, the reduction of the bolus at occlusion ending is performed according to the specific characteristics of each selectable syringe.

This principle allows to much reduce this bolus whoever the flow rate and the alarm threshold adjustment be.

| Syringe | Flow rate | Low limit 300 mmHg | Medium limit 500 mmHg | High limit 900 mmHg |
|--------------------|-----------|-----------------------|--------------------------|------------------------|
| Bd Plastipak 50 cc | 1 ml/h | ≤ 0,1 ml | ≤ 0,2 ml | ≤ 0,2 ml |
| | 5 ml/h | ≤ 0,1 ml | ≤ 0,2 ml | ≤ 0,2 ml |
| Bd Plastipak 20 cc | 1 ml/h | ≤ 0,1 ml | ≤ 0,1 ml | ≤ 0,1 ml |
| | 5 ml/h | ≤ 0,1 ml | ≤ 0,1 ml | ≤ 0,1 ml |

1 bar = 750 mmHg = 1000 hPa

Used tubing: SE 1400 S

- Note : The bolus reduction at occlusion ending is signalled by the alarm LED.

1.5.2. Biological specifications

Infusion liquid only comes into contact with the syringe and associated disposable.

1.5.3. Mechanical characteristics

The mechanical system consists of an aluminium support hinged to the front casing. The electrical connections and installation are effected by pushing the Pilot backwards into the support (see the Master PCA Operator's Guide). The casing is then lowered and locked onto the Pilot. On the button side two holes to insert fixing screws are available.

1.5.4. Dimensions / Weight

- Height : 135 mm
- Width : 370 mm
- Depth : 240 mm
- Weight : 1.9 kg approximately

1.5.5. Electrical characteristics

The Master PCA is powered from a Pilot syringe pump. The power supply characteristics are as follows :

- Power supply : 5.4V to 7.2V DC.
- Max. consumption : 180 mA.
- Max. power consumption : 1.3 W.

1.5.6. Electronic characteristics

The Master PCA comprises the following electronic subassemblies :

- CPU circuit board.
- LCD graphic display.
- Alarm LED circuit board.
- Infusion LED circuit board.

1.5.7. Master PCA Operator's Guide

The Master PCA Operator's Guide included operating cautions is available on request from our Customer Service.

1.5.8. Components used for manufacturing

- Case.....Polycarbonate/ABS alloy
- Control buttonPolyacetal
- Fixing handle.....Polyamide
- Protection hood.....Polycarbonate
- SupportPolyurethan painted aluminium
- Keyboard.....Polyester
- LabelsPolyester

1.5.9. Compliance

See Operator's guide

1.5.10. Registrations cards

Registration information is available upon request from our After Sales Service

2. ELECTRONIC BOARDS

2.1. CPU BOARD

2.1.1. Functional description

The Master PCA CPU circuit board is built around on the 80C320 microcontroller used in open mode. This micro is equivalent to the 80C32, but its slightly modified structure gives an overall speed improvement.

Characteristics

| | |
|-------------------------------|--|
| EPROM | 256 Ko |
| FLASH EPROM | 521 Ko |
| RAM | 128 Ko saved by a battery |
| EEPROM | 2 Kbytes for the instrument's configuration |
| WATCH DOG | MAX 691 resetable and safety RAM |
| INTERNAL CLOCK | PCF 8583P addressable by IC bus |
| LED Driver | 16 LED's max |
| Keyboard Driver | 16 keys max. |
| Display | 128 x 64 pixel screen graphics |
| T6963C integrated controller. | Manual adjustment of contrast |
| Attachments | RS TTL Master PCA / Pilot |
| RS 232 PC/ Master PCA. | RS 232 option |
| Power Supply | Continuous from 5.4 V to 7.2 V, includes regulator of slight voltage drops. Battery power safety system for the RAM and the internal clock 3.6V 60 mAh. |

The CPU circuit board has the following functions :

- power supply and regulation of the module.
- communication with the module.
- link module for the Master PCA/ Pilot .
- keyboard.
- CPU, memory.
-

2.1.2. Regulation - Power supply

The regulation/power supply module consists of a linear regulator with a low drop out voltage of for which the input voltage is supply by the Pilot battery and then regulated.

2.1.2.1. ON/OFF command

The ON/OFF key of the device is connected to the CD ON input line of the Pilot. A short press on this key switch on the power supply of the Pilot. Then, the Pilot battery supplies the necessary energy to the Master PCA (mains or battery).

When pressing ON/OFF button, the microprocessor knows the status of the button by reading TOFF by the input line from U24. It prepares the switch off mode and sends the command to turn off the voltage to the Pilot using RS 232 interface from the Pilot.

2.1.2.2. Protection of the power supply

In case of breakdown of power supply, the Pilot battery provides the necessary energy to the two devices. When the Master PCA is switched off, an internal battery saves the RAM in order to keep patient's history, and the internal clock continues to be supplied.

The battery is recharged when the device is in use. When the instrument is not in use, the minimum life expectation for a fully charged battery is 2 months.

The battery voltage of the Pilot is present on the connector subD 15 points J2, connecting the Pilot to the Master PCA.

- the minimum supplied voltage which does not provoke the resetting of the device is 5.1 V.
- the maximum supply voltage is 7.2 V

| J2 | Description |
|----|-------------|
| 6 | VBAT |
| 14 | VBAT |
| 5 | GND |
| 13 | GND |

2.1.2.3. Pre-alarms and alarms

Alarms and pre-alarms are identical to those of the Pilot. The alarm status parameters of the Pilot is transmitted via the RS 232 interface to the Master PCA.

2.1.2.4. Battery voltage levels

This voltage is measured at the battery terminals lead of the Pilot.

- pre-alarm : 5,8 V min / 6,0 V max.
- alarm : 5,6 V min / 5,8 V max.

The voltage drop between the input of the instrument's regulator and the battery is due to the internal fuse of the Pilot and various connections between. This voltage drop is 0.5 V max. for a current of 350 mA.

- the battery voltage of the Pilot can be measured inTP4

2.1.3. Watch dog

The watch dog manage the RESET line, the rising voltage of VP of the circuit voltage Vp to safety (RAM, INTERNAL CLOCK), the switching of the CE line of the RAM to keep the integrity of the data in the RAM in case of drop of the power supply is broken.


The voltage V.P. is supplied by the battery BT1.

The timer is activated by the signal WDOG generated by the output socket comprised of circuit U22 (

2.1.3.1. Voltage of Power Supply

the signal RESET is activated for a standard voltage value, after regulation, of : 4,65 V \pm 0,15 V.

2.1.3.2. Resetting



- maximum reset time of WDOG : 100 ms.
- activation of signal RESET : 200 ms.
- power supply voltage = 5.00 V \pm 0.25 V.

2.1.4. Communication modules

The CPU circuit board has 3 serial interface plug.
Connection socket for the Master PCA/ Pilot .

- RS 232-1 plug.
- RS 232-2 plug (optional).

2.1.4.1. Plug socket for the Master PCA/ Pilot

Series linking TTL of full duplex communication between the Pilot and the CPU.

2.1.4.2. Communication Master PCA/ Pilot

Two-directional series link TTL

| | |
|----------------------------|--------------------------------------|
| Transmission data speed : | 4800 Baud set value |
| | 31250 Baud by software configuration |
| Transmission data format : | 1 start bit |
| | 7 data bit |
| | 1 even parity bit |
| | 1 stop bit |

The signals are found on the male connector SubD 15 points J2 (connection Master PCA/ Pilot).

| J2 | Description |
|-----------|--------------------|
| 2 | RXD Pilot |
| 10 | TXD Pilot |
| 3 | GND power supply |

2.1.4.3. RS 232-1 connector

This connector is used for several types of communication.

2.1.4.3.1 Connection PC/ Master

Bi-directional series data RS 232-1 Bus

| | |
|----------------------------|---|
| Transmission data speed : | 4800, 9600, 19200 Baud selectable in configuration menu |
| Transmission data format : | 1 start bit |
| | 7 data bit |
| | 1 even parity bit |
| | 1 stop bit |

2.1.4.3.2 Communication Master/ Printer

This mode is permanently available when the instrument is being used to print histories.
Bi-directional link RS232 :

| | |
|-----------------------|--------------|
| Transmission speed : | 9600 Baud |
| Transmission format : | 1 start bit |
| | 8 data bit |
| | 0 parity bit |
| | 1 stop bit |

The control of flux is carried out by a hardware link between the DTR pin of the printer to the RTS pin 8 of the Master PCA.

- When the RTS pin is at + 12 V, the Master PCA considers that the printer is ready to receive some text.
 - When the printer is not connected, the RTS pin is set to - 12 V ; and the PCA therefore sends no text.
- The socket signals for the 3 modes of communication are found on the female connector SubD 9 points J7.

| J7 | Description |
|-----------|--------------------|
| 2 | RXD |
| 3 | TXD |
| 7 | CTS |
| 8 | RTS |
| 5 | GND |

2.1.4.4. RS 232-2 connector (Optional)

The unused series link RS 232 is reserved for future extensions.
The socket signals are found on the female connector sub D 9 points J6.

| J6 | Description |
|-----------|--------------------|
| 2 | RXD |
| 3 | TXD |
| 7 | CTS |
| 8 | RTS |
| 5 | GRD |

2.1.5. Error message

During the process of linking between the Master PCA and the Pilot . The following types of breakdown are detected and displayed on the LCD display screen.

Three types of messages are displayed:

❑ Alarms:

| Code | Description | Variable state |
|------|-------------------------------|----------------|
| 10 | Low battery | bit 10 |
| 11 | Wrong position of the syringe | bit 4 |
| 12 | Empty syringe | bit 9 |
| 14 | Disengagement | bit 7 |
| 15 | Syringe head | bit 5 |
| 16 | Occlusion | bit 6 |
| 18 | Dose limit reached | bit 14 |
| 19 | Wrong position of the syringe | bit 13 |
| 22 | Occlusion memory | bit 2 |
| 23 | Battery alarm | |
| 25 | Flange detection | |

❑ Recoverable errors

They allow to continue to use the device after their detection.

| Code | Description |
|------|---|
| 01 | Rotation control |
| 32 | Displacement control (on one segment) |
| 52 | Displacement control(during slack adjustment) |
| 72 | Displacement control (on total length) |
| 82 | Displacement control (versus flow rate) |

❑ **Locking errors**

The locking errors, worst, allow only to switch off the device by the ON/OFF to release.

| Code | Description |
|-------------|------------------------------------|
| 40 | Pilot E2prom |
| 50 | Pilot ADC |
| 60 | syringe parameters |
| 70 | Motor frequency |
| 03 | Pilot communication |
| 13 | Absent Pilot |
| 23 | Link to pilot closed |
| 33 | Reception |
| 43 | Transmission |
| 53 | No answer from Pilot |
| 63 | Bad Pilot answer |
| 73 | Bad Pilot type |
| 83 | Bad Pilot version |
| 93 | Master activation mode |
| 14 | period verification |
| 24 | Rotating direction check |
| 34 | Motor period check |
| 15 | Ram Master (internal) |
| 25 | Ram Master (external) |
| 35 | EPROM Master |
| 45 | E2prom Master |
| 55 | Internal clock |
| 65 | LCD Ram |
| 75 | Uarts |
| 16 | Date / time verification |
| 26 | Date / time comparison |
| 46 | Bad records history |
| 18 | Infusion maximum value reached |
| 28 | incorrect language file |
| 48 | Volume control on motor impulse |
| 58 | Flow rate control on motor impulse |
| 68 | Stop check |
| 78 | Maximum duration |

2.1.6. BUS I2C

The I2C bus is a series bus synchronous in communication with the internal clock and the EEPROM memory.

With regard to overlapping memory space between the two peripherals and to the fact that they are not addressable simultaneously, a single clock CLK I2C ensures the synchronicity of data transfer, two lines of transfer make sure of the exchange between each peripheral.

- EEPROM : 2 Kbytes 24C16 U13
- INTERNAL CLOCK : PCF 8583 P U21

| Line of communication | Socket |
|---|--------|
| CLK : clock generated by the microprocessor | U22.19 |
| SDA EEPROM : line Bi-directional exchange | P 1.0 |
| SDA HORO : line Bi-directional exchange | P 1.1 |

The clock is built around the circuit PCF 8583. It provides the day of the year, the month and hour. Voltage Vp of the battery BT1 ensures it is functioning via the RESET circuit, U11.

The system clock is generated by quartz X3 with frequency 32,768 kHz.

The signal frequency of HORO is 1 Hz.

A frequency control is carried out by an internal timer of the microprocessor and the value is then compared to the run time read by the bus CLK I2C.

Any differences are detected immediately. The frequency must be $1 \text{ Hz} \pm 10^{-6}$.

2.1.7. Master/ Pilot link

The linking connector Pilot / Master, J2, carries the logical information other than that described in paragraph 2.4. .

| J2 | Name | Function |
|----|------------|---|
| 1 | MAINS LED | The Pilot tells the Master PCA that the supply is connected. The battery is in use when the signal is at logical level 1. |
| 9 | ON KEY | The ON button on the front panel is connected in parallel with that of the Pilot. The Master's ON/OFF button is connected to the input CD ON of the Pilot. A short push of the button provokes the switching-on the Pilot power supply. The Pilot battery supplies the necessary energy to the Master PCA. |
| 11 | BUZZER | The Master PCA uses the Pilot's buzzer. The buzzer is activated by setting this line to zero. |
| 4 | OPTO MOTOR | This signal is sent from the Pilot towards the instrument and then redirected towards the input pin U1-16 . This line generates an impulse with each turn of the motor from the Pilot |
| 12 | OFF KEY | This signal tells the Master PCA that the Pilot has received a command to turn off the voltage. This line changes to logical level 1 to prepare for the voltage turn-off. |
| 6 | VBAT | This line provides the energy supply from the Pilot battery. |
| 14 | VBAT | This line provides the energy supply from the Pilot battery. |
| 5 | GND | — |
| 13 | GND | — |

2.1.8. Command inputs and visualisation

2.1.8.1. Keyboard

The keyboard is organised as a matrix of 4 lines and 4 columns. The role of each button is described the table below :

| | line 0 | line 1 | line 2 | line 3 |
|----------|---------------------|------------------|-------------|---------|
| column 0 | ENTER (Rotary knob) | | | |
| column 1 | | ENTER (keyboard) | NEW PATIENT | DISPLAY |
| column 2 | | ALARM SILENCE | START | HISTORY |
| column 3 | | PRIME | | STOP |

The two ENTER buttons are operated separately but have the same effect. The diodes D2 to D5 act to protect against short circuiting of the outputs of U26 when two buttons of the same line are pressed simultaneously.

The columns are activated for a logical level 0 by writing in U26 at address \$1000. The status of each button for lines 0 to 3 is read by a logical level 0 from the time of activation through reading the buffer U33 at address \$2000.

The ON/OFF button is connected to the input CD ON of the Pilot via the flexible plate link (see Flexible plate link).

The ON/OFF button is activated by a logical level 0. The status of the ON/OFF button is known by the CPU by reading one of the inputs of U5 set at the address \$2800.


2.1.8.1.1 Written register

| U26 | Address \$5000 | Set to zero at the RESET | |
|-----|----------------|--------------------------|-------------|
| D0 | COLUMN 0 | J11-5 | Active to 0 |
| D1 | COLUMN 1 | J11-6 | |
| D2 | COLUMN 2 | J11-7 | |
| D3 | COLUMN 3 | J11-8 | |

The data are reproduced at the outputs of U26 by a upright front on U 26-11. Only one column is activated at a time.

2.1.8.1.2 Writing cycle

- Time between two successive readings : 5 ms
- Time kept at logical level 0 : 5 μ s
- Refresh period : 20 ms.



2.1.8.1.3 Read register


| U24 | Address \$2800 | |
|-----|----------------|-------|
| D4 | LINE 0 | J11-1 |
| D5 | LINE 1 | J11-2 |
| D6 | LINE 2 | J11-3 |
| D7 | LINE 3 | J11-4 |

After the addressing of each column, the status of the lines 0 to 3 are read simultaneously through the buffer U5 on the data bus.

2.1.8.2. Coder


2.1.8.2.1 Rotation

The coder is a two phases (A, B), incremental type one It includes a push button having the same function as the ENTER key. The two phases are dry contacts, the common is earthen.



2.1.8.2.2 Coder characteristics

- 30 impulses/ 360°
- contacts normally open
- Max intensity per contact : 10 mA
- min intensity per contact : 1 mA



The signals from A and B are filtered to remove any erroneous coding due to rebounds of the contacts. A decoding consisting of a PAL U15 generates information for both decrementation and incrementation DEC, INC, which are reset to zero by the signal RAZ BOUTON. The information INC, DEC is read by the data bus every 5 ms.

Read and written register

| | |
|----|----------------|
| U | Address \$1800 |
| D0 | DEC |
| D1 | INC |

2.1.8.2.3 Push button

The coder includes a key-sensitive push button having the same function as the enter button on the keyboard.

Pressing the push button provokes the closure of the contact. The contact is connected to the keyboard matrix.

2.1.8.2.4 Description of the coder cable

| J14 | Description | |
|-----|-------------------|----------|
| 1 | CHANNEL B CODER | |
| 2 | CHANNEL A CODER | |
| 3 | PUSH BUTTON CODER | column 0 |
| 4 | PUSH BUTTON CODER | line 0 |
| 5 | GND | |

2.1.8.2.5 Indicators

All the right indicators are electro-luminescent diodes. The status of each diode is defined by the level of the corresponding output in the registers U23 and U26.

2.1.8.2.6 Written records

| | | |
|-----|----------------------|----------------------|
| U23 | Address \$800 | Set to zero at RESET |
| D0 | INFUSION LED 1 | J6.2 |
| D1 | INFUSION LED 2 | J6.3 |
| D2 | INFUSION LED 3 | J6.4 |
| D3 | INFUSION LED 4/START | J6.5 |
| D4 | STOP LED | J6.6 |
| D5 | OPEN | J6.7 |
| D6 | BATTERY LED | J4.2 |
| D7 | PRE ALARM LED | J4.5 |
| U26 | Address \$1000 | Set to zero at reset |
| D4 | ALARM LED | J4.4 |
| D5 | OPEN | J4.7 |
| D6 | OPEN | J4.6 |
| D7 | PATIENT LED | J2.12 |

2.1.8.2.7 Writing cycle

- Reset time : 10 ms
- Time kept at logical level 0 : 0.35 μ s



The electro-luminescent diodes ALARM and PRE ALARM consist of two pairs diodes set in parallel.

2.1.8.2.8 "MAINS SUPPLY PRESENT" signal

This signal is physically controlled by the signal Main LED (LED SECTEUR) given by the Pilot.

2.1.8.3. Additional outputs


The CPU circuit board has some unused peripheral outputs for its internal management, and others not in use.

| U22 | Address | Description |
|-----|-------------|---|
| D0 | PG-RAM0 | This signal allows the addressing of the RAM U6 to be extended in combination with the address bit A15. The decoded memory space is a block of 32 Kbytes of addresses between \$10000 and \$1FFFF. |
| D1 | PG-RAM1 | This signal allows the addressing of the RAM to be extended in combination with the address bit A15 and PG-RAM1. The supplementary decoded memory space is a block of 32 Kbytes of addresses between \$20000 et \$3FFFF. |
| D2 | RTS OPT | This signal is the "Request to Send" of the RS 232-2 link before conversion of the voltage level. |
| D3 | WDOG | This signal is connected to the input WDI of U11. It's the activation signal for the watch-dog. |
| D4 | CD CLEF | This signal allows to detect the key state. |
| D5 | CD CAPOT | This signal allows to detect the hood state. |
| D6 | CD PATIENT. | This signal allows to detect the patient hand set state. |
| D7 | CLK I2C . | This signal is the clock for activation of the bus I2C. |

2.1.8.4. CPU circuit board configuration

The CPU board is configured to be able to function with a program in EPROM or in flash EPROM.

- The configuration is carried out by contacts G1, G7 and G6.



2.1.8.4.1 Configuration EPROM

The circuit board can be fitted out with an EPROM 27C512 to 27C040, of 100 ms access time.

| | | |
|-----------|------|------------------------|
| G1 and G7 | open | EPROM type |
| G6 | 2.3 | EPROM 27C010, 020, 040 |
| G6 | 1.2 | EPROM 27C512 |

2.1.8.4.2 Configuration flash EPROM

The board is fitted with a flash EPROM 28F004.

- G1 and G7 closed using FLASH EPROM

2.1.9. Description of connectors

2.1.9.1. J2 : connection to the Pilot and external peripherals

This connector is attached to a flexible plate which redistributes signals to the Pilot and to the external connectors.

| Pin | Description | |
|-----|-----------------|----------------------------|
| 1 | GND | power supply |
| 2 | RTS | printer |
| 3 | TXD | printer |
| 4 | RIS | printer |
| 5 | RXD | printer |
| 6 | GND | power supply |
| 7 | CTS | option |
| 8 | TXD | option |
| 9 | RTS | option |
| 10 | RXD | option |
| 11 | V BAT | power supply |
| 12 | PATIENT LED | LED patient hand switch |
| 13 | PATIENT | patient hand switch input |
| 14 | CD PATIENT CALL | patient hand switch output |
| 15 | GRD | power supply |
| 16 | LOCK | hood state input |
| 17 | CD ILS LOCK | hood state output |
| 18 | VBAT | power supply |
| 19 | VBAT | power supply |
| 20 | GND | power supply |
| 21 | GND | power supply |
| 22 | TP8 | not in use |
| 23 | CTR MOT | motor rotation |
| 24 | CD BUZZER | BUZZER command |
| 25 | GND | power supply |
| 26 | TXD | Pilot |
| 27 | RXD | Pilot |
| 28 | TON | stop command |
| 29 | LED SUPPLY | supply presence |
| 30 | GND | power supply |

2.1.9.2. J3 : connector for LCD display

| Pin | Description | |
|-----|-------------|---------------------------|
| 1 | GND | power supply |
| 2 | GND | power supply |
| 3 | + 5V | power supply |
| 4 | BACK PLANE | polarisation display |
| 5 | WR* | transmission control line |
| 6 | RD* | reception control line |

| | | |
|----|--------|-------------------------------------|
| 7 | CD LCD | line for validation of memory space |
| 8 | A00 | line address A00 |
| 9 | RESET* | initialisation display |
| 10 | D0 | data line D0 |
| 11 | D1 | data line D1 |
| 12 | D2 | data line D2 |
| 13 | D3 | data line D3 |
| 14 | D4 | data line D4 |
| 15 | D5 | data line D5 |
| 16 | D6 | data line D6 |
| 17 | D7 | data line D7 |
| 18 | GND | power supply |

2.1.9.3. J4 : alarm display connector

| Pin | Description | |
|-----|---------------|---------------|
| 1 | GND | power supply |
| 2 | BATTERY LED | alarm circuit |
| 3 | SUPPLY LED | alarm circuit |
| 4 | ALARM LED | alarm circuit |
| 5 | PRE-ALARM LED | alarm circuit |
| 6 | NU | not in use |
| 7 | NU | not in use |
| 8 | GND | power supply |

2.1.9.4. J6 : infusion display connector

| Pin | Description | |
|-----|----------------------|--------------|
| 1 | GND | power supply |
| 2 | INFUSION 1 LED | power supply |
| 3 | INFUSION 2 LED | power supply |
| 4 | INFUSION 3 LED | power supply |
| 5 | INFUSION 4/START LED | On signal |
| 6 | STOP | Stop signal |
| 7 | NU | not in use |
| 8 | GND | power supply |

2.1.9.5. J8 : Volume control connector

| Pin | Description | |
|-----|-----------------|----------------|
| 1 | Motor impulse | Volume control |
| 2 | Motor direction | not in use |

2.1.9.6. J9 : connector back light of the display

| Pin | Description | |
|-----|----------------------|--------------------|
| 1 | BACK LIGHT ANODE +5V | back light anode |
| 2 | BACK LIGHT CATHODE | back light cathode |

2.1.9.7. J10 : debug connector

| Pin | Description | |
|-----|-------------|---------------------------|
| 1 | GND | power supply |
| 2 | WR* | transmission control line |
| 3 | RD* | reception control line |
| 4 | INT UART | Interruption line 0 |
| 5 | RESET* | reset line |
| 6 | DEVAL | peripheral validation |
| 7 | A15 | line of address A15 |

2.1.9.8. J11 : keyboard connector

| Pin | Description | |
|-----|-------------|------------------------|
| 1 | LINE 0 | keyboard matrix line |
| 2 | LINE 1 | keyboard matrix line |
| 3 | LINE 2 | keyboard matrix line |
| 4 | LINE 3 | keyboard matrix column |
| 5 | COLUMN 0 | keyboard matrix column |
| 6 | COLUMN 1 | keyboard matrix column |
| 7 | COLUMN 2 | keyboard matrix column |
| 8 | COLUMN 3 | |
| 9 | T ON OFF | ON/OFF |
| 10 | GND | power supply |

2.1.9.9. J12 : power supply monitor connector

| Pin | Description | |
|-----|--------------|-----------------------|
| 1 | VBAT FILTRED | filtered power supply |
| 2 | GND | power supply |

2.1.9.10. J14 : connector coder

| Pin | Description | | |
|-----|-------------------|-----------------------|----------|
| 1 | CHANNEL B CODER | decoding PAL | |
| 2 | CHANNEL A CODER | decoding PAL | |
| 3 | GND | power supply | |
| 4 | CODER PUSH BUTTON | keyboard matrix | column 0 |
| 5 | CODER PUSH BUTTON | keyboard matrix | line 0 |
| 6 | GND | power supply | |
| 7 | +5V | power supply | |
| 8 | CD KEY | key presence detector | |
| 9 | CD KEY | not in use | |
| 10 | GND | power supply | |

2.1.10. Electronic layout

See Annex 2.

2.1.11. Installation layout

See Annex 2.

2.2. ALARM LEDS BOARD

2.2.1. Description of the connector J1

This connector joins the LED's to the elevated current outputs from the CPU board. It is connected via J4 on the CPU board.

| Pin | Description | |
|-----|---------------|---------------|
| 1 | GND | power supply |
| 2 | MAIN LED | alarm circuit |
| 3 | MAIN LED | alarm circuit |
| 4 | ALARM LED | alarm circuit |
| 5 | PRE ALARM LED | alarm circuit |
| 6 | NU | not in use |
| 7 | NU | not in use |
| 8 | GND | power supply |

2.2.2. Electrical layout -

See Annex 2

2.2.3. Installation layout

See Annex 2.

2.3. INFUSION LEDS BOARD

2.3.1. Description of the connector J1

This connector joins the LED's to the outputs at an elevated current outputs from the CPU board. It is connected via J6 on the CPU board.

| Pin | Description | |
|-----|-------------------|------------------|
| 1 | GND | power supply |
| 2 | INFUSION LED 1/ON | infusion circuit |
| 3 | INFUSION LED 2 | infusion circuit |
| 4 | INFUSION LED 3 | infusion circuit |
| 5 | INFUSION LED 4 | infusion circuit |
| 6 | LED STOP | infusion circuit |
| 7 | NU | not used |
| 8 | GND | power supply |

2.3.2. Electronic layout .

See Annex 2.

2.3.3. Installation layout

See Annex 2.

2.4. FLEXIBLE LINK BOARD

2.4.1. Description of the connectors

2.4.1.1. J1 : link connector to the flexible plate

This connector joins the Pilot and the RS 232 input socket to the CPU board.

| Pin | Description | |
|-----|---------------------|------------------------------|
| 1 | GND | power supply |
| 2 | CTS | printer |
| 3 | TXD | printer |
| 4 | RTS | printer |
| 5 | RXD | printer |
| 6 | GND | power supply |
| 7 | CTS | option |
| 8 | TXD | option |
| 9 | RTS | option |
| 10 | RXD | option |
| 11 | VBAT | power supply |
| 12 | LED INF | hand switch Led anode |
| 13 | RETURN PATIENT CALL | hand switch LED cathode |
| 14 | CD PATIENT CALL | not in use |
| 15 | GND | power supply |
| 16 | RETURN ILS LOCK | hood contact |
| 17 | CD ILS LOCK | hood contact |
| 18 | VBAT | power supply |
| 19 | VBAT | power supply |
| 20 | GND | power supply |
| 21 | GND | power supply |
| 22 | OFF BUTTON | not in use |
| 23 | OPTO MOTOR | rotation motor |
| 24 | BUZZER | start command for the BUZZER |
| 25 | GND | power supply |
| 26 | TXD MASTER | Pilot |
| 27 | RXD MASTER | Pilot |
| 28 | ON BUTTON | ON/OFF switch |
| 29 | LED SUPPLY | supply presence |
| 30 | GND | power supply |

2.4.1.2. J7 : series link connector for printer

| Pin | Description | |
|-----|-------------|--------------|
| 1 | NC | not in use |
| 2 | RXD | printer |
| 3 | TXD | printer |
| 4 | NC | not in use |
| 5 | GND | power supply |
| 6 | NU | not in use |
| 7 | RTS U | printer |
| 8 | CTS | printer |
| 9 | NU | not in use |

2.4.1.3. J6 : optional series link connector

| Pin | Description | |
|-----|-------------|--------------|
| 1 | NC | not in use |
| 2 | RXD | option |
| 3 | TXD | option |
| 4 | NC | not in use |
| 5 | GND | power supply |
| 6 | NC | not in use |
| 7 | RTS | option |
| 8 | CTS | option |
| 9 | NC | not in use |

2.4.1.4. J2 : connector linking Pilot

| Pin | Description | |
|-----|-------------|-----------------|
| 1 | LED SUPPLY | supply presence |
| 2 | RXD MASTER | Pilot |
| 3 | GND | power supply |
| 4 | OPTO MOTOR | rotation motor |
| 5 | GND | power supply |
| 6 | VBAT | power supply |
| 7 | NU | not in use |
| 8 | NU | not in use |
| 9 | ON BUTTON | ON/OFF switch |
| 10 | TXD MASTER | Pilot |
| 11 | BUZZER | buzzer command |
| 12 | OFF BUTTON | stop button |
| 13 | GND | power supply |
| 14 | VBAT | power supply |
| 15 | NU | not in use |

2.4.1.5. J3 :

not installed

2.4.1.6. J5 : hand set connector

This connector is destined to be used in PCS mode.

| Pin | Description | |
|-----|------------------|------------------|
| 1 | CATHOD | Hand set LED |
| 2 | HAND SET CONTACT | hand set contact |
| 3 | HAND SET CONTACT | hand set contact |
| 4 | ANOD | Hand set LED |

2.4.2. Electronic layout .

See Annex : 2.

2.4.3. Implantation layout

See Annex : 2.

3. CONFIGURATIONS, CALIBRATIONS AND CHECK


3.1. Configurations

CPU board configuration

CPU board is configurable to work with an EPROM or FLASH EPROM program.
Configuration is carried out by the drops G1, G7 et G6.

3.1.1. Configuration EPROM

The board can be equipped with a 27C512 to 27C040 EPROM, access time : 100 ms.



| | | |
|----------|------|------------------------|
| G1 et G7 | open | EPROM use |
| G6 | 2.3 | EPROM 27C010, 020, 040 |
| G6 | 1.2 | EPROM 27C512 |

3.1.2. Flash EPROM configuration

The board is equipped of a 28F004 flash EPROM

| | | |
|----------|--------|-----------------|
| G1 et G7 | closed | FLASH EPROM use |
|----------|--------|-----------------|

For the other configurations see operator's guide.

3.2. Calibrations

Aimless

3.3. Check

3.3.1. Electrical safety tests

In compliance with EN 60 601.1 complies.

3.3.2. Integrated tests

The device has integrated auto-tests of the following components :



- Screen .
- LED's .
- Keyboard .
- Rotary knob .
- Serial links .
- Internal clock.


The tests can be perform with a Pilot CE 0459 or with an external power supply of 7 V via the 15 points sub D connector.


The following polarisation must be respected in this case. The pin 10 should be connected to the pin 2.

| Pin | Description | |
|-----|-------------|--------------|
| 6 | + 7 V | Power supply |
| 14 | + 7 V | Power supply |
| 5 | GND | Power supply |
| 13 | GND | Power supply |
| 2 | RXD | Pilot |
| 10 | TXD | Pilot |


3.3.3. Test mode


Press the two buttons simultaneously during the instrument auto-test SILENCE  + START .




The screen shown here invites the user to enter the test mode by pressing ENTER .

If the user doesn't quickly validate the entry into this mode, the instrument will return to the programming menu.



A turn of the Dial  allows the user to select the type of test or to display information.



Pressing STOP  allows the user to leave the test and return to the selection menu.

3.3.3.1. Functioning duration




This screen displays the following information :

- total instrument running time.
- total usage time since the last use.

The instrument assumes that the average length of a month is 30 days. Successive action on  and  allows one to change the maintenance date.


The time elapsed since the last maintenance is renewed after modification of the maintenance date.

3.3.3.2. Evolution dates



This screen displays the last evolution dates and time :
- date and time of loading the Master PCA application in flash memory.


3.3.3.3. Software version



This screen displays the numbers of the software versions :

- version, revision of the Master PCA application.
- date generated and the checksum of the software.



Pressing  displays the language screens :

- language.
- version, revision of the files as well as the creation date of the language.


All the different screens can be consecutively displayed using the rotary knob.

3.3.3.4. Screen LCD test



This test alternates between lighting all the pixels then every odd pixels, even pixels.


3.3.3.5. LED's test



This screen shows the status of the display LED's.

At the onset of this test, all the LED's are illuminated at the same time for 3 s.

3.3.3.6. Keyboard test




This screen shows the positions of the buttons according to the following order :

- squares 1 to 4 : line 0 column 0 to 3.
- squares 5 to 8 : line 1 column 0 to 3.
- squares 9 to 12 : line 2 column 0 to 3.
- squares 13 to 16 : line 3 column 0 to 3.

As soon as the button is pressed the corresponding symbol appears on the screen in the case according to how its matrix is organised. To stop the test, it is necessary to hold down the STOP button for more than 2 seconds.


3.3.3.7. Rotary knob test




This screen shows the number of impulses made by the rotary knob in rotating as well as an indication of the speed by means of a bar-graph. The sense of rotation is the same as that shown on the horizontal scale.

- a single bar corresponds to slow speeds.
- two bars correspond to greater speeds.

3.3.3.8. Patient switch test




This screen shows the state of the cover, programming key and patient switch
Important : to check Hand set functionality, the wiring layout must be respected.



- 1 : Patient switch LED cathode
- 2 : Patient switch contact
- 3 : Patient switch contact
- 4 : Patient switch LED anode


3.3.3.9. Serial link test



Before carrying out this test, it's necessary to fitted the device with test plugs made from subD 9 points with TxD and RxD links and a second link between CTS and RTS. This plug equip the RS232-1 links (PC link).

This screen shows the states of the different serial link lines. If one link is defective, the error message will appear.


3.3.3.10. Latest events



This test allows to check the last 10 events. Each event is display on one screen-page, each screen-page is selected with the rotary knob The events are numbered chronologically with the last as number 1.

The possible errors are out of 3 types : see §2.15 for details


3.3.3.11. Clock period



This screen displays the measurement of the internal period clock which is cyclically updated. For correct functioning of the instrument, the displayed value must equal 1.000 s.

3.3.3.12. Pilot information


This test extracts information on the latest 3 connected Pilots. Each Pilot number is displayed on one screen-page, and each page is selected using the rotary knob. The different screens are numbered chronologically with the most recent numbered 1.



This screen displays the following information :

- Pilot type.
- series number of the Pilot.
- software version of the Pilot.
- date and checksum of the Pilot EPROM.

3.3.3.13. Battery load



This is a timer to indicate the time spent during the menu enter after 96 hours charge, the bargraph is full and indicates that the RAM battery is fully charged.

To charge this battery:

- Connect the device to main.
- Go to test mode.
- Select : battery load
- Charge for at least 96 hours

Exit this menu reset and stop the timer.

4. REPLACING SUB-ASSEMBLIES

Important : a complete check of the instrument must be made after any internal investigation.

4.1. Replacing the electronic circuit boards

Important : be very careful with the flexible plate on opening the device.

1. Remove the 7 screws in the front face, and detach the fixing clamp.
2. Disconnect the LCD flat cable, taking care not to damage it.
3. Disconnect the connection cables from the keyboard connector, the display circuit connectors the coder connector and the flexible plate connector.
4. Take out the board, very carefully.
5. Replace the LED or the keyboard circuit if necessary.

4.2. Dismounting the support

Remove the two fixing screws on the hinge situated next to the CPU circuit board and pull.

4.3. Replacing the flexible circuit

Note : dismantling the flexible circuit systematically implies its replacement.

1. Dismantle the two connector support plates.
2. Detach the flexible circuit from the aluminium support.
3. Detach the connectors from their plate support.
4. Clean the glue from the surface with 95° methylated spirits.
5. Remove the protective sheet from the new flexible circuit, equipped with its connectors. Glue the flexible circuit making sure it is well positioned with regard to the 2 fixing screw-holes in the Pilot / Master PCA.
6. Glue on the protective film of the flexible circuit.

4.4. Replacing the handle

1. Unscrew the handle until it unclips from its socket and remove the screw from its thread.
2. Unscrew the plate enclosing the mechanism situated below the support.
3. Disengage the 2 wingnuts from their housing with pliers and pull out the set.
4. Insert a new handle into its socket.
5. Fix back the lock support plate or change the faulty parts of the lock.

5. MAINTENANCE

5.1. Cleaning and disinfecting

The Master PCA is part of the patient's immediate environment. It is advisable to clean and disinfect the device's external surfaces on a daily basis in order to protect patient and staff. Disconnect the device from its mains supply before starting to clean.

- Do not place in an AUTOCLAVE nor IMMERSE the device. Do not let liquids enter the device's casing.
- If the device is placed in a high contamination risk unit, it is advisable to leave it in the room during aerial disinfecting, after having disinfected it with a moist cloth.
- Use a cloth soaked in DETERGENT-DISINFECTANT, previously diluted with water if required, to destroy micro-organisms. Avoid abrasive scrubbing which could scratch the casing. Do not rinse or wipe surfaces.
- Do not use: TRICHLOROETHYLENE-DICHLOROETHYLENE - AMMONIA - AMMONIUM CHLORIDE - CHLORINE and AROMATIC HYDROCARBON - ETHYLENE DICHLORIDE-METHYLENE CHLORIDE - CETONE. These aggressive agents could damage the plastic parts and cause device malfunction.
- Take care also with ALCOHOL BASED SPRAYS (20% - 40% alcohol). They lead to tarnishing of and small cracks in the plastic, and do not provide the necessary cleaning prior to disinfecting. Using disinfecting applies by SPRAYS may be done, in accordance with the manufacturer recommendation, from a distance of 30 cm of the device, avoid the accumulation of the product in liquid form.

Please contact the appropriate service, handling suitable cleaning and disinfecting products, in your establishment for further details.

5.2. Storage

The device should be stored in a dry, cool place. In case of prolonged storage, the battery should be disconnected via the battery access flap situated underneath the device. This should be done by a qualified technician.

- Storage temperature: -10°C + 60°C.
- Permissive relative humidity: maxi 85%, no condensation.

5.3. Servicing

To ensure normal performance of the device, it is recommended to replace the internal battery each 3 years. This should be done by a qualified technician.

The qualified technicians in your establishment or our After-Sales Service should be informed if the device is dropped or if any of malfunction occurs. In this case, the device must not be used.

For further information concerning the pump servicing or its use, please contact our After-Sales Service or our Customer service.

If the device has to be returned to our After-Sales Department, proceed to its cleaning and disinfecting. Then, pack it very carefully, if possible in its original packaging, before sending it with a detailed description of the fault, to the official representative of Vial Medical.

Vial Medical is not liable for loss or damage to the device during transport to our After-Sales Department.

5.4. Regular inspections - Preventive maintenance

5.4.1. Before using checking

(See operator's guide of Master PCA : §14.2 page 28)

This check must be performed before every use of the Master

Note : In the framework of continuous improvement, this checklist may change any time. Please contact our after sales service for up-to-date version.(Addresses at the end of this document)

5.4.2. Preventive maintenance

(See next page : Technical check certificate)

In order to insure preventive maintenance, preventive technical check is recommended every 24 months. This technical check must be performed by qualified technician and is not covered by any contract from FRESNIUS VIAL.

For more information contact After Sale Service.(Addresses at the end of this document)

5.4.3. Internal historical saving battery

- Changing Internal historical saving battery is recommended every 2 years
- Changing this battery obliged to dismount CPU board following procedure described in § 4.1.
- Deweld battery, avoiding excess heating, short circuit and electrostatic charges when manipulate the board.
- The place of the battery is marked by a sticker.
- Prior placing a new battery make sure implant direction is correct and polarities are in accordance with serigraphy.
- Battery type is : 60 mA, CdNi, weldable Gb.
- Use procedure set up by battery maker to destroyed the removed battery.

Technical check certificate

--> See STK - Protokoll

| | |
|--|--|
| | |
| | |

①

| | |
|--|--|
| | |
|--|--|




⚡ ⏪

| | |
|--|--|
| | |
| | |

| | |
|--|--|
| | |
|--|--|

| | |
|--|--|
| | |
|--|--|

| |
|--|
| |
|--|



6. ANNEX 1 : ILLUSTRATED PARTS LIST

6.1. Traceability


6.1.1. Introduction

The aim of this chapter is to guide technicians looking for spare parts when servicing the device.

6.1.2. traceability table

| Serial number | From :....16230001 To :.....16470140 | From :.....1670141 To :.....1653270 | From.....16530271 |
|----------------|---|--|-------------------|
| UC Card | 182004 | 182004 | 182004 |
| Rotating coder | 182997 | 182997 | 182989 |
| Magnet | 182208 | 182232 | 182232 |
| Magnet support | 182408 | 182408 | 182408 |

6.2. General view



6.3. Mechanical parts list

6.3.1. Aluminium support

| Ref.: | Rep: | Quantity | Description |
|---------|------|----------|------------------------------|
| 182161 | 1 | 1 | External support plate |
| 182705 | 3 | 1 | Communication foam joint |
| 182106 | 4 | 1 | Moulded complete handle |
| 182705 | 5 | 1 | External support foam joint |
| 182157 | 6 | 1 | Pilot communication plate |
| 182205 | 7 | 1 | Aluminium PCA support |
| 182108/ | 8 | 1 | Handle screw end |
| 199593 | 9 | 2 | Inox M 5* 6 split screw |
| 199597 | | 2 | Sub d lg6 spacer |
| 199599 | | 10 | Black TF M 3*8 Taptite screw |
| 182779 | 12 | 1 | Flexible circuit protection |

6.3.2. Front panel

| Ref.: | Rep: | Quantity | Description |
|--------|------|----------|---------------------------|
| 182226 | 21 | 1 | Front panel hinge |
| 182703 | 22 | 1 | Hinge plate joint |
| 182702 | 23 | 1 | Plate band joint |
| 182124 | 24 | 1 | Master lever locker lever |
| 182155 | 25 | 1 | Front panel plate |
| 199587 | 26 | 5 | ff 10h3050-5 spacer |
| 199591 | 27 | 1 | Inox Z1 M 3*40 TCB screw |
| 182407 | 28 | 1 | PCA technical closer |
| 182128 | 29 | 1 | PCA front panel |
| 182012 | 30 | 1 | Anti-blink LCD screen |

6.3.3. cover

| Ref.: | Rep: | Quantity | Description |
|--------|------|----------|---------------------------------|
| 190096 | 43 | 2 | Nylon washer |
| 182401 | 44 | 2 | ACME 0180400 pillow-block |
| 199570 | 45 | 2 | CB M 4*16 screw |
| 199534 | 46 | 2 | Flat inox ZAN washer |
| 182408 | 47 | 1 | 57050 overmoulded magnet |
| 182223 | 48 | 1 | PCA magnet support |
| 199602 | 49 | 2 | Black TF 2.5 split screw |
| 182229 | 50 | 1 | Left pin |
| 182228 | 51 | 1 | Right pin |
| 199605 | 52 | 1 | Polyamide TF M2 * 6 Nylon screw |
| 182227 | 53 | 1 | PCA cover |
| 182207 | 54 | 1 | PCA locking plate |

6.3.4. Locker

| Ref.: | Rep: | Quantity | Description |
|--------|------|----------|-----------------------------|
| 182162 | 60 | 1 | PCA handle mechanism closer |
| 182222 | 61 | 1 | Lock bolt |
| 182224 | 67 | 1 | Cover hook |
| 182402 | 62 | 1 | Ronis 14700-22 locker |
| 182219 | 63 | 1 | Blade spring |
| 182215 | 64 | 1 | Locking spring |
| 182230 | 65 | 1 | Locker cam |
| 182214 | 66 | 1 | Ronis locker nut |
| 199559 | 72 | 2 | H5 8 * 4.2 M4 Nylstop nut |
| 199594 | 73 | 1 | Spacer |

6.4. Electronic parts list

6.4.1. Programmation pass

| Ref.: | Rep: | Quantity | Description |
|--------|------|----------|---------------------|
| 182225 | 68 | 1 | Pass socket hood |
| 182217 | 69 | 2 | Pass socket contact |
| 182204 | 70 | 1 | Programmation pass |
| 182403 | 71 | 1 | Bicolor LED |

6.4.2. Control and command

| Ref.: | Rep: | Quantity | Description |
|--------|------|----------|------------------------|
| 182834 | 31 | 1 | IEC / VIAL front panel |
| 182833 | 31 | 1 | DIN / VIAL front panel |
| 182830 | 31 | 1 | NEUTER front panel |
| 182023 | 32 | 1 | PCA key board |
| 182989 | 33 | 1 | HE 13 rotating coder |
| 182100 | 34 | 1 | Rotating knob |
| 182018 | 35 | 1 | Gelded LCD Flextrip |
| 182001 | 36 | 1 | MGLS graphical display |
| 182991 | | 1 | LCD / UC link |
| 182406 | 15 | 1 | Proximity sensor |

6.4.3. Electronic

| Ref.: | Rep: | Quantity | Description |
|--------|------|----------|-------------------------------------|
| 199598 | 37 | 5 | SSR 4-3-01 spacer |
| 182201 | 38 | 1 | PCA syringe LED integrated circuit |
| 182015 | 39 | 1 | Gelded alarm LED integrated circuit |
| 182004 | 40 | 1 | PCA CPU board |

6.4.4. Patient hand switch

| Ref.: | Rep: | Quantity | Description |
|--------|------|----------|----------------------------|
| 182051 | 11 | 1 | Patient hand switch socket |
| 182998 | 42 | 1 | Wired patient hand switch |

6.4.5. Connectors

| Ref.: | Rep: | Quantity | Description |
|--------|------|----------|------------------------------|
| 170689 | 10 | 1 | 9 pts sub d female connector |
| 182010 | 13 | 1 | Master flexible circuit |
| 170688 | 14 | 1 | 15 pts sub d male connector |
| 161761 | 16 | 1 | 6160.15.2D1200 connector |

7. ANNEX 2 : Installation and Electronic layouts

7.1. CPU BOARD

| Description | FV REF | REF | Rev. |
|----------------------------------|---------|-----------|------|
| Installation layout (# 2/2 x A3) | A301210 | D394GS004 | C2 |
| Electrical layout (# 7/7 x A3) | A301204 | D194GS000 | C2 |

7.2. ALARM LEDS BOARD

| Description | FV REF | REF | Rev. |
|---------------------------------|---------|-----------|------|
| Installation layout (#1/1 x A4) | A301334 | D395GU004 | B0 |
| Electrical layout (#1/1 x A4) | A301323 | D195GU000 | B0 |

7.3. INFUSION LEDS BOARD

| Description | FV REF | REF | Rev. |
|----------------------------------|---------|-----------|------|
| Installation layout (# 1/2 x A4) | A301400 | D395GU009 | C1 |
| Electrical layout (# 1/1 x A4) | A301327 | D195GU002 | C1 |

7.4. FLEXIBLE LINK BOARD

| Description | FV REF | REF | Rev. |
|--------------------------------|---------|-----------|------|
| Electrical layout (# 1/1 x A3) | A301199 | D195CR000 | B0 |
| Printed circuit (# 1/1 x A3) | A301201 | D195GU001 | B0 |

8. ADDENDA

9. Useful addresses

All request for information or documentation (technical file, tubing catalogue or commercial documentation) should be addressed to :

CUSTOMER SERVICE INTERNATIONAL

| | | |
|--|--|--|
| | Fresenius Vial Le Grand Chemin, 38590 Brézins FRANCE | Tel. : 33 (0)4 76 67 10 81 or 10 54 Fax : 33 (0)4 76 65 52 22 |
|--|--|--|

AFTER-SALES SERVICES

| | | |
|----------------------|--|---|
| INTERNATIONAL | Fresenius Vial Le Grand Chemin, 38590 Brézins FRANCE | Tel. : 33 (0)4 76 67 10 76 Fax : 33 (0)4 76 65 56 66 |
|----------------------|--|---|

| | | |
|----------------|--|---|
| BELGIUM | Fresenius NV/SA Belgique DIVISION VIAL MEDICAL Molenberglei 7 2627 Schelle BELGIUM | Tel. : 32/380 73 07 Fax : 32/880 50 07 |
|----------------|--|---|

| | | |
|----------------|--|---|
| GERMANY | FRESENIUS MCM AM-NEUNEN BERG 8 63749 ALZENAU GERMANY | Tel. : 49/60 23 97 22-0 Fax : 49/60 23 43 06 |
|----------------|--|---|



It is possible that this document contains typing errors or mistakes. Changes may occur at any time in subsequent editions.

COPYRIGHT © 1998, **Fresenius Vial S.A.**

This technical manual may not be reproduced in whole or in part without the written consent of **Fresenius Vial S.A.**

Fresenius Vial S.A. - siège social : Le Grand Chemin - 38590 Brézins (FRANCE)
S.A. à directoire et conseil de surveillance au capital de 90 128 000 FF - SIREN Grenoble B 408 720 282

 **Fresenius Vial**
Infusion Technology