

SUMMIT INDUSTRIES INC. 2901 W. Lawrence Avenue Chicago, IL 60625 1-800-972-9776

TABLE OF CONTENTS

1.0	INTRODUCTION	1
2.0	PRECAUTIONS	2
2	1 Installation	2
2	2 Operation	2
3.0	GENERAL PRODUCT OVERVIEW	4
3	1 Major Components	4
3	2 Specifications	7
3	3 Options	9
4.0	FUNCTIONS	10
4	1 Getting Started	10
4	2 Positioning the Tube	11
	3 Collimator Operation	
4	4 Setting Techniques	12
	5 Making Exposures	
4	6 Prep/Expose Tubestand Status Messages	15
4	7 Latching PREP	15
5.0	TECHNICAL DATA	16
	1 Radiation Leakage (Meets Federal Government Standards)	
	2 Rated line voltage: 240 VAC, 50 / 60 Hz., single phase	
	3 Range of line voltage regulation at max line current: 1% to 5%	
	4 Maximum line current (based on input voltage of 240 VAC): 140 AMP	
	5 Technique factors that constitute the max. line current: 300 mA @ 125 kVp	
5	6 Maximum Deviation of technique factors:	16
6.0	CLEANING THE InnoVet Select [™]	18
7.0	SCHEDULED DEALER MAINTENANCE	19
-	1 X-Ray Generator Maintenance	-
	2 X-Ray Tube Housing	
	3 Collimator	
7	4 Tubestand	20
7.	5 Other	20

8.0 T	ROUBLESHOOTING	
	Error Message	
	No Exposure	
8.3	Collimator light bulb burned out	
	No digital display	
	Circuit tripping.	
8.6	Full kVp range inaccessible	
	Poor film quality	
9.0 W	ARRANTY NOTICE	

1			
Revision	Pages Affected/Revision Description	Release Date	ECR #
A	RELEASED	-	-
В	CONVERTED TO MS WORD	MARCH 2002	-
С	UPDATED COMPONENT DIAGRAM	APRIL 2002	2564
D	ADDED "COLD" AND "READY" DISPLAY MESSAGES, EXPOSURE COUNT DISPLAY	OCTOBER 2002	2815
E	ADDED 50 HZ. REFERENCE AND TIME STATIONS	JULY 2006	5195

Installation and Service Manual Revisions History

1.0 INTRODUCTION

Congratulations! You've purchased an *InnoVet Select* Veterinary X-Ray System. The *InnoVet* line represents a new concept in veterinary radiography by being the first system to incorporate the x-ray generator into the table. This major innovation led to other significant design breakthroughs yielding the most practical, efficient and compact unit ever created. The result of this innovation and improvement is the *InnoVet Select* you have purchased.

Because the *InnoVet Select* is so unique, its features and functions are a bit different from those of traditional x-ray systems. The purpose of this manual is to familiarize you with the *InnoVet Select* and assist you with operating and maintaining your new system. The manual is divided into indexed sections for quick and easy access to information.

Should a situation arise which cannot be resolved through the troubleshooting procedures described in this manual, please contact the Summit dealer where you purchased your machine for assistance, or call Summit directly at 1-800-972-9776.

2.0 PRECAUTIONS

2.1 Installation

The *InnoVet Select* must be installed and serviced by an authorized Summit dealer. All warranties will be void if such installation or service is performed by persons not authorized by Summit.

- 2.2 Operation
 - 2.2.1 Radiation

X-rays are dangerous to both operator and others in the vicinity unless established safe exposure procedures are strictly observed.

The useful and scattered beams can produce serious, genetic or potentially fatal bodily injuries to any persons in the surrounding area if used by an unskilled operator. Adequate precautions must always be taken to avoid exposure to the useful beam, as well as to leakage radiation from within the source housing or to scattered radiation resulting from the passage of radiation through matter.

Those authorized to operate, test participate in or supervise the operation of the equipment must be thoroughly familiar and comply completely with the currently established safe exposure factors and procedures described in publications such as Sub-Chapter J of Title 21 of the Code of Federal Regulations "Diagnostic X-Ray Systems and their Major Components" and the National Council on Radiation Protection (NCRP) No.33, "Medical X-Ray and Gamma-Ray protection for energies up to 10 MeV-Equipment Design and Use", as revised or replaced in the future.

Failure to observe these warnings may cause serious, genetic or potentially fatal bodily injuries to the operator or those in the area.

2.2.2 Electrical

Only properly trained and qualified personnel should be permitted access to any internal parts. Live electrical terminals may be deadly. Be sure line disconnect switches are opened and other appropriate precautions are taken before opening access doors, removing enclosure panels, or attaching accessories. Do not remove the flexible high tension cables from the x-ray tube housing or high tension generator or the access covers from the generator until the main and auxiliary power supplies have been disconnected.

When disconnecting high-voltage cables, they must be grounded immediately in order to dissipate any electrical charge that may remain on the cables or the tube.

Failure to comply with the foregoing may result in serious or potentially fatal bodily injuries to the operator or those in the area.

2.2.3 Mechanical

All of the movable assemblies and parts of x-ray equipment should be operated with care. Only properly trained and qualified personnel should be permitted access to any internal parts.

Particular care should be taken when servicing the inside of the tubestand. There is an extreme threat of mechanical pinching between the vertical slide and counterweight due to their close proximity and opposite directions of motion.

3.0 GENERAL PRODUCT OVERVIEW

3.1 Major Components



AWK287-01

- a. Tubestand Supports collimator and x-ray tube. Moves both horizontally and vertically.
- b. X-Ray Tube
- c. Angulation Dial Displays degree of x-ray tube angulation.
- d. Angulation Lock When loosened, allows angulation of x-ray tube.
- e. Handle Bars Used to move tube vertically or to angle the tube.
- f. Light Switch Button Turns collimator bulb on and off.
- g. Source-Image Distance (SID) Indicator When marker aligns with label on tubestand, SID (the distance between the tube focal spot and the film) is 40". There is a position for tabletop work and one for grid cabinet applications.
- h. Grid Cabinet Lead-lined cabinet holds cassette tray and grid.
- i. Grid, mounted to upper carriage of grid cabinet, absorbs scatter radiation for improved film clarity.
- j. Cassette Tray Molded indentions hold cassettes in correct position.
- k. Cassette Box (optional) Opens to store cassettes. Lead lined to protect film.
- I. Rustproof Support Feet Adjust to level the table.
- m. Exposure Foot Switch Used to make exposures.
- n. Release Foot Pedal When depressed, opens front of table to access control panel.
- o. Table Base All steel, welded construction.
- p. Circuit Breaker Turns generator on/off (Located on the side of the tilt-out panel).
- q. kVp Dials One for major adjustment (15 kVp increments); One for minor (2 kVp increments).
- r. mA Dial Adjusts mA. Two stations; 300 and 100 when using small tube focal spot and 300 when using large. Optional third mA station: either 400 (large focal spot) or 30 (small focal spot).

s. Time Dial - Adjusts time. 18 stations available:

Line Frequency	Time Stations Available			
50 Hz.	1/100, 2/100, 3/100, 4/100, 5/100, 6/100, 8/100,			
	1/10, 12/100, 18/100, 24/100, 30/100, 36/100,			
	48/100, 6/10, 9/10, 1 2/10, 1 1/2			
60 Hz.	1/120, 1/60, 1/40, 1/30, 1/24, 1/20, 1/15, 1/12,			
	1/10, 1/8, 3/20, 1/5, 1/4, 2/5, 1/2, 5/8, 3/4, 1			

- t. Drain Trough Minimizes spillage behind table.
- u. Table Top (float table top) 24" x 53" of Formica (33" x 57" of Formica).
- v. Tubestand/Grid Cabinet Interlock When engaged, allows tubestand, grid cabinet and cassette tray to move together.
- w. Digital Displays Upper one for kVp; Center one for mAs, and Lower one for status & error messages.
- x. Collimator Controls size and rectangular shape of primary beam emitted from x-ray tube and provides coincident light field.
- y. Shutter Adjustment Knobs One each to determine width and depth of primary beam.
- z. Swivel Mounting Plate Allows collimator to be rotated.
- aa.Lamp Switches (optional) Turns Collimator bulb on and off. Located on the front and both ends of the table.
- bb. Float Top Unlock Button (optional). Allows transverse travel of the table top. Located on the front of the table at both corners.

3.2 Specifications

3.2.1 INTEGRATED GENERATOR

- Tilt-out control.
- Dual-digital timer with 18 time stations (1/100 to 1 1/2 second for 50 Hz. or 1/120 to 1 second for 60 Hz.).
- Adjustable solid state audio exposure signal.
- SCR contactor.
- Programmable tube protect system.
- Tube overload sensor.
- Independent kVp and mA adjustment circuits for each mA station.
- Two or three mA stations: 100S and 300L (30S, 400L Optional).
- mA overload protection.
- kVp selection, major 15 kVp increments, minor 2 kVp increments.
- mA, time and kVp knobs.
- Built-in kVp compensation.
- Easy-access exposure count display
- Filament stabilization.
- 16-foot line cable; 4-foot control-to-transformer cable.
- Integral 24-Volt A.C. power supply.
- Programmable line frequency 50 or 60 Hz.
- Programmable power saving mode.
- Latching prep option.

3.2.2 TABLE / TUBESTAND

- Enclosed, integrated table. Removable 24" x 53" top.
- Welded base construction.
- Full-travel, lead-lined grid cabinet.
- 103-line, 8:1 aluminum grid.
- Grid cabinet/tubestand interlock.
- Tubestand-mounted digital kVp, mAs and status displays.
- Attached, counterbalanced tubestand; travels full table length.
- Angulation dial and operator handles.
- Two position foot switch.

3.2.3 CASSETTE HOLDER

- Easy-fit, removable, ABS plastic holder.
- Accepts 8" x 10", 10" x 12" and 14" x 17" cassettes in either direction.
- 3.2.4 TUBE
 - ✤ 1.0-2.0 mm focal spots; 140,000 H.U.

3.2.5 CABLES

✤ 10' Max-Flex high voltage cables. Federal terminations.

3.2.6 COLLIMATOR

- Manual with light field.
- Swivel-mount.
- 90-second bulb "on" time.
- Gradual turn off warning (Soft Start Bulb Life feature).

3.3 Options

3.3.1 TABLE OPTIONS

- Exposure foot switch.
- Extra lamp switches.
- ✤ Float top 33" x 57" or Float top 33" x 67".
- 12" table-top extension. Specify right or left (fixed top only).
- Upgrade grid to 10:1, 103 line.
- Animal restraining device.
- Third mA Station 400 (large filament) or 30 (small filament).
- Cassette Storage Bin.

3.3.2 TUBE / COLLIMATOR OPTIONS

- ✤ 0.6-1.5 mm; 200,000 H.U.
- Programmable lamp times 30, 60, 90, and 120 seconds.

3.3.3 OTHER OPTIONS

- Remote exposure hand switch.
- Standard manual tray with sliding jaws; accepts all cassette sizes up to 14" x 17".

4.0 FUNCTIONS

- 4.1 Getting Started
 - 4.1.1 Opening Unit Depress the Release Foot Pedal to open the front access door.
 - 4.1.2 Power Switch To turn the generator on, move the circuit breaker switch to the up position. The switch is located at the left side of the hinged front door.

It is a good idea to turn the unit off at the end of each day. This ultimately will extend your system's life. If power saving mode is selected, the unit turns off automatically.

- 4.1.3 Generator Warm up When the generator is first turned on after being turned off for two hours or more, the tubestand display will indicate "COLD". After ten minutes this display will change to "READY" which will scroll across the tubestand status window. The generator will take exposures when "COLD", but the outputs will be most accurate when the display indicates "READY". For consistent films wait for the "READY" indicator prior to taking exposures.
- 4.1.4 Tube Warm up It is recommended that you warm up the unit prior to the first higher-powered exposure to prolong tube life. This is done by taking a series of three, 1/10 second exposures using the 300 mA station. Start with a kVp of 70 and increase to 80; then 90. Exposures should be spaced at least 20 seconds apart.
- 4.1.5 Storing Cassettes Cassettes can be stored in the optional storage box or boxes mounted to either side of the table.
- 4.1.6 Loading Cassettes Three sizes of cassettes can be used: 8" x 10", 10" x 12" and 14" x 17". Simply place the cassette into the specially-formed grooves of the tray. It will accept the cassette in either lengthwise or crosswise position.

IMPORTANT:

Although the storage box is lead-lined, it will not protect film from fogging after excessive exposure.

- 4.2 Positioning the Tube
 - 4.2.1 Moving the Grid Cabinet The grid cabinet can be moved in tandem with the tubestand, assuring alignment of the primary beam to the center of the film. The grid cabinet can be positioned anywhere along the entire length of the table. Grasp the tubestand handle bars and move to the desired location.
 - 4.2.2 Disengaging the Grid Cabinet Alternatively, the grid cabinet may be moved independent of the tubestand. Release the interlock by sliding the spring-loaded lever extending from the grid cabinet to the left. The grid cabinet may then move while the tubestand remains stationary. To re-engage the interlock, position grid cabinet in-line with tubestand. The interlock is engaged when you hear a metallic "click" sound and feel the interlock mechanism move into place.
 - 4.2.3 Vertical Movement The tube can be moved from 10-1/2" to 40" above the table top. For convenience, the SID markers indicate 40" positions from the table top or the grid cabinet. During vertical motion, you will note a pause at one height, which is the pre-set height for optimal exposure, set by adjusting the detents during installation.
 - 4.2.4 Float Top (option) Push in float top unlock button and move the table. Release the button to lock the top position.
 - 4.2.5 Angulation Release the lock, which is located on the right side of the tube arm. Grasp the handle bars and turn the tube to proper angle. Turn the lock handle clockwise to hold tube position. The degree of rotation is indicated by the dial located between the handle bars.
- 4.3 Collimator Operation
 - 4.3.1 Light Field Once the tube is positioned relative to the grid cabinet, press the light button on the front of the collimator or one of the optional "no hands" collimator switches on the front or ends of the table. A light field will appear. WARNING! DO NOT LOOK DIRECTLY INTO LIGHT SOURCE AS RETINA DAMAGE COULD RESULT.

- 4.3.2 Programmable Lamp Times Lamp gradually turns off in one of the following selectable times: 30, 60, 90, 120 seconds. The factory setting is 90 seconds. If you push the button during the lamp's power down ramp, a new "on" cycle starts.
- 4.3.3 Collimator Adjustment Use the sliding knobs to move the shutters. Scales indicate knob positions that correspond to common cassette sizes. Use the pointer to quickly move shutters into relative position, then fine tune to collimate the x-ray beam to the area of radiographic interest.
- 4.3.4 Swivel Movement When doing table top work, you may want to swivel the collimator to avoid awkward patient positioning. Simply swivel the collimator to the preferred position. WARNING: Do <u>not</u> touch the black lamp cover at the rear of the collimator as it gets quite hot. The collimator hesitates at the 0° and 90° positions for correct re-alignment. This is in order to help find these position accurately.
- 4.4 Setting Techniques
 - 4.4.1 mA Adjust mA by moving the pointer to one of two (or three) stations: 100S or 300L (30S, 400L).

When set at 100 mA, the small tube focal spot (1.0 mm in standard tube) is used, yielding a more finely detailed image. At 300 mA, the large focal spot (2.0 mm in standard tube) is used.

4.4.2 kVp - Two knobs control kVp - one for coarse adjustments of 15 kVp increments and one for fine tuning with 2 kVp adjustments.

Turn the knobs to select the desired kVp.

The kVp value is displayed on the tubestand-mounted display.

marked both in fractions of seconds and in fratural numbers.		
Line Frequency	Time Stations Available	
50 Hz.	1/100, 2/100, 3/100, 4/100, 5/100, 6/100, 8/100,	
	1/10, 12/100, 18/100, 24/100, 30/100, 36/100,	
	48/100, 6/10, 9/10, 1 2/10, 1 1/2	
60 Hz.	1/120, 1/60, 1/40, 1/30, 1/24, 1/20, 1/15, 1/12,	
	1/10, 1/8, 3/20, 1/5, 1/4, 2/5, 1/2, 5/8, 3/4, 1	

4.4.3 Time - The time dial can be set at any of 18 stations; these are marked both in fractions of seconds and in natural numbers:

4.5 Making an X-Ray Exposure.

Note: If it is desired, to view the total number of X-ray exposures, simply turn the KVp major, KVp minor & mA selector knobs fully counterclockwise, the bottom display will indicate "EXPO". The KVp display and the mAs display will indicate the total number of exposures taken on this x-ray control.

4.5.1 Making Exposures with the Standard Foot Switch

There are two ways to make a x-ray exposure with the 2-position footswitch.

Prep and Expose

When proper patient position is achieved, partially depress the footswitch to activate the x-ray tube prep cycle. When "PREP" is displayed on the Tubestand, watch the patient for respiration phase or proper restraint. At the desired moment, make the x-ray exposure by fully depressing the footswitch. An audible tone will indicate an exposure has been made.

Expose

When proper patient positioning is achieved, fully depress the exposure footswitch. The x-ray tube will begin its prep cycle and upon completion of the prep cycle, the x-ray tube will make an exposure. An audible tone will indicate an exposure has been made. Note that in this procedure an x-ray exposure will not be made immediately upon depressing the footswitch. A delay will be experienced due to the x-ray tube prep cycle.

4.5.2 Making Exposures with the Optional Foot Treadle

First depress the Exposure Foot Treadle until the display reads "WAIT." This starts the tube prep cycle. After 1 to 2 seconds the tubestand display will read "PREP." The tube is now ready for an exposure. The exposure will start when you depress the Exposure Foot Treadle. The generator will instantaneously activate and sound an audible "beep" to indicate that an exposure was made.

If the Exposure Foot Treadle is initially depressed and then held, the unit will automatically go through its two-second prep cycle and then expose. Hence, you will note a slight pause between pressing the foot switch and hearing the exposure signal.

- 4.6 Prep/Expose Tubestand Status Messages
 - WAIT Tube is not yet warmed up. Please wait about 2 seconds. Exposures are not allowed yet.
 - **PREP** Tube is ready. Exposure starts immediately after initiating with foot switch.
 - **X-RAY** X-Ray exposure is in progress. Radiation is present.
 - **DONE** Exposure is completed.
 - **COLD** Unit needs to warm up for about ten minutes to achieve the most accurate exposure outputs.
 - **READY** Unit is warmed up and ready for taking exposures.

EXPOSURE COUNT

With the kVp major, kVp minor and mA selector knobs turned fully counter-clockwise, the bottom display will indicate "EXPO". The KVp display and the mAs display will indicate the total number of exposures taken on this x-ray control.

- 4.7 Exposure Foot Treadle/Latching PREP (option)
 - 4.7.1 Press exposure foot switch or treadle momentarily and release.
 - 4.7.2 **WAIT** is displayed for 2 seconds while tube is warmed up.
 - 4.7.3 **PREP** is displayed for 20 seconds. Exposure can be made any time during this interval.
 - 4.7.4 For exposure push foot switch or treadle down and hold for duration of exposure.
 - 4.7.5 If after 20 seconds no exposure is made, **PREP** is released.

5.0 TECHNICAL DATA

- 5.1 Radiation Leakage (Meets Federal Government Standards)
 - 5.1.1 From the tube housing: less than 50 mR/hr at 1m from focal spot (@4mA, 125kVp)
 - 5.1.2 From collimator: less than 50 mR/hr at 1m (@4mA, 125 kVp)
- 5.2 Rated line voltage: 240 VAC, 50 / 60 Hz., single phase.
- 5.3 Range of line voltage regulation at max line current: 1% to 5%.
- 5.4 Maximum line current (based on input voltage of 240 VAC): 140 AMP
- 5.5 Technique factors that constitute the max. line current: 300 mA @ 125 kVp.
 - 5.5.1 Generator rating:
 - 5.5.1.1. Output current: 100-300 mA.
 - 5.5.1.2. Output voltage: 40-125 kVp.
 - 5.5.2 Generator duty cycle:
 - 5.5.2.1. 100 mA @ 125 kVp: 4%
 - 5.5.2.2. 300 mA @ 125 kVp: 1%
- 5.6 Maximum Deviation of technique factors:

The maximum deviation statements provided below are broad and considered "worst case" conditions.

Typical performance will normally be significantly better.

5.6.1 Peak Tube Potential (kVp)

The Maximum Deviation of the Peak Tube Potential during an exposure from its pre-indicated value is + 5 kVp.

5.6.2 Tube Current

The Maximum Deviation of the tube current during an exposure from its pre-indicated value is + 10% of selected value within a kVp range of 50-125 kVp (inclusive), and + 20% of selected value within a kVp range of 40-50 kVp (exclusive).

5.6.3 Exposure Time

The Maximum Deviation of the measured exposure time from its pre-indicated value:

Max Deviation	Time Range
± 1 pulse	1/100 – 12/100 sec. for 50Hz. 1/120 – 1/10 sec. for 60Hz.
	1/120 – 1/10 sec. for 60Hz.
± 15%	18/100 – 1 ½ sec. for 50 Hz.
	1/8 – 1 sec. for 60 Hz.

6.0 CLEANING THE INNOVET SELECT™

The *InnoVet Select* was designed to withstand the sometimes extreme conditions of a busy veterinary practice. In addition to the drip trough and rust-proof table legs, you'll note that there are no grooves or holes to collect debris.

To periodically clean the unit, use a sponge and a non-abrasive cleanser. Knobs on the control panel are removable for this purpose.

The system must be located in a different room than the processor. Failure to observe this could lead to premature corrosion not covered by warranty.

7.0 SCHEDULED DEALER MAINTENANCE

The following schedule of maintenance for the *InnoVet Select* is required for safety of operation, continued ease of use, and continued long life of the product.

The maintenance program should only be performed by service personnel authorized by Summit Industries. Frequency of the service should be 30 days after installation and annually thereafter.

- 7.1 X-Ray Generator Maintenance
 - 7.1.1 Verify that all technique selector knobs are mounted solidly.
 - 7.1.2 Verify proper operation of kVp and mAs digital displays.
 - 7.1.3 Verify proper operation of audio and visual x-ray termination.
 - 7.1.4 Verify time, mA, and kVp accuracy.
 - 7.1.5 Examine all electrical connections of the high voltage generator.
 - 7.1.6 Examine oil level of high voltage transformer and refill if fluid level is greater than 3/4" from the top of the lid.
 - 7.1.7 Examine the high voltage cable ends for carbon tracking.
 - 7.1.8 Examine inside of power supply box for any possible fluid leakage.
- 7.2 X-Ray Tube Housing
 - 7.2.1 Inspect for possible oil leakage.
 - 7.2.2 Assure that the housing is tightly fastened to the tube mount and collimator.
 - 7.2.3 Inspect stator cable for fraying or damage.
 - 7.2.4 Inspect high voltage cable ends for carbon tracking. Clean and Regrease HV Cable ends annually.

- 7.3 Collimator
 - 7.3.1 Verify accuracy of field size.
 - 7.3.2 Verify accuracy of light field to x-ray field alignment.
 - 7.3.3 Check lamp on-off switch for proper operation.
 - 7.3.4 Inspect collimator cable for fraying or damage.
- 7.4 Tubestand
 - 7.4.1 Verify set-screws are securely holding tubearm.
 - 7.4.2 Inspect counterweight cables for fraying or damage.
 - 7.4.3 Inspect all tubestand movements for binding or interferences.
- 7.5 Other

Inspect high voltage cable and all other cabling for damage.

8.0 TROUBLESHOOTING

8.1 Error Message

Upon error detection, the *InnoVet Select* will display one of the following error codes:

- **FOOT** Early release of the exposure foot switch or foot treadle during exposure. Films are likely to turn out light. Foot must be depressed for the entire selected exposure time.
- **TUBE** Single shot tube power limits were exceeded by the selection technique. Reduce kVp, mA and / or time (PREP and EXPOSE locked out). Message clears when the techniques change.
- **HEAT** Tube heat exceeded. PREP and EXPOSE locked out until technique reduced or tube cools down.
- **ROTR** Tube rotor error. Call for service.
- **ZCO** Equipment component (hardware) error. Call for service.
- **ZZZZ** (Power saver mode only.) Unit will go to sleep within 5 minutes.
- FIL1-FIL4 Hardware error. Call for service.
- **COLD** The unit needs to warm up for best output accuracy, but exposures can be taken if needed. This display should clear about ten minutes after first turn on of the day.
- 8.2 No Exposure

If, after fully depressing the foot switch, there is no audible tone, your unit may not be turned on. Use the circuit breaker to turn on power. If there is still no exposure, call a service representative.

8.3 Collimator light bulb burned out.

When the collimator light bulb burns out, you may either call a service person for replacement/realignment or may choose to replace the bulb yourself.

The replacement bulb must be FCS type; 150 Watt at 24 V with a pin base.

WARNING #1:

DO NOT LOOK DIRECTLY INTO THE COLLIMATOR LIGHT SOURCE AS RETINA DAMAGE COULD RESULT.

WARNING #2:

TURN THE CIRCUIT BREAKER OFF BEFORE ATTEMPTING BULB REPLACEMENT TO AVOID POSSIBLE ELECTRIC SHOCK.

First, remove the back cover of the collimator. The bulb is now visible.

WARNING #3:

THE BULB, WHEN LIT, GETS VERY HOT. ALLOW IT TO COOL SEVERAL MINUTES TO AVOID BURNING YOUR FINGERS.

Pull the bulb straight out from the ceramic base.

WARNING #4:

DO NOT MAKE DIRECT CONTACT WITH THE REPLACEMENT BULB AS OIL FROM YOUR FINGERS WILL CREATE HOT SPOTS ON THE BULB WHICH SUBSEQUENTLY MAY LEAD TO CRACKING.

To avoid direct contact between bulb and fingers, the replacement bulb is normally packaged in a cellophane wrapper. Tear open one end so that prongs are exposed; hold the glass portion, still in its wrapper, and insert the prongs into the ceramic base. Make sure that each prong is securely fastened into its socket so that the resulting light field is straight. Replace the back cover.

If the light field is out of alignment with the grid cabinet, call an authorized service representative for re-alignment.

8.4 No digital display

If, after turning on the circuit breaker, digital displays do not light up, this indicates lost power. The source may be a blown fuse within the unit or a problem with incoming line power. Call an authorized service representative.

8.5 Circuit tripping

If the circuit breaker shuts off on its own, this indicates mA overload. This could be due to miscalibration of the mA or another malfunction.

In the event that the unit is miscalibrated, you may still be able to operate the unit by moving to 100 mA. This will allow you to take X-Rays until an authorized service representative is called.

8.6 Full kVp range inaccessible.

If this occurs, it means that incoming line voltage is too low. Call an authorized service representative to correct.

8.7 Poor film quality

There are many interrelated factors in the radiographic process which impact the degree of darkening visible on a film. The following is a list of items to check whenever you notice a significant change in film quality:

- 8.7.1 Film Processing Has chemistry been changed recently? Is a new brand/type being used? Have developing temperatures varied? All of these affect film quality.
- 8.7.2 Technique Were proper technique factors used? Consult with a Summit dealer or call Summit for information on setting up a technique chart.

- 8.7.3 Screens/Cassettes Screens within the cassettes degrade with time. They typically lose 10% of their light-emitting capability per year. It is recommended that cassettes be re-screened every five to seven years. All screens should be replaced at one time to maintain uniform exposure.
- 8.7.4 Film Has a new brand/type/speed of film been used? Is the film beyond its expiration date? If so, new film may be purchased or technique factors may need adjustment.
- 8.7.5 Foot switch If thorax/chest films are blurring more than usual, it could be that the exposure switch is being used incorrectly. See section 4.5 to review how to make exposures.

9.0 WARRANTY NOTICE

The *InnoVet Select* is warranted as follows: any Summit manufactured part proving defective will be repaired or replaced free of charge, F.O.B. factory, if the defective part is returned for inspection, freight charges prepaid. This warranty covers parts only and does not include any on-sight labor costs.

This warranty is five years from the date of installation.

This warranty does not apply to collimator glassware or x-ray tubes or damage caused by accident, misuse, neglect, or shipment, and is void if service is performed by persons other than authorized Summit dealers or if equipment is interconnected with components not approved by Summit for compatibility.

The x-ray tube manufacturer's warranty is one year, pro-rated. Summit will extend this warranty to five years, pro-rated on manufacturer's list price at time of claim, provided the tube is purchased and used exclusively with an *InnoVet* **Select** system.

Please call your Summit dealer to coordinate the return of a warranty part or call Summit direct at 1-800-972-9776.

INNOVET SELECT[™] Veterinary Radiographic

System

Installation & Service Manual

K286

Rev J

July 2006

SUMMIT INDUSTRIES INC. 2901 W. Lawrence Avenue Chicago, IL 60625 773-588-2444 www.summitindustries.net

Table of Contents

MANUAL REVISION CHART1-1
UNPACKING2-1
RADIATION AND MECHANICAL/ELECTRICAL WARNING
EQUIPMENT SPECIFICATIONS
PRELIMINARY MECHANICAL ASSEMBLY
5.1RADIOGRAPHIC TABLE5-15.2HIGH VOLTAGE TRANSFORMER5-15.3TUBESTAND INSTALLATION5-25.4TUBE ARM, X-RAY TUBE, CABLES AND COLLIMATOR5-85.5FOOT TREADLE INSTALLATION5-95.6FILM BIN INSTALLATION5-95.7GENERATOR INSTALLATION5-10
ELECTRICAL CONNECTIONS
6.1 LOCATION OF ELECTRICAL CONNECTIONS6-16.2 MAKING CONNECTIONS6-3
CALIBRATION
7.1PRECALIBRATION CHECKS.7-17.2MA CALIBRATION7-37.3KVP CALIBRATION.7-67.4.CALIBRATION OF KVP WITH A MEANS TO MEASURE ACTUAL KVP7-77.5CALIBRATION OF KVP WITHOUT A MEANS TO MEASURE ACTUAL KVP7-77.6CONTROL BOARD DIP SWITCH SETTINGS7-117.7X-RAY TERMINATION BEEPER VOLUME CONTROL7-127.8TUBESTAND DISPLAY MESSAGES7-127.9DETAILED TROUBLE SHOOTING GUIDE7-13

TABLE OF CONTENTS

8.0	FINAL ADJUSTMENTS AND QUALITY ASSURANCE	8-1
	 8.1 LEVELING THE TUBEARM	8-1 8-1 8-1
9.0	OPTIONS FOR FIELD INSTALLATION	9-1
	 9.1 FILM BIN 9.2 TABLE TOP EXTENSION 9.3 ANIMAL RESTRAINTS	9-1 9-1
10.0	LOCATION OF FIGURES WITHIN THE TEXT SCHEMATICS	

1.0 MANUAL REVISION CHART

TABLE 1

InnoVet Select™ Installation and Service Manual Text Revisions (K286-T)			
Revision Date Description of Revision		<u>Ecr No.</u>	
А	OCT 1998	REFORMATTED AND EXPANDED TEXT/FIGURES	-
В	JUNE 1999	REVISED TEXT IN SECTION 9.1	-
С	JUNE 2000	ADDDED 4 TH POTENTIOMETER FOR KVP CALIBRATION	-
D	FEB 2002	CONVERTED TO MS WORD	-
E	APRIL 2002	REVISED TUBEARM INSTALLATION/LEVELING	2564
F	JULY 2002	REVISED INSTALLATION INSRUCTIONS FOR ADDING TRIM WEIGHT TO THE TUBE- STAND COUNTERBALANCE	2667
G	OCT 2002	ADDED RC TIMER BOARD, EXPOSURE COUNTER INFORMATION	2815
Н	OCT 2004	EXPANDED TROUBLESHOOTING GUIDE	4335
J	JUL 2006	ADDED 50 HZ. REFERENCE AND TIME STATIONS	5195

2.0 UNPACKING

The *InnoVet Select*[™] Radiographic System is delivered in three packages:

- High voltage transformer, control, x-ray tube (palletized).
- Table, collimator and accessories.
- Tubestand.

It is the installer's responsibility to inspect the shipment for damage and proper count. Upon receipt of the merchandise, any visible damage to the cartons should immediately be examined while the shipper is present. If the visible damage to the cartons also includes damage to the merchandise, the installer is responsible for making all claims with the shipping company.

If there is hidden damage to the merchandise, it is the installer's responsibility to discover that damage within a reasonable amount of time and contact the shipping company.

3.0 RADIATION AND MECHANICAL/ELECTRICAL WARNING

3.1 RADIATION

WARNING

X-rays are dangerous to both operator and others in the vicinity unless established safe exposure procedures are strictly observed.

The useful and scattered beams can produce serious, genetic or potentially fatal bodily injuries to any persons in the surrounding area if used by an unskilled operator. Adequate precautions must always be taken to avoid exposure to the useful beam, as well as to leakage radiation from within the source housing or to scattered radiation resulting from the passage of radiation through matter.

Those authorized to operate, test, participate in or supervise the operation of the equipment must be thoroughly familiar and comply completely with the current established safe exposure factors and procedures described in publications such as Sub-Chapter J of Title 21 of the <u>Code of Federal Regulations</u>, "Diagnostic x-ray Systems and their Major Components," and the National Council on Radiation Protection (NCRP) No. 33, "Medical x-ray and Gamma-Ray protection for energies up to 10 MeV-Equipment Design and Use," as revised or replaced in the future.

Failure to observe these warnings may cause serious, genetic or potentially fatal bodily injuries to the operator or those in the area.

3.2 ELECTRICAL

WARNING

Failure to comply with the following may result in serious or potentially fatal bodily injuries to the operator or those in the area.

Only properly trained and qualified personnel should be permitted access to any internal parts. Live electrical terminals may be deadly; be sure line disconnect switches are opened and other appropriate precautions are taken before opening access doors, removing enclosure panels, or attaching accessories.

Do not remove the flexible high voltage cables from the x-ray tube housing or high voltage transformer or the access covers from the generator until the main and auxiliary power supplies have been disconnected.

When disconnecting high voltage cables, they must be grounded immediately in order to dissipate any electrical charge that may remain on the cables or the tube.

3.3 MECHANICAL

WARNING

Particular care should be taken when servicing the inside of the tubestand. There is an <u>extreme</u> threat of mechanical pinching between the vertical slide and counterweight due to their close proximity and opposite directions of motion.

The control is secured to the table with an air cylinder, hinge pin and Truarc retaining ring only. Once the ring and hinge pin are removed, care must be taken to insure that the control does not swing forward in an uncontrolled manner.

All of the movable assemblies and parts of x-ray equipment should be operated with care. Only properly trained and qualified personnel should be permitted access to any internal parts.
4.0 SPECIFICATIONS

4.1 TECHNICAL RATINGS

TABLE 2 – TECHNICAL RATINGS

		[]
	300mA Station	400 mA Station (Optional)
<u>Rated Line Voltage</u> -	240 VAC, 50 / 60 Hz., Single Phase.	240 VAC, 50 / 60 Hz., Single Phase.
Acceptable Line Voltage Regulation at Maximum Line Current	Not to exceed 5%.	Not to exceed 5%.
Maximum Line Current (240 VAC input)	140A @ 300mA	187A @ 400mA
<u>Technique Factors that</u> <u>Constitute the Maximum</u> <u>Line Current</u>	300 mA @ 125 kVp.	400mA @ 100kVp
<u>Control Rating</u> – Output Current	100 mA and 300mA are standard. Three mA stations of 30, 100, and 300 mA are available as an option.	Three mA stations of 100, 300, and 400 mA are available as an option.
<u>Control Rating</u> – Output Voltage	40 to 125kVp	40 to 125 kVp
Control Duty Cycle –	100 mA @125 kVp-4% 300 mA @ 125 kVp-1%	100 mA @125 kVp-4% 300 mA @ 125 kVp-1% 400 mA @ 100 kVp-1%

4.2 SPACE REQUIRMENTS

1. For STANDARD Fixed-Top tables



2. For tables with Float-Top option



4.3 ELECTRICAL REQUIREMENTS

TABLE 3 – ELECTRICAL POWER SUPPLY REQUIREMENTS

	300mA Station	400 mA Station (Optional)
Equipment Category	300mA @ 125 kVp, Single Phase400mA @ 100 kVp. Single Phase	
Nominal Line Voltage	240 VAC, 50 / 60 Hz., Single Phase	240 VAC, 50 / 60 Hz., Single Phase
Line Voltage Range Allowed (Alternate line voltages which will provide normal operation.)	194-284 VAC. 194-284 VAC.	
<u>Maximum Momentary</u> Line Current	140 amperes (at 240 VAC)	187 amperes (at 240 VAC)
Note: Maximum momentary line current at alternate line voltages can be determined using the formula:	e current at V_2 V_2 V_2 oltages can V_2 = alternate line voltage V_2 = alternate	
Line Voltage Regulation under load	The line voltage drop under load is not to exceed 5% at maximum line current.	The line voltage drop under load is not to exceed 5% at maximum line current.

Calculating line voltage regulation	$\frac{(V_{NL} - V_{L})}{V_{L}} X100$ V_{L} $V_{NL} = \text{ line voltage under}$ "no load" conditions $V_{L} = \text{ line voltage under}$	$\frac{(V_{NL} - V_{L})}{V_{L}} X100$ V_{L} $V_{NL} = \text{ line voltage under "no load" conditions}$ $V_{L} = \text{ line voltage under "full}$
	"full load" conditions	load" conditions
Minimum Over Current Protection Rating (Service disconnect)	50% of maximum line current rating or greater. (100 Amps is recommended.)	50% of maximum line current rating or greater. (100 Amps is recommended.)
Distribution Transformer Requirements	Minimum 30 kVA dedicated to the x-ray control	Minimum 37.5 kVA dedicated to the x-ray control
Wire Size from Power Transformer to Disconnect Switch	for 50 feet use #2 AWG; for 100 feet use #00AWG; for 200 feet use 250 mcm	For 50 feet use #2 AWG; for 100 feet use #00AWG; for 200 feet use 250 mcm

The information provided in Table 3 above and text below is taken from NEMA standards. <u>Minimum Power Supply Requirements for x-ray Machines</u>, and <u>National</u> <u>Electric Code</u>.

Connection to Supply Circuit (taken from N.F.P.A. 70-1984)

A disconnecting means of adequate capacity for at least 50% of the input required for the momentary rating shall be provided in the supply circuit. The disconnecting means shall be operable from a location readily accessible from the x-ray control. Underwriters Laboratories also requires that this disconnecting means to be mounted on either the wall behind the x-ray control or the wall directly adjacent to it.

5.0 PRELIMINARY MECHANICAL ASSEMBLY

5.1 RADIOGRAPHIC TABLE

- 1. During shipment, the table top is banded to the table base. Cut loose the banding and set the table top aside for assembly later. Loosen the hold-down screws underneath the front rail of the table to release the grid cabinet. Remove packing material from rear of grid cabinet.
- 2. Install the five leveling feet by threading them into the holes at the two front corners of the table and at the ends and middle of the back rail. The two feet with a 4" long body should be installed at the front corners of the table. The three feet with a 1" long body should be installed at the back rail. (In the accessory bag there are six blue and six red "U" shaped plastic spacers. When leveling the table fully tighten the foot until the shoulder of the foot is firmly against either these spacers or the bottom of the table/back rail). Use of the spacers makes the table base/leveling feet union a solid "one piece" structure, minimizing any wobble or table base motion.
- 3. Position the table close to its final location, using the leveling feet and spacers to bring the table close to level at this time. Final leveling can be done at end of assembly.

5.2 HIGH VOLTAGE TRANSFORMER

1. Verify that the oil level in the HV transformer is approximately ½" to ¾" below the underside of the cover. Refill with "Shell Diala oil AX" or remove oil as necessary.

- 2. Using a piece of cardboard to protect the floor, push the transformer onto the transformer shelf in the table by tilting the transformer and sliding it up onto the shelf. Once the transformer is in place, loosen the screw on the vent plug to allow for oil expansion. Cover the vent plug and flange with one of the loose fitting plastic covers provided with the HV transformer to prevent contaminants from entering the oil.
 - **Note:** The high voltage transformer must be installed before the control is put in place. It is also recommended that the transformer be installed before the tubestand is mounted to the table for maximum stability during installation.

5.3 TUBESTAND INSTALLATION

The sequence of assembly is slightly different between Fixed-Top and Float-Top tables.

For Fixed-Top (standard) tables proceed directly to section 5.3.1 below.

For Float-Top (optional) tables the table top frame must first be removed before the tubestand can be mounted. Skip ahead to section 5.3.2 for instructions for this type of table.

5.3.1 For Standard Fixed-Top Table

1. Remove tubestand stop from the left end of the table's lower rear rail. See figure 3.



Figure 3: Tubestand stop at lower left rail



Figure 4 – Tubestand Positioning

2. Position the tubestand so that it is a few inches from the end of the back rail of the table with the lower tubestand bearings in line with the back rail. See Figure 4.

(The tubestand shown in Figures 4 and 5 is for a Float-Top table, but mounting the tubestand to the rail is done the same for either configuration.)



Figure 5 – Mounting Tubestand on Rail

 Carefully tilt the tubestand away from the table until the elevated tubestand bearing can be pushed into the lower bearing rail. Now straighten the tubestand and push it completely into the upper and lower bearing rails. See Figure 5.

4. Reinstall the tubestand stop on the lower rail. If leveled properly in step 5.1, the tubestand should glide easily from end to end and stay in position without drifting to either side.

- **5.3.2** For Optional Float–Top Table (Refer to photos on preceding pages for reference.)
- 1. Mount the tubestand bearing carriage to the tubestand using the four (4) ¹/₄-20" x 3/8" hex head screws provided. Figures 4 and 5 on page 5-3 show a tubestand with a mounted bearing carriage. Verify that the safety rollers are mounted on the tubestand bottom to prevent the table from tipping backwards from the weight of the tubestand. (This may occur if the high voltage transformer is not yet installed.)
- 2. To mount the tubestand, the table top frame will need to be moved forward for access to the upper tubestand bearing rail at the left end of the table. Refer to Figures 6 and 7 below for reference.
 - A) Remove the single, large Phillips head screw found at the front and rear interior corners of the table top frame. Their location is shown in Figure 6 (front only shown).



Figure 6 – Table Top Frame Mounting Screws

B) Now that the table top frame is free from the bearing tracks, you can raise the frame slightly to gain access to the four small Phillips head screws which secure the table top frame's bearing track at the exterior head end of the table. Remove these four screws, sliding the bearing track forward. See Figure 7.



Figure 7 – Table Top Raised; Bearing Track Moved Forward For Mounting Tubestand

C) Manually release the transverse lock solenoid pins and slide the table frame forward to expose the upper tubestand bearing rail.

- Remove the tubestand stop on the left end of the table's lower rear rail. (As viewed from the front.) <u>See Figure 3 Tubestand Stop</u> on page 5-3 for reference. Position the tubestand so that it is a few inches from the end of the back rail of the table with the lower tubestand bearings in line with the back rail.
 See Figure 4 Tubestand Positioning, Page 5-4 for reference.
- 4. Carefully tilt the tubestand away from the table until the lower tubestand bearing can be pushed into the bearing rail. Now straighten the tubestand and push it completely into the upper and lower bearing rails. See Figure 5 Mounting Tubestand on Bearing Rail, Page 5-4 for reference. Reinstall the tubestand stop on the lower rail. A well leveled table will allow smooth motion of the tubestand with no drift in either direction.
- 5. Reassemble the table frame by repeating step 2 above in reverse. Replace the tubestand stop on the lower rail. If not done previously, remove the phenolic table top panel at this time via the nuts located at the underside corners of the table top frame.

5.4 TUBE ARM, X-RAY TUBE AND COLLIMATOR INSTALLATION

- 1. Install the tube mount by inserting the 1-1/2" diameter tube arm of the tube mount into the collar on the face of the tubestand. Be sure to fully insert the tube arm. Secure the tube arm by tightening the two set screws on the outside diameter of the collar.
- **Note**: Do not remove the two hex head shipping bolts which lock the main counterweight and vertical carriage together until the tube, collimator and high voltage cables have been mounted.
- 2. Mount the x-ray tube and collimator. When using the standard tube and collimator, two 1/8" spacers will achieve the proper field size coincidence. The standard tube has a part number of C458, and the collimator G800 (refer to their respective manuals for detailed instructions as needed). Refer to Final adjustments Section 8 for central ray to grid cabinet adjustments after x-ray to light field verification and generator calibration is completed later in this manual.
- 3. Install the high voltage cables. Cable drape brackets are included with the hardware accessories bag. The cable drape brackets are installed above and below the handle on the right side of the tubestand. Loosely drape the HV cables with the collimator and stator cables within the brackets until the vertical slide has been freed to move. The cables will route directly into the right rear of the table base, not beneath the lower tubestand rail.

<u>WARNING</u>

Particular care should be taken when servicing the inside of the tubestand. There is an **extreme** threat of mechanical pinching between the vertical slide and counterweight due to their close proximity and opposite directions of motion. All of the movable assemblies and parts of x-ray equipment should be operated with care. Only properly trained and qualified personnel should be permitted access to any internal parts.

4. Remove the shipping bolts which lock the main counterweight and vertical carriage together, taking care to insure that the counterbalance is adequate to avoid uncontrolled motion of the vertical slide. The trim weights provided in the hardware kit should be used to achieve proper counterbalancing of the tubestand. Remove the back panel of the tubestand, and then insert trim weights as needed into the channels at the top of the main counterweight.

5.5 FOOT TREADLE INSTALLATION (skip ahead to 5.6 if using standard foot switch)

If the *InnoVet Select*[™] has the latching prep Foot Treadle Option, the treadle **must be installed prior to mounting the control into the table.** (If the Foot Treadle Option is being added to an installed system as a field upgrade, refer to Section 9 "Options for Field Installation".)

- 1. Put the treadle in position across the front of the table base.
- 2. Thread the spring guide pins located at either end of the treadle up through the table base and lock each guide pin in place with the cotter pin provided.
- 3. Bolt the foot treadle bracket to the table base using the hex head bolts provided.
- 4. Plug the foot treadle cable into the S7 location of the junction box at the right side of the table frame. Refer to Figure 14 as needed.

5.6 FILM BIN INSTALLATION

(skip ahead to 5.7 if no film bin is to be installed)

The film bin must be installed prior to mounting the control into the table base. The film bin comes in either a right-hand or left-hand configuration. Slide the PEM studs of the film bin through the holes in the front edge of the appropriate side of the table and secure it using the hardware provided. Bolt the back side of the film bin to the rear of the table using the bolts provided. Install the plastic organizer in the top drawer of the film bin. Note that for left-hand film bins the control serial tag may be partially obscured. Be certain to record the control's serial number in the user and installation manual if this is the case.

5.7 GENERATOR INSTALLATION

- 1. Lay the control face down on cardboard, centered in front of the table.
- 2. Adjust the table legs to measure 34-5/8" between the V grooves of the pivot rollers at the base of the table.
- 3. Engage the notches on the bottom of the control with the pivot rollers and lift the control so that it hinges up into the table base.
- 4. Connect the air cylinder to the control using the pin and locking ring as shown in Figure "9". Note that there is no other device holding the control within the table. If the pin is ever to be removed from the air cylinder, exercise caution to maintain balance of the control within the table.



Figure 9: Air Cylinder Attachment – AWW430-06

5. Bolt the foot loop to the bottom of the control as shown in Figure "10".



Figure 10: Mounting the Foot Loop – AWW430-07

6. Center the control within the table opening by adjusting the two front feet. If the left side of the control is too close to the left pedestal leg, raise the left front corner of the table slightly. If the right side of the control is too close to the right pedestal leg raise the right front corner slightly. Use the red or blue shims as required to achieve a snug fit between the top of the leveling feet and the bottom of the table for maximum stability.

6.0 ELECTRICAL CONNECTIONS

6.1 LOCATION OF ELECTRICAL CONNECTIONS

1. Remove the technique select knobs, each of which are secured to their respective switch shafts by two Allen set screws. Remove the two Phillips head screws with rubber gaskets at the top of the control cover. Remove the two shoulder bolts from the bottom corners of the control cover and lift cover away as shown in Figure 11 below.



Figure 11 Removing the Front Cover

Once the front cover has been removed, the interfacing terminal strips within the control will be accessible as shown in Figure 12 below.



Figure 12 – Interior of Control

2. If not already done, remove the phenolic table top panel for access to the rear of the control. The table top panel can be removed by a hex nut at each underside corner of the table top. The Line, Primary, Filament, Collimator, Stator, Foot Switch and Table Interface cables will all enter the control as shown in Figure 13 below.



Figure 13 – Rear of Control Showing Cable Entry Points

3. There is a "Table Interface Junction Box" where the cables for all Options will plug in. (Refer to Figure 14 below.) There will be one cable from this Junction Box which will enter the control and connect to TB - 5, acting as the interface between the control and all possible table options. The connections will be described in detail in section 6.2.



Figure 14 – Table Interface Junction Box

6.2 MAKING CONNECTIONS

- 1. Connect the Primary cable (labeled P1, P2 and GND) and the Filament cable (labeled XS, XL, XC and GND) to the High Voltage transformer terminals which carry the same markings.
- Route the collimator cable (labeled CSW1, CSW2, C1, C2 and GND) through one of the openings in the rear of the control and connect it to TB 4 terminals 21, 22, 23, 24 and 25 respectively, which are labeled the same as the cable. (Refer to Figure 15 for reference).

- 3. Route the tubestand display cable out of the tubestand and into the black connector on the back of the generator. The connectors are keyed for proper orientation. This cable, as the stator and collimator cables, must not limit nor be strained by tubestand motion.
- 4. Route the x-ray tube stator cable (labeled T5, T6, 07, 08 and 09) through one of the openings in the control and connect to TB -4 terminals 31, 32, 33, 34 and 35 respectively, which are labeled the same as the cable. Refer to Table 4 and Figure 15 for reference. Remove jumper between terminals T5 & T6 of TB4 when connecting tubes thermal cut-out SW.

DESCRIPTION		COLOR CODE FOR STATOR CABLE		
TERMINAL TB4 MARKINGS	CONNECTIONS	EUREKA TUBES	MACHLETT TUBES	SUMMIT & TOSHIBA TUBES
Τ5	THERMAL CUTOUT SWITCH	ORANGE/ BROWN		YELLOW
Т6	THERMAL CUTOUT SWITCH	BROWN/ ORANGE		BLUE
07	RUN WINDING OF TUBE STATOR	BLACK	BLACK	BLACK
08	START WINDING OF TUBE STATOR	RED	GREEN	RED
09	COM. LEAD OF TUBE STATOR	WHITE	WHITE	WHITE

Table 4 – Common Stator Cable Designations

5. If the system is using the standard two-position Foot Switch, route the cable (labeled S1, S2 and S3) through one of the openings in the control and connect to TB – 4 terminals 27, 28 and 29 respectively, which are labeled the same as the Foot Switch cable. (Refer to Figure 15 below for reference.)



Figure 15 - TB3 and TB4 Terminal Designations

- 6. If the latching prep Foot Treadle option has been installed, push the air hoses from the foot pads on the air cylinder nozzles protruding from the rear of the control. Plug the two conductor cable from the treadle into the Table Interface Junction Box at S4 (refer to Figure 14 above).
- At the Table Interface Junction Box plug the interface cable into S7. Route this cable through an opening in the rear of the control and connect it to TB - 5 as show in Table 5 below. (Note that there are no connections made to terminal 1 or 3 on TB - 5).

S7 Interface Cable Wire Color	Goes To Control
Brown	TB 5 - 2
Blue	TB 5 - 4
Red	TB 5 - 5
Green	TB 5 - 6
Orange	TB 5 - 7
Violet	TB 5 - 8
Black	TB 5 - 9
White	TB 5 - 10

Table 5 Interface Cable From Table to Control

- 8. All options (such as table locks, collimator hip switches, latching prep foot treadle) plug into the Table Junction Box. The Junction Box is designed so that each cable can only plug into the proper connector. If a plug will fit into more than one connector then it is a parallel circuit and can go in either connector. Refer to the Electrical Schematic Section 10 for detailed information about how the Junction Box is wired and how it fits into the overall electrical system as needed.
- 9. Connect the high voltage cables, taking note of the following important points:
 - A) Be careful of the three pins at each end of the cables as they can break if the cable is dropped.
 - B) Route the cables into the rear of the table so that they do not limit or restrict tubestand motion and they are not stressed by the tubestand motion.
 - C) The pins at the end of the HV cables should not be pinched tightly together but be slightly spread apart, about the thickness of a utility knife blade. (For old style split-pins only).
 - D) The HV cable terminals must be thoroughly cleaned and then coated with Vapor Proofing Compound (normally provided with x-

ray tube) prior to insertion into the x-ray tube or HV transformer receptacles. It is common to use dielectric oil in the HV transformer, and vapor proofing compound in the x-ray tube.

- E) Be sure that "Anode" of HV transformer connects to "Anode" of xray tube, and "Cathode" to "Cathode". Insert the HV cable terminals into the appropriate receptacle and screw the cable nut as tightly as possible by hand – Do not use tools for tightening.
- **NOTE:** Be certain to verify that the control's on/off switch is in the "OFF" position, and that the main incoming power at the service disconnect is also in the "OFF" position. It is strongly recommended that a meter be used to confirm no voltage is present before connecting the Line cable to the service disconnect switch.
- 10. Route the line cable under the lower tubestand rail and wire it into the safety disconnect. Connect "L1" to one leg of the incoming line power, "L2" to the other incoming leg of line power and "G" to the ground lug.

7.0 CALIBRATION

7.1 PRECALIBRATION CHECKS

- 1. It is recommended that the assembler read and understand the information provided with the x-ray tube prior to making any x-ray exposures. Particular attention should be given to:
 - A) Initial seasoning of the x-ray tube
 - B) Single exposure tube ratings
 - C) Accumulated heat and anode ratings of tube
- 2. Line voltage adjustments:

Terminals on TB2 are provided for coarse and fine adjustments of line voltage. Measure the line voltage at the disconnect switch and relocate wires marked "LVAC" and "LVAF" to two terminals where the sum of the two terminals markings equal the measured line voltage, ± 2 VAC. For example: If line voltage is 240 VAC, connect "LVAC" to "224" and "LVAF" to "+16".

3. Prior to turning "on" power, set each selector switch as follows:

Power "On-Off"	"Off"
kVp major and minor	Fully Counterclockwise
mA selector	300L
Time selector	1/120 sec. for 60 Hz. (or 1/100 sec. for 50 Hz.)

TABLE 6 – Selector Switch Settings

4. In order to prevent accidental production of x-rays during initial check-out, disconnect the leads marked P1 and P2 from TB 1 and install a jumper from TP – 13 to TP – 4 of J400 board (this will hold the SCR open). Switch the power safety disconnect switch "On" and then switch the control's "on-off" switch to "On." Observe the following:

mAs display indicates	2.5
kVp display indicates	<40

- 5. Rotate kVp knobs so that display indicates >40 kVp. While observing the tube filaments through the port of the x-ray tube, verify that the correct filament is lit for each mA station. The small focal spot should illuminate for 30 or 100 mA, and the large focal spot should illuminate for 300 or 400 mA.
- 6. Verification of Rotor circuit.
 - A) <u>Traditional Foot Switch</u> Depress the foot switch to the 1st stage (PREP), or the "PREP" button on the optional remote switch, and verify anode rotation and filament boost. The bottom status display should show "WAIT" momentarily, then "PREP".
 - B) <u>Foot Treadle with Latching Prep</u> Depress and release the foot switch. The bottom status display on the tube stand should show "WAIT", then "PREP". Verify anode rotation and filament boost. The unit will automatically remain in PREP mode for about 20 seconds, then return to idle condition.
- 7. The kVp meter circuit is factory adjusted, but display accuracy should be checked prior to other calibration steps.
 - 1. With power "Off" connect an AC voltmeter capable of reading 0-300 VAC across the wipers (center terminals of tap switches) of the kVp selector switches (major and minor) of TB3 terminals 9 and 10. This will allow measurement of the "no-load" primary voltage.
 - 2. Turn on power and adjust the kVp selector switches to achieve 220 VAC on the voltmeter. Refer to Table 8 for correct kVp display values
 - 3. Reset the kVp tap switches to achieve 160 VAC on the meter. Refer to Table 9 for the correct kVp display values.

TABLE 0 - RUP/IIIA Readings @ 220 VAC		
mA Station	kVp display on tubestand	
30 (if included)	116	
100	114	
300 87		
400 (if included)	70	

TABLE 8 – kVp/mA Readings @ 2	20 VAC	
-------------------------------	--------	--

TABLE 9 - KVp/ma Readings @ 160 VAC		
mA Station	kVp display on tubestand	
30 (if included)	82	
100	80	
300	53	
400 (if included)	-	

TABLE 9 –	kVn/mA	Readings	0	160 VAC
	N P/IIIA	Reduings	9	IUU IAU

If the kVp meter varies by more than 5 kVp to the above values, proceed to the kVp calibration procedure **7.4** before adjusting mA.

Switch the disconnect switch to "Off" and reconnect the leads P1 and P2 to TB 1 on the control. Be sure to remove the jumper between TP – 13 and TP – 4 on the J400 PCB.

7.2 mA CALIBRATION

Achieving maximum accuracy of tube current (mA) involves two types of adjustments:

- Overall mA level (bands of filament resistor RX), and
- mA balance throughout the useful kVp range (bands of space charge compensating resistor RSCC).

The leads connected to the resistor bands of resistor RX (the third resistor on the right of the three at the right corner of the Control Board) are marked with letters which correspond to the mA stations (small-SM, medium-MD, large-LG). Moving one of these bands upwards increases the mA for the corresponding mA station.

The leads connected to the bands of RSCC are also marked with SM, MD, LG which correspond to the mA stations. Moving one of these bands upwards increases the space charge compensation, which increases mA at low kVp, and reduces mA at high kVp. Moving one of these bands down reduces the space charge compensation, decreasing mA at low kVp and increasing mA at high kVp.

Normal mA tracking will result in values slightly high at 80 kVp and mA values roughly equal to one another but slightly low at the high and low ends of the kVp range (125 and 50). See Figure 16 below for a visual representation of this effect.

(2) terminals, of TB3, are marked "mA" (term #19 & 20). Remove the factory installed jumper from these 2 terminals & connect an mA_{AC} meter to these terminals to measure tube current.



Figure 16 – mA Tracking Example

For the following steps, use mA values and band designations in the table below:

Table 10 -	mA Stations/Band L	<u>_abels</u>

mA Station	<u>100/300</u> Configuration	<u>30/100/300</u> Configuration	<u>100/300/400</u> <u>Configuration</u>
30mA		SM	
100mA	SM	MD	SM
300mA	MD	LG	MD
400mA			LG

7.2 mA Calibration (continued)

In the following procedures, use the SM, MD or LG bands as appropriate for the mA stations available per this unit's configuration according to Table 10 above. Calibration is typically started with the 100mA station.

- 1. Connect a Dynalyzer or **AC** mAs meter to perform the following procedure. An AC mAs meter can be inserted into the circuit by removing the jumper between TB3-19 "MA~" and TB3-20 "MA~" and connecting the meter to those terminals (**AC readings only**). Be certain to replace the jumper after removing meter.
- 2. Select 100 mA and 80 kVp.
- 3. Adjust the appropriate band of filament resistor RX for the mA station chosen to produce the desired mA. Moving the band upward increases mA.
- 4. Select 50 kVp and note the mA produced.
- 5. Select 120 kVp and note the mA produced.
- 6. If the mA at 120 kVp is higher than the mA at 50 kVp, move the band of RSCC appropriate for the selected mA station upward. If lower, move the band downward. Achieve relative balance per Figure 16 above. The unit is capable, by fine tuning, to track within <u>+</u>5% of selected value.
- 7. Repeat steps (3) through (6) until no further adjustments are necessary.
- 8. Repeats steps 2 through 7 above for the one or two remaining mA options available on your station.
- 9. <u>Replace jumper</u> after removing the AC mAs meter.

7.3 kVp CALIBRATION

The *InnoVet Select*TM kVp compensation circuit has four offset potentiometers (one for each possible mA station), and four slope potentiometers. These adjustment pots are located on the J400 control board. See Figure 17 – Electrical Control Board Component Identification for component orientation below.



Figure 17 – Control Board Component Identification

kV POTS (8 TOTAL)

7.4 kVp CALIBRATION WITH A MEANS TO MEASURE ACTUAL kVp

Note: In order to make a valid comparison between the actual kVp output and the kVp meter's preread indication, the tube current must be calibrated accurately per section 7.2. Failure to do so will result in miscalibration of kVp.

For the following steps, use values in the table 11 below"

Adjustment Pot	mA Configuration
P1	30 mA offset
P2	100 mA offset
P3	300 mA offset
P10	400 mA offset
P4	30 mA slope
P5	100 mA slope
P6	300 mA slope
P11	400 mA slope

TABLE 11 – Adjustment Pot Identification

In the following procedure, fine tuning of the kVp display calibration depends on the assumption that mA has been calibrated.

- 1. Offset Adjustment:
 - a. Select 50 kVp and 100 mA. While monitoring mA and kVp output, compare the measured kVp to the kVp meter's readings.
 - b. Adjust the kVp offset pot (from Table 11 above) that corresponds to the mA of your Station (For a 100 mA station the pot is P2) until the kVp meter's preread indication matches the actual measured kVp output.
- 2. Slope Adjustment:
 - c. Select 120 kVp. Make an x-ray exposure and adjust the slope pot (See Table 11 above) until the kVp meter's preread indication matches the measured kVp output. Repeat steps 1 and 2 until no further adjustments are necessary.

- 3. Verify that the kVp display tracks well throughout the kVp range.
- 4. Repeat steps 1, 2 and 3 above for each mA Station, using Table 11 as a guide for identifying the appropriate adjustment potentiometer.

The offset and slope adjustments are interdependent; as an adjustment in the offset is made, it will affect the slope, and as an adjustment is made in the slope it will affect the offset. By anticipating the interactive response it will reduce the number of adjustments required to achieve the desired accuracy of kVp display at each end of the kVp range.

Under normal conditions the kVp meter's indication will not deviate from the measured kVp output by more than 5 kVp.

7.5 CALIBRATION OF kVp WITHOUT A MEANS TO MEASURE ACTUAL kVp

The preceding Tables 10 - (mA Stations/Band Labels) and 11 - (Adjustment Pot Identification) contain settings and readings referred to in the following section. In the following procedures, an mA will be chosen, a kVp tap switch voltage will be set, and the kVp display will be adjusted to indicate the value shown in the column on the left side of Table 12 below.

Typically if the 50 kVp and highest (120 or 100) kVp are set in this manner, all values in between will also be correct. The tap switch voltages for kVp outputs between 50 and 100 kVp shown below are not needed for this procedure, but are provided for reference.

kVp	30 mA Station Unloaded Tap Switch Voltage	100 mA Station Unloaded Tap Switch Voltage	300 mA Station Unloaded Tap Switch Voltage	400 mA Station Unloaded Tap Switch Voltage
50	104	108	155	183
60	121	125	173	201
70	139	143	190	218
80	156	160	208	236
90	174	178	225	253
100	191	195	243	271
110	209	213	260	-
120	226	230	278	-

TABLE 12 – TYPICAL KVP TAP SWITCH VOLTAGES

7.5 CALIBRATION OF kVp WITHOUT MEANS TO MEASURE ACTUAL kVp (continued)

- 1. Connect an AC voltmeter capable of measuring 0-300 VAC to the common terminals of the minor and major kVp tap switches, or to TB3 terminals 10 and 11.
- 2. Select an mA station. Adjust the kVp tap switches to achieve the voltage which should result in 50 kVp per Table 12 above. (For example, at the 100 mA station adjust the tap switches to achieve 108 VAC).
- 3. Adjust the kV offset pot shown in Table 11 for the mA station selected. (For 100 mA offset adjust P2), until the kVp meter reads 50 kVp
- 4. Readjust the kVp tap switches to achieve the voltage which should result in 120 kVp per Table 12. (For example, for the 100 mA station adjust the tap switches to achieve 230 VAC.)
- 5. Adjust the appropriate slope pot until the kVp meter reads 120 kVp.
- 6. Repeat steps 2 through 5 until no further adjustments are required.
- 7. Repeat these steps, adjusting tap switch voltages, then offset and slope pots, until the kVp display indicates the appropriate value across the operating range for each mA station.

7.6 CONTROL BOARD DIP SWITCH SETTINGS

The InnoVet SelectTM has a variety of dip switch settings to enable or disable various functions and to configure the system for the type of x-ray tube used in the system. These switches are typically set at the factory, but can be changed to accommodate the desires of the operator or the needs of the application. The dip switch combinations are shown in Figure 18 below.



DIP SWITCH SETTINGS

Figure 18 – Dip Switch Settings for Configuring Options and Features

7.7 X-RAY TERMINATION BEEPER VOLUME CONTROL

(See Figure 17 – Electrical Control Board Component Identification)

The volume of the x-ray termination beeper can be adjusted by P7, located on the control board. Turning P7 counter-clockwise increases the volume from quiet to very loud. Note: The front cover of the control will quiet the sound considerably so test the sound with the cover in place.

7.8 TUBESTAND DISPLAY MESSAGES

The following Messages are displayed by the bottom display on the InnoVet Select[™] tubestand. The COOL and READY features can be disabled via DIP2-4 (see FIG 18).

- **WAIT** PREP has been requested but the tube is not warmed up yet. The message clears as soon as the tube is ready for an exposure. This typically takes between one and two seconds.
- **PREP** The tube is warm and now ready for an exposure.
- **XRAY** Exposure is occurring. Ionizing radiation is present.
- **GOOD JOB** Exposure completed successfully. The message is cleared when the exposure switch is released.
- **TUBE** Indicates that the singe shot tube power limits are exceeded by the selected techniques. It is accompanied by an error beep and prep is not allowed. The message clears when allowed techniques are selected.
- **EXPOSURE COUNT** With the kVp major, kVp minor and mA selector knobs turned fully counter-clockwise, the kVp display the bottom display will indicate "EXPO". The mAs and the kVp display will indicate the total number of exposures taken on this x-ray control.
- **HEAT** Indicates tube heat interlock. The message is accompanied by an error beep and cleared as soon as the x-ray tube cools. PREP is locked out.

- **ZCO** Line voltage zero crossing signal not detected. Accompanied by error beep. PREP locked out. Message clears as soon as zero crossing signal detected.
- **FOOT** Foot switch released before exposure is finished. Accompanied by error beep. The message is cleared when any of the techniques change or PREP is requested.
- **ROTR** Tube rotor error. Accompanied by error beep. PREP is terminated. The message is cleared when any of the techniques change or PREP is requested.
- FIL1 Standby filament current too low, or not present. Accompanied by error beep. PREP locked out. Message clears when standby filament current is acceptable.
- FIL2 Standby filament current too high. Accompanied by error beep. PREP locked out. Filament circuit disabled. Message will not clear without rebooting system and resolving cause of high idle current.
- **FIL3** Boost filament current too low. Accompanied by error beep. PREP locked out. The message is cleared when any of the techniques change, PREP is requested, or cause of low boost current is found and resolved.
- FIL4 Boost filament current too high. Accompanied by error beep. PREP locked out. The message is cleared when any of the techniques change, PREP is requested, or cause of high boost current is found and resolved.
- **COLD** The control is cold. Exposures are possible, but the outputs will not be at their maximum accuracy. "COLD" will change to "READY" once the unit is on for about ten minutes.
- **READY** READY scrolls across the tubestand display once the unit has warmed up to a level that delivers the most accurate outputs, about ten minutes after initial turn-on.
- **ZZZZ** Sleep timer is 5 minutes away from turning off system. Accompanied by error beep. Active in power saving mode only.

7.9 Detailed Troubleshooting Guide for K200 Generator

I. CIRCUIT BREAKER TRIPS AT TURN ON

- A. Confirm mechanical latching of breaker
 - 1. Turn OFF main input power at the wall breaker
 - 2. Flip the control circuit breaker ON and OFF several times to confirm the breaker will mechanically latch in the ON position. If not, replace breaker.

B. Check for shorted or leaking SCR

- 1. With the control turned OFF, select DEMO mode by setting DIP2-5 of J400 board to the ON (up) position.
- 2. If the control turns on, look for voltage between P1 and P2. If more than 4 VAC is present, replace the SCR. If control does not turn on see C below.
- C. Confirm valid kVp input and J400 board
 - 1. If the breaker still trips when in DEMO mode, disconnect the wire from terminal W2 on the J400 board and turn the control ON. If the display is stuck in "EXPO" mode, the kVp tap switch input voltage is missing at the J400 board header H2. Check the wire-to-plug and plug-to-board connections at H2, and the kVp tap switches for open/resistive contacts or connections.
 - 2. If the kVp display is not stuck in "EXPO" mode, and shows a valid selection between 40 kVp and 125 kVp, replace the J400 board.

II. CIRCUIT BREAKER TRIPS AT EXPOSURE

The circuit breaker will trip due to either excessive mA or excessive line current. It may appear to be related to a specific mA station, or occur above a specific kVp range, or occur at longer time stations.

- A. Determine if the mA calibration is valid
 - 1. Remove the jumper between TB3-19 and TB3-20 and connect a mAs meter between these points, ensuring that the meter is set to read AC mA. *The range of the meter scaling must be set to maximum.*
 - 2. Typical filament voltage between XC/XS for 100 mA is 39 VAC during PREP (+/- 2 VAC). Typical filament voltage between XC/XL for 300 mA is 51 VAC during PREP (+/- 2 VAC). Under no circumstances does this voltage need to be more than 56 VAC.
 - 3. Take exposures and measure the mA or mAs output. Readjust the filament voltage as needed to bring mA into range.
- 4. If the filament voltage is correct but the output is ½ the expected mAs the SCR or high voltage transformer may be half-waving. See II.D below.
- B. Determine if the line voltage drop is excessive (this test is only valid after the mA calibration has been verified).
 - Place a voltmeter (set for min/max reading or to the hundreds of volts scale) between L1 and L2 of the control. Record the no-load line voltage. Measure and record the line voltage during an exposure of 90 kVp, 300 mA, 300 mAs. The drop should be less than 7%. If the lights dim, the supply wires rattle in the conduit or the wall breaker buzzes, the wire gauge or distribution transformer may be undersized.
- C. Determine if high voltage breakdown is present (hypot test)
 - 1. Inspect the high voltage cables, cable ends, and high voltage receptacles for carbon tracking and evidence of high voltage breakdown.
 - 2. Verify that the anode rotates and has no visible melt marks.
 - 3. With the high voltage cables removed from the A700 transformer, fill the receptacles with at least two inches of dielectric oil.
 - 4. Remove the jumper between TB3-19 and TB3-20 and insert a mAs meter set for AC and at maximum range. Connect a jumper between TP12 and TP15 (near kVp comp pots lower left of board).
 - 5. Remove U14 from the J400 board. Turn on the control and adjust P12 to achieve 1.8 VDC (+/- 0.1VDC) at TP15, referenced to TP1.
 - 6. Adjust tap switches for 170 VAC between kVp major and minor common. This will produce a no-load potential of 110 kVp on the transformer.
 - 7. Take exposures. Does the breaker trip or is mA present? If yes, the diode sticks in the transformer are shorting and must be replaced. If no, there may be a problem with the tube or cables not found with visual inspection.
- D. Identification of poor SCR gating or half-waving transformer
 - 1. Select 50 kVp, 100 mA, and 100 mAs. During exposure there should be about 108 VAC measured between P1 and P2 of the transformer. If the voltage at P1 and P2 is 108 VAC, replace the diode sticks. If the voltage is near 54 VAC, verify the SCR gate leads are tight, replace the SCR.
 - 2. A heavy groan from the transformer during exposure (that worsens as time and power are increased) is likely due to half-waving from the SCR.

3. Open diode sticks will appear as high current draw with one-half the expected mA (when the filament voltages are normal as described in II.A above, and the test equipment is used properly). Non-invasive kVp meters may indicate an 8 millisecond exposure regardless of selected time, due to their resetting during the missing pulse.

III. CIRCUIT BREAKER TRIPS WITH ROTATION OF KVP TAP SWITCH

- A. Look for a broken connection or poor solder joint on the noise suppression assembly between the common poles of the kVp major and minor tap switches. Replace if suspect, as this problem can damage rectifier sticks in the high voltage transformer.
- IV. "FIL 1" ERROR (low filament current in standby/idle mode)
 - A. Open connection at H8 or H9 of J400 board
 - 1. This could be a pin-to-plug or plug-to-board connection.
 - B. Open/loose terminal at either end of filament cable (XC, XL, XS)
 - C. Open F1 or F2 fuse in control
 - D. Open between resistor and resistor band (RSCC or RX)
 - E. Open or loose cathode cable, open filament of x-ray tube
 - 1. Remove the cathode cable from the tube end. Use a collimator bulb as a dummy filament load by connecting it between XC/XS or XC/XL pins and see if it will light use standby (idle) mode only. If the collimator bulb lights the tube filament may be open. If not, cathode cable may be open.
 - 2. Switch cables anode for cathode. If collimator bulb now lights replace the cable formerly in the cathode side.
 - F. Check the J400 board
 - Once it is confirmed there are no open primary connections or open secondary loads, check for about 1.7 VDC (+/- 0.2 V) between TP12 and TP1. If the voltage is missing, confirm LED4 is illuminated. If the light is OFF, jumper pin3 to pin4 of U18 optocoupler. Replace U18 optocoupler if function returns, or the J400 board.

- V. "FIL 2" ERROR (high filament current during standby/idle)
 - A. Confirm J400 operation New Install? Check if H.V. Cables crossed.
 - 1. Filament voltage during idle is typically about 24 VAC for small, 28 VAC for large spot as measured XC to XS or XL. This idle voltage typically produces 1.6 to 1.8 VDC at TP12. Turn pot P9 clockwise 1/2 turn. If error goes away, reduce the idling voltage by moving the bottom band on the boost resistor as needed.
 - 2. If FIL2 is still present, remove H8 header from the J400 board and cycle power. The error should change to FIL1, and voltage at TP12 should be near zero. If FIL2 error is still present replace the board.

VI. "FIL 3" ERROR (low filament current during prep)

- A. Confirm contact of boost resistor band to the resistor
- B. Confirm connection at H8 pins 7 and 8
 - 1. These contacts close to bypass most of the boost resistor. If contacts are good at H8 but closure does not occur at prep, replace A455 motor start board.
- VII. "FIL 4" ERROR (high filament current during prep)
 - A. Gross miscalibration
 - 1. Filament voltage is set too high.
 - a) 54 VAC at prep (XC to XS) is more than enough to achieve 300mA. A resistive cathode cable, filament cable or a miswire of XC/XL/XS (causing both filaments to be illuminated at one time) would require setting the filament voltage above this level and create the error.
 - b) 54 VAC at prep will produce about 2.68 VDC at TP12.
 - 2. Faulty or misused test equipment may give an improperly low reading for mA, causing an over-adjustment of the filament circuit and a "FIL4" error. Meters must be set to AC mA, and if equipped with an adjustable range scale (20, 200, 2000) *it must be set to maximum*.
 - 3. Functional test equipment does not see the mA produced.
 - a) Always remove the jumper between the two "AC mA" terminal locations when measuring mA, and always replace the jumper when finished.
 - b) Ensure spark gap on top of the HV transformer at "M1" is not shorted to ground.
 - 4. Failure of J400 board. This is by far the least likely possibility.

VIII. DIGITAL DISPLAY PROBLEMS

- A. Status display does not change or is at wrong state
 - 1. On J400 board set DIP2-4 to OFF (down), disabling "ready/cold" feature. If the display now updates, the "ready/cold" feature is not available with the display's software version. Update the software or leave switch OFF.
 - 2. If locked in "wait" check for shorted footswitch.
- B. Scrambled characters/missing segments
 - 1. Verify that each of the conductors on the cable going into H6 of J400 board is secured to the plug by looking beneath the plug cover.
 - 2. Confirm the Mate & Lock plug pins of the display cable at the rear of the control are seated and making contact.
 - 3. Confirm the plug and wires are well seated at J350 board in tubestand.
- C. Blank display
 - 1. Verify the Display cable connections at the H6 header of the J400 board, the pins in the Mate & Lock plug at the back of the control, and at the J350 display board.
 - 2. Confirm +20 VDC at TP3, +12VDC at TP5 of J400.
- D. kVp changes without tap switch adjustment
 - Select 100 mA and 143 VAC between the leads of MOV1 on the J400 board, located near H2. 70 kVp should be displayed. If the voltage varies check H2 and kVp tap switch connections.
 - 2. Confirm a stable 1 VAC between the junction of R35A/R36 and TP1. If this voltage is unstable replace the J400 board. If the voltage is stable and the kVp is still variable replace the J350 board.
- E. "HEAT" Error
 - 1. Verify TB4 "T5/T6" tube thermal switch connections are closed.
 - 2. Verify zero ohms on J400 board between U22 pin 1 and pin 2.
 - 3. Replace U22 or J400 board.
- F. "TUBE" Error
 - 1. Select 70 kVp and 10 mAs. The "Tube" Error should go away.
 - 2. If 70 kVp and 10 mAs cannot be selected, the J400 board may not be getting a valid kVp, time or mA input, or the board may be bad.
 - a) Adjust kV major/minor tap switches to achieve 180 VAC between the two leads of MOV1 on the J400 board. If no

voltage is present inspect H2 and kVp tap switch common poles for open connection.

- b) Remove the H10 plug from the J400 board. Place a jumper between pins 1 and 3 on the board. If the "Tube" error is cleared inspect the connections at the cable plug, replace time selector.
- c) Remove the H11 plug from the board, and jumper pin 1 to pin 2. If the error clears, inspect the switch and connections, replace the mA selector.
- d) If the error still exists replace the J400 board.

G. "ZCO" ERROR

1. Verify 14 VAC between TP1 and TP2 on J400 board. If voltage is present switch the wires in H3 pins 2 and 3. If error is still present, switch the wires back to their original locations and replace the J400 board.

IX.LONG "WAIT" DISPLAY, "ROTR" ERROR AT EXPOSURE ATTEMPT

- A. Verify rotor circuit
 - 1. The A455 motor start board controls rotor delay time. 240 VAC is applied between stator main (07) and common (09) for 1.5 seconds, and then the voltage drops to 50 VAC. The red LED on the motor start board turns on at that time. If the voltage stays at 240 VAC for longer than 1.5 seconds, adjust R17, replace C4, or replace the board.
 - If the red LED on the A455 board does not light, confirm 240/50 VAC at 07-09 stator connections of the control and at tube end. Verify stator resistance (30 ohms 07-black to 09-white, 60 ohms 08-red to 09white, 90 ohms 07-black to 08-red), replace A455 board if voltages and connections are correct.

X. EXPOSURE BEEP or PREP AT TURN ON

A. Check for shorted exposure switch.

- 1. If Latching Prep option is installed, verify open contacts between TB5-2 and TB5-4. If a short is measured, remove red and black wires from air switches at TB5-1 and TB5-4 and recheck measurement to determine if the short is in the air or strip switches.
- 2. If a two-position foot switch is installed, confirm open contacts between S1, S2 and S3 on TB3.
- 3. Remove H5 at J400 board. If problem continues replace J400.

XI.mAs DISPLAY IS INCORRECT FOR SETTINGS

- A. Verify mA switch inputs
 - 1. Confirm proper number of switch positions (two and three mA station units are available). If three positions are available on a two mA station unit, rotate switch fully CCW and remove the switch. With shaft pointing upwards, move the stop tab to the fully CCW position.
 - 2. Remove H11 plug on J400 to force in a good low mA station selection. Jumper pin 1 to pin 2 on J400 board H11 for medium mA station. Jumper pin 1 to pin 3 for high mA station (if third station exists). If mAs display responds to change of mA station replace mA switch, if not, replace J400 board.
- B. Verify time selector switch inputs
 - 1. Select 100 mA. Remove H10 plug from J400 board. Without H10 connected mAs display should be 0.83. Jumper pin 1 to pin 3 on the board, and display should change to 8.3. If it does, change time selector switch. If not, replace J400 board.

XII. NO mA MEASURED WITH RADIATION PRESENT

- A. Spark gap is shorted to top of transformer.
 - 1. The metal tab should not make contact with transformer lid. Normal spacing is equal to the thickness of radiographic film.
- B. Jumper is still present between TB3-19 and TB3-20.
- C. Faulty mA / mAs meter or leads.

XIII. UNIT WILL NOT GO INTO SLEEP MODE

- A. Ensure DIP2-1 on J400 is set to UP (ON) to enable feature.
- B. Poor connection between "W2" on J400 and mA overload coil of ON/OFF breaker.
- C. Shorted spark gap on top of transformer between M1 and GND.
 - 1. Proper spacing is the thickness of a piece of radiographic film.
- D. High line voltage fluctuation; new software may be less sensitive.

8.0 FINAL ADJUSTMENTS AND QUALITY ASSURANCE

8.1 LEVELING THE TUBE ARM

An adjustment screw for leveling the tube arm is provided on the front of tubestand, just below the tube arm base. Turn the Allen screw clockwise to raise the tube arm or counterclockwise to lower the tube arm.

8.2 40 INCH SID DETENTS

To position the tube 40" from either the table top or the film cassette, detents are provided. These detents are factory adjusted, however some adjustment is provided if needed. Access to the adjustments are from the rear of the tubestand.

8.3 VERIFICATION OF X-RAY TO LIGHT FIELD COINCIDENCE

The coincidence of x-ray to light field must be verified. Specific instructions can be found in the manual for the G800 collimator.

8.4 VERIFICATION OF X-RAY TO IMAGE RECEPTOR ALIGNMENT

Any left-to-right adjustment required to achieve coincidence between the central ray and the center of the image receptor can be achieved by adjusting the grid cabinet-to-tubestand interlock bracket, located on the lower front surface of the tubestand.

1. Any front-to-back adjustment of the central ray necessary (after the arm has been leveled from 8.1 above) can be achieved by the slotted holes in the tube mounting platform. Slightly loosen the collimator mounting bolts, shift the tube and collimator as required, and re-tighten the mounting bolts.

8.5 SECURING TABLE TO FLOOR

Once the table is in final position and all checks have been made, it can be secured t the floor using the four (4) clamps provided with the feet. Merely slide the clamp into the notch at the base of each foot, turn at 45° to the table and secure with a lag screw.

9.0 OPTIONS FOR FIELD INSTALLATION

Many options for the *InnoVet Select* TM will arrive at the jobsite pre-installed. However, all of the following features can be added after the installation is complete if the user so desires.

9.1 FILM BIN

The film bin comes in either a right hand or left hand configuration. See instructions that are provided with film bin kit.

9.2 TABLE TOP EXTENSION (For Fixed-Top tables only)

The table top extension adds 12" of additional table top surface and can be mounted on either end of the table. To attach the extension align notches in extension with buttons mounted on end of table support and slide over buttons as shown in Figure 19.



Figure 19: Table Top Extension – AWW430-08

9.3 ANIMAL RESTRAINTS

The animal restraints used are nautical style tie-down cleats. They are screwed to the underside of the table top frame at either end of the table. This prevents lead aprons from snagging on the cleats.

9.4 FOOT TREADLE

If the *InnoVet Select* TM has the latching prep Foot Treadle Option, the treadle must be installed prior to mounting the control into the table. If the Foot Treadle Option is being added to an existing system, the control will need to be removed from the table base and laid flat on the floor to install the treadle.

- 1. Remove the control from the table base by unclipping the truarc ring and extracting the hinge pin which hold the control to the air cylinder. Lay the control face down on cardboard and pull it slightly away from the table.
- 2. Put the treadle in position across the front of the table base.
- 3. Thread the spring guide pins located at either end of the treadle up through the table base until they protrude out the top of the base, securing each guide pin in place with cotter pin provided.
- 4. Bolt the foot treadle brackets to the underside of the table base using the hex head bolts provided.
- 5. Plug the foot treadle cable into the S4 location of the junction box at the right side of the table frame.
- 6. Slide the control so that the notches engage the pivot rollers, lift the control up into the table opening and resecure the control to the air cylinder using the pin and locking ring.
- 7. Remove the front panel of the control. Install the air switch packets using the screws and nuts provided so that the nozzle of each exits the rear of the control.
- 8. Connect the wires from the air switches to TB5-2 and TB5-3 in the bottom of the control. Connect the air hoses to the air switches by pushing them over the nozzles which exit the rear of the control.
- Place a jumper between TB5 1 and TB 5 2, and another jumper between TB5-3 and TB5-4. Verify on the J400 control board that DIP 1-5 is in the UP position. (To enable the Latching Prep mode as shown in Figure 18, page 7-12).
- 10. If any existing exposure switches are to be left in the system, alert the user that they will now operate like the Treadle; that is step on and release to start PREP, then step on and release to start EXPOSURE. Generally, removing these old two-position switches is recommended.

10.0 ELECTRICAL SCHEMATICS AND DIAGRAMS

Figure/Drawing No.

Description

Fig. 1	Floor Space for Fixed-Top Option
Fig. 2	Floor Space for Float-Top Options
Fig. 3	Tubestand Stop at Lower Left Rail
Fig. 4	Tubestand Positioning 5-4
Fig. 5	Mounting Tubestand on Rail 5-4
Fig. 6	Table Top Frame Screws & Lock Release Tabs
Fig. 7	Elevated Frame & Transverse Bearing Track 5-6
Fig. 9	Air Cylinder Attachment 5-10
Fig. 10	Mounting the Foot Loop 5-11
Fig. 11	Removing the Front Cover
Fig. 12	Interior of Control
Fig. 13	Rear of Control Showing Cable Entry Points
Fig. 14	Table Interface Junction Box 6-3
Fig. 15	TB3 and TB4 Connections
Fig. 16	mA Tracking Example7-4
Fig. 17	Control Board Component Identification
Fig. 18	DIP Switch Settings
Fig. 19	Table Top Extension 9-2
Fig. 20	System Cable Routing 10-2
Fig. 21	High Voltage Transformer Diagram 10-3
K158	Schematic, InnoVetSelect Table
J499	Schematic, K200 Control
J402	Schematic, Control Board (J400)
J352	Schematic, Display Board (J350)
01032	Schematic, RC Timer Board (01030-000)
	· · · · · ·

Summit Industries, Inc.

Innovet Select[™] Installation and Service Manual



Figure 20: System Cable Routing – AWW430-09



Figure 21: High Voltage Transformer Diagram – AWW430-10













COMPONENT LAYOUT for J350 DISPLAY PCBA



00130-000 RC Timing Control







Technical Support Dept Ph: (800) 729 - 9729 Fax: (773) 588 - 3404

A455 Interfacing to the Generator



***** PHASE SHIFT CAPACITOR: 25*uF*: Summit Part Number HAB18

W434 REV. A 30-Filament Boost 120VAC at PREP 50-Unregulated (Chassis F2 Via Switch) ĎC 12vdc 20 - 24v CR F1 (.5A) VR1 Fuse rating reduced to 0.04A not used CR9 CR10 K3 T1 C3 240 VAC > C6 QC(N.O.) F2 (10A) 09 K2 Anode Motor Common $15 \, \mathrm{c}$ ≷R16 $\leq R4$ 50 VAC > -12 O R11 CR11 7 Time Delay R6 Ŵ Adjust **₹**___] ₹ R17 R8 R9 IC2 5 87 70 Expose Interlock to Timer R5 CR32 F3 (20A) 110-IC1 R12 + R18 Expose ⋛ Ŵ *₹*R15 R7 R10≶ ______ C4 Ready 07 Motor Main CR2 LT1 *≩*_{R14} Indicator 140-C5 R20 $\leq R13$ CR6 K1 , + K2 R19 R2 CR5 \square C2 08 Motor Aux. Z1 CR1 Q3 130-Q1 + Q2 ____ ≶R3 TC1 - R1 Q4 0 VAC - 100-

Annotations by Roger F, Tech Support

Chassis F1.

SCHEMATIC, MOTOR START CIRCUIT



InnoVet Technique Chart 60Hz NOTE: 400 speed (except where noted) - automatic processing @ 94 degrees.

				processing e						
Head VD	Table To			Table w/Grid		-				
cm	4	6	8	10	12	14	16	18	20	22
kVp	56	60	64	68	72	76	80	84	88	92
mA	100	100	100	100	100	100	100	100	100	100
time	1/20	1/20	1/10	1/10	1/10	1/10	1/10	1/10	1/10	1/10
mAs	5	5	10	10	10	10	10	10	10	10
Head Lateral	Table To	p		Table w/Grid						
cm	4	6	8	10	12	14	16	18	20	22
kVp	56	60	64	68	72	76	80	84	88	92
mA	100	100	100	100	100	100	100	100	100	100
time	1/20	1/20	1/10	1/10	1/10	1/10	1/10	1/10	1/10	1/10
mAs	5	5	10	10	10	10	10	10	10	10
Thorax VD	Table To	p	Table w/Gri	d						
cm	4	7	10		16	19	22	25	28	31
kVp	52	58	64	70	76	82	88	94	100	106
mA	300	300	300	300	300	300	300	300	300	300
time	1/120	1/120	1/60	1/60	1/60	1/60	1/60	1/60	1/60	1/60
mAs	2.4	2.4	5	5	5	5	5	5	5	5
Thorax Lateral	Table To	p	Table w/Gri	d						
cm	4	7	10	13	16	19	22	25	28	31
kVp	54	60	66	72	78	84	90	96	102	108
mÅ	300	300	300	300	300	300	300	300	300	300
time	1/120	1/120	1/60	1/60	1/60	1/60	1/60	1/60	1/60	1/60
mAs	2.4	2.4	5	5	5	5	5	5	5	5
Thorax Decubitus	Table To	p	Table w/Gri	d						
cm	4	7	10	13	16	19	22	25	28	31
kVp	54	60	66	72	78	84	90	96	102	108
mA	300	300	300	300	300	300	300	300	300	300
time	1/120	1/120	1/60	1/60	1/60	1/60	1/60	1/60	1/60	1/60
mAs	2.4	2.4	5	5	5	5	5	5	5	5
Abdomen VD	Table To	p	Table w/Gri	d						
cm	4	7	10	13	16	19	22	25	28	31
kVp	52	58	62	68	74	80	86	92	98	104
mA	300	300	300	300	300	300	300	300	300	300
time	1/120	1/60	1/30	1/30	1/30	1/30	1/30	1/30	1/30	1/30
mAs	2.5	5	10	10	10	10	10	10	10	10
Abdomen Lateral	Table To	p	Table w/Gri	d						
cm	4	7	10		16	19	22	25	28	31
kVp	52	58	64	-	76	82	88	94	100	106
mA	300	300	300		300	300	300	300	300	300
time	1/120	1/60	1/30	1/30	1/30	1/30	1/30	1/30	1/30	1/30
mAs	2.5	5	10	10	10	10	10	10	10	10
			•							

InnoVet Technique Chart 60Hz - 400 speed (except where noted) - automatic processing @ 94 degrees.

NOTE: 400 speed (except wi	here noted) - automatic	processing	@ 94 degre	es. ¹			
Cervical VD	Table Top)		Table w/Grid			
cm	4	6	8	10	12	14	
kVp	56	60	64	68	72	76	
mÅ	300	300	300	300	300	300	
time	1/120	1/120	1/40	1/40	1/40	1/40	1/
							_

InnoVet Technique Chart 60Hz

)		Table w/Griu						
cm	4	6	8	10	12	14	16	18	20	22
kVp	56	60	64	68	72	76	80	84	88	92
mA	300	300	300	300	300	300	300	300	300	300
time	1/120	1/120	1/40	1/40	1/40	1/40	1/40	1/40	1/40	1/40
mAs	2.5	2.5	7.5	7.5	7.5	7.5	7.5	7.5	7.5	7.5
Cervical Lateral	Table Top)		Table w/Grid						
cm	4	6	8	10	12	14	16	18	20	22
kVp	56	60	64	68	72	76	80	84	88	92
mA	300	300	300	300	300	300	300	300	300	300
time	1/60	1/60	1/60	1/30	1/30	1/30	1/30	1/30	1/30	1/30
mAs	5	5	5	10	10	10	10	10	10	10
Cervical Soft Tissue Lateral	Table Top)		Table w/Grid						
cm	4	6	8	10	12	14	16	18	20	22
kVp	56	60	64	68	72	76	80	84	88	92
mA	300	300	300	300	300	300	300	300	300	300
time	1/120	1/120	1/120	1/60	1/60	1/60	1/60	1/60	1/60	1/60
mAs	2.5	2.5	2.5	5	5	5	5	5	5	5
Thoracic VD	Table Top)	Table w/Gri	id						
cm	4	7	10	-	16	19	22	25	28	31
kVp	52	58	64	-	76	82	88	94	100	106
mA	300	300	300		300	300	300	300	300	300
time	1/60	1/60	1/30	1/30	1/30	1/30	1/30	1/30	1/30	1/30
mAs	5	5	10	10	10	10	10	10	10	10
Thoracic Lateral	Table Top Table		Table w/Gri	id						
cm	4	7	10		16	19	22	25	28	31
kVp	52	58	64		76	82	88	94	100	106
mA	300	300	300		300	300	300	300	300	300
time	1/60	1/60	1/30	1/30	1/30	1/30	1/30	1/30	1/30	1/30
mAs										10
	5	5	10	10	10	10	10	10	10	10
Lumbar VD	5 Table Top		10 Table w/Gri	id	10	-		-	-	
cm	Table Top 4	7	Table w/Gri 10	id 13	16	19	22	25	28	31
cm kVp	Table Top 4 52	7 58	Table w/Gri 10 64	id 13 70	16 76	19 82	22 88	25 94	28 100	31 106
cm	Table Top 4	7	Table w/Gri 10	id 13 70 300	16	19	22	25	28	31 106 300
cm kVp	Table Top 4 52 300 1/120	7 58 300 1/120	Table w/Gri 10 64 300 1/40	id 13 70 300 1/40	16 76 300 1/40	19 82 300 1/40	22 88 300 1/40	25 94 300 1/40	28 100 300 1/40	31 106 300 1/40
cm kVp mA	Table Top 4 52 300	7 58 300	Table w/Gri 10 64 300	id 13 70 300	16 76 300	19 82 300	22 88 300	25 94 300	28 100 300	31 106 300
cm kVp mA time	Table Top 4 52 300 1/120	7 58 300 1/120 2.5	Table w/Gri 10 64 300 1/40	id 13 70 300 1/40 7.5	16 76 300 1/40	19 82 300 1/40	22 88 300 1/40	25 94 300 1/40	28 100 300 1/40	31 106 300 1/40 7.5
cm kVp mA time mAs Lumbar Lateral cm	Table Top 4 52 300 1/120 2.5 Table Top 4	7 58 300 1/120 2.5 7	Table w/Gri 10 64 300 1/40 7.5 Table w/Gri 10	id 13 70 300 1/40 7.5 id 13	16 76 300 1/40 7.5 16	19 82 300 1/40 7.5 19	22 88 300 1/40 7.5 22	25 94 300 1/40	28 100 300 1/40 7.5 28	31 106 300 1/40 7.5 31
cm kVp mA time mAs Lumbar Lateral cm kVp	Table Top 4 52 300 1/120 2.5 Table Top	7 58 300 1/120 2.5 7 58	Table w/Gri 10 64 300 1/40 7.5 Table w/Gri 10 64 300 1/40 7.5	id 13 70 300 1/40 7.5 id 13 70	16 76 300 1/40 7.5 16 76	19 82 300 1/40 7.5 19 82	22 88 300 1/40 7.5 22 88	25 94 300 1/40 7.5 25 94	28 100 300 1/40 7.5 28 100	31 106 300 1/40 7.5 31 106
cm kVp mA time mAs Lumbar Lateral cm	Table Top 4 52 300 1/120 2.5 Table Top 4 52 300 1/120 2.5 300 300	7 58 300 1/120 2.5 7	Table w/Gri 10 64 300 1/40 7.5 Table w/Gri 10 64 300 1/40 7.5	id 13 70 300 1/40 7.5 id 13 70 300	16 76 300 1/40 7.5 16 76 300	19 82 300 1/40 7.5 19	22 88 300 1/40 7.5 22 88 300	25 94 300 1/40 7.5 25	28 100 300 1/40 7.5 28 100 300	31 106 300 1/40 7.5 31 106 300
cm kVp mA time mAs Lumbar Lateral cm kVp	Table Top 4 52 300 1/120 2.5 Table Top 4 52	7 58 300 1/120 2.5 7 58	Table w/Gri 10 64 300 1/40 7.5 Table w/Gri 10 64 300 1/40 7.5	id 13 70 300 1/40 7.5 id 13 70	16 76 300 1/40 7.5 16 76	19 82 300 1/40 7.5 19 82	22 88 300 1/40 7.5 22 88	25 94 300 1/40 7.5 25 94	28 100 300 1/40 7.5 28 100	31 106 300 1/40 7.5 31 106

Shoulder CC	Table To		Table w/Gri							
cm	6	8	10	12		16	18	20	22	24
kVp	60	64	68	72		80	84	88	92	96
mA	100	100	100	100	100	100	100	100	100	100
time	1/20	1/20	1/10	1/10	1/10	1/10	1/10	1/10	1/10	1/10
mAs	5	5	10	10	10	10	10	10	10	10
Shoulder Lateral	Table To	2	Table w/Gri							
cm	6	8	10	12		16	18	20	22	24
kVp	60	64	68	72		80	84	88	92	96
mA	100	100	100	100	100	100	100	100	100	100
time	1/20	1/20	1/10	1/10	1/10	1/10	1/10	1/10	1/10	1/10
mAs	5	5	10	10	10	10	10	10	10	10
Pelvis VD	Table To	0	Table w/Gri							
cm	6	8	10	12		16	18	20	22	24
kVp	60	64	68	72		80	84	88	92	96
mA	300	300	300	300	300	300	300	300	300	300
time	1/60	1/60	1/30	1/30	1/30	1/30	1/30	1/30	1/30	1/30
mAs	5	5	10	10	10	10	10	10	10	10
Pelvis Lateral	Table To	ט	Table w/Gri							
cm	6	8	10	12	14	16	18	20	22	24
kVp	60	64	68	72		80	84	88	92	96
mA	300	300	300	300	300	300	300	300	300	300
time	1/40	1/40	1/20	1/20	1/20	1/20	1/20	1/20	1/20	1/20
mAs	7.5	7.5	15	15	15	15	15	15	15	15
Pelvis Distraction	Table To	ט	Table w/Gri							
cm	6	8	10	12		16	18	20	22	24
kVp	60	64	68	72		80	84	88	92	96
mA	300	300	300	300	300	300	300	300	300	300
time	1/60	1/60	1/30	1/30	1/30	1/30	1/30	1/30	1/30	1/30
mAs	5	5	10	10	10	10	10	10	10	10
Humerus CC - <u>100 speed</u>	Table To				Table w/Grid					
cm	2	4		8	10	12	14	16	18	20
kVp	54	58		66	70	74	78	82	86	90
mA	100	100		100	100	100	100	100	100	100
time	1/60	1/60	1/60	1/60	1/30	1/30	1/30	1/30	1/30	1/30
mAs	1.7	1.7	1.7	1.7	3.3	3.3	3.3	3.3	3.3	3.3
Humerus Lateral - <u>100 speed</u>	Table To				Table w/Grid					
cm	2	4		8	10	12	14	16	18	20
kVp	54	58		66	70	74	78	82	86	90
mA	100	100		100	100	100	100	100	100	100
time	1/60	1/60	1/60	1/60	1/30	1/30	1/30	1/30	1/30	1/30
mAs	1.7	1.7	1.7	1.7	3.3	3.3	3.3	3.3	3.3	3.3

Elbow CC - <u>100 speed</u>	Table Top)			Table w/Gr	id				
cm	2	4	6	8	10			-	-	20
kVp	52	56	60	64	70	74	78	82	86	90
mA	100	100	100	100	100	100	100	100	100	100
time	1/60	1/60	1/60	1/60	1/40	1/40	1/40	1/40	1/40	1/40
mAs	1.7	1.7	1.7	1.7	2.5	2.5	2.5	2.5	2.5	2.5
Elbow Lateral - <u>100 speed</u>	Table Top)			Table w/Gr	id				
cm	2	4	6	8	10			-		20
kVp	52	56	60	64	70	74	78	82	86	90
mA	100	100	100	100	100	100	100	100	100	100
time	1/60	1/60	1/60	1/60	1/40	1/40	1/40	1/40	1/40	1/40
mAs	1.7	1.7	1.7	1.7	2.5	2.5	2.5	2.5	2.5	2.5
Radius & Ulna CC - <u>100 speed</u>	Table Top)			Table w/Gr	id				
cm	2	4	6	8	10					20
kVp	54	58	62	66	70	74	78	82	86	90
mA	100	100	100	100	100	100	100	100	100	100
time	1/60	1/60	1/60	1/60	1/30	1/30	1/30	1/30	1/30	1/30
mAs	1.7	1.7	1.7	1.7	3.3	3.3	3.3	3.3	3.3	3.3
Radius & Ulna Lateral - <u>100 speed</u>	Table Top)			Table w/Gr					
cm	2	4	6	8	10			-	-	20
kVp	54	58	62	66	70					90
mA	100	100	100	100	100					100
time	1/60	1/60	1/60	1/60	1/30	1/30	1/30	1/30	1/30	1/30
mAs	1.7	1.7	1.7	1.7	3.3	3.3	3.3	3.3	3.3	3.3
Carpus & Metacarpus DP - <u>100 speed</u>	Table Top									
cm	1	2	3	4	5	6	7	8		10
kVp	60	62	64	66	68	70	72	74		78
mA	100	100	100	100	100	100	100	100		100
time	1/40	1/40	1/40	1/40	1/40	1/40	1/40	1/40	1/40	1/40
mAs	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5
Carpus & Metacarpus Lateral - 100 spe	<u>er</u> Table Top									
cm	1	2	3	4	5	6	7	8		10
kVp	60	62	64	66	68	70	72	74	76	78
mA	100	100	100	100	100	100	100	100	100	100
time	1/40	1/40	1/40	1/40	1/40	1/40	1/40	1/40	1/40	1/40
mAs	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5
Femur VD	Table Top		Table w/G							
cm	6	8		10 12						24
kVp	60	64		68 72						96
mA	100	100	10							100
time	1/30	1/30	1/15	1/15	1/15	1/15	1/15	1/15	1/15	1/15
mAs	3.3	3.3	6.7	6.7	6.7	6.7	6.7	6.7	6.7	6.7

$InnoVet\ Technique\ Chart\ 60Hz$ NOTE: 400 speed (except where noted) - automatic processing @ 94 degrees.

Femur Lateral	Table To	0	Table w/C	Grid						
cm	6	8	1	10 12	14	l 16	5 18	20	22	24
kVp	60	64	6	68 72		6 80) 84	. 88	92	96
mA	100	100	10	00 100	100) 100) 100	100	100	100
time	1/30	1/30	1/15	1/15	1/15	1/15	1/15	1/15	1/15	1/15
mAs	3.3	3.3	6.7	6.7	6.7	6.7	6.7	6.7	6.7	6.7
Stifle CC - <u>100 speed</u>	Table Top	0			Table w/Gr					
cm	2	4	6	8	10				18	20
kVp	52	56	60	64	68				84	88
mA	100	100	100	100	100				100	100
time	1/60	1/60	1/60	1/60	1/30	1/30	1/30	1/30	1/30	1/30
mAs	1.7	1.7	1.7	1.7	3.3	3.3	3.3	3.3	3.3	3.3
Stifle Lateral - <u>100 speed</u>	Table To				Table w/Gr					
cm	2	4	6	8	10				18	20
kVp	52	56	60	64	68				84	88
mA	100	100	100	100	100				100	100
time	1/60	1/60	1/60	1/60	1/30	1/30	1/30	1/30	1/30	1/30
mAs	1.7	1.7	1.7	1.7	3.3	3.3	3.3	3.3	3.3	3.3
Tibia & Fibula CC - <u>100 speed</u>	Table To	0			Table w/Gr					
cm	2	4	6	8	10				18	20
kVp	52	56	60	64	68				84	88
mA	100	100	100	100	100				100	100
time	1/60	1/60	1/60	1/60	1/30	1/30	1/30	1/30	1/30	1/30
mAs	1.7	1.7	1.7	1.7	3.3	3.3	3.3	3.3	3.3	3.3
Tibia & Fibula Lateral - <u>100 speed</u>	Table To	0			Table w/Gr					
cm	2	4	6	8	10				18	20
kVp	52	56	60	64	68		-		84	88
mA	100	100	100	100	100				100	100
time	1/60	1/60	1/60	1/60	1/30	1/30	1/30	1/30	1/30	1/30
mAs	1.7	1.7	1.7	1.7	3.3	3.3	3.3	3.3	3.3	3.3
Tarsus & Metatarsus DP - <u>100 speed</u>	Table Top									
cm	1	2	3	4	5	6	7	8		10
kVp	60	62	64	66	68	70	72	74		78
mA	100	100	100	100	100	100	100	100		100
time	1/40	1/40	1/40	1/40	1/40	1/40	1/40	1/40	1/40	1/40
mAs	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5
Tarsus & Metatarsus Lateral - <u>100 spee</u>	<u>ed</u> Table To									
cm	1	2	3	4	5	6	7	8		10
kVp	60	62	64	66	68	70	72	74		78
mA	100	100	100	100	100	100	100	100		100
time mAs	1/40 2.5	1/40 2.5	1/40 2.5	1/40 2.5	1/40 2.5	1/40 2.5	1/40 2.5	1/40 2.5	1/40 2.5	1/40 2.5

$InnoVet\ Technique\ Chart\ 60Hz$ NOTE: 400 speed (except where noted) - automatic processing @ 94 degrees.

Birds VD	Table Top)								
cm	1	2	3	4	5	6	7	8	9	10
kVp	50	51	52	52	53	53	54	54	55	56
mÅ	100	100	100	100	100	100	100	100	100	100
time	1/20	1/20	1/20	1/20	1/15	1/15	1/15	1/15	1/15	1/15
mAs	5	5	5	5	6.7	6.7	6.7	6.7	6.7	6.7
Birds Lateral	Table Top									
cm	1	2	3	4	5	6	7	8	9	10
kVp	50	51	52	52	53	53	54	54	55	56
mA	100	100	100	100	100	100	100	100	100	100
time	1/20	1/20	1/20	1/20	1/15	1/15	1/15	1/15	1/15	1/15
mAs	5	5	5	5	6.7	6.7	6.7	6.7	6.7	6.7
Pocket Pets DV	Table Top									
cm	1	2	3	4	5	6	7	8	9	10
kVp	53	54	55	56	57	58	59	60	61	63
mA	100	100	100	100	100	100	100	100	100	100
time	1/20	1/20	1/20	1/20	1/20	1/20	1/20	1/20	1/15	1/15
mAs	5	5	5	5	5	5	5	5	6.7	6.7
Pocket Pets Lateral	Table Top									
cm	1	2	3	4	5	6	7	8	9	10
kVp	53	54	55	56	57	58	59	60	61	63
mA	100	100	100	100	100	100	100	100	100	100
time	1/20	1/20	1/20	1/20	1/20	1/20	1/20	1/20	1/15	1/15
mAs	5	5	5	5	5	5	5	5	6.7	6.7
Reptiles DV	Table Top									
cm	1	2	3	4	5	6	7	8	9	10
kVp	53	54	55	56	57	58	59	60	61	63
mA	100	100	100	100	100	100	100	100	100	100
time	1/20	1/20	1/20	1/20	1/20	1/20	1/20	1/20	1/15	1/15
mAs	5	5	5	5	5	5	5	5	6.7	6.7
Reptiles Lateral	Table Top)								
cm	1	2	3	4	5	6	7	8	9	10
kVp	53	54	55	56	57	58	59	60	61	63
	100	100	100	100	100	100	100	100	100	100
mA	100	100	100	100	100	100				
mA time	1/20	1/20	1/20	1/20	1/20	1/20	1/20	1/20	1/15	1/15

InnoVet Technique Chart 60Hz NOTE: 400 speed (except where noted) - automatic processing @ 94 degrees.







No Exposure condition for InnoVets and 325 / 360 generators. A455 Motor Start Board Troubleshooting

The A455 has a quick visual indicator LED which lights when all is ok for exposure. This should light approx. 1.5 seconds after PREP begins. Does Anode energize? If yes, does LED light?

Included for you here are the collected documentation for the A455 Motor Start board, and we suggest you always consider these failure instances in the following logical sequence (slightly revised sequence from prior documents included in attachment).

- 1. Measure Mains voltage to generator and adjust Line Compensation if needed.
 - a. Confirmation measurement, Chassis F1 to F2 should measure 120VAC, when Generator ON.
 - b. Test DC volts measured on motor start board across C3 capacitor, typical 20 to 24volts (unregulated, only at boost/run).
- 2. Measure proper connection of anode motor windings with ohm Meter, generator OFF. You do not need to disconnect the cable unless ohms measured are incorrect.
 - a. COM.09 to MAIN.07 should be lowest measured value = 25 to 30 ohms.
 - b. COM.09 to AUX.08 should be approx 2x main = 50 to 60 ohms.
 - c. AUX.08 to MAIN.07 should be the sum, both windings = 75 to 90 ohms.
 - d. ALL windings must be Open to chassis ground, >1MEG ohm.

(Past cases a reverse wired tube had run for 5 or 6 years, and eventually current sensing reduced to below allowable limit. In some cases correcting the wiring recovers operation without parts.)

- 3. The anode motor Auxilliary Capacitor (also known as Start cap or Shift cap), location C1 on generator chassis, Summit part number HAB18, 25*u*F 450 VAC.
 - a. Measure AC Volts COM09 to MAIN07, note both boost (240v) and run (50v). Do not use Auto-Range option on voltmeter. Instead, pre-set meter range to 400V or higher.
 - b. C1: Rarely service personnel have a meter to accurately measure this as a capacitor, so we test this as AC Volts measured COM09 to AUX08, expected are 350v boost and 70v run.
 - c. In case of erroneous readings, double-check by measuring directly on C1 capacitor.
 - d. Consider measured values in relation to main motor winding (07 to 09) values. MAIN winding values need to be correct before expecting AUX winding values to be correct.
- 4. Define Failure mode after windings and capacitor have been confirmed.
 - a. Time delay for motor and filament boost:
 Generally this is confirmed with test 3a, but should be duration 1.5 seconds.
 Adjust R17 to fine tune. If R17 has no effect, check also R16, CR11, C4, IC2 (type LM358).
 - b. Current sensing protection to prevent exposure to a stationary anode.
 Confirm Reference to IC1 pin6 by measuring dc volts across R7, typical 1.6vdc.
 Confirm Main sensed current by measuring dc volts across R20, typical 5.4vdc.
 Confirm Aux. sensed current by measuring dc volts across R3, typical 5.9vdc
- 5. Bypass current sensing options:
 - c1. Short across R7 (setting reference to 0v).
- 6. SWAP or replace IC1 and IC2 (type LM325) and repeat 4a and 4b to test for chip failure.



NO Exposure & No LED light on A455 Motor Start Board

FIRST – Check your line voltage and re-tap auto transformer if needed. For the Spectra, this is simple operator task: hold the "Line Check" switch and turn the "Line Adjust" Tap switch so that the needle of the kV Meter is at "V" position.

Second, find which of the 2 general functions are failing, and then swap IC chips U1 and U2 to see if failure remains consistent, or if failure mode changes with change of chips then get new IC chips LM358N or equivalent.

Function 1 failure: Time Delay, Boost to Run. Function 2 failure: Current sensing of the Rotor windings. Failure 3: IC Chip failure.

1. Measure voltage to the Anode motor main winding Where?: AC volts 07 to 09, located at bottom of TB3) Normally 0V in idle, 240 VAC Boost (start of PREP) After 1.5 seconds 50 VAC Run (R17)

If correct, continue on to #2 current sensing. If NOT correct, try adjusting the potentiometer R17. Time delay is created by R16, R17 and C4, and charging "ramp" can be checked at anode of Diode CR11.

2. Check resistance Ω of Anode windings, to be sure nothing has changed. Main (07 to 09) measures 25 to 30 Ohms. Phase (08 to 09) measures 50 to 60 Ohms. SUM both (07 to 08) measures 75 to 90 Ohms. NOTE: ALL SHOULD BE OPEN TO GROUND.

2b. Use (-) negative side of C3 capacitor for common or ground reference. Measure the following voltages during Boost or Run mode (sometimes measure separate for both boost and run), because this board is not powered at other times.

Positive side of C3 = 20 to 24 vdc (unregulated DC supply).

check both sides of R7. measurements should be 0 vdc, and 1.6 vdc (approx), this is your reference.

Phase current sensing: Cathode of CR1 should be approx 5.9 vdc Main current sensing: Cathode of CR2 should be approx 5.4 vdc

If Main is OK but Phase is lower, then maybe you need to replace the C1 Rotor Capacitor.

Please call if you need any additional help.

Roger Frueh Technical Support Summit Industries Inc. Toll Free 1 (800) 729-9729 ext 4037



Troubleshooting the A455 Motor Start Board

FIRST – check your incoming line voltage and if needed re-tap your auto transformer.

Note that the K1 relay energizes at boost, putting 240 VAC between the "07" main and "09" common terminals on the tube stator cable. This voltage drops to 50 VAC during the run stage of PREP as K1 drops out and K2 is energized. Due to the stator capacitor, there should be about 360 VAC during boost and 70 VAC during run between "07" main and "08" phase terminals.

- 1. Throughout the PREP sequence (both boost and run) there must be 1.5 VDC across R7. Use the right side of R20 at top of board for ground reference. This sets the threshold voltage for minimum main and phase rotor current.
- 2. Throughout the PREP sequence there should be about 5 VDC on IC1 pin 3, and 8 VDC on IC1 pin 5. This is the sensed rotor current for main and phase windings respectively.
- 3. As a result of these voltages, IC1 pin 1 goes to 12 VDC, and IC1 pin 7 goes to 10 VDC.
- 4. Throughout the PREP sequence IC2 pins 2 and 6 each go to 6 VDC.
- 5. IC2 pin 5 charges to about 7 VDC during the PREP sequence. Once it charges to about 5 VDC, the output at IC2 pin 7 goes to 10 VDC. The time it takes to charge is controlled by R17 and C4, typically 1.5 seconds.
- 6. When IC2 pin 7 goes to 10 VDC, K1 drops out and K2 pulls in, the 240 VAC is removed and 50 VAC is applied to the stator as the board transitions from boost to run.
- 7. The 10 VDC output of IC2 pin 7 also changes the 3.4 VDC seen at IC2 pin 3 in boost to 7 VDC, making the output of IC1 pin 1 go to 10 VDC. This 10 VDC output turns on Q4, which energizes K4, turns on the LED, and closes the connection between J6 pins 7 and 11, passing 120 VAC hot to the timer to initiate the exposure.

Please call if you need any additional help.

Best Regards,

Technical SupportSummit Industries Inc.Toll Free1 (800) 729-9729Direct1 (773) 588-2444



Technical Support Dept Ph: (800) 729 - 9729 Fax: (773) 588 - 3404

A455 Interfacing to the Generator



***** PHASE SHIFT CAPACITOR: 25*uF*: Summit Part Number HAB18



A455 Troubleshooting

SCHEMATIC, MOTOR START CIRCUIT Page 5 of 6



PART REV Α 00425-000 SCR WIRING INSTRUCTIONS The gate terminal layout of SCR assemblies manufactured by SANREX, are slightly different from SCR assemblies manufactured by GE, WESTINGHOUSE, AEG, EUPEC, etc. Therefore, when replacing an SCR assembly with one of a different gate terminal layout, care must be taken to properly connect the corresponding gate leads. Incorrect wiring of the SCR gate terminals will not damage the SCR but will result in failure of the SCR to turn ON. It is strongly suggested that the serviceman make a wiring diagram of the original SCR, before any wires are removed. Then, using the terminal locations table indicated below, reconnect the new SCR to the proper terminals. EUPEC SCR SANREX SCR ALL LEADS CONNECTED TO: TO BE CONNECTED TO: × AK ์ 🗙 ์ AK TERMINAL "1" TERMINAL "3" TERMINAL "2" TERMINAL "2" × ĸ × Κ TERMINAL "1" TERMINAL "3" З TERMINAL "4" TERMINAL "K2" × Α \times А TERMINAL "2" TERMINAL "G2" 4 S L TERMINAL "3" TERMINAL "K1" TERMINAL "1" TERMINAL "G1" EUPEC SCR SANREX SCR A316 00194-000 TOLERANCE UNLESS OTHERWISE SPECIFIED X.X ±.030 X.XX ±.015 SUMMIT INDUSTRIES, IN X.XXX ±.005 $1/X \pm 1/64$ | 1/X ±1/64 2901 W. L ANGLES ±1/2* | SUMMIT TEL 773-509-6273 2901 W. LAWRENCE AVE, CHICAGO, IL FAX 773-588-6820 INSTRUCTIONS, SCR WIRING INSTRUCTION NEXT ASS'Y 00424-000 SEE ABOVE SCALE FINISH NONE DATE SIZE A SHEET 1 NO. 00425-000 A 2289 RELEASED JP 07/16/01 DWG BY DATE APPROVED REVA 06/25/01 LTR | ECR | REVISION BY DATE TΡ

ORIGINAL INNOVET, GENERATOR MODEL W300

SCR WIRING INSTRUCTIONS

The gate terminal layout of SCR assemblies manufactured by SANREX, are slightly different from SCR assemblies manufactured by GE, WESTINGHOUSE, AEG, EUPEC, etc. Therefore, when replacing an SCR assembly with one of a different gate terminal layout, care must be taken to properly connect the corresponding gate leads. Incorrect wiring of the SCR gate terminals will <u>not</u> damage the SCR but <u>will</u> result in failure of the SCR to turn ON. Each of the INNOVET x-ray controls, model numbers W300 & W400, contain two SCR assemblies, one as the MAIN SCR & the other as the BACK-UP SCR. The following wiring diagrams can be used to properly connect the SCR gate leads, if required:

SCR's manufactured by GE, WESTINGHOUSE, AEG, EUPEC, etc.





SCRs manufactured by SANREX





BACK-UP SCR ASS'Y

CCR PKG 00194-000

SCR (SILICON CONTROLLED RECTIFIER) BLOCK (ALSO KNOWN "PHASE CONTROL THYRISTOR MODULE") DUAL SCR PACKAGE RMS ON – STATE CURRENT – 86 AMPS © 81 DEGREES C. AVERAGE ON–STATE CURRENT–55 AMPS © 81 DEGREES C. SURGE ON–STATE CURRENT (1/2 CYCLE) – 1190 AMPS © 50 HZ – 1300 AMPS © 60 HZ REPTITIVE REVERSE VOLTAGE – 800 V NON–REPETITIVE REVERSE VOLTAGE – 960 V GATE TRIGGER CURRENT (MAXIMUM REQUIRED TO ASSURE CONDUCTION) – 50 mA GATE TRIGGER VOLTAGE (MAXIMUM REQUIRED TO ASSURE CONDUCTION) – 3 V

MFR: SANREX CORP. (PORT WASHINGTON, NY) MFR P/N: PK55FG80

NOTE: (SIMILAR TO A316)





^{PART}N□. 00194

B



THYRISTOR MODULE **PK(PD,PE)55FG**

Power Thyristor/Diode Module **PK55FG** series are designed for various rectifier circuits and power controls. For your circuit application, following internal connections and wide voltage ratings up to 1600V are available. and electrically isolated mounting base make your mechanical design easy.

- IT(AV) 55A, IT(RMS) 86A, ITSM 1300A
- ●di/dt 100A/µs
- dv/dt 1000V/µs

(Applications) Various rectifiers AC/DC motor drives Heater controls Light dimmers Static switches Internal Configurations



Maximum Ratings

		Ratings								
Symbol	Item	PK55FG40 PD55FG40 PE55FG40	PD55FG40 PD55FG80 PD55FG120 PD55FG		PK55FG160 PD55FG160 PE55FG160	Unit				
Vrrm	* Repetitive Peak Reverse Voltage	400	800	1200	1600	V				
Vrsm	* Non-Repetitive Peak Reverse Voltage	480	960	1300	1700	V				
Vdrm	* Repetitive Peak off-state Voltage	400	800	1200	1600	V				

Symbol	lt	em	Conditions	Ratings	Unit
T(AV)	*Average On-s	state Current	Single phase, half wave, 180° conduction, Tc=81°C	55	Α
T(RMS)	*R.M.S. On-st	ate Current	86	А	
Ітѕм	*Surge On-state Current		1/2 Cycle, 50/60Hz, Peak Value, non-repetitive	1190/1300	А
l²t	* I ² t		Value for one cycle surge current	7040	A ² S
Рсм	Peak Gate Power Dissipation			10	W
PG(AV)	Average Gate P	ower Dissipation		1	W
FGM	Peak Gate Current			3	Α
Vfgm	Peak Gate Voltage (Forward)			10	V
Vrgm	Peak Gate Vo	oltage (Reverse)		5	V
di/dt	Critical Rate of Rise	e of On-state Current	100	A/µs	
Viso	* Isolation Breakdo	wn Voltage (R.M.S)	A.C. 1minute	2500	V
Tj	* Operating Junc	tion Temperature		-40 to +125	°C
Tstg	* Storage Tem	perature		-40 to +125	°C
	Mounting	Mounting (M5)	Recommended Value 1.5-2.5 (15-25)	2.7 (28)	N∙m
	Torque	Terminal (M5)	Recommended Value 1.5-2.5 (15-25)	2.7 (28)	(kgf•cm)
	Mass		Typical Value	170	g

Electrical Characteristics

Item	Conditions	Ratings	Unit
Repetitive Peak off-state Current,max	Tj=125℃, Vd=Vdrm	15	mA
* Repetitive Peak Reverse Current,max	Tj=125℃, Vd=Vdrm	15	mA
*On-state Voltage,max	I⊤=165A	1.6	V
Gate Trigger Current, max	VD=6V, IT=1A	50	mA
Gate Trigger Voltage,max	VD=6V, IT=1A	3	V
Gate Trigger Voltage, min	Тј=125℃, VD=½VDRM	0.25	V
Critical Rate of Rise of off-state Voltage, min	Tj=125℃, VD=⅔VDRM	1000	V/µs
*Thermal Impedance,max	Junction to case	0.5	°C/W
	Repetitive Peak off-state Current,max * Repetitive Peak Reverse Current,max * On-state Voltage,max Gate Trigger Current,max Gate Trigger Voltage,max Gate Trigger Voltage,min Critical Rate of Rise of off-state Voltage,min	Repetitive Peak off-state Current,maxTj=125°C, VD=VDRM* Repetitive Peak Reverse Current,maxTj=125°C, VD=VDRM* On-state Voltage,maxIT=165AGate Trigger Current,maxVD=6V, IT=1AGate Trigger Voltage,maxVD=6V, IT=1AGate Trigger Voltage,maxTj=125°C, VD= $\frac{1}{2}$ VDRMCritical Rate of Rise of off-state Voltage,minTj=125°C, VD= $\frac{2}{3}$ VDRM	Repetitive Peak off-state Current,maxTj=125°C, VD=VDRM15* Repetitive Peak Reverse Current,maxTj=125°C, VD=VDRM15* On-state Voltage,maxIT=165A1.6Gate Trigger Current,maxVD=6V, IT=1A50Gate Trigger Voltage,maxVD=6V, IT=1A3Gate Trigger Voltage,minTj=125°C, VD=1/2VDRM0.25Critical Rate of Rise of off-state Voltage,minTj=125°C, VD=2/3VDRM1000

*mark: Thyristor and Diode part. No mark: Thyristor part

SanRex 50 Seaview Blvd. Port Washington, NY 11050-4618 PH.(516)625-1313 FAX(516)625-8845 E-mail: semi@sanrex.com

PK(PD,PE)55FG





