

ABI PRISM[®] 7700 Sequence Detection System and TaqMan[®] Card Upgrade

Installation Manual

© Copyright 2001. All rights reserved.

For Research Use Only. Not for use in diagnostic procedures.

For limited label license information, please refer to the ABI PRISM 7700 Sequence Detection System User's Manual, P/N 904989

ABI PRISM and its design, Applied Biosystems, Masterpiece, and MicroAmp are registered trademarks of Applied Biosystems Corporation or its subsidiaries in the U.S. and certain other countries.

ABI is a trademark of Applied Biosystems Corporation or its subsidiaries in the U.S. and certain other countries.

AmpliAmp Gold and TaqMan are registered trademarks of Applied Biosystems Corporation, Inc.

AppleScript and Macintosh are registered trademarks of Apple Computer, Inc.

All other trademarks are the sole property of their respective owners.

Applied Biosystems Corporation is committed to providing the world's leading technology and information for life scientists. Applied Biosystems Corporation consists of the Applied Biosystems and Applied Genomics businesses.

07/2001

4316192C

Contents

1 Start Here

Overview	1-1
Before You Begin	1-2
Installation Schedules/Checklists for Different Installations	1-3

2 Safety

Overview	2-1
Laboratory and Instrument Safety Requirements	2-2
Laser Safety	2-5
Laboratory Requirements	2-6
General Warnings	2-8

3 Setting Up and Testing the System

Overview	3-1
Getting Started	3-2
Electrical Configuration	3-3
Computer and Software Setup	3-6
System Testing At Startup	3-8
Thermal Cycler Tests	3-9
Laser Alignment Test	3-10
Fluorescence and Background Checks	3-11

4 Performing a Spectral Calibration

Overview	4-1
Spectral Calibration Overview	4-2
Materials	4-2
Archiving Current Component Files	4-3
Generating the Background Component File	4-4
Generating a Pure Dye File	4-6

5 Verifying the Instrument

Overview	5-1
The ABI PRISM 7700 Instrument Verification Run	5-2
Preparing the β -Actin Installation Plate	5-3
Setting Up and Starting the β -Actin Installation Run	5-6

Preparing and Running the RNase P Instrument Verification Plate	5-9
Analyzing Data	5-13
7700 Sequence Detection System Verification Calculations.....	5-15

6 *TaqMan Card Upgrade Protocol*

Overview	6-1
<i>Section: TaqMan Card Upgrade Overview</i>	6-2
Installation Procedure	6-3
Materials	6-4
<i>Section: Performing a TaqMan Card Upgrade</i>	6-5
Installing the Firmware and Software	6-6
Alignment Pin Replacement	6-9
Setting Up the Filling Station and Vacuum Pump Assembly	6-10
Preparing the PCR Reaction Mix and Filling the TaqMan Card.....	6-12
Loading the Card into the 7700 Instrument.....	6-19
Determining Exposure Time for Data Collection	6-22
Setting Up and Starting a TaqMan Card Run	6-25
<i>Section: Analyzing Card Data, Exporting Results, and Troubleshooting</i>	6-27
Analyzing the Card Data and Verifying System Performance	6-28
Exporting Results	6-31
TaqMan Card Verification Calculations	6-32
Troubleshooting.....	6-34

7 *Customer Training*

Overview	7-1
Customer Training for the 7700 Sequence Detection System.....	7-2
Customer Training for the TaqMan Card Upgrade	7-3
Completing the Installation	7-4

A *Abbreviations*

B *Material and Equipment*

ABI PRISM 7700 Instrument Installation Kit.....	B-2
TaqMan Card Upgrade Kit	B-4
Customer-Supplied Materials and Equipment.....	B-5

C *7700 Installation Specifications*

7700 Instrument Installation Specifications	C-2
TaqMan Card Installation Specifications.....	C-3

D Installation Report Checklists

Preinstallation ChecklistsD-2
Customer Training ChecklistsD-5
Laboratory Safety Checklist.....D-7
Installation Reports.....D-9

E Creating a TaqMan Card Template

F Software Upgrade Install Procedure

Index

Start Here

1

Overview

About This Chapter This chapter describes the purpose of this manual and installation schedules and checklists for different installations.

In This Chapter This chapter contains the following topics:

Topic	See Page
Before You Begin	1-2
Installation Schedules/Checklists for Different Installations	1-3

Before You Begin

Purpose of This Manual

The instructions presented in this manual are intended to assist a trained Applied Biosystems service engineer in the installation and performance verification of the:

- ◆ ABI PRISM® 7700 Sequence Detection System
- ◆ TaqMan card upgrade for the ABI PRISM 7700 instrument
- ◆ ABI PRISM 7700 instrument and TaqMan card upgrade
- ◆ Upgrade of software and firmware on the ABI PRISM 7700 instrument

This manual may serve as a reference during installation.

For a thorough understanding of the instrument and the installation procedure, it is essential to receive proper service training and refer to the *ABI PRISM 7700 SDS User's Manual*. If more detailed information pertaining to the service of the instrument is needed, refer to the *ABI PRISM 7700 Sequence Detection System Service Manual*.

Installation Procedures

Use the following table as a guide for chapters to reference for specific types of installations:

If you are installing the...	Refer to...	Install Time
ABI PRISM 7700 Sequence Detection System	Chapters 1–5, 7.	2 days
TaqMan card only	Chapters 1, 2, 4- 7.	1 day
ABI PRISM 7700 Sequence Detection System software and the TaqMan card upgrade	Chapters 1–7.	2.5 days
SDS software upgrade	Appendix F, "Software Upgrade Install Procedure."	10 minutes

Preinstallation Checklists

Before scheduling the installation, the customer should be contacted in order to review the preinstallation checklist (see "Preinstallation Checklists" on page D-2).

Review the completed checklist before starting the installation in order to ensure that everything needed to complete the installation is available. In some instances, the installation can be started even though some preparations have not been made; this decision should be made based on the experience of the service engineer and after consulting the local service management.

Installation Schedules/Checklists for Different Installations

Suggested Installation Schedules The following suggested installation schedules are to be used as guidelines. Experienced engineers may modify the schedules to complete installations more quickly.

ABI PRISM 7700 Installation The ABI PRISM 7700 SDS installation procedure will take approximately 2 days.
Check Off Each Item As It Is Completed

	Verify that all equipment, chemicals, and supplies have arrived. Check "ABI Prism 7700 Instrument Installation Kit" on page B-2 and "TaqMan Card Upgrade Kit" on page B-4.
	Verify that all customer-supplied equipment and supplies needed to run the chemical installation kit are readily available. Refer to "Customer-Supplied Materials and Equipment" on page B-5 and "Preinstallation Checklists" on page D-2.
	Unpack and Setup the Power Macintosh®. This should include reformatting the hard drive, reloading the Macintosh computer operating system, and installing the Sequence Detection System software.
	Unpack the ABI PRISM 7700 instrument.
	Check the internal electrical connections on the ABI PRISM 7700 instrument.
	Test the ABI PRISM 7700 instrument software/hardware functions.
	Create a new Background file and a new Spectral Components file.
	Prepare and run an installation tray using the chemical installation kit and accessories from the packing kit.
	Demonstrate to the customer the instrument operation using the Sequence Detection System software, including Sample Setup, Thermal Cycler Conditions, and Starting Runs.
	Analyze the run from the previous day and determine if it passes the installation specifications. The customer should do this as part of their training.
	Show the customer how to display Raw Data, Standard Curves, and Amplification Plots.
	Review data analysis procedures with the customer, including setting the threshold, reanalyzing data, and examining the Experimental Report for statistical information.
	Demonstrate system maintenance procedures.
	Demonstrate Macintosh computer maintenance and file backup.
	Fill out and send in all software registration cards and the installation postcard that are found with the system. These include the cards for the Macintosh computer software.
	Read the MSDSs provided by the chemical manufacturers.

**TaqMan Card
Upgrade Installation**

The TaqMan® Human Cytokine Card Upgrade installation will take approximately one day.

Check Off Each Item As It Is Completed

	<p>Verify that all equipment, chemicals, and supplies have arrived.</p> <p>Check “ABI Prism 7700 Instrument Installation Kit” on page B-2 and “TaqMan Card Upgrade Kit” on page B-4.</p>
	<p>Verify that all customer-supplied equipment and supplies needed to run the chemical installation kit are readily available.</p> <p>Refer to “Customer-Supplied Materials and Equipment” on page B-5 and “Preinstallation Checklists” on page D-2.</p>
	<p>Create a new Background file and a new Spectral Components file.</p>
	<p>Run an RNase P installation tray to verify system performance.</p>
	<p>Set up the ABI PRISM® Card Filling station and vacuum pump assembly and verify that the vacuum pump pulls the necessary vacuum for filling cards.</p>
	<p>Analyze the RNase P run and determine if it passes the installation specifications.</p>
	<p>Install the card upgrade firmware and software.</p>
	<p>Replace 7700 alignment pins if necessary.</p>
	<p>Fill and run a card using human genomic DNA and Universal PCR master mix. Before starting chemistry run, check 7700 exposure time for card data collection and adjust if necessary.</p> <p>Demonstrate to customer how to fill card and set up card run and how to determine the proper exposure time for a run.</p>
	<p>Analyze the run and determine if VIC dye layer passes the installation specifications.</p> <p>The customer should do this as part of their training.</p>
	<p>Review card troubleshooting procedures: checking for/fixing vacuum leaks in filling system, aligning card and adapters in the 7700 instrument, etc.</p>
	<p>Fill out and send in the installation post card that is included in the installation kit.</p>

Safety

2

Overview

About This Chapter This chapter contains general safety information and guidelines for working in a customer's laboratory.

In This Chapter This chapter contains the following topics:

Topic	See Page
Laboratory and Instrument Safety Requirements	2-2
Laser Safety	2-5
Laboratory Requirements	2-6
General Warnings	2-8

Laboratory and Instrument Safety Requirements

Documentation User Attention Words Five user attention words appear in the text of all Applied Biosystems user documentation. Each word implies a particular level of observation or action as described below.

Note Calls attention to useful information.

IMPORTANT Indicates information that is necessary for proper instrument operation.

⚠ CAUTION Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices.

⚠ WARNING Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

⚠ DANGER Indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury. This signal word is to be limited to the most extreme situations.

Chemical Hazard Warning **⚠ WARNING CHEMICAL HAZARD.** Some of the chemicals used with Applied Biosystems instruments and protocols are potentially hazardous and can cause injury, illness, or death.

- ◆ Read and understand the material safety data sheets (MSDSs) provided by the chemical manufacturer before you store, handle, or work with any chemicals or hazardous materials.
 - ◆ Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (*e.g.*, safety glasses, gloves, or protective clothing). For additional safety guidelines, consult the MSDS.
 - ◆ Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with adequate ventilation (*e.g.*, fume hood). For additional safety guidelines, consult the MSDS.
 - ◆ Check regularly for chemical leaks or spills. If a leak or spill occurs, follow the manufacturer's cleanup procedures as recommended on the MSDS.
 - ◆ Comply with all local, state/provincial, or national laws and regulations related to chemical storage, handling, and disposal.
-

Chemical Waste Hazard Warning **⚠ WARNING CHEMICAL WASTE HAZARD.** Wastes produced by Applied Biosystems instruments are potentially hazardous and can cause injury, illness, or death.

- ◆ Read and understand the material safety data sheets (MSDSs) provided by the manufacturers of the chemicals in the waste container before you store, handle, or dispose of chemical waste.
- ◆ Handle chemical wastes in a fume hood.
- ◆ Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (*e.g.*, safety glasses, gloves, or protective clothing). For additional safety guidelines, consult the MSDS.
- ◆ Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with adequate ventilation (*e.g.*, fume hood). For additional safety guidelines, consult the MSDS.
- ◆ After emptying the waste container, seal it with the cap provided.

- ◆ Dispose of the contents of the waste tray and waste bottle in accordance with good laboratory practices and local, state/provincial, or national environmental and health regulations.

Site Preparation and Safety Guide

A site preparation and safety guide is a separate document sent to all customers who have purchased an Applied Biosystems instrument. Refer to the guide written for your instrument for information on site preparation, instrument safety, chemical safety, and waste profiles.

Ordering MSDSs

You can order free additional copies of MSDSs for chemicals manufactured or distributed by Applied Biosystems using the contact information below.

To order MSDSs...	Then...							
Over the Internet	a. Go to our Web site at www.appliedbiosystems.com/techsupp b. Click MSDSs <table border="1" data-bbox="829 758 1463 1024"> <thead> <tr> <th>If you have...</th> <th>Then...</th> </tr> </thead> <tbody> <tr> <td>The MSDS document number or the Document on Demand index number</td> <td>Enter one of these numbers in the appropriate field on this page.</td> </tr> <tr> <td>The product part number</td> <td rowspan="2">Select Click Here, then enter the part number or keyword(s) in the field on this page.</td> </tr> <tr> <td>Keyword(s)</td> </tr> </tbody> </table> c. You can open and download a PDF (using Adobe® Acrobat® Reader™) of the document by selecting it, or you can choose to have the document sent to you by fax or email.	If you have...	Then...	The MSDS document number or the Document on Demand index number	Enter one of these numbers in the appropriate field on this page.	The product part number	Select Click Here , then enter the part number or keyword(s) in the field on this page.	Keyword(s)
If you have...	Then...							
The MSDS document number or the Document on Demand index number	Enter one of these numbers in the appropriate field on this page.							
The product part number	Select Click Here , then enter the part number or keyword(s) in the field on this page.							
Keyword(s)								
By automated telephone service	Use "To Obtain Documents on Demand" under "Technical Support."							
By telephone in the United States	Dial 1-800-327-3002 , then press 1 .							
By telephone from Canada	<table border="1" data-bbox="829 1318 1463 1444"> <thead> <tr> <th>To order in...</th> <th>Dial 1-800-668-6913 and...</th> </tr> </thead> <tbody> <tr> <td>English</td> <td>Press 1, then 2, then 1 again</td> </tr> <tr> <td>French</td> <td>Press 2, then 2, then 1</td> </tr> </tbody> </table>	To order in...	Dial 1-800-668-6913 and...	English	Press 1 , then 2 , then 1 again	French	Press 2 , then 2 , then 1	
To order in...	Dial 1-800-668-6913 and...							
English	Press 1 , then 2 , then 1 again							
French	Press 2 , then 2 , then 1							
By telephone from any other country	See the specific region under "To Contact Technical Support by Telephone or Fax" under "Technical Support."							

For chemicals not manufactured or distributed by Applied Biosystems, call the chemical manufacturer.

Instrument Safety Labels

Safety labels are located on the instrument. Each safety label has three parts:

- ◆ A signal word panel, which implies a particular level of observation or action (*e.g.*, CAUTION or WARNING). If a safety label encompasses multiple hazards, the signal word corresponding to the greatest hazard is used.
- ◆ A message panel, which explains the hazard and any user action required.

- ◆ A safety alert symbol, which indicates a potential personal safety hazard. See the *ABI PRISM® 7700 Site Preparation and Safety Guide* for an explanation of all the safety alert symbols provided in several languages.

About Waste Disposal

As the generator of potentially hazardous waste, it is your responsibility to perform the actions listed below.

- ◆ Characterize (by analysis if necessary) the waste generated by the particular applications, reagents, and substrates used in your laboratory.
- ◆ Ensure the health and safety of all personnel in your laboratory.
- ◆ Ensure that the instrument waste is stored, transferred, transported, and disposed of according to all local, state/provincial, or national regulations.

Note Radioactive or biohazardous materials may require special handling, and disposal limitations may apply.

Moving and Lifting the Instrument

⚠ CAUTION PHYSICAL INJURY HAZARD. Improper lifting can cause painful and sometimes permanent back injury.

Use proper lifting techniques when lifting or moving the instrument. Safety training for proper lifting techniques is recommended.

Do not attempt to lift or move the instrument without the assistance of others. Depending on the weight of the instrument, this action may require two or more people.

Before Operating the Instrument

Ensure that everyone involved with the operation of the instrument has:

- ◆ Received instruction in general safety practices for laboratories
- ◆ Received instruction in specific safety practices for the instrument
- ◆ Read and understood all related MSDSs

⚠ CAUTION Avoid using this instrument in a manner not specified by Applied Biosystems. Although the instrument has been designed to protect the user, this protection can be impaired if the instrument is used improperly.

Laser Safety

Contains Argon Gas Laser

In addition to the standard hazards in electromechanical instrumentation, the ABI PRISM® 7700 Sequence Detector contains an argon ion gas laser that emits up to 10 mW of electromagnetic radiation, principally at 488.0 nm. The instrument has been designed to comply with Title 21, U.S. Government DHEW/BRH Performance Standards, Chapter 1, Subchapter J, Section 1040, as applicable. This product falls into Class I category when all interlocks are in place.

Safety Features Incorporated

The following safety features have been incorporated in the ABI PRISM 7700 instrument:

- ◆ The cabinet is designed to prevent access to collateral laser radiation exceeding the accessible emission limits listed in Performance Standards for Laser Products, 21 CFR 1040.10.
 - ◆ The heated cover assembly, when it is not forward and locked down, activates interlock switches that cause the laser light path to be blocked.
 - ◆ Safety labels for Class I Standards have been affixed to the instrument to warn of the danger from laser radiation.
 - ◆ These interlocks and labels must not be removed or defeated by customers. The interlocks and labels ensure compliance with the above mentioned performance standards of the U.S. Code of Federal Regulations. Only authorized service personnel are allowed to defeat these interlocks during service, maintenance, or installation of the instrument.
-
-

WARNING ! WARNING ! VISION HAZARD: LASER SOURCE. When the laser beam is exposed, as it will be during some of the installation procedures, it is Class IIIb, which is extremely hazardous. Direct or reflected laser beam exposure at 10mW for 0.1 second can burn the eye's retina, leaving a permanent blind spot.

Safety Precautions When Interlock Is Defeated

When performing procedures that require the interlock to be defeated, you must follow these safety precautions:

- ◆ Post warning signs (P/N T-5350) on the laboratory door and in the work area to warn of the danger from laser radiation.
 - ◆ Make sure the laser power is at its minimum setting to decrease the hazard from the light intensity.
 - ◆ Never look directly into the laser beam or the laser head opening.
 - ◆ Always use the laser alignment tool when checking the laser light path during service procedures to prevent the reflection of laser light into the room.
 - ◆ Remove all rings, watches, jewelry, and metal frame eyeglasses that might reflect the laser beam into someone's eyes.
 - ◆ Wear laser safety goggles (P/N 100355) specially designed to protect eyes against laser light.
-
-

Laboratory Requirements

Preinstallation Checklist

Before scheduling the installation, the customer should be contacted in order to go through the preinstallation checklist (see page D-2). Review the completed checklist before starting the installation in order to ensure that everything needed to complete the installation is available. In some instances, the installation can be started even though some preparations have not been made; this decision should be made based on the experience of the service engineer and after consulting the local service management.

Laboratory Space Requirements

ABI PRISM 7700 Instrument Dimensions

Width	Depth	Height	Weight
94 cm (37 in.)	72.5 cm (28.5 in.)	61 cm (24 in.)	130 kg (286 lb.)

IMPORTANT A minimum of 15 cm (6.0 in.) clearance at the rear of the instrument and 69 cm (27 in) above the instrument is necessary to provide adequate ventilation and service access.

The Macintosh® computer should be on the same table as the ABI PRISM 7700 instrument. If this is impractical, the equipment should not be located further than 1.8 m (6 ft.) from the instrument. Both the ABI PRISM 7700 instrument and the Macintosh computer fit on a 1.8 meter (6 ft.) wide table or bench.

Environmental Conditions and Requirements

The laboratory temperature should be maintained between 15 °C (59 °F) and 30 °C (85 °F). The relative humidity in the laboratory should be no greater than 80%. The unit should not be placed near heaters or cooling ducts.

The thermal output of the instrument when it is operating under normal conditions is 8800 Btu/h (~2600 W).

This unit is for indoor use only and for altitudes not exceeding 2,000 meters.

The installation category (transient overvoltage category) for this instrument is II and it is classified as portable equipment. The instrument has a pollution degree rating of 2 and may be installed in an environment that has non-conductive pollutants only.

Electrical Requirements

The electrical receptacle should have a dedicated electrical line with a circuit breaker rated at 30A (with a power rating of 6.0 kVA). There should also be an isolated electrical ground between the instrument and the building main electrical service. The electrical receptacle must be located within 3 m (10 ft.) of the instrument rear panel. The following table specifies the electrical operating range for various locations:

Voltage (AC)	Frequency (Hz)	Location
200 +/- 10%	50/60 +/- 1%	Japan
208 +/- 10%	50/60 +/- 1%	USA/Canada
220 +/- 10%	50/60 +/- 1%	pre-1992 Europe
230 +/- 10%	50/60 +/- 1%	EC
240 +/- 10%	50/60 +/- 1%	Australia and pre-1992 UK

CAUTION Connecting the instrument to the wrong voltage source can damage the instrument. Do not assume that the laboratory's voltage supply is correct. Always measure the voltage at the wall receptacle before configuring the power supplies, connecting the power cord, and turning the instrument on.

In the U.S.A., Canada, and Japan, the instrument is supplied with a detachable electrical cord equipped with a NEMA L6-30P 30A/250V twistlock plug. The electrical receptacle which accepts this plug is the NEMA L6-30R twistlock receptacle (30A/250 V). This receptacle is not supplied by Applied Biosystems. It is available from Hubbell, Inc. (part number 2620A).

In Europe and Australia, the instrument is supplied with a detachable electrical cord equipped with a IEC-309 30A/240V plug. The electrical receptacle which accepts this plug is available in two types, panel mount (part number 58074) or surface mount (part number 58004) These receptacles are not supplied by Applied Biosystems. They are available from any LEGRAND distributor.

Voltage Quality

Line voltage must be within 10% of the nominal value. High or low voltages may have adverse effects on the electronic components of the ABI PRISM 7700 instrument.

A dedicated line and isolated ground between the instrument and building main electrical service are necessary to help prevent such "voltage spikes". If the laboratory environment contains devices which are electrically "noisy", or if you are in an area with frequent electrical storms, a line conditioner may enhance the system's reliability. A 5 kVA capacity is generally recommended to help ensure the proper operation of the instrument.

Certain types of electrical noise are greatly exaggerated by poor or improper electrical ground connections. To prevent these problems, it is very important to have a dedicated line and isolated ground between the instrument and building main electrical service.

General Warnings

Chemical Hazard Warning ! **WARNING ! CHEMICAL HAZARD.** Some of the chemicals used with Applied Biosystems instruments are potentially hazardous and can cause injury, illness or death.

- ◆ Read and understand the material safety data sheets (MSDSs) provided by the chemical manufacturer before you store, handle, or work with any chemicals or hazardous materials.
- ◆ Minimize contact with and inhalation of chemicals. Wear appropriate personal protective equipment when handling chemicals (e.g., safety glasses, gloves, or clothing). For additional safety guidelines consult the MSDS.
- ◆ Do not leave chemical containers open. Use only with adequate ventilation.
- ◆ Check regularly for chemical leaks or spills. If a leak or spill occurs, follow the manufacturer's cleanup procedures as recommended on the MSDS.
- ◆ Comply with all local, state/provincial, or national laws and regulations related to chemical storage, handling, and disposal.

Chemical Waste Hazard Warning ! **WARNING ! CHEMICAL WASTE HAZARD.** Wastes produced by Applied Biosystems instruments are potentially hazardous and can cause injury, illness, or death.

- ◆ Read and understand the material safety data sheets (MSDSs) provided by the manufacturers of the chemicals in the waste container before you store, handle, or dispose of chemical waste.
- ◆ Handle chemical wastes in a fume hood.
- ◆ Minimize contact with and inhalation of chemical waste. Wear appropriate personal protective equipment when handling chemicals (e.g., safety glasses, gloves, or clothing).
- ◆ After emptying the waste container, seal it with the cap provided.
- ◆ Dispose of the contents of the waste tray and waste bottle in accordance with good laboratory practices and local, state/provincial, or national environmental and health regulations.

Electrical Shock Hazard ! **WARNING ! ELECTRICAL SHOCK HAZARD.** Severe electrical shock, which could cause physical injury or death, can result from working on an instrument when the high voltage power supply is operating. To avoid electrical shock, disconnect the power supply to the instrument, unplug the power cord, and wait at least 1 minute before working on the instrument.

Setting Up and Testing the System

3

Overview

About This Chapter This chapter describes the instrument electrical configuration, and the setup of the computer and software. It also contains procedures that test: the electrical system, the thermal cyclers, the laser alignment, and the fluorescence and background.

In This Chapter This chapter contains the following topics:

Topic	See Page
Getting Started	3-2
Electrical Configuration	3-3
Computer and Software Setup	3-6
System Testing At Startup	3-8
Thermal Cycler Tests	3-9
Laser Alignment Test	3-10
Fluorescence and Background Checks	3-11

Getting Started

Checking Equipment

To make sure equipment and supplies are ready for installation and system testing:

1	Verify that the safety and space requirements have been met per the ABI PRISM 7700 <i>Site Preparation and Safety Guide</i> to ensure there is access to the instrument on all sides.
2	Verify that all equipment, chemicals, and supplies have arrived.
3	Verify that all customer-supplied equipment and supplies needed to run the chemical installation kit are readily available. Refer to the preinstallation checklist in "Preinstallation Checklists" on page D-2.
4	Unpack and setup the Power Macintosh® computer. This should include reformatting the hard drive, reloading the Macintosh® computer operating system, and installing the Sequence Detection System software.
5	Unpack the ABI PRISM® Sequence Detector.

Electrical Configuration

Checking the Electrical Receptacle

Inspect the electrical receptacle where the instrument is to be plugged in, to ensure that it matches the power cord plug. In the USA, Canada, and Japan, the receptacle should be a NEMA L6-30R (rated at 30A/250V). In Europe and Australia, the receptacle should accept an IEC-309 30A/250V plug. Use of receptacles other than these may void the instrument's safety certification and should not be used.

CAUTION Always check the laboratory voltage at the wall receptacle with a voltmeter before configuring the instrument; incorrect voltage can lead to instrument damage.

Plugging the Instrument In for the First Time

To connect the instrument to the power source:

Step	Action
1	Use a DVM to check voltage at the electrical receptacle.
2	After confirming the AC voltage, plug the power cord into the back of the instrument and push the instrument's electrical breaker down to the OFF position.
3	Plug the power cord into the wall receptacle.

Configuring the Internal Electrical Connections

Before powering up the instrument, check and configure electrical connections to ensure proper instrument operation.

To configure the internal electrical connections:

Step	Action						
1	Open the main cover on the instrument (secured by four screws).						
2	Inspect all electrical harnesses, especially those connected to the main CPU PCB, to ensure that the connections are secure.						
3	Locate the instrument voltage configuration plug in the packing kit, and check that it is the proper plug by comparing it to the illustration in "7700 Instrument Voltage Configuration Plugs" on page 3-4.						
4	Inspect the front of the Laser Power Supply and set the voltage configuration switch to the setting that corresponds to the laboratory electrical supply. <table border="1" data-bbox="581 1346 1161 1465"><thead><tr><th>Use this setting...</th><th>For these areas...</th></tr></thead><tbody><tr><td>LOW</td><td>200–208 Vac</td></tr><tr><td>HIGH</td><td>220–240 Vac</td></tr></tbody></table> <p>CAUTION Always check the laboratory voltage at the wall receptacle with a voltmeter before configuring the instrument; incorrect voltage can lead to instrument damage.</p>	Use this setting...	For these areas...	LOW	200–208 Vac	HIGH	220–240 Vac
Use this setting...	For these areas...						
LOW	200–208 Vac						
HIGH	220–240 Vac						

To configure the internal electrical connections: *(continued)*

Step	Action						
5	<p>The thermal cycler operates on 230 Vac (supplied from the Vac Distribution Assembly), but it must be configured for the proper frequency (either 50 or 60 Hz).</p> <p>Insert a properly configured 36-pin voltage select plug into J-13 on the AC-Control PCB in the 9600. The following table lists the plug and the pins that are connected:</p> <p>Note The plug has five jumpers with each wire connecting two pins.</p> <table border="1"> <thead> <tr> <th>This plug...</th> <th>Has these pins jumpered...</th> </tr> </thead> <tbody> <tr> <td>230 V/50 Hz plug</td> <td>2-5, 9-17, 28-34, 31-36, and 33-35 jumpered</td> </tr> <tr> <td>230 V/60 Hz plug</td> <td>2-5, 9-16, 28-34, 31-36, and 33-35 jumpered</td> </tr> </tbody> </table> <p>For configuration information, see “7700 Instrument Voltage Configuration Plugs,” below.</p> <p>CAUTION Damage to the 9600 may occur if the wrong voltage configuration plug is inserted into J-13. Always check the wiring of the 9600 voltage configuration plug before installing it in the instrument.</p>	This plug...	Has these pins jumpered...	230 V/50 Hz plug	2-5, 9-17, 28-34, 31-36, and 33-35 jumpered	230 V/60 Hz plug	2-5, 9-16, 28-34, 31-36, and 33-35 jumpered
This plug...	Has these pins jumpered...						
230 V/50 Hz plug	2-5, 9-17, 28-34, 31-36, and 33-35 jumpered						
230 V/60 Hz plug	2-5, 9-16, 28-34, 31-36, and 33-35 jumpered						
6	<p>The thermal cycler must also be configured for the proper electrical frequency (either 50 or 60 Hz.).</p> <p>See “Configuring the Thermal Cycler Frequency for 50 hertz Areas” on page 3-5.</p>						

7700 Instrument Voltage Configuration Plugs

CAUTION If the ABI PRISM 7700 instrument is not configured properly, the voltage supplied to the thermal cycler may be incorrect, resulting in electronics or firmware damage.

Listed from left to right, the voltage select plugs are for 200, 208, 220, and 240 Vac.

1	4	7	1	4	7	1	4	7	1	4	7
2	5	8	2	5	8	2	5	8	2	5	8
3	6	9	3	6	9	3	6	9	3	6	9

ABI PRISM 7700 Voltage Configuration Plugs

Pins	Description...
#1 and #3	are oriented to the left and have a flat side.
Jumpered pins (shown in shaded colors)	Are for this voltage...
2-5, 4-9, and 6,7	200 V
2-5,3-9, and 6,7	208 V
2-5, 1-9, 6-7	220 V
5-9 and 6-7	230/240 V.

Configuring the Thermal Cycler Frequency for 50 hertz Areas


After the thermal cycler voltage-frequency configuration plug has been installed and the instrument is plugged in, the thermal cycler firmware must also be configured if the instrument is to be run in 50 hertz areas.

To configure the 7700 instrument thermal cycler's electrical frequency:

Step	Action
1	<p>Start the thermal cycler.</p> <ul style="list-style-type: none">a. Plug a 9600 keyboard assembly into the thermal cycler CPU PCB.b. While holding the keyboard ENTER key, turn the 7700 instrument on. <p>The power switch is located at the lower left front of the instrument. The instrument circuit breakers are located at the back near the power cord receptacle.</p> <p>The 9600 display will appear when the key is released.</p>
2	<p>Calibrate the thermal cycler for 50 hertz areas.</p> <ul style="list-style-type: none">a. After the 9600 display comes on, press MORE 999 to access the 9600 service diagnostics.b. Press the OPTION key twice to move the cursor to the CALIB selection and then press ENTER.c. After the first calibration test appears, press STEP to go to Calibration Test #2.d. Press ENTER to view the present voltage and frequency configuration.e. Press the STEP or OPTION keys until the correct configuration is reached.f. Press STOP 5 times to reset the thermal cycler.g. Turn off the instrument and remove the 9600 keyboard when you are finished.

Computer and Software Setup

System Setup Cable connections of the ABI PRISM 7700 Sequence Detection System.



NOTE: Depending on the type of monitor used, its power cord MAY plug into the back of the PowerPC Macintosh.

Setting Up the Computer

Step	Action
1	Follow the basic setup procedure in the Macintosh computer documentation shipped with the computer. CAUTION Always check the voltage configuration at the back of the computer to ensure that it matches the site voltage. Failure to properly set the voltage switch can cause damage to the computer.
2	Connect the computer to the 7700 instrument with the 6 ft. cable found in the instrument packing kit. Attach the cable to the modem port at the back of the computer, and then to the computer port on the back of the instrument.
3	Turn on the monitor.
4	Turn on the computer by pressing the Power key on the keyboard. The power key is marked with a triangle.
5	Install the Macintosh Computer Operating System: a. Insert the Apple CD that contains the operating system. b. Restart the computer while pressing the c key. c. Reinitialize the hard drive (see the Disk Tools folder). d. Reload the Macintosh computer operating system: "INSTALL O.S." e. Restart the computer.

Installation of the Software

Perform the following steps to load the SDS application on the hard disk. This procedure should also be used if the application is corrupted and must be reloaded.

Step	Action
1	Insert the ABI PRISM 7700 CD-ROM into the CD-ROM drive.
2	Use the Installer program to load the application. After installation is complete, restart the computer. Note Always run Norton Utilities from the CD. Do not install Norton Disk Utilities or Norton Disk Doctor on the hard drive.
3	Microsoft Excel is shipped with the 7700 instrument, but it is not used for all service procedures. Install Excel only if performing a TaqMan card upgrade, or if the customer requests it.
4	Open the Utilities folder and drag the Techtool® application to the hard drive. Techtool should be used once a week to rebuild the desktop.
5	Go to the Energy Save Control Panel and set to NEVER for all settings.
6	Go to the Memory Control Panel and turn Virtual Memory off, and ensure that there is at least 32 MB of built-in RAM.
7	Open the Extensions Manager Control Panel , and check that the PowerPC Interrupt extension is turned on (restart the computer if necessary). Note The Power PC Interrupt extension is not needed with Macintosh OS 7.6 and newer.
8	Rebuild the Macintosh computer desktop when the installation is complete. Also, rebuild the Macintosh computer desktop as needed in the event that unexpected computer crashes cause operational problems.
9	Restart the computer before starting any run for service purposes.

System Testing At Startup

Power-Up Tests Ensuring that the electronics and electrical system of the 7700 instrument are functioning properly is important. The following procedure can be used to check the system electronics.

To perform the 7700 instrument Power-up Test:

Step	Action						
1	Turn on the instrument.						
2	Look at the three (3) LEDs on the front panel near the heated cover: <ul style="list-style-type: none"> ◆ If the READY LED is lit, the instrument has booted up properly. ◆ The COMM LED will blink when the 7700 instrument is communicating with the computer. 						
3	Open the instrument main cover, and check that the 683324 PCB LEDs are “marching”, which is the ready state. (This confirms that the 7700 instrument has booted-up properly.)						
4	After the 7700 instrument has booted-up, launch the Sequence Detection System (SDS) software on the computer and open a background plate. <p>Note When the SDS is launched a dialog boxes appears with the message, “Could not open document because needed Pure Spectra are missing (FAM, TAMRA). New Pure Spectra must be extracted.” Click OK.</p> <p>To open a background plate:</p> <ol style="list-style-type: none"> a. From the File menu, select New Plate. b. Select background plate and click OK. <p>Note The following steps require an open plate document, although you will not actually perform a run.</p>						
5	<p>Run Instrument Verification.</p> <ol style="list-style-type: none"> a. Click the Show Analysis button. b. From the Instrument menu, select Diagnostics, and 7700 instrument Verification from the submenu. The Instrument Tests dialog box appears. c. Click the Run Tests button. <table border="1" style="width: 100%;"> <thead> <tr> <th>If...</th> <th>Then...</th> </tr> </thead> <tbody> <tr> <td>the test passes</td> <td>proceed to the “Thermal Cycler Tests” on page 3-9.</td> </tr> <tr> <td>any test fails</td> <td>ensure that the ABI PRISM 7700 instrument booted from the proper firmware.</td> </tr> </tbody> </table>	If...	Then...	the test passes	proceed to the “Thermal Cycler Tests” on page 3-9.	any test fails	ensure that the ABI PRISM 7700 instrument booted from the proper firmware.
If...	Then...						
the test passes	proceed to the “Thermal Cycler Tests” on page 3-9.						
any test fails	ensure that the ABI PRISM 7700 instrument booted from the proper firmware.						

Thermal Cycler Tests

System Performance Test

Note Thermal Cycler tests, which are run with the use of a 9600 Service Keyboard Assembly, ensure that the internal ABI PRISM 7700 thermal cycler operation is not being affected by the 7700 CPU. The 9600 Service Keyboard Assembly must be used because the 7700 instrument contains a modified 9600 that does not have a keyboard and display.

Step	Action						
1	Make sure that the instrument power is off.						
2	Plug a 9600 keyboard assembly into the thermal cycler CPU PCB.						
3	<p>a. While holding the keyboard ENTER key down, turn the power switch on. The power switch is located at the lower left front of the instrument.</p> <p>b. Release the ENTER key after a few seconds.</p>						
4	After the 9600 display appears, let the 7700 thermal cycler warm up for at least 5 minutes before running any tests.						
5	Press the Option key to step to UTIL , and then to the DIAG menu to access the user diagnostics.						
6	Select the System Performance test , and then click OK .						
7	<p>Use the System Performance test result to determine the next step:</p> <table border="1"> <thead> <tr> <th>If...</th> <th>Then...</th> </tr> </thead> <tbody> <tr> <td>the System Performance test passes</td> <td>proceed to the "Laser Alignment Test" on page 3-10.</td> </tr> <tr> <td>the System Performance test fails</td> <td>proceed to step 8.</td> </tr> </tbody> </table> <p>Note System Performance test results can be unreliable on some newer thermal cycler compressors, which are more powerful than older compressors. If a failure is encountered, run the Chiller test and Heater test before repairing any thermal cycler.</p>	If...	Then...	the System Performance test passes	proceed to the "Laser Alignment Test" on page 3-10.	the System Performance test fails	proceed to step 8.
If...	Then...						
the System Performance test passes	proceed to the "Laser Alignment Test" on page 3-10.						
the System Performance test fails	proceed to step 8.						
8	<p>Run the Chiller test and the Heater test.</p> <p>Note The Chiller test and Heater test can also fail due to the new compressors and do not necessarily mean that the thermal cycler has a problem. Run the chemical installation kit to verify that there is, in fact, a problem with the thermal cycler before taking corrective action.</p>						
9	<p>If the results of the chemical installation run indicate that the thermal cycler has a problem, then take the following action:</p> <p>a. Unless a hard failure is observed, perform a temperature verification run before replacing components. See Chapter 5, "Verifying the Instrument."</p> <p>b. Determine what the problem is by inspecting the thermal cycler profile for the verification run from either the Raw Data view or Multicomponent view.</p>						
10	<p>Additional 9600 Service Tests can be run on the 7700 instrument using the 9600 Service Diagnostics.</p> <p>Note The 9600 Service Diagnostic Install Test will fail in many cases, because the heated cover is different in the 7700 instrument. If the test fails, press MORE 999 to bypass the test and automatically set the install flag.</p>						

Laser Alignment Test

About Checking the Alignment

The ABI PRISM 7700 instrument's laser and optics were aligned during the manufacturing process and probably do not need to be realigned. However, the alignment should be checked to ensure that the instrument was not damaged during shipment.

Use the following two procedures to check the alignment of the laser. The first procedure prepares the instrument for the alignment check. The second procedure checks the alignment of the laser.

Preparing the 7700 Instrument for Alignment



Step	Action
1	Make sure the instrument is off.
2	Remove any articles from the top of the 7700 instrument and raise the cover.
3	Remove the beam path cover.
4	Attach the laser alignment tool to the front of the optical cube.
5	Post a laser warning sign and put on laser safety goggles. ! WARNING ! VISION HAZARD: LASER SOURCE. When the laser beam is exposed, it is Class IIIb, which is extremely hazardous. Direct or reflected laser beam exposure at 10mW for 0.1 second can burn the eye's retina, leaving a permanent blind spot. Eye Hazard.
6	By-pass the laser interlock switches by either <ul style="list-style-type: none"> ◆ Pulling the heated cover forward or ◆ Using the Laser Interlock Defeat Switch.
7	Turn the instrument on.

Checking Alignment

1	Start the Sequence Detection System software on the Power Macintosh® computer.
2	Go to the Functions Test menu.
3	Open the Shutters, and check the position of the laser beam on the front and back targets. IMPORTANT Do not change the laser alignment if the laser beam is not centered exactly on the middle of the laser targets. The final laser alignment is never verified or adjusted using the laser alignment tool.
4	Close the shutters, remove the laser alignment tool, and put a Fluorescent Test tray in the sample block, and close the heated cover.
5	Open the shutters and read ("poll") the CCD at various MUX positions. Signal should be seen on all wells that contain a fluorescent sample. To check the laser alignment, use the 7700 Service Diagnostics Application Plate Read Tests.
6	When the alignment has been checked, replace the optical path cover, and close the 7700 instrument's main cover. Note If the laser requires alignment, refer to Service Bulletin #12.

Fluorescence and Background Checks


Fluorescent Tests (Signal Throughput)

Test	Description
Fluorescent Test tray (P/N T-6222)	<p>Signal throughput can be checked using a Fluorescent Test tray.</p> <p>The example in the figure below was obtained by running the T-6222 at 25 °C for 1 minute (open a background plate to collect data).</p> <p>Note Signal heights of the fibers vary depending on the well position, the specific 7700 instrument being used, and the temperature of the test tray.</p>  <p>The figure is a line graph titled "Raw Spectra - T-6222/Bkgrd/25C/1min/002". The y-axis is labeled "Fluorescent Intensity" and ranges from 0 to 22000 in increments of 2000. The x-axis is labeled "Wavelength (nm)" and ranges from 500 to 660 in increments of 10. Multiple colored lines represent different fiber measurements. All lines show a broad peak between 520 and 550 nm, a smaller peak around 590 nm, and a very sharp, high-intensity peak at approximately 615 nm. The peak intensity at 615 nm varies significantly between the different fiber measurements, ranging from about 8000 to 19000.</p>
Fluorescent Tests (Red Fluorescent Source)	<p>Signal throughput results can be obtained by scanning any Fluorescent Test tray for one minute at 25 °C. Use a background plate document to collect data.</p> <p>The figure below shows that the results vary depending on the type of tray being scanned.</p> <p>IMPORTANT Blue test trays or paper can be used to confirm signal throughput, but do not use them to confirm laser alignment.</p>  <p>The figure is a line graph titled "Raw Spectra - red tray 10ms". The y-axis is labeled "Fluorescent Intensity" and ranges from 0 to 16000 in increments of 2000. The x-axis is labeled "Wavelength (nm)" and ranges from 500 to 660 in increments of 10. Multiple colored lines represent different fiber measurements. The curves are broad and relatively flat, starting to rise around 550 nm and peaking between 610 and 640 nm. The peak intensity varies between approximately 8000 and 14000.</p>

**Blackboard Tests
(System Noise)**

The Blackboard test checks the background signal of the system (see the figure below). Use a background plate document in the Sequence Detection System software to scan a Blackboard for one minute at 25 °C. This will provide a good indication of the background signal.

Note The maximum height of the highest peak should be less than 1200 fluorescent units. If the system noise results are poor, refer to the *ABI PRISM 7700 SEQUENCE DETECTION SYSTEM SERVICE MANUAL* for troubleshooting.



Note Signal heights vary depending on the well position, the specific ABI PRISM 7700 instrument being used, and the temperature of the test tray. Results obtained from using a blackboard do not take into account background signal from the sample block or real PCR samples.

Performing a Spectral Calibration

4

Overview

About This Chapter This chapter describes how to archive the current component files, generate a Background Component file, and generate a Pure Dye file.

In This Chapter This chapter contains the following topics:

Topic	See Page
Spectral Calibration Overview	4-2
Archiving Current Component Files	4-3
Materials	4-2
Generating the Background Component File	4-4
Generating a Pure Dye File	4-6

Spectral Calibration Overview

Spectral Calibration Files The Spectral Calibration files are used by the ABI PRISM® Sequence Detection System software for data analysis.

When a spectral calibration is performed, two files are created:

- ◆ Background Component
- ◆ Pure Dye

These files are stored in a folder called Spectra Components, which is automatically created when the ABI PRISM 7700 software is launched for the first time.

The Spectra Components folder is in the following location on the hard drive:

System Folder > Preferences > SDS

Spectral Calibration Process To perform a spectral calibration:

- ◆ Archive or delete the current Spectra Components folder (if one already exists)
- ◆ Generate a Background Component file
- ◆ Generate a Pure Dye file

Materials

Materials Required Materials Required To Perform a Spectral Calibration

Materials	Source
SDS Spectral Calibration Kit	Applied Biosystems (P/N 4305822)
MicroAmp® Optical 96-Well Reaction Plate/Optical Caps	Applied Biosystems (P/N 403012)
Gloves, disposable, powder-free	MLS
Pipettors, positive-displacement or air-displacement	MLS
Pipette tips, aerosol resistant	MLS
Water, RNase-free, distilled, deionized	MLS

Archiving Current Component Files

Archiving To archive an existing Spectral Components folder:

Step	Action
1	Make sure that the SDS software is not running.
2	Navigate to the Spectra Components folder: System Folder > Preferences > SDS > Spectra Components
3	Click on the Spectra Components icon text and enter a new name for the folder.

Generating the Background Component File

Preparing Instrument and Background Tray

Step	Action
1	Turn on the ABI PRISM 7700 instrument and allow it to warm up for at least 30 minutes before starting the first run.
2	Prepare the background tray. <ol style="list-style-type: none">Pipette 50μL water into each well of a 96-well optical tray.Cover wells with optical caps.Make sure there are no bubbles in the wells.

Run the Background Tray


Step	Action
1	Launch the SDS software. Note A dialog box will appear with the message "Could not open document because needed Pure Spectra are missing (FAM, TAMRA). New Pure Spectra must be extracted." Click OK and proceed.
2	Open a new background plate. <ol style="list-style-type: none">From the File menu, choose New Plate...Select Background from the Plate Type pop-up menu.Click OK.
3	Place the tray in the sample block.
4	Toggle to the Analysis View.
5	Click the Run button.

Analyze the Run

To analyze the run:

Step	Action
1	Make sure the Analysis View is showing.
2	Highlight all 96 wells.
3	From the Analysis menu, select Raw Spectra .

To analyze the run: *(continued)*

Step	Action
4	<p>Verify that all 96 curves are of similar shape and acceptable value.</p> <p>IMPORTANT What is the spec? (≤ 1200 according to Jake F.) The background for all 96 wells should be similar in shape or the instrument may not analyze data properly. The fluorescent intensity of the background must be no more than 1200 for all wells. If there is a problem with these results, refer to the <i>ABI PRISM 7700 Sequence Detection System Service Manual</i>.</p> 

Extract and Save the Background Component Data


Step	Action
1	<p>From the Instrument menu, select Calibrate > Extract Background Component.</p> <p>Note At this point, the Background Component file is generated and automatically saved to the Spectra Components folder.</p>
2	To save the background plate document, select Save As... from the File menu.
3	Type a descriptive name for the file.
4	Click Save .
5	From the File menu, select Close to close the document.
6	Quit the SDS software.

Generating a Pure Dye File

Preparing a Pure Dye Tray Prepare a pure dye tray using the SDS Spectral Calibration Kit (P/N 4305822) that is shipped with the chemical installation kit.

Step	Action
1	Pipette 50 μ L of each of the seven dye standards into 4 separate wells. The result will be 28 wells containing 50 μ L each. Refer to the Plate Document Setup diagram (page 4-6, step 5) for placement suggestion. Note Do not dilute the dye solutions.
2	Cap the tray and make sure there are no bubbles in the wells.
3	Place the tray into the sample block.

Preparing a Pure Spectra Plate Document To prepare a Pure Spectra plate document:


Step	Action
1	Launch the Sequence Detection System software. Note A dialog box appears with the message “Could not open document because needed Pure Spectra are missing (FAM, TAMRA). New Pure Spectra must be extracted.” Click OK and proceed.
2	Open a new Pure Spectra plate. a. From the File menu, choose New Plate . b. Select Pure Spectra from the Plate Type pop-up menu. c. Click OK .
3	Open the Sample Type pop-up menu and verify that all seven of the following dyes are listed: FAM, JOE, ROX, SYBR, TAMRA, TET, VIC Note If all dyes are present, skip to “Assigning Dye Standards To the Plate Document” on page 4-7. If any of the dyes are missing, proceed to the next step.
4	To add dyes to the list, select Sample Type Setup... from the Sample Type pop-up menu. The following dialog box appears: 
5	Click Add . A new row appears at the bottom of the dye list. 

To prepare a Pure Spectra plate document: *(continued)*



Step	Action
6	<p>Click the acronym text field, and enter a name for the new dye that is no more than 5 characters long (<i>i.e.</i>, VIC or SYBR).</p> <p>Acronym text field</p> 
7	<p>Click the name text field and enter Pure Dye.</p> <p>Name text field</p> 
8	<p>Click the color field.</p> <p>Color field</p>  <p>The color palette dialog box appears.</p>
9	<p>Select a color for the new dye and click OK. The color field for the new dye fills with the new color.</p> <p>Note The color field must be changed, because the default color is gray, which is the same as the background.</p> <p>For more information on selecting a dye color, refer to the <i>ABI PRISM 7700 Sequence Detection Systems User's Manual</i>.</p>
10	<p>Repeat step 5–9 to add other dyes to the list. Click OK when finished.</p>

Assigning Dye Standards To the Plate Document


To assign dye standards to the plate document:

Step	Action
1	<p>From the Setup menu, select Sample Type Palette.</p> <p>The Sample Type Palette dialog box appears:</p> 

To assign dye standards to the plate document: *(continued)*

Step	Action
2	<p>Select the four wells from the plate document that correspond to the wells on the MicroAmp® Optical 96-Well Reaction Plate containing the FAM standard.</p> 
3	<p>Click the FAM checkbox in the palette dialog box to label the wells.</p> <p>Click Update from the palette box. The Sequence Detection System software designates the selected row with the FAM dye.</p> 

To assign dye standards to the plate document: *(continued)*

Step	Action
4	<p>Repeat steps 1–3 for each additional dye standard (JOE, ROX, SYBR Green, TAMRA, TET, and VIC) until your plate document resembles the figure below.</p> 

Running the Samples and Creating the Pure Dye File

To run the samples and create the Pure Dye file:

Step	Action
1	Click the Show Analysis button and click the Run button.
2	<p>Extract pure dye.</p> <ol style="list-style-type: none"> Highlight the first set of four wells that contain the FAM dyes. Select Calibrate from the Instrument menu, and select Extract Pure Dye from the submenu.
3	<p>When the spectra for the wells are displayed, inspect the quality of the four curves. The spectra of the four wells must be smooth and similar in shape, or the instrument will not analyze data correctly.</p> <p>Note Contaminated wells may produce spectra that deviate from average peaks. If a spectrum is unacceptable, move it to the discard box as described in “Removing an Outlier Spectrum From the Averaged Box” on page 4-10.</p> <p>IMPORTANT You may remove only one of the four spectra from each dye set.</p>
4	<p>Repeat steps 2–3 for the other six dyes.</p> <p>When finished, the Pure Dye file is automatically saved in the Spectra Components folder.</p>
5	<p>Save the Pure Dye plate document.</p> <ol style="list-style-type: none"> From the File menu, select Save As. A dialog box appears. Enter Pure Dye in the Save this document as text box, and click Save.

To run the samples and create the Pure Dye file: *(continued)*


Step	Action
6	Quit the SDS software when finished. A new Spectra Components folder is now in the SDS Preferences folder. Note The Pure Dye file should be replaced every 60 days to six months.

Removing an Outlier Spectrum From the Averaged Box

Step	Action
1	Click the outlier spectrum in the Averaged box.
2	Click the <input type="button" value="=>"/> button to move the spectrum to the discard box. You may remove only one of the four dyes from any particular set. Note Discarding a spectrum does not remove it from the graph, but the average is automatically recalculated.

Pure Dye Spectra and Peak Wavelength

Dyes	Peak Wavelength
FAM	~520 nm
JOE	~550 nm
ROX	~610 nm
SYBR Green	~520 nm
TAMRA	~580 nm
TET	~540 nm
VIC	~550 nm



**Wavelength Axis
Purpose**

The wavelength axis should not be used as an exact measure of wavelength but is intended as a guide to where an individual spectrum should appear.

When the ABI PRISM 7700 instrument is spectrally calibrated, the wavelength positions of different spectral emissions do not change. The spectrograph and CCD camera are bolted in place, and the pixel positions of the CCD camera cannot be changed by software commands.

Verifying the Instrument

5

Overview

About This Chapter This chapter describes how to prepare a β -actin installation plate, set up and start the β -actin installation run, and prepare and run the RNase P installation plate. It also includes how to analyze the data and verification calculations.

In This Chapter This chapter contains the following topics:

Topic	See Page
The ABI Prism 7700 Instrument Verification Run	5-2
Preparing the β -Actin Installation Plate	5-3
Setting Up and Starting the β -Actin Installation Run	5-6
Preparing and Running the RNase P Instrument Verification Plate	5-9
Analyzing Data	5-13
7700 Sequence Detection System Verification Calculations	5-15

The ABI PRISM 7700 Instrument Verification Run

Chemical Installation Kits Available

There are two chemical installation kits available for use in verification of the ABI PRISM® 7700 Sequence Detector. Either one may be used. The following tables list the steps to completing each type of verification run:

◆ β -Actin

Method	See Page
Preparing the β -Actin Installation Plate	5-3
Setting Up and Starting the β -Actin Installation Run	5-6
Analyzing Data	5-13

◆ RNase P

Method	See Page
Preparing and Running the RNase P Instrument Verification Plate	5-9
Analyzing Data	5-13

Preparing the β -Actin Installation Plate

Chemistry Preparation Guidelines

Guidelines When Preparing the ABI PRISM 7700 β -Actin Chemical Installation Plate

Guideline	Description
Have experience with PCR samples	The person preparing an installation tray should have experience in making up PCR samples.
Use sterile supplies and techniques	Sterile laboratory supplies and techniques are needed to help ensure the success of the TaqMan PCR reactions.
Have adequate master mix on hand	Adequate master mix should be made up to ensure sufficient solution is available to complete the tray.
Keep individual reagents on ice	When making up the standards and Unknowns, keep the individual reagents on ice until they are added to the tray. IMPORTANT Put the master mix on ice and any unused enzymes in the freezer in order to maximize their activity.
To prevent denaturation	Add the enzymes last and mix gently to help prevent denaturation.
Cap wells after filling	Cap all 96 wells with optical caps after the tray is filled with the appropriate reaction solutions.
Use dimpled-end of installation tool	Use the dimpled-end of the cap-installation tube to secure the caps after capping the PCR tubes.
Be careful not to splash reaction solutions	! WARNING ! CHEMICAL HAZARD. All chemicals on the instrument, including liquid in the lines, are potentially hazardous. Always determine what chemicals have been used on the instrument before changing reagents or instrument components. Wear appropriate eyewear, clothing, and gloves when working on the instrument. When the installation tray is removed from the base, be careful not to splash the reaction solutions on the inside of the optical caps. Trays that have extensive splashing should be centrifuged before they are run.

Preparing Reagents

Heat the DNA templates at 37°C for 20 minutes, then vortex and spin to ensure that the DNA goes into solution.

Making Master Mix Make the master mix according to the table below. You will have enough master mix for 112 reactions, to allow for loss during pipetting.

Note Reagents identified by a single asterisk (*) are found in the TaqMan[®] DNA Template Reagents Kit. Reagents identified by a double asterisk (**) are found in the TaqMan[®] Core Reagents Kit.

! WARNING ! CHEMICAL HAZARD. Be sure to familiarize yourself with the MSDSs before using reagents or solvents.

Component	Amount (μL)	Action
MgCl ₂ **	784	Vortex all tubes, except those containing enzymes, to mix contents thoroughly.
Water	1484	
TaqMan [®] Buffer A**	560	Spin all tubes to collect contents at the bottom before pipetting.
Forward Primer*	560	
Reverse Primer*	560	When delivering reagents to the master mix, pipette them up and down several times to mix, then evacuate all from the pipette.
β-Actin Probe*	560	
dATP/dCTP/dGTP/dUTP**	112 of each	
UNG (Enzyme) **	56	
AmpliTaq Gold [®] (Enzyme)**	28	
Total Volume	4480	

IMPORTANT Add the enzymes last and mix gently to help prevent denaturation. Put the master mix on ice and any unused enzymes in the freezer to maximize their activity.

Preparing Samples Mix Reagents

Sample	Add...	To...
Unknown 1	200-μL Unknown Template	1800-μL master mix
Unknown 2	100-μL Unknown Template and 100-μL TE buffer or water	1800-μL master mix
NTC	24-μL TE Buffer or water	216-μL master mix
1K Standard	24-μL Template 1	216-μL master mix
2K Standard	24-μL Template 2	216-μL master mix
5K Standard	24-μL Template 3	216-μL master mix
10K Standard	24-μL Template 4	216-μL master mix
20K Standard	24-μL Template 5	216-μL master mix

Making up the PCR Plate

Pipette 50 μ L of each sample into the appropriate wells of a 96-well optical plate. Use the setup as indicated in Table 5-1 and described below.

Item	Description	Plate Location	Volume to Pipette Into Each Well
Unknown 1	10K copy number	A 1 - 12 B 1 - 12 C 1 - 12	50 μ L
Unknown 2	5K copy number	F 1 - 12 G 1 - 12 H 1 - 12	50 μ L
NTC	No Template Control	D 1 - 4	50 μ L
Template 1	1K Standard	D 5 - 8	50 μ L
Template 2	2K Standard	D 9 - 12	50 μ L
Template 3	5K Standard	E 1 - 4	50 μ L
Template 4	10K Standard	E 5 - 8	50 μ L
Template 5	20K Standard	E 9 - 12	50 μ L

Table 5-1 Sample Setup in 96-Well Optical Tray

	1	2	3	4	5	6	7	8	9	10	11	12	
A													
B				Unknown #1(36 PCR tubes of 10K)									
C													
D	NTC			Template 1 (1K)				Template 2 (2K)					
E	Template 3 (5K)				Template 4 (10K)				Template 5 (20K)				
F													
G				Unknown #2 (36 PCR tubes of 5K)									
H													

IMPORTANT Make sure that there are no bubbles trapped in the bottom of the wells. These will affect final Ct values and must be removed prior to thermal cycling.

Setting Up and Starting the β -Actin Installation Run

Setting Up the Plate Document


To set up the plate document:

Step	Action																																
1	Launch the SDS application software.																																
2	In the Sample Type pop-up menu, select Sample Type Setup... to verify that TAMRA is the Quencher dye and ROX is the Reference dye. <div data-bbox="548 506 984 852" data-label="Image"> <table border="1"> <thead> <tr> <th>Acronym</th> <th>Name</th> <th>Color</th> <th>Reporter</th> </tr> </thead> <tbody> <tr> <td>IPC+</td> <td>Internal Positive</td> <td>Yellow</td> <td>JOE</td> </tr> <tr> <td>IPC-</td> <td>Internal Positive</td> <td>Blue</td> <td>JOE</td> </tr> <tr> <td>TARG</td> <td>RelQ Target</td> <td>Orange</td> <td>FAM</td> </tr> <tr> <td>ENDO</td> <td>RelQ Endogenous</td> <td>Brown</td> <td>None</td> </tr> <tr> <td>STND</td> <td>Standard</td> <td>Pink</td> <td>FAM</td> </tr> <tr> <td>UNKN</td> <td>Unknown</td> <td>Light Blue</td> <td>FAM</td> </tr> <tr> <td>NTC</td> <td>No Template Control</td> <td>Purple</td> <td>FAM</td> </tr> </tbody> </table> <p>Reference ROX</p> <p>Quencher <input checked="" type="checkbox"/> Quencher TAMRA</p> <p>Add Remove Cancel OK</p> </div>	Acronym	Name	Color	Reporter	IPC+	Internal Positive	Yellow	JOE	IPC-	Internal Positive	Blue	JOE	TARG	RelQ Target	Orange	FAM	ENDO	RelQ Endogenous	Brown	None	STND	Standard	Pink	FAM	UNKN	Unknown	Light Blue	FAM	NTC	No Template Control	Purple	FAM
Acronym	Name	Color	Reporter																														
IPC+	Internal Positive	Yellow	JOE																														
IPC-	Internal Positive	Blue	JOE																														
TARG	RelQ Target	Orange	FAM																														
ENDO	RelQ Endogenous	Brown	None																														
STND	Standard	Pink	FAM																														
UNKN	Unknown	Light Blue	FAM																														
NTC	No Template Control	Purple	FAM																														
3	Change the Quencher dye. <ol style="list-style-type: none"> Select the correct dye from the pop-up menu in the Sample Type Setup window. 																																
4	Change the Reference dye. <ol style="list-style-type: none"> Close the Sample Type Setup Window. Switch to Analysis View. From the Instrument menu, select Diagnostics and then Advanced Options. Select the correct dye from the Reference Dye pop-up menu. Click OK. 																																
5	Check that FAM is selected in the Dye Layer pop-up menu of the plate document.																																
6	In Setup View, enter the sample information into the plate document as shown below, using the following steps.																																

To set up the plate document: (continued)

Step	Action												
7	<p>Unknown 1 (10K replicate)</p> <ul style="list-style-type: none"> ◆ Highlight rows A through C. ◆ Select UNKN-Unknown from the Sample Type pop-up menu. ◆ Enter “A” into the Replicate field. 												
8	<p>Unknown 2 (5K replicate)</p> <ul style="list-style-type: none"> ◆ Highlight rows F through H. ◆ Select UNKN-Unknown from the Sample Type pop-up menu. ◆ Enter “B” into the Replicate field. 												
9	<p>NTC (No Template Control)</p> <ul style="list-style-type: none"> ◆ Highlight cells D 1 through 4. ◆ Select NTC-No Template Control from the Sample Type pop-up menu. 												
10	<p>Standards</p> <ul style="list-style-type: none"> ◆ Highlight the appropriate wells. ◆ Select STND-Standard from the Sample Type pop-up menu. ◆ Enter the appropriate copy number into the Quantity field for wells: <table border="1"> <thead> <tr> <th>Well Positions</th> <th>Copy Number</th> </tr> </thead> <tbody> <tr> <td>D 5 – 8</td> <td>1000</td> </tr> <tr> <td>D 9 – 12</td> <td>2000</td> </tr> <tr> <td>E 1 – 4</td> <td>5000</td> </tr> <tr> <td>E 5 – 8</td> <td>10 000</td> </tr> <tr> <td>E 9 – 12</td> <td>20 000</td> </tr> </tbody> </table>	Well Positions	Copy Number	D 5 – 8	1000	D 9 – 12	2000	E 1 – 4	5000	E 5 – 8	10 000	E 9 – 12	20 000
Well Positions	Copy Number												
D 5 – 8	1000												
D 9 – 12	2000												
E 1 – 4	5000												
E 5 – 8	10 000												
E 9 – 12	20 000												

To set up the plate document: *(continued)*

Step	Action
11	<p>From the Setup View, click the Thermal Cycler Conditions button. Set up the conditions as illustrated below.</p> 
12	Check that data will be collected during the PCR extension phase.
13	Check that the value in the Sample Volume text box is set to 50 μL , and click OK .
14	<p>Load the installation plate into the sample block.</p> <p>IMPORTANT Start PCR reactions as soon as possible to prevent degradation of the reaction components.</p>
15	Pull the heated cover forward and secure it in place by turning the knob clockwise to lower and secure it over the sample tray.
16	<p>Start the run.</p> <ol style="list-style-type: none"> Click the Show Analysis button to toggle to the Analysis View. Click the Run button. Wait for the 9600 to beep and the shutters to click before leaving. <p>Note Run will last approximately 2 hours.</p>
17	When run is finished, proceed to “Analyzing Data” on page 5-13.

Preparing and Running the RNase P Instrument Verification Plate

Preparing the Plate **IMPORTANT** The RNase P plate should be stored at 15–20 °C. Do not freeze.

Step	Action
1	Remove the RNase P plate from the freezer, thaw, and mix gently by inversion.
2	Spin the plate briefly in a table top centrifuge with a plate adapter to collect the contents at the bottom of the wells and to remove bubbles.
3	Load the installation tray into the sample block. Note The proper orientation of the plate in the ABI PRISM 7700 instrument is to have the A1 tube in the upper-left corner. Start PCR reactions as soon as possible to prevent degradation of the reaction components.
4	Pull the heated cover forward and secure it in place by turning the knob clockwise to lower and secure it over the sample tray.

Setting Up the Plate Document

To set up the plate document:

Step	Action
1	Launch the SDS application software.
2	Open a new plate. a. Close the plate document that opens automatically. b. Select New Plate from the File menu. The New Plate dialog box appears.
3	Select a plate type and run. ◆ Single Reporter from the Plate Type pop-up menu, ◆ Real Time from the Run pop-up menu, and click OK .

New Plate


Plate Type:

Data Acquisition


Instrument:

Run:

To set up the plate document: *(continued)*

Step	Action
4	<p>From the Sample Type pop-up menu, select Sample Type Setup... to verify that TAMRA is the Quencher dye and ROX is the Reference dye.</p> 
5	<p>Change the Quencher dye.</p> <p>a. Select the correct dye from the pop-up menu in the Sample Type Setup window.</p>
6	<p>Change the Reference dye.</p> <p>a. Close the Sample Type Setup Window.</p> <p>b. Switch to Analysis View.</p> <p>c. From the Instrument menu, select Diagnostics and then Advanced Options.</p> <p>d. Select the correct dye from the Reference Dye pop-up menu.</p> <p>e. Click OK.</p>
7	<p>Check that FAM is selected in the Dye Layer pop-up menu of the plate document.</p>
8	<p>In Setup View, enter the sample information into the plate document as shown below, using the steps on the next page.</p> 

To set up the plate document: *(continued)*

Step	Action												
9	<p>Unknown 2 (5K replicate)</p> <ol style="list-style-type: none"> Highlight rows A through C. Select UNKN-Unknown from the Sample Type pop-up menu. Enter B into the Replicate field. 												
10	<p>Unknown 1 (10K replicate)</p> <ol style="list-style-type: none"> Highlight rows F through H. Select UNKN-Unknown from the Sample Type pop-up menu. Enter A into the Replicate field. 												
11	<p>NTCs (No Template Control)</p> <ol style="list-style-type: none"> Highlight cells D 1 through 4. Select NTC-No Template Control from the Sample Type pop-up menu. 												
12	<p>Standards</p> <ol style="list-style-type: none"> Highlight the appropriate wells. Select STND-Standard from the Sample Type pop-up menu. Enter the appropriate copy number into the Quantity field for wells: <table border="1" data-bbox="581 835 1154 1079"> <thead> <tr> <th>Well Positions</th> <th>Copy Number</th> </tr> </thead> <tbody> <tr> <td>D 5 - 8</td> <td>1250</td> </tr> <tr> <td>D 9 - 12</td> <td>2500</td> </tr> <tr> <td>E 1 - 4</td> <td>5000</td> </tr> <tr> <td>E 5 - 8</td> <td>10 000</td> </tr> <tr> <td>E 9 - 12</td> <td>20 000</td> </tr> </tbody> </table>	Well Positions	Copy Number	D 5 - 8	1250	D 9 - 12	2500	E 1 - 4	5000	E 5 - 8	10 000	E 9 - 12	20 000
Well Positions	Copy Number												
D 5 - 8	1250												
D 9 - 12	2500												
E 1 - 4	5000												
E 5 - 8	10 000												
E 9 - 12	20 000												
13	<p>Check that the thermal cycler conditions are set up as illustrated below.</p> 												
14	Check that data will be collected during all PCR phases.												
15	Check that the value in the Sample Volume text box is set to 50 μ L, and click OK .												

To set up the plate document: *(continued)*

Step	Action
16	Start the run. a. Click the Show Analysis button to toggle to the Analysis View. b. Click the Run button. c. Wait for the 9600 to beep and the shutters to click before leaving. Note Run will last approximately 2 hours.
17	When run is finished, proceed to “Analyzing Data” on page 5-13.

Analyzing Data

Overview When the installation run is complete, analyze the data so that system verification calculations can be performed.


CAUTION Before removing a TaqMan PCR tray from the ABI PRISM 7700 instrument, ensure that the sample block is not holding at 4 °C. Failure to do so may cause the optical caps to stick to the heated cover from condensation, and when the cover is moved back, the caps can come off and spill the fluorescent reaction solutions inside of the instrument. Do not insert a 4 °C hold after stage 3.

Analyzing the Run To analyze the run:

Step	Action
1	When the installation run is complete, save the run.
2	Select Analyze from the Analysis menu. The Amplification Plot window appears. If the Amplification Plot window does not appear, select Amplification Plot (⌘G) from the Analysis menu.
3	Set the threshold. Note The proper threshold is 0.15 – 0.25 There are two ways to do this: <ul style="list-style-type: none"> ◆ Either enter a value in the Use Threshold text box and click the Update Calculations button to reanalyze with the new threshold <p>or</p> <ul style="list-style-type: none"> ◆ Click and drag the horizontal drag bar (see figure below). An update automatically occurs when the bar is moved.

Horizontal drag bar

To analyze the run: (continued)

Step	Action
	<p>IMPORTANT The threshold needs to be adjusted to make the C_T values as accurate as possible. Adjust the threshold so that it is within the region of exponential growth (that is, where the growth curves are linear). Setting it too close to either the baseline or the plateau will cause high standard deviations and loss of discrimination between populations. This can be chosen from the Amplification Plot view of the 36 “5,00 template” population unknowns, at the point where the signals converge.</p>
4	<p>From the Viewer pop-up menu, select Ct vs Well Position, and determine if there are any outliers.</p> <p>Outliers are reactions with “count values” significantly different from others in the replicate set.</p> <div style="text-align: right; margin-right: 100px;">Viewer pop-up menu</div>  <p>IMPORTANT At least 30 of each set of 36 Unknowns must be used in the analysis for ABI PRISM 7700 instrument performance verification.</p>
5	<p>If outliers are present, take the following action:</p> <ol style="list-style-type: none"> Note the well position of each outlier. Click OK and return to the Setup view of the plate document. Highlight identified wells and select Not In Use from the Sample Type pop-up menu. Return to the Analysis view and reanalyze.
6	<p>After the analysis is complete, view the results in the Experiment Report from the Window menu.</p>
7	<p>Record the Average Copy Number and the Standard Deviation of each of the replicate populations.</p> <p>Note If the standard deviations are too high, there may be a problem with either the standard curve or the threshold setting.</p>

7700 Sequence Detection System Verification Calculations

Experimental Report Contents For each replicate population, the Experimental Report lists the:

- ◆ Copy number for each well
- ◆ Average copy number for each replicate set
- ◆ Standard deviation for each replicate set

Required Number of Reactions At least 30 of each set of 36 Unknown reactions must be used in the calculation to verify instrument performance.

Note Any position not used in the calculations must be checked during a subsequent PCR run to ensure that the thermal cycler, optics, and data collection are all functioning properly.

Validation Specification Calculation After determining that at least 30 of the Unknown reactions are successful, determine whether the run passed the validation specifications.

If the calculation below is true, then the instrument passes the validation specification (99.7% confidence level when determining the difference between the 5K and 10K populations).

$$[(\text{Copy.Unk.1}) - 3(\text{STDev.Unk.1})] > [(\text{Copy.Unk.2}) + 3(\text{STDev.Unk.2})]$$

Calculation Term	Value	Replicate Population
Copy.Unk.1	Average copy number of Unknown #1	10K
STDev.Unk.1	Standard deviation of Unknown #1	
Copy.Unk.2	Average copy number of Unknown #2	5K
STDev.Unk.2	Standard deviation of Unknown #2	

In simple terms, the calculation states that the average copy number of the 10K replicate population minus three standard deviations is greater than the average copy number of the 5K replicate population plus three standard deviations.

Calculating Validation Specification

Step	Action	Example										
1	Determine the values for each calculation term. Note You may choose to omit up to six reactions.	The calculation terms in this example are defined below. <table border="1"> <thead> <tr> <th>Equation Term</th> <th>Value</th> </tr> </thead> <tbody> <tr> <td>Copy.Unk.1</td> <td>9900</td> </tr> <tr> <td>STDev.Unk.1</td> <td>870</td> </tr> <tr> <td>Copy.Unk.2</td> <td>5100</td> </tr> <tr> <td>STDev.Unk.2</td> <td>520</td> </tr> </tbody> </table>	Equation Term	Value	Copy.Unk.1	9900	STDev.Unk.1	870	Copy.Unk.2	5100	STDev.Unk.2	520
Equation Term	Value											
Copy.Unk.1	9900											
STDev.Unk.1	870											
Copy.Unk.2	5100											
STDev.Unk.2	520											
2	Insert the values into the specification calculation: $[(\text{Copy.Unk.1}) - 3(\text{STDev.Unk.1})] > [(\text{Copy.Unk.2}) + 3(\text{STDev.Unk.2})]$	$9900 - 3(870) > 5100 + 3(520)$ $7200 > 6600$										
3	Determine whether the calculation is true or false. <table border="1"> <thead> <tr> <th>If it is...</th> <th>Then the 7700 instrument...</th> </tr> </thead> <tbody> <tr> <td>true</td> <td>passes the validation specification.</td> </tr> <tr> <td>false</td> <td>fails the validation specification.</td> </tr> </tbody> </table>	If it is...	Then the 7700 instrument...	true	passes the validation specification.	false	fails the validation specification.	In this example, the calculation is true; therefore the 7700 instrument SDS passes the validation specifications.				
If it is...	Then the 7700 instrument...											
true	passes the validation specification.											
false	fails the validation specification.											

IMPORTANT After the instrument passes the installation specifications, record the passing run information on the installation postcard (included in the packing kit) and send to Applied Biosystems in Foster City.

After Verifying Instrument Performance

After verifying the instrument performance, take the following action:

If...	Then...
performing a TaqMan card upgrade installation	proceed to Chapter 6, "TaqMan Card Upgrade Protocol."
performing an ABI PRISM 7700 instrument installation only	proceed to Chapter 7, "Customer Training."

TaqMan Card Upgrade Protocol

6

Overview

About This Chapter This chapter describes performing a TaqMan® Human Cytokine Card Upgrade, and analyzing card data and verifying the system performance. It also includes information on exporting data, verifying calculations, and troubleshooting.

In This Chapter For limited license information, please refer to the *TaqMan Human Cytokine Card Protocol* (P/N 4307577).

This chapter contains the following topics:

Topic	See Page
Section: TaqMan Card Upgrade Overview	6-2
Installation Procedure	6-3
Materials	6-4
Section: Performing a TaqMan Card Upgrade	6-5
Installing the Firmware and Software	6-6
Alignment Pin Replacement	6-9
Setting Up the Filling Station and Vacuum Pump Assembly	6-10
Preparing the PCR Reaction Mix and Filling the TaqMan Card	6-12
Loading the Card into the 7700 Instrument	6-19
Determining Exposure Time for Data Collection	6-22
Setting Up and Starting a TaqMan Card Run	6-25
Section: Analyzing Card Data, Exporting Results, and Troubleshooting	6-27
Analyzing the Card Data and Verifying System Performance	6-28
Exporting Results	6-31
TaqMan Card Verification Calculations	6-32
Troubleshooting	6-34

IMPORTANT Sequence Detection System software version 1.6.3 or newer **MUST BE INSTALLED** in order to perform the upgrade. This is because the 1.6.3 XILINX chip, which is not included in the upgrade package, is needed to run SDS software version 1.7.1.

Section: TaqMan Card Upgrade Overview

In This Section

Topic	See Page
Installation Procedure	6-3
Materials	6-4

Installation Procedure

Installation Process The following table provides a description of the install process:

Step	Action	Approximate time	Refer to...
1	Extract Background and Spectral Components	30 minutes	Chapter 4, "Performing a Spectral Calibration."
2	Verify the 7700 instrument performance (Use either RNase P Plate or β -actin Plate for this purpose)	For RNase P: ◆ 5 minutes setup ◆ 2 hours run time	"Preparing and Running the RNase P Instrument Verification Plate" on page 5-9.
		For β -actin: ◆ 1hr setup ◆ 2 hours run time	"Preparing the β -Actin Installation Plate" on page 5-3.
	Ensure that the instrument passes install specifications	20 minutes for analysis	"Analyzing Data" on page 5-13.
3	Install the card upgrade Hardware and Software	20 minutes	"Installing the Firmware and Software" on page 6-6.
4	Replace alignment pins (if necessary)	5 minutes	"Alignment Pin Replacement" on page 6-9.
5	Set up the ABI PRISM® Filling Station and vacuum pump assembly	20 minutes during 7700 instrument run	"Setting Up the Filling Station and Vacuum Pump Assembly" on page 6-10.
6	Prepare reagents and fill the card	10 minutes	"Preparing the PCR Reaction Mix and Filling the TaqMan Card" on page 6-12.
7	Load the TaqMan card into the 7700 instrument	2 minutes	"Loading the Card into the 7700 Instrument" on page 6-19.
8	Verify the 7700 instrument Exposure time for the card data collection	2 minutes	"Determining Exposure Time for Data Collection" on page 6-22.
9	Run the card and analyze to ensure that the system passes install specifications	1 hour 40 minutes for ABI PRISM 7700 instrument run	"Setting Up and Starting a TaqMan Card Run" on page 6-25.
		20 minutes for analysis	"Analyzing the Card Data and Verifying System Performance" on page 6-28.
10	Train customer	30 minutes	"Customer Training for the TaqMan Card Upgrade" on page 7-3.
Total Time		7-8 hours	

Materials

Verify Material Verify the TaqMan Card Upgrade Kit (P/N 4311899) contents against the packing lists. For the kit contents, see “TaqMan Card Upgrade Kit” on page B-4. Ensure That the following materials are on-site.

Instruments	Recommended Source
ABI PRISM® 7700 Sequence Detection System	Applied Biosystems
Centrifuge with 96-well plate adapter	MLS ^a
Microcentrifuge	MLS
Vacuum pump, oil-based, that can pull a vacuum down to 25 mTorr	Welch Two-Stage Belt-Drive Vacuum Pump (P/N 1400B-01)
Vacuum trap, Kontes®	VWR Catalog (P/N KT9266300-0021)
Materials	Recommended Source
MicroAmp® Optical 96-Well Reaction Plate/Optical Caps	Applied Biosystems (P/N 403012)
Gloves, disposable, powder-free	MLS
Microcentrifuge tubes, sterile 1.5-mL	MLS
Pipettors, positive-displacement or air-displacement	MLS
Pipette tips, aerosol resistant	MLS
Polypropylene tubes	MLS
Water, RNase-free, distilled, deionized	MLS

a. MLS is Major Lab Suppliers.

Section: Performing a TaqMan Card Upgrade

In This Section

Topic	See Page
Installing the Firmware and Software	6-6
Alignment Pin Replacement	6-9
Setting Up the Filling Station and Vacuum Pump Assembly	6-10
Preparing the PCR Reaction Mix and Filling the TaqMan Card	6-12
Loading the Card into the 7700 Instrument	6-19
Determining Exposure Time for Data Collection	6-22
Setting Up and Starting a TaqMan Card Run	6-25
Analyzing the Card Data and Verifying System Performance	6-28
Exporting Results	6-31
TaqMan Card Verification Calculations	6-32
Troubleshooting	6-34

Installing the Firmware and Software

Before Starting the Installation

IMPORTANT Sequence Detection System software version 1.6.3 must be currently installed before performing the following procedure. This is because the 1.6.3 XILINX chip, which is not included in the TaqMan card upgrade package, is needed to run SDS software version 1.7.1.

Electronic Handling Procedures

CAUTION Ground yourself before touching any components. Failure to do so could damage electrical components.

Replacing Thermal Cycler EPROMS

Step	Action
1	Turn off the 7700 instrument.
2	Replace the thermal cycler EPROMS (circuit reference U4 and U5 on the 9600 CPU PCB) with the parts labeled "B8 U5 4EA2" (P/N 4311885) and "B8 U4 419A" (P/N 4311884). IMPORTANT Please note the orientation of the notch on each EPROM when placing it on the 9600 CPU PCB.
3	Turn on the 7700 instrument and wait for the ready light to appear.
4	Launch the current 7700 software, open a plate document, click Show Analysis and check that sample and cover temperatures of the thermal cycler are being displayed by the Sequence Detection System. This confirms that there is a communication link between the Macintosh® computer and the thermal cycler.
5	Quit the 7700 software and turn off the 7700 instrument.

Setting Dip Switches and Replacing 7700 EPROM

To set the dip switches and replace the 7700 EPROM:

Step	Action																		
1	Set the dip switches (circuit reference SW1) on the 7700 PCA Controller Board (P/N 683324) to the following configuration in order to temporarily disable firmware downloading. <table border="1" data-bbox="537 1304 945 1663"> <thead> <tr> <th>Dip Switch</th> <th>Set to</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>ON</td> </tr> <tr> <td>2</td> <td>OFF</td> </tr> <tr> <td>3</td> <td>OFF</td> </tr> <tr> <td>4</td> <td>OFF</td> </tr> <tr> <td>5</td> <td>OFF</td> </tr> <tr> <td>6</td> <td>OFF</td> </tr> <tr> <td>7</td> <td>OFF</td> </tr> <tr> <td>8</td> <td>OFF</td> </tr> </tbody> </table>	Dip Switch	Set to	1	ON	2	OFF	3	OFF	4	OFF	5	OFF	6	OFF	7	OFF	8	OFF
Dip Switch	Set to																		
1	ON																		
2	OFF																		
3	OFF																		
4	OFF																		
5	OFF																		
6	OFF																		
7	OFF																		
8	OFF																		
	Note The first five (5) 7700 instrument's were built with a MUX that requires dip switch #2 to be in the ON position. These instruments have a note on the top of the MUX that indicates the proper dip switch position.																		

To set the dip switches and replace the 7700 EPROM: *(continued)*

Step	Action
2	Replace the firmware EPROM (circuit reference U41) on the 7700 instrument's PCA Controller Board (P/N 683324) with the new part labeled "M24 U41" (P/N 4311878). Note Sequence Detection System software version 1.6.3 must be currently installed, because the 1.6.3 XILINX chip is needed to run SDS software version 1.7.1.
3	Turn on the 7700 instrument and wait for the ready light to appear.
4	Turn off the 7700 instrument.

Setting Dip Switches To Re-Enable Firmware Downloading

Step	Action																		
1	Set the dip switches as follows to re-enable firmware downloading: <table border="1" data-bbox="581 705 989 1066"> <thead> <tr> <th>Dip Switch</th> <th>Set to</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>OFF</td> </tr> <tr> <td>2</td> <td>OFF</td> </tr> <tr> <td>3</td> <td>OFF</td> </tr> <tr> <td>4</td> <td>OFF</td> </tr> <tr> <td>5</td> <td>OFF</td> </tr> <tr> <td>6</td> <td>OFF</td> </tr> <tr> <td>7</td> <td>OFF</td> </tr> <tr> <td>8</td> <td>OFF</td> </tr> </tbody> </table> Note The first five (5) 7700 instrument's were built with a MUX that requires dip switch #2 to be in the on position. These instruments have a note on the top of the MUX that indicates the proper dip switch position.	Dip Switch	Set to	1	OFF	2	OFF	3	OFF	4	OFF	5	OFF	6	OFF	7	OFF	8	OFF
Dip Switch	Set to																		
1	OFF																		
2	OFF																		
3	OFF																		
4	OFF																		
5	OFF																		
6	OFF																		
7	OFF																		
8	OFF																		
2	Turn on the 7700 instrument.																		

Installing New Sequence Detection System Software

To install the new Sequence Detection System software:

Step	Action
1	Delete all Sequence Detection System software on the Power Macintosh® computer. Verify that all of the application programs have been removed by trying to open a SDS software run document.
2	Install the new Sequence Detection System software v.1.7.1 from the CD-ROM.
3	Launch the new Sequence Detection System software.
4	Open a plate document, click Show Analysis and check that sample and cover temperatures of the thermal cyclers are being displayed. This confirms that there is a communication link between the computer and the thermal cyclers.

To install the new Sequence Detection System software: *(continued)*

Step	Action
5	Verify that the EPROM version is listed as 24.0.b24. a. Select Diagnostics from the Instrument menu, and 7700 Instrument Verification from the submenu. The Instrument Tests dialog box appears. b. Click Run Tests . c. When the tests are completed, verify that the EPROM version listed is 24.0.b24 and click Done . Note This test may show a false failure for the shutters. Listen for the sound of the shutters opening and closing during the test to verify that they are working.
6	From the File menu, select Quit .
7	Install the Relative Quantification software into the SDS Applications folder on the hard drive.

Alignment Pin Replacement


Why Replace Alignment Pins?

The alignment pins of some older ABI PRISM 7700 Sequence Detection Systems must be replaced to run TaqMan® Human Cytokine cards. The alignment pins of older ABI PRISM 7700 instruments can interfere with the fit of the TaqMan card sandwich (lens plate+TaqMan card+reflective plate, see “About the Card Design” on page 6-12). The older pins contain a lip that prevents the lens plate from aligning with row A of the TaqMan card. This misalignment may cause amplification errors for wells A1–A12, and may also cause the cards to leak during cycling.


What Pins to Replace

Replace the alignment pins if they appear as in the following diagram:

Old alignment pin



New alignment pin



Note These pins are included in the TaqMan card Installation Kit, with the card adapter (P/N 4311059).


Replacing the Alignment Pins

Step	Action
1	Remove the alignment pins by unscrewing them counterclockwise.
2	Screw the new alignment pins into the sample block.
3	Tighten the pins as needed.

Setting Up the Filling Station and Vacuum Pump Assembly

Introduction Samples are loaded into TaqMan cards using a specialized tool called the ABI PRISM® Card Filling Station. The station is the focal point of the vacuum-assisted loading. The station works in combination with a vacuum pump, a vacuum trap, and a gauge to establish the vacuum necessary for filling cards.

The ABI PRISM Card Filling Station is shown in the figure below.



Setting Up Vacuum Pump Assembly

To set up the vacuum pump, gauge, and hoses:

Step	Action
1	Ensure that all fittings are tight on the vacuum tube/gauge assembly. Also ensure that the teflon thread seal tape has been applied to all threads.
2	If vacuum pump is new, fill with vacuum oil to specified level on pump, and turn the pump on and let it run for a few minutes to lubricate parts. Note Watch fill level to verify that there is enough oil. Turn pump off and add more if necessary.
3	Check that pump is pulling a vacuum by disconnecting the hose from the fill station.
4	Turn vacuum pump off.
5	Fit the open end of the hose onto the intake hose fitting on the vacuum pump. Note You may need to remove the section of large hose (and adapter), depending on the size of the intake hose fitting on the specific vacuum pump that the customer has supplied.

To set up the vacuum pump, gauge, and hoses: *(continued)*

Step	Action
6	<p>If a vacuum trap is available, cut the hose between the gauge and the quick connect fitting as shown above, and attach the vacuum trap.</p> <p>The vacuum trap should be oriented so that the higher section is toward the quick connect fitting.</p> <p>IMPORTANT It is okay to proceed with installation if the customer has not supplied a vacuum trap, but customer should be warned that damage to the filling station, vacuum pump and/or the vacuum gauge could occur without the use of a vacuum trap.</p>
7	<p>Verify that there are no leaks in the system.</p> <ol style="list-style-type: none"><li data-bbox="586 590 1149 617">a. Disconnect the hose assembly from the fill station.<li data-bbox="586 627 1442 684">b. Turn the vacuum pump and gauge on, and verify that there are no leaks in the system. <p>Note With the hose assembly disconnected, the gauge reading should go down to at least 300 mTorr within a few minutes. When the vacuum hose assembly is connected to the fill station, and a card is in the fill station, the gauge should read lower than 600 mTorr.</p>
8	<p>If necessary, tighten fittings and hose connections to prevent leaks.</p>
9	<p>After verifying vacuum system performance, release the vacuum in the assembly by pressing the quick disconnect to allow air into the hoses.</p> <p>IMPORTANT Do not leave the assembly under vacuum for extended periods of time without allowing it to return to atmospheric pressure. This can cause vacuum oil to be pulled into the hoses and gauge, which can ruin components and contaminate the system.</p>

Preparing the PCR Reaction Mix and Filling the TaqMan Card

Preparing PCR Reaction Mix


Note This protocol is optimized for TaqMan Universal PCR Master Mix.

To prepare the PCR reaction mix:

Step	Action
1	Remove the 18S and Human cDNA from the freezer to allow it to thaw.
2	Prepare the following mixture in a microcentrifuge tube: <ul style="list-style-type: none"> ◆ 150 μL 2X TaqMan PCR Universal Master Mix ◆ 30 μL 20X 18S Primer and TaqMan Probe Mix ◆ 2 μL Human cDNA ◆ 118 μL water Total volume = 300 μ L
3	Cap the microcentrifuge tube, and mix the solution thoroughly.
4	Use a table top centrifuge to spin the tube briefly to eliminate air bubbles from the mixture. The control sample can now be loaded into a TaqMan card and run on an ABI PRISM 7700 instrument.

About the Card Design

The following shows an exploded view of a TaqMan card to illustrate the important components:



Number	Component	Description
1	Reaction card	Acts as the vessel for the PCR
2	Adhesive flap	Used to seal the reaction card after it has been filled with sample and master mix
3	Adhesive backing	
4	Fill hole	Connects to the fill consumable
5	Alignment dimples	Guides the attachment of the fill consumable
6	Fill port	Connects to the card fill hole
7	Alignment pins	Guides the attachment of the fill consumable

Number	Component	Description
8	Fill reservoir	The reservoir for the cDNA sample
9	Fill consumable	A disposable component that channels the fluid from the fill port into the reaction card
10	Alignment holes	Aid in aligning the card within the ABI PRISM Filling Station

Guidelines for Loading TaqMan Cards

Follow the guidelines below to ensure proper filling of the card.

- ◆ Do not remove a TaqMan card from its packaging until you are ready to load it with reaction mix.


Excessive exposure to light damages the fluorescent probes.

- ◆ Do not twist or bend the soft fill consumable.



The seal between the reaction card and the fill consumable is crucial to the loading procedure. If broken, the seal may leak and result in an inadequately filled card.

Preparing the Card


To prepare the TaqMan card for filling:

Step	Action
1	Remove a TaqMan card and fill consumable from the refrigerator, and allow it to warm to room temperature for 5-10 minutes. Note Leave card in its protective bag while it warms up.
2	Remove card from the sealed bag and assemble the fill consumable and card: a. Remove the adhesive backing from the fill consumable. b. Align the fill port and two pins to the holes in the card and press them together (see below). c. Make sure the fill consumable is flush against the card d. Realign the consumable if necessary. 

To prepare the TaqMan card for filling: *(continued)*



Step	Action
3	<p>Carefully load the card into the ABI PRISM Card Filling Station.</p> <p>Orient the card so that the pins on the station align with the holes in the soft fill consumable as indicated in the figure below. Once the pins are correctly aligned, press down firmly on the top of the fill consumable to ensure a good fit.</p> <p>IMPORTANT Do not press down on the junction between the fill consumable and cytokine card.</p> <p>IMPORTANT Do not twist or bend the soft fill consumable attached to the TaqMan cards.</p> 
4	<p>Fold the adhesive flap on the card assembly backward onto itself, so that it will not interfere with the fill port when the fill station is closed.</p>
5	<p>Close the filling station lid, pressing firmly on top plate to ensure that the fill station is closed completely.</p> <p>IMPORTANT Make sure that the adhesive flap is folded backwards, so that it will not interfere with the fill hole on the card.</p> 

To prepare the TaqMan card for filling: *(continued)*




Step	Action
6	<p>Attach the vacuum hose to the diaphragm on the filling station lid, if it is not already attached.</p> <p>The end of the vacuum hose contains a quick-release valve that “clicks” when locked into place.</p>  <p>The diagram illustrates the attachment of a vacuum hose to a filling station lid. A vacuum hose with a quick-release valve is being inserted into a port on the lid. An arrow points to the connection point. Labels 'Vacuum hose' and 'Diaphragm' are present. The diaphragm is a rectangular plate with a grid of small circular ports. The lid is shown in a perspective view, and the vacuum hose is being inserted into a port on the lid. The diagram is labeled 'GRT1510' in the bottom right corner.</p>
7	Turn on the vacuum pump.
8	<p>Allow the vacuum pump to evacuate the card until the digital gauge on the hose stabilizes at or below 600 mTorr.</p> <p>IMPORTANT Do not fill the card above 600 mTorr vacuum. Above that reading, the pump may not create a vacuum strong enough to adequately fill a card.</p>

Filling and Sealing the Card

To fill and seal the card for thermal cycling:

Step	Action
1	<p>After the vacuum reaches 600 mTorr, pipette 300 μL of the PCR reaction mix (see "Preparing PCR Reaction Mix" on page 6-12) into the fill port of the filling station.</p> 
2	<p>In one motion, firmly pull the filling station actuator to its maximum extension. The TaqMan card fills with sample.</p> <p>IMPORTANT Do not force the actuator into the closed position after filling is finished. It will retract automatically when the filling station lid is opened.</p> 
3	<p>Complete the following as quickly as possible:</p> <ol style="list-style-type: none"> Disconnect the vacuum hose from the filling station. Open the filling station lid. Turn off the vacuum pump. <p>IMPORTANT Remember to bleed the remaining vacuum in the hose when finished by briefly reconnecting the hose to the fill station.</p>
4	<p>Remove the card from the filling station.</p>

To fill and seal the card for thermal cycling: *(continued)*

Step	Action
5	<p data-bbox="537 275 1382 331">Detach and discard the fill consumable. Make sure to also remove the adhesive strip that attaches the fill consumable to the card</p> 
6	<p data-bbox="537 701 1187 730">Bend back the adhesive flap, and peel off the plastic backing.</p> 
7	<p data-bbox="537 1031 1365 1087">Fold the adhesive flap over the front edge of the card, making sure to align the holes in the flap to the wells of the card.</p> 
8	<p data-bbox="537 1360 1084 1390">Press firmly on the flap to ensure an adequate seal.</p> <p data-bbox="537 1409 1333 1465">The card is now filled and ready to load into the ABI PRISM 7700 Sequence Detection System.</p>


Loading the Card into the 7700 Instrument

Overview This section contains the following:

- ◆ Description of the ABI PRISM® Card Adapter
- ◆ Procedure for loading a TaqMan card into an ABI PRISM 7700 Sequence Detection System

ABI PRISM Card Adapter Design



Because of the unique properties of the TaqMan card design, an ABI PRISM Card Adapter is needed to run the card's on the 7700 instrument. The Card Adapter is a unique device that ensures that adequate heat transfer and fluorescent data collection occur during PCR. The components of the adapter are shown in the following figure:





Component	Description
Lens plate	Contains lenses that direct the focal point of the argon ion laser into the wells of the card.
TaqMan card	Contains a sample and the necessary reagents for the PCR.
Reflective plate	Ensures efficient conduction and heat transfer to all wells of the TaqMan card.

Loading the TaqMan Card

To load the card into the ABI PRISM 7700 instrument:

Step	Action
1	Slide the sample block cover back, exposing the sample block.
2	Place the reflective plate on top of the sample block so that the notched corner is located in the upper-right corner.
	
3	Place the filled TaqMan card on top of the reflective plate so that the keyed corner of the card aligns with the notched corner of the plate.
	

To load the card into the ABI PRISM 7700 instrument: *(continued)*

Step	Action
4	<p data-bbox="586 275 1422 331">Place the Lens plate on top of the card so that the holes in the plate fit over the alignment pins on the sample block.</p> 
5	<p data-bbox="586 974 1187 1003">Slide the cover over the sample block and tighten the lid.</p>  <p data-bbox="586 1268 1403 1297">IMPORTANT Tighten the cover firmly to ensure uniform heating of the card.</p>

Determining Exposure Time for Data Collection

Introduction The data collection exposure time can be set in the software and determines the amount of time that fluorescent data is collected for each data point. The proper setting is instrument-dependent and can be determined by performing a plate read.

IMPORTANT The exposure time **MUST BE DETERMINED** for each 7700 instrument separately before the first TaqMan card is run, in order to ensure that data is collected properly. This is because the correct exposure time is dependent on the sensitivity of the 7700 instrument and will vary from instrument to instrument.

IMPORTANT The exposure time is determined separately for plate runs and card runs.

Determining Exposure Time To determine if the default exposure time is correct for TaqMan card runs:


Step	Action								
1	Launch the Sequence Detection System version 1.7.1. The software displays a new plate document.								
2	From the File menu, select Close .								
3	From the File menu, select New Plate .								
4	Configure a new plate document with the following attributes: <table border="1" data-bbox="537 915 1302 1079"> <thead> <tr> <th>From Menu...</th> <th>Select...</th> </tr> </thead> <tbody> <tr> <td>Plate Type</td> <td>Single Reporter</td> </tr> <tr> <td>Plate Format</td> <td>The Card</td> </tr> <tr> <td>Run</td> <td>Plate Read</td> </tr> </tbody> </table>	From Menu...	Select...	Plate Type	Single Reporter	Plate Format	The Card	Run	Plate Read
From Menu...	Select...								
Plate Type	Single Reporter								
Plate Format	The Card								
Run	Plate Read								
5	Select all wells of the plate document.								
6	Select UNKN-Unknown from the Sample Type pop-up menu. The SDS software labels all selected wells as UNKN.								
7	Click the Show Analysis button.								



To determine if the default exposure time is correct for TaqMan card runs: *(continued)*

Step	Action
8	<p>Click the Pre-PCR Read button.</p> <p>The plate read will take a few seconds.</p>
9	<p>When the plate read is finished, select Raw Spectra from the Analysis menu to verify that the tops of peaks are visible for all 96 wells. Visible peaks indicate that the CCD camera is not being saturated.</p> <p>The raw data should look like the figure below.</p> <div data-bbox="586 520 1455 1020" data-label="Figure"> <p>The graph titled "Raw Spectra - Card Pre-read" displays fluorescent intensity on the y-axis (0 to 7000) against wavelength in nanometers on the x-axis (500 to 660). It contains numerous overlapping colored lines, each representing a different well. All lines show a similar broad peak centered around 550-600 nm, with the highest peak reaching approximately 6000 intensity units.</p> </div> <p>If the exposure time is too long, the signal will saturate the CCD camera and the tops of some peaks will be cut off like the figure below.</p> <div data-bbox="586 1142 1455 1642" data-label="Figure"> <p>The graph titled "Raw Spectra - CardPlateReadTooHighExposure" displays fluorescent intensity on the y-axis (0 to 40000) against wavelength in nanometers on the x-axis (500 to 660). It contains numerous overlapping colored lines. The peaks are much higher, reaching up to 30000-35000 intensity units. Many of the peaks are truncated at the top, indicating that the signal has saturated the CCD camera.</p> </div> <p>Note If the CCD camera is being saturated, proceed to the next step. If the raw data looks okay, proceed to "Setting Up and Starting a TaqMan Card Run" on page 6-25.</p>

To determine if the default exposure time is correct for TaqMan card runs: *(continued)*

Step	Action
10	<p>To adjust the card exposure time, select Diagnostics and then Advanced Options from the Instrument menu.</p> 
11	<p>Click on the Set 7700 Exposure Time checkbox as shown above, and change the exposure time in the for Card text box.</p> <p>Note The amount by which the exposure time should be reduced depends upon the severity of saturation and is a judgment call.</p>
12	<p>Repeat steps 2-9 to verify that the new exposure time setting is acceptable.</p>

Setting Up and Starting a TaqMan Card Run

Overview This section describes the procedures for creating a plate document from the template file and starting a TaqMan card run.


Where to Find the Template File The Sequence Detection System version 1.7.1 contains a template file configured specifically for TaqMan card runs. The template file is located in the Templates folder that is in the SDS 1.7.1 folder.

Template files are identical to plate documents, however they do not contain fluorescence data from a previous run.


Note If you do not find a Cytokine Card template file, refer to Appendix E, "Creating a TaqMan Card Template."

Creating Plate Documents From the Template File

To create a plate document from the Sequence Detection System template file:

Step	Action
1	Open the Sequence Detection System version 1.7.1. The software opens and displays a new plate document.
2	From the File menu, select Close .
3	From the File menu, select Open Plate . A directory dialog box appears.
4	Navigate to the location of the Cytokine Card template file: SDS 1.7.1 folder > Templates folder > Cytokine Card Template file Select the template file and click Open . The SDS software creates a plate document with attributes identical to that of the template file.
5	Click the Thermal Cycler Conditions button and verify that the conditions are as shown below. 
6	Click the Show Analysis button to toggle to the Analysis View.

To create a plate document from the Sequence Detection System template

Step	Action
7	<p>From the Instrument menu, select Diagnostics and then Advanced Options to verify that the settings are as shown below.</p> <p>IMPORTANT The card exposure time will vary depending upon the instrument and should be set according to procedure described in “Determining Exposure Time” on page 6-22.</p> 
8	<p>Click Run to begin thermal cycling.</p> <p>Note The ABI PRISM 7700 instrument may pause momentarily before initiating thermal cycling, to allow the heated cover to cool. Normal operating temperature for the heated cover is 60°C for a card run, due to the difference in thermokinetic properties of the TaqMan card from those of a MicroAmpOptical 96-Well Reaction Plate.</p>

Section: Analyzing Card Data, Exporting Results, and Troubleshooting

In This Section

Topic	See Page
Analyzing the Card Data and Verifying System Performance	6-28
Exporting Results	6-31
TaqMan Card Verification Calculations	6-32
Troubleshooting	6-34

Analyzing the Card Data and Verifying System Performance

General When the card installation run is complete, verify that the system is performing to specifications. Use the following information as a guide for data analysis.

Analyzing the Run Use the following as a guide when analyzing installation runs:

Step	Action																		
1	When the installation run is complete, save the run.																		
2	<p>Select Analyze from the Analysis menu.</p> <p>The Amplification Plot window appears.</p> <p>If the Amplification Plot window does not appear, select Amplification Plot from the Analysis menu.</p> <div data-bbox="542 701 1360 1306" data-label="Figure"> <table border="1" data-bbox="820 1102 1128 1239"> <thead> <tr> <th></th> <th>Ct</th> <th>Std Dev</th> </tr> </thead> <tbody> <tr> <td>FAM - A1</td> <td>35,000</td> <td>0.009</td> </tr> <tr> <td>FAM - A2</td> <td>35,000</td> <td>0.011</td> </tr> <tr> <td>FAM - A3</td> <td>35,000</td> <td>0.005</td> </tr> <tr> <td>FAM - A4</td> <td>35,000</td> <td>0.012</td> </tr> <tr> <td>FAM - A5</td> <td>35,000</td> <td>0.011</td> </tr> </tbody> </table> </div> <p>IMPORTANT The picture shown above is what you should expect to see when you first open the Amplification Plot window. For the Installation run, no amplification occurs in the FAM dye layer. Therefore, you will be concerned with the VIC dye layer only.</p>		Ct	Std Dev	FAM - A1	35,000	0.009	FAM - A2	35,000	0.011	FAM - A3	35,000	0.005	FAM - A4	35,000	0.012	FAM - A5	35,000	0.011
	Ct	Std Dev																	
FAM - A1	35,000	0.009																	
FAM - A2	35,000	0.011																	
FAM - A3	35,000	0.005																	
FAM - A4	35,000	0.012																	
FAM - A5	35,000	0.011																	
3	<p>From the Reporter pop-up menu, select VIC.</p> <p>You should now see a graph similar to the picture below.</p>																		

Use the following as a guide when analyzing installation runs: *(continued)*

Step	Action																		
4	<p>Set the threshold. There are two ways to do this:</p> <ul style="list-style-type: none"> ◆ Either enter a value in the Use Threshold text box and click the Update Calculations button to reanalyze with the new threshold <p>or</p> <ul style="list-style-type: none"> ◆ Click and drag the horizontal drag bar (see figure below). An update automatically occurs when the bar is moved. 																		
<p>Horizontal drag bar</p> <p>The screenshot shows the 'Amplification Plot' window. The plot title is 'Amplification - Card Install 10-29-99'. The y-axis is labeled 'ΔRn' and ranges from 10⁻² to 10¹. The x-axis is labeled 'Cycle' and ranges from 0 to 35. Multiple colored lines represent different samples (VIC-A1 to VIC-A9). A horizontal line with a drag bar is positioned at approximately 10⁻¹ on the y-axis. Below the plot is a 'Threshold Cycle Calculation' panel. It has a 'Threshold' section with 'Use Threshold' set to 0.12 and a 'Suggest' button. Below that is 'Mult. * Stddev' set to 10.0 and 'Omit Threshold' set to 2.0. There is also a 'Baseline' section with 'Start' set to 3 and 'Stop' set to 15. An 'Update Calculations' button is at the bottom of this panel. To the right of the panel is a table with columns 'VIC', 'Ct', and 'Std Dev'. The table contains data for VIC-A1 through VIC-A5. Below the table is a 'Please Set the Threshold Value on All Reporter Layers. Click OK to Continue.' message and an 'OK' button.</p> <table border="1" data-bbox="868 966 1169 1102"> <thead> <tr> <th>VIC</th> <th>Ct</th> <th>Std Dev</th> </tr> </thead> <tbody> <tr> <td>VIC - A1</td> <td>19.156</td> <td>0.011</td> </tr> <tr> <td>VIC - A2</td> <td>18.246</td> <td>0.009</td> </tr> <tr> <td>VIC - A3</td> <td>18.043</td> <td>0.009</td> </tr> <tr> <td>VIC - A4</td> <td>19.089</td> <td>0.011</td> </tr> <tr> <td>VIC - A5</td> <td>18.578</td> <td>0.012</td> </tr> </tbody> </table>		VIC	Ct	Std Dev	VIC - A1	19.156	0.011	VIC - A2	18.246	0.009	VIC - A3	18.043	0.009	VIC - A4	19.089	0.011	VIC - A5	18.578	0.012
VIC	Ct	Std Dev																	
VIC - A1	19.156	0.011																	
VIC - A2	18.246	0.009																	
VIC - A3	18.043	0.009																	
VIC - A4	19.089	0.011																	
VIC - A5	18.578	0.012																	
<p>IMPORTANT The threshold needs to be adjusted to make the C_T values as accurate as possible. Adjust the threshold so that it is within the region of exponential growth (= where the growth curves are linear). Setting it too close to either the baseline or the plateau will cause high standard deviations and loss of discrimination between populations.</p>																			

Use the following as a guide when analyzing installation runs: *(continued)*



Step	Action
5	<p data-bbox="537 275 1409 331">From the Viewer pop-up menu, select Ct vs Well Position, and determine if there are any outliers.</p> <p data-bbox="537 348 1369 380">Outliers are reactions with “count values” significantly different from the others.</p> <div data-bbox="537 411 1409 982" style="text-align: center;"> <p data-bbox="537 485 1409 982">The figure is a scatter plot titled "Amplification - Card install run 11-1-99". The vertical axis is labeled "Threshold Cycle (Ct)" and ranges from -5,000 to 45,000 in increments of 5,000. The horizontal axis is labeled "Well Position" and ranges from 0 to 100 in increments of 5. The plot shows a series of data points for 96 wells. Most points are clustered between 15,000 and 20,000 Ct, showing a slight wave pattern. Two points are significantly higher than the rest: one at well position 8 with a Ct of approximately 33,000, and another at well position 95 with a Ct of approximately 35,000. A bracket labeled "Outliers" spans these two points.</p> </div> <p data-bbox="537 1014 1409 1071">Note A slight wave pattern in “count values” over the 96 wells is characteristic of TaqMan card data and is considered to be normal.</p> <p data-bbox="537 1094 1349 1150">IMPORTANT At least 90 of the 96 reactions must be used for TaqMan card performance verification.</p>
6	<p data-bbox="537 1163 1040 1194">If outliers are present, take the following action:</p> <ol data-bbox="537 1199 1369 1367" style="list-style-type: none"> <li data-bbox="537 1199 959 1230">a. Note the well position of each outlier. <li data-bbox="537 1234 1211 1266">b. Click OK and return to the Setup view of the plate document. <li data-bbox="537 1270 1369 1327">c. Highlight identified wells and select Not In Use from the Sample Type pop-up menu. <li data-bbox="537 1331 1024 1362">d. Return to the Analysis view and reanalyze.
7	<p data-bbox="537 1379 1341 1409">After the analysis is complete, proceed to “Exporting Results” on page 6-31.</p>

Exporting Results

About Exporting Results

To analyze data from the TaqMan card, export the results to a data file. The Sequence Detection System version 1.7.1 can export raw data from a sequence detection run in a tab-delimited format compatible with most spreadsheet applications.

Exporting Run Results To a Data File

Step	Action
1	<p>Select Export from the File menu, and Results (⌘H) from the submenu.</p> 
2	<p>Enter a name in the Export result data as text box.</p> 
3	Click the All Wells radio button.
4	Navigate to the location on the hard drive you would like to store the file.
5	<p>Click Export.</p> <p>The Sequence Detection System exports the data. The 7700 data file icon is shown below.</p> 
6	Quit the Sequence Detection System.

TaqMan Card Verification Calculations

Installation Specifications There are two specifications that need to be met in order to pass the TaqMan card upgrade installation. Both of them involve the Ct values for the VIC dye layer only.

Passing Criteria for the VIC dye layer:

- ◆ Average Ct \leq 20
- ◆ Standard Deviation for Ct's \leq 1

Note For each set of 96 reactions, at least 90 must be used in the system performance verification calculations.

Results File The exported results file will list the Ct's for both the FAM dye layer and the VIC dye layer in column F.

IMPORTANT For the specification calculations, use data from the VIC dye layer only.

Validation Specification Calculation You will need Microsoft Excel to perform the following specification calculations. To verify the TaqMan card upgrade installation:

Step	Action																																																																																																																																																																																																	
1	<p>Navigate to the location on the hard drive containing the card run exported results file, and double click on the icon to open the file.</p> <p>This should automatically launch Microsoft Excel. You will see data similar to the following:</p> <table border="1"> <thead> <tr> <th colspan="12">Card Install Results 10-29-99</th> </tr> <tr> <th></th> <th>A</th> <th>B</th> <th>C</th> <th>D</th> <th>E</th> <th>F</th> <th>G</th> <th>H</th> <th>I</th> <th>J</th> <th></th> </tr> <tr> <th>1</th> <th>Well</th> <th>Reporter</th> <th>Type</th> <th>Baseline</th> <th>StdDev</th> <th>deltaRn</th> <th>Ct</th> <th>Quantity</th> <th>Replicate</th> <th>Qty Mean</th> <th>Qty StdDev</th> <th>IF</th> </tr> </thead> <tbody> <tr><td>2</td><td>1</td><td>FAM1</td><td>UNKN</td><td>3.35E-03</td><td>2.15E-02</td><td></td><td>35</td><td>-1.00E+00</td><td></td><td>0.00E+00</td><td>0.00E+00</td><td>-</td></tr> <tr><td>3</td><td>2</td><td>FAM1</td><td>UNKN</td><td>8.04E-03</td><td>2.95E-02</td><td></td><td>35</td><td>-1.00E+00</td><td></td><td>0.00E+00</td><td>0.00E+00</td><td>-</td></tr> <tr><td>4</td><td>3</td><td>FAM1</td><td>UNKN</td><td>7.22E-03</td><td>-1.97E-02</td><td></td><td>35</td><td>-1.00E+00</td><td></td><td>0.00E+00</td><td>0.00E+00</td><td>-</td></tr> <tr><td>5</td><td>4</td><td>FAM1</td><td>UNKN</td><td>9.73E-03</td><td>-2.00E-02</td><td></td><td>35</td><td>-1.00E+00</td><td></td><td>0.00E+00</td><td>0.00E+00</td><td>-</td></tr> <tr><td>6</td><td>5</td><td>FAM1</td><td>UNKN</td><td>8.50E-03</td><td>2.34E-02</td><td></td><td>35</td><td>-1.00E+00</td><td></td><td>0.00E+00</td><td>0.00E+00</td><td>-</td></tr> <tr><td>7</td><td>6</td><td>FAM1</td><td>UNKN</td><td>7.52E-03</td><td>1.65E-02</td><td></td><td>35</td><td>-1.00E+00</td><td></td><td>0.00E+00</td><td>0.00E+00</td><td>-</td></tr> <tr><td>8</td><td>7</td><td>FAM1</td><td>UNKN</td><td>7.42E-03</td><td>-1.43E-03</td><td></td><td>35</td><td>-1.00E+00</td><td></td><td>0.00E+00</td><td>0.00E+00</td><td>-</td></tr> <tr><td>9</td><td>8</td><td>FAM1</td><td>UNKN</td><td>2.51E-03</td><td>2.76E-03</td><td></td><td>35</td><td>-1.00E+00</td><td></td><td>0.00E+00</td><td>0.00E+00</td><td>-</td></tr> <tr><td>10</td><td>9</td><td>FAM1</td><td>UNKN</td><td>5.47E-03</td><td>8.63E-04</td><td></td><td>35</td><td>-1.00E+00</td><td></td><td>0.00E+00</td><td>0.00E+00</td><td>-</td></tr> <tr><td>11</td><td>10</td><td>FAM1</td><td>UNKN</td><td>5.51E-03</td><td>-3.86E-03</td><td></td><td>35</td><td>-1.00E+00</td><td></td><td>0.00E+00</td><td>0.00E+00</td><td>-</td></tr> <tr><td>12</td><td>11</td><td>FAM1</td><td>UNKN</td><td>7.68E-03</td><td>1.46E-02</td><td></td><td>35</td><td>-1.00E+00</td><td></td><td>0.00E+00</td><td>0.00E+00</td><td>-</td></tr> <tr><td>13</td><td>12</td><td>FAM1</td><td>UNKN</td><td>7.40E-03</td><td>-1.85E-02</td><td></td><td>35</td><td>-1.00E+00</td><td></td><td>0.00E+00</td><td>0.00E+00</td><td>-</td></tr> </tbody> </table> <p>Note The specification calculations are based on the 96-wells of Ct values for the VIC dye layer, which can be found in column F, Rows 112 to 207.</p>	Card Install Results 10-29-99													A	B	C	D	E	F	G	H	I	J		1	Well	Reporter	Type	Baseline	StdDev	deltaRn	Ct	Quantity	Replicate	Qty Mean	Qty StdDev	IF	2	1	FAM1	UNKN	3.35E-03	2.15E-02		35	-1.00E+00		0.00E+00	0.00E+00	-	3	2	FAM1	UNKN	8.04E-03	2.95E-02		35	-1.00E+00		0.00E+00	0.00E+00	-	4	3	FAM1	UNKN	7.22E-03	-1.97E-02		35	-1.00E+00		0.00E+00	0.00E+00	-	5	4	FAM1	UNKN	9.73E-03	-2.00E-02		35	-1.00E+00		0.00E+00	0.00E+00	-	6	5	FAM1	UNKN	8.50E-03	2.34E-02		35	-1.00E+00		0.00E+00	0.00E+00	-	7	6	FAM1	UNKN	7.52E-03	1.65E-02		35	-1.00E+00		0.00E+00	0.00E+00	-	8	7	FAM1	UNKN	7.42E-03	-1.43E-03		35	-1.00E+00		0.00E+00	0.00E+00	-	9	8	FAM1	UNKN	2.51E-03	2.76E-03		35	-1.00E+00		0.00E+00	0.00E+00	-	10	9	FAM1	UNKN	5.47E-03	8.63E-04		35	-1.00E+00		0.00E+00	0.00E+00	-	11	10	FAM1	UNKN	5.51E-03	-3.86E-03		35	-1.00E+00		0.00E+00	0.00E+00	-	12	11	FAM1	UNKN	7.68E-03	1.46E-02		35	-1.00E+00		0.00E+00	0.00E+00	-	13	12	FAM1	UNKN	7.40E-03	-1.85E-02		35	-1.00E+00		0.00E+00	0.00E+00	-
Card Install Results 10-29-99																																																																																																																																																																																																		
	A	B	C	D	E	F	G	H	I	J																																																																																																																																																																																								
1	Well	Reporter	Type	Baseline	StdDev	deltaRn	Ct	Quantity	Replicate	Qty Mean	Qty StdDev	IF																																																																																																																																																																																						
2	1	FAM1	UNKN	3.35E-03	2.15E-02		35	-1.00E+00		0.00E+00	0.00E+00	-																																																																																																																																																																																						
3	2	FAM1	UNKN	8.04E-03	2.95E-02		35	-1.00E+00		0.00E+00	0.00E+00	-																																																																																																																																																																																						
4	3	FAM1	UNKN	7.22E-03	-1.97E-02		35	-1.00E+00		0.00E+00	0.00E+00	-																																																																																																																																																																																						
5	4	FAM1	UNKN	9.73E-03	-2.00E-02		35	-1.00E+00		0.00E+00	0.00E+00	-																																																																																																																																																																																						
6	5	FAM1	UNKN	8.50E-03	2.34E-02		35	-1.00E+00		0.00E+00	0.00E+00	-																																																																																																																																																																																						
7	6	FAM1	UNKN	7.52E-03	1.65E-02		35	-1.00E+00		0.00E+00	0.00E+00	-																																																																																																																																																																																						
8	7	FAM1	UNKN	7.42E-03	-1.43E-03		35	-1.00E+00		0.00E+00	0.00E+00	-																																																																																																																																																																																						
9	8	FAM1	UNKN	2.51E-03	2.76E-03		35	-1.00E+00		0.00E+00	0.00E+00	-																																																																																																																																																																																						
10	9	FAM1	UNKN	5.47E-03	8.63E-04		35	-1.00E+00		0.00E+00	0.00E+00	-																																																																																																																																																																																						
11	10	FAM1	UNKN	5.51E-03	-3.86E-03		35	-1.00E+00		0.00E+00	0.00E+00	-																																																																																																																																																																																						
12	11	FAM1	UNKN	7.68E-03	1.46E-02		35	-1.00E+00		0.00E+00	0.00E+00	-																																																																																																																																																																																						
13	12	FAM1	UNKN	7.40E-03	-1.85E-02		35	-1.00E+00		0.00E+00	0.00E+00	-																																																																																																																																																																																						
2	Scroll down to the second set of data, which is for the VIC layer (Rows 112 to 207).																																																																																																																																																																																																	
3	<p>Click on any empty cell to select it, and enter the following formula EXACTLY as shown below.</p> <p>= Average (F112 : F207)</p>																																																																																																																																																																																																	
4	<p>Press Enter on the keyboard and the average value of the 96 boxes from F112 to F207 is displayed in the selected cell.</p> <p>This is the average Ct value for the VIC dye layer.</p>																																																																																																																																																																																																	
5	<p>Click on another empty cell to select it, and enter the following formula EXACTLY as shown below.</p> <p>=stdev (F112 : F207)</p>																																																																																																																																																																																																	

To verify the TaqMan card upgrade installation: *(continued)*

Step	Action
6	Press Enter on the keyboard and the standard deviation for the Ct values in the 96 boxes from F112 to F207 is displayed in the selected cell. This is the standard deviation of Ct's for the VIC dye layer.
7	Record these two values, save the document and quit Microsoft Excel.

Troubleshooting

Troubleshooting This section describes troubleshooting for three common tasks:

- Topics**
- ◆ Loading the card
 - ◆ Analyzing data
 - ◆ Interpreting results
-

Loading the Card

Troubleshooting Tips for Loading the Card

Observation	Possible Cause(s)	Recommended Action
Vacuum is not reaching the proper level (600 microns or below)	<ul style="list-style-type: none">◆ Loose connection in the plumbing between the pump and the fill station◆ Battery in the vacuum gauge is low◆ Fill fixture is not closed properly	<ol style="list-style-type: none">a. Check all connections.b. Change the battery in the vacuum gauge.c. Check the fill fixture to ensure it is fully closed.d. Adjust/align card assembly in fill fixture.
The sample does not completely enter the card when the plunger is pulled	<ul style="list-style-type: none">◆ The card and fill consumable are not aligned properly inside the fill station◆ Adhesive flap is covering the fill hole on the card◆ The fill station is out of adjustment	<ol style="list-style-type: none">a. Using a pipette, remove the sample from the reservoir, and store it temporarily.b. Unlock and open the fill station.c. Reposition the card and fill consumable and load the card again.d. Fold the adhesive flap on the card backwards and out of the way of the card fill hole.e. If the failure persists, the problem may be the fill station.

Analyzing Data

Troubleshooting Tips for Data Analysis

Observation	Possible Cause	Recommended Action							
<p>All well signals slowly degrade to background by the end of the run (as viewed in the Multi-component view)</p>	<ul style="list-style-type: none"> ◆ Card does not contain fluid ◆ The heated cover was not pressed down completely ◆ The rubber portion of the lens plate sealing the card fill port is damaged ◆ Instrument malfunction 	<p>Check the card for fluid.</p> <table border="1" data-bbox="870 411 1461 772"> <thead> <tr> <th data-bbox="870 411 1101 478">Does the card contain fluid?</th> <th data-bbox="1101 411 1461 478">Then...</th> </tr> </thead> <tbody> <tr> <td data-bbox="870 478 1101 709">No</td> <td data-bbox="1101 478 1461 709"> <p>Remember to press down heated cover completely for future runs.</p> <p>If the problem persists, replace the rubber on the top plate of the card adapter with the part included with the kit.</p> </td> </tr> <tr> <td data-bbox="870 709 1101 772">Yes</td> <td data-bbox="1101 709 1461 772"> <p>The error may be an instrument malfunction.</p> </td> </tr> </tbody> </table>		Does the card contain fluid?	Then...	No	<p>Remember to press down heated cover completely for future runs.</p> <p>If the problem persists, replace the rubber on the top plate of the card adapter with the part included with the kit.</p>	Yes	<p>The error may be an instrument malfunction.</p>
Does the card contain fluid?	Then...								
No	<p>Remember to press down heated cover completely for future runs.</p> <p>If the problem persists, replace the rubber on the top plate of the card adapter with the part included with the kit.</p>								
Yes	<p>The error may be an instrument malfunction.</p>								
<p>Growth curves slope down toward the later cycles of large growth curves</p>	<p>The CCD camera is becoming saturated.</p> <p>Sometimes the combination of a new 7700 instrument and a new card adapter can cause the CCD camera to become saturated for a few wells.</p>	<p>Correcting the problem is optional. The saturation will subside with use and does not affect the C_T calculations.</p> <p>To correct the problem for subsequent runs (the problem run cannot be corrected for the current run):</p> <ol style="list-style-type: none"> a. Select Diagnostics from the Instrument menu, and Advanced Options from the submenu. The Advanced Options dialog box appears. b. In the Miscellaneous section, select the checkbox labeled Set 7700 Exposure Time for Plate. c. Enter the appropriate value as determined using the procedure in "Determining Exposure Time" on page 6-22 in the text box, and click OK. <p>IMPORTANT Reset the exposure time before the next plate run.</p>							

Troubleshooting Tips for Data Analysis *(continued)*

Observation	Possible Cause	Recommended Action						
No amplification in the wells of the top row (wells A1–A12)	The alignment pins on the 7700 instrument sample block may be interfering with heat transfer	Replace the alignment pins.						
Growth curves slope up at the later cycles and/or the earlier cycles correlating with a loss in precision of the data	Plate file is in Standard Plate format	Check the top right corner of the plate file in the analysis view for the plate file type.						
		<table border="1"> <thead> <tr> <th>Plate File Type</th> <th>Then...</th> </tr> </thead> <tbody> <tr> <td>Standard Plate</td> <td>The problem run cannot be corrected after the run is complete. Rerun the sample.</td> </tr> <tr> <td>Card</td> <td>The following may be the cause: <ul style="list-style-type: none"> ◆ A faulty heated cover sensor, if the temperature of the heated cover is 95 °C or greater. ◆ A faulty card, possibly air leaks. </td> </tr> </tbody> </table>	Plate File Type	Then...	Standard Plate	The problem run cannot be corrected after the run is complete. Rerun the sample.	Card	The following may be the cause: <ul style="list-style-type: none"> ◆ A faulty heated cover sensor, if the temperature of the heated cover is 95 °C or greater. ◆ A faulty card, possibly air leaks.
		Plate File Type	Then...					
Standard Plate	The problem run cannot be corrected after the run is complete. Rerun the sample.							
Card	The following may be the cause: <ul style="list-style-type: none"> ◆ A faulty heated cover sensor, if the temperature of the heated cover is 95 °C or greater. ◆ A faulty card, possibly air leaks. 							

Interpreting Results

Troubleshooting Tips for Interpreting the Results

Observation	Possible Cause	Recommended Action
Average endogenous control C _T s are over 23	<ul style="list-style-type: none"> ◆ Not enough sample in the card ◆ Not enough cDNA from total RNA ◆ cDNA was reverse transcribed using poly-T primers 	<ul style="list-style-type: none"> ◆ Load more sample into the card for the next run. ◆ Ensure that the cDNA is reverse transcribed from random hexamers and unfiltered.

Customer Training

7

Overview

About This Chapter This chapter describes customer training for the ABI PRISM® 7700 Sequence Detector and training for the TaqMan® Human Cytokine Card Upgrade. It also includes information on verifying the ABI PRISM 7700 instrument and the TaqMan card system performance.

In This Chapter This chapter contains the following topics:

Topic	See Page
Customer Training for the 7700 Sequence Detection System	7-2
Customer Training for the TaqMan Card Upgrade	7-3
Completing the Installation	7-4

Customer Training for the 7700 Sequence Detection System

Basic Training The customer training checklist on page D-5 should be used as a guideline for training customers in using the 7700 Sequence Detection System.

Discuss the following with the customer:

- ◆ Instrument theory of operation (optics and thermal cycler)
- ◆ TaqMan PCR theory (if needed)
- ◆ Sample preparation
- ◆ Description of laser safety features

Instrument Operation Train the customer on the following subjects:

- ◆ Setting up runs (Sample Type info, Document set up, and thermal cycler conditions)
- ◆ Analyzing data (the applications specialist should participate if possible)
- ◆ Using the Sequence Detection System software, including how to display graphical results and pointing out the importance of the key menus (Preferences, Advanced Options, etc.)

Instrument Maintenance Discuss these items with the customer:

Subject	Instruction
Radioactive samples	Radioactive samples should not be used on the instrument.
Noncapped PCR tubes	Provide warning about running noncapped PCR tubes/wells. CAUTION Do not open PCR tubes while they are still in the 7700 instrument. This will lead to fluorescent contamination of the instrument, a problem that could require extensive service.
Securing heated cover	Describe how to properly secure the heated cover in place.
Cleaning sample block	Describe how to clean the sample block with alcohol and cotton swabs.
Fluorescent-test and Blackboard storage	Importance of storing the Fluorescent-test and Blackboard trays in a clean area.
Cleaning sample trays	Clean contaminated sample trays with 10% bleach and water.
Cleaning solutions	Use of cleaning solutions to clean the plexiglass window. CAUTION Using organic solvents will damage the finish of the window. Mild cleaning solutions such as normal glass cleaner should be used to clean a dirty window.
Test programs	How to run test programs.
Backing up and storing data	The importance of regularly backing up and storing data.
Maintaining computer	Importance of maintaining the computer to give optimal performance.

Customer Training for the TaqMan Card Upgrade

Basic Training The customer training checklist on page D-6 should be used as a guideline for training customers in upgrading the TaqMan card.

Discuss the following with the customer:

- ◆ ABI PRISM® Card Filling Station theory of operation (wheels and channels, vacuum)
- ◆ Purpose of ABI PRISM® Card Adapter components
- ◆ How to fill cards
- ◆ Sample preparation

Instrument Operation Train the customer on the following subjects:

- ◆ Card software and run differences: setting up runs (run type, thermal cycler conditions)
- ◆ Using a template file to create a run document
- ◆ Determining exposure time before a run (Plate Read)
- ◆ Analyzing data (the applications specialist should participate)
- ◆ Using the Sequence Detection System software, including how to display graphical results and pointing out the importance of the key menus (Preferences, Advanced Options, etc.)

Instrument Maintenance Discuss These Items With the Customer

Subject	Instruction
Radioactive samples	Radioactive samples should not be used on the instrument.
Securing heated cover	Describe how to properly secure the heated cover in place. Emphasize importance of tightening heated cover firmly to avoid leaks of the card.
Cleaning card adapter components	Describe how to clean the optical and reflective plates with water, alcohol and cotton swabs.
Backing up and storing data	The importance of regularly backing up and storing data.
Maintaining computer	Importance of maintaining the computer to give optimal performance.

Completing the Installation

Verifying 7700 Instrument System Performance

To...	Then...
prove that the system does distinguish between 5 000 and 10,000 populations when using the installation kit	<p>use the data analysis features of the software as described under “Analyzing Data” on page 5-13.</p> <p>At least 30 of the 36 PCR replicates from each of the two unknowns must be used to do these calculations.</p> <p>Tubes with known problems (empty, damaged cap, pipetting problems) are the ONLY reactions that can be removed from the analysis.</p> <p>The position of any failed PCR reactions must be verified in a subsequent run (either an installation run or a user run).</p>

Verifying TaqMan Card System Performance

To...	Then...
prove that the card system is getting adequate heat transfer and that reactions are uniform over the entire plate	<p>use the data analysis features of the software as described under “Analyzing the Card Data and Verifying System Performance” on page 6-28.</p> <p>At least 90 of the 96 endogenous control reactions (VIC dye layer) must be used to do these calculations.</p> <p>The position of any failed PCR reactions must be verified in a subsequent run (either an installation run or a user run).</p>

Completing the Training Checklist

When the installation and customer training are complete, obtain the customer’s signature on all reports. This indicates that they are satisfied with the installation and have been trained in the proper use and maintenance of the system. Make copies of the appropriate completed training checklist for both the customer and Applied Biosystems.

Customer Support Education

Discuss the warranty with the customer, including any and all customer responsibilities. Also, inform the customer about how to obtain technical support and how to contact the local Applied Biosystems service department.

Completing the Service Report

Fill out a service report. An accurate and complete service report helps ensure the timely reporting of service data.

Installation Report

Complete and return the installation report. This provides additional information to the product group and manufacturing.

Quality Postcards

Complete and return the quality postcard. This is an important part of the company's continuing efforts to improve product quality.

**Other Pertinent
Data or Information**

Communicate any information that you believe to be especially important about the product/installation to the service product specialist.

Abbreviations



Abbreviation	Meaning
A	Ampere
Btu/h	British thermal units per hour
cm	centimeter
cpm	cubic feet per minute
ft	foot
g	gram
Hz	hertz
in.	inch
kVA	kilovolt-ampere
L	liter
m	meter
min	minute
mL	milliliter
mm	millimeter
MSDS	material safety data sheet
mTorr	unit of pressure
NEMA	National Electrical Manufacturers Assoc.
No.	number
oz	ounce
par.	paragraph
pH	measure of acidity or alkalinity
Rev.	revision
Vac	volts, alternating current
Vdc	volts, direct current
W	watt
°C	degrees Celsius
°F	degrees Fahrenheit
μL	microliter

Material and Equipment

B

In This Appendix This appendix contains the following topics:

Topic	See Page
ABI PRISM 7700 Instrument Installation Kit	B-2
TaqMan Card Upgrade Kit	B-4
Customer-Supplied Materials and Equipment	B-5

ABI PRISM 7700 Instrument Installation Kit

Overview The ABI PRISM® 7700 Sequence Detector includes the following kits, which are used during installation.

Chemical Installation Kit **Note** The kit should be unpacked and stored as soon as it arrives at the laboratory.

The following table lists the contents of the Chemical Installation Kit (P/N 402065):

Chemical Installation Kit

Part Number	Description	Qty (boxes)	Storage Conditions
N8080228	TaqMan® PCR Core Reagents w/Gold	2	-20 °C
4305822	SDS Spectral Calibration Kit	1	-20 °C
401970	TaqMan® DNA Template Reagents	2	-20 °C
401846	β-Actin Control Reagents Kit	1	-20 °C
402823	Protocol, TaqMan® Reagent Kit w/Gold	1	N/A

7700 Packing Kit The following table lists the contents of the 7700 Packing Kit (P/N 604065):

7700 Packing Kit

Part Number	Description
006899	MicroAmp Base/Holder
100383	ABI Mouse Pad
200800	Mac Serial Cable (2: one to connect to the 7700 instrument, one for optional printer)
604094	Blackboard Tray
N8010438	PCR Cap Installation Tool
N8010560	MicroAmp® Optical 96-well Reaction Plate
T-6273	Plate Removal Tool
N8010935	Optical Caps

**Parts Shipped With
the System**

The following parts are also shipped with the system and they are needed to complete the installation:

Shipped with the System

Part Number	Description
4311876	ABI PRISM® 7700 Sequence Detection System software version 1.7 CD-ROM
4313106	Microsoft Excel software
254313	Power Cord
100803	15 inch Monitor
4313103	Power Macintosh® computer
604294	Primer Express® software Version 1.0
604064	Voltage Kit (shipped in the United States, Canada, and Japan)
604216	Voltage Kit (shipped to Europe, Australia, ECO, and Latin America)
4304907	ABI PRISM 7700 User's Manual Set

TaqMan Card Upgrade Kit

Overview The TaqMan® Human Cytokine Card Upgrade Package (P/N 4311899) includes the following the ABI PRISM 7700 and TaqMan Card Upgrade Installation Manual (P/N 4316192) and the following pieces.

Card Hardware Package The following table lists the contents of the Card Hardware Package (P/N 4311060):

Card Hardware Package

Part Number	Description	Qty
4311059	ABI PRISM Card® Adapter	1
4311061	ABI PRISM Card® Filling Station	1

Card H/W and F/W Upgrade Kit The following table lists the contents of the Card H/W and F/W Kit (P/N 4311877):

Card H/W and F/W Kit

Part Number	Description	Qty
4313011	7700 Software Kit, Version 1.7.1	1
4313010	ABI PRISM® 7700 Relative Quantification Software Kit	1
4312990	TaqMan Card Upgrade EPROM Kit	1

Card Reagent Package The following table lists the contents of the Card Reagent Package (P/N 4313072):

Card Reagent Package

Part Number	Description	Qty	Storage Conditions
4305822	SDS Spectral Calibration Kit	1 box	-15 to -25 °C
4329609	TaqMan® Human Cytokine Cards (10)	1 box	2-8 °C, dark
4307577	TaqMan Human Cytokine Card User Protocol	1	N/A
4304437	TaqMan Universal PCR Master Mix	1 box	2-8 °C, dark
N8080234	TaqMan RT Reagents	1 box	2-8 °C, dark
4310982	RNAse P Plate	1 box	2-8 °C, dark
4307281	TaqMan® Control Total RNA (Human)	1 bag	-15 to -25 °C
4318839	20X 18S Primer and TaqMan® Probe Mix	1 bag	-15 to -25 °C
4330147	Human cDNA (Raji)	1 bag	-15 to -25 °C

IMPORTANT The 7700 Software Update Version 1.6.3 Service Kit (P/N 4305070) is not supplied as part of the TaqMan Card Upgrade Kit, but may be required for installation if the customer does not currently have SDS software version 1.6.3 or newer installed on their instrument. Check with the customer when scheduling the card installation and order P/N 4305070 for the installation if necessary.

Customer-Supplied Materials and Equipment

Overview The following sections contain lists of items that are required to complete specific installations but are not included in the Installation Kits. Make sure that the customer has these items on hand before scheduling the installation.

7700 Installation Materials

Instruments	Recommended Source
Centrifuge with 96-well plate adapter	MLS ^a
Microcentrifuge	MLS

a. MLS is Major Lab Suppliers.

Materials	Source
Gloves, disposable, powder-free	MLS
Microcentrifuge tubes, sterile 1.5-mL	MLS
Pipettors, positive-displacement or air-displacement	MLS
Pipette tips, aerosol resistant	MLS
Water, RNase-free, distilled, deionized	MLS

TaqMan Card Upgrade Installation Materials

Instruments	Recommended Source
ABI PRISM 7700 Sequence Detection System	Applied Biosystems
Centrifuge with 96-well plate adapter	MLS
Microcentrifuge	MLS
Vacuum pump, oil-based, that can pull a vacuum down to 25 mTorr (microns)	Welch Two-Stage Belt-Drive Vacuum Pump (P/N 1400B-01)

Materials	Recommended Source
MicroAmp [®] Optical 96-Well Reaction Plate/Optical Caps	Applied Biosystems (P/N 403012)
Gloves, disposable, powder-free	MLS
Microcentrifuge tubes, sterile 1.5-mL	MLS
Pipettors, positive-displacement or air-displacement	MLS
Pipette tips, aerosol resistant	MLS
Polypropylene tubes	MLS
Water, RNase-free, distilled, deionized	MLS

7700 Installation Specifications

C

In This Appendix This appendix contains the following topics:

Topic	See Page
7700 Instrument Installation Specifications	C-2
TaqMan Card Installation Specifications	C-3

7700 Instrument Installation Specifications

Guarantee of Performance When the installation is complete, the ABI PRISM® 7700 Sequence Detection System is guaranteed to distinguish between 5 000 and 10 000 starting molecules at a 99.7% confidence level by analyzing at least thirty 50- μ L replicates per population.

How the Performance Is Verified To support this guarantee, the ABI PRISM 7700 instrument performance is verified during installation by running the chemical installation kit and analyzing the resulting data in a manner designated by Applied Biosystems.

What the Calculations Are Based On The calculations used to determine the confidence level are based on the average copy number and the standard deviation of each unknown population (at least 30 of the 36 PCR reactions of each unknown population must be used in the calculation), as compared to a standard curve that is generated using the PCR TaqMan standards found in the chemical installation kit.

The No Template Controls (NTC) that are also run with the chemical installation kit are used to give both qualitative information about the PCR reactions and quantitative baseline data used during the data analysis.

Calculation To Achieve Confidence Level To achieve the 99.7% confidence level, the following calculation must be true:
$$[(\text{Copy.Unk.1}) - 3(\text{STDev.Unk.1})] > [(\text{Copy.Unk.2}) + 3(\text{STDev.Unk.2})]$$

where

Copy.Unk.1 = Average Copy Number of Unknown #1 (10K replicate population)

STDev.Unk.1 = Standard Deviation of Unknown #1 (10K replicate population)

Copy.Unk.2 = Average Copy Number of Unknown #2 (5K replicate population)

STDev.Unk.2 = Standard Deviation of Unknown #2 (5K replicate population)

TaqMan Card Installation Specifications

Performance Verification	The TaqMan® Human Cytokine Card Upgrade performance is verified during installation by running the chemical installation kit and analyzing the resulting data in a manner designated by Applied Biosystems.
---------------------------------	---

Confidence Levels	The calculations used to determine the confidence levels are the average count and the standard deviation of Ct's for at least 90 of the 96 endogenous control PCR reactions (that is, the VIC dye layer).
--------------------------	--

Installation Calculations	To pass installation, the following must be true using at least 90 of the 96 endogenous control reactions: <ul style="list-style-type: none">◆ Average Ct \leq 20.0◆ Standard Deviation of Ct's \leq 1.0
----------------------------------	---

Installation Report Checklists

D

**Topics in This
Appendix**

This appendix contains the following checklists:

Topic	See Page
Preinstallation Checklists	D-2
Customer Training Checklists	D-5
Laboratory Safety Checklist	D-7
Installation Reports	D-9

Preinstallation Checklists

About these Checklists These checklists are for customer use to ensure that all preparations have been made for installation. Review the appropriate checklist with the customer before scheduling the installation, to make sure that everything has been checked off and is ready.

ABI 7700 Instrument Checklist

Checklist for the ABI PRISM® 7700 Sequence Detector

✓ if ready	Date Confirmed	Components
General		
		Received instrument(s) and verified that cartons are intact.
		Read the MSDSs provided by the chemical manufacturers.
		Read the Site Preparation and Safety Manual.
		Verified that instrument(s), serial number(s), and system configuration as shown on the packing list are the same as ordered.
		Reserved two uninterrupted days for in-lab training during installation.
		Unpacked and stored contents of chemical installation kit.
Electrical		
		A dedicated 5.0 kVA power line and ground, or a 5.0 kVA power line with a line conditioner or UPS, is in place.
		A power line circuit breaker rated at 30 A is in place.
		A NEMA Twistlock or similar power outlet receptacle is within 2.5 m (8 ft.) of the instrument location.
		Instrument voltage, if specified on the packing list, matches the voltage available in the laboratory.
Laboratory		
		Laboratory bench is of correct dimensions to accommodate the system and is situated so the instrument is accessible to the installer on all four sides.
		Instrument location accommodates a total vertical clearance of 130 cm (51 in.).
		Laboratory safety requirements as specified in the Site Preparation and Safety Manual are met.
		Laboratory environmental requirements as specified in the Site Preparation and Safety Manual are met.
		Room ventilation accommodates instrument heat output.
		Deionized water is on site.
		Proper waste disposal method for hazardous chemical waste is established.
Equipment		
		Macintosh®-compatible printer is available if not ordered with instrument.
		Jouan 422 or other large centrifuge to accommodate 96 plastic micro-Amp tube racks and generate a minimum of 1400 gs is available.

Checklist for the ABI PRISM® 7700 Sequence Detector (continued)

✓ if ready	Date Confirmed	Components
		Ice bucket is available.
		Autoclave is available.
		–20 °C laboratory freezer is available.
Consumables		
		Additional Applied Biosystems reagent kits ordered.
		Additional computer supplies (paper, disks, etc.) ordered.
		Chemically resistant disposable gloves available.
		Micro pipettes and tips (Pipetman models P20, P200, P1000; or Eppendorf 1–10 µL, 10–100 µL, and 100–1000 µL and tips) ordered.
		Disposable test tubes available.
		Capped tubes: 2 mL, and 5 mL, 10 mL, or 15 mL available.
		Tube racks available.
		Wet ice available.
		TE Buffer (1X or 100X) available.

TaqMan Card Upgrade Checklist

Checklist for the TaqMan® Human Cytokine Card Upgrade

✓ if ready	Date Confirmed	Components
General		
		Received (1) hardware and (2) chemical installation boxes (3 total) and verified that cartons are intact.
		Reserved one uninterrupted day for the installation.
		Unpacked and stored contents of chemical installation kit.
		Determined the current version of SDS software installed on the ABI PRISM 7700. It should be 1.6.3 or newer in order to perform the TaqMan card upgrade. IMPORTANT If the current software installed is 1.6 or older, CSE will need to order the 7700 Software Update Version 1.6.3 Service Kit (P/N 4305070) for installation during upgrade service visit.
Equipment		
		Welch or other oil-based vacuum pump has been connected to vacuum gauge/hose assembly that is included in hardware kit and verified to pull a vacuum down to 600 mTorr or better.
		Jouan 422 or other large centrifuge to accommodate 96 plastic micro-Amp tube-racks and generate a minimum of 1400 gs is available.
		Ice bucket is available.
		–20 °C laboratory freezer is available.
Consumables		
		Chemically resistant disposable gloves available.

Checklist for the TaqMan® Human Cytokine Card Upgrade *(continued)*

✓ if ready	Date Confirmed	Components
		Micro pipettes and tips (Pipetman models P20, P200, P1000; or Eppendorf 1–10 μL , 10–100 μL , and 100–1000 μL and tips) ordered.
		Disposable test tubes available.
		Capped tubes: 2 mL, and 5 mL, 10 mL, or 15 mL available.
		Tube racks available.
		Wet ice available.

Customer Training Checklists

ABI Prism 7700 Instrument Training Checklist

ABI PRISM 7700 Instrument Customer Training Checklist

✓ after training or discussion	
General Topics (User's Manual)	
	TaqMan theory
	Overview of the optical system
	Description of laser safety features
Electrical Information	
	Location of power switch and circuit breakers
	Connection of cables to the Macintosh computer
	Location of reset switch
Instrument Operation	
	How the chemical installation kit is used
	Proper storage of chemicals
	Loading TaqMan PCR Trays into the instrument
	Launching the Sequence Detection System software
	Setting up a run: Selecting wells, TC conditions, and Sample Type
	Starting, monitoring, and ending runs
	Use of the SDS Spectral Calibration Kit (P/N 4305822)
	Purpose and use of the Blackboard Tray
	Use of the Instrument's test programs
	Use of On-line help
System Maintenance	
	Use and care of heated cover's knob to prevent instrument damage
	Assembly of reaction trays to prevent contamination of the instrument
	Cleaning the sample block and post-PCR reaction trays and retainers
	Bringing the sample block to 25 degrees C before removing trays
	Cleaning the heated-cover's plexiglass door (do not use organic solvents)
	Warning about the use of radioactive samples and non-capped reaction tubes which can damage the instrument and can void the warranty
	Rebuilding the Macintosh computer desktop and routinely backing-up data files
Data Analysis	
	Description of analysis options
	Graphical and tabulated results of analyzed data
	Data back-up and removing files from the Macintosh computer
Instrument Performance Verification	
	Verify instrument performance (chemical installation kit results)
	Inform customer of technical support and service department phone numbers
	Discuss the customer's responsibilities and what is covered under the warranty

**TaqMan Card
Training Checklist**

TaqMan Human Cytokine Card Customer Training Checklist

✓ after training or discussion	
General Topics	
	TaqMan theory
	Sample preparation (importance of using random hexamers in RT reaction)
	Overview of the card filling system (vacuum filling technique)
	Importance of good vacuum for proper filling of cards
	Purpose of ABI PRISM® Card Adapter components
	How the chemical installation kit is used (VIC dye layer/endogenous control amplification ONLY)
	Proper storage of chemicals
	Use of the SDS Spectral Calibration Kit (P/N 4305822)
TaqMan Card Hardware Operation	
	Aligning card in ABI PRISM® Filling Station (folding adhesive flap back)
	Filling card with reagents
Instrument Operation	
	Loading TaqMan card and adapter components into the 7700 instrument
	Launching the Sequence Detection System software
	Changes to 7700 instrument software for TaqMan card upgrade
	Determining correct exposure time for card runs vs. plate runs
	Setting up a card run: Selecting wells, TC conditions, and Sample Type
	Creating a run document from a template file
	Starting, monitoring, and ending runs
	Proper tightening of heated cover's knob to prevent card reagent leakage
System Maintenance	
	Cleaning the card adapters
	Warning about the use of radioactive samples and noncapped reaction tubes which can damage the instrument and can void the warranty
	Rebuilding the Macintosh computer desktop and routinely backing-up data files
Data Analysis	
	Description of analysis options
	Graphical and tabulated results of analyzed data
	Exporting results files for downstream analysis (tab-delimited format)
Instrument Performance Verification	
	Verify instrument performance (chemical installation kit results)
	Inform customer of technical support and service department phone numbers
	Discuss the customer's responsibilities and what is covered under the warranty

Laboratory Safety Checklist

Laboratory Safety For instrument-specific safety information, refer to the ABI PRISM 7700 *Site Preparation and Safety Guide* (P/N 4304908).

Before installing or servicing any instrument, take a moment to review the following conditions in the vicinity of your work area:

Laboratory Safety Conditions

Condition	What you should do...
Who is responsible for the lab you are working in?	Identify a person who knows what is going on in the laboratory and is readily available to help or warn you in an emergency.
Do you see any hazard warning signs indicating the presence of radioactive or pathogenic materials, carcinogens, poisonous gases, exposed high-voltage sources, lasers, etc.?	Check with the person responsible for the laboratory, and satisfy yourself that the area is safe to work in. Satisfy yourself that there is no exposure hazard.
Where are the safety shower, eye wash, first-aid kit, fire extinguisher, spill collection materials, and other safety equipment located?	Locate emergency safety equipment in the lab. Some Applied Biosystems reagents are hazardous, and there may be other hazardous chemicals nearby. If something were to spill or leak, know where to wash yourself, how to put out a fire, and what you will use to clean up a spill.
Does the room have good ventilation?	Make sure that any chemical handling hoods are in good working order. Check that the ABI PRISM 7700 instrument is properly vented. If you notice any chemical odors, find the source and ask to see the material safety data sheet on the chemical causing the odor. Do not inhale any chemicals, and do not leave any containers uncapped. Consider every chemical to be potentially harmful.
Are there chemical wastes?	Do not handle or dispose of chemical waste. Have the customer do so. Tell the customer to refer to the waste profile for complete information on handling and disposing of waste. Advise customers to wear gloves, protective eye wear, and a laboratory coat when handling or disposing of chemical waste, which should be in accord with all local, regional, and federal regulations and laws.

Laboratory Safety Conditions *(continued)*

Condition	What you should do...						
<p>What is the evacuation plan for the building you are in?</p>	<p>Know what the posted routes are and where your assembly point is in an evacuation emergency.</p> <p>Check that your escape route is not blocked or obstructed by locked doors, boxes, chemicals on the floor, empty crates, etc.</p>						
<p>What special safety precautions are required in the lab?</p>	<p>Determine what safety precautions are followed by other laboratory personnel, and adhere to the same (or better) practices yourself.</p> <p>Be sure to wear safety glasses, gloves, a laboratory coat, and any other attire listed in the material safety data sheets.</p> <table border="1" data-bbox="943 709 1421 1060"> <thead> <tr> <th data-bbox="948 716 1182 751">If...</th> <th data-bbox="1182 716 1416 751">Then...</th> </tr> </thead> <tbody> <tr> <td data-bbox="948 751 1182 934"> <p>you have concerns about your personal safety when working in a customer's laboratory</p> </td> <td data-bbox="1182 751 1416 934"> <p>voice your concerns politely.</p> </td> </tr> <tr> <td data-bbox="948 934 1182 1054"> <p>you cannot resolve the situation to your satisfaction</p> </td> <td data-bbox="1182 934 1416 1054"> <p>telephone the regional Field Service office for support.</p> </td> </tr> </tbody> </table>	If...	Then...	<p>you have concerns about your personal safety when working in a customer's laboratory</p>	<p>voice your concerns politely.</p>	<p>you cannot resolve the situation to your satisfaction</p>	<p>telephone the regional Field Service office for support.</p>
If...	Then...						
<p>you have concerns about your personal safety when working in a customer's laboratory</p>	<p>voice your concerns politely.</p>						
<p>you cannot resolve the situation to your satisfaction</p>	<p>telephone the regional Field Service office for support.</p>						

Installation Reports

ABI Prism 7700 Instrument Installation Report

Note This information should be filled in at the completion of the installation and given to the customer. Forward a copy of this report to the ABI PRISM 7700 instrument service product specialist. The installation report that is shipped with the 7700 instrument includes the preceding ABI PRISM 7700 instrument training checklist.

7700 instrument serial number _____

Board revision _____

Firmware revision _____

9600 Performance Test Pass? _____

Test results (Fluor Test) _____

Laser current at 10 mW (in volts) _____

Dark current (Max) (in amps) _____

Peak bin intensity _____

(Intensity of lowest well) (approx.) _____

(Ratio of highest/lowest well) (approx.) _____

Background (Blackboard tray-intensity of highest well) _____

Noise (Filter tray-Std. Dev. of highest well) _____

Well-to-well variability (Filter tray-Std. Dev. of intensities) _____

7700 instrument specification calculation results _____

Standard curve correlation coefficient _____

ABI PRISM 7700 instrument customer training checklist _____

Comments _____

Customer Name: _____

Customer Signature: _____

Service Engineer: _____

Date: _____

**TaqMan Card
Installation Report**

Note This information should be filled in at the completion of the upgrade installation and given to the customer. Forward a copy of this report to the ABI PRISM 7700 instrument service product specialist.

7700 instrument serial number _____

ABI PRISM Card Filling Station serial number _____

Board revision _____

Firmware revision _____

7700 instrument specification calculation results _____

Standard curve correlation coefficient _____

TaqMan card specification verification: average 18S Ct _____

TaqMan card specification verification: standard deviation 18S Ct _____

TaqMan card customer training checklist _____

Comments _____

Customer Name: _____

Customer Signature: _____

Service Engineer: _____

Date: _____

Creating a TaqMan Card Template



Why Create a Template File?

A template file alleviates the need for repetitious construction of TaqMan® Human Cytokine Card plate documents. An unlimited number of documents can be created from a single template file.

Creating a Template for a Card Run

Step	Action								
1	Open the Sequence Detection Systems (SDS) software.								
2	From the File menu, select New .								
3	Configure a new plate document with the following attributes. <table border="1"><thead><tr><th>From menu...</th><th>Select...</th></tr></thead><tbody><tr><td>Plate Type</td><td>Single Reporter</td></tr><tr><td>Plate Format</td><td>The Card</td></tr><tr><td>Run</td><td>Real Time</td></tr></tbody></table> <p>Note TaqMan Human Cytokine Cards can be used only with version 1.7.1 or later of the SDS software. If the options above do not appear in the New Plate dialog box, update your instrument with a newer version of the software.</p>	From menu...	Select...	Plate Type	Single Reporter	Plate Format	The Card	Run	Real Time
From menu...	Select...								
Plate Type	Single Reporter								
Plate Format	The Card								
Run	Real Time								

Configuring the PCR Thermal Profile

Step	Action																														
1	Click the Thermal Cycler Conditions button.																														
2	Configure the thermal cycling profile. <div style="text-align: center; margin: 10px 0;"> </div> <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 10px;"> <thead> <tr> <th>Stage</th> <th>1</th> <th>2</th> <th colspan="2">3</th> </tr> </thead> <tbody> <tr> <td>Description</td> <td>UNG Activation^a</td> <td>AmpliTaq Gold Activation^b</td> <td colspan="2">PCR</td> </tr> <tr> <td></td> <td>HOLD</td> <td>HOLD</td> <td colspan="2">CYCLE (35 cycles)</td> </tr> <tr> <td></td> <td></td> <td></td> <td>Denature</td> <td>Anneal/Extend</td> </tr> <tr> <td>Temperature</td> <td>50 °C</td> <td>99 °C^c</td> <td>99 °C^c</td> <td>60 °C</td> </tr> <tr> <td>Time</td> <td>2 min</td> <td>10 min</td> <td>15 sec</td> <td>1:00 min</td> </tr> </tbody> </table> <p style="margin-top: 10px;"> ^a. Required for optimal AmpErase® UNG enzyme activity. ^b. Required for AmpliTaq Gold® DNA Polymerase activation. ^c. Because of the unique thermokinetic properties of the TaqMan Human Cytokine Card, a setting of 99 °C is necessary for this step. </p>	Stage	1	2	3		Description	UNG Activation ^a	AmpliTaq Gold Activation ^b	PCR			HOLD	HOLD	CYCLE (35 cycles)					Denature	Anneal/Extend	Temperature	50 °C	99 °C ^c	99 °C ^c	60 °C	Time	2 min	10 min	15 sec	1:00 min
Stage	1	2	3																												
Description	UNG Activation ^a	AmpliTaq Gold Activation ^b	PCR																												
	HOLD	HOLD	CYCLE (35 cycles)																												
			Denature	Anneal/Extend																											
Temperature	50 °C	99 °C ^c	99 °C ^c	60 °C																											
Time	2 min	10 min	15 sec	1:00 min																											
3	Click the Sample Volume text field and type 100 µL .																														

Configuring the FAM Dye Layer

Step	Action
1	Select FAM from the Dye Layer pop-up menu.
2	Select all wells of the plate document.
3	Select TARG - RelQ Target from the Sample Type pop-up menu. The SDS software labels all selected wells as TARG.

Sample Type: TARG - RelQ Target Thermal Cycler Conditions: 7700 Single Reporter The Card
 Sample Name: _____
 Replicate: _____
 Quantity: _____
 Comment: _____

Show Analysis Dye Layer: FAM

	1	2	3	4	5	6	7	8	9	10	11	12
A	TARG A1	TARG A2	TARG A3	TARG A4	TARG A5	TARG A6	TARG A7	TARG A8	TARG A9	TARG A10	TARG A11	TARG A12
B	TARG B1	TARG B2	TARG B3	TARG B4	TARG B5	TARG B6	TARG B7	TARG B8	TARG B9	TARG B10	TARG B11	TARG B12
C	TARG C1	TARG C2	TARG C3	TARG C4	TARG C5	TARG C6	TARG C7	TARG C8	TARG C9	TARG C10	TARG C11	TARG C12
D	TARG D1	TARG D2	TARG D3	TARG D4	TARG D5	TARG D6	TARG D7	TARG D8	TARG D9	TARG D10	TARG D11	TARG D12
E	TARG E1	TARG E2	TARG E3	TARG E4	TARG E5	TARG E6	TARG E7	TARG E8	TARG E9	TARG E10	TARG E11	TARG E12
F	TARG F1	TARG F2	TARG F3	TARG F4	TARG F5	TARG F6	TARG F7	TARG F8	TARG F9	TARG F10	TARG F11	TARG F12
G	TARG G1	TARG G2	TARG G3	TARG G4	TARG G5	TARG G6	TARG G7	TARG G8	TARG G9	TARG G10	TARG G11	TARG G12
H	TARG H1	TARG H2	TARG H3	TARG H4	TARG H5	TARG H6	TARG H7	TARG H8	TARG H9	TARG H10	TARG H11	TARG H12

Configuring the VIC Dye Layer

To configure the VIC dye layer:




Step	Action
1	Select VIC from the Dye Layer pop-up menu.
2	Select all wells of the plate document.
3	Select ENDO - RelQ Endogenous Control from the Sample Type pop-up menu.
4	Select cells A1–A4.

Sample Type: ENDO - RelQ Endogenous...
 Sample Name: _____
 Replicate: _____
 Quantity: _____

Show Analysis Dye Layer: VIC


	1	2	3	4	5
A	ENDO A1	ENDO A2	ENDO A3	ENDO A4	ENDO A5
B	ENDO B1	ENDO B2	ENDO B3	ENDO B4	ENDO B5

To configure the VIC dye layer: (continued)

Step	Action																																																																																																																					
5	<p>Click the Sample Name text field and type IL-1alpha.</p>  <p>Sample Type: ENDO - RelQ Endogenous... Sample Name: IL-1alpha Replicate: Quantity: Show Analysis Dye Layer: VIC</p> <table border="1"> <thead> <tr> <th></th> <th>1</th> <th>2</th> <th>3</th> <th>4</th> <th>5</th> </tr> </thead> <tbody> <tr> <th>A</th> <td>ENDO IL-1 alpha</td> <td>ENDO IL-1 alpha</td> <td>ENDO IL-1 alpha</td> <td>ENDO IL-1 alpha</td> <td>ENDO A5</td> </tr> <tr> <th>B</th> <td>ENDO B1</td> <td>ENDO B2</td> <td>ENDO B3</td> <td>ENDO B4</td> <td>ENDO B5</td> </tr> </tbody> </table>		1	2	3	4	5	A	ENDO IL-1 alpha	ENDO IL-1 alpha	ENDO IL-1 alpha	ENDO IL-1 alpha	ENDO A5	B	ENDO B1	ENDO B2	ENDO B3	ENDO B4	ENDO B5																																																																																																			
	1	2	3	4	5																																																																																																																	
A	ENDO IL-1 alpha	ENDO IL-1 alpha	ENDO IL-1 alpha	ENDO IL-1 alpha	ENDO A5																																																																																																																	
B	ENDO B1	ENDO B2	ENDO B3	ENDO B4	ENDO B5																																																																																																																	
6	<p>Repeat steps 4–5 for each target cytokine so that the plate document mirrors the assay configuration of the TaqMan Human Cytokine Card, as shown below.</p>  <p>1 2 3 4 5 6 7 8 9 10 11 12 A IL-1α IL-1β IL-2 B IL-3 IL-4 IL-5 C IL-6 IL-7 IL-8 D IL-10 IL-12p35 IL-12p40 E IL-13 IL-15 IL-17 F IL-18 G-CSF GM-CSF G M-CSF IFN-γ LT-β H TGF-β TNF-α TNF-β</p> <p style="text-align: right;">GRT540</p>  <p>Sample Type: ENDO - RelQ Endogenous... Sample Name: Replicate: Quantity: Show Analysis Dye Layer: VIC</p> <table border="1"> <thead> <tr> <th></th> <th>1</th> <th>2</th> <th>3</th> <th>4</th> <th>5</th> <th>6</th> <th>7</th> <th>8</th> <th>9</th> <th>10</th> <th>11</th> <th>12</th> </tr> </thead> <tbody> <tr> <th>A</th> <td>ENDO IL-1 ALPHA</td> <td>ENDO IL-1 ALPHA</td> <td>ENDO IL-1 ALPHA</td> <td>ENDO IL-1 ALPHA</td> <td>ENDO IL-1 BETA</td> <td>ENDO IL-1 BETA</td> <td>ENDO IL-1 BETA</td> <td>ENDO IL-1 BETA</td> <td>ENDO IL-2</td> <td>ENDO IL-2</td> <td>ENDO IL-2</td> <td>ENDO IL-2</td> </tr> <tr> <th>B</th> <td>ENDO IL-3</td> <td>ENDO IL-3</td> <td>ENDO IL-3</td> <td>ENDO IL-3</td> <td>ENDO IL-4</td> <td>ENDO IL-4</td> <td>ENDO IL-4</td> <td>ENDO IL-4</td> <td>ENDO IL-5</td> <td>ENDO IL-5</td> <td>ENDO IL-5</td> <td>ENDO IL-5</td> </tr> <tr> <th>C</th> <td>ENDO IL-6</td> <td>ENDO IL-6</td> <td>ENDO IL-6</td> <td>ENDO IL-6</td> <td>ENDO IL-7</td> <td>ENDO IL-7</td> <td>ENDO IL-7</td> <td>ENDO IL-7</td> <td>ENDO IL-8</td> <td>ENDO IL-8</td> <td>ENDO IL-8</td> <td>ENDO IL-8</td> </tr> <tr> <th>D</th> <td>ENDO IL-10</td> <td>ENDO IL-10</td> <td>ENDO IL-10</td> <td>ENDO IL-10</td> <td>ENDO IL-12p35</td> <td>ENDO IL-12p35</td> <td>ENDO IL-12p35</td> <td>ENDO IL-12p35</td> <td>ENDO IL-12p40</td> <td>ENDO IL-12p40</td> <td>ENDO IL-12p40</td> <td>ENDO IL-12p40</td> </tr> <tr> <th>E</th> <td>ENDO IL-13</td> <td>ENDO IL-13</td> <td>ENDO IL-13</td> <td>ENDO IL-13</td> <td>ENDO IL-15</td> <td>ENDO IL-15</td> <td>ENDO IL-15</td> <td>ENDO IL-15</td> <td>ENDO IL-17</td> <td>ENDO IL-17</td> <td>ENDO IL-17</td> <td>ENDO IL-17</td> </tr> <tr> <th>F</th> <td>ENDO IL-18</td> <td>ENDO IL-18</td> <td>ENDO IL-18</td> <td>ENDO IL-18</td> <td>ENDO G-CSF</td> <td>ENDO G-CSF</td> <td>ENDO G-CSF</td> <td>ENDO G-CSF</td> <td>ENDO GM-CSF</td> <td>ENDO GM-CSF</td> <td>ENDO GM-CSF</td> <td>ENDO GM-CSF</td> </tr> <tr> <th>G</th> <td>ENDO M-CSF</td> <td>ENDO M-CSF</td> <td>ENDO M-CSF</td> <td>ENDO M-CSF</td> <td>ENDO IFN-GAM</td> <td>ENDO IFN-GAM</td> <td>ENDO IFN-GAM</td> <td>ENDO IFN-GAM</td> <td>ENDO LT-BETA</td> <td>ENDO LT-BETA</td> <td>ENDO LT-BETA</td> <td>ENDO LT-BETA</td> </tr> <tr> <th>H</th> <td>ENDO TGF-BETA</td> <td>ENDO TGF-BETA</td> <td>ENDO TGF-BETA</td> <td>ENDO TGF-BETA</td> <td>ENDO TNF-ALP</td> <td>ENDO TNF-ALP</td> <td>ENDO TNF-ALP</td> <td>ENDO TNF-ALP</td> <td>ENDO TNF-BETA</td> <td>ENDO TNF-BETA</td> <td>ENDO TNF-BETA</td> <td>ENDO TNF-BETA</td> </tr> </tbody> </table>		1	2	3	4	5	6	7	8	9	10	11	12	A	ENDO IL-1 ALPHA	ENDO IL-1 ALPHA	ENDO IL-1 ALPHA	ENDO IL-1 ALPHA	ENDO IL-1 BETA	ENDO IL-1 BETA	ENDO IL-1 BETA	ENDO IL-1 BETA	ENDO IL-2	ENDO IL-2	ENDO IL-2	ENDO IL-2	B	ENDO IL-3	ENDO IL-3	ENDO IL-3	ENDO IL-3	ENDO IL-4	ENDO IL-4	ENDO IL-4	ENDO IL-4	ENDO IL-5	ENDO IL-5	ENDO IL-5	ENDO IL-5	C	ENDO IL-6	ENDO IL-6	ENDO IL-6	ENDO IL-6	ENDO IL-7	ENDO IL-7	ENDO IL-7	ENDO IL-7	ENDO IL-8	ENDO IL-8	ENDO IL-8	ENDO IL-8	D	ENDO IL-10	ENDO IL-10	ENDO IL-10	ENDO IL-10	ENDO IL-12p35	ENDO IL-12p35	ENDO IL-12p35	ENDO IL-12p35	ENDO IL-12p40	ENDO IL-12p40	ENDO IL-12p40	ENDO IL-12p40	E	ENDO IL-13	ENDO IL-13	ENDO IL-13	ENDO IL-13	ENDO IL-15	ENDO IL-15	ENDO IL-15	ENDO IL-15	ENDO IL-17	ENDO IL-17	ENDO IL-17	ENDO IL-17	F	ENDO IL-18	ENDO IL-18	ENDO IL-18	ENDO IL-18	ENDO G-CSF	ENDO G-CSF	ENDO G-CSF	ENDO G-CSF	ENDO GM-CSF	ENDO GM-CSF	ENDO GM-CSF	ENDO GM-CSF	G	ENDO M-CSF	ENDO M-CSF	ENDO M-CSF	ENDO M-CSF	ENDO IFN-GAM	ENDO IFN-GAM	ENDO IFN-GAM	ENDO IFN-GAM	ENDO LT-BETA	ENDO LT-BETA	ENDO LT-BETA	ENDO LT-BETA	H	ENDO TGF-BETA	ENDO TGF-BETA	ENDO TGF-BETA	ENDO TGF-BETA	ENDO TNF-ALP	ENDO TNF-ALP	ENDO TNF-ALP	ENDO TNF-ALP	ENDO TNF-BETA	ENDO TNF-BETA	ENDO TNF-BETA	ENDO TNF-BETA
	1	2	3	4	5	6	7	8	9	10	11	12																																																																																																										
A	ENDO IL-1 ALPHA	ENDO IL-1 ALPHA	ENDO IL-1 ALPHA	ENDO IL-1 ALPHA	ENDO IL-1 BETA	ENDO IL-1 BETA	ENDO IL-1 BETA	ENDO IL-1 BETA	ENDO IL-2	ENDO IL-2	ENDO IL-2	ENDO IL-2																																																																																																										
B	ENDO IL-3	ENDO IL-3	ENDO IL-3	ENDO IL-3	ENDO IL-4	ENDO IL-4	ENDO IL-4	ENDO IL-4	ENDO IL-5	ENDO IL-5	ENDO IL-5	ENDO IL-5																																																																																																										
C	ENDO IL-6	ENDO IL-6	ENDO IL-6	ENDO IL-6	ENDO IL-7	ENDO IL-7	ENDO IL-7	ENDO IL-7	ENDO IL-8	ENDO IL-8	ENDO IL-8	ENDO IL-8																																																																																																										
D	ENDO IL-10	ENDO IL-10	ENDO IL-10	ENDO IL-10	ENDO IL-12p35	ENDO IL-12p35	ENDO IL-12p35	ENDO IL-12p35	ENDO IL-12p40	ENDO IL-12p40	ENDO IL-12p40	ENDO IL-12p40																																																																																																										
E	ENDO IL-13	ENDO IL-13	ENDO IL-13	ENDO IL-13	ENDO IL-15	ENDO IL-15	ENDO IL-15	ENDO IL-15	ENDO IL-17	ENDO IL-17	ENDO IL-17	ENDO IL-17																																																																																																										
F	ENDO IL-18	ENDO IL-18	ENDO IL-18	ENDO IL-18	ENDO G-CSF	ENDO G-CSF	ENDO G-CSF	ENDO G-CSF	ENDO GM-CSF	ENDO GM-CSF	ENDO GM-CSF	ENDO GM-CSF																																																																																																										
G	ENDO M-CSF	ENDO M-CSF	ENDO M-CSF	ENDO M-CSF	ENDO IFN-GAM	ENDO IFN-GAM	ENDO IFN-GAM	ENDO IFN-GAM	ENDO LT-BETA	ENDO LT-BETA	ENDO LT-BETA	ENDO LT-BETA																																																																																																										
H	ENDO TGF-BETA	ENDO TGF-BETA	ENDO TGF-BETA	ENDO TGF-BETA	ENDO TNF-ALP	ENDO TNF-ALP	ENDO TNF-ALP	ENDO TNF-ALP	ENDO TNF-BETA	ENDO TNF-BETA	ENDO TNF-BETA	ENDO TNF-BETA																																																																																																										

Saving As a Template File

To save the plate document as a template file:

Step	Action
1	Select Save As from the File menu. The Save As dialog box appears.
2	Save as a template file. a. Click the Stationery Pad radio button. b. Click the Save this document as text field and type a name for the template file.
	 <p>Click here</p> <p>Click here and type a file name</p>
3	Select a location for the software to place the template file.
4	Click OK . The software saves the template file.


Software Upgrade Install Procedure

F

Overview ABI PRISM® 7700 Sequence Detection System software is continually being revised and improved to meet the needs of our customers. There have been several software updates since the instrument was first released, some of which required a service installation.

The procedure below describes how to perform a software update when hardware (chip) changes are involved. It is a generic procedure and may vary slightly depending on the specific version of software being updated. Ensure that the 7700 instrument is functioning properly.

Ensuring the Instrument Is Functioning Properly

Step	Action
1	<p>Ensure that the 7700 instrument is functioning properly before starting the upgrade.</p> <ol style="list-style-type: none"> Use a Fluorescent Test tray (ABD P/N 4305178 or T-6222) to perform a plate read and verify that the raw spectra for all 96 wells look as expected. Open a Real Time, New Plate, Single Reporter with all 96 wells designated as Unknown and the following thermal cycler program:  <p>While running this quick test, verify that the sample and cover temperatures of the thermal cycler are being displayed and that the raw data can be viewed upon completion of the run.</p>

**Replacing Thermal
Cycler EPROMS**

Step	Action
1	Turn off the 7700 instrument.
2	<p>Replace the thermal cycler EPROMS (circuit reference U4 and U5 on the 9600 CPU PCB) with the parts labeled “U4” and “U5”.</p> <p>IMPORTANT Please note the orientation of the notch on each EPROM when placing it on the 9600 CPU PCB.</p> <p>CAUTION Ground yourself before touching any components. Failure to do so could damage electrical components.</p>
3	Turn on the 7700 instrument and wait for the ready light to appear.
4	<p>Launch the current 7700 instrument software, open a plate document, click Show Analysis and check that sample and cover temperatures of the thermal cycler are being displayed by the Sequence Detection System software.</p> <p>This confirms that there is a communication link between the Macintosh® computer and the thermal cycler.</p>
5	Quit the 7700 instrument software and turn off the 7700 instrument.

**Setting Dip Switches
and Replacing 7700
EPROM**

Step	Action																		
1	<p>Set the dip switches (circuit reference SW1) on the 7700 instrument PCA Controller Board (P/N 683324) to the following configuration in order to temporarily disable firmware downloading.</p> <table border="1" data-bbox="532 997 943 1356"> <thead> <tr> <th>Dip Switch</th> <th>Set to</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>ON</td> </tr> <tr> <td>2</td> <td>OFF</td> </tr> <tr> <td>3</td> <td>OFF</td> </tr> <tr> <td>4</td> <td>OFF</td> </tr> <tr> <td>5</td> <td>OFF</td> </tr> <tr> <td>6</td> <td>OFF</td> </tr> <tr> <td>7</td> <td>OFF</td> </tr> <tr> <td>8</td> <td>OFF</td> </tr> </tbody> </table> <p>Note The first five 7700 instrument’s were built with a MUX that requires dip switch #2 to be in the ON position. These instruments have a note on the top of the MUX that indicates the proper dip switch position.</p>	Dip Switch	Set to	1	ON	2	OFF	3	OFF	4	OFF	5	OFF	6	OFF	7	OFF	8	OFF
Dip Switch	Set to																		
1	ON																		
2	OFF																		
3	OFF																		
4	OFF																		
5	OFF																		
6	OFF																		
7	OFF																		
8	OFF																		
2	<p>Replace the XILINX PROM (circuit reference U41) on the 7700 instrument PCA controller board (P/N 683324) with the chip labeled “U41”.</p> <p>Note Some software upgrades do NOT include replacement of the XILINX PROM.</p>																		
3	Replace the firmware EPROM (circuit reference U41) on the 7700 instrument’s PCA Controller Board (P/N 683324) with the new part labeled “U41”.																		
4	Turn on the 7700 instrument and wait for the ready light to appear.																		
5	Turn off the 7700 instrument.																		

**Setting Dip Switches
To Re-Enable
Firmware
Downloading**

Step	Action																		
1	<p>Set the dip switches as follows to re-enable firmware downloading:</p> <table border="1"> <thead> <tr> <th>Dip Switch</th> <th>Set to</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>OFF</td> </tr> <tr> <td>2</td> <td>OFF</td> </tr> <tr> <td>3</td> <td>OFF</td> </tr> <tr> <td>4</td> <td>OFF</td> </tr> <tr> <td>5</td> <td>OFF</td> </tr> <tr> <td>6</td> <td>OFF</td> </tr> <tr> <td>7</td> <td>OFF</td> </tr> <tr> <td>8</td> <td>OFF</td> </tr> </tbody> </table> <p>Note The first five 7700 instrument's were built with a MUX that requires dip switch #2 to be in the ON position. These instruments have a note on the top of the MUX that indicates the proper dip switch position.</p>	Dip Switch	Set to	1	OFF	2	OFF	3	OFF	4	OFF	5	OFF	6	OFF	7	OFF	8	OFF
Dip Switch	Set to																		
1	OFF																		
2	OFF																		
3	OFF																		
4	OFF																		
5	OFF																		
6	OFF																		
7	OFF																		
8	OFF																		
2	Turn on the 7700 instrument.																		

**Installing the New
Sequence Detection
System Software**

Step	Action
1	<p>Delete all Sequence Detection System software on the Macintosh computer.</p> <p>Verify that all of the application programs have been removed by trying to open a SDS software run document.</p>
2	Install the new Sequence Detection System software from the CD-ROM.
3	Launch the new Sequence Detection System software.
4	<p>Open a plate document, click Show Analysis and check that sample and cover temperatures of the thermal cycler are being displayed.</p> <p>This confirms that there is a communication link between the Macintosh computer and the thermal cycler.</p>
5	<p>Verify that the new EPROM version is listed.</p> <ol style="list-style-type: none"> Select Diagnostics from the Instrument menu, and 7700 Instrument Verification from the submenu. The Instrument Tests dialog box appears. Click Run Tests. When the tests are completed, verify that the new EPROM version is displayed and click Done. <p>Note This test may show a false failure for the shutters. Listen for the sound of the shutters opening and closing during the test to verify that they are working.</p>
6	Using a Fluorescent Test tray, perform a real-time run for all 96 wells and check the raw spectra to verify that it is possible to collect data.

Index

Numerics

- 50 hertz areas, configuring the Thermal Cycler 3-5
- 7700 instrument
 - completing the installation 7-4 to 7-5
 - customer training checklist D-5
 - installation kit B-2 to B-3
 - installation report D-9
 - installation specifications C-2
 - instrument verification
 - chemical installation kit guidelines 5-3
 - preparing and running RNase P install plate 5-9 to 5-12
 - preparing β -actin install plate 5-3 to 5-5
 - setting up and starting β -actin install run 5-6 to 5-8
 - laboratory requirements 2-6 to 2-7
 - environment conditions and requirements 2-6
 - laboratory space 2-6
 - pre-installation checklist 2-6
 - loading the card into the instrument 6-19 to 6-21
 - preinstallation checklist D-2 to D-3
 - software upgrade procedure F-1 to F-3
 - testing at startup 3-8
- 9600
 - configuring 3-5
 - system performance test 3-9

A

- abbreviations A-1
- alignment pin replacement 6-9
- analyzing data, instrument verification 5-13 to 5-14
- archiving, current spectra component files 4-2

B

- Background Component file, generating 4-4 to 4-5
- Blackboard test 3-12
- Blue test trays or paper, use to confirm signal throughput 3-11
- buffer hazards 5-3

C

- checklist
 - 7700 instrument preinstallation checklist D-2 to D-3
 - 7700 instrument, customer training D-5
 - TaqMan card upgrade D-3 to D-4
 - TaqMan card, customer training D-6
- chemical hazard warning 2-8
- chemical installation kit
 - guidelines 5-3
 - restart computer before using 3-7
- Chiller and Heater Tests, 9600 3-9
- computer
 - setup 3-6 to 3-7
- CT value criteria, adjusting for accuracy 5-14, 6-29

- customer supplied materials and equipment B-5
- customer training
 - 7700 instrument checklist D-5
 - for the 7700 SDS instrument 7-2
 - for the TaqMan card upgrade 7-3
 - TaqMan card checklist D-6
- Cytokine Cards. *See* TaqMan Human Cytokine Cards

E

- electrical configuration 3-3 to 3-5
 - checking electrical receptacle and operating voltage 3-3
 - configuring the thermal cycler in 50 hz areas 3-5
 - confirming internal electrical conditions 3-3 to 3-4
- exporting results 6-31
- exposure time, determining for data collection 6-22 to 6-24

F

- failed positions, further checking 5-15
- filling station, and vacuum pump assembly, setting up 6-10 to 6-11
- Fluorescent Test, signal throughput 3-11

I

- installation
 - before starting 1-2
 - completing the installation 7-4 to 7-5
 - computer and software setup 3-6 to 3-7
 - electrical configuration 3-3 to 3-5
 - checking electrical receptacle and operating voltage 3-3
 - confirming internal electrical conditions 3-3 to 3-4
 - confirming the thermal cycler in 50 hz areas 3-5
 - installation schedule 1-3
 - instrument installation specifications C-2
 - preparing and running RNase P install plate 5-9 to 5-12
 - setting up and starting β -actin install run 5-6 to 5-8
 - TaqMan card installation specifications C-3
 - testing at startup 3-8
- installation kit, 7700 instrument B-2 to B-3
- installation report
 - 7700 instrument D-9
- installation specifications, use to determine successful installation 5-15
- instrument performance, number of successful reactions needed 5-15, 6-32
- Instrument Verification test 3-8
- internal electrical connections, how to configure 3-3 to 3-4

L

- laboratory requirements 2-6 to 2-7
 - environmental conditions and requirements 2-6

laboratory space 2-6
pre-installation checklist 2-6
laboratory safety checklist
checklist
laboratory safety checklist D-7 to D-8
laser
alignment test 3-10
laser safety 2-5

M

manual, purpose of 1-2
materials and equipment
7700 instrument install kit B-2 to B-3
customer supplied materials and equipment B-5
TaqMan card upgrade kit B-4
TaqMan Human Cytokine Cards 6-4

N

No Template Controls, purpose of C-2

P

PCR reaction mix, preparing, and filling the card 6-12 to 6-18
performance, number of successful reactions
needed 5-15, 6-32
Power-Up tests, performing 3-8
preinstallation checklists
7700 instrument checklist D-2 to D-3
TaqMan card upgrade checklist D-3 to D-4
pure dye file, generating 4-6 to 4-11

Q

Quality Postcards, completing and returning 7-5

R

removing a TaqMan PCR tray, precaution to observe 5-13
report, installation
7700 instrument D-9
RNase P install plate, preparing and running 5-9 to 5-12

S

safety
general 2-2 to 2-8
general warnings 2-8
laser safety 2-5
Signal heights
criteria affecting 3-11
of the fibers, blackboard test 3-12
software
7700 instrument upgrade procedure F-1 to F-3
installation 3-6 to 3-7
spectra component files
archiving 4-2
spectral calibration
archiving current component files 4-3
generating Background Component file 4-4 to 4-5

generating pure dye file 4-6 to 4-11
materials 4-2
 β -actin install plate, preparing 5-3 to 5-5
Standard Deviation, recording 5-14
system noise test 3-12
system performance
distinguishing between 5,000 and 10,000
populations 7-4
test 3-9
verifying
analyzing the run 5-13 to 5-14

T

TaqMan Human Cytokine Cards
alignment pin replacement 6-9
analyzing the card and verifying system
performance 6-28 to 6-30
archiving current spectra component files 4-2
checklist D-3 to D-4
creating template E-1 to E-5
customer training checklist D-6
determining exposure time for data collection 6-22 to 6-24
exporting results 6-31
installation procedure described 6-2
installation specifications C-3
installing hardware and software 6-6 to 6-8
loading the card into the 7700 instrument 6-19 to 6-21
preparing PCR reaction mix and filling the card 6-12 to 6-18
setting up and starting a TaqMan card run 6-25 to 6-26
setting up filling station and pump assembly 6-10 to 6-11
TaqMan card verification calculations 6-32 to 6-33
troubleshooting 6-34 to 6-36
upgrade kit B-4
upgrade package contents 6-4
template, creating TaqMan Card template E-1 to E-5
tests
system
fluorescence and background checks 3-11 to 3-12
laser alignment test 3-10
testing at startup 3-8
thermal cycler tests 3-9
thermal cycler
configuring in 50 hz areas 3-5
system performance test 3-9
threshold, criteria for setting 5-14, 6-29

U

Unknown tubes with known problems, only reactions to be removed from analysis 7-4

V

vacuum pump assembly, setting up 6-10 to 6-11
vision hazard, laser source 3-10

voltage quality, about 2-7

W

warning, chemical hazard warnings 2-8

warranty, discussing with the customer 7-4

wavelength axis, use only as a guide 4-11

Headquarters

850 Lincoln Centre Drive
Foster City, CA 94404 USA
Phone: +1 650.638.5800
Toll Free: +1 800.345.5224
Fax: +1 650.638.5884

Worldwide Sales Offices

Applied Biosystems vast distribution and service network, composed of highly trained support and applications personnel, reaches into 150 countries on six continents. For international office locations, please call our local office or refer to our web site at www.appliedbiosystems.com.

www.appliedbiosystems.com



Applera Corporation is committed to providing the world's leading technology and information for life scientists. Applera Corporation consists of the Applied Biosystems and Celera Genomics businesses.

Printed in the USA, 06/2001
Part Number 4316192C

an **Applera** business