
Service Manual

TEM-Kombinationssystem
TEM Combination System

2232.201
2232.211 (USA)

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Important general instructions for use

Ensure that this product is only used as intended and described in the instruction manual by adequately trained and qualified personnel, and that maintenance and repair is only carried out by authorized specialized technicians.

Operate this product only in the combinations and with the accessories and spare parts listed in the instruction manual. Use other combinations, accessories and spare parts only if they are expressly intended for this use and if the performance and safety requirements are met.

Reprocess the products before every application and before returning them for repair as required by the instruction manual in order to protect the patient, user or third parties.


Subject to technical changes!

Due to continuous development of our products, illustrations and technical data may deviate slightly from the data in this manual.

CAUTION - USA only:

Federal law restricts this unit to be used or sold, except under the supervision of a medical doctor.

Safety instructions and levels of danger

Symbol	Level of danger
	WARNING! Failure to observe can result in death or severe injury.
	CAUTION! Failure to observe can result in slight injury or damage to the product.
	IMPORTANT! Failure to observe can result in damage to the product or surrounding.
	NOTE! Tips for optimum use and other useful information.

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ENGLISH

Contents

1	General information	1
1.1	Ordering spare parts	1
1.2	Disassembly of housing cover	2
2	Maintenance	3
2.1	Important notes	3
2.2	Maintenance of device and accessories	3
2.3	Visual check	4
2.4	Electrical safety test	4
2.5	Measuring device and auxiliary equipment	5
2.6	Function test	6
2.7	Troubleshooting	12
2.7.1	Error messages for Laparo CO2 Pneu automatic insufflator	12
2.7.2	Troubleshooting on Laparo CO2 Pneu	13
2.7.3	Error messages for TEM-PUMP	14
2.7.4	Troubleshooting on TEM-PUMP	15
3	Assemblies	16
3.1	Switching power pack 2232.041/.641 (item 150)	16
3.2	Control PCB 2232.041/.641 (item 200)	16
3.2.1	Disassembly and assembly	16
3.2.2	Functional checks	16
3.3	Control PCB 2232.041/.641 (item 210)	16
3.3.1	Measuring points and supply voltages	16
3.3.2	Disassembly and assembly of PCB	16
3.3.3	Function tests	16
3.4	Valve assembly 2232.041/.641 (item 160)	17
3.4.1	Disassembly and assembly	17
3.4.2	Function tests	17
3.5	Pressure reducing valve DM2 2232.041/.641 (item 170)	17
3.5.1	Function tests	18
3.6	Pressure reducing valve DM1 2232.244/.644 TEM-PUMP (item 230)	18
3.6.1	Function tests	18
4	Anhang	19
4.1	Reparaturteile / Spare parts Laparo-CO2-Pneu 2232.041 / 2232 .641 (USA)	19
4.1.1	Explosionszeichnung / Exploded View Laparo-CO2-Pneu	20
4.2	Reparaturteile / Spare parts TEM-Pump 2232.244 / 2232.644 (USA)	21
4.2.1	Explosionszeichnung / Exploded View TEM-Pump	22
4.3	Pläne / Diagrams Laparo CO2-Pneu 2232.041 / 2232.641 (USA)	23
4.3.1	Lageplan des Schaltnetzteils (Pos. 150) / Component drawing of switching power supply (item 150)	23
4.3.2	Lageplan der Bedien-E-Karte (Bestückungsseite) (Pos. 200) / Component drawing of Operation PCB (component side) (item 200)	24
4.3.3	Lageplan der Bedien-E-Karte (Lötseite) (Pos. 200) / Component drawing of Operation PCB (solder side) (item 200)	25
4.3.4	Lageplan der Steuer-E-Karte (Bestückungsseite) (Pos. 210) / Component drawing of Control PCB (component side) (item 210)	26
4.3.5	Lageplan Steuer-E-Karte (Lötseite) (Pos. 210) / Component drawing Control PCB (solder side) (item 210)	27

4.3.6	Lageplan der Schnittstellen-E-Karte (Pos. 140) / Component drawing of Interface PCB (item 140)	28
4.3.7	Blockschaltplan 2232.041 / Block diagram 2232.041	29
4.3.8	Blockschaltplan 2232.641 (USA) / Block diagram 2232.641 (USA)	30
4.3.9	Verdrahtungsplan 2232.041 /.641 / Wiring diagram 2232.041 /.641	31
4.3.10	Pneumatikplan 2232.041 / Pneumatic diagram 2232.041	32
4.3.11	Pneumatikplan 2232.641 (USA) / Pneumatic diagram 2232.641 (USA)	33
4.4	Pläne / Diagrams TEM-PUMP 2232.244 / 2232.644 (USA)	34
4.4.1	Lageplan der Steuer-E-Karte (Pos. 070) / Component drawing of Control PCB (item 070)	34
4.4.2	Lageplan der Schnittstellen-E-Karte (Pos. 030) / Component drawing of Interface PCB (item 030)	35
4.4.3	Blockschaltplan 2232.244 / Block diagram 2232.244	36
4.4.4	Blockschaltplan 2232.644 (USA) / Block diagram 2232.644 (USA)	37
4.4.5	Verdrahtungsplan 2232.244 /.644 / Wiring diagram 2232.244 /.644	38
4.4.6	Pneumatikplan 2232.244 / Pneumatic diagram 2232.244	39
4.4.7	Pneumatikplan 2232.644 (USA) / Pneumatic diagram 2232.644 (USA)	40
5	Protokolle / Reports	41
5.1	Prüfprotokoll / Test Report	41
5.2	Wartungsprotokoll / Maintenance report	42

1 General information

1.1 Ordering spare parts

The part number required for ordering a spare part is listed in the repair parts list under the item number used in the service manual.

 **IMPORTANT!**

Specify the following numbers when ordering spare parts:


- Part number of the spare part.
- Model/type number of the device.
- Serial number of the device.

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
Exploded view

1.2 Disassembly of housing cover

◆ Required tools: Phillips screwdriver.

- ◇ Unscrew the 4 housing feet (1.2), **fig.1**.
- ◇ Slide trim plates (1.3) from the rails (1.4), **fig.2**.
- ◇ Remove the 4 bottom Phillips screws (1.5), **fig.2**.
- ◇ Pull the housing cover (1.0) as indicated by the arrow (a) until the tongues (1.8) are disengaged from the corresponding grooves in the front plate (1.1), **fig. 3**.
- ◇ Lift off the housing cover (1.0) as indicated by the arrow (b), **fig. 3**.
- ◇ Disconnect the earthing/ground cable (1.9), **fig. 3**.

2 Maintenance

2.1 Important notes

This Service Manual describes the external service measures defined for the product.

IMPORTANT!

For carrying out the service measures the product instruction manual must be observed.

2.2 Maintenance of device and accessories

To protect the technician who carries out the test, the maintenance jobs and tests must be carried out in the order specified below.

- ◇ Visual check
- ◇ Electrical safety test
- ◇ Functional check

The test and measuring procedures specified in the chapter "Electrical Safety Test" relate to the test in accordance with EN / IEC 60601-1.

Alternatively, the electrical safety test can be carried out in accordance with DIN 751.

The operator shall bear responsibility for the limit values and collection / documentation of output values measured.

After repair all measurement values (specified output values) of the device must be checked as described in the Service Manual, and readjusted if deviations are observed.

CAUTION!

Do not use this product if the specified measuring values or functions are not fulfilled.

NOTE!

Any maintenance or testing job on the device or accessories must be documented.

2.3 Visual check

Item to be checked	Check for
Device and accessories	<ul style="list-style-type: none"> ◆ Safety-relevant soiling / contamination and general cleanliness. ◆ Mechanical damage ◆ Loose or missing parts
Controls	◆ Mechanical function and easy operation
Lettering / Symbols	◆ Completeness and legibility
Safety-relevant labeling (e.g warnings)	◆ Completeness and legibility
Device fuses	◆ The values specified by the manufacturer on the identification plate (Nominal current and melting characteristics)
Cables	◆ for perfect condition (position, insulation and brittleness)
PCBs/cards	◆ Corrosion or other damage
Protective earth connection	◆ Between LAPARO-CO ₂ PNEU and TEM-PUMP
	◆ Replace damaged parts immediately

2.4 Electrical safety test

Item to be checked	Check for
Protective earth connection to EN / IEC 60601-1	<p>Measured values:</p> <ul style="list-style-type: none"> ◆ without power/mains cable: $\leq 0.1 \text{ Ohm}$ ◆ with firmly attached power/mains cable: $\leq 0.2 \text{ Ohm}$ <p>Test conditions: $I_{\text{meas}} 25 \text{ A} \pm 10\%$, $V_0 \leq 6 \text{ V}$, $t_{\text{test}} 5 \text{ s to } 10 \text{ s}$, 50 Hz / 60 Hz</p> <p>To be tested: The resistance between protective earth contact or protective earth pin in the power plug and any other exposed (bare) metal part connected to the protective earth conductor.</p>
Leakage current to EN / IEC 60601-1	<p>Measured values:</p> <ul style="list-style-type: none"> ◆ earth leakage current: $\leq 170 \mu\text{A}$ limit value to EN / IEC 60601-1, normal condition: $\leq 500 \mu\text{A}$ ◆ patient leakage current: $\leq 10 \mu\text{A}$ limit value to EN / IEC 60601-1: normal condition for BF type: $\leq 100 \mu\text{A}$ normal condition for CF type: $\leq 10 \mu\text{A}$ <p>Test conditions: Measurement arrangement (MD) and measurement setup to EN/IEC 60601-1</p> <p>Measuring points: Insufflation connector, pump head; pressure measurement connector and irrigation bottle</p> <p>To be tested: The leakage current which can flow from each pole of the power supply unit through or via the insulation through the protective earth conductor or from the applied part via the patient to earth.</p>
Additional information	Disconnect the video connection before making the measurement.

2.5 Measuring device and auxiliary equipment

To carry out the maintenance, repair and test routines described in this service manual, the test devices and auxiliary equipment listed below are indispensable. Before you start a repair job, make sure that the measuring devices are properly calibrated and in perfect condition.

Designation	Additional information
DC power source	24V DC , no special current requirements
RMS multimeter	
CCD camera	with Y/C(S-Video)- and FBAS (BNC) output
Video monitor	with Y/C(S-Video)- and FBAS (BNC) input
Hygiene filter	4171.111
High-Flow insufflation tube	8170.232 (comprises connector 64202.147)
CO ₂ flow tube	Measuring range: 4-40 l /min, e.g. made by Platon Flow Management Jays Close, Viables, Basingstoke RG22 4 BS UK internet: www. platon.co.uk
CO ₂ flow tube	Measuring range: 1-10 l/min
Digital pressure gauge	Required resolution: 1 mm Hg (in the range below 0.5 bar) e.g.: "VDM 100" made by ASF, Munich, or "GDH 14 AN" made by Greisinger electronic, D-93128 Regenstauf, Germany
CO ₂ cylinder	0.75 kg (2046.961) For the tests the cylinder pressure must be in the "green range". The green LED of the cylinder pressure gauge is lit. If during the test the cylinder pressure indicator switches from green to yellow, you may still complete the test without changing the cylinder.
High pressure connection tube	8170.801
Test bag	must be able to withstand 30 mm Hg, such as an empty 3 litre infu- sion bag (should be as elastic as possible)
Syringe	50 ml
Y connector	for silicone tubes (see page 11)
Silicone tube	1.0 m, 4 x 2 mm
TEM tube set	8170.841
Pump tube, complete	8170.861
Hygiene filter	4171.112
Container	with 1 litre of water
Insufflation tube 2.5 m	8170.101
Connection tube for irrigation bottle	8170.862
Connection cable	2232.981
2 x connection tube	2232.851
Foot switch	2030.231

2.6 Function test

NOTE!

The section marked * relate to instruction manual GA-A 164.

Item concerned	Required checks
Preparation	Measuring setup as in section 2.5. Connect tube 2232.851 and connection cable 2232.981. Connect a sufficiently filled gas cylinder and open the valve. Do not connect the hygiene filter and pressure measuring tube (A) to the TEM-PUMP.
State after switch-on (Self test) Display software version	The self test takes less than 2 minutes. In case of a defect, 2 minutes are exceeded and all LEDs of the bar are triggered. The digital display for the actual flow indicates the current software version of the insufflator (2 digits). On the actual pressure digital display the current software of the TEM-PUMP is displayed (2 digits). A brief signal is sounded to indicate the end of the self test. TEM-PUMP: All LED indicators of the TEM-PUMP light up briefly. The stricture pump starts operating briefly and then stops. All buttons and the pump symbol are illuminated in night mode, all other indicators and displays are off. Upon successful self test: Laparo-CO₂ Pneu: Upon successful completion of the self test, the digital displays for gas volume consumed, flow and intraabdominal pressure read "0.0". There must be no error message on the display for gas volume consumed. The level indicators for pressure preselection and flow preselection must read the value preselected last. TEM-PUMP: Upon successful completion of the self test, all pushbuttons and the pump symbol are illuminated in night mode, the malfunction symbol is not illuminated. If a filter is plugged into the pressure measurement input, the level indicator reads the speed preselected last. If the self test could not be completed successfully, error messages E01 to E15 for the Laparo CO ₂ Pneu automatic insufflator or E30 to E41 for the TEM-PUMP appear on the digital display for the gas volume consumed.
Mode change	Note: The TEM combination system can now be configured for LAPAROSCOPY mode or TEM mode. To change the mode, the devices must be switched off for more than 30 seconds. To select LAPAROSCOPY mode, actuate the insufflation button on the laparo CO ₂ Pneu after a successful self test. After switching on insufflation, all LEDs on the TEM-PUMP must go out and remain in this state, until the TEM combination is switched off for more than 30 seconds and then switched on again.
LAPAROSCOPY mode	Connect hygiene filter with insufflation tube (B) to the patient output connector.


Item concerned	Required checks
Flow control	<p>Connect hygiene filter with insufflation tube to patient output connector. Preselect a nominal pressure of 25 mm Hg on the Laparo CO₂ Pneu. Actuate the "INSUFFLATION" button. The button must be lit. Read off the current flow on the digital display.</p> <p>Measurement values:</p> <ul style="list-style-type: none"> ◆ 1.1 l/min < Current flow < 1.5 l/min <p>Connect the CO₂ flow tube (Description in section 2.5) to the insufflation tube. Push the "HIGH FLOW" button. The button must be lit.</p> <p>Preselect a nominal flow of 15 l/min. Read off the current flow on the digital display and on the flow tube.</p> <p>Measurement values:</p> <ul style="list-style-type: none"> ◆ 14 l / min ≤ Current flow in acc. with digital display ≤ 16 l / min ◆ 12 l / min ≤ Current flow in acc. with flow tube ≤ 18 l / min <p>Preselect a nominal flow of 30 l/min. Read off the current flow on the digital display and on the flow tube.</p> <p>Measurement values:</p> <ul style="list-style-type: none"> ◆ 28 l/min ≤ Current flow in acc. with digital display ≤ 32 l/min ◆ 24 l/min ≤ Current flow in acc. with flow tube ≤ 33 l/min <p>Switch off insufflation.</p>
Consumption indicator	<p>Actuate the "RESET" button. The digital display for gas volume consumed reads "00.0".</p>
Pressure control	<p>Select a nominal pressure of 15 mm Hg and a nominal flow of 10 l/min. Connect the pressure measuring device/gauge (see section 2.5) and the insufflation tube to the bag representing the abdomen, using the Y-connector.</p> <p>First press the "INSUFFLATION" button and then the "HIGH FLOW" button. After a short time the preselected pressure of 15 mm Hg prevails in the bag. It is admissible that the pressure is exceeded briefly and then released.</p> <p>Tolerances:</p> <ul style="list-style-type: none"> ◆ Level indicators for nominal and actual pressure: ± 2 segments ◆ Digital display: ± 3 mm Hg <p>Switch off insufflation briefly.</p> <p>Preselect a nominal pressure of 20 mm Hg and push the "INSUFFLATION" button. After a short time a preselected pressure of 20 mm Hg prevails in the bag. It is admissible that the pressure is exceeded briefly and then released.</p> <p>Tolerances:</p> <ul style="list-style-type: none"> ◆ Level indicators for nominal and actual pressure: ± 2 segments ◆ Digital display: ± 3 mm Hg
Gas supply monitoring	<p>Press the "RESET" button. Close the gas cylinder valve.</p> <p>Select a nominal pressure of 20 mm Hg and a nominal flow of 5 l/min. Switch on "HIGH FLOW" insufflation.</p> <p>After a short time the gas supply indicator switches from continuous green light to continuous yellow light. The "Low Alarm"* is sounded. Then the gas supply indicator switches from continuous yellow light to a flashing yellow light. "Medium Alarm"* is sounded.</p>

Item concerned	Required checks
<p>Alarm, overpressure, overpressure switch</p>	<p>Plug the filter and the connector (64202.147) of the high flow insufflation tube 8170.232 onto the pressure connector (insufflation opening) on the front panel as shown in the figure below. Instead of the insufflation tube (which you should pull off the connector), connect the Y-connector described in section 2.5 via a 4 x 2 mm silicone tube. Connect the 50 ml syringe and the digital pressure measuring device/gauge to the other ends of the Y-connector.</p> <p>Preselect a nominal pressure of 15 mm Hg. Close the pressure release hole in the bottom plate of the housing and push the "INSUFFLATION" button. The pressure in the measuring setup rises to 15 mm Hg. If you exceed the nominal pressure by 4 mm Hg by actuating the syringe, "Medium Alarm"* is sounded.</p> <p>If the pressure exceeds 30 mm Hg in this measuring setup, "overpressure > 30 mm Hg" alarm* is triggered with a 5 second delay. In this case "High Alarm"* is sounded. The alarm indicator remains lit. On the digital display for the gas volume consumed error message E11 appears and the digital display for pressure blinks.</p> <p>If the pressure exceeds 45 mm Hg in this measurement setup, after a delay period the overpressure switch is triggered. This causes the "High Alarm" to switch to continuous sound, the alarm indicator remains lit and valve 2 switches off (switching noise).</p>
<p>Brief interruption</p>	<p>Switch off the device for approx. 10 seconds. After switching on all functions and settings previously selected must still be active.</p>
<p>Setup for overpressure simulation</p>	<p>see end of this section (function check)</p>

Item concerned	Required checks
TEM mode	<p>Switch off the device for 30 seconds. After switch-on, wait until the self-test is completed. Connect the pressure measuring tube (A)* with hygiene filter to the measuring input connector of the TEM-PUMP. Use the "PUMP ON/OFF" button of the TEM-PUMP to select TEM mode (HIGH FLOW button on laparo CO₂ pneu blinks). The insufflation gas flow in TEM mode is no longer cycled but continuous.</p>
Flow control	<p>On the laparo CO₂ pneu preselect a nominal pressure of 25 mm Hg and a flow of 1 l/min. Start insufflation. The "HIGH FLOW" button is lit. The digital display must read a current flow of 0.8 - 1.2 l/min. Switch off insufflation. Connect the CO₂ flow tube to the insufflation tube (B)*, with the outlet of the flow tube being left open. On the Laparo CO₂ pneu, preselect a nominal flow of 8 l/min. A flow of more than 8 l/min must not be preselectable.</p> <ul style="list-style-type: none"> ◆ ≥ 7 l/min current flow in acc. with digital display ≤ 9 l/min ◆ ≥ 6 l/min current flow in acc. with CO₂ flow tube ≤ 10 l/min <p>A constant flow must be visible in the flow tube. A starting of the roller pump can be tolerated in spite of the open system as it is energized at as little as 2 mm Hg.</p>
Pressure control	<p>Connect the pressure measuring device/gauge, insufflation tube (B)*, pressure measuring tube (A)* and the pump tube (suction D and Y)* to the bag representing the abdomen and to the devices. We recommend using a T-CONNECTOR to connect the tube set to the bag (type and shape depending on the structure of the abdomen model). On the Laparo CO₂ pneu preselect a nominal pressure of 15 mm Hg and a nominal flow of 4 l/min. Set the TEM-PUMP to minimum output and start. Start insufflation. In TEM mode the pressure measurement is effected via measuring tube (A)* of the TEM-PUMP, as the insufflator generates a constant flow (in laparoscopy mode the pressure measurement is effected directly in the insufflation flow). At 2-5 mm Hg the displacement pump must start operating. After a transient response phase, the preselected pressure of 15 mm Hg is reached.</p> <ul style="list-style-type: none"> ◆ Tolerance of the nominal and actual pressure level indicators: 15 mm Hg +/- 2 segments ◆ Tolerance of actual pressure digital display: 15 mm Hg +/- 3mmHg ◆ Tolerance of the pressure measuring device/gauge display: 15 mm Hg +/- 3 mmHg <p>Preselect a nominal pressure of 20 mm Hg and a nominal flow of 8 l/min. Set the TEM-PUMP to maximum output. Start insufflation. After a phase of transient response, the preselected pressure of 20 mm Hg is reached.</p> <ul style="list-style-type: none"> ◆ Tolerance of the nominal and actual pressure level indicators: 20 mm Hg +/- 2 segments ◆ Tolerance of actual pressure digital display: 20 mm Hg +/- 3 mmHg ◆ Tolerance of digital multimeter display: 20 mm Hg +/- 3 mmHg <p>Stop insufflation.</p>

Item concerned	Required checks
Displacement pump	<p>Use the same measuring setup as for pressure control (pressure gauge not required). On the TEM-PUMP, select minimum pump speed. Start displacement pump with the "PUMP ON/OFF" button. Start insufflation. At approx. 2-5 mm Hg the pump starts operating. Close the pump drain tube (Y)* approx. 30 cm after the displacement pump by kinking it. The displacement pump must now work against the overpressure in the tube; the pump must not stop. Let go of tube. Briefly actuate the foot switch. The pump speed must increase to its maximum value. When you press the foot switch again or actuate either the UP or the DOWN button, the pump must rotate at the preselected speed. Open the cover of the roller pump. The malfunction/fault indicator must light up. The pump stops and must not start again. Close the cover. The fault indicator goes out. Remove hose from roller pump. Use a tool to block the pump head through the hole in the front, then close the cover. Start the roller pump with the "PUMP ON/OFF" button. The fault indicator must blink. After approx. 5 seconds: The pump symbol also starts blinking, the roller pump switches off automatically, the fault indicator goes out, an alarm is sounded. It must not be possible to start the roller pump before the pump cover has been opened and closed again.</p>
Filter switch	<p>Measuring setup as above. Select a nominal pressure of 15 mm Hg and a nominal flow of 8 l/min on the insufflator. Press the "START INSUFFLATION" button. On the TEM-PUMP remove the hygiene filter from the measuring input connector (21). Error message "E 38" is displayed and an alarm is sounded. The insufflation pressure must drop to 8 mm Hg - 15 mm Hg. Switch off insufflation.</p>
Ventilation	<p>Use the same measuring setup as above. Preselect a nominal pressure of 20 mm Hg and start insufflation. When the nominal pressure is reached, exert an overpressure between 24 mm Hg and 30 mm Hg on the bag. The Laparo CO₂ pneu automatic insufflator must release the overpressure without issuing an alarm.</p>
Irrigation function	<p>The irrigation function is performed with an empty irrigation bottle. TEM-Pump: Connect the CO₂ 1-10 l/min flow tube to the irrigation output connector via irrigation tube (C)* (outlet of flow tube remains open). Actuate foot switch and keep it depressed (continuously). There must be a flow of 4 +/- 0.5 l/min. After approx. 5 seconds the flow is interrupted and an alarm is triggered.</p>

Setup for overpressure simulation:

2.7 Troubleshooting

The error messages are displayed on the “gas volume consumed” display.

2.7.1 Error messages for Laparo CO₂ Pneu automatic insufflator

Fault/error	Possible cause	Remedy
Error message E01	- Gas supply interrupted - Input pressure switch defective.	◆ Open gas cylinder or connect full gas cylinder ◆ In case of repeated error message with full, open gas cylinder: Contact service department
Error message E02	- Pushbutton depressed when switching on the device - If there is no operator's error: keyboard defective	◆ Repeat switch-on test ◆ Replace control PCB
Error message E03	- Internal power supply out of tolerance	
Error message E04	- Valve V1 defective; Flow measurement defective; V3 defective	
Error message E05	- Insufflation tube or hygiene filter connected when switching on the laparo CO ₂ pneu - Valve V2 defective	◆ Remove tube (hygiene filter) and switch off the device and switch on again after 20 seconds. ◆ Replace valve assembly
Error message E06	- Pressure measurement defective	
Error message E07	- Flow measurement defective	
Error message E08	- Valve V3 defective; V2 defective; Pressure measurement defective; Flow measurement defective	◆ Replace valve assembly
Error message E09	- Nominal value invalid	◆ Replace control PCB
Error message E10	- Helium	◆ Use CO ₂ gas
Error message E11	- Overpressure alarm: P _{act} greater than 30 mm Hg for more than 5 seconds	◆ Reduce intraabdominal pressure
Error message E12	- General hardware error; Pressure measurement defective; V1 defective; Flow measurement defective	
Error message E13	- “3 bar” or “Input pressure, 30 bar” pressure switch defective	◆ Replace both pressure switches as a precaution
Error message E14	- Pressure monitoring defective	
Error message E15	- DIL switch coding (position of electronic pressure monitoring fast/slow) incorrect; pressure switch defective	◆ Change DIL switch coding to suit the software and control PCB versions.

2.7.2 Troubleshooting on Laparo CO₂ Pneu

Fault/error	Possible cause	Remedy
Device without function	Power switch not switched on Power cable not connected Fuse in device defective No mains/line voltage	<ul style="list-style-type: none"> ◆ Switch on power switch ◆ Connect power cable ◆ Replace fuse ◆ Check in-house power supply
No gas flow	<ul style="list-style-type: none"> - Insufflation not switched on - Gas cylinder empty or closed 	<ul style="list-style-type: none"> ◆ Switch on insufflation ◆ Change or open the gas cylinder
Preselected pressure is exceeded	<ul style="list-style-type: none"> - Manipulations of the user increase the abdominal pressure - Laparo CO₂ Pneu defective - Other device causes pressure increase 	<ul style="list-style-type: none"> ◆ Eliminate cause ◆ Contact the service department ◆ Eliminate cause
Pressure won't build up	<ul style="list-style-type: none"> - Tube system leaking - Flow preselection too low - Instrument stopcock closed 	<ul style="list-style-type: none"> ◆ Check tubes and connectors ◆ Increase flow preselection ◆ Open instrument stopcock

2.7.3 Error messages for TEM-PUMP

Fault/error	Possible cause	Remedy
Error message E30	- Hygiene filter connected to the measuring input connector during switch-on - Microswitch at measuring input jammed	◆ Remove hygiene filter and switch off the power switch of the Laparo CO ₂ Pneu after 20 seconds; then switch on again ◆ Make sure pin moves easily or replace socket
Error message E31	- Pressure sensor for pressure measurement defective	◆ Contact the service department
Error message E32	- Pressure sensor for pressure measurement defective - Valve V2 defective	◆ Check pneumatic connections between Laparo CO ₂ Pneu and TEM-PUMP ◆ If error is not eliminated, contact the service department ◆ Adjust V2, replace valve assembly if necessary
Error message E33	- Current rectal value has exceeded limit value	◆ Reduce rectal pressure ◆ If fault is not eliminated, contact the service department
Error message E34	- Pushbutton depressed during switch-on	◆ Switch off the TEM-PUMP and then on again
Error message E35	- CO ₂ supply TEM-PUMP faulty	◆ Contact the service department
Error message E36	- not assigned	
Error message E37	- wrong pump tube - Roller pump defective	◆ Use original pump tube ◆ Contact the service department
Error message E38	- Emergency insufflation, measuring line with hygiene filter hasn't been connected to TEM-PUMP	◆ Connect measuring line with hygiene filter to TEM-PUMP ◆ If error is not eliminated, contact the service department
Error message E39	- Pressure switch defective	◆ Contact the service department
Error message E40	- Data transmission defective	◆ Switch off system and switch on again after 20 seconds ◆ If error is not eliminated, contact the service department ◆ Check connection of data lines
Error message E41	- Pneumatic error - TEM-PUMP leaking - Throttle in pneu defective	◆ Check tightness of measuring line for rectal pressure on rear panel of device ◆ If error cannot be eliminated, contact the service department

2.7.4 Troubleshooting on TEM-PUMP

Fault/error	Possible cause	Remedy
Device without function	- Power switch not switched on	◆ Actuate the power switch
Indicator lights and background illumination remain dark after switch-on	- Power cable not connected - Device fuse defective - No power supply/mains	◆ Connect power supply/mains connection ◆ Replace fuse ◆ Check mains/power supply
Suction power cannot be preselected	- Device defective	◆ Contact the service department or return the device for repair
Roller pump cannot be activated or switched on by means of "pump ON/OFF" button	- "Cover of pump housing open" indicator lamp is lit	◆ Close pump cover
Level indicator for "measured vacuum" does not show any reading although the pump is running	- Vacuum tube or hygiene filter incorrectly connected or defective - Other tubes incorrectly connected or defective	◆ Check tube and filter for tightness and correct connection, replace if necessary ◆ Check tube connections for tightness and correct connection, replace if necessary
No or poor suction	- Tubes incorrectly connected or defective - Tubes clogged - Instrument clogged	◆ Check tube connections for clogging, replace if necessary ◆ Check if tube connections clogged, replace if necessary ◆ Check the instrument for patency of the channel, replace if necessary
Alarm message of flow monitoring circuit despite of flow	- Flow sensor cable incorrectly connected or defective - Flow sensor defective	◆ Check cable for correct connection, replace if necessary ◆ Replace the flow sensor
No alarm message from flow monitoring with insufficient flow	- Flow sensor soiled from the inside - Flow sensor defective	◆ Reprocess or replace the flow sensor ◆ Replace flow sensor

3 Assemblies

3.1 Switching power pack 2232.041/.641 (item 150)

- ◇ Install preadjusted PCB and check for function in acc. with GA-A 164 and section 2.6 of this service manual.

3.2 Control PCB 2232.041/.641 (item 200)

3.2.1 Disassembly and assembly

- ◇ Remove the complete front frame and replace the PCB.

3.2.2 Functional checks

- ◇ Install the preset card and check for function (in acc. with GA-A 164 and section 2.6 of this service manual).

3.3 Control PCB 2232.041/.641 (item 210)

3.3.1 Measuring points and supply voltages

- ◇ Plug X1:
 - Pin 1 = PFD
 - Pin 2 = - 12 V \pm 0.4 V
 - Pin 3 = GND
 - Pin 4 = + 24 V
 - Pin 5 = GND
 - Pin 6 = + 12 V \pm 0.4 V
 - Pin 7 = GND
 - Pin 8 = + 5 V
 - Pin 9 = NC
 - Pin 10 = + 5 V
- ◇ Plug X3:
 - Pin 1 = DS 2 (low)
 - Pin 2 = DS 1 (high)
 - Pin 3 = 12 V \pm 0.4 V
- ◇ Plug X7:
 - Pin 1 = PFD
 - Pin 2 = - 12 V
 - Pin 3 = GND
 - Pin 4 = + 24 V
 - Pin 5 = GND
 - Pin 6 = + 12 V
 - Pin 7 = GND
 - Pin 8 = + 5 V
 - Pin 9 = NC
 - Pin 10 = + 5 V

3.3.2 Disassembly and assembly of PCB

- ◇ Remove the complete front frame and replace the PCB.

3.3.3 Function tests

- ◇ Install preadjusted PCB and check for function (in acc. with GA-A 164 and section 2.6 of this service manual).

3.4 Valve assembly 2232.041/.641 (item 160)

3.4.1 Disassembly and assembly

- ◇ Unscrew union nut and pull off output tube. Unscrew additional screw connections in the bottom plate from the bottom of the device. The valves cannot be replaced individually.


3.4.2 Function tests

- ◇ After replacement of assembly check all functions (in acc. with GA-A 164 and section 2.6 of this service manual).

3.5 Pressure reducing valve DM2 2232.041/.641 (item 170)

Test and adjustment:

- ◇ Laparo CO₂ pneu 2232 is switched off. Connect the device to the CO₂ cylinder using a high pressure tube.
- ◇ Open main valve of gas cylinder, to allow pressurized CO₂ to flow to the laparo CO₂ pneu.
- ◇ The valve assembly contains 3 valves (in the following order, from the front side of the device towards the rear: V2, V3, V1). Disconnect the two upper plug connectors on V1, the cable at the bottom (protective earth) remains connected.
- ◇ Disconnect the two 2 x 1.5 mm silicone tubes for the pressure sensors from V3 (connectors [a]) and connect the digital pressure gauge (see 2.5) to the two free connections using the Y-connector. The Y-connector is shown below:

- ◇ Disconnect the 4x2 mm silicone tube in front from connector [b] and close securely. Caution: if you do not close the tube the flow sensor may be damaged in the subsequent test.
- ◇ Apply 24-V-DC from a DC power supply to the free connectors of V1 (polarity is not important).
- ◇ Close the connector [b] in front with your finger.
- ◇ Read the pressure gauge. Turn the plastic screw on top of DM 2 until a stable pressure of 70 ± 1.5 mmHg ($\triangleq 93 \pm 2.0$ mbar) is achieved.
- ◇ Open and close the connector [b] in front several times. Check that the pressure is controlled correctly.
- ◇ Disconnect DC voltage from V1. Only then (!) should the closed tube to the flow sensor be opened again.
- ◇ Connect all connectors to their initial positions.
For connections [a] and [b] see section 4.1.1

3.5.1 Function tests

- ◇ After changing the assembly perform a complete check of the function in acc. GA-A 164 and section 2.6 of this service manual.

3.6 Pressure reducing valve DM1 2232.244/.644 TEM-PUMP (item 230)

Test and adjustment:

- ◇ Remove the cover of the TEM-PUMP. Place the pump at a right angle on the left side next to the insufflator.
 Establish TEM pressure and data connections, CO₂ and mains/power supply.
- ◇ Switch on the combination.
 Upon successful completion of the self test, remove the silicone tube from IC 27 and connect to the pressure measuring device/gauge.
- ◇ Check pressure and adjust to 1.0 bar if necessary.
- ◇ Reconnect silicone tube to IC 27.

3.6.1 Function tests



- ◇ Perform a complete check of the function in acc. GA-A 164 and section 2.6 of this service manual.

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RICHARD WOLF Belgium S.A.
RICHARD WOLF France S.A.R.L.