
Service Manual

1 CCD ENDOCAM

5512

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Important general instructions for use

Ensure that this product is only used as intended and described in the instruction manual by adequately trained and qualified personnel, and that maintenance and repair is only carried out by authorized specialized technicians.

Operate this product only in the combinations and with the accessories and spare parts listed in the instruction manual. Use other combinations, accessories and spare parts only if they are expressly intended for this use and if the performance and safety requirements are met.

Reprocess the products before every application and before returning them for repair as required by the instruction manual in order to protect the patient, user or third parties.


Subject to technical changes!

Due to continuous development of our products, illustrations and technical data may deviate slightly from the data in this manual.

CAUTION - USA only:

Federal law restricts this unit to be used or sold, except under the supervision of a medical doctor.

Safety instructions and levels of danger

Symbol	Level of danger
	WARNING! Failure to observe can result in death or severe injury.
	CAUTION! Failure to observe can result in slight injury or damage to the product.
	IMPORTANT! Failure to observe can result in damage to the product or surrounding.
	NOTE! Tips for optimum use and other useful information.

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ENGLISH

Contents

1	General information	1
1.1	Ordering spare parts	1
1.2	Disassembly of housing cover	2
2	Maintenance	3
2.1	Important notes	3
2.2	Maintenance of device and accessories	3
2.3	Visual check	4
2.4	Electrical safety test	4
2.5	Measuring devices and other equipment	5
2.6	Function tests	6
3	Assemblies	7
3.1	Control PCB (Item 0070)	7
3.1.1	Disassembly and assembly of control PCB	7
3.1.2	Adjustments	7
3.1.3	Function tests	7
3.2	Switching power pack (Item 0170)	7
3.2.1	Disassembly and assembly of assembly	7
3.2.2	Adjustments	7
3.2.3	Function tests	7
3.3	Main board	7
3.4	Camera head	7
4	Anhang / Annex	8
4.1	Reparaturteile / Repair Parts	8
4.2	Explosionszeichnung / Exploded View	9
4.3	Blockschaltplan / Block diagram	10
4.4	Verdrahtungsplan Kamera Controller/ Wiring Diagram Camera Controller	11
4.5	Verdrahtungsplan Kamerakabel / Wiring Diagram Camera Cable	12
5	Protokolle / Reports	13
5.1	Wartungsprotokoll / Maintenance report	13
5.2	Prüfprotokoll / Test report	14

1 General information

1.1 Ordering spare parts

The part number required for ordering a spare part is listed in the repair parts list under the position number used in the service manual.

 **IMPORTANT!**

Specify the following numbers when ordering spare parts:


- Part number of the spare part.
- Model/type number of the device.
- Serial number of the device.

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
Exploded view

1.2 Disassembly of housing cover

► Tools required: Phillips screwdriver

- ◇ 4 Unscrew device feet (1.2) **fig.1**.
- ◇ Slide trim plates (1.3) from holding rail (1.4), **fig.2**.
- ◇ Remove the lower 4 phillips screws (1.5), **fig.2**.
- ◇ Pull housing cover (1.0) in direction of arrow (a), until springs (1.8) project freely from front plate (1.1), **fig. 3**.
- ◇ Remove housing cover (1.0) in direction of arrow (b), **fig. 3**.
- ◇ Disconnect earthing cable (1.9), **fig. 3**.

2 Maintenance

2.1 Important notes

This Service Manual describes the external service measures defined for the product.

IMPORTANT!

For carrying out the service measures the product instruction manual must be observed.

2.2 Maintenance of device and accessories

To protect the technician who carries out the test, the maintenance jobs and tests must be carried out in the order specified below.

- ◇ Visual check
- ◇ Electrical safety check
- ◇ Functional check

The test and measuring procedures specified in the chapter "Electrical Safety Test" relate to the test in accordance with EN / IEC 60601-1.

Alternatively, the electrical safety test can be carried out in accordance with DIN 751.

The operator shall bear responsibility for the limit values and collection / documentation of initial values measured.

After repair all measurement values (specified output values) of the device must be tested as described in the Service Manual, and readjusted if deviations are observed.

CAUTION!

Do not use this product if the specified measuring values or functions are not fulfilled.

NOTE!

Any maintenance or testing job on the device or accessories must be documented.

2.3 Visual check

Item to be checked	Check for
Device and accessories	<ul style="list-style-type: none"> ◆ Safety-relevant soiling / contamination and general cleanliness. ◆ Mechanical damage ◆ Loose or missing parts
Controls	◆ Mechanical function and easy operation
Labeling / Symbols	◆ Completeness and legibility
Safety-relevant labeling (e.g warnings)	◆ Completeness and legibility
Device fuses	◆ the values specified by the manufacturer on the identification plate (Nominal current and melting characteristics)
Cables	◆ for perfect condition (position, insulation and brittleness)
PCBs/cards	◆ Corrosion or other damage
	◆ Replace damaged parts immediately


2.4 Electrical safety test

Item to be checked	Check for
Protective earth connection to EN / IEC 60601-1	<p>Measured values:</p> <ul style="list-style-type: none"> ◆ without power/mains cable: $\leq 0.1 \text{ Ohm}$ ◆ with firmly attached power/mains cable: $\leq 0.2 \text{ Ohm}$ <p>Test conditions: $I_{\text{meas}} 25 \text{ A} \pm 10\%$, $V_0 \leq 6 \text{ V}$, $t_{\text{test}} 5 \text{ s to } 10 \text{ s}$, $50 \text{ Hz / } 60 \text{ Hz}$</p> <p>To be tested: The resistance between protective earth contact or protective earth pin in the power plug and any other exposed (bare) metal part connected to the protective earth conductor.</p>
Leakage current to EN / IEC 60601-1	<p>Measured values:</p> <ul style="list-style-type: none"> ◆ earth leakage current: $\leq 100 \mu\text{A}$ limit value to EN / IEC 60601-1, normal condition: $\leq 500 \mu\text{A}$ ◆ patient leakage current: $\leq 10 \mu\text{A}$ limit value to EN / IEC 60601-1: normal condition for BF type: $\leq 100 \mu\text{A}$ normal condition for CF type: $\leq 10 \mu\text{A}$ <p>Test conditions: Measurement arrangement (MD) and measurement setup to EN/IEC 60601-1</p> <p>To be tested: The leakage current which can flow from each pole of the power supply unit through or via the insulation through the protective earth conductor or from the applied part via the patient to earth.</p>
Additional information	Disconnect the video connection cables between monitor, light projector and camera controller before making the measurement.

2.5 Measuring devices and other equipment

The measuring devices and other equipment listed below are compulsory for performing the maintenance, repair and test routines described in this Service Manual. Before any repair job, ensure that the measuring devices are properly calibrated and in perfect condition.

Designation	Additional data
Multimeter	
Monitor	with S-Video and video composite ("FBAS") connector (BNC).
Light source	with video control and light cable
Endoskop	∅ 10 mm
Remote control cable	5502.991
Objective lens	C-mount
BNC video cable	
S-VHS cable	

2.6 Function tests

Functions	Test procedure
	Switch on monitor. Connect 5512 with video and S-video cable to monitor (see manual). Do not connect the camera head.
LEDs on the device.	Switch on the 5512. All LEDs and buttons on the front panel will light up briefly. Then the following states are possible. The LED of the power/mains switch is lit. The remaining LEDs and buttons can be illuminated or dimmed (which is only visible in darkened environment).
Monitor	A color bar test image appears on the monitor. Check all video and S-video output signals (monitor image)
Video functions	Connect the camera head. Image white object (e.g. white piece of paper) with objective lens and telescope. Press the „white balance“ button on the controller. During „White balance“ the button is dimmed, then it remains lit. A brief acoustic signal is sounded to confirm that automatic white balance has started, a long acoustic signal is sounded to indicate the end. The object now appears white on the monitor. Use your cupped hand as test image. The result must be an image rendering the natural color of the skin. Press the button on the camera head for more than 1 second. A white balance is performed after the second sound signal. Successful white balance is acknowledged by a sound signal of approximately 1 second.
Shutter mode	Switch light source to manual mode, set to medium light intensity. Press the „Shutter function“ button. The button is illuminated. Image various scenes with different brightnesses. The camera should control brightness as a function of object brightness. Set the „Shutter“ button to SHUTTER FIX (button illumination is dimmed, shutter control deactivated). In the case of bright scenes the video image will bloom.
Image gain increase (GAIN)	Set light source to manual mode, medium light intensits. Set controller to shutter mode (button is illuminated, shutter control activated). Push the „GAIN“ button. The button is illuminated, a brief signal is sounded. Check the single-stage gain increase by example of a dark image. The increase in brightness must be clearly visible. A well-illuminated image must not show any brightness/gain increase.
Fixed aperture	This setting is only possible in conjunction with a light source with automatic video control. Switch light source to automatic mode and spot, at medium light intensity. Switch controller to fixed aperture (SHUTTER FIX). Image various scenes with different brightnesses. The control should compensate various illuminations/brightnesses.
Remote control	Connect remote control cable 5502.991 to „Remote“ socket. Use the multimeter to check continuity of the switch contact on the jack plug between the tip and the second ring, while the button on the camera head is depressed briefly (nominal value: < 100 Ohm)

3 Assemblies

3.1 Control PCB (Item 0070)

3.1.1 Disassembly and assembly of control PCB

- ◇ Remove the complete front frame and replace PCB.

3.1.2 Adjustments

- ◇ The PCB is delivered in preset condition. Not adjustments or settings are required.

3.1.3 Function tests

- ◇ See function check (section 2.6).

3.2 Switching power pack (Item 0170)

3.2.1 Disassembly and assembly of assembly

- ◇ Remove the complete switching power pack and replace.

3.2.2 Adjustments

- ◇ Voltage check
 - ◆ Input voltage 100 – 230 V AC
 - ◆ Output voltage 12V ± 0,2V DC

No adjustments are required if the assembly is replaced.

3.2.3 Function tests

- ◇ See function check (section 2.6).

3.3 Main board

 **NOTE!**

If there is a problem on the main board, please send in the complete camera with camera head to R.WOLF for repair.

3.4 Camera head

 **NOTE!**

In the case of problems /damage on the camera head we recommend sending in the complete unit (camera + camera head) to R.WOLF for repair and have the necessary adjustments carried out..

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RICHARD WOLF UK Ltd.
RICHARD WOLF Belgium S.A.
RICHARD WOLF France S.A.R.L.