

QA-ST User and Service Manual

QA-ST Safety Tester

METRON



North American Version 1.1

Copyright © 2000 by METRON. All rights reserved.

METRON:

USA

1345 Monroe NW, Suite 255A

Grand Rapids, MI 49505

Phone: (+1) 888 863-8766

Fax: (+1) 616 454-3350

E-mail: support.us@metron-biomed.com

FRANCE

30, rue Paul Claudel

91000 Evry, France

Phone: (+33) 1 6078 8899

Fax: (+33) 1 6078 6839

E-mail: info@metron.fr

NORWAY

Travbaneveien 1

N-7044 Trondheim, Norway

Phone: (+47) 7382 8500

Fax: (+47) 7391 7009

E-mail: support@metron.no

Disclaimer

METRON provides this publication as is without warranty of any kind, either express or implied, including but not limited to the implied warranties of merchantability or fitness for any particular purpose. Further, METRON reserves the right to revise this publication and to make changes from time to time to the content hereof, without obligation to METRON or its local representatives to notify any person of such revision or changes. Some jurisdictions do not allow disclaimers of expressed or implied warranties in certain transactions; therefore, this statement may not apply to you.

Limited Warranty

METRON warrants that the QA-ST Safety Tester will substantially conform to published specifications and to the documentation, provided that it is used for the purpose for which it was designed. METRON will, for a period of twelve (12) months from date of purchase, replace or repair any defective simulator, if the fault is due to a manufacturing defect. In no event will METRON or its local representatives be liable for direct, indirect, special, incidental, or consequential damages arising out of the use of or inability to use the QA-ST Safety Tester, even if advised of the possibility of such damages. METRON or its local representatives are not responsible for any costs, loss of profits, loss of data, or claims by third parties due to use of, or inability to use the QA-ST Safety Tester. Neither METRON nor its local representatives will accept, nor be bound by any other form of guarantee concerning the QA-ST Safety Tester other than this guarantee. Some jurisdictions do not allow disclaimers of expressed or implied warranties in certain transactions; therefore, this statement may not apply to you.

Table of Contents

1. INTRODUCTION	1-1
1.1 Features.....	1-1
1.2 Specifications.....	1-1
1.3 General Information.....	1-2
2. INSTALLATION	2-1
2.1 Receipt, Inspection and Return.....	2-1
3. OPERATING QA-ST	3-1
4. QA-ST SAFETY TESTS	4-1
4.1 Linearity Test of Voltmeter (Lead to Lead).....	4-1
4.2 Mains Supply Voltage.....	4-2
4.3 Ground Wire Resistance DC.....	4-3
4.4 Leakage Current – Internal Case.....	4-4
4.5 Leakage Current – External Case.....	4-5
4.6 Patient Leakage Current – Lead to Lead.....	4-6
4.7 Patient Leakage Current – Lead to Ground.....	4-7
4.8 Leakage Current – Lead to ISO.....	4-8
5. CALIBRATION	5-1
5.1 Linearity Test of Voltmeter (Lead to Lead).....	5-1
5.2 Mains Supply Voltage.....	5-2
5.3 Ground Wire Resistance DC.....	5-3
5.4 Leakage Current – Internal Case.....	5-4
5.5 Leakage Current – External Case.....	5-5
5.6 Patient Leakage Current – Lead to Lead.....	5-6
5.7 Patient Leakage Current – Lead to Ground.....	5-7
5.8 Leakage Current – Lead to ISO.....	5-8
APPENDIX A – CIRCUIT BOARD SCHEMATIC (NORTH AMERICAN VERSION)	A-1
APPENDIX B – ERROR REPORT FORM	B-1
APPENDIX C - IMPROVEMENT SUGGESTION FORM	C-1

This page intentionally left blank.

This page intentionally left blank.

1. Introduction

This chapter describes the METRON QA-ST Safety Tester's features and specifications.

1 1 F e a t u r e s

The QA-ST is designed to quickly check electrical safety on medical and industrial devices. It is palm-sized and very portable. The tests can be performed manually or in step-by-step mode. The results are displayed in a bright 4-digit LCD display.

The QA-ST performs all electrical safety tests necessary to satisfy AAMI ESI, NFPA-99, VDE 701, VDE 751 and HEI 95. It also features differential current measurements obtained by an internal current clamp. This presents a new leakage current measurement technique for convenient testing of medical equipment with a fixed power cord (hardwired to wall, not detachable).

1 2 S p e c i f i c a t i o n s 1

1. Line Voltage Measurement

Range: 0 - 300 V

Accuracy: $\pm 1\%$ reading, ± 1 digit

¹ For North American QA-ST Version 1.1

Resolution: 1 V L-N

2. Current Capacity

15 Amps @ 100 Volts

15 Amps @ 120 Volts

10 Amps @ 230 Volts

3. Resistance Measurement

DC-Source:

Range: 0 - 19.99 Ohms

Accuracy: $\pm 2\%$ reading, ± 1 digit

Resolution: 0.01 ohms

Current Source: 100 mA DC

AC-Source:

Range: 0 - 1999 mOhms

Accuracy: $\pm 2\%$ reading, ± 1 digit

Resolution: 1 mOhm

Current Source: 200 mA AC

4. Leakage Current

Range: 0 - 1999 μ A RMS

Accuracy: DC and 25 to 1.0 kHz, $\pm 2\%$ reading, ± 1 digit

1.0 kHz to 100 kHz, $\pm 2.5\%R$, ± 1 digit

100 kHz to 1 MHz, $\pm 5.0\%R$, ± 1 digit

Resolution: 1 μ A

5. Input Impedance

1000 ohm, AAMI ES1-1993 load or

IEC 60601.1 load

6. Power Consumption

110 V (50/60) Hz = 200 mA

115 V (50/60) Hz = 200 mA

240 V (50/60) Hz = 100 mA

.
1
. 3
G
e
n
e
r
a
l

I
n
f
o
r
m
a
t
i
o
n

1. Temperature Requirements:
 - +15/59 to +35/95 °C/F while operating
 - 0 to +50/122 °C/F for storage
2. Display:
 - Type: LCD
 - Numeric format: 4 characters
 - Display control: Keypad
3. Power: 100V/110V or 230V/240V, 50/60 Hz
4. Housing: Plastic case
5. Physical Dimensions (D x W x H):
 - 191mm x 135mm x 45mm
 - 7.5 in x 5.3 in x 1.8 in
6. Weight: 930 g / 2.0 lbs.

1
·
4
Q
A
-
S
T
O
r-
d
e
r-
i
n
g
l
i
n
e
f
o
r-
m
a
t
i
o
n

Order no:

11040: QA-ST Electrical Safety Analyzer
(Specify power supply socket.)

Standard Accessories:

11045: User/Service Manual QA-ST

17290: Test Lead (Kelvin Cable)

Various power supply socket including:

European Schuco, French Schuco, UK, Swiss,
Australian, US.

Optional Accessories:

11041: Carrying Case

2. Installation

This chapter explains unpacking, receipt inspection, claims, and the general procedures for QA-ST setup.

2

1

**R
e
c
e
i
p
t
,**

**I
n
s
p
e
c
t
i
o
n**

**a
n
d**

**R
e
t
u
r
n**

1. Inspect the outer box for physical damage.
2. Carefully unpack all items from the box and check to see that you have the following items:
 - QA-ST Electrical Safety Analyzer (P.N. 11040)
 - User/Service Manual QA-ST (P.N. 11045)
 - Test Lead (Kelvin Cable) (P.N. 17290)
 - Specified power supply socket
3. If you note physical damage, or if the unit fails to function according to specification, inform the supplier immediately. When METRON or the company's Sales Agent, is informed, measures

will be taken to either repair the unit or dispatch a replacement. The customer will not have to wait for a claim to be investigated by the Sales Agent. The Customer should issue a new purchase order to ensure delivery.

4. When returning an instrument to METRON, or the Sales Agent, fill out the address label, describe what is wrong with the instrument, and provide the model and serial numbers. If possible, use the original packaging material for return shipping. Otherwise, repack the unit using:
 - a reinforced cardboard box, strong enough to carry the weight of the unit.
 - at least 5 cm of shock-absorbing material around the unit.
 - nonabrasive dust-free material for the other parts.


Repack the unit in a manner to ensure that it cannot shift in the box during shipment.

METRON's product warranty is on page ii of this manual. The warranty does not cover freight charges. C.O.D. will not be accepted without authorization from METRON or its Sales Agent.

This page intentionally left blank.

3. Operating the QA-ST

This chapter describes the controls and terminals of the QA-ST Safety Tester.



NOTE
Remove the batteries and disconnect the AC Adapter / Battery Eliminator if you do not intend to use the PS-410 for an extended period of time.

1.

Inlet supply condition indicators:

Input Power Condition			Inlet Power Condition LEDs Lit		
Live	Neutral	Ground	L (Red)	G (Yellow)	N (Green)
L	N	G	X	X	X
L	N	NC	O	X	X
L	NC	G	X	O	X
L	G	N	X	O	X
N	L	G	O	X	X
N	G	L	X	O	X
N	L	NC	O	X	X
NC	L	G	X	O	X
N	NC	L	X	O	X
NC	G	L	X	O	X
G	L	N	X	O	X
G	N	L	X	O	X
G	L	NC	O	O	O
NC	L	N	X	O	X
G	NC	L	X	O	X
NC	N	L	X	O	X

L=Live N=Neutral G=Ground NC=Not connected X=ON O=Off

Binary display format:

LEAD NUMBER	DISPLAY X=ON O=OFF
1	OOOX
2	OOXO
3	OOXX
4	OXOO
5	OXOX
6	OXXO
7	OXXX
8	XOOO
9	XOOX
10	XOXO

4. QA-ST Safety Tests

This chapter describes QA-ST safety tests by function. If you are unfamiliar with the simulators' basic operation, refer to *paragraph 3.3, Operating QA-ST*.

NOTE: When the QA-ST is initially powered up, the unit will be in a default condition. At this time no measurements are being performed.

—
4
.
1
L
e
a
k
a
g
e
C
u
r-
r
e
n
t

(
L
e
a
d
t
o
L
e
a
d
)



The purpose of this test is to test for any leakage current between patient leads (applied parts).

1. Connect patient leads to the QA-ST starting with port 1, then 2 and so forth.

2. Press **"Patient Current Lead to Lead"**.




3. The result for Lead 1 will be displayed on the LCD display. To test Lead 2, press **"Patient Current Lead to Lead"**. Observe that the lead under test number is indicated in binary format by the LEDs in the **"Measurements"** area located below the LCD display.



4
2
M
a
i
n
s

S
u
p
p
l




y
V
o
l
t
a
g
e

Arrythmias: Premature Beats


- 4.
- 5.

The purpose of this test is to have the QA-ST display the Mains Supply Voltage coming into the QA-ST. This is measured between Hot and Neutral.

1. Press **“Supply V”**.
Supply
V Square Wave.
2. The value will be displayed on the LCD. The LED “ON” in the “Measurements” box will illuminate.



4
.
3
G
r
o
u
n
d
W
i
r
e
R
e
s
i
s
t
a



n
c
e
D
C
.

The purpose of this test is to determine if there is any DC resistance between the cabinet of the Device Under Test (DUT) and the Mains Supply Ground on the DUT.


1. Press "Ground Wire Res. DC". QA-ST will display over range (1).



2. Connect the Kelvin Test Lead to the cabinet of the DUT.
3. The resistance value will be displayed on the LCD.

4.4 Leakage Current – Internal Case

The purpose of this test is to determine if there is any leakage current between DUT Ground and Mains Supply Ground.



1. Plug in the DUT into the “Outlet Supply”.

2. Press “Leakage current Int. Case”.




3. Next, press “Outlet Supply RS” (Reversed Supply). The displayed QA-ST value may change.



4. Then, press “Outlet Supply OG” (Open Ground). The displayed QA-ST value may change.



4.5 Leakage Current – External Case



The purpose of this test is to determine if there is any leakage current between DUT Enclosure (chassis) and Mains Supply Ground.

1. Connect the Kelvin Clamp to the chassis of the DUT.

2. Press “Leakage current Ext. Case”.



3. The value will be displayed on the LCD. Note the value.
4. Next, press “Outlet Supply RS” (Reversed Supply). The displayed QA-ST value may change.



5. Then, press “Outlet Supply OG” (Open Ground). The displayed QA-ST value may change.



- 6.

4

6
P
a
t
i
e
n
t

L
e
a
k
a
g
e

C
u
r
r
e
n
t

L
e
a
d


t
o

L
e
a
d

The purpose of this test is to determine if there is any patient leakage current between DUT Enclosure (chassis) and Mains Supply Ground.

1. Connect all patient leads to the QA-ST starting with lead 1, then lead 2, and so forth.

2. Press "Patient leakage current – lead to lead".



3. The result will be displayed on the LCD. Note the value.
4. Press "lead to lead" to change the patient lead number. The diode panel will display, in binary format, the number of the actual patient lead that is under test.

**.
4
. 7
P
a
t
i
e
n
t

L
e
a
k
a
g
e

C
u
r
r
e
n
t

L
e
a
d


t
o

G
r
o
u
n
d**

The purpose of this test is to determine if there is any leakage current between all the patient leads in parallel, and DUT Ground.

1. Connect all patient leads to the QA-ST starting with lead 1, then lead 2, and so forth.

2. Press "Patient leakage current lead to GND" (Ground).




QA-ST

3. The value will be displayed on the LCD. Note the value

4. .

4.8 Leakage Current - Lead to ISO

The purpose of this test is to determine the quality of isolation between the DUT chassis and all of the patient leads in series.



QA-ST

1. Connect all patient leads to the QA-ST starting with lead 1, then lead 2, and so forth.
2. Connect the Kelvin Clamp to the chassis of the DUT.
NOTE: All patient leads are connected internally within the QA-ST to the same potential.

3. Press "lead ISO".



CAUTION: During the lead isolation test high voltage is present at the patient lead connections. While there is a resistor in the test circuit, making the risk of dangerous shock minimal, caution should be taken to avoid contact with the patient lead inputs or conductive applied parts during this test.

4. The results will be displayed on the LCD.

5. Calibration


This chapter describes the QA-ST calibration by function. If you are unfamiliar with the simulators' basic operation, refer to *paragraph 3.3, Operating QA-ST*.

**5
.
1
L
i
n
e
a
r-
it
y
T
e
s
t
o
f
V
o
l
t
m
e
t
e
r

(
L
e
a
d
t
o
L
e
a
d
)**

The purpose of this test is to calibrate the display/voltmeter.

2. Connect a voltage source in series with an ammeter and 10K Ω resistance as shown below.



6. Press “Patient Current Lead to Lead”.



7. Adjust the voltage source until the ammeter displays 1000μA. Then, adjust the QA-ST potentiometer R63 (see Appendix A for location) until the QA-ST LCD displays the same 1000μA value.

5
.
2
M
a
i
n
s

S
u
p
p
l
y

V
o
l
t
a
g
e

The purpose of this test is to calibrate the QA-ST to display the correct mains voltage.

3. Press “**Supply V**”.



4. Adjust the mains voltage with an external variable AC voltage (VARIAC) source. Then, adjust the QA-ST potentiometer R67 (see *Appendix A for location*) until the QA-ST LCD displays the same value.
5. Perform the test for all the given values in the test protocol under “Mains Supply Voltage: Reference”.


5
.
3
G
r
o
u
n
d
W
i
r
e
R
e
s
i
s
t
a
n
c
e
D
C

The purpose of this test is to determine if there is any DC resistance between the cabinet of the device under test, or DUT, and the Mains Earth on the DUT.

4. Press "Ground Wire Res. DC". QA-ST will display over range (1).



5. Connect a Kelvin Test Lead in series with an ammeter directly to DUT Earth. QA-ST will display $102\text{mA} \pm 2\%$.
6. Next, connect as shown in the below figure, without the ammeter and test leads.




7. Perform the test for tall the given values in the test protocol under "Ground Wire Resistance DC: Reference".

5.4 Leakage Current – Internal Case

The purpose of this test is to determine if there is any leakage current between DUT Earth and Mains Earth.

5. Connect a voltage source in series with an ammeter and 10K Ω resistance as shown below.



6. Press “Leakage current Int. Case”.



7. Adjust the voltage until the ammeter displays the first reference value in the test protocol. QA-ST will display the same value as the reference ammeter. Note the value.
8. Next, press “Outlet Supply RS” (Reversed Supply). The displayed QA-ST value will not change.



9. Then, press “Outlet Supply OG” (Open Ground). The displayed QA-ST value will not change.




- Perform the test for all the given values in the test protocol under “Leakage Current Int. Case: Reference”.

5.5 Leakage Current – External Case

The purpose of this test is to determine if there is any leakage current between DUT Enclosure (chassis) and Mains Earth.

- Connect a voltage source in series with an ammeter and 10KΩ resistance as shown below.



- Press “Leakage current Ext. Case”.



- Adjust the voltage until the ammeter displays the first reference value in the test protocol. QA-ST will display the same value as the reference ammeter. Note the value.
- Next, press “Outlet Supply RS” (Reversed Supply). The displayed QA-ST value will not change.



10. Then, press “Outlet Supply OG” (Open Ground). The displayed QA-ST value will not change.



7. Perform the test for tall the given values in the test protocol under “Leakage Current Ext. Case: Reference”.

5

6
P
a
t
i
e
n
t

L
e
a
k
a
g
e

C
u
r
r
e
n
t


L
e
a
d

t
o

L
e
a
d

The purpose of this test is to determine if there is any patient leakage current between DUT Enclosure (chassis) and Mains Earth.

5. Connect a voltage source in series with an ammeter and 10KΩ resistance as shown below.



6. Press “Patient leakage current – lead to lead”.



7. Adjust the voltage until the ammeter displays the first reference value in the test protocol. QA-ST will display the same value as the reference ammeter. Note the value.
8. Press “lead to lead” to change the patient lead number. The diode panel will display, in binary format, the number of the actual patient lead that is under test.
9. Perform the test for all the given values in the test protocol under “Patient Leakage Current Lead to Lead: Reference”.

**5
.
7
P
a
t
i
e
n
t

L
e
a
k
a
g
e

C
u
r
r
e
n
t


L
e
a
d

t
o

G
r
o
u
n
d**

The purpose of this test is to determine if there is any leakage current between all the patient leads in parallel, and DUT Earth.

5. Connect a voltage source in series with an ammeter and 10K Ω resistance as shown below.



6. Press "Patient leakage current lead to GND" (Ground).




7. Adjust the voltage source until the ammeter displays the first reference value in the test protocol. QA-ST will display the same value as the reference ammeter. Note the value.
8. Perform the test for all the given values in the test protocol under "Patient Leakage Current Lead to GND: Reference".

5.8 Leakage Current - Lead to ISO

The purpose of this test is to determine the quality of isolation between the DUT chassis and all of the patient leads in series.

5. Connect Kelvin test leads in series with an ammeter and a reference resistance to one of the patient leads (*see below*). Note that all patient leads are connected internally within the QA-ST to the same potential.



6. Press "lead ISO".



7. Adjust the voltage source until the ammeter displays the first reference value in the test protocol. QA-ST should display the same value as the reference ammeter. Note the value. If the ammeter value is not the same as the reference value, adjust the mains using potentiometer R67. (*See Appendix A*)

8. Perform the test for tall the given values in the test protocol under "Leakage Current to ISO: Reference".

8.

This page intentionally left blank.

USA1345 Monroe NW, Suite 255A
Grand Rapids, MI 49505

Phone: (+1) 888 863-8766

Fax: (+1) 616 454-3350

E-mail: support.us@metron-biomed.com**FRANCE**30, rue Paul Claudel
91000 Evry, France

Phone: (+33) 1 6078 8899

Fax: (+33) 1 6078 6839

E-mail: info@metron.fr**NORWAY**Travbaneveien 1
N-7044 Trondheim, Norway

Phone: (+47) 7382 8500

Fax: (+47) 7391 7009

E-mail: support@metron.no

From: (name) _____

Address: _____

Phone: _____

Fax: _____

E-mail: _____

Date: _____

Error Report

Product: _____

Version: _____

Serial no.: _____

Description of the situation prior to the error:

Description of the error:

(METRON AS internally)

Comments:

Received date:

Correction date:

Ref No.

 Critical Normal Minor

USA1345 Monroe NW, Suite 255A
Grand Rapids, MI 49505

Phone: (+1) 888 863-8766

Fax: (+1) 616 454-3350

E-mail: support.us@metron-biomed.com**FRANCE**30, rue Paul Claudel
91000 Evry, France

Phone: (+33) 1 6078 8899

Fax: (+33) 1 6078 6839

E-mail: info@metron.fr**NORWAY**Travbaneveien 1
N-7044 Trondheim, Norway

Phone: (+47) 7382 8500

Fax: (+47) 7391 7009

E-mail: support@metron.no

From: (name) _____

Address: _____

Phone: _____

Fax: _____

E-mail: _____

Date: _____

Improvement Suggestion

Product: _____

Version: _____

Description of the suggested improvement:

(METRON AS internally)

Comments:

Received date:

Correction date:

Ref No.

 Critical Normal Minor

USA

1345 Monroe NW, Suite 255A
Grand Rapids, MI 49505

Phone: (+1) 888 863-8766
Fax: (+1) 616 454-3350

E-mail: support.us@metron-biomed.com

FRANCE

30, rue Paul Claudel
91000 Evry, France

Phone: (+33) 1 6078 8899
Fax: (+33) 1 6078 6839

E-mail: info@metron.fr

NORWAY

Travbaneveien 1
N-7044 Trondheim, Norway

Phone: (+47) 7382 8500
Fax: (+47) 7391 7009

E-mail: support@metron.no

