

Instructions For Use

OASIS

Oxygen and Suction Integrated
System



702-0012.3

© Copyright

No part of this manual may be copied for any reason without the written permission of BPR Medical Ltd.

Manufacturer

BPR Medical Ltd, 22 Hamilton Way
Mansfield, NG18 5BU, United Kingdom

Tel: +44 (0) 1623 628281

Fax: +44 (0) 1623 628289

Approvals

OASIS is a CE Marked device in accordance with the Medical Device Directive (93/42/EEC) and meets the requirements of the following Harmonised Standards:

EN738-1Pressure Regulators and Pressure Regulators with Flowmetering Devices

EN10079-3Suction Equipment Powered from Vacuum or Pressure Source

Index

1.0 General Safety Warnings	3
2.0 Specific Safety Warnings Concerning Oxygen	4
3.0 Description of Symbols	5
4.0 Product Description	5
5.0 User Instructions	6
5.1 Connection of the Gas Pressure Regulator to a Cylinder	
5.2 Oxygen Delivery	
5.3 Suction Delivery	
5.4 Removing the Regulator from Cylinder	
5.5 Suction Jar Systems	
5.6 Suction Jar Compatibility	
6.0 Maintenance	9
6.1 Frequency of Maintenance	
6.2 Inspection	
6.3 Cleaning	
6.4 Leak Tests	
6.5 Flow Rate Test	
6.6 Suction Test	
7.0 Servicing	10
7.1 Frequency of Servicing	
7.2 Service Manual	
7.3 Decontamination Certificate	
8.0 Specification	11

8.0 Specification

General			
Operating, Storage and Transport Temperature Range		-20°C to 60°C	
Warranty		1 Year	
Flowmeter			
Nominal Flow Rate		0, 0.5, 1, 1.5, 2, 3, 4, 5, 6, 8,10, 15	
Flow Accuracy		±10% of setting at 1 lpm and above, ±20% of setting below 1 lpm across the full input range	
Effect on accuracy as a result of:			
i) varying inlet pressure		Typically less than 15% change across the input pressure range	
ii) varying outlet resistance		Less than 1% of reading up to 5kPa (50cmH ₂ O) back pressure	
iii) varying temperature		Less than 0.2% of reading per degree Celsius	
Delivery Hose		6.0 mm inside diameter	
Suction	Medium Suction Option (Blue Venturi Cartridge)	Medium/High Suction Option (Gold Venturi Cartridge)	High Suction Option (Silver Venturi Cartridge)
Vacuum Range	0 to –60 kPa (+5kPa/-30kPa)	0 to –55 kPa (+5kPa/-30kPa)	0 to –50 kPa (+5kPa/-30kPa)
Vacuum Flow	10 – 16 lpm	15-25 lpm	38 – 60 lpm
Gas consumption (typical)	20 lpm	30 lpm	85 lpm max
Recommended Suction Hose	7.5mm inside diameter	7.5mm inside diameter	7.5mm inside diameter
ISO EN 10079-3 Rating	Medium Vacuum Equipment	Medium Vacuum equipment	Medium Vacuum Equipment
Regulator			
(P ₁) Maximum Input Pressure [kPa]	20 000 kPa		
(P ₃) Minimum Input Pressure [kPa]	900 kPa		
Outlet Pressure [kPa] Stated at 10 000 kPa Input Pressure and flow of 1 lpm	400 kPa		
Minimum Flow Capacity [lpm]	100 lpm		
(Q ₁) Standard Discharge [lpm]	30		
(P ₂) Outlet Pressure [kPa] Stated at Input Pressure P ₃ and flow Q ₁	400		

7.0 Servicing

7.1 Frequency of Servicing

OASIS and its oxygen pressure regulator require servicing every 4 years. A Service Due date is given on the OASIS and pressure regulator rating label as a reminder.

7.2 Service Manual

A service manual is available for the OASIS Plus, which provides further information on spare parts and service procedures.

7.3 Decontamination Certificate

A valid decontamination certificate must accompany equipment returned for servicing. For reasons of Health & Safety, equipment returned for service without a valid decontamination certificate will be returned without further processing.

1.0 General Safety Warnings

IMPORTANT: Please ensure that all personnel who may use the OASIS have read this User Instruction Manual, are fully conversant with the Safety Warnings included within it and have received adequate general instruction on the clinical application of oxygen and suction equipment.

- W1 Before using this equipment read through the entire instruction manual. Attempting to use this equipment without an appreciation of its correct operation, its limitations and the general safety warnings associated with compressed oxygen may result in patient or user injury.
- W2 Oxygen therapy may be a critical treatment. A flowmeter should be used in strict accordance with the prescription and instructions of a qualified clinician. The effectiveness of supplemental oxygen therapy can only be determined by continuous monitoring of blood oxygen levels. It is essential that PaO2 or SpO2 monitoring is carried out when supplemental oxygen is prescribed.
- W3 Ensure that the oxygen cylinder is securely fastened to the trolley base.
- W4 When the flowmeter is turned off isolate the patient from the equipment by disconnecting the delivery tube from either the patient or the flowmeter.
- W5 Only appropriately trained service personnel working in controlled conditions must perform disassembly, assembly and testing of this equipment.
- W6 Do not submerge in water or allow any fluid to enter the equipment. If you have reason to suspect that fluid or other ingress has occurred remove the device from use and contact the manufacturer.
- W7 The performance of the device may be affected if it is stored or transported in temperatures outside of the range -20°C to +60°C.
- W8 Check the cylinder contents before use and at regular intervals during use as low cylinder pressure may result in poor or non-performance of the equipment. Always change the oxygen cylinder when the oxygen cylinder contents gauge is showing red.
- W9 Hydrophobic filters for re-useable suction jars and disposable liner systems are designed for single use only and must be replaced after every patient use. Single Use products should be disposed of in accordance with hospital protocols. Do not attempt to clean and re-use disposable filters or liners.
- W10 OASIS has been tested with the suction jar systems indicated in this manual. The use of suction jar systems other than those listed and/or approved by the manufacturer may result in poor or non-performance of the equipment.
- W11 The accuracy of the flowmeter and suction device may be affected if the input pressure is other than that stated in the specifications.

2.0 Specific Safety Warnings Concerning Oxygen

- W12 Oxygen is not flammable, but the presence of oxygen will drastically increase the rate and severity of combustion. Hydrocarbons such as oil or grease become highly combustible in the presence of oxygen. Oxygen must never be allowed to contact oil, grease or other petroleum-based substances.
- W13 Never administer oxygen while smoking or when near an open flame.
- W14 Oxygen cylinders have fill pressures up to 20 000 kPa (200 Bar). Never use medical oxygen from a cylinder without first reducing the pressure to a safe level through a suitable pressure regulator.
- W15 Ensure that the threaded fittings or indexing pins on the regulator are properly mated and tightened before opening the cylinder valve.
- W16 Never install a pin index regulator with more than one Yoke Seal between the cylinder and the regulator. Before attaching the regulator verify that the post valve is not already fitted with a Yoke Seal. Never add a second Yoke Seal to the regulator inlet while one is in place.
- W17 Make sure that the pin-indexed or threaded fittings on the regulator inlet are compatible with the gas cylinder on which it is to be fitted. Never attempt to force an incompatible connection.
- W18 The threaded bullnose connection on a UK medical oxygen cylinder may be the same as that for a UK medical air cylinder. Ensure that you have the correct cylinder for the application. OASIS and its regulator are intended for use with oxygen only and must not be connected to an air cylinder.
- W19 Before removing a regulator from a cylinder fully close the cylinder valve and release all gas from the regulator.
- W20 Never use oxygen as a pressure medium to purge obstructed pipelines or equipment, to operate pneumatic tools, or to build up any pressure in tanks.
- W21 Never permit compressed medical gases to enter a regulator suddenly. Always open the cylinder valve slowly. Do not stand in front of a regulator outlet when opening the cylinder valve.
- W22 Do not use or store oxygen equipment near excessive heat (>50°C or 125°F) or an open flame.
- W23 Use only lubricants recommended by the manufacturer when servicing the OASIS. The use of lubricants other than those recommended by the manufacturer may result in fire or explosion.
- W24 OASIS must always be used in an upright position. When fitting OASIS to a patient trolley ensure that the front membrane panel is vertical.

6.0 Maintenance

The following safety warnings apply. Please read them carefully before you proceed: W1, W5, W6, W7, W10, W11, W12, W14, W15, W16, W17, W18, W19, W20, W21, W22, W23 & W24

6.1 Frequency of Maintenance

Undertake Inspection every 6 months or when there is reason to believe that the device has sustained damage or may not be working correctly.

6.2 Inspection

Regulator Checklist:-

Input seal integrity (o-ring or dowty seal)

Cylinder contents gauge in good condition. Needle registers in the green area of the gauge on a full cylinder and zero when not connected to a cylinder

Sintered bronze filter in place and not blocked by debris

Labels intact and legible

Service due in the next 6 months?

OASIS Manifold Checklist:

Vacuum gauge in good condition

Low Pressure indicator shows when system pressure is removed and disappears when system pressure is restored

High pressure hose between regulator and manifold is in good condition.

OASIS is securely fixed to the trolley, no movement

Service due in the next 6 months?

6.3 Cleaning

Use a mild disinfectant soap solution and hot water to clean external parts in accordance with hospital procedures.

6.4 Leak Tests

Flowmeter Leak Through - Connect and charge the system to a full cylinder in accordance with Section 5. Check that the oxygen and suction control knobs are at zero (fully clockwise) and the OASIS low-pressure indicator is not showing red. Connect one end of a length of hose to the flowmeter outlet barb and place the other end into a container of water. Check that there are no bubbles coming from the end of the tube, which would indicate leak through.

External Leak – Following on from the test above, turn off the cylinder and leave the system charged for 1 hour. At the end of this time ensure that the system pressure hasn't dropped such that the low-pressure indicator has come forward and is now showing red.

Note:- If the low pressure indicator does come forward, re-charge the system and check each joint for leaks using a proprietary leak detector spray until the leak is found.

6.5 Flow Rate Test

Verify oxygen flow rates at all flow settings against specification (i.e. within $\pm 10\%$ at 1 l/min and above or $\pm 20\%$ below 1 l/min).

6.6 Suction Test

Verify maximum suction capability by occluding the suction outlet, turning the control knob fully anticlockwise and checking the suction gauge.

3.0 Description of Symbols



Attention, see
Instructions for Use



Use no oil



Ascending diagonal line on the decal
label indicates increase in flow.



CE Marking

4.0 Product Description

OASIS is an oxygen flowmeter and a suction regulator integrated into one device, driven by a single oxygen cylinder.

It is designed to provide controllable delivery of supplemental oxygen and suction in emergency situations or when patients are being moved on patient trolleys.

OASIS is intended to be fitted to the head end of patient trolleys to provide easy access to both oxygen and suction where and when it is needed most. The integrated, sturdy design provides protection from damage in a harsh environment.

The oxygen flowmeter is of fixed orifice type, providing eleven pre-set flow rates selectable by turning the control knob to the required position. The flow rate at each position is shown through a window in the body of the flowmeter and is given in units of litres per minute (lpm). The patient or outlet connection barb is situated underneath the OASIS directly below the flowmeter control knob.

The suction device comprises a variable vacuum source to which a proprietary suction jar system can be connected. Vacuum is generated using the venturi principle. On turning the control knob, oxygen is fed to the venturi generator, which in turn provides a vacuum at the suction barb. The amount of vacuum is indicated by the vacuum gauge indicated in units of minus Kilo Pascals (-kPa) and minus Bar (-Bar). Vacuum is only noticeable on the gauge when an obstruction is perceived by the system, occluding the output port or hose with your finger is a simple way of testing system performance.

To provide early indication that a cylinder is empty or has been turned off, a low pressure indicator is included. The indicator, situated on the front panel, shows red when there is no gas present. A pressure regulator is provided with the OASIS to reduce the high cylinder pressure to a safe working level. The regulator is permanently connected to the OASIS manifold.

Suction Jar Systems

	Suction Jar Size		
	0.8 / 1.0 Litre	1.2 / 1.5 Litre	2.0 Litre
Abbott			
Receptal	Re-useable suction jar for disposable liner Disposable liner Atcupal cradle type	W064 L212 A	E423 F806 F806
Allegiance			
Medi-Vac Guardian	Disposable canister Atcupal cradle type	1212CE B	4220CE
Medi-Vac Flex Advantage	Re-useable suction jar for disposable liner Disposable liner Atcupal cradle type	65652-611 65651-910 B	65652-616 65651-920 B
Medi-Vac CRD	Re-useable suction jar for disposable liner Disposable liner Atcupal cradle type	65652-511 1510CE B	65652-516 1515CE B
Tyco			
Sep-T-Vac	Disposable canister Atcupal cradle type	8888-310003 B	8888-310062 B
	Re-useable suction jar for disposable liner Disposable liner Atcupal cradle type	8888-310839 8888-310821 B	8888-310843 8888-310805
VacSax			
	Re-useable suction jar for disposable liner Disposable liner Atcupal cradle type	3833-132 9910-340 B	3833-090 9910-208 B
Flowmeter Spa			
Monokit Autoclavable	Autoclavable suction jar Atcupal cradle type	MAK1000 A	
Monokit Disposable	Re-useable suction jar for disposable liner Disposable liner Atcupal cradle type	5113080 B	

5.0 User Instructions

The following safety warnings apply. Please read them carefully before you proceed: W1, W2, W3, W4, W6, W7, W8, W9, W10, W11, W12, W13, W14, W15, W16, W17, W18, W19, W20, W21, W22, W23 & W24

5.1 Connection of the Gas Pressure Regulator to a Cylinder

Check that the cylinder type and regulator inlet connector are compatible.

Check the presence and condition of the input connector seal on the regulator.

Note: For pin index regulators this is a washer type seal (Yoke Seal) and for bull-nose regulators an o-ring.

Fit the regulator to the cylinder ensuring that it is properly mated and the connectors are done up. Hand tight is sufficient.

Open the cylinder valve with a suitable cylinder key. Check the regulator cylinder contents gauge to verify that there is sufficient gas available.

Check that the 'Low Pressure' indicator on the front face of the OASIS manifold does not show red (i.e. indicator fully forwards).

5.2 Oxygen Delivery

Connect a sufficient length of oxygen tubing (not supplied) to the outlet barb. Connect the other end of the oxygen tubing to the patient or patient enclosure using the appropriate connector (not supplied).

Determine the appropriate oxygen flow rate in l/min required for the patient. Turn the control knob anticlockwise until the required rate is visible in the viewing window.

Continue to monitor the patient as advised by the clinician. Continue to check the cylinder contents gauge on the regulator.

Upon completion of the therapy turn the control knob to zero and disconnect the delivery hose from the patient.

5.3 Suction Delivery

Prepare the suction jar system in accordance with the manufacturer's instructions.

If you are intending to use an autoclavable suction jar system ensure that a hydrophobic filter and overflow protection device is fitted.

Note:- Autoclavable suction jar systems have single use hydrophobic filter/overflow protection devices, which should be replaced with each new use/patient. The disposable liner systems recommended for OASIS (refer appropriate table) have an integral hydrophobic filter and overflow protection device.

Connect the suction hose (supplied) between the suction barb on the Acupal manifold and the barb on the suction jar marked 'Vacuum'.

Connect suction hosing and wand (not supplied) to the barb on the suction jar marked 'Patient'. This tubing is used for removing fluids from the patient and will deposit into the suction jar.

Note:- If these hoses are reversed the suction device will not operate correctly.

Fit the suction jar into the cradle provided on the bottom of the trolley.

To activate the suction unit turn the control knob anticlockwise. The level of suction (vacuum pressure) can be adjusted by turning the control knob as required.

5.4 Removing the Regulator from the Cylinder

To remove the regulator from the cylinder, first ensure the cylinder valve is closed. Bleed off any residual gas pressure in the system by turning on either the oxygen or suction delivery systems until flow ceases. Undo and remove the regulator from the cylinder.

5.5 Suction Jar Systems

Suction jar systems can generally be broken down into autoclavable and disposable types.

Autoclavable systems have a rigid suction jar with a sealable lid. Following cleaning and sterilisation they may be reused. A disposable hydrophobic filter is fitted to such systems and should be replaced after every use.

Disposable systems comprise either a rigid suction jar into which a disposable liner is fitted (the rigid jar being reusable) or a rigid or semi-rigid canister requiring no other components and being wholly disposable. The disposable systems recommended for use with the OASIS (refer appropriate table) have integrated hydrophobic filters and overflow protection devices and no other provision need be made in this regard.

Components from single use or disposable systems should never be re-used and no attempt should be made to clean or sterilise them.

5.6 Suction Jar Compatibility

OASIS has been tested and found to be compatible with the suction jar systems detailed in the table provided.

Some 2 Litre size jars shown in the table, whilst compatible with the OASIS suction regulator, are too tall and are liable to be broken should the trolley be lowered or put into the Trendelenberg position. These suction jar systems are therefore considered incompatible with OASIS and are shown with the manufacturer's part number crossed out (e.g. ~~XXX XXXX~~). For the same reason all 3 Litre and larger suction jar systems are incompatible with OASIS.