

APC ErgoTouch Pro

Airborne Particle Counter

ETP-1.1 en
April 2010

Operation and Service Manual
(English)



APC ErgoTouch Pro

Airborne Particle Counter

*ETP-1.1 en
April 2010*

Operation and Service Manual (English)

SHIP TO/MAIL TO:

Biotest AG
Technology Center Microbiology
Wernher von Braun Straße 9
D-69214 Eppelheim
Germany

Technical Support:

+49 (0) 6221 726 5130

Fax:

+49 (0) 6221 726 5139

E-mail address:

Hycon@biotest.de

Website:

<http://www.biotest.com>

Local Technical Support:

USA

Biotest Diagnostics Corporation
400 Commons Way, Suite F
Rockaway, NJ 07866 USA
Tel: +1 973 625 1300
Toll Free: +1 800 522 0090
Fax: +1 973 625 5882

Email: customerservice@biotestusa.com

Japan

Biotest K.K.
Daini Yasuda Building, 2F
3-32-13, Tsuruya-cho, Kanagawa-ku,
Yokohama-shi Kanagawa-ken
221-0835 Japan
Tel: +81 45 317 5111
Fax: +81 45 31 5110

Benelux

Biotest Seralc° NV
Generaal De Wittelaan 17B bus 14
2800 Mechelen
Belgium
Tel: +32 (0) 15 28 50 50
Fax: +32 (0) 15 28 50 55
Email: info@biotest.be



Manual History

The following is a manual history of the  **APC** ErgoTouch Pro Airborne Particle Counter:

Revision	Date
A	November 2008
ETP-1.0 en	January 2009
ETP-1.0 de	January 2009
ETP-1.1 en	April 2010

Warranty

Revision	ETP-1.1 en / April 2010
Copyright	©Biotest AG / 2010 / All rights reserved.
Address	Biotest AG, Technology Center Microbiology, Wernher von Braun Straße 9, 69214 Eppelheim, Germany
E-mail Address	HYCON@biotest.de
Limitation of Warranty and Liability	<p>Seller warrants the goods sold hereunder, under normal use and service as described in the operator's manual, shall be free from defects in workmanship and material for (24) months, or the length of time specified in the operator's manual, from the date of shipment to the customer. This warranty period is inclusive of any statutory warranty. This limited warranty is subject to the following exclusions:</p> <p>a. Parts repaired or replaced as a result of repair services are warranted to be free from defects in workmanship and material, under normal use, for 90 days from the date of shipment.</p> <p>b. Seller does not provide any warranty on finished goods manufactured by others or on any fuses, batteries or other consumable materials. Only the original manufacturer's warranty applies.</p> <p>c. Unless specifically authorized in a separate writing by Seller, Seller makes no warranty with respect to, and shall have no liability in connection with, goods which are incorporated into other products or equipment, or which are modified by any person other than Seller.</p> <p>The foregoing is IN LIEU OF all other warranties and is subject to the LIMITATIONS stated herein. NO OTHER EXPRESS OR IMPLIED WARRANTY OF FITNESS FOR PARTICULAR PURPOSE OR MERCHANTABILITY IS MADE.</p> <p>TO THE EXTENT PERMITTED BY LAW, THE EXCLUSIVE REMEDY OF THE USER OR BUYER, AND THE LIMIT OF SELLER'S LIABILITY FOR ANY AND ALL LOSSES, INJURIES, OR DAMAGES CONCERNING THE GOODS (INCLUDING CLAIMS BASED ON CONTRACT, NEGLIGENCE, TORT, STRICT LIABILITY OR OTHERWISE) SHALL BE THE RETURN OF GOODS TO SELLER AND THE REFUND OF THE PURCHASE PRICE, OR, AT THE OPTION OF SELLER, THE REPAIR OR REPLACEMENT OF THE GOODS. IN NO EVENT SHALL SELLER BE LIABLE FOR ANY SPECIAL, CONSEQUENTIAL OR INCIDENTAL DAMAGES. SELLER SHALL NOT BE RESPONSIBLE FOR INSTALLATION, DISMANTLING OR REINSTALLATION COSTS OR CHARGES.</p> <p>No action, regardless of form, may be brought against Seller more than 12 months after a cause of action has accrued. The goods returned under warranty to Seller's factory shall be at Buyer's risk of loss, and will be returned, if at all, at Seller's risk of loss.</p> <p>Buyer and all users are deemed to have accepted this LIMITATION OF WARRANTY AND LIABILITY, which contains the complete and exclusive limited warranty of Seller. This LIMITATION OF WARRANTY AND LIABILITY may not be amended, modified or its terms waived, except by writing signed by an Officer of Seller.</p>
Service Policy	<p>Knowing that inoperative or defective instruments are as detrimental to Biotest as they are to our customers, our service policy is designed to give prompt attention to any problems. If any malfunction is discovered, please contact your nearest sales office or representative, or call Biotest's Customer Service department at +49 (0) 6103 801 496</p>
Trademarks	<p>APC, ERGOTOUCH, ERGOTOUCH PRO and ERGOTRAK are trademarks of Biotest AG. Biotest and the Biotest logo are trademarks of Biotest AG. Microsoft and Excel are registered trademarks of Microsoft Corporation.</p>


Contents

Manual History	ii
Warranty	iii
Safety Information	vii
Laser Safety	vii
Labels.....	viii
Description of Caution / Warning Symbols	viii
Caution and Warning Symbols.....	ix
Chapter 1 Introduction and Unpacking	1
Optional Accessories.....	2
Chapter 2 Getting Started	4
Using the Instrument Stand and Stylus	4
Providing Power	5
To Install the Lithium-ion Battery.....	5
To Use AC Power	6
Performing a Daily Zero Check	6
To Perform a Zero Check.....	6
Installing an Isokinetic Inlet.....	7
Installing a Temperature / Relative Humidity Probe	7
Chapter 3 Operation	8
Screen Layout and Functionality	8
Main Tab	9
Zoomed Data Screen	11
Setup Tab.....	12
Software Input Panel (Keyboard).....	12
Data Setup Screen.....	13
Count Units Screen.....	13
Clear Samples Screen	13
System Setup Screen	14
Change Power On Password Screen	14
Change Setup Password Screen.....	15
System Configuration Screen	15
Print Settings Screen	16
Print Schedule Screen	17
Device Setup Screen	17
Date and Time Screen	18
Display Screen.....	18
Diagnostics Screen.....	19
Sampling Setup Screen	20
Channels Screen	20
Sample Timing Screen.....	21
Count Mode Screen	22

Environment Screen	23
Locations Screen	23
Alarms Screen	24
Recipe Screen	25
Data Tab	26
Reports Tab	28
Room Definition Screen	29
Generate Screen.....	30
Chapter 4 Data Handling	31
System Requirements	31
Installation	31
ERGOTRAK Lite Software Configuration	35
Data Presentation	36
Data Export	37
File Management	38
Device Connection	39
Working with Connected Devices.....	40
Displaying Records	41
Device Related Application Menus.....	42
Main Application Menu	42
Delete All Records	42
Export all Records.....	43
Export Selected Records	43
Recipes	44
Locations.....	46
File Menu.....	47
Device Menu	47
Help Menu	48
Troubleshooting.....	48
Chapter 5 Maintenance	49
Maintenance Schedule	49
Recommended Field Maintenance Schedule	49
Daily Zero Check.....	49
Cleaning the Instrument Enclosure	49
Chapter 6 Troubleshooting	50
Chapter 7 Contacting Customer Service	52
Technical Contacts.....	52
Returning the APC ERGOTOUCH PRO Airborne Particle Counter for Service	52
APPENDIX A Specifications	53
Temperature / RH Probe (optional accessory)	54

Safety Information

This section gives instructions to promote safe and proper handling of the

 handheld airborne particle counters.

IMPORTANT

There are no user-serviceable parts inside the instrument. Refer all repair and maintenance to a qualified factory-authorized technician. All maintenance and repair information in this manual is included for use by a qualified factory-authorized technician.

Any attempt to open or service the instrument by non-authorized agencies or personnel will void the warranty.

Laser Safety

The APC ERGOTOUCH PRO Airborne Particle Counter is a handheld Class I laser- based instrument.

- During normal operation, the user will **not** be exposed to laser radiation.
- Precaution should be taken to avoid exposure to hazardous radiation in the form of intense, focused, visible light.
- Exposure to this light may cause blindness.

Take these precautions:

- **DO NOT** remove any parts from the particle counter unless you are specifically told to do so in this manual.
- **DO NOT** remove the housing or covers. There are no user-serviceable components inside the housing.








WARNING

The use of controls, adjustments, or procedures other than those specified in this manual may result in exposure to hazardous optical radiation.

Labels

Advisory labels and identification labels are attached to the outside of the particle counter housing and to the optics housing on the inside of the instrument:

1. Serial Number Label (back panel)	 9306-01BT S/N: 93060849001    COMPLIES WITH 21 CFR 1040.10 AND 1040.11 U.S. Patent 6167107 Wernher von Braun Strabe 9, D-69214 Eppelheim, Germany HYCONService@biotest.de
2. Laser Radiation Label (internal)	<div style="border: 1px solid black; padding: 10px; text-align: center;"> <p>DANGER!</p> <p>VISIBLE LASER RADIATION WHEN OPEN. AVOID DIRECT EXPOSURE TO BEAM</p> <p>WARNING: NO USER SERVICABLE PARTS INSIDE. REFER SERVICING TO QUALIFIED PERSONNEL</p> </div>
3. European symbol for non-disposable item. Item must be recycled.	

Description of Caution / Warning Symbols

Appropriate caution/warning statements are used throughout the manual and on the instrument that require you to take cautionary measures when working with the instrument:

Caution



C a u t i o n

Failure to follow the procedures prescribed in this manual might result in irreparable equipment damage. Important information about the operation and maintenance of this instrument is included in this manual.

Warning







W A R N I N G

Warning means that unsafe use of the instrument could result in serious injury to you or cause damage to the instrument. Follow the procedures prescribed.

Caution and Warning Symbols

The following symbols may accompany cautions and warnings to indicate the nature and consequences of hazards:

	Warns that uninsulated voltage within the instrument may have sufficient magnitude to cause electric shock. Therefore, it is dangerous to make contact with any part inside the instrument.
	Warns that the instrument contains a laser and that important information about its safe operation and maintenance is included in the manual.
	Warns that the instrument is susceptible to electro-static dissipation (ESD) and ESD protection procedures should be followed to avoid damage.
	Indicates the connector is connected to earth ground and cabinet ground.

Getting Help

To obtain assistance with this product or to submit suggestions, please contact Biotest Customer Service, Biotest Technical Center Microbiology, or your local Biotest representative or service center:

Biotest AG
Technology Center Microbiology
Wernher von Braun Straße 9
69214 Eppelheim
Germany
Telephone: +49 (0) 6221 726 5130

Fax: +49 (0) 6221 726 5139

E-mail Address: HYCON@biotest.de
Web site: www.biotest.com

Biotest AG
Customer Service Microbiology
Landsteinerstraße 5
63303 Dreieich
Germany
Telephone: +49 (0) 6103 801 496

Fax: +49 (0) 6103 801 505

Chapter 1

Introduction and Unpacking

The APC ERGOTOUCH PRO Airborne Particle Counter is a lightweight, handheld particle counter with a touchscreen interface. It operates on the included lithium-ion battery or AC power.

The device has a 0.1 ft³ / min (2.83 L / min) flow rate and counts bin sizes from 0.3 to 10 µm. Up to 10,000 data sets can be downloaded for analysis and reporting using the ERGOTRAK LITE Data Download Software included with the device.




Typical applications for this particle counter include clean room monitoring, research, exposure assessment, indoor air quality, filter testing, clearance testing, quality assurance, and contaminant migration studies. All APC ERGOTOUCH airborne particle counters meet the JIS B 9921: 1997 standard.

Unpacking the **APC** ErgoTouch Pro Airborne Particle Counter

Carefully unpack the APC ERGOTOUCH PRO airborne particle counter from the shipping container and verify that all the items shown in the photos below and listed in the following tables are present. Contact your local Biotest representative immediately if items are missing or broken:

APC ERGOTOUCH PRO Airborne Particle Counter Standard Included Items





Qty.	Item Description	Part/Model	Reference Picture
1	APC ERGOTOUCH PRO Airborne Particle Counter	942 550	
1	Power Cord	US, UK, EURO (included)	
1	AC Power Adapter	942 562 (included)	
1	Isokinetic Inlet	942 602 (included)	
1	Battery Pack	942 561 (included)	






Qty.	Item Description	Part/Model	Reference Picture
1	Computer Cable (2 m), USB A to B	191 551 (included)	
1	Stylus	(included)	
1	Purge Filter Assembly (HEPA)	942 601 (included)	
1	ERGOTRAK LITE data download utility CD (includes Operation and Service Manual)	(included)	(not shown)
1	Operation and Service Manual	(included)	(installed on CD)
1	Calibration Certificate	(included)	(not shown)
1	Quick Start Guide	(included)	(installed on CD)

Optional Accessories

The following photos and table list optional accessories. If you ordered optional accessories, make certain they have been received and are in working order:

APC ERGOTOUCH PRO Airborne Particle Counter Optional Accessories

Item Description	Part/Model	Reference Picture
External Battery Charger w/ AC Adapter and Power Cord	942 563	
External Printer	942 610	
Custom Carrying Case	942 570	
Printer Paper (10 rolls)	942 611	(not shown)
Temperature / Relative Humidity Probe	942 565	

Item Description	Part/Model	Reference Picture
Isokinetic Inlet SS	942 602	
Isokinetic Probe SS (remote sampling)	942 603	
Barbed Inlet (remote sampling)	191 824	
Tubing, 30 m / 100 ft (remote sampling)	191 822	
APC Compressed Gas Adapter	942 675	
APC Compressed Gas Adapter Nozzle Set (ErgoTouch)	942 680	(not shown)

Chapter 2

Getting Started

This chapter provides information to help you use the APC ERGOTOUCH PRO Airborne Particle Counter including:

- using the instrument stand and stylus
- providing power
- performing a daily check
- installing an isokinetic inlet
- installing a Temperature / Relative Humidity Probe (Item # 942 565)

Using the Instrument Stand and Stylus

The APC ERGOTOUCH PRO is equipped with an integral instrument support stand. To open the stand, grasp it by the large finger hole and pull it out until it locks into place. Be careful not to overextend the stand. To store the stand out of the way when not in use, simply push the stand back until it snaps into place.



The APC ERGOTOUCH PRO is also equipped with a plastic stylus for use with the touchscreen interface. The stylus locks into place in the case near the top of the unit when not in use.



Providing Power

The APC ERGO TOUCH PRO may be powered using a removable rechargeable lithium-ion battery, or through an AC power cord.

Notes:

1. When using AC power, the battery (if installed) charges when the instrument is on, but not while actively sampling.
2. Removing / changing the lithium-ion battery or disconnecting the AC power does not cause loss of data.
3. A new battery is installed as part of routine service, if necessary.

To Install the Lithium-ion Battery

1. Remove the battery cover from the back of the instrument by lightly depressing the textured tab on the cover located on the lower left.
2. Slide the lithium-ion battery into the slot, press down lightly and slide it forward (toward the top of the unit) until it locks into place.
3. Replace the battery cover and slide it in place until you hear a click.



W A R N I N G


The battery supplied by Biotest has built in protection against explosion and fire hazard. Do **not** use a substitute.



W A R N I N G

Do **not** use non-rechargeable batteries in this instrument. Fire, explosions, or other hazards may result.

To Use AC Power

1. Connect the AC power adapter to the country-specific power cord.
2. Insert the AC power adapter into the bottom of the APC ERGOTOUCH PRO.
3. Connect the power cord to an outlet.
4. Press the on/off button  (located in the center of the front of the instrument).
5. After a splash screen displays the Biotest logo, a brief start-up sequence begins as the Windows® CE operating system boots up.




Performing a Daily Zero Check

A zero check should be performed at least once a day. It should also be performed before conducting any important testing or certification.

To Perform a Zero Check

Turn on the instrument and wait until the main menu appears.

1. Remove the isokinetic inlet if attached. The daily zero check cannot be performed when the isokinetic inlet is attached to the instrument.
2. Attach the zero filter to the inlet nozzle located on the top of the instrument.
3. Press the **Start** button  or the blue triangular key above the **On/Off** key and allow the instrument to purge for 2 minutes.
4. After the 2-minute purge, continue to sample. In accordance with JIS standards, there should be no more than 1 particle counted at any size in 5 minutes.



Note: If the instrument does not go to zero (1 particle in 5 min is considered zero), refer to Chapter 6, [Troubleshooting](#), for additional information.

5. Remove the zero filter and put the isokinetic inlet back on; the instrument is now ready for operation.

Installing an Isokinetic Inlet

The isokinetic inlet smoothly accelerates air into the inlet of the instrument. To install, simply thread the inlet directly onto the inlet nozzle until finger tight. The inlet seals over an o-ring so it doesn't have to be very tight to seal.



Installing a Temperature / Relative Humidity Probe

To install the optional Temperature / Relative Humidity Probe (Item # 942 565):

1. Align the small red dot at the base of the probe to the corresponding red dot on the socket.
2. Press the probe into the socket until it clicks.
3. Temperature and relative humidity are automatically displayed in the upper left corner.
4. Remove the probe by pulling straight up.





Chapter 3

Operation

The APC ERGOTOUCH PRO Airborne Particle Counter is controlled using a touchscreen display. Use the plastic stylus or your finger tip. **DO NOT** use sharp objects (such as a pen point) that may damage the screen overlay.

To turn on the instrument, press the **on/off** button (located in the center of the front of the instrument). After a splash screen displays the Biotest logo, a brief start-up sequence begins as the Windows® CE operating system boots up.

The instrument is ready for operation when the main tab (shown right) appears. If an optional Temperature / Relative Humidity Probe (Item # 942 565) is attached, those values will be shown in the upper-left corner also.

		11/20/08 11:43:03 AM	
		Automatic	
µm	Σ	Location	
0.3	0		
0.5	0	Delay: 00:00:05	
1.0	0	Count: 0/10	
3.0	0	Time: 00:01:00	
5.0	0	Hold: 00:00:10	
10.0	0	Recs: 123/10000	
Main		Setup	Data Reports

Screen Layout and Functionality

There are four main screens (tabs): Main, Setup, Data, and Reports. The operation of each of these screens, the information displayed on them, and the operations you can perform from each are described in the remainder of this chapter.

Some screens require or allow you to enter information. To enter information, tap on the touchscreen and an on-screen keyboard appears.

Main Tab



The **Main Tab** is the default screen. The left side of the screen summarizes the concentrations for the currently selected location. Tap on the size and count portion of the screen to enable **Zoomed Data** screen (see [Setup Tab](#)).









The display shows:

- Temperature*
- Relative humidity*
- Bin sizes
- Particle count/concentration



*Temperature and humidity are displayed only if the optional T/RH probe is installed.



The status bar at the top of the screen shows the current time and date settings (see the [Setup Tab](#)) and indicates:

		 11/20/08 11:43:03 AM
		Automatic 
µm	Σ	Location
0.3	0	<div>▼</div> Delay: 00:00:05 Count: 0/10 Time: 00:01:00 Hold: 00:00:10 Recs: 123/10000
0.5	0	
1.0	0	
3.0	0	
5.0	0	
10.0	0	
Main	Setup	Data Reports

Icon	Description
	Laser requires service
	Sufficient flow through the APC ERGOTOUCH PRO
	Insufficient flow through the APC ERGOTOUCH PRO
	Operating on AC power, no battery installed
	Operating on AC power, battery is installed and charging
	Battery charged
	Low battery
	Battery must be charged

The right side of the **Main Tab** shows locations and other information (delay, cycles, and so on). These can be configured using the Setup Tab.

Sampling		11/20/08 11:45:57 AM	
23.9 °C 52 %RH		 Sampling Automatic 	
µm	Σ	Location	
0.3	107623	Loc001	
0.5	8300	Delay: 00:00:00	
1.0	613	Count: 1/10	
3.0	48	Time: 00:00:57	
5.0	14	Hold: 00:00:10	
10.0	6	Recs: 1/10000	
Main		Setup	Data Reports

Field	Description
Location	Use this dropdown box to display information about any of the available locations.
Delay	The initial delay between the time the Start/Stop button is pressed and the instrument begins sampling. Valid only when Automatic mode is selected.
Count	The number of samples that have been taken/the total number of samples. Valid only when Automatic mode is selected.
Time	The time for each sample. Valid only when Automatic mode is selected.
Hold	The time between samples. Valid only when Automatic mode is selected.
Recs	The total number of records in the database /10,000 (maximum number of records).
Manual/Automatic/Beep	Mode indicator refers to the “Data Count Mode” (see section below).
	Press the Start/Stop button to begin sampling in the configured mode. Start/Stop may also be entered using the triangle-shaped blue button above the power button on the front of the instrument.
	Prints the current sample to the optional printer.



Zoomed Data Screen



The **Zoomed Data** screen is entered by touching in the size and count part of the **Main Tab** display. The bottom portion of the screen summarizes the concentrations for the currently selected location. Tap the size and count portion of the display to switch back to the **Main Tab** display.

The display shows:

- Temperature *
- Relative humidity *
- Bin sizes
- Particle count/concentration

* Temperature and Humidity are displayed only if the optional T/H probe is installed.


Sampling			11/18/08 02:25:07 PM
25.0 °C		Unknown	
45 %RH			
µm	Δ	Σ	
0.5	5597	5671	
1.0	17	74	
3.0	37	57	
5.0	12	20	
10.0	8	8	
Main	Setup	Data	Reports

Field	Description
	Prints the current sample to the optional printer.
Location	Label that indicates the location selected when the sample was taken.
	Press the Start/Stop button the begin sampling in the configured mode.

Setup Tab


The setup tab provides access to the following:

Data Setup	View Count Units and Clear Samples .
System Setup	Change Power On Password, Setup Password, System, Print Settings, and Print Schedule.
Device Setup	Set Date and Time, Display and Diagnostics.
Sampling Setup	Set up particle Channels, Sample Timing, Alarms, Count Mode, Locations, and Environment.
Recipes	Save a group of settings (recipes) that you use over and over so you don't have to reset individual settings.




Software Input Panel (Keyboard)

Throughout the setup screens, a keyboard will appear on the screen when text may be entered. Data may be entered using this keyboard. When the entry is complete, press either the ↵ (**Enter**) or **Esc** keys. The keyboard will then be hidden until another text entry box is selected.



Data Setup Screen


This screen lets you access the **Count Units** screen and the **Clear Samples** screen:



Count Units Screen


This screen lets you set the way in which particle concentration information is displayed:

Field	Description
Differential	Select to display particle count data as a differential Δ (count data indicates the number of particles within a size range from the current channel to the next larger one).
Cumulative	Select to display particle count data as cumulative Σ (the total number of counts includes all particles larger than the indicated channel size).
Concentration	Display concentration in ft^3 or m^3 . If Beep mode is selected, display of concentration values is not allowed.




Clear Samples Screen

The **Clear Samples** screen lets you clear all samples from the internal database. Select **Yes** to clear all samples. Select **No** to return to the **Data Setup** screen:



System Setup Screen

From the **System Setup** screen you can select (or change) a password for using the unit (**Power On Password**), set up a password (**Setup Password**), select system configuration parameters (**System**), configure the printer (**Print Settings**), and schedule print jobs (**Print Schedule**):



Change Power On Password Screen

If a **Power On Password** has been previously set, you must enter that password before being allowed to change the **Power On Password**. If a **Power On Password** is set, then on instrument startup a password screen will ask for the password before the instrument can be used. A blank password is regarded as no password and if set as the new password, will not prompt you for a password on system startup.

Tap on the screen to display the on-screen keyboard and enter the required information.

Note: Keep the password in a safe place. It is very difficult to reset the password and requires contacting the factory.

Field	Description
Old Password	Enter your existing password (if one has already been set).
New Password	Enter a new password. The password can be any length and use any characters.
Confirm New Password	Retype the new password then press OK. A confirmation message appears if the password is changed.

Change Setup Password Screen

If a **Setup Password** has been previously set, you must enter that password before being allowed to change the **Setup Password**. If a **Setup Password** is set, clicking on the **Setup** tab at the bottom of the main screen brings up a password screen. That password must be entered in order to change instrument settings.


Tap on the screen to display the on-screen keyboard and enter the required information.

Change Setup
10/03/08
6:14:29 AM

Old Password

New Password

Confirm New Password



OK

Cancel

Note: Entering a blank password will turn off password protection.

Field	Description
Old Password	Enter your existing password. (if one has already been set).
New Password	Enter a new password The password can be any length and use any characters.
Confirm New Password	Retype the new password then press OK. A confirmation message appears if the password is changed.


System Configuration Screen

Use this screen to set system configuration parameters. Press **OK** when finished:

Field	Description
Δ and Σ on Zoom	Select to zoom in on both cumulative (Σ) and differential (Δ) counts on the Main Tab . To zoom the Main Tab , click on the left side of the Main Tab . (It takes a moment for the screen to update.) Click on the screen again to return to normal view.

System
10/24/08
7:40:33 AM

☐ **Δ and Σ on Zoom**





Defaults

OK

Cancel

Print Settings Screen

A hard copy of a sample set or statistics can be printed from the instrument using an optional thermal printer. Use this screen to set print parameters. Press **OK** when finished.

Print Settings			10/24/08 7:41:42 AM
<input checked="" type="checkbox"/> Serial Number			
<input checked="" type="checkbox"/> Model Name			
<input checked="" type="checkbox"/> Separator			
<input type="checkbox"/> Differential			
<input checked="" type="checkbox"/> Cumulative			
		Defaults	
		OK	
		Cancel	


Field	Description
Serial Number	Indicates that the serial number of the particle counter used to collect the data will be printed.
Model Name	Indicates that the model number of the particle counter used to collect the data will be printed.
Separator	Indicates a line separator will be printed after the Model Name and Serial Number in the header of all printouts
Differential	Indicates that the differential value of the data will be printed.
Cumulative	Indicates that the cumulative value of the data will be printed.

Note: Printer paper has a colored strip printed on the last meter of each roll to indicate time to change the paper roll.


Print Schedule Screen

Use this screen to schedule automatic printing. You can choose to either print when an alarm occurs or print whenever a sample is complete:

Field	Description
Automatic Printing	Enables automatic printing
On Alarm	Print data when an alarm condition occurs.
On Sample	Print data whenever a sample completes.

Print Schedule  10/24/08 7:42:27 AM


☒ **Automatic Printing**
☒ **On Sample**
☐ **On Alarm**







Defaults
OK
Cancel

Device Setup Screen

Use this screen to access screens that let you set or change the date and time, make audio visual selections, and run diagnostics:

Device Setup  10/22/08 1:17:41 PM

 **Date and Time**  **Display**  **Diagnostics**

 **Back**


Date and Time Screen

This screen lets you set the current date and time and set the date format. Press **OK** when finished. You can select options using the arrows or tapping on the screen.

Field	Description
Date	Press the down arrow to display a calendar then select the date from the calendar.
Time	Select the time component you want to change (hours; minutes; seconds) and then use the left and right arrows to adjust to the current time.
Date Format	Highlight the date format you want to use from the list.
24 Hour	Time display is in 24 hour format.


Display Screen

This screen lets you set or change visual parameters

Field	Description
Screen Alignment 	Press this item to reset the screen alignment, and follow the directions on the alignment screen.


Diagnostics Screen

This screen lets you view the system's model, serial number, copyright, manufacture date, calibration date, next calibration date, firmware version, and USB IP address. Press **Close** when finished.




Sampling Setup Screen

Use this screen to access screens that let you set up how sampling is displayed and handled. You can select which channels to use, the sample timing, the count mode, environment, sampling locations, and alarm thresholds.



Channels Screen

This screen lets you choose the channels that are enabled and set their particle size. Press **OK** when finished.



Field	Description
Enable	Select the channels you want to view on the main display.
Size	This box allows for changing the default size for any channel. Highlight the size information and use the on-screen keyboard to change its value. Channels cannot be set below 0.3 or above 10.0 µm and they may not overlap one another.

Sample Timing Screen

This screen lets you select parameters for sampling. Use the up and down arrows or the on-screen keyboard to change or enter information. These parameters are only valid when the APC ERGOTOUCH PRO is running in **Automatic** mode. Press **OK** when finished.

Sample Timing 10/20/08 2:50:04 PM

Count: 10

Delay: 00 : 00 : 04

Hold: 00 : 00 : 10

☒ **Time:** 00 : 01 : 00

☐ **Volume:** 0

Cubic Feet

Cubic Meter

Defaults


OK

Cancel

Field	Description
Count	Count is the total number of samples you want collected. In Automatic mode, a Count value of 0 will cause the instrument to count continuously using the settings for Delay , Time , and Hold until the Start/Stop button is pressed again. Use the up and down arrows or the on-screen keyboard to set the count.
Delay	Delay indicates how long it will be before the first sample is taken. It takes approximately 6 seconds for the pump to reach the flow set point; taking a measurement before the pump is functioning properly may result in a data error. Highlight the time component you want to change (hours, minutes, seconds) and use the up and down arrows or the on-screen keyboard to change the value.
Hold	Hold indicates how long the instrument pauses between samples. Highlight the time component you want to change (hours, minutes, seconds) and use the up and down arrows or the on-screen keyboard to change the value.
Time	Time indicates the duration of each sample run (count particles). Highlight the time component you want to change (hours, minutes, seconds) and use the up and down arrows or the on-screen keyboard to change the value.
Volume	Volume sets the volume of air that will pass through the instrument for each sample. If you select volume, you must select Cubic Feet , Cubic Meters or Liters for measurement using the arrows.

Count Mode Screen


Use this screen to set the sample count mode. Press **OK** when finished.



Field	Description
Automatic	If you select this mode, the APC ERGOTOUCH PRO starts counting in Automatic mode when you press the Start/Stop button according to the setting on the Sample Timing screen.
Manual	If you select this mode, the APC ERGOTOUCH PRO starts sampling immediately when you press the Start/Stop button and stops at the end of the sample time, which is configured on the Sample Timing screen.
Beep	If you select this mode, the APC ERGOTOUCH PRO starts sampling data immediately and beeps whenever the threshold for the smallest bin is reached, as specified in Alarms screen. This can be very useful when searching for leaks, especially around filters. When this mode is selected the particle count data is displayed in total counts (not concentration).

Environment Screen

Use this screen to set the units for temperature, which is displayed on the **Main Tab**, and the printouts when a humidity and temperature probe is hooked up to the instrument.




The screenshot shows the 'Environment' screen with a blue header bar. On the right, it displays the date '10/03/08' and time '9:03:06 AM' next to a battery icon. The main title is 'Temperature Units'. Below it is a list box containing '°F' and '°C', with '°C' selected and highlighted in blue. To the right of the list box are icons for a thermometer and a water drop. At the bottom right are three buttons: 'Defaults', 'OK', and 'Cancel'.

Field	Description
°F	Display temperature in degrees Fahrenheit.
°C	Display temperature in degrees Celsius.


Locations Screen

Associating collected samples with labeled locations can help keep your data organized. The APC ERGOTOUCH PRO allows you to create up to 250 labeled locations (up to 10 characters in length). Use this screen to add, remove, or modify a location names to the list of locations.



The screenshot shows the 'Locations' screen with a blue header bar. On the right, it displays the date '11/21/08' and time '11:43:43 AM' next to a battery icon. The main area is a list box containing 'Room121', 'Room156', 'Room177', and 'Room203'. To the right of the list box is a compass rose icon. At the bottom are five buttons: 'Add...', 'Edit...', 'Remove', 'OK', and 'Cancel'.

To modify a location name, highlight the name in the list, then click the **Edit...** button. In the **Enter Location** screen click the edit box in the middle and use the on-screen keyboard to modify a location name. (You cannot edit an empty location). Click **OK** when finished.




The screenshot shows the 'Enter Location' screen with a blue header bar. On the right, it displays the date '10/03/08' and time '9:12:24 AM' next to a battery icon. The main text says 'Enter new name for location.' Below this is a text input box containing 'Room121'. To the right of the input box is a compass rose icon. At the bottom right are two buttons: 'OK' and 'Cancel'.

To add a location, click on the **Add...** button . In the **Add Location** screen click in the edit box in the middle and use the on-screen keyboard to add a location name. Click **OK** when finished.


To remove a location, click on location to be removed and click the **Remove** button.

Back in the main **Locations** screen, after all editing has been completed, press **OK** when finished.



Alarms Screen

Use this screen to set the alarm threshold for each channel. Press **OK** when finished.



Field	Description
Enable	Select the channels on which you want to enable alarms.
Threshold	To change the threshold for any channel, click the up and down arrows for that channel or use the on-screen keyboard to change its value. The threshold value units use the currently selected display units (see Count Units screen).

When a channel value exceeds the threshold value you set, the channel data is highlighted in red on the **Main Tab**, an audible alarm sounds, and the alarm icon appears on the **Main Tab**.

To clear the alarm, click the alarm icon . In addition, the record is printed if you have selected that option on the **Print Schedule** screen.

Recipe Screen

Use this screen to load and save recipes. **Recipes** let you save a group of settings (recipe) that you use over and over so you don't have to reset individual settings. There may be up to 100 recipes stored in the unit.

Recipe

10/03/08
9:23:21 AM

Save
Save As...

Fed Std

Load
Delete
Close

Field	Description
Save	<p>When you select Save, a new window opens that lets you enter a name for the recipe you want to save. You can also save the current settings in an existing recipe by selecting the recipe and using Save. The settings / parameters that are saved include:</p> <p>For each channel (1-6):</p> <ul style="list-style-type: none"> • Alarm setting (on/off) • Alarm threshold (value) • Channel setting (enabled/disabled) • Channel size (value) <p>Sample Timing settings</p> <ul style="list-style-type: none"> • Count mode • Count total • Start delay (in seconds) • Hold delay (in seconds) • Sample time (in seconds) <p>Count Mode/Units Settings</p> <ul style="list-style-type: none"> • Display normalized • Units (count, ft³ or m³) • Cumulative/Differential • Volume units <p>Printing settings</p> <ul style="list-style-type: none"> • Auto print and mode • Print cumulative/differential • Print model, separator, serial number
Save As	When you select Save As , a new window opens that lets you enter a name for the recipe you want to save.
Load	Highlight the recipe you want to load and press Load . The settings/parameters are reset to the values of that recipe.
Delete	Highlight the recipe you want to delete and press Delete . The recipe is deleted.

Data Tab


The **Data Tab** lets you preview data that has been collected. Use the elevator (slide) on the right to scroll through the records. The record number is displayed at the bottom of the tab. As each record displays, its data and relevant parameters are displayed.

				11/18/08 02:07:02 PM
m³	Size μm	Δ m³	Σ m³	
	0.3	275695	286863	
	0.5	9733	11168	
	1.0	1214	1435	
	3.0	53	221	
	5.0	145	168	
Print	10.0	23	23	
Location: Sample 003 Laser: OK Alarm: YES Sample: 00:01:00 Vol: 131.0 L Flow: ALRM Date: 11/18/08 Temperature: 24.4 °C Time: 11:59:56 AM Humidity: 53 %				
Record:		20	Records:	123 / 10000
Main	Setup	Data	Reports	

Field	Description
#, ft ³ , m ³	Button used to change between counts and concentration displays.
Size μ m	Channel size.
Δ	Differential concentration.
Σ	Cumulative concentration.
Location	Location where the data was collected.
Sample	Duration of the sampling period.
Date	Date on which the data was collected.
Time	Time at which data was collected.
Temperature	Temperature at the end of the time the data was collected (if probe connected during sampling).
Humidity	Relative Humidity level at the end of the time the data was collected (if probe connected during sampling).
Flow	Status of the flow.
Alarm	Alarm threshold was triggered (YES) or not (NONE).
Laser	Status of the laser.

The print button will allow a range of sample data to be printed using the optional external printer.

Print

 10/20/08 2:55:30 PM

Select range to print.

☒ All

☐ Records

1

▲▼


–

107

▲▼


Print

Close



The print data screen will show progress on the current selected range of sample data to be printed. Press the **Cancel Printing** button to cancel the rest of the print job.


Print

 10/22/08 3:19:26 PM

Printing Sample
16


15% Completed

Cancel Printing



Reports Tab

Use this screen to select various standard reports for viewing and printing



The standard reports are shown below:

Fed Std 209E FT 11/20/08 11:55:25 AM

Room Area:	96.88 ft ²	
Class Level:	100	
Air Flow:	Unidirectional	
Room Status:	Operational	
Min Locations:	4	
Min Samples:	5	

Allowable Sizes	Min Volume	
0.2 um	1.00E-001	f3
0.3 um	1.00E-001	f3
0.5 um	2.00E-001	f3

Generate

Close

Fed Std 209E M 11/20/08 11:56:35 AM

Room Area:	9.00 m ²	
Class Level:	M2	
Air Flow:	Unidirectional	
Room Status:	Operational	
Min Locations:	4	
Min Samples:	5	

Allowable Sizes	Min Volume	
0.1 um	5.71E+000	L
0.2 um	2.64E+001	L
0.3 um	6.47E+001	L
0.5 um	2.00E+002	L

Generate

Close

EC-GMP 11/20/08 09:01:59 AM

Room Area:	9.00 m ²	
Class Level:	C	
Air Flow:	Unidirectional	
Room Status:	At-Rest	
Min Locations:	3	
Min Samples:	3	

Allowable Sizes	Min Volume	
0.5 um	2.00E+000	L
5.0 um	6.90E+000	L

Generate

Close

ISO 14644-1 11/20/08 09:00:51 AM

Room Area:	9.00 m ²	
Class Level:	3	
Air Flow:	Unidirectional	
Room Status:	At-Rest	
Min Locations:	3	
Min Samples:	3	

Allowable Sizes	Min Volume	
0.1 um	2.00E+001	L
0.2 um	8.44E+001	L
0.3 um	1.96E+002	L
0.5 um	5.71E+002	L
1.0 um	2.50E+003	L

Generate

Close

Field	Description
Room Area	Displays the area of the room in ft ² or m ² .
Class Level	Depends on the Room Definition , see below.
Air Flow	Displays the airflow characteristics of the room.
Room Status	Displays the status of the room. See Room Definition screen below.
Min Locations	Displays the minimum number of locations that must be sampled in the room.
Min Samples	Displays the minimum number of samples that must be taken at each location.
Allowable Sizes	Allowable channel sizes for the selected Class Level for that Standard.
Minimum Volume	Displays the minimum volume (in cubic feet or meters) that must be sampled on each channel.
Room Definition	Press to set definitions for the room. See Room Definition screen below.
Generate	Select to print a single record or a range of records. See Print screen below.

Room Definition Screen

Use this screen to define specific values for the room. Press **OK** when finished:

Room Definition

10/03/08
9:45:27 AM

Room Status
Operational


Air Flow
Unidirectional

Class
M2

Area
0.83612

☐ ft²
☒ m²

OK



Field	Description
Room Status	Select the room status: As Built, At Rest, or Operational.
Air Flow	Select the air flow: Unidirectional or Multidirectional.
Class	Select the class of the room: The class is dependent on the standard: FED FT3: 1, 20, 100, 1000, 10000, 100000 FED M3: M1.0, M1.5, M2.0, M3.0, M3.5, M4.0, M4.5, M5.0, M5.5, M6.0, M6.5, M7.0 ISO14644-1: 1, 2, 3, 4, 5, 6, 7, 8, 9 EC GMP: A, B, C, D
Area	Use the on-screen keyboard to enter the area of the room in ft ² or m ³ .

Generate Screen

This screen lets you generate the report using either a single record or a range of records. Press the **Generate** button to generate the selected report.

Generate... 10/03/08 9:47:47 AM

Select range used for reports.

☒ All
☐ Records

1 107

Generate...
Close

The generated report will be displayed on the screen and may be viewed on the screen or printed (optional printer must be attached) by pressing the **Print** button.

Fed Std 209E M 10/03/08 9:52:23 AM

Federal Standard 209E_M

** Model 9306 **

Serial Number: xxxxxxxxx

Targeted Class: M7.0

Room Area: 100m^2

Room Status: Operational

Air Flow: Unidirectional

Min Locations: 3

Min Samples/Room: 5

Print Close

Chapter 4

Data Handling

USB Computer Communication

The APC ERGOTOUCH PRO Airborne Particle Counter is equipped with a USB compatible cable for uploading and downloading information to a PC. The cable plugs into the right side of the instrument.



System Requirements

Hardware (minimum requirements):

- Pentium 90 MHz, 32 MB RAM, 20 MB of hard drive space, or minimum requirements for the operating system, whichever is higher. Free hard drive space depends on the amount of records to be stored.
- A free USB port.

Software (minimum requirements):

- Windows® XP operating system with Service Pack 2 or higher, Windows® Vista® operating system with Service Pack 1 or higher.
- .NET 2.0 or higher.

Installation

The ERGOTRAK Lite version 2 software is supplied on a CD that loads software and communications drivers for ErgoTouch or ErgoTouch Pro Handheld Airborne Particle Counters. Hardware manuals are also supplied on the CD along with this software manual.

Installation consists of two parts:

- **Installation of ERGOTRAK Lite software.**
Run "setup.exe" from the provided CD and follow on-screen instructions.
- **Installation of device drivers.**
Device driver installation is executed in the background during the setup process and does not require user input. Once the installation is finished, drivers are ready for use.

The computer must be restarted to complete the installation. Once the ErgoTouch or ErgoTouch Pro Airborne Particle Counter is connected and powered on, the system will automatically detect the device and will proceed with installing the necessary drivers.

- **ErgoTouch**

The system will detect the device and install all necessary drivers without requiring user input.

- **ErgoTouch Pro**

When a particle counter is connected for the first time, the system will automatically detect the device and will start driver installation process.



When asked "Can Windows connect to Windows Update search for software?", select "No, not this time" and click Next.



Select Install the software automatically (Recommended) and click Next. Hardware Wizard will search for the drivers and locate them in \System32\drivers directory.

Once that is done, the following screen will appear.



Depending on your system setup you may see a warning message.



Click <Continue Anyway> and the installation will proceed.




Once everything is completed, click Finish.

This procedure is required only if connecting an ErgoTouch or ErgoTouch Pro device for the first time. For all subsequent devices Windows® operating system will automatically locate necessary drivers and install them without requiring user input.


ERGOTRAK Lite Software Configuration

It is recommended that ERGOTRAK Lite software be configured prior to connecting to the ErgoTouch or ErgoTouch Pro Airborne Particle Counter. The software can be configured at any time, but configuring it before connecting any devices will ensure smoother operation. ErgoTrak Lite software can be started either by clicking on its icon on the desktop or by selecting it through the Start menu. Shortcut to the application is located at **Start->All Programs->Biotech->ERGOTRAK Lite v2.0**.

If ERGOTRAK Lite software detects that the data delimiter for comma-separated values file is the same as the decimal point symbol on your machine, the following warning will be presented:



Clicking <Yes> will bring you to the [preferences](#) page. If you click <No> or <Cancel>, you can select preferences page by clicking on the Preferences tab inside main application window:




The main parameters related to the ErgoTrak Lite software application are configured on the Preferences Page.

Note: All settings on Preferences page are for ERGOTRAK Lite data storage and presentation only. The ErgoTouch or ErgoTouch Pro Airborne Particle Counter device configuration is not changed.

Data Presentation

This section allows you to select specific parameters for displaying reports on data page.

Option	Description
Format	Is used to specify whether columns containing number of particles for any given channel are displayed as differential or cumulative values.
Units	This subsection specifies if data is displayed as the total number of particle counts during sampling period or as the particle concentration (Particle Counts/Sampled Volume of Air). When concentration is selected, the units of volume are determined based on the selection made in Volume subsection.
Volume	Allows you to select which units to use for concentration calculations.
Temperature	Specifies temperature readings are displayed in units of °F or °C. This is only applicable to the Airborne Particle Counter which has external temperature probe capability with the probe installed.



These settings can be changed at any point when the software is operating, and the displayed reports automatically refresh calculated values to reflect new changes.


Data Export

In addition to importing data from ErgoTouch or ErgoTouch Pro Airborne Particle Counters, ERGOTRAK Lite software allows you to select additional data which can be contained in the saved export file:




Option	Description
Prefix serial number to records	If selected, each data record will be prefixed with the serial number of the device from which it was retrieved.
Include statistics	Includes basic statistics such as minimum, maximum, average, standard deviation for all records being exported.
Include device info	Includes device information and provides information that contains data about device and units used in calculations.
Include data headings	Provides headings for data records.

Option	Description
Delimiter	Specifies which character should be used to delimit data in the comma separated value file. It is crucial that this character differs from the character used to display decimal point in numbers in order to avoid data misalignment.



File Management

Save location is a default place where data files should be saved. Every time you choose to save information to a file, the Save File dialog is still presented allowing you to choose a different location if necessary.




All preference settings are automatically saved once ERGOTRAK Lite software is closed and will be recalled next time ERGOTRAK Lite software is started. The settings are saved on the Windows® user level. This allows a single installation of ERGOTRAK Lite software for all users on a single computer while still allowing a user to save individual preferences.

Device Connection

ERGOTRAK Lite software is capable of connecting multiple ErgoTouch or ErgoTouch Pro Airborne Particle Counters. It will recognize device type and will download data accordingly.

The ERGOTRAK Lite application automatically detects connected devices and initiates records download as well as download of location names and recipes (if supported by the specific ErgoTouch or ErgoTouch Pro Airborne Particle Counters). Although ERGOTRAK Lite software searches for new devices every 5 seconds, it is recommended that devices are fully powered on before connecting to the computer. This will ensure faster response time and fewer connection errors.

Once devices are detected, they will be displayed on left application pane.




Working with Connected Devices

Clicking on the + sign next to the ErgoTouch or ErgoTouch Pro Airborne Particle Counter icon expands the information supported by that model of ErgoTouch or ErgoTouch Pro particle counter potentially including records, locations, and recipes.



Expanding “Records”, “Locations”, or “Recipes” in same fashion allows access to corresponding data.



Displaying Records

Clicking on the check box next to **Records** icon will select all the available records for displaying in the report tab page. Each record is marked with date and time it was recorded.

The screenshot shows the software interface with a file tree on the left and a report window on the right. The file tree shows a folder named '9303-0193030903003' containing a 'Records' folder with 14 records, each with a date and time. The report window displays the following data:

Record Number	Location Number	Date & Time	Sample Time	Hold Time	Ch1 Size (um)	Cumulative	Ch2 Size (um)	Cumulative	Ch3 Size (um)	Cumulative
1	250	5/14/2009 11:41:07	00:00:10	00:00:01	0.3	11215	0.5	1221	5.0	12
2	250	5/14/2009 11:41:18	00:00:10	00:00:01	0.3	10859	0.5	1110	5.0	23
3	250	5/14/2009 11:41:29	00:00:10	00:00:01	0.3	10263	0.5	1132	5.0	29
4	250	5/14/2009 11:41:40	00:00:10	00:00:01	0.3	10028	0.5	1049	5.0	16
5	250	5/14/2009 11:41:51	00:00:10	00:00:01	0.3	10373	0.5	1110	5.0	16
6	250	5/14/2009 11:42:02	00:00:10	00:00:01	0.3	9733	0.5	998	5.0	21
7	250	5/14/2009 11:42:13	00:00:10	00:00:01	0.3	10142	0.5	1084	5.0	21
8	250	5/14/2009 11:42:24	00:00:10	00:00:01	0.3	10577	0.5	1115	5.0	14
9	250	5/14/2009 11:42:35	00:00:10	00:00:01	0.3	10085	0.5	1024	5.0	21
10	250	5/14/2009 11:42:46	00:00:10	00:00:01	0.3	10392	0.5	1043	5.0	33
11	250	5/14/2009 11:42:57	00:00:10	00:00:01	0.3	10053	0.5	1053	5.0	33
12	250	5/14/2009 11:43:08	00:00:10	00:00:01	0.3	10049	0.5	971	5.0	33
13	250	5/14/2009 11:43:19	00:00:10	00:00:01	0.3	9881	0.5	975	5.0	25
14	250	5/14/2009 11:43:30	00:00:10	00:00:01	0.3	9901	0.5	1011	5.0	40

Data in the report is displayed corresponding to the settings that were chosen on Preferences page Chapter 2, "[Data Presentation](#)". Data is only displayed from the records which are checked. The report is automatically refreshed every time the application preferences are changed or new records are (un)checked.

In order to view records from different ErgoTouch or ErgoTouch Pro Airborne Particle Counters, select the records of interest and the corresponding report will be generated and displayed.

The screenshot shows the software interface with a file tree on the left and a report window on the right. The file tree shows a folder named '9306-V 70806175' containing a 'Records' folder with 14 records, each with a date and time. The report window displays the following data:

Record Number	Location Number	Date & Time	Sample Time	Hold Time	Ch1 Size (um)	Cumulative	Ch2 Size (um)	Cumulative	Ch3 Size (um)	Cumulative
1	250	5/14/2009 11:41:07	00:00:10	00:00:01	0.3	11215	0.5	1221	5.0	12
2	250	5/14/2009 11:41:18	00:00:10	00:00:01	0.3	10859	0.5	1110	5.0	23
3	250	5/14/2009 11:41:29	00:00:10	00:00:01	0.3	10263	0.5	1132	5.0	29
4	250	5/14/2009 11:41:40	00:00:10	00:00:01	0.3	10028	0.5	1049	5.0	16
5	250	5/14/2009 11:41:51	00:00:10	00:00:01	0.3	10373	0.5	1110	5.0	16
6	250	5/14/2009 11:42:02	00:00:10	00:00:01	0.3	9733	0.5	998	5.0	21
7	250	5/14/2009 11:42:13	00:00:10	00:00:01	0.3	10142	0.5	1084	5.0	21
8	250	5/14/2009 11:42:24	00:00:10	00:00:01	0.3	10577	0.5	1115	5.0	14
9	250	5/14/2009 11:42:35	00:00:10	00:00:01	0.3	10085	0.5	1024	5.0	21
10	250	5/14/2009 11:42:46	00:00:10	00:00:01	0.3	10392	0.5	1043	5.0	33
11	250	5/14/2009 11:42:57	00:00:10	00:00:01	0.3	10053	0.5	1053	5.0	33
12	250	5/14/2009 11:43:08	00:00:10	00:00:01	0.3	10049	0.5	971	5.0	33
13	250	5/14/2009 11:43:19	00:00:10	00:00:01	0.3	9881	0.5	975	5.0	25
14	250	5/14/2009 11:43:30	00:00:10	00:00:01	0.3	9901	0.5	1011	5.0	40

Data from different ErgoTouch or ErgoTouch Pro Airborne Particle Counters can be viewed by selecting (clicking) the record in the file display belonging to a desired ErgoTouch or ErgoTouch Pro instrument and the corresponding report will be displayed.

The report viewer tab page allows some basic operations such as printing, moving between pages, etc. Those options are accessible via the toolbar above the report.



Device Related Application Menus

All allowable operations for a given ErgoTouch or ErgoTouch Pro instrument are accessible via the right-click menu or by the “Device” menu at the top of application. Right-clicking on the device of interest or selecting the device by left-clicking and then selecting **Device** menu, potentially provides the following options (certain ErgoTouch or ErgoTouch Pro Airborne Particle Counters do not support all of these options.)



Main Application Menu

This is the menu that provides access to all ERGOTRAK Lite software functionality for the specified device. If device does not support certain options, those options will be grayed out and inaccessible.

Delete All Records


Immediately deletes all the records from the device. This operation is irreversible and you are prompted with the following warning to avoid selecting this item inadvertently.



Once you confirm decision by clicking **Yes** button, records are deleted on the device as well as from ERGOTRAK Lite software application. It is always recommended to save data before deleting it from the device.

Export all Records

This allows you to export all the records from the device in the following formats: XML, HTML, Excel, CSV, and Raw Data. [Data Presentation](#) settings affect how the data from the records is saved and exported.




File Format	Description
XML (extensible Markup Language)	Saving data in this format allows for transfer between ERGOTRAK LITE software and other third party applications that need access to data without any formatting.
HTML	Creates a web page file that contains data in tabular format that is convenient for displaying with Web browsers.
Excel	This file saves data in Microsoft® Office XML format (not to be confused with regular XML) which can be opened by Microsoft® Excel® XP or newer. This file format allows for data formatting such as different fonts and cell background colors.
CSV	A comma separated value file that stores all data in plain text format delimited by the character selected under Data Export settings on Preferences tab. page.
Raw Data	This option allows you to save records for future ERGOTRAK Lite software use. Data is saved in binary format and can be opened through File->Open Data File menu. Only data <i>records</i> are saved if this option is selected. Recipes and locations have to be saved separately utilizing the corresponding menu selections. This is done in order to allow you to save only recipes or location in order to import them to other devices.

Export Selected Records

Same functionality as “Export All Records” but will only export selected records instead of all the records from the device.

Recipes


If the device selected supports recipes, this menu item allows you to perform certain manipulations on recipes that are loaded on the device.



Options	Description
Export All	Allows you to export all the recipes to a binary file that can be opened by ERGOTrak Lite software. This option is used for making backups of recipes on the device and enabling you to move recipes to different devices without the need to create them one by one.
Load from File	Loads user-specified recipe file that was created by using "Export All" option. Recipes will be added to the device in ERGOTrak Lite software but will not be transferred to the device until you select Upload to Device option.
Upload to Device	Upload all the changes made to recipes to the device immediately.
Delete from Device	Deletes all the recipes from the device immediately.
New Recipe	Creates new default recipe.
Delete Recipe	Deletes selected recipe. Note: Only "Delete from Device" and "Upload to Device" operations are executed immediately, for all other operations, you have to select "Upload to Device" in order for changes to take effect.

If the connected device supports recipes, ERGOTrak Lite software allows you to modify, create, delete and export recipes.

Selecting the recipe of interest in the left pane opens Recipe tab page in the right pane:




This page allows you to modify existing recipe's parameters.

Options	Description
Recipe Title	Used to change the name of recipe.
Number of Runs	The number of sample runs to be executed.
Start Delay	Amount of time to wait before starting taking samples.
Sample Delay	Time delay between sampling runs.
Volume/Time per Run	This setting allows you to set up the length of sampling time for each run. You can either set up runtime by entering the volume of air to be sampled or by choosing sampling time values.
Count Mode	Is used to set up device report display on the unit. It does not affect actual data stored in the record.
Bin (1...x)	Sets up bin boundary for each channel the instrument supports. If instrument does not support variable binning, this area is grayed out and unavailable. This menu option also enables setting of alarm thresholds, enable/disable alarms, and enable or disable particular bin from displaying in the report on the device. These settings do not affect values in downloaded records. The data from all channels is stored in the data record, even if the channel is not displayed on the ErgoTouch or ErgoTouch Pro screen.

Locations


If the device selected supports location, naming this item allows you to manipulate downloaded locations.



You can save existing labels, create new ones, or rename existing ones. All the changes except selecting “Delete from Device” menu option are not applied to the device until you select **Upload to Device** from the right-click menu under Locations or from the Device menu at the top of application’s window.

Options	Description
Export All	Allows you to export all location labels to a binary file that can be opened by ErgoTrak Lite software. This option is used for making backups of location labels on the device and enabling you to move location labels to different devices without the need to create them one by one.
Load from File	Loads user-specified location labels file that was created by using “Export All” option. Locations will be added to the device in ErgoTrak Lite software but will not be transferred to the device until you select Upload to Device option.
Upload to Device	Upload all the changes made to location labels to the device immediately.
Delete from Device	Deletes all the location labels from the device immediately.
New Location	Creates new location label.
Delete Location	Deletes selected location label.

If the device connected supports location labels, they will be displayed under the device icon on the left side pane.




Note:

In order to rename existing location, left-click on the location label of interest and then left-click on it again. This will change the label to edit mode. You can now type in a new name:




File Menu

This menu enables you to load existing ERGOTRAK Lite software records file and exit the application.



Selecting **Open Data File** presents you with a dialog window to choose which file to open.



Once the selected file is loaded, ERGOTRAK Lite software will display device information of the device form which records have been downloaded initially in the left pane. This option is only for records' data files. Recipe and Location data files can only be opened through the right-click menu or Device menu when the device that supports them is selected in left pane.

Exit allows you to exit the application. If there are any changes that have not been saved, you will be prompted to update devices with the changes.

Device Menu

This menu is the same as right-click menu.

Help Menu

Allows you to display this help file and information about the software.

Troubleshooting

If connection is broken or device stops responding during data transfer, the ERGOTRAK Lite software will display only data that has been received up to that point. The ErgoTouch will not generate an application error. However, the ErgoTouch Pro will display the application error. This is due to the nature of communication between devices and ERGOTRAK Lite software.

Solution: Power down the device and then turn it back on. Close the application and start it again.

If you have any difficulty setting up or have any technical or application questions about this system, contact us at the address provided in Chapter 7.

Chapter 5

Maintenance

Note

There are no user-serviceable parts inside this instrument. Opening the instrument case may void the warranty. Biotest recommends that you return the APC ERGOTOUCH PRO Airborne Particle Counter to the factory for any required maintenance or service not described in this manual.

Maintenance Schedule

Biotest recommends annual factory cleaning and calibration for the APC ERGOTOUCH PRO Airborne Particle Counter. See [Chapter 7, "Contacting Customer Service"](#) for service/calibration.

Recommended Field Maintenance Schedule

Item	Frequency
Daily zero check	Daily (or before each use).
Factory cleaning and calibration	Annually.
Cleaning the instrument enclosure	As needed

Daily Zero Check

The daily zero check ensures that the instrument is properly assembled and free from leaks, residual particles and electronic noise. Please see Chapter 2, "[Getting Started](#)" for detailed instructions on performing the zero check.




Cleaning the Instrument Enclosure

To clean the enclosure, dampen a lint-free cloth and gently wipe the surface until surface contamination is removed.

Chapter 6

Troubleshooting

Symptom	Possible Cause	Corrective Action
Counts are too low	<ul style="list-style-type: none"> Instrument is being operated outside temperature or relative humidity specifications. Internal parts have been damaged because instrument was stored at a temperature greater than 122 °F (50 °C). Instrument has contamination on the optics due to condensation or excessive loading Laser or pump control is damaged Unit is due for calibration 	<ul style="list-style-type: none"> Operate instrument within specifications Return to factory for service Return to factory for service Return to factory for service Return to factory for service
Instrument does not turn on	<ul style="list-style-type: none"> Battery is dead AC cord is not plugged into unit 	<ul style="list-style-type: none"> Recharge battery or connect to AC power Connect AC cord
Instrument does not meet zero count specification (≤ 1 particle/5 mins)	<ul style="list-style-type: none"> HEPA filter is not connected properly and room air is leaking into the HEPA filter assembly Residual particles from previous samples are shedding off internal parts and into the optics An internal component has been damaged due to operation outside of temperature specifications or one or more excessive bumps or jolts, and electronic noise is inducing false counts A leak has developed in the aerosol flow path Internal optics have become dirty 	<ul style="list-style-type: none"> Check that the HEPA filter has been tightly connected to the inlet. Check that rubber o-ring (black) on the inlet is in place Purge instrument by running the instrument for 10–15 minutes before attempting zero count test Return to factory for service Return to factory for service Return to factory for service
Battery does not charge	<ul style="list-style-type: none"> The unit must be turned on but not in sampling mode for the battery to charge 	<ul style="list-style-type: none"> Turn on unit.

Symptom	Possible Cause	Corrective Action
LOW BATTERY ERROR 	<ul style="list-style-type: none"> Low battery 	<ul style="list-style-type: none"> Recharge battery or connect AC cord
PHOTODETECTOR ERROR	<ul style="list-style-type: none"> Direct light is entering the aerosol inlet Laser has become misaligned due to excessive bumps or jolts Internal optics have become dirty 	<ul style="list-style-type: none"> Remove instrument from direct light Return to factory for service Return to factory for service
SYSTEM ERROR	<ul style="list-style-type: none"> Information is not being read properly by microprocessor 	<ul style="list-style-type: none"> Restart instrument. If problem persists, contact Biotest technical support
TEMPERATURE HUMIDITY PROBE ERROR	<ul style="list-style-type: none"> Temperature/RH probe was not recognized 	<ul style="list-style-type: none"> Detach and reconnect probe. If problem persists, contact Biotest technical support
FLOW ERROR 	<ul style="list-style-type: none"> Instrument was unable to control flow rate (if any tubing is connected to particle counter) Pressure drop across inlet may be too large Inlet not at ambient pressure 	<ul style="list-style-type: none"> Restart measurement Lessen pressure drop across inlet by using larger diameter tubing, less tubing, and/or adding a bleed valve Do not subject the unit to other than ambient pressure conditions
LASER POWER WARNING 	<ul style="list-style-type: none"> Laser power has fallen outside of specification 	<ul style="list-style-type: none"> Return to factory for service

Chapter 7

Contacting Customer Service

This chapter gives directions for contacting people at Biotest AG for technical information and directions for returning the APC ERGO TOUCH PRO Airborne Particle Counter for service.

Technical Contacts

- If you have any difficulty setting up or operating the APC ERGO TOUCH PRO, or if you have technical or application questions about this system, contact an applications engineer at Biotest AG, +49 (0) 6221 726 5130 or e-mail hycon@biotest.de.
- If the APC ERGO TOUCH PRO, does not operate properly, or if you are returning the instrument for service, visit our website at <http://www.biotest.com>, or contact Biotest Customer Service at +49 (0) 6103 801 496.
- Call your local Biotest Subsidiary or Representative.

Returning the APC ERGO TOUCH PRO Airborne Particle Counter for Service

First contact your local Biotest representative to make arrangements. If you don't know or don't have a Biotest representative, you can contact Biotest directly:

Visit our website at <http://www.biotest.com> or call Biotest Technical Support at +49 (0) 6221 726 5130 or Biotest Customer Service at +49 (0) 6103 801 496 for specific return instructions. Customer Service will need this information when you call:

- The instrument model number
- The instrument serial number
- A purchase order number (unless under warranty)
- A billing address
- A shipping address

All instruments must be properly decontaminated prior to packing and return to Biotest. Use the original packing material to return the instrument to Biotest. If you no longer have the original packing material, seal off any ports to prevent debris from entering the instrument and ensure that the display and the connectors on the instrument front and back panels are protected.

APPENDIX A

Specifications

All specifications meet or exceed JIS B 9921: 1997. They are subject to change without notice.

Size Range	0.3–10 µm
Channel Sizes	0.3, 0.5, 1.0, 3.0, 5.0, 10.0 µm (user configurable in 0.1 µm units)
Counting Efficiency	50% @ 0.3 µm; 100% for particles > 0.45 µm (per JIS)
Concentration Limits	2,000,000 particles/ft ³ @ 5% coincidence loss
Light Source	Laser diode
Zero Count Level	≤1 count/5 minutes (<2 particles/ft ³) Meets JIS B 9921:1997
Flow Rate	0.1 CFM (2.83 L/min) with ±5% accuracy
Flow Control	Automatic
Calibration	NIST traceable
Sample Probe/Tubing	Isokinetic sampling probe
Sampling Modes	Manual, automatic, beep, concentration and cumulative/differential
Sampling Time	1 second to 99 min 59 sec
Sampling Frequency	1 to 999 cycles or continuous
Sample Output	Internal HEPA filter
Vacuum Source	Internal pump
Communication Mode	RS-232, USB
Data Storage	10,000 samples
Data Security	Password protected
External Alarm	Internal audible alarm. Indicators for counts, low battery, and sensor failure
Environmental Sensors	Optional temperature/RH probe supported
Display	QVGA 3.7-inch touch screen with Windows CE
Languages	English
Reports	ISO-14644-1, FS-209E & EC GMP
Printer	Optional external printer supported
External Surface	High impact injection molded plastic
Power	110 to 240 VAC 50 to 60 Hz Universal in-line power supply
Battery	Removable Li-Ion
Battery Life	Up to 3 hours of continuous use
Recharge Time	3 hours
Dimensions (L x W x H)	25.4 x 11.4 x 7.6 cm (10 x 4.5 x 3 in.)
Weight	1.0 kg (2.2 lbs) with battery
Warranty	2 years

Operating Conditions	5 °C to 35 °C; 20% to 95% non-condensing relative humidity
Storage Conditions	0 °C to 50 °C; Up to 98% non-condensing relative humidity
Included Accessories	Power supply, power cord, battery, isokinetic inlet, stylus, purge filter, ERGOTRAK Lite data download software, operational manual on CD, computer cable, calibration certificate, and Quick Start Guide.
Optional Accessories	Temp R/H probe, stainless steel isokinetic inlet and probe, tubing, barbed inlet fitting, printer, printer paper, carrying case and external battery charger

Temperature / RH Probe (optional accessory)

Temperature Range	32 to 115°F (0 to 45°C)
Accuracy	±2°F (±1°C)
Relative Humidity Range	10 to 90% RH
Accuracy	±5% RH