

Technical Documentation

EvitaXL Intensive Care Ventilator



Revision 5.0 5664.590 9036054

Because you care

General

1	Symbols and Definitions		
2	Note	S	3
Fun	ction	Description	
1	Gene	eral	7
2	Basi	c principle	7
	2.1	Control panel	8
	2.2	Electronics	8
	2.3	Pneumatics	8
	2.4	Simplified block diagram	9
3	Elect	ronics	10
	3.1	CPU PCB 68332	10
	3.2	CO2 Carrier PCB	12
	3.3	CO2 measurement	.15
	3.4	CO2 sensor	16
	3.5	Processor Board PCB	. 18
	3.6	Power Supply PCB	19
3.7 Power pack		Power pack	20
	3.8	Graphic Controller 8 PCB	21
	3.9	Communication PCB (optional)	24
	3.10	Pediatric Flow PCB (optional)	26
	3.11	IFCO Carrier PCB for additional optional features	.27
	3.12	SmartCare PCB (option)	. 28
	3.13	Pneumatic Controller PCB	29
	3.14	HPSV Controller PCB	33

Contents

4	Pneu	Imatics	34
	4.1	Gas connection block	. 35
	4.2	Parallel mixer	. 38
	4.3	Inspiratory block	.42
	4.4	Patient system	.44
	4.5	PEEP/PIP valve	.46
5	Theo	ory of operation	47
	5.1	AIR supply	.47
	5.2	O2 supply	.48
	5.3	Inspiratory phase	.49
	5.4	Paw high alarm limit	. 50
	5.5	Emergency pressure relief mechanism (safety valve)	. 50
	5.6	Expiration	. 50
6	Gas	mixture	52
	6.1	Correcting the oxygen cartridge valve	. 52
	6.2	Dependence on the supply pressure	. 52
7	Fund	tion diagram	53
Mair	ntena	ance Procedures	
1	Cool	ing-air filter	57
	1.1	Replacing the cooling-air filter	.57
2	Amb	ient-air filter	59
	2.1	Replacing the ambient-air filter	. 59
3	Lithi	um battery (CPU 68332 PCB)	61
	3.1	Replacing the lithium battery (CPU 68332 PCB)	.61

Contents

4	Real-time clock		
	4.1	Replacing the real-time clock	. 63
Ser	vice	Diagnosis Mode	
1	Gen	eral	69
2	Call	up of service diagnosis mode	71
3	Diag	nosis 'Front'	72
	3.1	Diagnosis page 'Processor'	. 72
4	Diag	nosis 'Electronic'	73
	4.1	Diagnosis page 'Processor'	. 73
	4.2	Diagnosis page 'Sensors'	. 76
	4.3	Diagnosis page 'add. HW'	. 78
5	Diag	nosis 'Pneumatic'	79
	5.1	Diagnosis page 'Processor'	. 79
	5.2	Diagnosis page 'Sensors'	. 80
	5.3	Diagnosis page 'Valves'	. 83
6	Diag	nosis 'Other'	85
	6.1	Diagnosis page 'Lokbook' (error list)	. 85
	6.2	Diagnosis page 'Option Release'	. 88
	6.3	Diagnosis page 'SW Info'	. 88

Contents

Schematics and Diagrams

1 Piping and function diagram

91

Annex

Test List

Technical Information

General

1	Symbols and Defini-		
	tions	WARNING	
		A WARNING statement per tially hazardous situation or serious injury.	rovides important information about a poten- a which, if not avoided, could result in death
		CAUTION A CAUTION statement pro hazardous situation which, injury to the user or patient erty.	vides important information about a potentially if not avoided, may result in minor or moderate or in damage to the equipment or other prop-
		NOTE A NOTE provides additiona during operation.	al information intended to avoid inconvenience
		Definitions according to Ge	rman standard DIN 31051:
		Inspection	= examination of actual condition
		Maintenance	 measures to maintain specified condition
		Repair	 measures to restore specified condition
		Servicing	= inspection, maintenance, and repair
2	Notes	This Technical Documentat	ion conforms to the IEC 60601-1 standard.
		Read each step in every pro Always use the proper tools the instructions and/or reco the equipment may operate damaged.	ocedure thoroughly before beginning any test. and specified test equipment. If you deviate from mmendations in this Technical Documentation, improperly or unsafely, or the equipment could be
		It is our recommendation to	use only Dräger parts and supplies.
		The maintenance procedure be performed by qualified s dures do not replace inspec	es described in this Technical Documentation may ervice personnel only. These maintenance proce- ctions and servicing by the manufacturer.
		The information in this Tech be disclosed to third parties turer.	nnical Documentation is confidential and may not without the prior written consent of the manufac-
		This Technical Documentat descriptions found in this Te for reading and studying the with the product at the time	ion is for the purpose of information only. Product echnical Documentation are in no way a substitute e Instructions for Use/Operating Manual enclosed of delivery.
		Know-how contained in this change through research as right to make changes to th	Technical Documentation is subject to ongoing nd development and Dräger Medical reserves the is Technical Documentation without notice.

NOTE

Unless otherwise stated, reference is made to laws, regulations or standards (as amended) applicable in the Federal Republic of Germany for equipment used or serviced in Germany. Users or technicians in all other countries must verify compliance with local laws or applicable international standards.

Function Description

2

- 1 General The machine is a time-controlled, constant-volume long-term ventilator for adults and children.
 - **Basic principle** The machine communicates via a serial interface (CAN) and consists of the following assemblies:
 - Control panel
 - Electronic assembly
 - Pneumatic assembly



Control unit (front)



2.1	Control panel	The control unit is the interface between the machine and the operator. The control unit is used to set parameters, it displays measured values, and generates warnings. The control unit comprises the following subassemblies:
		 15" TFT display
		 Membrane keypad
		 Touchscreen with resistive touch
		 Rotary encoder
		 Graphic Controller 8 PCB
		– Connector PCB
2.2	Electronics	The electronic assembly is the central control unit. The electronic assembly includes the CPU 68332 PCB, the CO2 Carrier PCB with the Processor PCB and the Power Supply PCB, and the power pack (Communication PCB, Pediatric Flow PCB, IFCO Carrier PCB, and the SpO2 PCB as option).
2.3	Pneumatics	The pneumatic assembly controls the pneumatic valves following preset ven- tilation parameters. The pneumatic assembly includes an independent micro- processor system and the valve control. The pneumatic assembly comprises the Pneumatic Controller PcB, the HPSV Controller AIR/O2 PCB, the PEEP valve, the mixer, the compressed-gas connection, the flow sensor, and the O2 sensor.

2.4 Simplified block diagram



Figure 2 Basic principle

1	Keys	13	Supply voltages
2	Rotary encoder	14	Power switch
3	Touchscreen	15	Second inspiratory Paw
4	TFT display	16	Reset pneumatics processor and venting
5	Information LEDs and Alarm LEDs	17	Electronic processor reset and second loudspeaker alarm
6	CAN bus	18	Inspiratory Paw
7	Graphics processor reset	19	O2 sensor
8	Not applicable	20	FiO2 (HPSV mixer)
9	Loudspeaker with sound chip	21	AIR (HPSV mixer)
10	Second loudspeaker (piezo)	22	flow sensor
11	Voltage monitoring (activates reset of the processors and the piezo)	23	Expiratory valve with PEEP
12	Rechargeable battery (Goldcap capacitor)	24	Expiratory Paw

- 3 Electronics
- **3.1 CPU PCB 68332** The CPU 68332 PCB is integrated in the electronic assembly of the machine. The CPU 68332 PCB includes an independent microprocessor system, two external interfaces, three internal interfaces, the loudspeaker control and a serial EEPROM.



Figure 3 Block diagram of the 68332 PCB

- **3.1.1 EEPROM** The EEPROM is connected to the synchronized, serial interface of the 68332. The EEPROM characterizes the machine (enabled options, serial number, etc.).
- **3.1.2 Microprocessor system** The microprocessor systems consists of a 68332 CPU, a 512 kB randomaccess memory (RAM), and a 1 MB electrically programmable and electrically erasable read-only memory (flash EPROM). The RAM is battery-buffered. When the battery is being replaced a Goldcap capacitor ensures voltage supply of the RAM. Programming of the flash EPROMs is only possible if the system identifies the "SERVICE-Q" signal.

Evita	XL	Function Description
3.1.3	RS232 interface	The CPU 68332 PCB provides an RS232 interface in the Evita. The interface is labeled COM1. Optocouplers electrically isolate the RS232 interface from the machine.
3.1.4	ILV port	The ILV interface allows for independent lung ventilation with two machines. The ILV interface is not electrically isolated. Pin 3 of the ILV interface is pro- vided with a filler plug. This filler plug prevents confusion with the RS232 interface.
3.1.5	Driver	The driver adjusts the access times between the 68332 CPU, the clock and the DUART.
3.1.6	Clock	The clock makes sure that the current time is displayed. The clock is battery- buffered. This is to ensure that the clock is supplied with the required operat- ing voltage after the device is switched off.
3.1.7	DUART	The DUART (Dual Universal Asynchronous Receiver/Transmitter) has two serial interfaces and digital inputs and outputs. The SpO2 module and the CO2 module are connected to the serial interfaces.
3.1.8	DC/DC converter	The DC/DC converter generates the voltage "+5 V ISO" required for the inter- face. The input voltage of the DC/DC converter is +5V.
3.1.9	CAN interface	The CAN interface is a fast, serial interface. Via the CAN interface the control unit can communicate with the electronic and pneumatic assemblies. The transmission rate is 800 kbit/s.
3.1.10	Bus driver	Via the bus driver, the signals from the address bus, the data bus, and the control bus are transmitted to the motherboard. The 68332 CPU uses the bus driver to communicate with the optional printed circuit boards installed on the motherboard (currently the Pediatric flow PCB (optional Neoflow feature)).
3.1.11	Sound generator	The sound generator controls the loudspeaker in the control unit. The sound generator includes the sound volume control and the tone generation of the loudspeaker. The sound volume is controlled by the DUART.
3.1.12	Reset logic	The CPU 68332 PCB uses a reset signal to reset (restart) the control unit and the pneumatic assembly. A reset is also triggered when the +5 V operating voltage is not reached (too low or too high).
		The pneumatic assembly can also send a reset signal to the CPU 68332 PCB, thus resetting the CPU 68332 PCB. The reset logic controls and displays the resets.

3.2 CO2 Carrier PCB

The CO2 Carrier PCB is integrated in the electronic assembly of the Evita. The printed circuit board includes the mount and the electrical isolation of the CO2 module and the SpO2 module, the mains voltage failure logic, the temperature measurement and the voltage monitoring.



Figure 4 Block diagram of the CO2 Carrier PCB







Figure 6 Block diagram of the CO2 Carrier PCB, part 2

3.2.1 Electrical isolation Optocouplers electrically isolate the printed circuit boards from the interfaces. The printed circuit board is provided with plug-in contacts. The SpO2 and CO2 modules are plugged into these plug-in contacts. The X3 connector is part of the optional "SpO2" module and is not equipped.

3.2.2 Mains voltage failure logic The mains voltage failure logic monitors the mains voltage supply. In the event that a mains voltage failure occurs while the machine is operating an audible alarm will sound.

3.2.3	Voltage monitoring	The microprocessor monitors the following voltages:
		– 15 V
		– +10 V
		– +24 V
		– +12 V
		– +5 V
		All voltages are present at a voltage divider. An A/D converter reads out the respective voltages. The CPU 68332 PCB reads out the AD converter.
3.2.4	Measurement of temper- ature and pressure	A temperature sensor measures the current temperature. The temperature sensor is a thermistor (NTC). A temperature hybrid outputs the respective analog voltage value. The output signal from the temperature hybrid is transmitted to an A/D converter which converts the analog voltage value in a digital value. The CPU 68332 PCB reads out the digital value.
		A pressure sensor measures the current ambient pressure. The output signal from the pressure sensor is transmitted to an A/D converter which converts the analog voltage value in a digital value. The CPU 68332 PCB reads out the digital value.
3.3	CO2 measurement	The CO2 measuring system comprises three modules: – CO2 sensor

- Processor Board PCB
- Power Supply PCB





3.4 CO2 sensor

The CO2 sensor comprises the CO2 measuring unit and a microprocessor system. A lamp generates a light spectrum of up to 4.5 m. This light travels through the cuvette and two sapphire discs and reaches the detectors. The detectors emit electrical signals depending on the CO2 concentration. The microprocessor analyzes these signals and transmits them to the Processor Board PCB via an RS232 interface.

The CO2 measuring unit is kept at a constant temperature to avoid condensation.



Figure 8 Sectional view of the CO2 sensor



Figure 9 Block diagram of the CO2 sensor

3.5 Processor Board PCB

The Processor Board PCB controls the heating (CO2 measuring unit) and the measured-data transfer of the CO2 sensor.

EvitaXL





3.6 Power Supply PCB

The Power Supply PCB provides the supply voltages for the Processor Board PCB and the CO2 sensor. The supply voltages are electrically isolated from the machine. The CO2 measured values are transmitted to the CPU 68332 PCB via the Power Supply PCB.





Function Description

3.7 Power pack

The switched-mode power pack provides the following output voltages:

- +24 V
- +15 V
- -15 V
- +12 V
- +5 V

The output voltages are short-circuit-proof.

3.7.1 DC module (option in converted Evita 4 or Evita 2 dura units)

The DC module makes sure that the device is powered in the event of a mains voltage failure.



Figure 12 Schematic circuit diagram of the switched-mode power pack

EvitaXL

3.8 Graphic Controller 8 PCB

The Graphic Controller 8 PCB is fitted to the control unit. Connection to the electronic assembly is made via a 15-pin SUB-D connector.

The printed circuit board includes an independent processor system, the voltage supply of the control unit, the interface to the electronic assembly and the display control.





3.8.1	MPC 850	The power PC "MPC 850" has a maximum clock frequency of 80 MHz. The bus frequency rate is 40 MHz.
		The MPC 850 contains the following components:
		– 32-bit RISC microprocessor
		 2 kB data cache, 1 kB program cache
		 MMU (Memory Management Unit)
		 On-chip emulation debug mode
		 Timer, bus monitor, watchdog, interrupt controller, closk synthesizer, JTAG interface
		 Memory controller
		 CPM (Communications Processor Module)
		 4 baud rate generators
		 2 SCC (Serial Communication Controller)
		 2 SMC (Serial Management Controller)
		– 1 SPI
		– 1 I2C
		 3.3 V operation with 5 V-tolerance inputs
3.8.2	SDRAM	The RAM has a memory area of 16 MB. This memory area is provided by two SDRAMs with 4 Mbit x 16 bit each.
3.8.3	Flash EPROM	The program of the control unit is stored in the flash EPROM. It has a total capacity of 8 MB. The flash EPROMs are connected to the MPC 850 through a logic circuit in order to avoid unintentional changing of the contents. The logic circuit avoids that data is written without a connector being connected to the "COM1" interface.
3.8.4	LED Control	The seven keys on the membrane keypad have a status LED each. In order to activate a status LED, the MPC 850 sends signals to a parallel port. This parallel port triggers the respective driver of the status LED.
3.8.5	Key inputs	All keys are connected to ground. For interference suppression purposes, the key inputs are provided with ferrite cores and suppressor capacitors.
		The control knob key and all other key are connected to the MPC 850.
3.8.6	Loudspeaker monitoring	The MPC 850 monitors the current across the loudspeaker. To this end, a 5- ohm resistor is looped in the loudspeaker lead as a combination of eight 10- ohm resistors. Capacitors tap the voltage present at the 5-ohm resistor and supply it to a differential amplifier. The amplified signal is available at a com- parator.
3.8.7	Reset generator	The 3.3-V reset generator DS1819A delivers the reset signals for the MPC 850.
3.8.8	CAN controller	The CAN controller 82527 is connected to the MPC 850 through level converters. The CAN controller is supplied with 5 V.

3.8.9	ADuC 812	The microcontroller ADuC 812 controls the control knob and the touchscreen.
		Rotary transducer input
		The rotary transducer supplies two signals that are phase-shifted by 90 degrees. The microcontroller ADuC 812 scans both signals periodically (at approx. 1 ms period duration) thus determining the position of the rotary transducer.
		Touchscreen control
		The glass pane used on the control unit is coated and conductive. In front of it, there is a foil which is also coated and conductive. The foil's coating resistance is approx. 1 kohm. Electrodes are attached to the foil (at the top and bottom) and to the glass pane (at the left and right).
		When a finger touches the glass pane, the front foil is pressed against the glass pane thus creating a conductive connection.
		The position of the finger is calculated as follows:
		A current flows from top to bottom. The glass pane is switched to high resis- tance. The voltages are measured from top to bottom. The conductive, coated foil functions as a voltage divider the tap point of which is formed by the contact point of the finger. The Y coordinate is determined by measuring the voltage.
		Now the foil becomes high-resistant and a current flows across the glass pane from the left to the right instead. The X coordinate is determined by measuring the voltage.
		The glass pane has a 8-wire design. That means that two electrodes each are attached to the edges. Changes in resistance on the glass pane itself or in the leads have no influence on the measurement result. A 4-wire measurement allows a current to flow across two electrodes, the reference voltage is then measured at the other two electrodes.
3.8.10	Graphic controller	The graphic controller MB86290A generates graphic signals and includes character operations. The graphic controller is connected to the MPC 850, it occupies an address area of 64 MB. The output signals are analog. The video memory has 16 MB.
3.8.11	Voltage supply	The control units uses the 24 V supplied from the power pack to generate the following voltages: - 12 V (15" TET display), can be switched off with MPC 850
		- 5 V (ADuC 812)
		– 3.3 V (MPC 850)
		 2.5 V (grafic controller)
3.8.12	Switched-mode regula- tor	The switched-mode regulator (DC/DC converter) converts the 24 V into 12 V. A second switched-mode regulator (DC/DC converter) converts the 24 V into 5 V.
3.8.13	Series regulators	Two series regulators filter any interference from the DC/DC converter. As a result, the voltages (3.3 V and 2.5 V) do not receive any clock signals from the DC/DC converter.

3.8.14 Comparators

The comparators make sure that the 5-V converter only starts operating at an operating voltage of approx. 17 V or higher.

3.9 Communication PCB (optional)

The CommunicationPCB is integrated in the electronic assembly of the machine. The printed circuit board includes an independent microprocessor system, the voltage supply of the interfaces, an internal CAN interface, an external CAN interface, two RS232 interfaces and two analog outputs.



Figure 14 Block diagram of the Communication PCB

3.9.1	Microprocessor system	The microprocessor system comprises the 68332 CPU, one 512-KB RAM, and one 1-MB flash EPROM (electrically programmable and erasable read-only memory).
3.9.2	Uprog	Uprog generates the voltage required for programming the flash EPROMs. An enable logic prevents unintentional supply of the programming voltage to the flash EPROMs.
3.9.3	Reset	The reset logic generates a defined reset after power-up. The CPU 68332 can be reset by the CPU 68332 PCB.

EvitaXL		Function Description
3.9.4	CAN interfaces	The Communication PCB is not connected to the data bus of the CPU 68332 PCB. The data are transmitted via an internal CAN interface (Controller Area Network – fast, serial interface). The external CAN interface is electrically isolated from the machine. Electrical isolation is made by means of optocouplers.
3.9.5	RS232 Interfaces	The Communication PCB provides an RS232 interface in the machine. The interfaces are labeled COM2 and COM2. The interfaces are electrically iso- lated from the machine. Electrical isolation is made by means of optocouplers.
3.9.6	Analog outputs	The analog outputs supply voltages between 0 V and 4.095 V. The assignment of analog outputs is freely selectable. The resolution of the output voltage is 1m V per bit.
3.9.7	DC/DC converter	The DC/DC converters generate +5V ISO each for the voltage supply of the interfaces. The input voltages of the DC/DC converters are +5 V.
3.9.8	EEPROM	The EEPROM stores internal data of the interface. The EEPROM has a 2 kB capacity.

EvitaXL

3.10 Pediatric Flow PCB (optional)

The Pediatric Flow PCB is integrated in the electronic assembly. The printed circuit board has two flow measuring channels for connection of the Babylog flow sensor, one four-channel multiplexer, one 12-bit A/D converter, and one interface to the CPU 68332 PCB.



Figure 15 Block diagram of the Pediatric Flow PCB

3.10.1 Flow Measurement The Babylog flow sensor measures the flow. The Babylog flow sensor is equipped with two measuring wires. One measuring wire is covered by a plastic bar, the Babylog flow sensor recognizes the direction of the flow. A measuring bridge analyzes the flow.

- **3.10.2 Multiplexer** The multiplexer consists of four analog-value selectors. The software controls the analog-value selectors. The multiplexer transmits the flow sensor measurement signals "CGND and Uref LIMIT" to a buffer (AMP). Then the measurement signals are transmitted to an A/D converter (ADC).
- **3.10.3** A/D converter The input voltage of the A/D converter ranges from 0 V to 10 V. The A/D converter converts the flow measurement signals into digital data. The CPU PCB controls the A/D converter and the multiplexer via an interface (DATA OUT, DATA IN). The voltage drop across the multiplexer, the buffer and the A/D converter is measured using the "UREF LIMIT" reference voltage and can be taken into account when measuring the flow.

3.11 IFCO Carrier PCB for additional optional features The IFCO Carrier P supports the nurse on one additional ambi

The IFCO Carrier PCB functions as carrier board for other optional features; it supports the nurse call and Remote Pad (cable remote control) functions and one additional ambient pressure sensor (optional).



Figure 16 Block diagram of the IFCO Carrier PCB

3.11.1 Nurse Call The nurse call transmits top-priority alarms (!!!) displayed on the screen to, e.g., a central station. An alarm is also transmitted if the internal loudspeaker for audible alarms fails. The alarm is reset automatically as soon as the cause of the alarm disappears. The alarms are suppressed during the boot phase (start phase of the device).

The central station alarm signaling is carried out by relay contacts. The alarm status can be determined by scanning these relay contacts. An alarm is considered to be triggered if contacts 3 and 5 of the nurse call connectors are closed by the relay contacts. These relay contacts are electrically isolated from the rest of the electronics. A bistable relay is used in order to keep this alarm signaling function active even when the power supply fails.

The nurse call function monitors the contacts of the power switch thus being able to detect whether the device is switched on or off. A power failure in the device can thus be detected. This alarm (power failure) can be cancelled by switching the device off. However, the device must be switched off within a defined time window. This time window is dependent on the charge status of specific capacitors. The time window is typically 2 minutes.

3.11.2 Remote PadThe Remote Pad is a cable remote control that is connected to the IFCO Carrier PCB. The Remote Pad has 6 key to control the ventilator. The Remote Pad is also provided with alarm LEDs. The inputs and outputs of the Remote Pads are short-circuit protected and protected against inadvertent wrong connections.

Pressure sensorThe measured value supplied by the ambient pressure sensor on the IFCO(optional)Carrier PCB is evaluated by the safety software. This measured value is used
to monitor the ambient pressure sensor on the CO2 Carrier PCB.

3.11.3

3.12 SmartCare PCB (option)

The SmartCare PCB allows the patient to be weaned much faster from the ventilator.

EvitaXL

The PCB has a microcontroller, an EEPROM, an SDRAM, a flash EPROM, a reset generator and "CAN" and "Ethernet" interfaces.



Figure 17 Block diagram of the SmartCare PCB

3.12.1 Microcontroller The microcontroller controls the hardware of the SmartCare PCB. 3.12.2 EEPROM The serial, rewritable memory module "EEPROM" stores the SmartCare PCB's data. The EEPROM has a 16 kByte capacity. 3.12.3 SDRAM The SDRAM is a synchronous RAM with a capacity of 32 MBytes. 3.12.4 Flash EPROM The flash EPROM (electrically programmable and electrically erasable ROM) has a capacity of 16 MBytes. 3.12.5 Reset generator The reset generator generates reset signals for the microcontroller. 3.12.6 Ethernet LTX 905 The "Ethernet" interface allows the communication with external devices. The "Ethernet" interface is electrically isolated from the machine. Electrical isolation is made by means of an optocoupler. 3.12.7 CAN Intel 82527 The "CAN Intel 82527" interface allows the communication with other microcontrollers in the device. The "CAN Intel 82527" interface is connected directly to the Motherboard PCB.

3.13 Pneumatic Controller PCB

- The Pneumatic Controller PCB is located in the pneumatic assembly. The printed circuit board provides the following functions:
 - Supply pressure measurement
 - Inspiratory and expiratory airway pressure measurement
 - Esophageal pressure measurement
 - Measurement of the flow
 - Measurement of the O2 concentration in the respiratory gas
 - Fan monitoring
 - Solenoid valve control
 - PEEP valve control
 - Interface to the HPSV Controller PCBs
 - CAN Interface
 - 68332 CPU with RAM and flash memory
 - Serial EEPROM for storage of device configuration



Figure 18 Flow measurement, Pneumatic Controller PCB



Figure 19 Pneumatic Controller PCB block diagram

loop control system. If gas flows past this wire, heat is dissipated. The larger

3.13.1 Pressure measurement Two pressure sensors measure the pressure. The respective pressure sensors measure the pressure during the inspiratory phase and the expiratory phase. The airway pressure signals are transmitted to a multiplexer. A zero calibration of the pressure sensors is carried out automatically every three minutes.
 3.13.2 Flow Measurement The flow sensor works according to the principle of a constant-temperature hot-wire flowmeter. Respiratory gas flows along a very thin, electrically heated platinum wire in a measuring tube. The wire is heated to a temperature of 180 degrees centigrade and kept at this temperature with a closed-
EvitaXL		Function Description
		the gas volume per time unit flowing past, the higher the heat dissipation. The heating power required to keep the wire at a constant temperature is a measure of the gas flow.
3.13.3	Cleaning of the flow sen- sor	The control for cleaning (glowing) of the flow sensor supplies a defined cur- rent to both measuring wires. During calibration the measuring wires begin to glow and burn any impurities. The sensor is cleaned when pressing the "flow calibration" key or automatically after drug nebulization. The cleaning process occurs during the next inspiratory phase or after 15 seconds at the latest.
3.13.4	Oxygen measurement	The O2 sensor works according to the galvanic cell principle. Oxygen mole- cules contained in the gas mixture to be measured diffuse through a plastic diaphragm into the electro-chemical cell and are reduced at the noble metal electrodes. At the same time a base electrode is being oxidized. The base electrode is spent by the oxidation process and thus determines the life of a sensor. The current flowing through the cell is proportional to the oxygen par- tial pressure in the gas mixture to be measured.
		Provided the pressure and temperature of the gas mixture to be measured are kept constant, the measured value will be directly proportional to the oxy- gen partial pressure. The O2 amplifier on the O2 Top PCB is mounted exter- nally on the inspiratory block The output signal is transmitted to the O2 Contact PCB via spring contacts. From there the output signal is transferred to the Pneumatic Controller PCB. The O2 cell is also connected to the O2 Top PCB via spring contacts.
3.13.5	Fan monitoring	At the front panel of the machine a fan is mounted to limit the temperature and the O2 concentration in the electronic unit of the pneumatic assembly in case of failure. The electronics monitors the fan.
3.13.6	Multiplexer	The multiplexer consists of 16 analog-value selectors. The software controls the analog-value selectors. The multiplexer routes the measurement signals from the pressure sensors, the O2 amplifier, the FAN UREF and the flow sensor to a buffer. The output signals of the buffer are then available at an A/D converter.
3.13.7	Solenoid valves	Two power drivers control the solenoid valves. A comparator monitors the power driver outputs for the nebulizer and the O2/Air switchover. The power drivers can be switched off from the electronic unit.
3.13.8	PEEP valve control	A voltage-controlled current source with power MOSFET controls the PEEP valve.
		A quad operational amplifier serves to adapt the D/A converter output signal to the current range of the PEEP valve. The CPU controls the D/A converter.
		The PEEP valve control is calibrated to the PEEP valve. The calibration data are stored in the serial EEPROM.
3.13.9	HPSV interface	The status lines of the HPSV Controller PCB are led to the Pneumatic Con- troller PCB via the pneumatics motherboard. Two bus drivers transmit the data to the data bus of the CPU. The data are transmitted to the HPSV Con- troller PCB by two power drivers. The data are accepted by "power swing" of the respective chip-select pin.

3.13.10 CAN Interface	The CAN interface comprises a CAN controller and a series-connected driver. The CAN controller is directly connected to the data bus of the CPU. The control unit, the electronics and the pneumatics communicate via a CAN interface. The transmission rate is 1 Mbit/s.
3.13.11 Microprocessor system	The microprocessor system on the Pneumatic Controller PCB consists of a 68332 CPU, a 256 kByte flash EPROM (electrically programmable and eras- able read-only memory) and a 256 kByte RAM.
3.13.12 Serial EEPROM	The serial EEPROM stores the data of the pneumatics. The EEPROM has a 128 Byte capacity.

3.14 HPSV Controller PCB

The pneumatic assembly contains two identical HPSV Controller PCBs. The board slot determines which of the boards is assigned to O2 and AIR. The HPSV Controller PCB comprises the following functions:

- Microcontroller with EPROM and RAM
- A/D converter for measurement of supply pressure
- D/A converter for current set-point specification
- Closed-loop circuit for current control
- Power transistor (power source)

Note: The characteristic of the HPSV cartridge is stored in the cartridge itself. The HPSV Controller PCB reads out this chracteristic.



Figure 20 HPSV Controller PCB with HPSV

- 4 Pneumatics
- The machine needs a driving gas pressure (AIR and O2) of 2.7 to 6 bar.

The pneumatic assembly consists of the following subassemblies:

- Dräger gas connection block/FAS gas connection block
- Parallel mixer or mixer block
- Pressure sensors
- PEEP/PIP valve
- Inspiratory block
- patient system





Figure 21 Layout of the pneumatic subassemblies

4.1 Gas connection block The gas connection block comprises the "O2" gas connection (M12x1 female) and the "compressed air" connection (M 20x1.5 male). The following connections are available: NIST, DISS (USA) and DIN. The connections are fitted with filters F1.1 and F1.2 (metal fiber web). The check valves D1.1 (AIR) and D1.2 (O2) prevent the gas from flowing back into the central gas supply system.

The pressure regulators DR1.1 and DR1.2 are set to 2 bar. The control gas flows past the DR1.1 to the 3/2-way valve Y1.1, from there to the emergency valve Y1.3, to the PEEP/PIP valve Y4.1 and finally to the emergency valve Y3.1.

The gas also flows to the expiratory pressure sensor S6.2 (purge flow) via the restrictor R1.1 (0.08 L/min).

Gas flows to the nebulizer via the 3/2-way valve Y1.4, if appropriately adjusted.

In the event of "AIR" supply failure, the machine will switch over to O2 supply (O2 switchover function)





Function Description

4.1.1 Characteristics of various gas connection blocks

- Characteristics of the FAS gas connection block:
- Connector plate
- Pressure regulator





Figure 23 Dräger gas connection block (left) and FAS gas connection block (right)

Legend

DR1.1	AIR pressure regulator
DR1.2	O2 pressure regulator
Y1.1	3/2-way solenoid valve, O2/AIR
Y1.2	3/2-way solenoid valve pressure, sensor calibration inspiration
Y1.3	3/2-way solenoid valve, venting
Y1.4	3/2-way solenoid valve, nebulizer

The gas connection blocks are interchangeable since their mounts for attaching to the machine are identical.

External distinguishing feature of the "Dräger gas connection block":

- four fixing screws

External distinguishing feature of the FAS gas connection block:

two fixing screws

EvitaXL



Figure 24 Dräger gas connection block (left), FAS gas connection block (right)

4.2 Parallel mixer The parallel mixer is a fast, electrically controllable proportional valve for gas flows between 5 and 180 L/min at supply pressures of 3 to 6 bar. Partial flows of less than 5 L/min are pulsed at a constant flow of 5 L/min. The supply gases compressed air (AIR) and oxygen (O2) available at the parallel mixer have a supply pressure of 2.7 bar to 6 bar. The parallel mixer mixes the two gases in accordance with the set parameters. The parallel mixer supplies the inspiratory gas to the patient.

The parallel mixer consists of the following components:

- Mixer connection block
- 1 cartridge valve with displacement sensor system for compressed air (AIR)
- 1 cartridge valve with displacement sensor system for oxygen (O2)
- 2 supply pressure sensors measuring the inlet pressure of the supply gases

The HPSV Controller PCBs control the parallel mixer electrically. The control signals are transmitted to the parallel mixer via the Pneumatic Motherboard PCB and the Pneumatic Controller PCB.



Figure 25 Parallel mixer

4.2.1 Mixer connection block

The two cartridge valves are mounted to the mixer connection block. The inspiratory gases in the mixer connection block are supplied to the respective cartridge valve. The respiratory gas available at the outlet of the cartridge valves is mixed in the mixer connection block and supplied to the inspiratory unit.

4.2.2 Cartridge valves (HPS valves O2 and AIR) The cartridge valve (HPS valve = HPSV = high-pressure servo-valve) supplies a defined amount of gas to the patient in accordance with the preset adjustment parameters for e.g. inspiration, trigger pressure, leak flow compensation.



Figure 26 Cartridge valve (HPS valve O2; HPS valve AIR)

During expiration the supply gas is available at the cartridge valve and at the supply pressure sensor In the cartridge valve the ball "A" is pressed into the valve seat "B"; this action closes the cartridge valve (see the following Figure).

During inspiration, the drive system applies a current to the cartridge valve. The drive system is equipped with a coil working according to the principle of a moving coil of the type used e.g. in loudspeakers. The plunger is deflected in proportion to the supplied current and pushes ball "A" out of valve seat "B". This causes an annular gap. The cartridge valve opens and gas flows to the patient.



Figure 27 Sectional view of the cartridge valve

Consequently, the size of the annular gap determines the gas flow. The annular gap between the ball and the valve seat is determined by the deflection of the plunger (i.e., by the displacement of the ball from the valve seat). The deflection is measured and controlled by a displacement measuring system.

A supply pressure sensor at the cartridge valve measures the supply pres- sure of the gas to be dosed. The supply pressure at the cartridge valve is also used in the calculation of the gas flow. Twice the pressure for the same width of gap will produce twice the flow.
A larger width of gap (annular gap) or a higher supply pressure results in a higher flow.
The flow is determined indirectly by the displacement signal and the supply pressure. The cartridge valves supply a flow of 5 to 180 L/min. The plunger displacement depends on the supply pressure. The supply pressure is between 3 and 6 bar absolute.
For a flow requirement of 120 L/min and a supply pressure of 5 bar the dis- placement will be approx. 0.2 mm. For a flow requirement of 120 L/min and a supply pressure of 1.5 bar the displacement will be approx. 0,6 mm.
As the parallel mixer must permanently operate with high precision over the full flow range, it is measured and the data (non-linearity) are stored in an EEPROM on the Displacement Sensor PCB.

Therefore no calibration is required when the cartridge valves are replaced. The cartridge valves must not be interchanged since they are specifically dimensioned and fitted for each individual gas.

4.2.3 Displacement sensor system The displacement sensor system measuring the position of the plunger in the valve is integrated in the cartridge valve. The displacement sensor system consists of the displacement measuring system and the Displacement Sensor PCB.

The displacement measuring system is a differential transformer. The a.c. voltage applied to the primary winding of the transformer has a frequency of approx. 1 MHz. The two secondary windings are switched such that their output voltages balance out. If the ferrite core (plunger of the cartridge valves) moves in the differential transformer, the output voltage of the displacement sensor system will change.

As the displacement output signal is not linear to the gas flow, the characteristic of the cartridge valve is measured and stored in the EEPROM. The micro-controller on the HPSV PCB thus balances the non-linearity of the cartridge valve.

The two circuits of the cartridge valves of parallel mixers are operated asynchronously in parallel (AIR, O2). To avoid beat interferences, the frequencies of the two oscillators must differ by a minimum of 200 kHz. Therefore the two cartridge valves have two different frequencies. The cartridge valves are measured at a special test stand. **4.2.4 Supply pressure sensors** The supply pressure sensors are calibrated to absolute pressure (0 bar). They measure the inlet pressure of the supplied gas. The supply pressure sensor is fitted with a P/U converter generating a pressure-dependent output voltage.

Measuring range	0 -7 bar
Sensitivity	1.58 V/bar 8 mV/bar
Offset voltage	300 mV 30 mV

The supply pressure sensors are linked via a flex-strip to the Displacement Sensor PCB. The Displacement Sensor PCB is installed in the cartridge valves.

4.2.5 Airway pressure sensors The pressure sensor mount comprises the airway pressure sensors S6.1 for the inspiratory side and S6.2 for the expiratory side. During inspiration, S6.1 monitors the airway pressure "high" (Paw high) and the airway pressure "low" (Paw low).

Measuring range	140 mbar
Sensitivity	36.5 mV/mbar 0.3 mV/mbar
Offset voltage	1.74 V 0.04 V

Calibration of the airway pressure sensors

The solenoid valves Y6.1 and Y6.2 expose the relevant airway pressure sensors to atmospheric pressure at specific time intervals, consequently, the airway pressure sensors are automatically calibrated. The airway pressure sensors S6.1 and S6.2 are zero-calibrated every 3 minutes. (Calibration of the possible electronic zero drift). To do so, the solenoid valves Y6.1 and Y6.2 are subsequently exposed to atmospheric pressure and the airway pressure sensors automatically calibrated.





4.3	Inspiratory block	The safety valve D3.3 limits the pressure in the inspiratory line to 100 mbar max.
		In the event of compressed air failure or mains voltage failure the pneumati- cally controlled emergency air valve Y3.1 will open so that the patient can breathe ambient air passing the filter F3.1. The non-return valve D3.1 pre- vents rebreathing of the air through the inspiratory line. The spring-loaded non-return valve D3.2 allows pressure to drop if valve Y3.1 opens.
		In the case of emergency air spontaneous breathing the patient can exhale through the expiratory valve Y5.1 on account of the spring loading (5 mbar) thus preventing rebreathing.
		The inspiration block is provided with the plug-in connection for the oxygen sensor.
		The restrictor R1.2 limit the drug nebulizer flow to 9 L/min.
4.3.1	Emergency air	If the gas supply or the voltage supply fails, the emergency air valve Y3.1 will no longer be controlled. The patient can breathe spontaneously through filter F3.1, non-return valve D3.1 and emergency air valve Y3.1.



Figure 29 Emergency air function diagram





4a

3a



Legend

1-1a	Patient connection (inspiration)
2-2a	Nebulizer port
3-3a	Esophageal pressure P1
4-4a	Esophageal pressure P2
5	10 mbar pressure relief valve D3.2
6	100 mbar pressure relief valve D3.3 (adjustable; NOTE: until middle of '96 valve slot milled in the block, after '96 separate piece screwed onto the block)
7	Emergency pressure relief mechanism control
8	Pressure measurement (inspiration)
9	O2 calibration control
10	O2 sensor chamber (behind the O2 amplifier)
11	O2 calibration purge flow outlet
12	Emergency air non-return valve D3.1 in O2 amplifier
13	O2 calibration diaphragm lattice Y3.3 with R3.1
14	Valve Y3.1
15	O2 amplifier

4.4 Patient system

The expiratory gas flows from the patient directly to the expiratory valve Y5.1. The copper measuring line at the 8a connection has a germicidal effect and connects the expiratory side to the pressure sensor S6.2.

The expiratory valve has a transmission ratio of approx. 1:1. The non-return valve D5.1 allows flow in one direction only and makes sure that gases do not travel backwards. The flow sensor S5.1 measures the expiratory flow.



Figure 31 Patient system function diagram



Figure 32 Sectional view of the patient system

The ratio between the control pressure at the 7a connection of the PEEP/PIP valve and the resulting pressure at the expiratory port is linear to the following values:

Control pressure of 3 mbar => expiratory pressure of 0 mbar

Control pressure of 33 mbar => expiratory pressure of 33 mbar

4.4.1 Patient system with RS water trap

The water trap avoids flow measurement faults caused by water droplets. Such faults may occur if the water traps on the ventilation hoses are not positioned at the lowest possible point. In this case the condensation water is collected in the water trap of the patient system.

The collector jar of the water trap can be removed during operation. The opening to the patient system is sealed automatically.









4.5 PEEP/PIP valve The PEEP/PIP valve Y4.1 consists of a diaphragm valve acting as a flowcontrol device and the linear drive whose plunger closes the diaphragm valve. A coil drives the PEEP/PIP valve Y4.1. The preset ventilation defines the settings. A computer program processes these settings, and the coil is driven by an appropriate current. The PEEP valve opens and adjusts a pressure proportional to the adjusted electric current (Note: 0 mA will correspond to -1 mbar, 500 mA to 120 mbar).

The PEEP/PIP valve Y4.1 controls the expiratory valve Y5.1 in the patient system via a servo-line. The solenoid valve Y1.3 supplies control gas to the restrictor R4.1 and to the PEEP/PIP valve Y4.1. The non-adjustable restrictor R4.1 is set to a flow of 3.5 L/min.

Depending on the PEEP setting the plunger coil is activated causing an appropriate servo-pressure to be applied to the diaphragm of the expiratory valve.

The software compares the preset and measured airway pressures. This comparison is a measure of the Pneumatic Controller PCB's control action on the PEEP/PIP valve. The PEEP/PIP valve is calibrated to the electronics. The calibration data are stored on the Pneumatic Controller PCB.



Figure 35 PEEP valve function diagram



Figure 36 Section view of the PEEP valve

5 Theory of operation

5.1 AIR supply

AIR flows through the filter F1.1 passing the non-return valve D1.1 reaching the mixer and flow control unit (pressure sensor S2.1 and HPSV Y2.1). At the same time, the gas flows to the 3/2-way solenoid valve Y1.1 via the pressure regulator DR1.1 which is set to 2 bar. From here the gas flows through the 3/2-way solenoid valve Y1.3 to the emergency air valve Y3.1 which closes. Furthermore, the gas passes the restrictor R4.1 to reach the PEEP/PIP valve Y4.1 and from there – depending on the setting – to the expiratory valve Y5.1. Finally, the gas passes the restrictor R1.1 to flow to the expiratory pressure sensor S6.2 connecting line on the patient side. At this point, expiratory humidity is prevented from reaching the pressure sensor S6.2.



Figure 37 AIR supply function diagram

5.2 O2 supply

O2 flows through the filter F1.2 passing the non-return valve D1.2 reaching the mixer and flow control unit (pressure sensor S2.2 and HPSV Y2.2). At the same time, the gas flows to the 3/2-way solenoid valve Y1.1 via the pressure regulator DR1.2 which is set to 2 bar.





5.2.1 O2/AIR changeover valve

The 3/2-way solenoid valve Y1.1 switches under the following circumstances:

- If the "AIR" supply fails
- When the O2 sensor is calibrated
- When the nebulizer function (depending on FiO2 setpoint) is triggered

If these conditions are given, the servo-system will be supplied with O2.

5.3 Inspiratory phase

Depending on the settings (O2 concentration, inspiratory volume, frequency, Ti, inspiratory flow, inspiratory pressure) the cartridge valves (HPSVs Y2.1 and Y2.2) open. The gas flows to the patient through the inspiratory port. At the same time, gas flows to the O2 sensor S3.1 and to the safety valve D3.3; from there, it flows through the 3/2-way solenoid valve Y6.1 to the inspiratory pressure sensor S6.1.

The safety valve D3.3 is fixed to 100 mbar and serves as an additional safety device in the event of a complete failure of the electronic control.

When calibrating the O2 sensor S3.1, valve Y3.3 disconnects the sensor from the inspiratory gas. The O2 sensor S3.1 is supplied with 100% O2 from the valve Y1.2, the restrictor R1.3, the restrictor R3.1 and the valve Y3.2. The O2 concentration and the flow will not be influenced by the inspiratory gas.

The pressure sensors S6.1 and S6.2 monitor the inspiratory pressure. During the entire inspiratory time the PEEP/PIP valve Y4.1 provides pressure to the expiratory valve Y5.1.



Figure 39 Inspiratory phase function diagram

Function Description

- **5.4 Paw high alarm limit** In case the "Paw high" alarm limit is exceeded during inspiration, the HPSV Y2.1 and Y2.2 interrupt the gas flow. The PEEP/PIP valve Y4.1 is switched to expiration and the patient can exhale.
- 5.5 Emergency pressure relief mechanism (safety valve) In case the "Paw high" alarm limit is exceeded by 5 mbar, an additional safety valve, the so-called "emergency pressure relief mechanism" Y1.3 will open. As a result, the emergency air valve Y3.1 opens and the pressure is vented through the non-return valve D3.2.
- **5.6 Expiration** At the start of expiration the cartridge valves (HPSV Y2.1 and Y2.2) are closed. No gas flows to the patient. The PEEP/PIP valve Y4.1 is switched to the set PEEP value. The expiratory valve Y5.1 is also relieved and the patient can exhale through the non-return valve D5.1 and the flow sensor S5.1. The flow sensor S5.1 measures the expiratory flow.



Figure 40 Expiratory phase function diagram

EvitaXL

5.6.1 Drug nebulizer

After pressing the "drug nebulizer" button, the drug nebulizer is switched on for 30 minutes. At the same time, the solenoid valve Y1.4 switches through in the flow-active inspiratory phase. The restrictor R1.2 supplies the drug nebulizer with drive gas. At the end of the inspiratory gas supply phase, the solenoid valve Y1.4 also switches back. The minute volume remains constant while the flow setting is being corrected. After completion of the drug nebulization the flow sensor S5.1 is automatically glowed clean.

Note: the drug nebulizer needs a minimum inspiratory flow of 16 L/min.



Figure 41 Nebulizer function diagram

6	Gas mixture	One cartridge valve (HPSV) each controls the air flow or the oxygen flow, directly from the supply lines. Both sub-flows are joined in the parallel mixer and then supplied to the inspiration block. Depending on the set oxygen concentration (FiO2: 0.21 to 1.00), the total flow to be metered is split into an "AIR" sub-flow and an "oxygen" sub-flow. Sub-flows of less than 5 L/min are no longer metered continuously, but in pulses lasting at least 8 ms and with a constant flow value of 5 L/min. This results in a pulse/pause ratio that corresponds to the sub-flow value.
6.1	Correcting the oxy- gen cartridge valve	Corresponding to its operating principle, the cartridge valve (HPSV) meters a mass flow. Owing to the differing gas density values of the compressed air and the oxygen, different volume flows would therefore be applied if these values were not corrected.
		Owing to the different gas density values, the portion of the inspiration flow demanded of the oxygen cartridge valve is therefore increased by 5%.
6.2	Dependence on the supply pressure	With regard to the gas supply, the operating range of the machine is specified from 2.7 bar to 6 bar gauge pressure. The machine monitors this operating range with the aid of the absolute pressure sensors on the cartridge valves on the basis of the following criteria:
		 Supply pressure always higher than 1.2 relative/absolute
		 At no flow, supply pressure higher than 2.5 bar relative
		At low supply pressures (below 3 bar gauge pressure), the cartridge valve can no longer apply high flows without errors, i.e. the actually supplied flow is less than the demanded flow. The inspiratory flow must be limited for reasons of flow accuracy and thus also for oxygen concentration accuracy reasons.
		Prior to the start of an inspiration, the maximum inspiratory flow for spontane- ous inspiration of the cartridge valve is defined on the basis of the measured supplied gauge pressure:
		Psupply > 4 bar absolute> max. inspiratory flow = 180 L/min
		Psupply < 4 bar absolute> max. inspiratory flow = 150 L/min
		For mandatory inspirations, the inspiration flow is limited by the setting to 120 L/min.

7 Function diagram



Figure 42 Function diagram

Legend

AIR	Air connection	Y3.1	Emergency air valve
02	Oxygen connection	Y3.3	Inspiratory valve
		Y4.1	PEEP/PIP valve
F1.1	Filter	Y5.1	Expiratory valve
F1.2	Filter	Y6.1	3/2-way solenoid valve, insp.
F3.2	Filter	Y6.2	3/2-way solenoid valve, exp.
D1.1	Non-return valve	S2.1	AIR pressure sensor (HPSV)
D1.2	Non-return valve	S2.2	O2 pressure sensor (HPSV)
D3.1	Non-return valve	S6.1	Inspiratory pressure sensor
D3.2	Non-return valve 10 mbar	S6.2	Expiratory pressure sensor

D3.3	Non-return valve 100 mbar	S3.1	O2 sensor
D5.1	Non-return valve	S5.1	Flow sensor
DR1.1	AIR pressure regulator	R1.1	Restrictor 0.08 L/min/2bar
DR1.2	O2 pressure regulator	R1.2	Restrictor 9 L/min/2bar
		R1.3	Restrictor 0.4 L/min/2bar
		R3.1	Restrictor (hole in the dia- phragm in Y3.3) 0.25 L/min/1.4 bar
		R4.1	Restrictor 3.5 L/min/2bar
Y1.1	3/2-way solenoid valve, O2/AIR		
Y1.2	3/2-way solenoid valve c	alibration	O2 sensor
Y1.3	3/2-way solenoid valve, venting		
Y1.4	3/2-way solenoid valve, nebulizer		
Y2.1	AIR HPSV (high-pressur	e servo-v	alve) parallel mixer
Y2.2	O2 HPSV (high-pressure	e servo-va	alve) parallel mixer

Maintenance Procedures

- 1 Cooling-air filter
- 1.1 Replacing the cooling-air filter
- 1. Remove cooling-air filter (Figure 1 /1) from its slot on the back of machine.



Figure 1 Rear view of the machine

- 2. If cooling-air filter is dirty, clean it in warm water with detergent added; allow to dry well.
- 3. If cooling-air filter is faulty, dispose of it according to local waste disposal regulations.
- 4. Insert cooling-air filter in slot (white side inside), taking care not to crease it, see Figure 1/1.

- 2 Ambient-air filter
- 2.1 Replacing the ambient-air filter
- 1. Swivel inspiratory port down (Figure 2/1).



Figure 2 Front view of the machine, protective cover

2. Loosen screw (Figure 3/1) with an appropriate screwdriver and remove the protective cover.



Figure 3 Protective cover



Figure 4 Ambient-air filter

cover.

If ambient-air filter is dirty, clean it in warm water with detergent added; allow to dry well. Continue with step 4.

3. Remove the ambient-air filter (Figure 4/1) from the tabs on the protective

If ambient-air filter is faulty, dispose of it according to local waste disposal regulations.

- 4. Slide new or cleaned ambient-air filter under the tabs on the protective cover, making sure the fitting position is correct, see Figure 4/1.
- 5. Replace the protective cover on the machine, and tighten the screw (Figure 3/1) with an appropriate screwdriver.

- 3 Lithium battery (CPU 68332 PCB)
- 3.1 Replacing the lithium battery (CPU 68332 PCB)
- 1. Switch on the machine and wait 2 minutes (Note: the capacitor for power supply to the memory modules (RAM) needs to be charged).
- 2. Write down the customer-specific settings (alarm limits, screen set-up, etc.) so that you can re-enter them after replacing the lithium battery.
- 3. Switch off the machine.

CAUTION

Electrostatic discharge may damage the electronic components. When handling electrostatic sensitive devices, use a static-dissipative mat and a static dissipative wrist strap.

- 4. Create ESD conditions.
- 5. Remove the screws (Figure 5/1) on the CPU 68332 PCB (Note: The CPU 68332 PCB might also be in another slot than that shown here).



Figure 5 Rear view of the machine

CAUTION

A data loss will occur if the RAM power supply is missing. Replace the lithium battery within one minute! When installing the lithium battery make sure that it is properly seated in its holder.

6. Remove the CPU 68332 from the machine.



Figure 6 Detail of the CPU PCB

8. Push the lithium battery into the battery holder (see Figure 6), making sure the polarity is correct.

7. Remove the spent lithium battery (Figure 6/1) from the battery holder and

dispose of according to local waste disposal regulations.

- 9. Push the CPU 68332 PCB into the machine.
- 10. Secure the CPU 68332 PCB to the machine using the screws (Figure 7/1).



Figure 7 Rear view of the machine

11. Check the electrical safety and the correct functioning of the device as described in the PMS procedure.

4 Real-time clock

4.1 Replacing the realtime clock

- 1. Switch off the machine.
- 2. Unplug the power plug of the machine from the socket-outlet.

CAUTION

Electrostatic discharge may damage the electronic components. When handling electrostatic sensitive devices, use a static-dissipative mat and a static dissipative wrist strap.

- 3. Create ESD conditions.
- 4. Remove the screws (Figure 8/1) (Note: The CPU 68332 PCB might also be in another slot than that shown here).

		A T
COMT RS232C		
		-
COM2 RS232C CC243 CAN	CAN Analog 🛁 🦉	
Sync.	Temp.	4
© Sp02 ▲ CO		

Figure 8 Rear view of the machine

5. Remove the CPU 68332 from the machine, and place it on the static-dissipative mat. 6. Carefully remove the real-time clock (Figure 9/1) from its IC socket.



Figure 9 Detail of the CPU PCB

7. Dispose of the spent real-time clock according to local waste disposal regulations.

NOTE

If, in the next work step, the real-time clock is not fitted correctly it will not display the time and date. Make sure the fitting position is correct when fitting the real-time clock in the next work step.

8. Firmly push the new real-time clock into its IC socket, making sure the fitting position is correct (coding), see Figure 10.



Figure 10 Real-time clock

9. Push the CPU 68332 PCB into the machine.



10. Secure the CPU 68332 PCB to the machine using the screws (Figure 11/1).

Figure 11 Rear view of the machine

- 11. Plug the power plug of the machine into the socket-outlet.
- 12. Switch on the machine.
- 13. Using the rotary knob, set the current time and date in the "Configuration" and "Basic Settings" menu.
- 14. Check the electrical safety and the correct functioning of the machine as described in the PMS procedure.
Service Diagnosis Mode

1	General			
		NOTE The complete	Service Diagr	osis Mode is only available in SW 6.n or higher
		The service dia time during op in English, no	agnosis mode eration. The S matter which l	enables on-line information to be read out at any rervice Diagnosis Mode user interface is always anguage version is used in the unit.
		The diagnosis	pages are spl	it up as follows:
		Front:	Readout of c	control panel data:
			Processor:	Keypad test
				Loudspeaker monitoring
				RAM Test
		Electronic:	Readout of e	electronics section data:
			Processor:	Voltage supply
				Internal unit temperature
				AWT sensor
				inspiratory pressure
				Ambient pressure sensor
				RAM/ROM test results
				Service connector
				CAN
				SW options
			Sensors:	Ambient pressure measurement
				Pediatric flow measurement
				CO2 measurement
				SpO2 measurement (not provided for the time being)
			add. HW:	IFCO-Carrier PCB

Pneumatic:	Readout of p	oneumatics section data:
	Processor:	RAM/ROM test results
		Fan monitoring
		Reference voltages
	Sensors:	4 airway pressure sensors
		Flow measurement
		Basic flow
		Pulse threshold
		O2 sampling
	Valves:	Setting of switching valves
		HPSV (supply pressure, status, O2 measure- ment, flow)
		PEEP/PIP valve (calibration values)
		Nebulizer flow
Other:	Special func	tions – error list and SW options:
	Logbook:	Readout of error list
	Option Release:	Enabling of options
	SW Info:	Software information, e.g. SmartCare.

2 Call-up of service diagnosis mode

Service diagnosis mode can be called up at any time in the course of operation. This mode only involves data readout; settings are not altered.

- Press the 'System Setup' key.
- Select the 'Service' menu.
- Enter the number code '4655'.
- Select the required menu.

System	Setup					X
System	Therapy					
Please en	ter access o	ode:				
						Overview
Code			1	2	3	Sound, Day/Night
EvitaXL Part-No.	8416923		4	5	6	Screen
06.10 ID 16448255	08.09.2005 5		7	8	9	户 Country
11336 h				0		Interface
						Service

Figure 1 Call-up of service diagnosis mode

3 Diagnosis 'Front'

3.1 Diagnosis page 'Processor'



Figure 2 Diagnosis page 'Processor'

Keys pressed:	Keypad test, determined by Graphics Controller PCB. 0 = key not pressed 1 = key pressed
	Important: Keys which affect screen content switch Evita out of diagnosis mode.
Sound detected:	Result of current measurement through loudspeaker in control panel. Detected flow of current through loud-speaker (= yes) continues to be displayed for 10 seconds.
RAM:	Output of result of RAM test, Graphics Controller PCB.

Evita XL

4 Diagnosis 'Electronic'

4.1 Diagnosis page 'Processor'

Service						
Front Electron	nic Pneur	natic Other				
5 V	5.137V	3806	GoldCap	8.783V	3206	
10 V	9.999V	3933	Piezo	0.000V	1	F
12 V	12.006V	2729	PowerSW	0.000V	0	F
15 V	14.832V	3102	Batt ext.	0.256V	844	S
-15 V -	15.120V	1528	AD1	2.499V	4094	
24 V	23. 947 V	3566	AD2	2.499V	4094	á
Тсри [С]	30.0	1.726V	EEPROM	ok		
Tair [C]	0.2	2.399V	Loss of data	no		
Pinsp [mbar]	10.5	0.573V	RAM	ok		
Pambient [mbar]	1032	2.085V	RAM-Comm			
ILV in / out		00/11				
Service connecto	or	no				
CAN component	s 1	1000000				
SW options	1_1	0000010_1	1100011_111111	00_0000000	0	
GoldCap (mF; LS	B)	225	???; ???	225 5	??	
GoldCap (Ohm; I	LSB)		???; ???	???		

Figure 3 Diagnosis page 'Processor'

5 V / 12 V / 15 V / -	Power pack supply voltages, measured on the CO2
15 V / 24 V:	Carrier PCB. Output of voltages and decimal values of
	AD conversion. There are no device error messages
	with the following values:

5 V = 4.5 V to 5.5 V	
------------------------	--

12 V = 9.0	3 V	to	14.4	V
------------	-----	----	------	---

15 V = 14 V to 18 V

-15 V = -12 V to -18 V

24 V = 19.2 V to 28.8 V

Important: The power pack is specified as follows:

5 V	=	5.07 V to 5.23 V
12 V	=	11.76 V to 12.24 V
15 V	=	4.55 V to 15.45 V
-15 V	=	-14.55 V to -15.45 V
24 V	=	22.32 V to 24.72 V

10 V:

10V reference voltage, generated from the 15V supply voltage on the CO2 Carrier PCB. Output of voltage and decimal value of AD conversion. Error message is given in the event of deviation greater than 4 %. Set value = 9.65 V to 10.35 V.

GoldCap:	Output of GoldCap voltage. The GoldCap capacitor supplies the substitute horn with current. Output of capacity and resistance value and decimal values of the AD conversion. After startup of the Evita '???' are displayed first. The first measured values will appear after approx. 10 minutes.
Piezo:	The current through the substitute horn is only mea- sured on the old CO2 Carrier PCB. Piezo detected, >30 decimals (0.018 V). Piezo off, <30 decimals.
PowerSW:	Output of voltage at mains-switch auxiliary contact. Voltage is important for detecting mains failure. Output of voltage and decimal value of AD conversion. Mains switch on, <164 decimals (0.100 V). Mains switch off, >512 decimals (0.312 V).
Batt ext.:	Output of voltage of externally connected DC voltage supply (option DC module). Output of voltage and dec- imal value of AD conversion.
AD1:	Output of AD converter reference voltage 1. Output of voltage and decimal value of AD conversion. Set value approx. 2.500 V.
AD2:	Output of AD converter reference voltage 2. Output of voltage and decimal value of AD conversion. Set value approx. 2.500 V.
TCPU [C]:	Temperature inside unit; measured on CO2 Carrier PCB. Output of temperature in °C and voltage. Permissible temperature range <65 °C. The message 'Fan malfunction !' is displayed as of 65 °C. (removed at 60 °C) and as of 70 °C the alarm 'Fan failure? !!!'
Tair [C]:	AWT sensor respiratory-gas temperature. Output of
	temperature in °C and voltage.
Pinsp [mbar]:	Second channel of the inspiratory pressure sensor.
Pambient [mbar]:	Air pressure currently being used in Evita (filtered). Output of air pressure in mbar and voltage. The permissible measuring range is = 600 mbar to 1100 mbar.
EEPROM:	Result of EEPROM test. EEPROM D22 is located on the CPU 68332 PCB.

Loss of data:	Output of data loss of RAM and EEPROM D22 on CPU 68332 PCB: no = no data loss yes = data loss.
	A data loss in the RAM may due to a flat battery on the CPU 68332 PCB.
	Following data loss, the Evita attempts to restore the data. If this proves possible, the red alarm 'Loss of data !!!' disappears from the screen. The Evita is OK if this alarm is no longer present when the unit is next switched on.
RAM:	Result of RAM test on the CPU 68332 PCB.
RAM-Comm:	Result of RAM test on Communication PCB. Valid values are only output if the PCB fitted (option).
ILV in/out:	Outputs only defined with ILV ventilation.
Service connector:	Output of whether or not RS 232 cable with service encoding is connected to COM1 on CPU 68332 PCB. no - service encoding not recognized, normal opera- tion yes - service encoding recognized, download of oper-
	ating software or external service mode via PC possible.
CAN components:	Output of PCBs connected to internal CAN at CPU 68332 PCB (0 = not fitted, 1 = fitted). Sequence is as follows from left to right: Graphics Controller PCB / Pneumatic Controller PCB / Communication PCB / 5 x undefined
SW options:	The released options are output in the following form: Output form: '1'_(option 0 - 7)_(option 8 - 15)_(option 16 - 23)_(option 24 - 31) Meaning: '1' = block 1 with 32 options (1 = option fitted,
	0 = option not fitted) The following options are currently provided in block 1: Option 7 = SpO2 measurement Option 8 = CO2 measurement Option 9 = DC power pack Option 10 = Comfort Breath (PPS and tube compensation)

4.2 Diagnosis page 'Sensors'

Front	Electronic	Pneumati	c Other					
Ambient Pr Pambient Raw P_Referenz Pair Poisp	essure z	[mbar] 1032 1032 1032 1032 1018 1020 6.4	[V] 2.083V	Paediatric f Flow/MV Flow NTPD/ Cal D/A ZeroOffs U_FLimit/U	low [L/min] 'S [L/min] [V] [LSB_DA] [LSB_AD] _Gnd [V]	measured 0.00 0.00 4.379 1272 -1792 3.833	shadow 0.00 4.384 851 -1792 0.000	Proces Sensor add. H
CO2 SW rev Head SN rev cal	DW H 6870	HL/ANDR 277 1.11		SpO2 SpO2 [%] Puls [bpm] Pleth Bargraph	0 0 0 0			
Status Selftest Føet Fwin	8800 FF 0 0	0000						

Figure 4 Diagnosis page 'Sensors'

Ambient Pressure:	Output of ambient pressure measurement values:				
	Pambient:	Air pressure used in Evita.			
	Raw:	Current raw voltage from air pressure sensor with calculated air pressure.			
	P_Reference:	Air pressure established as air-pressure reference value when Evita is switched on. If the deviation from PAir or PO2 is more than 8 %, the message 'Execute device check' is given.			
	Pair / PO2:	Air pressure of supply pressure sensors PAIR and PO2, determined in last device check.			
	Pinsp:	Second channel of the inspiratory pressure sensor.			

Paediatric	Pediatric flow measurement values. Valid values are only pro-
flow:	vided if option has been fitted.

CO2:	CO2 measurement values. A list of important CO2 measurement data is given below:				
	CO2:	Output of PCB type of the Signal Processor PCB of CO2 measurement on the CO2 Carrier PCB, e.g. 'DW HL' or 'Andros 4210'			
	SW rev:	Output of SW revision on Signal Processor PCB, e.g. '6870277 1.11' for Dräger software revision 1.11.			
	Tdet / Twin:	Output of decimal value of detector (Tdet) and window (Twin) temperature in CO2 sensor. Approx. 43 °C should have been attained in each case after 3 minutes. The temperatures may vary between 40 °C and 50 °C, depending on the operating conditions. 40 °C = 5326 41 °C = 5116 42 °C = 4916 43 °C = 4725 44 °C = 4543 45 °C = 4368 46 °C = 4201 48 °C = 3888 50 °C = 3602			

SpO2: Not provided for the time being.

4.3 Diagnosis page 'add. HW'

Servic	e			
Front	Electronic	Pneumatic	Other	
IFCO-C	Carrier PCB		[mbar]	[m\/]
PCB de	etected	yes	[mbar]	[iuv]
Remote Remote	e number e keys	1 000000		
Pambie SW-Ver	ent	1.00	0	0
311-76	151011	1.00		

Figure 5 Diagnosis page 'add. HW'

IFCO Carrier PCB: Output whether IFCO Carrier PCB (option) is fitted or not ('PCB detected').

Evita XL

5 Diagnosis 'Pneumatic'

5.1 Diagnosis page 'Processor'

Service								×
Front	Electronic	Pneumatic	Other					
Peripherie RAM EEPROM Fan	e of process ok ok ??1	- sor ???		ADC Uref5V ADCref DACref	5.014V 3.000V 3.001V	ok ok ok		Processor Sensors
								Valves

Figure 6 Diagnosis page 'Processor'

Periphery of processor:	Status of processor system of Pneumatic Controller PCB and fan for cooling electronics in pneumatics sec- tion. Error message is given in the event of deviation. Fan must be detected within 1 minute.
ADC:	Output of reference voltages of Pneumatic Controller PCB. Error message is given in the event of deviation. Target values: Uref5V = 4.9 V to 5.1 V ADCref = 2.9 V to 3.1 V DACref = 2.9 V to 3.1 V

5.2 **Diagnosis page** 'Sensors'

Service	9							×
Front	Electronic	Pneumati	c Other					
Pressure insp exp oes	e Sensors [mbar] 12.1 8.7 -47.6	[V] 2.100 2.421 0.000	Cal [V] 1.687 1.754 1.740	Flow-Sensor FlowExp	[L/min] 0.0 Uinop 0.795	[L/min] F0 4.037	NTPD/S 0.0 Fsig 2.018	Processo Sensors
aux Paw [mbar] (PCB CF R [mbar. C [mL/m	-47.6 PU) /L/s] ibar]	0.000 insp 12.2 12.3 0.0 0.0 0.35	1.740 exp 11.8 0.6 0.0	Cal DA Switch NTPS -> BTPS Quotient BTPS -> NTPD MVinsp [L/min]	2514 on 1.046 1.017 0.914 5.88			Valves
O2 Sens FiO2	or [Vol.%] 5397.6	[V] 2.5567 2.556	Cal [V] 2.729	Pmean O2 Cal P	[mbar] 6.7 1032			

Figure 7 Diagnosis page 'Sensors'

Pressure Sensors:	Output of measured values from the 4 airway pressure sensors insp., exp., oes and aux (oes and aux are envisaged for future options). Three values are output for each sensor:				
	[mbar]:	Pressure in mbar of calibrated sensor.			
	[V]	Output voltage of sensor; read in by Pneu- matic Controller PCB			

	Voltage = Calibration voltage + sensitivity
	Sensitivity = 36.5 ±0.3 mV/mbar
	Measuring range = 140 mbar
Cal [V]	Calibration voltage of sensor (zero at ambi- ent pressure) Set value = 1.74 ±0.50 V.
The pressure at	the Y-piece (Paw) is also calculated:
Paw [mbar]:	Pressure Paw calculated in each case by

w [mbar]: aw calculated in each ca insp. and exp. pressure sensor with allow-

difference between the two values is 5 mbar. (PCB CPU): Comparison value in mbar. Output voltage of inspiratory pressure sensor is additionally read in on CO2 Carrier PCB (second channel). The max. difference with respect to Paw is 5 mbar.

ance for tubing resistance. The permissible

	R [mbar/L/s]:	Tubing resistance: insp: From inspiratory port to the Y-piece.
		exp: From the Y-piece to the expiratory port.
		Top row: Measured in operation by safety software.
		Bottom row: Measured in device check.
		These values differ as they are measured at different flow rates.
	C [mL/mbar]:	Inspiratory tubing compliance.
Flow Sensor:	Measured values values are display stated. Explanation:	of the expiratory flow measurement. Note: All red under BTPS conditions unless otherwise
	BTPS	Based on 37 °C, ambient pressure + inspira- tory pressure, 100 % relative humidity. All Evita measured values and settings are based on BTPS.
	NTPS	Based on 20 °C, 1013 mbar, 100 % relative humidity. Corresponds to raw value of measured expiratory flow.
	NTPD	Based on 20 °C, 1013 mbar, dry. The mixer delivers the flow under the following condi- tions. The Evita setting under BTPS is con- verted to NTPD for the mixer.
	[L/min]:	Measured expiratory flow converted to BTPS. Corresponds to displayed value.
	[L/min] NTPD/S:	Raw flow value measured under NTPS or NTPD conditions
	[V]:	Flow measurement voltages in V:
		Uinop: Voltage set value >0.1 V, sensor probably OK if F0 and Fsig OK. 'Flow measurement inop. !!!' is displayed if voltage is <0.1 V.
		F0: 4-fold amplified bridge voltage of the flow measurement. Set value in no-flow condition = 4.04 V.
	Cal DA:	Fsig: 2-fold amplified bridge voltage of the flow measurement. Set value in no-flow condition = 2.02 V. Decimal DA converter value for calibration of flow sensor
		Set value = 2200 to 3200.

	Switch:	Microswitch for detecting position of flow sensor: on - Flow sensor in right-hand operating position. off - Flow sensor not ready for operation in replacement position
	NTPS -> BTPS:	Factor for converting raw value of measured flow NTPS to displayed flow BTPS value. Permissible deviation from 'quotient' = 5 %.
	Quotient:	NTPS/BTPS conversion factor calculated in a different way.
	BTPS -> NTPD:	Factor for converting flow BTPS to NTPD
	MVinsp [L/min]:	Inspiratory minute volume calculated from actuation signals of mixer. The measured expiratory minute volume may be max. 20 % more than this value. 'Flow measurement inop' is displayed if the difference is greater. The problem may be due to the following: - Flow measurement is faulty
		- Impermissible feed-in of external flow (e.g. by the nebulizer)
		- Mixer supplying excessive flow rate.
O2 Sensor:	Measured values surement. Output the sensor.	and calibration values for inspiratory O2 mea- voltage of O2 sensor is amplified directly at
	FiO2 [Vol%]:	Measured O2 concentration, pressure-compensated.
	[V]:	Amplified sensor voltage at Pneumatic Con- troller PCB input. This voltage is read in twice on this PCB. Refer to Cal [V] for permissible voltage range.
	Cal [V]:	Amplified sensor voltage during calibration to 100 vol.% O2. Permissible range: 1.257 to 5.644 V. Voltage too low: Sensor spent. Voltage too high: O2 measurement is faulty (O2 amplifier or Pneumatic Controller PCB)
	Pmean O2:	Current mean pressure in tubing system in mbar. Required for pressure compensation of O2 measurement.
	Cal P:	Air pressure allowance in mbar on calibrat- ing O2 measurement.

5.3 Diagnosis page 'Valves'

Service			×
Front	Electronic Pneuma	Other	
Valves Y1.1 Y1.2 Y1.3 Y1.4 Y6.1 Y6.2 Yx.x Yx.x Yx.x Yx.x Yx.x Yx.x	on on on off off off off off off off off	HPSVPEEPPressureStateO25.90bar0Air6.00bar0SettingSettingFiO230.0Vol.%Flow0.0L/minMeasuredFiO2FiO25397.6Vol.%	Processor Sensors Valves
Nebuliser Flow	r 6.0L/min		

Figure 8 Diagnosis page 'Valves'

Valves:

Actuation of switching valves in pneumatics section. For the valves Y1.1 and Y1.4 the actuation voltage is read back and output in the second value. Actuation and feedback must be identical.

on: Valve actuated.

off: Valve de-energized.

- Y1.1: O2/Air switching valve, off = Air.
- Y1.2: O2 measurement calibration valve, on = O2 measurement calibration.
- Y1.3: Safety valve, on = ventilation, off = safety shutdown active.
- Y1.4: Nebulizer valve, off = nebulizer off.
- Y6.1: Calibration valve for inspiration pressure sensor S6.1,
 - on = calibration.
- Y6.2: Calibration valve for expiration pressure sensor S6.2, on = calibration.
- Yx.x: Future option; a total of 10 valves can be connected to the Pneumatic Controller PCB.

HPSV:	Supply press set values for output, as we	ures and status of cartridge for O2 and Air. The flow and FiO2 and the measured FiO2 value is II.
	Pressure:	The absolute pressures for Air and O2 are output. Absolute pressure = rel. supply pres- sure + ambient pressure. Measuring range of supply-pressure sen- sors = 0 to 7 bar. Sensitivity = 1.58 V/bar ±8 mV/bar. Offset voltage = 300 mV ±30 mV.
	State:	Status messages from the respective HPSV Controller PCB for O2 and Air: 0 = no error.
		 2 = supply pressure less than 2 bar absolute, pressure supply failure. 1, 3 to 15 = errors on HPSV Controller PCB or in HPSV cartridge.
PEEP:	Values for PE	EP/PIP valve Y4.1
	Cal. values:	Calibration values for actuation of PEEP/PIP valve. Valve or PCB replacement must be followed by calibration.
	Setting PEEP:	PEEP setpoint input in mbar.
	Measured PEEP:	PEEP measured by airway pressure mea- surement in mbar.
Nebuliser	Flow value se	et on the nebulizer.

Flow:

6 Diagnosis 'Other'

6.1 Diagnosis page 'Lokbook' (error list)

Service				×
Front Electronic	Pneumatic Other			
First occurrence 06.12. 16:40 06.12. 16:41	Last occurrence 06.12. 16:44 06.12. 16:44	Code 07.08.082 17.99.911	Count 5 2	Logbook Option Release SW Info



6.1.1 Structure of the error list

All system errors (software and hardware) are stored in the error list. The system errors are stored as follows:

First occurrence	Last occurrence	Code	Count
24.01. 16:45	24.01. 16:45	LL.NN.MMM	1
24.01. 16:55	15.05. 07:25	LL.NN.MMM	4

The error messages are stored in chronological order. The most recent error message will appear at the bottom of the error list. If a specific error message is already on the list, only the most recent event will be updated and the error number will be incremented. In this case, no line will be added to the error list.

CAUTION

After each replacement of a processor board the software must be updated by DrägerService. Otherwise, different software versions would impair the function of the Evita.

6.1.2 Error code key

LL	=	Allocation to a specific component, PCB, or software.
NN	=	Type of error.
MMM	=	More detailed description or enumeration.

LL	NN	МММ	Description of message
00	XX	XXX	Normal monitoring. Stored in the user log only.
01	XX	XXX	Settings, etc. Stored in the user log only.
02	01 up to 61	001 up to 999	Errors that are detected by the safety soft- ware.
	71	001	Loudspeaker not detected.
		002	Flow measurement is faulty.
		003 up to 006	GoldCap capacitor.
		007	BOOT test is faulty.
		008	Generation of auxiliary alarm.
		009	The current through the piezo alarm generator is too high.
		010	The current through the piezo alarm generator is too low.
		011	The nebulizer valve monitoring is faulty.
		012	The loudspeaker monitoring circuit has detected an error.
		013	The loudspeaker monitoring circuit has detected an error.
		014	+15 V is too low.
		015	+10 V is too low.
		016	AD conversion for O2 measurement is faulty.
		017	Activation of O2/Air switching valve is faulty.
		018	Wrong nebulizer gas.
		019	Cold start detection is faulty.
		020	Hardware initialization is faulty.
		021	Different quartz times.
	72	000 up to 006	Errors that are detected by the safety soft- ware.
		007 up to 070	Errors that are detected by the ventilation soft- ware.
03	XX	XXX	Control unit components.

Service Diagnosis Mode

11	NN	ммм	Description of message
			Description of message
04	XX	XXX	Electronic assembly components. Components which are not located on the printed circuit boards.
	01	XXX	Power pack.
05	XX	XXX	Pneumatics components.
	01	XXX	Fan
	02	001	Reset-up line
		002	Reset-down line
		003	Disable line
	04	001 up to 004	O2/Air switching valve
06	XX	XXX	Extension box components
07	XX	XXX	Software error (ROSI = operating system soft- ware)
08	XX	XXX	Future expansion
09	XX	XXX	Miscellaneous
10	XX	XXX	Pneumatics Controller PCB
11	XX	XXX	HPSV Controller Air PCB
12	XX	XXX	HPSV Controller O2 PCB
13	XX	XXX	CPU PCB 68332
	98	001	BOOT error
14	XX	XXX	CO2 Carrier PCB
15	XX	XXX	Communication PCB
16	XX	XXX	Pediatric Flow PCB
17	XX	XXX	Graphic Controller PCB
18	XX	XXX	Pneumatics Motherboard PCB
19	XX	XXX	Electronics Motherboard PCB

6.2 Diagnosis page 'Option Release'



Figure 10 Diagnosis page 'Option Release'

Code: Input of release codes for options.

6.3 Diagnosis page 'SW Info'

ront	Electronic	Pneumatic	Other
CPU Pneu Disp	2102 2102 2102	624 624 625	
SmartCare Build ID IP Addres: Subnet ma Gateway	2004 10.55 sk 255.2 10.55	08041740 5.121.247 255.240.0 5.112.1	



Encoded software information about the installed software, e.g. SmartCare.

Schematics and Diagrams

4 x 1.5 Si nat. col. 30 lg 2 x 1 Si rd, or, br, gn, bl (Elastosil E 50 Wacker) 4 x 1 PAE wh, blck/wh Glued with Si glue 4 x 1.5 Si nat. col. 2 x 1.5 Si nat. col. 3b P2a P1a P1a 3c 6alt ്ര 49 e Parallel mixer with connection block P 25 88 10 ±2 I I 88 AIR Detail Z Hoses 29 a I $^{\circ}$ I 220 bl 20 Gas connection block 1 220 200 blck/wh 150 rd AIR AIR R1.3 DR1.2 190 <u>wh</u> 02 I 1.17 ī TOP R1.1 Ģ 1 Y1.3 I Y1.4 þ 50 bl 350 22 OL JO 52 I I ę ł 7a C I 2.1Y I 8a(si I ĝ R1.2 <u>6</u> Expiratory valve gn S 140 Ъ 20 **PEEP/PIP valve** Pressure measuring block Sockets have identical function Plnsp. 110 (since 1997 connections are reversed, hose length to rear R4.1 a panel has been reduced to S6.1 128 mm with hose 7 x 2.5 Si nat. col.) 170 1.8Y N Ν Y6.2 Y4.1 30 S6.2 1 PExp. 250

Figure 1 Piping diagram (Dräger Gas supply block)

Version 1.2_Released_Printed on_09.03.06_S5664590_Schematics_and_diagrams.fm Schutzvermerk DIN 34 beachten. Copyright reserved.

Piping and function diagram



EvitaXL



Version 1.2_Released_Printed on_09.03.06_S5664590_Schematics_and_diagrams.fm Schutzvermerk DIN 34 beachten. Copyright reserved.

EvitaXL

Annex

Test List

Technical Information

Test list (TL)

Dräger medical

A Dräger and Siemens Company

Evita XL

This test list can be processed with standard commercially available test aids and tools, but does not replace the inspections and maintenance work carried out by the manufacturer.

Notes on field of application: Tests marked with the "check" symbol are listed in the test report. The test results are to be documented in the test report. These test instructions apply to software version 5.n.



Table of contents

1	Unit configuration
1. 1	General
1. 2	Serial number (SN)
1. 3	Software version
1. 4	Operating hours
2	Electrical safety
2. 1	Visual check
2. 2	General note
2. 3	Electrical safety to VDE 0751
2. 4	Electrical safety to IEC 60-601/UL 60601-1
3	Function and condition test
3. 1	Accompanying documents
3. 2	Visual check
3. 3	Non-return valve in expiratory valve
3. 4	Power-on test
3. 5	Checking as per internal unit checklist
3. 6	Testing of safety-related valves
3. 7	Temperature measurement (option)
3. 8	CO2 measurement (option)
3. 9	Unit prior to handover
3.10	Unit handover
4	Test equipment
4. 1	Standard commercially available test equipment
5	Tubing diagrams and pneumatic diagram
5. 1	Tubing diagram of pneumatic system with Dräger gas connection block
5. 2	Tubing diagram of pneumatic system with FAS gas connection block
5. 3	Pneumatic diagram

1 Unit configuration

1.1 General

The basic unit configuration covers the following components:

- EvitaXL
- TrolleyVentilation tubing system
- Expiratory valve
- Temperature sensor (option)
- CO2 sensor (option)

Serial number (SN) 1.2

Ø	1.2.1	EvitaXL Serial number (SN) of EvitaXL (Note: The serial number is located on the rating plate	[_txt]
Ø	1.2.2	Expiratory valve 1 Serial number (SN) of expiratory valve 1	[_txt]
Ø	1.2.3	Control panel Serial number (SN) of control panel	[_txt]
	1.3	Software version		
Ø	1.3.1	Unit software version Connect the power plug of the unit to AC power supply.	[_txt]
		Switch the device on. The unit completes its self-test.		
		Set unit to operating mode.		
		Press "System Setup" button. The "System Setup" window is shown on the display.		
		Press "Service" softkey.		
		Read and note down the software version from the display.		
	1.4	Operating hours		
V	1.4.1	Unit operating hours Connect the power plug of the unit to AC power supply. Switch the device on. The unit completes its self-test.	[h]

Set unit to operating mode.

Press "System Setup" button. The "System Setup" window is shown on the display.

Press "Service" softkey.

Read and note down the operating hours from the display.

2 Electrical safety

☑ 2.1 Visual check

- The unit's power supply cord is undamaged.

- The screw-type terminal of the power cable is undamaged.

- The power switch including power-switch protective cover is undamaged. The power-switch protective cover automatically covers the power switch.

[OK]

- The CO2 sensor, the housing and the connecting cable are undamaged.

- The fitted fuse links match the ratings indicated on the power supply unit. [OK]

2.2 General note

NOTICE:

Devices that are connected to a network must be disconnected from this network first before performing electrical safety tests. Otherwise other devices that are connected to the network could be damaged.

The following section details the additional tests of electrical safety to VDE 0751 and IEC 60601-1/UL 60601-1. Which standard is applied depends on the relevant national regulations. If necessary, make the following settings for the unit on the measuring devices: Protection class "SK 1", safety class "BF".

2.3 Electrical safety to VDE 0751

2.3.1 Medical systems

When testing according to VDE 0751, test the system, not the individual devices.

Systems must be handled as devices.

A medical system is a combination of several devices of which at least one is a medical electrical device connected by functional connections or a movable multiple socket-outlet.

Additional ground connections of medical systems such as data cables or similar should be disconnected during testing.

Medical systems must be subjected to a visual check to find out whether individual units of the system have been replaced, added or removed in comparison to the previous configuration used for determination of the initial value. Any changes must be recorded. Measured values must be recorded as initial values after a change in the system.

The sum of the leakage currents of the individual devices must not exceed the specified limit value.

Note

A multiple socket outlet, if applicable, should be included in the test.

Ø	2.3.2	Protective earth conductor test according to VDE 0751 OF The value of the protective earth conductor resistance must not exceed 0.3 ohms in each case.	<]
	2.3.3	Device leakage current test to VDE 0751	
Ø	2.3.3.1	Initial value[μAThe initial measured value must not exceed 500 μA.	¥]
		INFO: Each initial measured value must be transferred to a new test report.	
V	2.3.3.2	Recurrent measurement[μAThe value of the repeat measurement may exceed the initial measured value by max. 50%, but must not exceed 500 μA.	¥]
	2.3.4	Patient leakage current to VDE 0751, ILV port or temperature senso	
V	2.3.4.1	Initial measured value, patient leakage current DC [μ/4 INFO: Each initial measured value must be transferred to a new test report. The initial measured value must not exceed 10 μA.	7]
V	2.3.4.2	Initial measured value, patient leakage current AC [μA INFO: Each initial measured value must be transferred to a new test report. The initial measured value must not exceed 100 μA.	7]
V	2.3.4.3	Recurrent measurement, patient leakage current DCµAThe recurrent measurement must not exceed 10 μA.	¥]
V	2.3.4.4	Recurrent measurement, patient leakage current ACµAThe recurrent measurement must not exceed 100 μA.	¥]
	2.3.5	Trolley multiple socket strip (option) INFO: The electrical safety test is performed within the system test, see note in test step 2.3.1.	

V	2.3.5.1	Fuse links The fuse links match the ratings on the rating plate.	[OK]
	2.4	Electrical safety to IEC 60-601/UL 60601-1		
V	2.4.1	Protective earth conductor resistance WARNING: Test the PE resistance with the power cable connected!	[_Ohm]
	2.4.2	The protective earth conductor resistance must not exceed 0.2 ohms in each case. Earth leakage current		
V	2.4.2.1	Normal condition N.C. The earth leakage current must not exceed 500 μA.	[µA]
V	2.4.2.2	Single fault condition S.F.C. Open circuit in neutral conductor. The earth leakage current must not exceed 1000 μA.	[µA]
		Subsequently the earth leakage current test is repeated with the power plug rotated condition can be established internally on some testers.	d. This	
V	2.4.2.3	Normal condition N.C. The earth leakage current must not exceed 500 μA.	[µA]
V	2.4.2.4	Single fault condition S.F.C. Open circuit in neutral conductor. The earth leakage current must not exceed 1000 μA.	[µA]
	2.4.3	Patient leakage current, ILV port or temperature sensor		
V	2.4.3.1	Normal condition N.C. The patient leakage current must not exceed 100 µA.	[µA]
Ø	2.4.3.2	Single fault condition S.F.C.Open circuit in neutral conductor.The patient leakage current must not exceed 500 μA.Subsequently the patient leakage current test is repeated with the power plug rotat	[ed. This	μA]
		condition can be established internally on some testers.		
Ø	2.4.3.3	Normal condition N.C. The patient leakage current must not exceed 100 μA.	[µA]
---	---------	--	---------------------------------------	---------------------
Ø	2.4.3.4	Single fault condition S.F.C. Open circuit in power cable conductor. The patient leakage current must not exceed 500 μA.	[µA]
	2.4.4	Trolley multiple socket strip (option)		
Ø	2.4.4.1	Fuse links The fuse links match the ratings on the rating plate.	[OK]
V	2.4.4.2	Protective earth conductor resistance The protective earth conductor resistance of the socket outlets must not exceed 0, case.	[2 ohms	_Ohm] in each
Ø	2.4.4.3	Earth leakage current (N.C.) The earth leakage current must not exceed 50 µA.	[µA]
V	2.4.4.4	 Single fault condition (S.F.C.) Open circuit in neutral conductor. The earth leakage current must not exceed 100 μA. Subsequently the earth leakage current test is repeated with the power plug rotate condition can be established internally on some testers. 	[µA]
Ø	2.4.4.5	Earth leakage current (N.C.) The earth leakage current must not exceed 50 μA.	[µA]
Ø	2.4.4.6	Single fault condition (S.F.C.) Open circuit in neutral conductor. The earth leakage current must not exceed 100 μA.	[µA]
	3	Function and condition test		
V	3.1	Accompanying documents - Instructions for Use are available according to user (Note: Comply with national I standards!). - Instructions for Use of options are available according to the user (Note: Comply laws and standards!). - The Medical Products Logbook (applies to Germany only) is available according (Note: Comply with national laws and standards!). [OK]	[aws and with nat to the us	OK] ional ser

	3.2	Visual check			
		Check condition of unit and essential accessories.			
V	3.2.1	Housing The housing is not damaged or dirty. The housing has no corrosion damage. Also paintwork damage. [OK]	[repair	Oł any	<]
V	3.2.2	Control and display elements The operating and display elements are not damaged or dirty. [OK]	[0	〈]
V	3.2.3	AIR and O2 compressed gas ports The AIR and O2 compressed gas ports are undamaged. [OK]	[Oł	<]
Ø	3.2.3	Labelling The labeling, including on adhesive option labels, is complete and legible. [OK]	[Oł	<]
Ø	3.2.5	PCBs The PCBs are all affixed by two screws to the unit. [OK]	[Oł	<]
V	3.2.6	Sockets (sensor, communication) The connectors on the PCBs are undamaged. [OK]	[Oł	〈]
Ø	3.2.7	Power supply unit The power supply unit is not damaged or dirty. [OK]	[Oł	<]
Ø	3.2.8	DC connection for external battery (option) The DC connection is undamaged. [OK]	[Oł	<]
Ø	3.2.9	Remote socket (optional) INFO: If a remote socket is installed on the Evita 4 with the "XL" option, it must be closed There must also be a warning label affixed next to the remote socket.	[off by	Oł a cap.	<]
		The remote socket is closed off by a cap. The cap is undamaged.			
		The warning notice is legible, and not dirty or damaged. [OK]			

☑ 3.2.10 Cooling-air fan on power pack

The cooling air fan is not damaged or dirty. The cooling air fan rotates smoothly. Clean out the cooling air fan with a vacuum cleaner if it is dirty. [OK]

☑ 3.2.11 Cooling-air filter (on rear, right)

The cooling air filter is not damaged, dirty or squashed. Replace the cooling air filter as necessary. [OK]

Remove the front panel (control panel) to perform the following steps.

☑ 3.2.12 Cover

The cover to the right of the inspiratory block is undamaged. [OK]

Remove the screw from the cover and lay the cover aside.

- The room air filter F3.1 in the cover is not damaged, dirty or squashed.

- The fan behind the cover is not damaged or dirty. The fan rotates smoothly and is not dirt clogged. Vacuum the fan to clean it as necessary.

- Unscrew O2 amplifier with non-return valve. The O2 amplifier with non-return valve D3.1 is not damaged or dirty.

- Check expiry date of O2 sensor capsule. Replace the O2 sensor capsule as necessary.

- The lip seal of the O2 sensor mount in the inspiratory block is undamaged.

- The seal between the inspiratory block and the O2 amplifier is undamaged. Insert O2 sensor capsule into inspiratory block. Fit O2 amplifier with non-return valve on unit.

- The inspiratory socket is not dirty or damaged.

- The nebulizer port is identified by a nebulizer icon. The nebulizer port is not dirty or damaged. [OK]

Fit cover on unit.

☑ 3.2.13 Expiratory block

The flow sensor is undamaged.

The flow sensor mount is not dirty or damaged.

The joint seal between the flow sensor and the expiratory valve is undamaged and the connector seal is correctly seated.

[OK]

☑ 3.2.13.1 Expiratory valve without water trap

WARNING:

When fitting, keep to the preferential direction of the plug with copper pipe!

Visual check of individual components - see Fig.. The individual components of the expiratory valve are undamaged.

2005-07 Released



OK]

ſ

OK]

OK]

[OK]

ſ

ſ

[OK]

Assemble the expiratory valve(s) in functional condition.



☑ 3.2.13.2 Expiratory valve with water trap

WARNING:

When fitting, keep to the preferential direction of the plug with copper pipe!

Visual check of individual components as described in previous test step plus additional visual check of individual components – see Figure. The individual components of the expiratory valve and the water trap are undamaged. [OK]

Assemble the expiratory valve(s) in functional condition.



☑ 3.2.13.3 Expiratory valve mount

____OK]

ſ

Γ

OK]

The expiratory valve mount is not damaged or dirty. The expiratory valve locks safely into the mount. [OK]

☑ 3.2.13.4 Patient system heater ventilation grille (without old filter)

____OK]

The patient system heater ventilation grille is not dirty. Clean the ventilation grille, or suction it clean with a vacuum cleaner, as necessary (Note: The ventilation grille is fitted on the right under the expiratory valve. Block the fan during suction.). [OK]

The dust protective filter is not damaged, dirty or squashed. [OK] ☑ 3.2.13.6 Lip seals The lip seals between the expiratory valve and the unit are undamaged. Replace lip seals if necessary. [OK] ☑ 3.2.14 Control panel with attachment - The housing and touchscreen are not dirty or damaged. - The keypad and the control knob are undamaged. The labelling is complete and legible. - The control panel can be swiveled and locked in place. The control panel can be removed without problem. - The rail mounting with interlock is undamaged. - The cable winder is undamaged. - The control panel connecting cable is undamaged. The connecting cable is secured to the control panel and to the front with 2 screws each. - The interlock for mounting of the control panel is undamaged. - The lid of the base unit is undamaged. If no tray or Dräger monitor is fitted to the unit, the four openings in the lid must be closed off with sealing plugs.

Dust protective filter of the patient system heater (new)

- The tray is undamaged.

[OK]

☑ 3.2.13.5

☑ 3.2.15 CO2 sensor (option)

- The housing with anti-kink sleeve is undamaged. The housing is bonded tight.

- The cable with plug is undamaged. The cable insulation is undamaged.

- The window in the CO2 sensor is clean and undamaged.

- The mounting of the test filter on the connecting cable is undamaged. The test filter is clean and not scratched.

- The windows of the cuvettes are bonded in flat, clean and undamaged.

- The windows in the park bracket for the CO2 sensor are clean and undamaged. [OK]

☑ 3.2.16 Two-column trolley (Evita 4 with "XL" option)

The castors are firmly fixed to the trolley and are undamaged. The castors run smoothly and swivel freely. The braking function of castors with parking brakes is safeguarded. [OK]

☑ 3.2.16.1 Accessories

Humidifier bracket The humidifier bracket is undamaged.

Cabinet The cabinet is undamaged.

Multiple socket strip The multiple socket strip is undamaged.

O2 distributor The O2 distributor is undamaged.



[

OK]

OK1

OK]

OK1

OK]

[

[OK]

Ø	3.2.17	Single-column trolley The castors are firmly fixed to the single-column trolley and are undamaged. The c smoothly and swivel freely. The braking function of castors with parking brakes is s [OK]	[castors run safeguardeo	_OK] d.
V	3.2.17.1	Accessories Respiratory gas humidifier bracket (option) The respiratory gas humidifier bracket is undamaged.	[OK]
		Cylinder bracket (option) The cylinder bracket is undamaged.		
		Multiple socket outlet (option) The multiple socket outlet is undamaged.		
		External rechargeable batteries (option) The external rechargeable batteries in the pedestal of the trolley are undamaged.		
		Connecting cable of external battery (option) The connecting cable of the external battery is undamaged. [OK]		
V	3.2.18	Wall rail bracket (option) The wall rail bracket is undamaged. [OK]	[_OK]
V	3.2.19	Hinged arm (option) The hinged arm is moving smoothly and is undamaged. [OK]	[<u>_</u> OK]
V	3.2.20	Temperature sensor (option) The temperature sensor is undamaged. [OK]	[OK]
V	3.2.21	Compressed gas connecting hoses INFO: Comply with national laws and standards for compressed gas connecting hoses!	[OK]
		The compressed gas connecting hoses are undamaged. There are no leaks in the gas connecting hoses. [OK]	compresse	ed

"Battery in cabinet" conversion kit with connecting cable (option) The "Battery in cabinet" conversion kit with connecting cable shows no signs of damage.

[OK]

		The adult test lung comprises the mask tube, ISO size 7 catheter connection socket and 2 L breathing bag. The mask tube, the ISO size 7 catheter connection socket and the 2 L breathing bag are complete and undamaged. (Note: The breathing bag must not be over-inflated or torn.) [OK]
V	3.2.24	SIEMENS test lung (adults) [OK] The SIEMENS test lung is not cracked or porous and is undamaged. [OK]
V	3.2.25	Bellows K (neonatal) [OK] The bellows K are undamaged. [OK]
V	3.2.26	Accessories for neonatal flow (option) [OK] - Flow sensor cable The flow sensor cable is undamaged. The socket on the flow sensor cable locks safely into the flow sensor insert.
		- Flow sensor ISO 15 The ISO 15 flow sensor and the flow sensor insert are undamaged.
		 Water trap for expiratory valve (option) The water trap for the expiratory valve is undamaged. The water trap for the expiratory valve must be used with the optional "neonate flow" feature. [OK]
	3.2.27	Drug nebuliser (white housing) WARNING: If other pneumatic drug nebulisers are used, significant variances in the breathing volume and the inspiratory oxygen concentration may result!
		INFO: Do not connect the drug nebulizer with the "black" housing to the unit!
Ø	3.2.27.1	Components of the drug nebuliser [OK] Check the following components for damage and dirt contamination:
		1 = Patient connection 2 = Drug nebulizer housing, white 3 = Flat seal 4 = O-ring 5 = Nozzle 6 = O-ring 7 = Atomizer 8 = Container

Tubing systems as per Instructions for Use

The tubing systems as per the Instructions for Use are undamaged.

9 = Drug nebulizer tube

☑ 3.2.22

☑ 3.2.23

[OK]

Dräger test lung (adults)

[_____OK]

_OK]

[___]

The listed drug nebuliser components are not damaged or dirty. [OK]



Ø	3.2.27.2	Drug nebulizer function test Fill drug nebulizer with water up to the "3" mark.	[_OK]
		Supply drug nebulizer with 1.8 bar to 2.2 bar.		
		A mist is produced at the outlet of the drug nebuliser. [OK]		
V	3.2.28	Special accessories (option) The special accessories listed below are undamaged:	[_OK]
		- Resutator 2000 - Child Resutator 2000 - Baby resuscitator		
V	3.3	Non-return valve in expiratory valve INFO: Test all expiratory valves listed on the Inspection Report/Log.	[_OK]
		Precondition: The unit is off and the expiratory valves are functional.		
		Connect the non-return valve in the expiratory valve to the switched-off unit.		
		Connect the test set-up below to the expiratory socket.		
		Use the syringe to generate a negative pressure of -7 mbar (cmH2O) and then re mbar (cmH2O).	duce it to -4	
		After approx. 5 seconds, the vacuum for each expiratory valve tested is still at lea (cmH2O).	st -1 mbar	





☑ 3.6.2 Non-return valve D3.2

____mbar]

Precondition: The unit is connected to the compressed gas supply and switched off.

Prepare the following test set-up:

Use the syringe to feed a volume of 50 mL into the inspiratory socket within 2 seconds to 4 seconds.

The measured value on the reference pressure gauge should be 5 mbar (cmH2O) to 10 mbar (cmH2O).

Withdraw the syringe.



☑ 3.6.3 Safety valve D3.3

mbar]

INFO:

The emergency air valve vents automatically when the "Standby" mode is activated (After approx. 3 minutes up to and including SW 6.0. With higher SW versions, immediately)! In order to activate the safety for another 3 minutes, the leak test in the "Check" menu must be started once using the "Check" softkey and then stopped immediately using the "Abort" softkey.

INFO:

The alarm "Pressure meas. inop." may be given in this test; it can be ignored. This alarm must disappear again when the test set-up is removed. There may also be a rise in pressure in the test set-up. The reason for this lies in a leak in the mixer, which has not been zeroed after power-on. The mixer is only zeroed when ventilation starts or after 3 minutes at Standby. The test can be performed immediately, however.

Precondition: The unit is connected to the compressed gas and mains power supply and switched on.

Prepare the following test set-up:

Set EvitaXL to "Standby" mode and confirm. From the "Check" menu, start the leak test once using the "Check" softkey and then stop it immediately using the "Abort" softkey.

Use the syringe to feed a volume of 50 mL into the inspiratory socket within 2 seconds to 4 seconds.

The measured value on the reference pressure gauge should be 101 mbar (cmH2O) to 110 mbar (cmH2O).

Remove the test set-up.



Evita XL

3.7 **Temperature measurement (option)**

Precondition: The unit is connected to the compressed gas and mains power supply and switched on.

☑ 3.8 CO2 measurement (option)

Connect the test gas cylinder and the cuvette to the tube from the calibration set. Slot the CO2 sensor onto the cuvette and enter the test gas cylinder concentration values into the unit. Calibrate the CO2 Sensor.

The displayed concentration corresponds to the value of the test gas.

☑ 3.9 Unit prior to handover

WARNING:

If the unit includes the external rechargeable battery option, check that the red cable connector is connected to the rear of the power supply unit.

Charge the internal and external batteries. Otherwise it will not be possible to run the unit on the fitted batteries in the event of a mains power failure!

☑ 3.10 Unit handover

Place fully functional unit at the user's/owner's disposal.

4 **Test equipment**

4.1 Standard commercially available test equipment

Pressure gauge. -10 mbar to 10 mbar, class 1.6 Pressure gauge, 0 mbar to 160 mbar, class 1.6 Pressure regulator, adjustable up to 2 bar minimum Flowmeter, 0.01 L/min to 14 L/min, class 1.6 Flowmeter, 10 L/min to 120 L/min, class 1.6 Syringe, 50 mL T-piece or Y-piece Thermometer

5 Tubing diagrams and pneumatic diagram

5.1 Tubing diagram of pneumatic system with Dräger gas connection b

[____OK]

OK]

[OK]



5.2 Tubing diagram of pneumatic system with FAS gas connection bloc



5.3 Pneumatic diagram





Test Report TL

Institution:

Delivery date:

Serial no.:

Other:

OK Para Name	Result
1 Unit configuration	
1.1 General	
1. 2 Serial number (SN)	
□ 1. 2. 1 EvitaXL	
□ 1. 2. 2 Expiratory valve 1	
□ 1.2.3 Control panel	
1. 3 Software version	
1. 3. 1 Unit software version	
1. 4 Operating hours	
1.4.1 Unit operating hours	h
2 Electrical safety	
2.1 Visual check	
2. 2 General note	
2. 3 Electrical safety to VDE 0751	
2. 3. 1 Medical systems	
2.3.2 Protective earth conductor test according to VDE 07	
2. 3. 3 Device leakage current test to VDE 0751	
□ 2. 3. 3. 1 Initial value	μA
□ 2.3.3.2 Recurrent measurement	μA
2. 3. 4 Patient leakage current to VDE 0751, ILV port or temperate	ire sensor
2.3.4.1 Initial measured value, patient leakage current DC	μA
2.3.4.2 Initial measured value, patient leakage current AC	μA
2.3.4.3 Recurrent measurement, patient leakage current	μA
2.3.4.4 Recurrent measurement, patient leakage current.	μA
2. 3. 5 Trolley multiple socket strip (option)	
□ 2. 3. 5. 1 Fuse links	
2. 4 Electrical safety to IEC 60-601/UL 60601-1	
□ 2.4.1 Protective earth conductor resistance	Ohm
2. 4. 2 Earth leakage current	
□ 2. 4. 2. 1 Normal condition N.C.	μA
□ 2. 4. 2. 2 Single fault condition S.F.C.	μA
□ 2. 4. 2. 3 Normal condition N.C.	μA
□ 2. 4. 2. 4 Single fault condition S.F.C.	μA
2. 4. 3 Patient leakage current, ILV port or temperature sensor	
□ 2.4.3.1 Normal condition N.C.	μA
□ 2.4.3.2 Single fault condition S.F.C.	μA
2.4.3.3 Normal condition N.C.	μA
□ 2.4.3.4 Single fault condition S.F.C.	μA
2. 4. 4 Trolley multiple socket strip (option)	
□ 2.4.4.1 Fuse links	
□ 2.4.4.2 Protective earth conductor resistance	Ohm
2.4.4.3 Earth leakage current (N.C.)	μA
□ 2.4.4.4 Single fault condition (S.F.C.)	μA
2.4.4.5 Earth leakage current (N.C.)	μA
□ 2. 4. 4. 6 Single fault condition (S.F.C.)	μA
3 Function and condition test	
3.1 Accompanying documents	

ОК	Para Name	Result
3. 2	2 Visual check	
	3. 2. 1 Housing	
	3. 2. 2 Control and display elements	
	3. 2. 3 AIR and O2 compressed gas ports	
	3. 2. 3 Labelling	
	3. 2. 5 PCBs	
	3. 2. 6 Sockets (sensor, communication)	
	3. 2. 7 Power supply unit	
	3. 2. 8 DC connection for external battery (option)	
	3. 2. 9 Remote socket (optional)	
	3. 2.10 Cooling-air fan on power pack	
	3. 2.11 Cooling-air filter (on rear, right)	
	3. 2.12 Cover	
	3. 2.13 Expiratory block	
	3. 2.13. 1 Expiratory valve without water trap	
	3. 2.13. 2 Expiratory valve with water trap	
	3. 2.13. 3 Expiratory valve mount	
	3. 2.13. 4 Patient system heater ventilation grille (without	
	3. 2.13. 5 Dust protective filter of the patient system heat	
	3. 2.13. 6 Lip seals	
	3. 2.14 Control panel with attachment	
	3. 2.15 CO2 sensor (option)	
	3. 2.16 Two-column trolley (Evita 4 with "XL" option)	
	3. 2.16. 1 Accessories	
	3. 2.17 Single-column trolley	
	3. 2.17. 1 Accessories	
	3. 2.18 Wall rail bracket (option)	
	3. 2.19 Hinged arm (option)	
	3. 2.20 Temperature sensor (option)	
	3. 2.21 Compressed gas connecting hoses	
	3. 2.22 Tubing systems as per Instructions for Use	
	3. 2.23 Dräger test lung (adults)	
	3. 2.24 SIEMENS test lung (adults)	
	3. 2.25 Bellows K (neonatal)	
	3. 2.26 Accessories for neonatal flow (option)	
	3. 2.27. 1 Components of the drug nebuliser	
	3. 2.27. 2 Drug nebulizer function test	
	3. 2.28 Special accessories (option)	
	3. 3 Non-return valve in expiratory valve	
	3.4 Power-on test	
	3.5 Checking as per internal unit checklist	
3.6	iesting of safety-related valves	
Ц	3. 6. 1 Emergency air valve Y3.1/non-return valve D3.1	mbar
	3. 6. 2 Non-return valve D3.2	mbar
	3. b. 3 Safety valve D3.3	mbar
3.7	(Iemperature measurement (option)	
IП	3. 8 CO2 measurement (option)	

□ 3.9 Unit prior to handover	4 Test equipment
□ 3.10 Unit handover	5 Tubing diagrams and pneumatic diagram

Report:

Tested according to test specifications.

Name:

Date/signature:

2005-07 Released



2005-12-06

Technical Documentation for EvitaXL according to EMC standard IEC/EN 60601-1-2: 2001

General Information

The EMC conformity includes the use of following external cables, transducers and accessories (see the following table):

Designation	Order no.
Flow sensor for neonates (set of 5)	8410179
Flow sensor connecting cable	8409626
Medibus cable	8306488
ILV cable, Evita 4/2 dura	8411794
CO2 mainstream sensor	8470300 or 6871500
External lead-gel rechargeable batteries 12 V/17 Ah	1843303
Battery cable conversion kit	8411822

EvitaXL should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is inevitable, EvitaXL should be observed to verify normal use in the configuration in which it will be used.

Other equipment which can be used adjacent to or stacked with the EvitaXL are listed in the Instructions for Use manual, in the Order List chapter or in the following table.

Designation	Order no.		
Aquapor EL	8414698, various		
Fisher & Paykel humidifier - accessories	8414144, various		
ECG monitoring	Various, on demand		
Compressor	8414350, various		



A Dräger and Siemens Company

Electromagnetic Emissions

Electromagnetic Emissions							
EvitaXL is intended for use in that is used in such an enviro	EvitaXL is intended for use in the electromagnetic environment specified below. The operator should assure that is used in such an environment.						
Emissions Compliance Electromagnetic environment according to							
RF emissions (CISPR 11)	Group 1	EvitaXL uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.					
	Class A	EvitaXL is suitable for use in all establishments other than domestic and those directly connected to the public low- voltage power supply network that supplies buildings used for domestic purposes.					
Harmonic emissions (IEC 61000-3-2)	Not applicable	Not applicable because RF emissions are class A.					
Voltage fluctuations / flicker (IEC 61000-3-3)	Not applicable	Not applicable because RF emissions are class A.					

Information re electromagnetic emissions (IEC 60101-1-2: 2001, table 201)

Dräger medical

A Dräger and Siemens Company

Electromagnetic Immunity

Electromagnetic Immunity						
EvitaXL is intended for use in the electromagnetic environment specified below. The operator should assure that is used in such an environment.						
Immunity against	IEC 60601-1-2 test level	Compliance level (EvitaXL)	Electromagnetic environment			
electrostatic discharge, ESD (IEC 61000-4-2)	contact discharge: ± 6 kV	± 2, 4, 6 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.			
electrostatic discharge, ESD (IEC 61000-4-2)	air discharge: ± 8 kV	\pm 2, 4, 8 kV, except for the interfaces marked with the \clubsuit ESD symbol.	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.			
electrical fast transients / bursts (IEC 61000-4-4)	power supply lines: ± 2 kV longer input / output lines: ± 1 kV	± 2 kV ± 1 kV	Mains power quality should be that of a typical commercial or hospital environment.			
surges on AC mains lines (IEC 61000-4-5)	common mode: ± 2 kV differential mode: ± 1 kV	± 2 kV ± 1 kV	Mains power quality should be that of a typical commercial or hospital environment.			
power frequency magnetic field 50/60 Hz (IEC 61000-4-8)	3 A/m	3 A/m	In close vicinity to EvitaXL, no equipment with extraordinary power frequency magnetic fields (power transformers, etc.) should be operated.			
voltage dips and short interruptions on AC mains input lines (IEC 61000-4-11)	dip >95%, 0.5 periods dip 60%, 5 periods dip 30%, 25 periods dip >95%, 5 seconds	>95%, 0.5 per. 60%, 5 per. 30%, 25 per. >95%, 5 sec.	Mains power should be that of a typical commercial or hospital environment. If operator requires continued operation during power mains interruptions, it is recommended to power EvitaXL from an uninterruptible supply or a battery.			
radiated RF (IEC 61000-4-3)	80 MHz – 2.5 GHz: 10 V/m	10 V/m	Recommended separation distance from portable and mobile RF transmitters with transmission power P_{EIRP} to EvitaXL including its lines: 1.84 m * $\sqrt{P_{EIRP}}^{X1}$			
RF coupled into lines (IEC 61000-4-6)	150 kHz – 80 MHz: 10 V within ISM bands, 3 V outside ISM bands ^{X2}	10 V 3 V	Recommended separation distance from portable and mobile RF transmitters with transmission power P_{EIRP} to EvitaXL including its lines: 1.84 m * $\sqrt{P_{EIRP}}^{X1}$			

Information re electromagnetic immunity (IEC 60601-1-2: 2001, tables 202, 203, 204)

X1: For P_{EIRP} the highest possible "equivalent isotropic radiated power" of the adjacent RF transmitter has to be inserted (value in Watt). Also in the vicinity of equipment marked with the symbol (w)

(interference may occur. Field strengths from fixed, portable or mobile RF transmitters at the location of EvitaXL should be less than 3 V/m in the frequency range from 150 kHz to 2.5 GHz and less than 1 V/m above 2.5 GHz.

^{X2}:ISM bands in this frequency range are: 6.765 MHz - 6.795 MHz, 13.553 MHz - 13.567 MHz, 26.957 MHz - 27.283 MHz, 40.66 MHz - 40.70 MHz.



A Dräger and Siemens Company

Recommended separation distances

Recommended separation distances between portable and mobile RF telecommunication devices and the EvitaXL			
max. P _{EIRP} (W)	3 V/m distance* (m)	1 V/m distance* (m)	Note
0.001	0.06	0.17	
0.003	0.10	0.30	
0.010	0.18	0.55	
0.030	0.32	0.95	e.g. WLAN 5250 / 5775 (Europe)
0.100	0.58	1.73	e.g. WLAN 2440 (Europe), Bluetooth
0.200	0.82	2.46	e.g. WLAN 5250 (not in Europe)
0.250	0.91	2.75	e.g. DECT devices
1.000	1.83	5.48	e.g. GSM 1800- / GSM 1900- / UMTS- mobiles, WLAN 5600 (not in Europe)
2.000	2.60	7.78	e.g. GSM 900 mobiles
3.000	3.16	9.49	

Information re separation distances (IEC 60601-1-2: 2001, tables 205 and 206)

* 3 V/m distance to transmitters with frequencies from 150 kHz to 2.5 GHz, otherwise 1 V/m distance.

Manufacturer: Dräger Medical AG & Co. KG Moislinger Allee 53 – 55 D-23542 Lübeck Germany

Phone: (+49) (0) 1805-3723437 Fax: (+49) 451/882 - 3779

C E ²²³

Directive 93/42/EEC concerning Medical Devices

Subject to change without notice Will not be replaced in the event of modifications.

© Copyright by Dräger Medical AG & Co. KG, Lübeck, Germany.

The warranty and liability conditions of the general terms and conditions for business transactions of Dräger Medical AG & Co. KG are not extended by this Technical Documentation.