

# Digital Imaging System

## MAXXray

### Service Manual

Vieworks Co., Ltd.

2009/03/17

**Manufacturer; vieworks Co., Ltd.**

#604, SuntechcityII, 307-2, Sangdaewon-dong Jungwon-gu

Seongnam-city Gyeonggi-do, 462-806, South Korea

TEL : +82-70-7011-6161 FAX : +82-31-737-4954

[www.vieworks.com](http://www.vieworks.com)

## Revision History

Rev 0.0    17 Mar. 2009    First release

---

## Accident Reporting

The FDA Medical Device Reporting regulation, 21, CFR 803 and the CE Council Directive 93/42/EEC concerning Medical Devices require that “the manufacturer of medical devices submit a report to the FDA or Local competent authorities whenever he becomes aware of information that reasonably suggests that one of his installed devices:

- may have caused or contributed to a death or serious injury, or
- has malfunctioned and, if the malfunction recurs, is likely to cause or contribute to a death or serious injury.

In order for Vieworks to comply with these requirements, all users of this equipment, operators and service technicians, are required to provide the Quality Assurance Manager at Vieworks with the following information regarding all reportable events as soon as possible:

- 1) Identification of the model and serial number.
- 2) Description of the event. Including whether any serious injury or death has been occurred.
- 3) Identification of the person who is submitting the information including phone number and fax number if available.

## Reference to standards

EN60601-1	Medical electrical equipment Part 1: general requirements for safety
EN60601-2-32	Medical electrical equipment Part 2: Particular requirement for the safety of associated equipment of X-ray equipment
EN60601-1-2	Medical electrical equipment Part 2: electromagnetic compatibility–requirements and tests

## Authorized representatives

If you have any accident, please contact the authorized representatives

## Table of Contents

<b>Revision History</b> .....	2
<b>Accident Reporting</b> .....	3
<b>Reference to standards</b> .....	3
<b>Authorized representatives</b> .....	3
<b>Table of Contents</b> .....	4
<b>1. Safety</b> .....	8
<b>1.1 Safety Guidelines</b> .....	8
<b>1.2 General Hazards</b> .....	9
<b>1.2.1 Radiation hazards</b> .....	9
<b>1.2.2 Electric shock hazard</b> .....	9
<b>1.2.3 Explosion Hazard</b> .....	9
<b>1.2.4 Implosion Hazard</b> .....	9
<b>1.3 Owner’s Responsibility</b> .....	9
<b>1.4 System Diagnostic</b> .....	10
<b>1.5 Calibration</b> .....	10
<b>1.6 Distances measurements</b> .....	10
<b>1.7 Left/Right Marker</b> .....	10
<b>1.8 Images Back-up</b> .....	11
<b>1.9 User Limitations</b> .....	11
<b>1.10 Cleaning the system</b> .....	11
<b>1.11 Overheating</b> .....	11
<b>1.12 Electrical fire</b> .....	11
<b>1.13 EMI/EMC Precaution</b> .....	11
<b>1.15 Maintenance precautions</b> .....	18
<b>1.16 Disposal</b> .....	18
<b>1.17 Changing Fuse</b> .....	18
<b>1.18 Others</b> .....	18
<b>1.19 Appropriation</b> .....	19
<b>1.20 Using together with other equipment</b> .....	19
<b>1.21 Classification (UL)</b> .....	19
<b>1.22 Installation and Maintenance</b> .....	19
<b>2. System Description</b> .....	20
<b>2.1 Intended use</b> .....	20
<b>2.2 System components</b> .....	20
<b>2.3 Component description</b> .....	21
<b>3. System Specification</b> .....	22

---

<b>3.1 CCD</b> .....	22
<b>3.3 Scintillator</b> .....	22
<b>3.4 Time of capture and transmission</b> .....	22
<b>3.5 CCD Cooling System</b> .....	22
<b>3.6 Image Specification</b> .....	22
<b>3.7 Image Tilt</b> .....	22
<b>3.8 Image Distortion</b> .....	22
<b>3.9 Dynamic Range</b> .....	23
<b>3.10 Black Level</b> .....	23
<b>3.11 Resolution</b> .....	23
<b>3.12 Electrical Specification</b> .....	23
<b>3.12.1 AC Power Input(PSU)</b> .....	23
<b>3.12.2 Power Consumption</b> .....	23
<b>3.13 Interface</b> .....	24
<b>4. System Installation</b> .....	25
<b>4.1 Matrox Meter II Digital Image Grabber Installation</b> .....	25
<b>4.1.1 Install Matrox Meter II Digital</b> .....	25
<b>4.1.2 Install Matrox Meter II Digital Driver</b> .....	25
<b>4.2 RXDN-USB2N Installation</b> .....	37
<b>4.2.1 Hardware Installation</b> .....	37
<b>4.2.2 Software Installation</b> .....	37
<b>4.3 RXDN-USB2M Installation</b> .....	40
<b>4.3.1 Hardware Installation</b> .....	40
<b>4.3.2 Software Installation</b> .....	41
<b>4.4 MAXXvue Installation</b> .....	43
<b>4.5 Hardware Installation</b> .....	47
<b>4.5.1 Connect DC Power cable to Detector</b> .....	48
<b>4.5.2 Connect Camera Interface cable to Detector</b> .....	48
<b>4.5.3 Connect RS232C cable to Detector</b> .....	49
<b>4.5.4 Connect DC Power cable to Power Supply Unit</b> .....	49
<b>4.5.5 Connect Data cable to capture board in workstation</b> .....	50
<b>4.5.6 Connect Data cable to RXDN-USB2N</b> .....	50
<b>4.5.7 Connect Data cable to RXDN-USB2M</b> .....	51
<b>4.5.8 Connect RS232C cable to COM port of the workstation</b> .....	51
<b>4.5.9 Connect RS232C cable to COM port of the RXDN-USB2M</b> .....	52
<b>4.5.10 Plug in AC Power cable to power supply unit</b> .....	52
<b>5. Preparation for operating MAXXRAY system</b> .....	53
<b>5.1 Select Image Grabber and MAXXRAY type</b> .....	53
<b>5.2 Diagnosis and calibration</b> .....	54

---

---

5.3 Configuration of MAXXvue.....	54
6. Diagnosis.....	54
6.1 Procedures of Diagnostic.....	54
6.2 Preparation for Diagnostic.....	54
6.3 Password Input.....	54
6.4 Serial Communication Test.....	55
6.5 Set parameter of Detector.....	55
6.6 CCD Cooling Test.....	56
6.7 Black Level Test.....	56
6.8 Flat Field Test.....	57
6.9 Save Result of Diagnosis.....	58
7.1 Procedures of Calibration.....	58
7.2 Preparation for Calibration.....	58
7.3 Acquire Black image.....	59
7.4 Acquire 20% flat field image.....	60
7.5 Acquire 40% flat field image.....	61
7.6 Acquire 60% flat field image.....	62
7.7 Distortion Correction.....	63
7.8 Adjust FOV(Field of view).....	64
7.9 Generate Calibration Data.....	65
7.10 Test Calibration Data.....	66
7.10.1 Resolution test.....	67
8. MAXXvue Configuration.....	68
8.1 What is Configuration.....	68
8.2 General parameter setting.....	68
8.2.1 Hospital information.....	69
8.2.2 USB Image Grabber Selection.....	69
8.2.3 MAXXRAY type.....	69
8.2.4 Detector direction compensation.....	69
8.2.5 Thumbnail tab information.....	69
8.2.6 Expansion setting.....	70
8.2.7 Dicom header info.....	70
8.2.8 Dummy fields display name.....	70
8.2.9 Patient comparison condition.....	70
8.2.10 File-worklist location.....	71
8.3 DICOM parameter setting.....	71
8.3.1 Worklist & PACS.....	72
8.3.2 File-worklist.....	72
9. Trouble Shooting.....	73

---

---

<b>9.1 Failure Mode</b> .....	73
<b>9.2 Repairing Power Failure</b> .....	74
<b>9.3 Repairing RS232 communication Failure</b> .....	74
<b>9.4 CCD Cooling Test Failure</b> .....	74
<b>9.5 Black Level Test Failure</b> .....	74
<b>9.6 Flat Field Test Failure</b> .....	74
<b>9.7 Resolution Test Failure</b> .....	74
<b>10. WARRANTY</b> .....	75
<b>Appendix A MAXXRAY Interconnection Diagram</b> .....	76
<b>Appendix B Symbols</b> .....	77
<b>Appendix C How to use line trigger</b> .....	78

# 1. Safety

## 1.1 Safety Guidelines

 **Caution**

Always be alert when operating this equipment. If a malfunction occurs, do not use this equipment until qualified personnel correct the problem.

This Product was designed and manufactured to ensure maximum safety of operation and to meet all the safety requirements applicable to electronic medical equipment. However, anyone attempting to operate the system must be fully aware of potential safety hazards. It should be operated and maintained in strict compliance with the following safety precautions and operating instruments contained herein:

- 1) The product should be installed, maintained and serviced according to Vieworks maintenance procedures and by Vieworks personnel or other qualified maintenance personnel approved in writing by Vieworks. Operation and maintenance should be done in strict compliance with the operation instructions contained in the maintenance manuals.
- 2) The system, in whole or in part, cannot be modified in any way without written approval from Vieworks.
- 3) Before authorizing any person to operate the system, verify that the person has read and fully understand the Service Manual. The owner should make certain that only properly trained and fully qualified personnel are authorized to operate the equipment. An authorized operators list should be maintained.
- 4) Prevent unauthorized personnel from access to the system.
- 5) It is important that this Service Manual be kept at hand, studied carefully and reviewed periodically by the authorized operators.
- 6) The owner should ensure continuous power supply to the system, with voltage and current according to the product specifications. If power failures are not infrequent, a UPS(Uninterrupted Power Supply) should be installed to avoid loss of data.
- 7) If the product does not operate properly or if it fails to respond to the controls described in this manual, the operator should immediately contact Vieworks field service representative, report the incident and await further instructions.
- 8) The images and calculations provided by this system are intended to be used as tools for the competent user. They are explicitly not to be regarded as a sole incontrovertible basis for clinical diagnosis. Users are encouraged to study the literature and reach their own professional conclusions regarding the clinical utility of the system.
- 9). The user should be aware of the product specifications and of the system's accuracy and stability limitations. These limitations must be considered before making any decision based on quantitative values, in case of doubt, please consult a Vieworks representative.

## 1.2 General Hazards

### 1.2.1 Radiation hazards

This system can be interfaced to x-ray generating equipment. Be certain to follow the safety instructions and specifications for wearing proper lead shielding when in the presence of x-ray generating equipment.

All personnel must wear dosimeters during all phases of installation, operation and maintenance of the system and the equipment to which it is interfaced.

### 1.2.2 Electric shock hazard

A three conductor AC power is supplied with this system to provide the proper electrical grounding. To minimize the shock hazard, the power cable must be plugged into a UL-approved three-contact electrical outlet.

Do not remove or open system covers or plugs. The internal circuits of the system use high voltages that can cause serious injury or death from electrical shock. The operator should never be allowed to open the panels of the system.

### 1.2.3 Explosion Hazard

Do not operate the equipment in the presence of flammable or explosive liquids, vapors or gases. Do not plug in or turn on the system when hazardous substances are detected in the environment.

If flammable substances are detected after the system has been turned on, do not attempt to turn off the system or unplug it. Evacuate and ventilate the area before turning the system off.

### 1.2.4 Implosion Hazard

Do not subject the system to serious mechanical shocks, as the cathode ray tube(CRT) can explode if struck or jarred. This may result in flying pieces of glass and coating that can cause serious injury.

## 1.3 Owner's Responsibility

 **Caution**

Do not use the system if unsafe conditions are known to exist. In case of hardware failure that could cause hazardous conditions (smoke, fire and etc), turn the power OFF and unplug the power cords of all subsystems.

The owner is responsible for ensuring that anyone using the system reads and understand the Service Manual and other relevant literature, and fully understands them. Vieworks makes no representation, however, that the act of reading this manual renders the reader qualified to operate, test and calibrate

the system.

## 1.4 System Diagnostic

The *MAXXvueCalibration* software runs a system diagnostic. Run *MAXXvueCalibration* software when install system or every 1 year after installation.

If an error is detected, report detailed error to Vieworks field service representative.

 **Caution**

The owner is responsible for ensuring that diagnostic of system is performed every year. Do not try to use the system if system diagnostic is fail.

## 1.5 Calibration

To ensure the optimal performances of the system it is important to verify that system is calibrated.

 **Caution**

The owner is responsible for ensuring that the system calibration is performed at installation time or if the system is repaired. Do not try to use the system if system calibration is not performed.

## 1.6 Distances measurements

Distances measurements in millimeters are possible only after distance calibration has been performed using a reference object (see operation manual).

 **Caution**

The operator is responsible for performing distance calibration with a reference object and verifying the results of the distance calibration before taking any distance measurements on an image.

## 1.7 Left/Right Marker

The operator is responsible for the correct and clear marking on the left or right side of the image to eliminate possible errors.

The software includes an option to mark the image with L (left) or R (right) indicator from acquisition phase through printing and archiving. If the operator chose, for any reason, not to use L/R markers, he must use an alternative way to eliminate any possible mistake.

## 1.8 Images Back-up

To avoid the possibility of losing images, which might result in patient being exposed to additional doses of radiation, it is important to back-up the images by filming or by using the CD or DVD option. This should be done as a routine operation for every patient.

If the hard disk of workstation is about to full, the operator should back-up images and delete the images to make room on hard disk for new patient.

### Caution

The operator is responsible for backing-up images of each patient. Do not accumulate images in the system without having a back-up.

## 1.9 User Limitations

The MAXXvue software has service mode, this mode could only be operated with the inputting PASSWORD. The service mode should be operated by the personnel who is qualified by Vieworks.

## 1.10 Cleaning the system

Use only isopropyl alcohol to clean surfaces of the system. Do not use detergents or organic solvents to clean the system. Strong detergent, and organic cleaners may damage the finish and cause structural weakening. Do not clean the system with turning the power on.

## 1.11 Overheating

Do not block the ventilation ports of the detector to prevent overheating of the detector. Overheating can cause system malfunction and damages.

## 1.12 Electrical fire

- This equipment is not suitable for use in the presence of a flammable an aesthetic mixture with air or with oxygen or nitrous oxide.
- Conductive fluids that drain into the active circuit components of the system may cause short circuits that can result in electrical fire. Therefore, do not place fluids or food on any part of the system.
- To avoid electric shocks and burns caused by use of the wrong type of fire extinguisher, make sure that the fire extinguisher at the site has been approved for use on electrical fires.

## 1.13 EMI/EMC Precaution

During installation of the system, care must be taken to prevent the potential risk of electromagnetic interference between this equipment and other devices. The device has been tested for EMI/EMC compliance, but interference can still occur in an electromagnetically noisy environment. Attempt to maintain a suitable distance between electrical devices to prevent cross-interference.

## 1.14 EMC Information

### 1.14.1 Guidance and manufacturer’s declaration - electromagnetic emissions

<p>The EUT is intended for use in the electromagnetic environment specified below.                  The customer or the user of the EUT should assure that it is used in such an environment.</p>		
Immunity test	Compliance	Electromagnetic environment -guidance
RF Emissions CISPR 11	Group 1	The EUT uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF Emissions CISPR 11	Class B	The EUT is suitable for use in ail establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	

### 1.14.2 Guidance and manufacturer’s declaration - electromagnetic immunity

<p>The EUT is intended for use in the electromagnetic environment specified below.                  The customer or the user of the EUT should assure that it is used in such an environment.</p>			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment -guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV Contact ±8kV air	±6kV Contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ± 1kV for input/output lines	±2kV for power supply lines ± 1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.

Surge IEC 61000-4-5	±1kV differential mode ±2kV common mode	±1kV differential mode ±2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycle 70% $U_T$ (30% dip in $U_T$ ) for 25 cycle <5% $U_T$ (<95% dip in $U_T$ ) for 5 s	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycle 70% $U_T$ (30% dip in $U_T$ ) for 25 cycle <5% $U_T$ (<95% dip in $U_T$ ) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EUT image intensifier requires continued operation during power mains interruptions, it is recommended that the EUT image intensifier be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE  $U_T$  is the a.c. mains voltage prior to application of the test level.

### 1.14.3 Guidance and manufacturer’s declaration – electromagnetic immunity

The EUT is intended for use in the electromagnetic environment specified below.

The customer or the user of the EUT should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz	3 Vrms 150 kHz to 80MHz	Portable and mobile RF communications equipment should be used no closer to any part of the EUT, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5GHz	3 V/m 80MHz to 2.5GHz	<p><b>Recommended separation distance</b></p> $d = \left[ \frac{3,5}{V_1} \right] \sqrt{P}$ $d = \left[ \frac{3,5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol :</p> 

NOTE 1) At 80MHz and 800MHz, the higher frequency range applies.

NOTE 2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an

electromagnetic site survey should be considered. If the measured field strength in the location in which the EUT is used exceeds the applicable RF compliance level above, the EUT should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the EUT.

<sup>b</sup> Over the frequency range 150kHz to 80MHz, field strengths should be less than  $[V_1]$  V/m.

### 1.14.4 Recommended separation distances between portable and mobile RF communications equipment and the EUT

The is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the EUT can help Prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the EUT as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150kHz to 80MHz	80MHz to 800MHz	800MHz to 2.5GHz
	$d = \left[\frac{3,5}{V_1}\right]\sqrt{P}$	$d = \left[\frac{3,5}{E_1}\right]\sqrt{P}$	$d = \left[\frac{7}{E_1}\right]\sqrt{P}$
	$V_1=3V_{rms}$	$E_1=3V/m$	$E_1=3V/m$
0.01	0.116	0.1166	0.2333
0.1	0.368	0.3687	0.7378
1	1.166	1.1660	2.3333
10	3.687	3.6872	7.3785
100	11.660	11.6600	23.333

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1) At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

### 1.14.5 Immunity and Compliance Level

Immunity test	IEC 60601 Test Level	Actual Immunity Level	Compliance Level
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80MHz	3Vrms	3Vrms
Radiated RF IEC 61000-4-3	3Vrms 80MHz to 2.5GHz	3V/m	3V/m

### 1.14.6 Guidance and manufacturer’s declaration - electromagnetic immunity

The EUT is intended for use in the electromagnetic environment specified below.

The customer or the user of the EUT should assure that it is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz	3 Vrms 150 kHz to 80MHz	The EUT must be used only in a shielded location with a minimum RF shielding effectiveness and, for each cable that enters the shielded location with a minimum RF shielding effectiveness and, for each cable that enters the shielded location
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5GHz	3 V/m 80MHz to 2.5GHz	Field strengths outside the shielded location from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than 3V/m. <sup>a</sup>  Interference may occur in the vicinity of equipment marked with the following symbol:  

NOTE 1) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

NOTE 2) It is essential that the actual shielding effectiveness and filter attenuation of the shielded location be verified to assure that they meet the minimum specification.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an

electromagnetic site survey should be considered. If the measured field strength outside the shielded location in which the EUT is used exceeds 3V/m, the EUT should be observed to verify normal operation.

If abnormal performance is observed, additional measures may be necessary, such as relocating the EUT or using a shielded location with a higher RF shielding effectiveness and filter attenuation.

## 1.15 Maintenance precautions

Do not open enclosures, disconnect or connect any cables or accessories. Only qualified personnel by Vieworks can do the maintenance.

## 1.16 Disposal

This product contains harmful materials such as lead. Improper disposal of this product may result in environmental contamination.

When disposing of this equipment, contact Vieworks representative. Do not dispose of any part of this equipment without consulting a Vieworks representative first.

Vieworks does not assume any responsibility for damage resulting from disposal of this equipment without consulting Vieworks.

## 1.17 Changing Fuse

### **Caution**

For Continued Protection Against Risk of Fire, Replace Only with Same Type and rating of Fuse. Disconnect Power Before Changing Fuse.

Use only fuse to meet the specification of the system when you replace fuse with another one.

## 1.18 Others

### **Caution**

No User- Serviceable Parts Inside.

## 1.19 Appropriation

 **Caution**

Don't make operation except for the intended purpose

The system, in whole or in part, cannot be modified in any way without written approval from Vieworks.

## 1.20 Using together with other equipment

**Warning:** When the unit is used together with other equipment in the patient area, the equipment shall be connected according to Standard UL 60601-1 and IEC 60601-1.

## 1.21 Classification (UL)

- 1) CLASS I EQUIPMENT
- 2) NO APPLIED
- 3) NO protection against ingress of water
- 4) NOT suitable for use in the presence of a flammable an aesthetic mixture with air or with Qxide
- 5) Continuous operation

## 1.22 Installation and Maintenance

 **Caution**

Only qualified service personnel, who have received training from Vieworks should perform this installation and troubling shooting.

Only qualified service personnel, who have received training from Vieworks should perform this installation and trouble shooting. Calibration procedures should be performed at the system installation time or if the x-ray generator is changed otherwise the system quality is decreased.

## 2. System Description

### 2.1 Intended use

MAXXRAY system is indicated for digital imaging solution designed for general radiographic system for human or veterinary anatomy. It is intended to replace film or screen based radiographic systems in all general purpose diagnostic procedures.

It controls x-ray exposure and x-ray dosage by means of interfacing with x-ray generator.

Various features of this system enable the operator to diagnose easier and faster than conventional non-digital techniques.

Computerized window, image inversion, image processing, zooming, panning, window level adjustment, contrast adjustment, and various features enable the operator to view diagnostic details difficult to see using conventional non-digital techniques.

### 2.2 System components

MAXXray system consists of detector, image transfer device, power supply unit, software and its accessories.

And there are two kind of detectors

#### MAXXRAY system components

Detector	RXDN-6000D RXDN-6500D
Image Transfer Device	Matrox Meter II dig RXDN-USB2N RXDN-USB2M
Power Supply Unit	RXDN-6000P-15
Software	
Viewer	<i>MAXXvue</i>
Configuration	<i>MAXXvueConfigure</i>
Calibration and Diagnostic	<i>MAXXvueCalibration</i>
Accessories	Camera Interface Cable (25M, P/N : 1110-3517-01A) DC Power Cable (10M, P/N : 1200-3407-01A) RS232 Cable (25M, P/N : 1170-3414-01A) Generator Interface Cable (15M, P/N:1170-3417-01A)
Workstation(option)	
OS	Windows XP professional
CPU	Minimum Pentinum 4, 3.0 GHz
Memory	Minimum 2G Byte
Hard Disk	Minimum 80G Byte
Ethernet	Minimum 100 Mbit/s
Monitor	1600x1200, Color
CD Rom	CD or DVD R/W

X-ray grid(option)

12:1 or 13:1, 200lines/inch

## 2.3 Component description

1) Detector

Create X-ray image by using CCD camera and output the result in RS-644 format.

2) Power supply unit

Supply DC-power to detector.

3) Image transfer device

Transfer the result from detector to workstation in USB2.0 specification.

4) MAXXvue

Software to view X-ray image. Get image from detector, process it to ease the diagnostic, save it in database and manage it.

5) MAXXvueCalibration

Diagnostic detector and report the result. Calibrate the system.

6) MAXXvueConfigure

Configure parameters for MAXXRAY Digital Imaging System

## 3. System Specification

### 3.1 CCD

CCD	FFT(Full Frame Transfer) CCD
Effective Pixel Number	2040 × 2478 (RXDN-6000D) 2040 X 2780 (RXDN-6500D)
Cell pitch	9 μm × 9 μm
Fill factor	100%

### 3.2 Area of Image

14" X 17" (RXDN-6000D)

12" X 16" (RXDN-6500D)

### 3.3 Scintillator

Gadolinium (545nm peak) or CsI

### 3.4 Time of capture and transmission

Shorter than 3.5 second

### 3.5 CCD Cooling System

CCD cooling system cools CCD to reduce thermal noise of the image.

1) Cooling Element

Thermoelectric Cooler(TEC)

2) Operation Temperature(CCD)

Maintain  $-5\pm 1^{\circ}\text{C}$  at  $25^{\circ}\text{C}$  of surrounding temperature.

3) Cooling Speed(CCD)

To be kept  $-5\pm 1^{\circ}\text{C}$  within 15 minutes after power on at  $25^{\circ}\text{C}$  of surrounding temperature.

### 3.6 Image Specification

1) Image Format

2040 X 2478 (RXDN-6000D)

2040 X 2780 (RXDN-6500D)

2) Field of View


Captured image should be reflected  $100\pm 2\%$  of scintillator size

### 3.7 Image Tilt

The amount of the rotation of the captured image should not exceed  $2^{\circ}$  comparing with the subject on the screen.

### 3.8 Image Distortion

Distortion(%) =  $100 \cdot (A+A') / (B+C)$  , It should be within 4.5% under this formula



### 3.9 Dynamic Range

#### 1) Definition

In order to measure dynamic range properly, the below procedures is necessary.

- Wait until CCD cooled enough after turning the detector power on.
- Take a dark image using Diagnostic software.

Suppose standard deviation(A) as an average pixel value of the pixels located within diameter of 10 pixels from the center of the image, then dynamic range is defined as followings

$$\text{Dynamic Range} = 20 \log(16383/A) \text{ (dB)}$$

#### 2) Specification

It should be more than 60dB

### 3.10 Black Level

#### 1) Definition

It is the average pixel value of the pixels located within diameter of 20 pixels from the center of the dark image. Dark image should be captured using Diagnostic software after cooling CCD enough.

#### 2) Specification

It should be 2500 +/- 1500(ADU)

### 3.11 Resolution

It should be more than 2.8 lp/mm(RXDN-6000D), 3.3 lp/mm(RXDN-6500D) at the center of the detector screen.

The measure resolution, follow the below procedure.

- Locate resolution chart (Nuclear Associates Model :07-523 or the equivalent) at the center of the detector screen with the diagonal position.
- Line pairs that could be separated by adjusting window level is the resolution.

### 3.12 Electrical Specification

#### 3.12.1 AC Power Input(PSU)

100 - 240VAC +/- 10%, 50/60Hz

#### 3.12.2 Power Consumption

Less than 95W (Detector itself)

### 3.13 Interface

- |                   |   |
|-------------------|---|
| 1) Image data     | RS-644(LVDS)  |
| 2) System Control | RS-232C   |
| 3) Image transfer | USB2.0  |
| 4) RS-232         | Baud rate : 19,200 bps<br>Parity bit : Even<br>Data bit : 8 bit<br>Stop bit : 1 bit |

## 4. System Installation

MAXXray system use Matrox Meter II Digital frame grabber, RXDN-USB2N, RXDN-USB2M as an image grabber

- Refer to 4.1 for Matrox Meter II Digital Image grabber Installation
- Refer to 4.2 for RXDN-USB2N
- Refer to 4.3 for RXDN-USB2M

### 4.1 Matrox Meter II Digital Image Grabber Installation


#### 4.1.1 Install Matrox Meter II Digital

- 1) Install 'Image grabber board(Meteor II Digital-4L)' into one of the PCI slots



#### 4.1.2 Install Matrox Meter II Digital Driver


- 1) Insert the MilLite 7.5 CD into the CD-ROM drive
- 2) Select 'Matrox Imaging Library' and 'Intellicam'. Then Click Continue




#### Caution

If autorun is not enabled in your computer, browse the contents of the MilLite 7.5 and double click the Setup.exe file to run the CD

3) Click "Next" button




4) Click "I Agree" button




5) Click "Next" button



6) Click "Next" button




7) Click "No" button to continue




8) Check "Meteor-II/Digital" and "VGA", then click "Next" button




9) Click "No" button




10) Click "OK" button




11) Click "Yes" button




12) Click "Next" button




13) Click "Next" button




14) Click "Next" button




15) Select "Custom" and input "128" into the edit box and Click "Next" button




16) Click "Next" button




17) Click "Finish" button




18) Click "Next" button




19) Click "Next" button




20) Click "Ok" button



21) Click "Ok" button



22) Click "Ok" button




23) Check " Yes, I want restart my computer now" and click "Finish" button



Then the computer will restart automatically.

24) Run the 'Start->Program->Matrox Imaging Products->Intellicam->Intellicam'. When the following window display, click the 'Close' Button



25) When the following window display, confirm that 'Meteor\_II\_Dig Device 0' is activated (It means that Meteor II Digital has been successfully installed)



Note) If you can not see the window, you can open it running 'Options->System Selection'

Click "Close" button

Now the installation of the frame grabber driver is finished

## 4.2 RXDN-USB2N Installation

※ If you have RXDN-USB2M, skip this section

### 4.2.1 Hardware Installation

- 1) Connect the power-cable and USB-cable to RXDN-USB2N Grabber



- 2) Turn the RXDN-USB2N power on
- 3) Connect the USB-cable to the PC




### 4.2.2 Software Installation

If USB cable is connected with the RXDN-USB2N powered on, the following installation wizard will be executed automatically.


- 1) Select "Install from a list or specific location (Advanced)" then click "Next" button




- 2) Select "Search for the best driver in these location" and check "Include this location in the search" and click "Browse" button.



3) Select folder in which driver file USB2.0 Grabber is located, then click "OK" button.



Then following window will be displayed.  
Click "Next" Button.



4) Click "Continue Anyway" button



5) Click "Finish" button



## 4.3 RXDN-USB2M Installation

※ If you have RXDN-USB2N, skip this section

### 4.3.1 Hardware Installation

- 1) Connect the power-cable and USB-cable to RXDN-USB2M



- 2) Turn on the power switch of RXDN-USB2M
- 3) Connect the other end of the USB-cable to the PC




### 4.3.2 Software Installation

If you connect the USB cable with RXDN-USB2M powered on, the following installation wizard will be executed automatically.


- 1) Select "Install from a list or specific location (Advanced)" then click "Next" button.




- 2) Select "Search for the best driver in these location" and check "Include this location in the search" and click "Browse" button.



3) Select folder in which driver file of RXDN-USB2M is located, then click "OK" button.




If the following window will be displayed, click "Next" Button.



4) Click "Continue Anyway" button



- 5) Click "Finish" button




### 4.4 MAXXvue Installation


- 1) Insert MAXXvue CD in CD Drive
- 2) Run "Setup.exe" program then InstallShield Wizard will be displayed, click "Next" button




- 3) Input User Name and Organization , and select user, then click "Next" button




4) Select setup type as Complete then click "Next" button




5) "Ready to Install the Program" window will be displayed, press "Install" Button




6) "Installing MAXXvue" will be displayed, Now the InstallShield is installing MAXXvue.



7) After finishing installation of MAXXvue, "InstallShield Wizard Completed" window will be displayed, Then press "Finish" button.



8) The following message box will be displayed, press "Yes" button to restart computer.



Now the *MAXXvue* software is successfully installed to the following directory.

Software                   c:\#program files\#*MAXXvue*

Data                        d:\#*MAXXvueData*

Executable File List

MAXXvue.exe	Viewer program
<i>MAXXray</i> Calibration.exe	Diagnosis and Calibration program
MAXXvueConfigure	Configure program


 To run MAXXvue, "East Asian Language" option in windows should be selected, if not , "Selected collating sequence not supported by the operating system" or "the selected order is not supported by this operating system" message will be displayed.

To select this option

1. Start->Control Panel-> Regional and Language...
2. Check "Install files for East Asian languages" option and click "OK" button to install

### 4.5 Hardware Installation


Keep the following step to connect cables between Detector and PC, If the system is used with High Frequency x-ray generator, the interface between x-ray generator and Detector is not needed, but this system is used with single phase x-ray generator x-ray generator should be interfaced with the detector.



### 4.5.1 Connect DC Power cable to Detector



### 4.5.2 Connect Camera Interface cable to Detector



### 4.5.3 Connect RS232C cable to Detector



### 4.5.4 Connect DC Power cable to Power Supply Unit



### 4.5.5 Connect Data cable to capture board in workstation

※ If you are provided with RXDN-USB2N or RXDN-USB2M, skip this section



### 4.5.6 Connect Data cable to RXDN-USB2N

※ If you are provided with Matrox Meteor II dig Grabber or RXDN-USB2M, skip this section



### 4.5.7 Connect Data cable to RXDN-USB2M

※ If you are provided with Matrox Meteor II dig Grabber or RXDN-USB2N, skip this section



### 4.5.8 Connect RS232C cable to COM port of the workstation

※ If you are provided with RXDN-USB2M, skip this section



#### 4.5.9 Connect RS232C cable to COM port of the RXDN-USB2M.

※ If you are provided with Matrox Meteor II dig Grabber or RXDN-USB2N, skip this section



#### 4.5.10 Plug in AC Power cable to power supply unit



## 5. Preparation for operating MAXXRAY system

To operate MAXXRAY system, following step should be executed

- 1) Run MAXXvueConfigure
- 2) Select image grabber type
- 3) Exit MAXXvueConfigure
- 4) Run MAXXvuecalibration
- 5) Diagnosis and calibration
- 6) Exit MAXXvuecalibration
- 7) Run MAXXvueConfigure
- 8) Do Configuration
- 9) Exit MAXXvueconfigure

### 5.1 Select Image Grabber and MAXXRAY type

- 1) Run MAXXvueConfigure in "C:\#program files\#MAXXvue" folder
- 2) Select image grabber

3) Exit MAXXvueConfigure

## 5.2 Diagnosis and calibration

Refer to clause 6 for diagnosis and clause 7 for calibration

## 5.3 Configuration of MAXXvue

Refer to clause 8

## 6. Diagnosis

After the installation of the Detector, check if the detector is operating normally.

This diagnostic procedure will check operating state of detector and will generate diagnostic file for reporting about installed detector.

Diagnostic will be performed in following sequence with the diagnostic software.

### 6.1 Procedures of Diagnostic

- 1) Serial Communication Test
- 2) Set Parameter of Detector
- 3) CCD Cooling Test
- 4) Black Level Test
- 5) Flat Field Test


### 6.2 Preparation for Diagnostic

- 1) Pull out protect cover of Power Supply Unit
- 2) Turn on the power of power supply unit
- 3) Close protect cover of Power Supply Unit
- 4) Run MAXXvueCalibration Software.

### 6.3 Password Input

Input Password, only after inputting password this software will operate.

The default password is "vieworks".



**⚠ Caution**

The password of the following program is same and if you change the password in the one of those program, the password of those program will be changed simultaneously

- MAXXvue.exe
- MAXXvueCalibration.exe
- MAXXvueConfigure.exe

### 6.4 Serial Communication Test

Select a communication port and click "Port Open" button to establish RS232C connection. The selection of a communication port varies depending on the selection of image grabber.


1) If you chose "RXDN-USB2M"



2) If you chose Matrox Meter II Digital or "RXDN-USB2N"




Select an appropriate port and device and click "Serial Communication Test" button. Then the test result will be displayed on the text box



### 6.5 Set parameter of Detector

Select "System Parameter" tab.

Click "Refresh System Control Value" button to get parameter of detector, then the parameters of detector will be displayed on the parameter field




Select "Line Trigger Polarity" and "Target Temperature" for proper one.  
 Set "Trigger Time Expand", this field is used to delay starting of CCD readout after trigger.

## 6.6 CCD Cooling Test


To test CCD cooling

- 1) Select "Cooling monitor" tab.
- 2) Click "Start Temp Monitor" button, then the temperature of CCD and ambient temperature will be displayed every 5 second.
- 3) Check the CCD temperature is under -4°C within 15 minutes
- 4) Click "Stop Temp Monitor" to finish CCD cooling test



## 6.7 Black Level Test

Select "Pre-resolution Test" tab.  
 Click "Get Test Dark Frame" button then dark frame image and the center dark level will be displayed.  
 If "center pixels avg level" is not within 1000 to 4000, then the detector is abnormal




By clicking the "Save Current Image" button, you can save current image.

### 6.8 Flat Field Test

Set SID is 1m, then set x-ray generator to 60KV 5mas

Make an x-ray expose without any target, then the flat field frame image and the center flat level will be displayed.

If the center flat level is not within 500 - 16000, then the detector is abnormal.



By clicking the "Save Current Image" button, you can save current image.

## 6.9 Save Result of Diagnosis

By clicking the "Save Diag Result" button, you can save diagnostic result.

## 7. Calibration

Calibration procedure compensates defect pixels and calibrates pixel gain using the installed x-ray generator and x-ray tube.

The calibration should be performed on the following case

Detector installation

X-ray generator replacement

X-ray tube replacement

### 7.1 Procedures of Calibration

- 1) Acquire Black image
- 2) Acquire 20% flat filed image
- 3) Acquire 40% flat filed image
- 4) Acquire 60% flat filed image
- 5) Distortion correction
- 6) Adjust FOV
- 7) Generate calibration data
- 8) Test calibration data

### 7.2 Preparation for Calibration


- Run Calibration software
- Do step1 to step 4 of the menu
- Set SID to 1.8m or 1m.

 **Caution**

Do not use x-ray grid when do calibration.

## 7.3 Acquire Black image

Select "Dark Frame Acquisition" tab.



Click "Test Shot" button in Dark Data Acquisition area to get black image and make sure black center value is within 1000 to 4500

Input Dark Accumulation Time data, the data means the number of dark frames used in calibration. The recommended data is 10.


Click "Get Dark" button to get dark frame data automatically, then the calibration software will get dark frame data automatically.

During the test of the Calibration, you can control the Window/Level by dragging the mouse. And you can change the level rapidly by dragging and pressing the SHIFT key at the same time.

## 7.4 Acquire 20% flat field image

Set x-ray generator to 60KV, 5mas

Select "20% Saturation Flat Image Acquisition" tab.



Click "Test Shot" button and make an x-ray exposure to get a flat field image and then, you can see center value in the "flat Average value" text box

Make sure the value in "flat Average value" text box should exist within the difference of 10% of the value in "Flat Average Target value" text box, if not, adjust the x-ray dose and repeat getting a flat field image until you get the flat level of the flat field image within the recommended level.


Input Flat Accumulation Time data, the data means the number of flat field images used in calibration. The recommended data is 10.

Wait until Ready of Nth expose message box appear, then make an x-ray exposure to get another flat field image.

Repeat getting the number of flat field image you selected

## 7.5 Acquire 40% flat field image

Select "40% Saturation Flat Image Acquisition" tab.



Adjust x-ray generator dose

Click "Test Shot" button and make an x-ray exposure to get a flat field image and then, you can see center value in the "flat Average value" text box

Make sure the value in "flat Average value" text box should exist within the difference of 10% of the value in "Flat Average Target value" text box, if not, adjust the x-ray dose and repeat getting a flat field image until you get the flat level of the flat field image within the recommended level.

Input Flat Accumulation Time data, the data means the number of flat field images used in calibration. The recommended data is 10.


Click "Get 40% Flat" button to get flat field images then Ready of Nth expose message box will appear and then make an x-ray exposure, then calibration software will get a flat field image and will display Nth expose made.

Wait until Ready of Nth expose message box appear, then make an x-ray exposure to get another flat field image.

Repeat getting the number of flat field image you selected

## 7.6 Acquire 60% flat field image

Select "60% Saturation Flat Image Acquisition" menu.



Adjust x-ray generator dose

Click "Test Shot" button and make an x-ray exposure to get a flat field image and then, you can see center value in the "flat Average value" text box

Make sure the value in "flat Average value" text box should exist within the difference of 10% of the value in "Flat Average Target value" text box, if not, adjust the x-ray dose and repeat getting the flat field image until you get the flat level of the flat field image within the recommended level.

Input Flat Accumulation Time data, the data means the number of flat field images used in calibration. The recommended data is 10.

Click "Get 40% Flat" button to get flat field images then Ready of Nth expose message box will appear and then make an x-ray exposure, then calibration software will get a flat field image and will display Nth expose made.


Wait until Ready of Nth expose message box appear, then make an x-ray exposure to get another flat field image.

Repeat getting the number of flat field image you selected

## 7.7 Distortion Correction

Select "Distortion correction" tab

Put object that is straight on the sides, and expose x-ray




Distortion corrected image will be displayed

If the distortion is not corrected then change "Lens Type" and click "Correct Distortion" button, then distortion corrected image will be displayed.

If disable "Enable Distortion Correction", then distortion will not be corrected.

## 7.8 Adjust FOV(Field of view)

Select "FOV Adjust" tab




Adjust Window/Level to see image clearly.

Adjust right or left line of the red square to locate vertical red line just inside of the white part of the image by dragging the mouse.

Adjust top or bottom line of the red square to locate horizontal red line just inside of the white part of the image by dragging the mouse.

## 7.9 Generate Calibration Data

Select "Gain generation" tab



Select "Flat for gain" and click "Save calibration information" button, then calibration software will generate calibration data file and save.

Flat images selected in "Flat for gain" option is used when generating gain data.

If click "Dark\_Labeled", then detected defect pixel and defect line will be indicated with ►◀ marker in dark image.

If click "flat20\_Labeled", then detected defect pixel and defect line will be indicated with ►◀ marker in Flat20% image

If click "flat40\_Labeled", then detected defect pixel and defect line will be indicated with ►◀ marker in Flat40% image


If click "flat60\_Labeled", then detected defect pixel and defect line will be indicated with ►◀ marker in Flat60% image

## 7.10 Test Calibration Data

Select "Post-resolution test" tab.

Check "Apply dead data" and "Apply gain data".

If you make an x-ray exposure, image will be displayed after being applied calibration data.



"Apply dead data" is an option applying dead pixel compensation.

"Apply gain data" is an option applying CCD pixel gain compensation.

"Apply Distortion correction" is an option applying distortion correction.

After changing options, click "Reprocess" button then you will get image processed with selected options.

By clicking the "Save Current Image" button, you can save current image.


"Open Dark", "Open flat20", "Open flat40", "Open flat60" button is used to open images saved during step 5 to step 8.

"Open information" button is used to check pixel data

Apply image processing option is used to do resolution test

### 7.10.1 Resolution test

1. Set SID to 1m and make sure grid is not installed
2. Attach resolution chart (Model : 07-523) on the center of the front panel of the detector with diagonal direction.
3. Set x-ray focal spot to small size
4. Set x-ray condition to 50KVP and 10mAs(in case of CsI screen) or 20mAs(in case of Gadolinium type screen)
5. Select "Post test" menu
6. Enable "apply dead data" check box  
Enable "apply gain data" check box.  
Enable "apply distortion correction" check box  
Enable " image processing" check box  
Enable "Auto ROI" check box  
Set "Sharpness" to 5  
Set "Contrast" to 2  
Set "Denoise" to 0  
Set "Freq layer" to 2  
Set "Lookup" to 1



7. Shot X-ray
8. Click "Reprocess button"
9. Adjust window width and level to make image looks clear
10. Click "200%" button to zoom in.
11. Check resolution

## 8. MAXXvue Configuration

### 8.1 What is Configuration

Using configuration software we can set the various parameters that is used in the MAXXvue, so the proper setting of parameter is important.

Configuration of MAXXvue is needed when the MAXXRAY is installed, and this job should be performed before using MAXXvue.


Configuration of MAXXvue is consist of two parts

General parameter setting : Basic information used by MAXXvue

Dicom parameter setting : Worklist related DICOM header setting

The configuration software is located at "C:\program files\MAXXvue" folder

### 8.2 General parameter setting



### 8.2.1 Hospital information

- They are the name, address and telephone number of the hospital in which the MAXXRAY system is installed.
- The values of the name field and address field is put into the 'institution name' tag and 'institution address' tag respectively when creating DCM files.

### 8.2.2 USB Image Grabber Selection


Select image grabber

### 8.2.3 MAXXRAY type


Select an appropriate detector type

### 8.2.4 Detector direction compensation


This is used to set the image direction  
Set direction according to installed state of detector




CW : 0  
CCW : 0  
H-Mirror : 1  
V-Mirror : 0



CW : 0  
CCW : 0  
H-Mirror : 0  
V-Mirror : 1



CW : 1  
CCW : 0  
H-Mirror : 0  
V-Mirror : 1



CW : 1  
CCW : 0  
H-Mirror : 1  
V-Mirror : 0

### 8.2.5 Thumbnail tab information



- Select the item(s) to be displayed in the thumbnail tab.
- Patient ID and Patient name is default and the additional information can be displayed in the thumbnail tab according to the setting

## 8.2.6 Expansion setting



If the Minimize Button check box is set, you will have the Minimize button above the 'Exit' button.

## 8.2.7 Dicom header info

### 1) Detector Unique Number

- The 12-digit number unique for each detector which constitutes 'Study Instance UID', 'Series Instance UID' and 'SOP Instance UID'
- If you run MAXXvueCalibration, the field is automatically filled in.

### 2) TAG\_STATION\_NAME

- A string for 'station name' tag
- This is used to distinguish the detectors if two or more detectors are installed in the same hospital

## 8.2.8 Dummy fields display name

- Dummy field is used for the connection with the worklist server
- If the worklist server send some information that is not mentioned in the DICOM standard, then MAXXvue can receive that information using dummy fields
- Dummy field will be displayed as set name in the order list


## 8.2.9 Patient comparison condition

- MAXXvue distinguishes among patients by comparing the selected fields of patients when inserting the patients from the worklist server into the local patient list database
- When we query study order from the worklist server, if all selected patient comparison field data is equal to already registered order, MAXXvue ignore that study order
- If none of comparison field is select, then Patient ID will be used as comparison field
- If any comparison field is selected than only the selected field will be used as comparison field

### 8.2.10 File-worklist location

- If worklist server is not available and the PACS system is providing order using text file MAXXvue can receive patients list from that file
- MAXXvue check assigned folder and if there is new order file MAXXvue will load study order from that file and after loading, MAXXvue will delete that file

## 8.3 DICOM parameter setting



MAXXvue use internally defined tag name related patient information for DICOM, these information will be filled when we register patient or receive study order from the worklist server.

After x-ray exposure these information will be stored to the internal database and will be put to the DICOM file as DICOM tag when we make DICOM file.

When MAXXvue receive order from the worklist server that might not be the standard DICOM tag, so MAXXvue has DICOM tag converting function to make standard DICOM file.

---

### 8.3.1 Worklist & PACS

#### 1) Definition

- Incoming : Receiving tag from the worklist server
- Internal : Internally defined field name of the patient information and study order in MAXXvue
- Outgoing : DICOM tag that will be stored to DICOM file

#### 2) Worklist Matching

- Select internal field name then currently matched tag will be displayed in the items table.
- Select incoming tag from the list and register it by clicking "<<" button
- Select outgoing tag from the list and register it by clicking ">>" button

### 8.3.2 File-worklist

- The format of file worklist should be text file and the field should be separated by delimiter
- The delimiter is user-defined and should be registered in the "delimiter" table
- Internally defined fields for patient information and study order is listed in the "internal field" table
- Set the matching internal field name according to the text order in the worklist file
- Internal field will be translated to the DICOM tag set at 3.1.2
- The delimiter could be one or more characters.
- The date and time format used in worklist file can be selected, and all characters located in the place of where star-closure(\*) exists is ignored.

<Date format>

YYYY : Year represented by full four digits

MM : Month as digits with leading zero for single-digit months

DD : Day of month as digits with leading zero for single-digit days

<Time Format>

HH : Hours with leading zero for single-digit hours; 24-hour clock

MM : Minutes with leading zero for single-digit minutes

SS : Seconds with leading zero for single-digit seconds



## 9.2 Repairing Power Failure

Check AC cable of Power Supply Unit is plugged securely, if not plug securely  
Check power switch in turned on, if not turn on power switch  
If not repaired, then replace Power Supply Unit

## 9.3 Repairing RS232 communication Failure

Check RS232 Cable is installed securely, if not secure screw of cable tightly  
Run calibration program and perform RS232 Communication Test procedure

If the failure is not repaired, Check RS232 port of workstation.  
If RS232 port of workstation id OK, replace RS232 Cable  
If still the failure is not repaired, replace main PCB.  
If still the failure is not repaired, replace Detector.

## 9.4 CCD Cooling Test Failure

Check fan is operating.  
if fan is not operating, check connector of fan is plugged correctly.  
If still the failure is not repaired, Check connector of TEC fan is plugged correctly.  
If still the failure is not repaired, , replace main PCB.  
If still the failure is not repaired, replace Detector.

## 9.5 Black Level Test Failure

Check cable between main PCB and CCD Block, if not plug connector and lock screw tightly  
If not repaired, replace main PCB.  
If still the failure is not repaired, replace Detector.

## 9.6 Flat Field Test Failure

Check x-ray is exposed properly  
If x-ray is exposed properly  
Check Generator Interface Cable is connected properly  
If not repaired, replace main PCB.  
If not repaired, replace Detector.

## 9.7 Resolution Test Failure

If resolution is less than specification, then the detector cannot be repaired in the field, so in this case please contact Vieworks Co., Ltd

---

## 10. WARRANTY

Vieworks Co. warrants that this product will be free from defects in materials and workmanship for a period of twelve (12) months from the date of delivery. If any such product proves defective during this warranty period, Vieworks Co., at its option, either will repair the defective product without charge for parts and labor, or will provide a replacement in exchange for the defective product. In order to obtain service under this warranty, Customer must notify Vieworks Co. of the defect before the expiration of the warranty period and make suitable arrangements for the performance of service. Customer shall be responsible for packaging and shipping the defective product to the service center designated by Vieworks Co. with shipping charges prepaid. Vieworks Co. shall pay for the return of the product to Customer if the shipment is to a location within the country in which the Vieworks Co. designated service center is located. Customer shall be responsible for paying all shipping charges, duties, taxes, and any other charges for products returned to any other locations.


This warranty shall not apply to any defect, failure, or damage caused by improper or inadequate maintenance and care. Vieworks shall not be obligated to furnish service under this warranty to repair damage resulting from attempts by personnel other than Vieworks Co.; or its representatives to install, repair, or service this product, to repair damage resulting from improper use or connection to incompatible equipment or power source; or to service a product that has been modified or integrated with other products when the effect of such modification or integration increases the time or difficulty of servicing the product.

THIS WARRANTY IS GIVEN BY VIEWORKS CO. WITH RESPECT TO THIS PRODUCT IN LIEU OF ANY OTHER WARRANTIES, EXPRESS OR IMPLIED. VIEWORKS CO. AND ITS VENDOR DISCLAIM ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. VIEWORKS CO. RESPONSIBILITY TO REPAIR OR REPLACE DEFECTIVE PRODUCTS IS THE SOLE REMEDY PROVIDED TO THE CUSTOMER FOR BREACH OF THIS WARRANTY. VIEWORKS AND ITS VENDORS WILL NOT BE LIABLE FOR ANY INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES IRRESPECTIVE OF WHETHER VIEWORKS CO. OR THE VENDOR HAS ADVANCE NOTICE OF THE POSSIBILITY OF SUCH DAMAGES.


There are no warranties which extend beyond the description mentioned in this document.

# Appendix A MAXXRAY Interconnection Diagram

## Using RXDN-USB2M Image Grabber



## Using RXDN-USB2N Image Grabber



## Appendix B Symbols



DIRECT CURRENT



ALTERNATING CURRENT



PROTECTIVE EARTH (GROUND)



EQUIPOTENTIALITY



OFF (POWER : DISCONNECTION FROM THE MAINS)



ATTENTION, CONSULT ACCOMPANYING DOCUMENTS



ON (POWER : CONNECTION FROM THE MAINS)



TYPE BF APPLIED PART

## Appendix C How to use line trigger

### 1. Scope


This document explain how to operate MAXXRAY using external trigger.

To take a picture of patient the x-ray will be exposed to the patient and at the same time the CCD in the MAXXRAY should be operated to receive the light from the screen.

For the CCD to receive light and generate an image of patient, MAXXRAY should know the starting time, and the ending time of the x-ray exposure.

The external trigger function is used for the MAXXRAY to know the starting time, and the ending time of the x-ray exposure.

### 2. Interface Circuit



Note : To activate trigger function, the current through LED of TLP281 should be in the range of 5.5mA to 16mA, so the voltage between VCC-EXT and Expose pin should be in the range of 12V to 36V .

### 3. Operation mode of trigger function

#### 1) Active Low Mode

- A. Standby : x-ray is not exposing, so the detector is waiting x-ray exposure
  - Do not flow current through the LED(TLP281).
  - To make above condition, do not apply voltage between *VCC\_EXT* and *Expose*
- B. Exposure : X-ray is exposing

- Flow current through the LED(TLP281) while x-ray is exposing
- To make above condition, apply voltage between *VCC\_EXT* and *Expose* while x-ray is exposing

## 2) Active High Mode


A. Standby : x-ray is not exposing, so the detector is waiting x-ray exposure

- Flow current through the LED(TLP281) while x-ray is exposing
- To make above condition, apply voltage between *VCC\_EXT* and *Expose* while x-ray is exposing

B. Exposure : X-ray is exposing

- Do not flow current through the LED(TLP281).
- To make above condition, do not apply voltage between *VCC\_EXT* and *Expose*

## 4. External trigger timing



<Fig. 5-1>

A. *A* should begin before or concurrently with *B* and *A'* should be terminated after or concurrently with *B'*

## 5. Software setting

Run MAXXvueCalibration and select line trigger polarity.