


combi[®] scan 100



Service Manual



**Analyzer for
Combi Screen[®] 11 SYS
urine test strips**

Contents

1. General information	4
2. Check of operation.	5
3. Check and repair	7
3.1 Check and repair of the power supply	7
3.2 Replacement and repair of the LCD	8
3.3 Replacement and repair of the printer.	9
3.4 Handling of the driving unit defect.	10
3.5 Handling of the optics defect	11
3.6 Repair of the push button panel and cable defect	12
3.7 Handling of SW defect.	12
4. Spare parts.	13
5. Preparation and required instruments	15
6. Error messages, possible reasons and trouble shooting.	16
7. Assembling instructions.	17
7.1 Replacement of the assembled upper housing.	17
7.2 LCD replacement.	20
7.3 Replacement of the push button PCB and the connecting cable.	22
7.4 Replacement of printer and panel	24
7.5 Feeding the printer paper	26

Contact:

Analyticon Biotechnologies AG
Am Mühlenberg 10
35104 Lichtenfels - Germany
Phone: +49(0)6454 7991-0
Fax: +49(0)6454 7991-71
e-mail: info@analyticon.de
www.analyticon.de

1. General information


This service manual describes the control and repair of the Combi Scan 100. It is written for use by a well trained service engineer, who is experienced in repairing and maintaining instruments for diagnostic investigations. The operations, which are described in this manual, should never be carried out by the enduser of the instrument.


The Combi Scan 100 is a precise and costly calibrated optical measurement system. The listed operations below must be accomplished with special attention and precision. Please take special care to disconnect the device from the power supply if this is required by the instruction. If spare parts or accessories are needed, it is absolutely necessary to use only original spare parts for this instrument.

Nevertheless, even a well trained service engineer will not be able to do all operations to get a defective Combi Scan 100 working again, because for several operations special equipment is needed. Depending on the problem, it might be necessary to send the instrument back to Analyticon for repair. More detailed information can be found in this manual. If a return of the instrument to Analyticon is recommended by the manual, please contact Analyticon first to agree on the further steps. Never send the meter back to Analyticon without any written confirmation of Analyticon to do so.

If repair operations are carried out, they are done in the responsibility of the company who does the repair. If the meter is opened for repair, Analyticon will not take any further warranty for this instrument.

2. Check of operation







3. Check and repair

3.1 Check and repair of the power supply


Warning! Never open power supply! If the power supply damaged replace it! Also replace power cord if any damage is observable!




3.2. Replacement and repair of the LCD



3.3. Replacement and repair of the printer



3.4. Handling of the defect driving unit



3.5. Handling of the optics defect

Optics unit defects.

- The LEDs do not light.
- The CCD does not detect.
- Optocoupler does not operate.

CHECK:

3. Insert a strip holder into the equipment and connect the power supply, push „START“ button. Watching it from sideways there must be a white light flashing up according to the rhythm of the movement of the strip holder!


REPAIR:

2. If the defects arisen during the check or the operation of the equipment do not cease

- by the replacement of firmware key (section 11.)
- then by a restart,

the equipment must be returned for repair!

3.6. Repair of the push button panel and cable defect




3.7. Handling of SW defect

If defects, arisen during check or operation, do not cease by the replacement of the firmware key and restart.

The equipment must be returned for repair!

4. Spare parts

Order Article No.	Denomination	Content	Figure
UCP-4451-1	Assembled LCD, with frame and plastic shim. The plastic shim is only at later models.	Display, LCD bear frame Assembly: section 7.2	
UCP-4452-1	Assembled printer, with panel	Printer, panel, cable Assembly: section 7.4	
UCP-4453-1	Assembled upper housing	Upper housing, LCD front panel, printer cover, push buttons Assembly: section 7.1	
UCP-4454-1	Spare assembled start- and push button PCB with cable	Start- and push button panel, cables, plugs Assembly: section 7.3	
1ASA7V1A	Power supply type SA125A-0735U-S		
612EPL19	Printer paper	Feeding: section 7.5	

Order Article No.	Denomination	Content	Figure
352USBAB	USB cable_type_a-b/1.8 m		
35200303	Power cord		
UCP-4404-1	Assembled strip holder for CombiScan100		
UCP-4405-1	Firmware		
UCP-4414-1	Assembled check-strip holder	Assembled strip holder with stuck test strip	
UCP-9201-1	User's manual		
UCP-9202-1	CD		

5. Preparation, required instruments, auxiliary material, labour safety regulations

Cross-grooved screwdriver PH1
Cross-grooved screwdriver PZ1
Straight grooved screwdriver size 4

There are no other labour safety regulations different from the general ones concerning the checking, assembly processes listed in this present Service Manual.


No dangerous material is used.


6. Error messages, possible reasons and trouble shooting methods


Error code	Denomination of the error	Reasons for the error	Trouble shooting
E1	Error in the driving of the strip holder ("0" position detection error).	The optical sensor does not detect the existence of the strip holder. Maybe no strip holder is inserted.	Replacement of the strip holder, if the error is still existing the device should be returned for service.
E3	Firmware (key-chip) error.	There is no firmware inserted or firmware damaged.	Insert firmware of the right version into the equipment to its determined place. If the error is still existing the device should be returned for service.
E4	Communication error (defect of the serial port or cable).	The serial port/plug or cable are defective or the connected unit (PC) doesn't work.	All connections should be fitted. Try to connect to an other device (PC). If the error is still existing the device should be returned for service.
E5	"Worklist" download error	Block, longer than 13 characters, or patient ID block not loaded up to 13 characters	Load the patient ID block up to 13 characters. Check communication with the connecting device. If the error is still existing the device should be returned for service.
E6	"Worklist full". The memory is fully loaded.	Work list memory is fully loaded.	Data down load from the memory, clear memory.
E7	Printer paper run out or is missing.		Feed the equipment with paper. section 7.5.
W1	No data in the memory to be transferred.		Fill up the memory by some measurements and try again. If the error is still existing the device should be returned for service.

7. Assembling instructions




7.1 Replacement of the assembled upper housing (UCP-4453-1)


	Operation	Figure
7.1.1.	<p>Disconnect power supply. Remove firmware (under the printer cover).</p> <p>Unscrew the screw (1 pc PT KA 22x8).</p>	 A close-up photograph of the printer's upper housing. A screw is being removed from the front panel. Two black arrows point from the text in the 'Operation' column to the screw and the front panel area.
7.1.2.	<p>Turn the equipment upside-down.</p>	 A photograph showing the printer's upper housing turned upside-down. The bottom of the housing is visible, showing several screws and a label with the number '120002'.
7.1.3.	<p>Unscrew the screws (4 or 5 pcs PT KA 30x8).</p>	 A photograph showing a person's hands using a red-handled screwdriver to remove a screw from the printer's upper housing. The housing is held in a way that allows for easy access to the screws.

	Operation	Figure
7.1.4.	Turn back the equipment.	
7.1.5.	Lift the upper case (upright), tilt it backward in 180 degree. Note the position of the cable.	
7.1.6.	Unscrew the two screws (PT KA 25x6) which fix the start button PCB, pull out the LED assembled onto the PCB from the housing.	
7.1.7.	Put aside the upper case, prepare the new one.	
7.1.8.	Fit LED into the housing (pay attention to the PCB position), screw the two screw in (2 pcs PT KA 25x6).	



	Operation	Figure
7.1.9.	Turn the upper housing onto the bottom case and holding it at right angles arrange the cables (they must not jam between the two parts). 	
7.1.10.	Turn the equipment upside down.	Figure according to 7.1.2.
7.1.11.	Screw the screws in (4 or 5 pcs PT KA 30x8, 1 pc PT KA 22x8).	Figure according to point 7.1.3.
7.1.12.	Turn the equipment back.	
7.1.13.	Insert the firmware into its place.	Figure according to point 7.1.1.


7.2. Replacement of the LCD (UCP-4451-1)

	<i>Operation</i>	<i>Figure</i>
7.2.1.	Disassemble the equipment.	According to points 7.1.1 – 7.1.5.
7.2.2.	<p>Release connector fixture on both end by pressing in direction as the red arrow shows.</p> <p>Pull out the flexible cable from the connector.</p>	
7.2.3.	<p>Unscrew the 2 screws which fix the button PCB to the frame. Remove plastic shim.</p>	
7.2.4.	<p>Unscrew the screws of the LCD frame (3 pcs PT KA 25x8.)</p>	




	Operation	Figure
7.2.5.	Replace the frame, assembled with LCD, fit it to its position.	
7.2.6.	Screw in the 3 pcs screws which fix the LCD bear frame (PTKA 25x8).	Figure according to point 7.2.4.
7.2.7.	Screw in the 2 pcs screws which fix the push button panel to the LCD bear frame (PTKA 25x8)	Figure according to point 7.2.3.
7.2.8.	Connect flexible cable. First insert cable into the slot of the connector. Then fix the connection by pressing both end of connector in direction shown by the red arrow.	
7.2.9.	Assemble the equipment.	According to 7.1.9 – 7.1.13



7.3 Replacement of the push button PCB and the connecting cable (UCP-4454-1)

	<i>Operation</i>	<i>Figure</i>
7.3.1.	Disassemble the equipment.	According to points 7.1.1 – 7.1.6.
7.3.2.	Put aside the upper housing.	
7.3.3.	Unscrew the 2 pcs screw which fix the push button panel to the LCD frame. Remove plastic shim.	
7.3.4.	Note the position of the cables. Disconnect the cable (10p). (Put aside the panels and the cables.)	
7.3.5.	Connect the connector of the cable of the spare assembled push button panels to the connector of the “main” panel.	Figure according to point 7.3.4
7.3.6.	Screw in the 2 pcs screws which fixes the printer panel to the LCD frame. (2pcs PT KA 25x8)	According to the point 7.3.3


	Operation	Figure
7.3.7.	<p>Insert the LED of the Start button into the housing (take care of the panel position), screw in the two screws fixing the panel (2 pcs PT KA 25x6)</p> 	
7.3.8.	<p>Assemble the equipment.</p>	<p>According to 7.1.9. - 7.1.13.</p>

7.4. Replacement of printer and panel (UCP-4452-1)

	<i>Operation</i>	<i>Figure</i>
7.4.1.	Disassemble the equipment	According to 7.1.1. – 7.1.5.
7.4.2.	Place the upper housing behind the lower housing.	
7.4.3.	Unscrew the 4 screws which fix the printer and the PCB, lift the panel.	
7.4.4.	Disconnect both connectors from the printer PCB. Note the directions and routing of the cables! (2p and 16p)	

	Operation	Figure
7.4.5.	Connect the cables to the spare assembled printer panel. (2p and 16p)	Figure according to point 7.4.4.
7.4.6.	Place the printer panel onto its location. Screw in the two short screws (PTKA 22x6) into the front holes, then screw in the two longer screws (PTKA 22x8) onto the rear. 	
7.4.7.	Assemble the equipment	According to the 7.1.9. – 7.1.13.

7.5. Feeding the printer paper (612EPL19)

	<i>Operation</i>	<i>Figure</i>
7.5.1	Open printer cover, insert paper roll into the holder.	
7.5.2	<p>Insert the end of the printer paper into the slot located at the rear of the printer.</p> <p>After it paper feeding goes on automatically if the device is turned on.</p> <p>The edge of the paper must not be torn or uneven.</p>	
7.5.3	Close the printer cover, stitch the end of the paper through the slot of the cover.	