# A-V Impulse™

## 6060 Controller and Accessories

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
-Mode d'emploi

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-Manual de funcionamiento

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-Manual de funcionamento

## TABLE OF CONTENTS


<b>1. A-V Impulse System</b> .....	<b>EN-1</b>
<b>2. Features of the A-V Impulse 6000 Series Controller</b> .....	<b>EN-1</b>
<b>3. Installation</b> .....	<b>EN-1</b>
<b>4. Accessories</b> .....	<b>EN-2</b>
<b>5. Product Overview</b> .....	<b>EN-2</b>
<b>6. Operating Precautions</b> .....	<b>EN-3</b>
<b>7. Contraindications</b> .....	<b>EN-3</b>
<b>8. Patient and Skin Care</b> .....	<b>EN-4</b>
<b>9. Indications For Use</b> .....	<b>EN-4</b>
<b>10. Fitting the Single Patient Use Impad Cover(s)</b> .....	<b>EN-5</b>
<b>11. Getting Started</b> .....	<b>EN-6</b>
11.1 Sequence of Operation.....	EN-6
11.2 Pausing Treatment.....	EN-6
11.3 Display.....	EN-6
11.4 Night Mode.....	EN-7
<b>12. Preset Parameters</b> .....	<b>EN-7</b>
12.1 Presets and their Uses.....	EN-7
12.2 Customizing Parameters.....	EN-8
<b>13. What to do if the Controller Detects a Fault</b> .....	<b>EN-9</b>
13.1 Fault Codes.....	EN-9
13.2 Hardware Error Codes.....	EN-9
<b>14. Maintenance</b> .....	<b>EN-9</b>
14.1 Routine Maintenance.....	EN-10
14.2 Repair.....	EN-10
14.3 AC Fuse Replacement.....	EN-10
14.4 Cleaning.....	EN-10
14.5 Disposal.....	EN-11
<b>15. Warranty and Factory Service</b> .....	<b>EN-11</b>
<b>16. Specifications</b> .....	<b>EN-12</b>
16.1 Description of Symbols.....	EN-12
16.2 Specifications.....	EN-14
16.3 Impulse Pressure Conversion Table mmHg/KPa.....	EN-14
<b>17. Guidance and Manufacturer's Declaration - Electromagnetic Emissions</b> .....	<b>EN-14</b>

## 1. A-V Impulse System

Clinical studies report that the veins in the sole of the foot act as a very powerful natural blood pump.\* During walking, blood is pumped from these veins in a high velocity surge, up the leg to the heart. The A-V Impulse foot compression system was developed to mimic this natural physiological process, maintaining blood circulation in the legs.

The system consists of a controller connected by tubing assemblies to specially designed inflation pads - Impad™ foot sole cover.

The Impad covers are rapidly inflated by an impulse of air from the controller. Following each impulse the controller automatically allows the inflation pad to deflate.

\*references available upon request

## 2. Features of The A-V Impulse 6000 Series Controller

- Liquid Crystal Display
- Simple set-up
- Microprocessor controller technology
- Alarms and displays to alert user to help with rapid troubleshooting.

## 3. Installation

**Warning: Ensure adequate clearance is provided around the controller to allow for free air circulation. Check that there is reasonable freedom from dust and lint. Do NOT use on wet surfaces. Do NOT cover the controller.**

Check the controller is complete, undamaged and that the voltage rating specified on the rating label is compatible with your AC power supply.

Items supplied with the A-V Impulse foot compression system are:

- A-V Impulse foot compression system controller
- Operators Manual
- Fitted AC Power Cable (as requested)
- Protective Packaging
- Tubing Assembly
- Service Manual


1. Plug the controller power cord into a properly grounded receptacle. The power cord must be situated in a way that makes disconnection easy and convenient. Do not block access to the power cord plug.
2. The A-V Impulse foot compression system is designed to fit on the footboard of a hospital bed using the bed hook located on the rear of the device. If a suitable footboard location is not available then the device can be located on a side table or on the floor at the foot of the bed. Regardless of the placement of the device, it must be situated so that the nurse or caregiver can see and have easy access to the user interface panel.

**Note: If located on the floor take care not to crush the controller when altering the bed height.**


## 4. Accessories

5065	Regular Impad Covers
5075	Large Impad Covers
5107	Extra Large Impad Covers
6067	Sterile Impad Covers, Left
6066	Sterile Impad Covers, Right
5007	Replacement tube set each

## 5. Product Overview

## 6. Operating Precautions

If this controller is dropped, subjected to excessive moisture, fluid spillage, humidity or high temperature, or suspected to have been damaged, remove it from service for inspection by qualified service personnel. When transporting or storing the controller use the original packaging, where possible, and adhere to the transport and storage conditions stated in the specifications section and on the outer packaging.

An electrical shock hazard exists if the casing is opened, removed, or if the unit is in contact with water. This system should not be used on wet surfaces, while the patient is bathing or otherwise in contact with water.

The A-V Impulse 6000 series controller is class 1 Medical Equipment. Warning: To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

An explosion hazard exists if the system is used in the presence of flammable anesthetics. Exercise care to locate the system away from hazardous sources.

When connected to AC power supply a three-wire (Live, Neutral and Earth ground) supply must be used. If the integrity of the external protective conductor is in doubt do not use the controller.

The power cord is the only reliable means to disconnect and isolate from the mains voltage. If required, unplug the device from the electrical receptacle.

A service manual containing servicing and testing information is available for the A-V Impulse foot compression system.

Illustrations used in this manual show typical values that may be used in the function of the controller. These settings and values are for illustrative purposes only. The complete range of settings and values are stated in the specification section.

Refer all servicing to authorized service personnel.

Operator Profile: the operator of this device must be capable of seeing the user interface. The operator should be able to move the device if needed and plug and unplug the air hoses. Although training on the use of the device is recommended, no special skills are required.

## 7. Contraindications

The A-V Impulse foot compression system is contraindicated for patients with conditions where an increase of fluid to the heart may be detrimental, including some patients with congestive heart failure, and those with pre-existing deep vein thrombosis, thrombophlebitis or pulmonary embolism. The device should be used with caution on the infected or insensitive extremity.

Special attention, additional padding and additional checks should be given to patients with poor circulation, fragile skin, insensitive extremities, diabetes and those who may be predisposed to tissue viability problems, including those receiving anti-coagulation therapy. To minimize pressure effects, reduce the impulse pressure and set the impulse duration to 1 second. Check for skin reddening and any early signs that may lead to tissue viability problems. Use additional padding, or discontinue treatment, according to clinical judgement.

## 8. Patient and Skin Care

For optimum results, good priming of the veins is required. This is assisted by a slight degree (approximately 15 degrees) of leg dependency down. Avoid leg elevation unless using an T.E.D.™ anti-embolism stocking, and do not let the limb get cool.

For optimum DVT prophylaxis, it is recommended the A-V Impulse foot compression system be used with T.E.D. anti-embolism stockings on moderate and high risk patients.

As with any treatment technique it is important to check regularly according to hospital protocol, for patient comfort and compliance and to pay particular attention to skin care and hygiene.

### Check Regularly:

- Impulse is felt directly under the arch of the foot.
- Impad cover fits snugly and comfortably.
- Skin integrity (remove stocking/stockinette).
- Skin (color) redness.

## 9. Indications For Use

The proper duration of use for each indication is subject to clinical judgement.

Indication	Recommended Guidelines
<b>Circulation Enhancement</b> Relieves circulatory disorders secondary to diminished blood flow, such as ischaemia secondary to peripheral vascular disease.	For temporary impairments such as temporary trauma or disease conditions, continuous use until the condition is resolved. For chronic impairments, daily use depending on the severity of the patient's condition and activity level.
<b>Deep Vein Thrombosis Prophylaxis and Pulmonary Embolism Prophylaxis</b> Use for effective prophylaxis for patients at risk of deep vein thrombosis (DVT) and pulmonary embolism (PE), including providing pre-intra and post-operative prophylaxis for DVT and PE.	Continuous use until the patient is fully ambulatory and weight bearing (not just mobilized).
<b>Acute Edema</b> Reduces acute edema, such as elevated compartment pressures, edema secondary to trauma and/or surgical procedures, post-bypass graft edema, post operative edema secondary to venous ligation or venous stripping and edema secondary to sprains, strains and sports related injuries of the extremity.	Continuous use until edema is reduced.
<b>Chronic Edema</b> Reduces chronic edema.	As required, but at least 4 hours per day.
<b>Extremity Pain Associated with Incident, Trauma or Surgery</b> Relieves pain, increase range of motion and limb mobility, and expedites return of function following trauma or surgery.	Continuous use until severity of pain is reduced or physician recommends alternative therapy.
<b>Venous Stasis/Venous Insufficiency</b> Treats venous stasis, venous insufficiency and varicose veins.	For temporary impairments such as temporary trauma or disease conditions, continuous use until condition is resolved. For chronic impairments, daily use depending on the severity of the patient's condition and activity level.

Indication	Recommended Guidelines
<b>Lymphedema</b> Reduces lymphedema, including lymphedema secondary to trauma and/or surgery and reduces or controls chronic lymphedema, including post-paralytic lymphedema due to stroke or spinal cord injury.	As required, but at least 4 hours per day.
<b>Leg Ulcers</b> Assists healing of cutaneous ulcers.	Continuous use until ulcer severity is reduced or physician recommends alternative therapy.

## 10. Fitting the Single Patient Use Impad Cover(s)

**Warning:** *The consumables designed for use with the A-V Impulse foot compression system are single patient use. The transfer of consumables between patients may result in cross infection or unexpected consumable failure due to extended use.*

**Caution:** *Check the skin for irritation and use additional padding according to clinical judgement.*

*The Impad cover should be fitted following the instructions supplied with the pad.*



*It is recommended that Impad covers are fitted to the patient in pairs.*

**Caution:** *For full Impad cover life do not walk or weight bear on the Impad cover when the controller is impulsing and only inflate when fitted to the foot. Ensure that the Impad cover is placed directly under the arch of the foot.*

1. For improved efficacy, apply T.E.D. anti-embolism stocking or stockinette over the foot and ankle as required. Avoid wrinkles.
2. Select the appropriate Impad foot cover. Red = LEFT. Blue = RIGHT. Place foot centrally on top of inflation pad as shown by the graphics on the Impad cover. The inlet tube should be on the medial side (inside) of the foot.

**CAUTION:** *The inflation pad must be placed directly under the arch of the foot.*

3. Wrap the medial side (inside) of the foot cover over the top of the foot and then overlap the lateral side (outside) of the cover and secure with the fastener strap. Next, wrap the rear strap around the back of the heel and secure with the fastener tab. Check that the foot cover is fitted snugly and the patient is comfortable.
4. For best results make sure the foot is BELOW the level of the heart.


## 11. Getting Started

**Caution:** Before you start to use The A-V Impulse foot compression system controller read the Operating Precautions, Contraindications and Symbol Definition sections.

**Note:** When the Impad cover(s) are fitted correctly and the controller is switched on the patient will feel an impulse under the arch of their foot.

### 11.1 Sequence of Operation

- Fit the Impad cover(s) to the patient.
- Connect the tubing assembly(s) to the appropriate Air Output Socket(s) on the rear of the controller (A).
- To switch on the controller press the I/O switch to the on (I) position. The display will illuminate and the controller will go through a self-test routine (B).
- The controller will display service and model data.
- Following the self-test a 5 second countdown is displayed in the top left of the display. The controller will display the current patient hours and total use hours.
- Press the / button(s) to start the treatment once the countdown is complete.

**Note:** During the countdown pressing the , , , or will access the buttons secondary function.

It is recommended that the "Patient Hours" function is reset to zero when the controller is used on a new patient.


Patient Hours reset is done by pressing either preset button during the 5 second self-test countdown.

### 11.2 Pausing Treatment

- To pause the treatment, press the / button(s) once.
- To continue the treatment, press the / button(s).

**Note:** If the controller remains paused for 2 minutes it will automatically go into Standby mode. To commence the treatment from Standby mode press the / button(s) twice.

### 11.3 Display

Controller Status

The bottom section of the display shows the controller's status. When functioning correctly, an Impad cover and tubing assembly will represent each selected channel and the smiling controller icon will be displayed.

Shortly before an impulse is due an arrow head flashes at the controller

end of the tubing assembly. The arrow head moves down the tubing assembly flashing one second before it is delivered. The impulse is shown on the Impad icon as it is delivered.


In the event of an alarm condition the frowning controller icon is displayed together with a fault code. For more information on fault codes refer to the label on the hanging handle or Fault Code section of this manual (Section 13.1).

### 11.4 Night Mode

In Night Mode the audible alarm sounds less frequently and the display flashes to draw attention to the controller in dimly lit conditions, minimising patient disturbance.

- To select Night Mode promptly double press .
- To cancel Night Mode promptly double press .

**Note:** the pressure status icon ✓ appears when set pressure is achieved.


## 12. Preset Parameters

The controller is pre-programmed with three commonly used impulse pressure/duration settings. These settings are available through use of the buttons.

### 12.1 Presets and their Uses

#### Preset 1


130 mmHg Impulse Pressure, 3 second Impulse Duration, 20 second Impulse Interval

## Preset 2

130 mmHg Impulse Pressure, 1 second Impulse Duration, 20 second Impulse Interval

**Preset/Action**  
Press (1-2-3) once


Recommended for patients with sensitive limb to maximize comfort

## Preset 3

80 mmHg Impulse Pressure, 1 second Impulse Duration, 20 second Impulse Interval

**Preset/Action**  
Press (1-2-3) twice


Recommended for patients with very sensitive extremities and allows patients to become accustomed to treatment.

## 12.2 Customizing Parameters

*Note: The impulse pressure, duration and interval should be specified by the prescribing physician, according to clinical judgement.*

The following parameters can be adjusted: impulse pressure, impulse duration and interval time. Customizing can be carried out while the controller is operating.

Parameters	Action
Impulse Pressure	Press (+) to increase in 10 mmHg steps. Press (-) to decrease in 10 mmHg steps.
Impulse Duration	Press (1-2-3) to choose 3 seconds or 1 second pulse duration.
Interval Time*	Press (A). Press (+) to increase in increments. Press (-) to decrease in increments. For optimal venous priming a 20 second minimum interval time is recommended.

\*The interval icon is displayed if the default is not used.

## 13. What to do if the Controller Detects a Fault

### 13.1 Fault Codes

The A-V Impulse foot compression system controller is capable of automatic adjustment to correct many situations that may lead to an alarm indication. If the controller is unable to make the necessary adjustment, then an audible alarm will beep and the controller graphic will show where an adjustment is needed.

The fault codes and corrective actions are shown below.

Fault Code	Description	Display	Actions
1	Tubing Not Connected		Check the tubing connections.
2	Tubing Kinked		Check Tubing for kinks or occlusion.
3	Low Pressure		Check Impad cover is not loose, has been disconnected from patient, or the Impad cover has been removed by the patient during use. Refit Impad cover.
4	High Pressure		Check Impad cover is not too tight and the tubing assembly is not kinked or occluded.

*Note: Explanation of the fault codes can also be found on the hanging handle label.*

### 13.2 Hardware Error Codes

When the A-V Impulse foot compression system controller detects a hardware error an audible alarm will beep and the controller graphic will show where an adjustment is needed.

The hardware error codes and corrective actions are shown below.

Fault Code	Description	Display	Actions
E01-E11	Hardware Error Codes	E01-E11	See Service manual
E12	Over Temperature Protection	E12	Check that the controller's air vents are free from obstruction and that there is free air circulation around the controller. Ensure that the environment is not excessively hot. Allowing the controller to cool will reset the alarm. If the problem re-occurs see the Maintenance section.

## 14. Maintenance

**Warning:** Before performing any servicing, maintenance or cleaning procedures, the equipment MUST be switched OFF (O) and disconnected from the AC supply. Servicing should only be performed by authorized service personnel with reference to the Service Manual (AV 6927-02). The 6-pin communication port on the rear of the controller is for use by authorized service personnel ONLY and there must be no patient connection.

**Warning:** The A-V Impulse foot compression system must not be modified in any way. Modification can present a risk of injury to the patient or user.

To ensure that the A-V Impulse foot compression system remains in good operating condition, it is important to keep it clean and carry out the routine maintenance procedures described below.

## 14.1 Routine Maintenance

Interval	Routine Maintenance Procedure
As Required	Thoroughly clean external surfaces of the controller before transfer between patients and after prolonged periods of storage. Attention should be paid to the fan grills, ensure these are clean at all times.
Annually / 3000 hours	<ol style="list-style-type: none"> <li>Inspect the AC power supply plug and cable for damage.</li> <li>Perform the functional tests outlined in the Service Manual.</li> </ol>

## 14.2 Repair

In the event of any problems occurring with this equipment, contact your local representative for repair and the returns procedure.

## 14.3 AC Fuse Replacement

If the controller does not operate when AC power is applied check the fuses on the rear panel and replace if necessary with fuses of the correct value and type; refer to the specification section of this manual (Section 16.2). If the fuses blow again or the controller fails to power up remove it from use and have it examined by authorized service personnel.

## 14.4 Cleaning

**Warning: Disinfectants that are known to be corrosive to metals must not be used.**

Use of Iodine can cause controller surface discoloration.

The controller enclosure can be cleaned with a soft cloth dampened with water or a mild detergent. To sanitize the device, apply cleaning agents with a cloth or wipe. Avoid excessive spraying, especially in the areas of the fan grills and connection ports on the back of the device. If any liquid enters the grill or ports, then internal component damage will likely result. The controller should be wiped with a clean, dry cloth afterward.

Do not immerse in any liquid. Do not use cleaning products containing ammonium chloride, acetone or other aromatic solvents, as those chemicals will degrade the integrity of the case and cause it to become embrittled and possibly crack.

The A-V Impulse foot compression system cannot be effectively sterilized by liquid immersion, autoclaving, or ETO sterilization, as irreparable damage to the system will occur.

The table below provides recommended detergent information and their chemical components.

RECOMMENDED CLEANERS	
Chemical component (with approximate concentrations)	Commercial Example
Dodecylbenzene Sulfonate, Coconut Diethanolamide diluted per instructions	Manu-klenz™**
0.5% bleach solution	Dispatch™**

## 14.5 Disposal

The controller should be disposed of as outlined by local country regulations taking environmental factors into consideration. Do not send the controller back to the manufacturer.

## 15. Warranty and Factory Service

Covidien warrants that your A-V Impulse foot compression system controller does not contain defective material or workmanship. This warranty is for one year from the delivery of the controller to the original purchaser. If within the first year there is a fault on the unit, the controller should be returned to the Covidien Service Department at the address shown below, and Covidien will free of charge, inspect the unit and replace any part which upon their examination, appears to be defective provided that there is no evidence that the unit has been tampered with, or mistreated. This warranty does not apply to the tubing assemblies or the disposable Impad inflation pads or equipment damaged through shipping, tampering, negligence, or misuse, including liquid immersion, autoclaving, or ETO sterilization.

If, within the first year of purchase the controller is serviced by personnel other than those expressly authorized by Covidien, the warranty becomes void, and the company cannot be held liable for any consequential damage caused as a result.

Some countries do not allow the exclusion or limitation of accidental or consequential damages, so the foregoing limitation or exclusion regarding damages may not apply. Also, this Limited Warranty gives you specific legal rights, and you may also have other rights, which vary, from country to country.

The Service Manual is intended as a guide to technically qualified personnel when evaluating instrument malfunctions, it is not to be construed as authorization to perform warranty repairs. Unauthorized service will void the warranty.

If you have any service maintenance problems, contact Customer Services at Covidien.

If possible, the original packaging should be used to ensure safe arrival. Before shipping your unit please call the number below to obtain a return material authorization number.

Covidien maintains a service facility with the capability to promptly repair an A-V Impulse foot compression system controller. Controllers requiring repair must be shipped fully assembled, prepaid and insured to:

**CANADA** Covidien Canada  
7500 Trans Canada Highway  
Pointe-Claire, Qc H9R 5H8  
877-644-8926

**UNITED STATES** Covidien  
5920 Longbow Drive  
Boulder, CO 80301  
1- (800) 255-8522

**AUSTRALIA** Australian Service Centre  
52A Huntingwood Drive  
(Corner Liberty Road)  
Huntingwood NSW 2148  
Toll Free: 1800 350 702  
Tel: (02) 9678 2256  
Fax: (02) 9671 8118

## 16. Specifications

### 16.1 Description of Symbols

<b>secs</b>	Impulse Duration in seconds
<b>mmHg</b>	Impulse Pressure in millimeters of mercury
<b>I</b>	ON - Power connection to the electrical supply
<b>O</b>	OFF - Power disconnection from electrical supply
	Run / Stop - Patient Channel Left and Right
	Increase Impulse Pressure Setting (secondary function - select configuration code. Patient Right Only)
	Decrease Impulse Pressure Setting (secondary function - select compressor pressure. Patient Right Only)
	Interval Adjustment - Display and select interval. (secondary function - select / deselect night mode )
	Preset Setting - Display and select preset parameter settings (secondary function - resets patient hours to zero)
	Dangerous Voltage
	Caution, Consult accompanying documents
	Consult Instructions for use
	Keep away from sunlight
	Keep dry
	Controller status — normal condition
	Controller status — alarm condition
	Channel indicator (not pressurized)
	Channel indicator (pressurized)
	Flow indicator
	Set pressure has been achieved
	Night mode
<b>Rx</b> ONLY	By prescription only

(CONTINUED ON NEXT PAGE)

### 16.1 Description of Symbols (cont'd.)

	WEEE (Waste from electrical and electronic equipment)
	Identification of a substance that is not contained or present within the product or packaging.
	Serial number
	Manufacturer
	Date of manufacture
	Type B applied part
	Fuse

## 16.2 Specifications

All values stated are nominal.

Height: 16.5 cm overall

Depth: 18.6 cm

Width: 32 cm

Weight: 3.57 kg

Electrical Supply: 120V 0.3A 60Hz

Fuses: T1AH 250V x 20 mm (Antisurge)

Standards: EN/IEC 60601-1, EN/IEC 60601-1-2, UL 60601-1,  
CAN/CSA C22.2 No. 601.1-M90, FCC Part 15  
Class 1 Type B device

Output Impulse Pressure: 60 - 200 mmHg

Impulse Duration: 1 second / 3 second (3 second factory set default)

Impulse Interval: 20 - 50 seconds (20 seconds factory set default)

Operating conditions: Temperature 0°C - +45°C

Transportation and

storage conditions: Temperature -40°C to +70°C

Relative humidity: 10%-95% condensing

Atmospheric pressure: 500mbar-1060mbar

## 16.3 Impulse Pressure Conversion Table mmHg/KPa

mmHG	KPa	mmHG	KPa	mmHG	KPa
60	8.0	110	14.7	160	21.3
70	9.3	120	16.0	170	22.7
80	10.7	130	17.3	180	24.0
90	12.0	140	18.7	190	25.3
100	13.3	150	20.0	200	26.7

## 17. Guidance and Manufacturer's Declaration - Electromagnetic Emission

**Warning:** Medical electrical equipment needs special precautions regarding EMC and needs to be installed according to the EMC information provided. Careful consideration of this information is essential when stacking or collocating equipment and when routing cables and accessories.

**Warning:** RF mobile communications equipment can effect medical electrical equipment.

Guidance and manufacturer's declaration - electromagnetic emissions		
The A-V Impulse 6000 Series Controller is intended for use in the electromagnetic environment specified below. The customer or the user of the A-V Impulse 6000 Series Controller should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The A-V Impulse 6000 Series Controller uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The A-V Impulse 6000 Series Controller is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity			
The A-V Impulse 6000 Series Controller is intended for use in the electromagnet environment specified below. The customer or the user of the A-V Impulse 6000 Series Controller should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U <sub>r</sub> (>95 % dip in U <sub>r</sub> ) for 0.5 cycle 40 % U <sub>r</sub> (60 % dip in U <sub>r</sub> ) for 5 cycles 70 % U <sub>r</sub> (30 % dip in U <sub>r</sub> ) for 25 cycles <5 % U <sub>r</sub> (>95 % dip in U <sub>r</sub> ) for 5 sec	<5 % U <sub>T</sub> (>95 % dip in U <sub>T</sub> ) for 0.5 cycle 40 % U <sub>T</sub> (60 % dip in U <sub>T</sub> ) for 5 cycles 70 % U <sub>T</sub> (30 % dip in U <sub>T</sub> ) for 25 cycles <5 % U <sub>T</sub> (>95 % dip in U <sub>T</sub> ) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the A-V Impulse 6000 Series Controller requires continued operation during power mains interruptions, it is recommended that the A-V Impulse 6000 Series Controller be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<b>NOTE:</b> U <sub>r</sub> is the a.c. mains voltage prior to application of the test level.			

## Guidance and manufacturer's declaration - electromagnetic emissions

The A-V Impulse 6000 Series Controller is intended for use in the electromagnetic environment specified below. The customer or the user of the A-V Impulse 6000 Series Controller should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance Level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the A-V Impulse 6000 Series Controller, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = [1.17] \sqrt{P}$ $d = [1.17] \sqrt{P}$ 80 MHz to 800 MHz $d = [2.33] \sqrt{P}$ 800 MHz to 2.5 GHz where <b>P</b> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <b>d</b> is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

**NOTE 1:** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup>Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the A-V Impulse 6000 Series Controller is used exceeds the applicable RF compliance level above, the A-V Impulse 6000 Series Controller should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the A-V Impulse 6000 Series Controller.

<sup>b</sup>Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

### Recommended separation distance between portable and mobile RF communications equipment and the A-V Impulse 6000 Series Controller @ 3Vrms

The A-V Impulse 6000 Series Controller is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the A-V Impulse 6000 Series Controller can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the A-V Impulse 6000 Series Controller as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[ \frac{3.5}{v_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{3.5}{E_T} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[ \frac{7}{E_T} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where **P** is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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**COVIDIEN**

IFU No. VVO100051

Covered by U.S. Patent 7,076,993 when used with an inflatable bag.

COVIDIEN, COVIDIEN with logo, and Covidien logo are U.S. and internationally registered trademarks of Covidien AG.

™\*Trademark of its respective owner.

Other brands are trademarks of a Covidien company.

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
### ADDITIONAL JOB INFORMATION (if available)

<b>Job #:</b>	VT620	<b>Dims:</b>	5.50" W x 8.50" H
<b>Artwork/Part #:</b>	VVO100051	<b>Metric Dims:</b>	13.97 cm W x 21.59 cm H
<b>Print Process:</b>	Offset Lithography	<b>Package Type:</b>	OPML (Instructions)
<b>Barcode Format:</b>	N/A	<b>File Name:</b>	VVO100051_GNRC_OPML
<b>Dieline:</b>	N/A	<b>Temp ID:</b>	N/A

<b>Check here if:</b>	<input type="checkbox"/> US Only <input checked="" type="checkbox"/> Generic Artwork	<b>Linked to Generic Artwork:</b>	AV6060-USA
		(if applicable)	

### NUMBER OF COLORS

### TINTS

### TECHNICAL INFORMATION

<b>Application/Version:</b>	InDesign/CS5
<b>Placed Images:</b>	N/A
<b>Fonts Used:</b>	Myriad Pro

### Comments:

---

---

---

### Printer Notes:

---

---

**COVIDIEN**

15 Hampshire Street,  
Mansfield, MA 02048

### GBU-SPECIFIC INFORMATION

<b>Graphic Designer:</b>	Noel Gouveia

### DESIGN AGENCY INFORMATION (for external use only)

<b>Agency Name:</b>	N/A
<b>Agency Contact:</b>	N/A

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©2011 Covidien.
**ECIREP** Covidien Ireland Limited,