Compex® runner

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













Art no. 4518164 v4

Compex® is a registered trademark of DJO, LLC.

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○ **Mi-SENSOR** location only for stimulators with muscle intelligence™ technology.

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Approximate placing of the electrodes makes for a less efficacious session, but is still risk-free. Compex declines all responsibility for the placing of electrodes on areas other than those recommended or in case of failure to respect the counter-indications.

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página 104

INTRODUCTION

Please read this manually carefully prior to using your Compex stimulator. You are strongly advised to read chapter I "Warnings" of this manual. The Compex Runner is a stimulator intended for muscular training and pain relief. Any person excepted those mentioned in the chapter I "Warnings" can use the Compex Runner.

CONTENTS

I. WARNINGS

1. Counter-indications	1
2. Safety measures	2

II. PRESENTATION

1. Reception of equipment and accessories	7
2. Guarantee	7
3. Maintenance	7
4. Storage and transport conditions	8
5. Conditions of use	8
6. Elimination	8
7. Standards	8
8. Patents	8
9. Normalised symbols	8
10. Technical characteristics	9

III. HOW DOES ELECTROSTIMULATION WORK? 10

IV. USAGE GUIDELINES

1. Electrode positions	12
2. Stimulation positions	12
3. Adjusting stimulation energies	12
4. Progression in the levels	13
5. Alternation of stimulation sessions/voluntary training	13

V. THE mi TECHNOLOGY

1. Practical rules of use	13
---------------------------	----

VI. HOW IT WORKS

1. Description of the stimulator	15
2. Connections	16
3. Preliminary settings	16
4. Selecting a programme category	17
5. Selecting a programme	17
6. Personalising a programme	18
7. During the stimulation session	18
8. Consumption and recharging	21
9. Problems and solutions	22

VII. PROGRAMMES 24

VIII. EMC TABLE (Electromagnetic Compatibility) 29

I WARNINGS

1. COUNTER-INDICATIONS

Important counter-indications

- Cardiac stimulator (pacemaker)
- Epilepsy
- Pregnancy (do not use on abdominal region)
- Serious arterial circulation problems in lower limbs
- Abdominal or inguinal hernia

Precautions when using the Compex

- After trauma or a recent operation (less than 6 months)
- Muscular atrophy
- Persistent pains
- Need for muscular rehabilitation

Osteosynthesis equipment

The presence of osteosynthesis equipment (metallic equipment in contact with the bone: pins, screws, plates, prostheses, etc.) is not a counterindication for the use of Compex programmes. The electrical currents of the Compex are specially designed to have no harmful effect on osteosynthesis equipment.

Important!

- Do not use Compex stimulator programmes if you have sensitivity problems.
- Never use the Compex for prolonged periods without medical advice.
- Consult your doctor if you are in any doubt whatsoever.
- Read this manual carefully, in particular chapter VII, which provides information concerning the effects and indications of each stimulation programme.

2. SAFETY MEASURES

What you should not do with the Compex and the *mi* system

- Do not use the Compex or the *mi*-SENSOR system in water or in a humid atmosphere (sauna, hydrotherapy, etc.).
- Do not use the Compex or *mi*-SENSOR system in oxygen-rich environments.
- Never carry out an initial stimulation session on a person who is standing. The first five minutes of stimulation must always be performed on a person who is sitting or lying down. In rare instances, people of a nervous disposition may experience a vasovagal reaction. This is of psychological origin and is connected with a fear of the muscle stimulation as well as surprise at seeing one of their muscles contract

without having intentionally contracted it themselves. A vasovagal reaction causes heart to slow down and blood pressure to drop, which can make you feel weak and faint. If this does occur, all that is required is to stop the stimulation and lie down with the legs raised until the feeling of weakness disappears (5 to 10 minutes).

- Never allow movement resulting from muscular contraction during a stimulation session. You should always stimulate isometrically; this means that the extremities of the limb in which a muscle is being stimulated must be firmly fixed, so as to prevent the movement that results from contraction.

- Do not use the Compex or **mi**-SENSOR system if you are connected to a high-frequency surgical instrument as this could cause skin irritation or burns under the electrodes.
- Do not use the Compex or **mi**-SENSOR system within X metre (see EMC table) of short wave or microwave devices as this could alter the currents generated by the stimulator. If you are in any doubt as to the use of the Compex close to another medical device, seek advice from the manufacturer of the latter or from your doctor.
- Do not use the Compex stimulator or **mi**-SENSOR system in areas where unprotected devices are used to emit electromagnetic radiation. Portable communications equipment can interfere with electrical medical equipment.
- Only use stimulation cables supplied by Compex.
- Do not disconnect any stimulation cables during a session while the stimulator is switched on. Switch the stimulator off first.
- Never connect stimulation cables to an external power supply. There is a risk of electric shock.
- Only use rechargeable battery units supplied by Compex.
- Never recharge the stimulator without first disconnecting the stimulation cables.
- Always use the charger provided by Compex to recharge the batteries.
- Never use the Compex or the charger if it is damaged (case, cables, etc.) or if the battery compartment is open. There is a risk of electric shock.
- Disconnect the charger immediately if the Compex "beeps" continuously, if there is abnormal heating or smell, or if smoke comes from the charger or the Compex.
- Do not recharge the battery in a confined space (carrying case, etc.). There is a risk of fire or electric shock.
- Keep the Compex and its accessories out of reach of children.
- Do not allow any foreign bodies (soil, water, metal, etc.) to penetrate the Compex, the battery compartment and the charger.
- Sudden temperature changes can cause condensation to build up inside the stimulator. To prevent this, allow it to reach ambient temperature before use.
- Do not use the Compex while driving or operating machinery.
- Do not use the stimulator at altitudes of over 3,000 metres.

Where never to apply the electrodes

- In the vicinity of the head.
- Counter-laterally, i.e. do not use two poles connected to the same channel on opposite sides of the body.
- On or near skin lesions of any kind (wounds, swelling, burns, irritation, eczema, etc.).

Precautions when using the **mi**-SENSOR system

- To access the Compex's **mi** technology functions, it is essential to connect the **mi**-SENSOR cable before switching on the device.
- Avoid connecting the **mi**-SENSOR cable when the Compex is switched on.
- Do not disconnect the **mi**-SENSOR cable while the sensor is being used.
- To work correctly, the **mi**-SENSOR system must not be blocked or subjected to pressure.

Precautions when using the electrodes

- Only use electrodes supplied by Compex. Other electrodes may have electrical properties that are unsuitable for the Compex stimulator.
- Always turn off the stimulator before moving or removing any electrodes during a session.
- Do not place the electrodes in water.
- Do not apply solvents of any kind to the electrodes.
- For best results, wash and clean the skin of any oil and dry it before attaching the electrodes.
- Attach the electrodes in such a way that their entire surface is in contact with the skin.
- For obvious reasons of hygiene, each user must have his/her own electrode set. Do not use the same electrodes on different people.
- Never use a set of adhesive electrodes for more than 15 sessions as their bonding power deteriorates over time and optimal contact is very important for both user comfort and efficacy.
- Some people with very sensitive skin may experience redness under the electrodes after a session. Generally, this redness is completely harmless and disappears after 10 to 20 minutes. Never start another stimulation session in the same area, however, if the redness is still visible.

II PRESENTATION

1. Reception of equipment and accessories

Your kit contains:

A 1 stimulator

B 1 charger

C 1 set of black Snap connectors for electrodes, colour-coded (blue, green, yellow, red)

D 2 bags of small electrodes (5x5 cm)

E 2 bags of large electrodes (5x10 cm)

F 2 user manuals

G 1 belt clip

H 1 travel pouch

2. Guarantee

Refer to the flyer attached.

3. Maintenance

Do not sterilise the stimulator.

Clean your stimulator using a soft cloth, and an alcohol-based, solvent-free cleaning product.

Use only a minimum amount of liquid when cleaning the Compex.

Do not attempt to repair the stimulator or any of its accessories.

Never dismantle the Compex or the charger containing high-voltage parts because of risk of electric shock.

Compex Médical SA declines all responsibility for any damages or consequences resulting from unauthorised attempts to open, modify, or repair the

stimulator. This may only be done by persons or repair services authorised by Compex Médical SA.

Your stimulator does not require calibration. Each Compex stimulator is always tested and validated prior to distribution. These are stable and do not vary under normal conditions.

Nonetheless, as the Compex is a high-quality electrical instrument, its lifespan depends on the use that is made of it and the care and maintenance it receives during its lifetime. If your stimulator contains parts that seem worn or defective, please contact your nearest Compex Médical SA customer service centre regarding an upgrade.

Medical and health professionals are obliged to service the unit in compliance with relevant national laws and regulations. This involves testing performance and safety parameters at regular intervals.

4. Storage and transportation conditions

The Compex contains rechargeable batteries. This is why it must be stored and transported in accordance with the following conditions:

Storage and transportation temperature:
- 20°C to 45°C

Maximum relative humidity: 75%

Atmospheric pressure: 700 hPa to 1,060 hPa

5. Conditions of use

Usage temperature: 0°C to 40°C

Relative humidity: 30% to 75%

Atmospheric pressure: 700 hPa to 1,060 hPa

Do not use in areas where there is a risk of explosion.

6. Elimination

The main purpose of the 2002/96/EEC Directive is to prevent the creation of waste electrical and electronic equipment (WEEE) and to reduce the amount of waste to be disposed of by encouraging reuse, recycling, and other forms of recovery. The wheelie bin pictogram with the bar through it means that the equipment cannot be discarded with household refuse, but that it must be collected selectively. The equipment must be delivered to a suitable collection point for treatment. By doing so, you will be contributing to the safeguarding of natural resources and health. Batteries must be disposed of in compliance with relevant national regulatory requirements.

7. Standards

The Compex complies with current medical standards.

To guarantee your safety, the Compex has been designed, manufactured, and distributed in compliance with the requirements of European Directive 93/42/EC on medical devices.

The Compex also complies with the IEC 60601-1 standard on general safety requirements for electro-medical devices, the IEC 60601-1-2 standard on electromagnetic compatibility, and the IEC 60601-2-10 standard on particular safety requirements for nerve and muscle stimulators.

Current international standards require that a warning be given concerning the application of electrodes to the thorax (increased risk of cardiac fibrillation).

The Compex also complies with Directive 2002/96/EEC on waste electrical and electronic equipment (WEEE).

8. Patents

The Compex incorporates several innovations with patents pending.

9. Normalised symbols

Caution: Read the user manual or operating instructions (symbol no. 0434 IEC 60878).

The Compex is a class II device with internal electric power and type BF applied parts (symbol no. 5333 IEC 60878).

Waste electrical and electronic equipment (WEEE) marking according to EN 50419 standard.

The On/Off button is a multifunction button (symbol no. 5009 IEC 60878).

Functions

On/Off (two stable positions)

Waiting or preparation for part of the unit

Stop (turns system off)

10. Technical characteristics

General information

941210 Battery unit

Nickel metal hybrid (NiMH) rechargeable (4.8V / $\geq 1,200$ mA/h).

68301x Chargers

Only battery chargers bearing the following information can be used to recharge the batteries supplied with Compex stimulators:

Europe

683010

Type TR1509-06-E-133A03

Input 100-240 VAC / 47-63 Hz / 0.5 A max.

Output 9V / 1.4 A / 15 W

UK

683012

Type TR1509-06-U-133A03

Input 100-240 VAC / 47-63 Hz / 0.5 A max.

Output 9V / 1.4 A / 15 W

601131 Black Snap-connection stimulation cables

Device connector: 6-pin

Electrode connector: Female Snap

Length: 1,500 mm

601160 Stimulation cable fitted with *Mi*-SENSOR

(accessory sold separately)

Device connector: 6-pin

Electrode connector: Female Snap

Length: 1,500 mm

Protection index

IPXO (IEC 60529)

Neuro-stimulation

All electrical specifications are given for an impedance of 500-1,000 ohms per channel.

Channels: Four independent and individually adjustable channels that are electrically isolated from each other and earthed.

Pulse shape: constant rectangular current with pulse compensation to eliminate any direct current component to prevent residual polarisation at skin level.

Maximum pulse intensity: 120 mA.

Pulse intensity increments: manual adjustment of stimulation intensity from 0 to 999 (energy) in minimum increments of 0.5 mA.

Pulse width: 60 to 400 μ s.

Maximum electrical charge per pulse: 96 microcoulombs (2 x 48 μ C, compensated).

Standard pulse ramp-up time: 3 μ s (20% 80% of maximum current).

Pulse frequency: 1 to 150 Hz.

III HOW DOES ELECTROSTIMULATION WORK?

The principle of electrostimulation is to stimulate nerve fibres by means of electrical impulses transmitted by electrodes.

The electrical pulses generated by Compex stimulators are high quality pulses - offering safety, comfort and efficiency - which can stimulate different types of nerve fibres:

1. the motor nerves, to stimulate a muscular response. The quantity and the benefits obtained depend on the stimulation parameters and this is known as electro-muscular stimulation (EMS).
2. certain types of sensitive nerve fibres to obtain analgesic or pain-relieving effects.


1. Motor nerve stimulation (EMS)

In voluntary activity, the order for muscular work comes from the brain, which sends a command to the nerve fibres in the form of an electrical signal. This signal is then transmitted to the muscular fibres, which contract.

The principle of electrostimulation accurately reproduces the process observed during a voluntary contraction. The stimulator sends an electrical current impulse to the nerve fibres, exciting them. This excitation is then transmitted to the muscular fibres causing a basic mechanical response (= muscular twitch). The latter constitutes the basic requirement for muscular contraction. This muscular response is completely identical to muscular work controlled by the brain. In other words, the muscle cannot distinguish whether the command comes from the brain or from the stimulator.

The parameters of the Compex programmes (number of impulses per second, contraction time, rest time, total programme time) subject the muscles to different types of work, according to muscular fibres. In fact, different types of muscular fibres may be distinguished according to their respective contraction speed: slow, intermediate and fast fibres. Fast fibres will obviously predominate in a sprinter, while a marathon runner will have more slow fibres.

With a good knowledge of human physiology and a perfect mastery of the stimulation parameters of the various programmes, muscular work can be directed very precisely towards the desired goal (muscular reinforcement, increased blood flow, firming up, etc.).

2. Stimulation of the sensitive nerves

The electrical impulses can also excite the sensitive nerve fibres to obtain an analgesic or pain-relieving effect.

The stimulation of the tactile sensitive nerve fibres blocks the transmission of pain by the nervous system. The stimulation of another type of sensitive fibres creates an increase in the production of endorphins and, therefore, a reduction of pain.

With pain relief programmes, electrostimulation can be used to treat localised sharp or chronic pains as well as muscular pains.

Caution: Do not use the pain relief programmes for a long period without medical advice.

Benefits of electrostimulation

Electrostimulation is a very effective way to make your muscles work:

- with significant improvement of different muscular qualities
- without cardio-vascular or mental fatigue
- with limited stress on the joints and tendons. Electrostimulation thus allows a greater quantity of work by the muscles compared with voluntary activity.

To be effective, this work must involve the greatest possible number of muscular fibres. The number of fibres working depends on the stimulation energy. The maximum tolerable energy should therefore be used. The user controls this aspect of stimulation. The higher the stimulation energy, the greater the number of muscular fibres that are working and, therefore, the more significant the progress achieved. To maximise results, Compex recommends that you complement your electrostimulation sessions with other efforts, such as:

- regular physical exercise
- proper and healthy nutrition
- balanced lifestyle.

IV USAGE GUIDELINES

The usage guidelines presented in this section should be considered as general rules. For all programmes, it is recommended that you read carefully the usage information and advice presented in chapter VII "Programmes".

1. Electrode positions

For optimal results, use the electrode positions recommended by Compex. To do this, refer to the pictures and pictograms shown on the inside cover of the manual.

Each stimulation cable has two poles:

A positive pole (+) = red connection

A negative pole (–) = black connection

A different electrode must be connected to each pole.

Note: It is possible and normal to have an electrode arrangement that leaves one electrode connection free from a cable.

Depending on the characteristics of the current, efficacy can be optimised in certain programmes by placing the electrode connected to the positive pole (red connection) "strategically".

When working with a muscle stimulation programme (programme involving muscle contractions), it is important to place the positive electrode on the motor point of the muscle.

It is crucial to choose the right size electrodes (large or small) and correctly position these on the muscle group you want to stimulate to ensure the efficacy of the treatment. Therefore, always use the size of electrodes shown in the pictures. Unless you have other specific medical instructions, always follow the placement directions in the pictures.

Where necessary, look for the best possible position by slowly moving the positive electrode over the muscle until you find the point that will produce the best contraction or the most comfort for you.

Compex disclaims all responsibility for consequences arising from electrodes placed in other positions.

2. Stimulation positions

This position will vary depending on the position of the electrodes, the muscle group you wish to stimulate, and the programme you are using. For programmes involving powerful muscular contractions, the muscle should always be stimulated in an isometric fashion. You must therefore fix the extremities of your limbs securely. In this way, you provide maximum resistance to the movement and prevent any shortening of the muscle during the contraction, which could create cramp pains and serious stiffness after the session. For example, when stimulating the quadriceps, the user should be in a seated position with the ankles fixed with straps to prevent extension of the knees.

For the other types of programme (for instance, the programmes in the **Pain relief** and **Recovery** categories and the programmes **Cramp prevention**, **Long run optimisation** and **Overcompensation**), which do not induce powerful muscular contractions, position yourself as comfortably as possible.

To determine the stimulation position to be adopted depending on how the electrodes are placed and the programme chosen, refer to the specific applications on our website www.compexrunner.com.

3. Adjusting stimulation energies

In a stimulated muscle, the number of recruited fibres depends on the stimulation energy. For programmes involving powerful muscular contractions, you must therefore use maximum stimulation energies (up to 999), always at the limit that you can endure, in order to recruit the maximum number of fibres.

4. Progression in the levels

In general, it is not advisable to go through the different levels quickly with the intention of reaching level 3 as fast as possible. In fact, the different levels correspond to progress with electrostimulation.

The simplest and most usual procedure is to start with level 1 and raise the level when changing to a new stimulation cycle.

At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session per week at the last level used.

5. Alternation of stimulation sessions / voluntary training

The stimulation sessions can be done outside or during voluntary training.

When voluntary training and stimulation are done during the same session, it is generally recommended that the voluntary training be done first followed by the stimulation. This means that the voluntary training is not done on muscle fibres which are already tired. This is particularly important for strength and explosive strength training. However, in resistance training, it can be very useful to proceed in the reverse order. Before the voluntary training, by means of stimulation for resistance, a "specific prefatigue session" is carried out on the muscle fibres without general and cardiovascular fatigue. In this way, the voluntary effort on the "prepared" fibres will push the glycolytic metabolism faster and further.

V THE **mi** TECHNOLOGY

mi stands for *muscle intelligence*™ (all items referring to this technology are preceded by the symbol **mi**).

This technology takes into account the specific features of each of our muscles and thus offers stimulation adapted to their characteristics.

It's simple, as the data is transmitted to the stimulator automatically!

It's personalised, as each of our muscles is unique!

The Runner stimulator bears the label **mi-READY**, as it can use certain functions of the **mi** technology if a stimulation cable fitted with the **mi-SENSOR** system (sold separately) is connected.

1. Practical rules of use

*To access the **mi** functions make sure that the stimulation cable fitted with the **mi-SENSOR** system (sold separately) is connected before you turn on the stimulator.*

*Avoid connecting the stimulation cable fitted with the **mi-SENSOR** system when the Compex is switched on.*

*To work correctly, the **mi-SENSOR** system must not be blocked or subjected to pressure.*

*During the stimulation session, the stimulation cable fitted with the **mi-SENSOR** system must always be connected to an electrode.*

mi-SENSOR (sold separately)

- This is a small sensor that links the stimulator to the electrodes. **mi-SENSOR** is the key which can be used to measure certain physiological characteristics of the muscle, analyse them and adapt stimulation parameters. This tailoring of each programme to the changing condition of the muscle clearly enhances your comfort and the therapeutic efficacy.

Mi-SCAN

- This function adapts the electrostimulation session to the physiology of each user. Just before starting the work session, **Mi-SCAN** probes the chosen muscle group and automatically adjusts the stimulator parameters to the excitability of this area of the body. It is a truly personalised measurement.
N.B.: The **Mi-SCAN** function operates as soon as a stimulation cable fitted with the **Mi-SENSOR** system (sold separately) is connected to the stimulator.
- This function is implemented at the start of the programme by a short sequence during which measurements are taken (a horizontal bar scans the figurine at the left of the screen).
- Throughout the entire measurement test, you must remain completely still and be perfectly relaxed. The **Mi-SENSOR** system is very sensitive: the smallest contraction or the slightest movement could disturb the measurement test.
- Note that certain people might feel an unpleasant tingling sensation during the test.
- Once the test is complete, the symbol  will appear, indicating that the programme can be started.

Mi-TENS

- The **Mi-TENS** function can be used to considerably reduce the number of unwanted muscular contractions, enhancing both user comfort and treatment efficacy.
- The **Mi-TENS** function can only be accessed for the **Acute tendinitis** and **Ankle twist** programmes.
- For these programmes, the **Mi-TENS** function allows you to control the stimulation energies while maintaining them at an effective level and considerably limiting the onset of muscular contractions.
- Short measurement tests (lasting 2 to 3 seconds) are carried out regularly throughout the entire programme.

- There is a test phase after each increase of the stimulation energies. To ensure the correct performance of these tests, you must remain perfectly still during this time.
- Depending on the results of the measurement tests recorded by the device, the level of the stimulation energies may be slightly decreased automatically.
- It is important always to adopt the most comfortable possible stimulation position. Also, you should try to remain still and not contract the muscles of the stimulated area.

Mi-RANGE


- This function indicates the ideal energy adjustment range for programmes whose efficiency requires vigorous muscular twitches.
- The **Mi-RANGE** function is therefore only available in programmes that use low stimulation frequencies (less than 10 Hz).
- For programmes allowing the use of the **Mi-RANGE** function, the stimulator checks whether you are in the ideal energy range. If you are below the range, the stimulator prompts you to increase the energy by displaying  signs.
- Once the stimulator has detected your ideal adjustment area, a square bracket will appear to the right of the bar graph of the channel to which the **Mi-SENSOR** system is connected. The figures in this bracket indicate the energy range you must work in order to achieve optimal results.
- If you set the energy level below this range, the stimulator will prompt you to increase it by flashing a series of  signs on the screen.

VI HOW IT WORKS

You are strongly advised to carefully read the counter-indications and safety measures described at the start of this manual (chapter 1: "Warnings") prior to using your stimulator.

1. Description of the stimulator

A On/Off button

B "+" button used to:
Increase stimulation energies in several channels simultaneously.
Access the LAST menu (the last programme used).

C Sockets for the 4 stimulation cables

D Stimulation cables

Channel 1 = blue Channel 2 = green
Channel 3 = yellow Channel 4 = red

E +/- buttons for the 4 stimulation channels

F Socket for the charger (*Slide the red cover to the right to free the charger connector*)

G Rechargeable battery compartment

H Belt clip socket

2. Connections

Connecting the cables

The stimulator cables plug into the 4 sockets on the front of the device.

Four cables can be connected simultaneously to the four channels of the device. Both the sockets and the cables are colour-coded to simplify use and facilitate identification of the different channels:
blue = channel 1 green = channel 2
yellow = channel 3 red = channel 4

The Runner stimulator bears the label **Mi-ready**, in other words it can use certain functions of the **Mi** technology if a stimulation cable fitted with the **Mi-SENSOR** system is connected. This cable, sold separately, gives access to the **Mi-SCAN**, **Mi-TENS** and **Mi-RANGE** functions. It is connected to any one of the four stimulation channel sockets located on the front of the device and two electrodes are "clipped" on to it in the same way as for a conventional cable.

Connecting the charger

The Compex has considerable operating autonomy, as it uses rechargeable batteries.

To recharge them, use the charger supplied with your device and connect it to the front of the device, then plug the charger in a socket.

Remove any stimulation cables connected to the stimulator before recharging it.

We strongly recommend you fully charge the battery before using it for the first time as this will improve its autonomy and life span.

3. Preliminary settings

Before using the unit for the first time, you should select the working language of the device which is displayed on the options screen. Follow the instructions below. Afterwards, for the greatest comfort, Compex offers you a number of setting options (choice of language, screen contrast adjustment, adjusting backlighting and sound volume setting). To change any of these settings, bring up the options screen by holding the On/Off button on the left of the stimulator for a few seconds.


fig.1

1b Use the channel 1 +/- button to choose the language you wish to use.

1c Use the channel 2 +/- button to adjust the contrast of the screen.

1d Use the channel 3 +/- button to adjust the volume.

1e Use the channel 3 +/- button to adjust the volume.

On : backlighting always on.

Off : backlighting always off.


Auto : backlighting activated whenever a button is pressed.

1a Press the On/Off button to validate the selected parameters. Your stimulator saves your options. It is now ready for use with the settings you selected.

4. Selecting a programme category

To switch on your stimulator, briefly press the On/Off button on the left of the Compex. A musical signal can be heard and a screen showing the different programme categories is displayed.

Before choosing a programme, you should select the desired category.

2a Press the On/Off button to switch off the unit.


2b Use the channel 1 +/- button to choose the category you wish to use.

2e Press the channel 4 +/- button to confirm your choice and access the programme selection screen.

N.B.: Press the "i" button to access the LAST menu.

5. Selecting a programme

To choose a programme, it is particularly useful to consult chapter VII "Programmes".

3a Press the On/Off button to return to the previous screen.

3b Press the channel 1 +/- button to choose a programme.

3e Press the channel 4 +/- button to validate your selection and, depending on whether you have the symbol **START** or *********


a) **START** = the stimulation session starts immediately

b) ********* = a parameter setting screen appears.

N.B.: Press the i button to access the LAST menu.

LAST

To make your Compex more pleasant and easier to use, the LAST menu allows you to directly access the last used programme. To access the menu, press the "i" button before selecting a programme. This function can be accessed from the programme category screen (fig. 2), or the programme list screen (fig. 3).

4a Press the On/Off button to return to the previous screen.

4d Press the channel 3 +/- button to personalize the programme.

4e Press the channel 4 +/- button to execute the programme.

6. Personalising a programme

The programme personalisation screen is not available for all programmes!


fig.5

5a Press the On/Off button to return to the previous screen.

5b Some programmes require manual selection of the muscle group you want to stimulate. This muscle group is shown in black on a small figurine above channel 1. Use the channel 1 +/- button to select your chosen group. The seven muscle groups proposed are shown successively in black on the small figurine:

-  Abdomen and lower back
-  Buttocks
-  Thighs
-  Legs and feet
-  Forearms and hands
-  Shoulders and arms
-  Thorax and back

Current international standards require that a warning be given concerning the application of electrodes to the thorax (increased risk of cardiac fibrillation).

N.B.: If a cable fitted with the **mi**-SENSOR system (sold separately) is connected to the device, the muscle group is selected automatically.

5c Press the channel 2 +/- button to stop the warming-up sequence (closing the small animated rising convection symbol above the radiator).

5d Press the channel 3 +/- button to choose the programme difficulty level.

5e Press the channel 4 +/- button to confirm your choices and launch the programme.

7. During the stimulation session

Preliminary **mi**-SCAN test

*The **mi**-SCAN test is launched only if a cable fitted with the **mi**-SENSOR system (sold separately) is already connected to the stimulator when this is turned on!*

It is essential that you remain completely still and relaxed during the entire test.

If the **mi**-SENSOR system is connected, the test starts immediately after programme selection and personalisation.


fig.6


6a Press the On/Off button to interrupt the test at any moment. To relaunch the test from the start, press any of the 4 +/- buttons.

6e The **mi** logo above channel 4 indicates that the **mi**-SENSOR cable is connected to the stimulator. The figures in channel 4 will automatically change during the test as the sensor records different measurements. The 4 +/- buttons will be disabled for the duration of the test.

Once the test is complete, you will automatically be taken to the default programme start screen and prompted to increase the stimulation energies.

Adjusting stimulation energies

When you start a programme, you will be prompted to increase the stimulation energies. This is critical to the success of any treatment. To find out what energy level needs to be reached for each programme, refer to the specific applications.


7a Press the On/Off button to place the unit in Pause mode.

7b c d e The Compex beeps and the symbols of the four channels flash, changing from + to 000: the four channels are at 000 energy. You must increase the stimulation energy so that the stimulation can start. To do this, press the + buttons for the relevant channels until the desired setting is reached. If you want to increase the energy level of all four channels simultaneously, press the "i" button, located below the On/Off button. Press the "i" button twice to increase the levels in the first 3 channels, and 3 times to increase the levels in the first 2 channels. When you press the "i" button, the associated channels are highlighted in white on a black background.

Programme progression

Stimulation actually starts when the stimulation energy has been increased. The examples reproduced below are sufficient to understand the general rules.

8a Press the On/Off button to interrupt the programme momentarily. To restart it, simply press the channel 4 +/- button. The session will resume at 80% of the energy levels that were in use prior to the interruption.

8 b c d e The different energies reached during the contraction phase are shown by a series of black bar graphs, that of the rest phase energies are shown by hatched bar graphs.

N.B.: Active rest phase stimulation energies are automatically set at 50% of contraction energies. These can be modified during the rest phase. Once modified, they will be totally independent of the contraction energies.

N.B.: If your Compex emits a beeping sound and the symbols **+** under the active channels begin to flash the stimulator is suggesting you increase the level of the stimulation energies. If you are working at the maximum tolerance level, simply ignore this message.

Statistics

Your stimulator includes a statistics menu that allows you to view in real time the important information for a programme. To access the statistics screen, you must put the stimulator into the "Pause" mode or wait for the end of the programme.


fig.9


fig.10

Fig. 9 For programmes inducing muscular contractions, the indication "MAX" shows the maximum energy level reached per channel during contraction phases. "AVG" indicates the average stimulation energy level used by the different channels during the contraction phase.

Fig. 10 In low frequency programmes, the **Mi-RANGE** function, symbolised by a single square bracket shows what percentage of stimulation time was spent in the optimal energy range.

9e-10e Press the channel 4 +/- button to resume the programme at the point where it was interrupted.

End of programme

At the end of each session, a small flag will be displayed on the screen and a short melody will be played. To switch off the stimulator, press the On/Off button. Depending on the programmes, usage statistics may be displayed (cf. **fig. 9** and **10**). It should be noted, however, that statistics are not available for all programmes.


fig.11

8. Consumption and recharging

Never recharge the stimulator without first disconnecting the stimulation cables. Always use the charger supplied by Compex to recharge the batteries.

The Compex runs on rechargeable batteries. The battery's life depends on the programmes and stimulation energy used. We strongly recommend you fully charge the battery before using your Compex for the first time as this will improve its autonomy and life span. If you do not use your device for a long period of time, please recharge the battery regularly.

Electricity consumption

The charge state of the battery is indicated by a small battery icon on the screen.


fig.12

If the battery icon contains just two lines, this means that power is running low. Stop the session and recharge the unit.

If the symbol **START** normally displayed above the channel 4 +/- button is not visible and if the small battery icon is flashing, the batteries are flat. It is no longer possible to use the device. Recharge immediately.

Recharging

Remove all stimulation cables from the stimulator before recharging it. Connect the charger to the mains and then connect the stimulator to the charger. The charge menu illustrated below appears automatically.


fig.13

The duration of the charging operation is shown on the screen (a complete charge may take 1.5 hours). The battery icon is animated while the battery is recharging. When fully charged, the icon will be full and the total time taken to recharge the battery will flash on the screen. Simply disconnect the charger: the Compex will turn off automatically.

9. Problems and solutions

Electrode fault


fig.14

The Compex beeps and alternatively displays the symbol of a pair of electrodes and an arrow pointing to the channel where a problem has been detected. In the example above, the stimulator has detected an error in channel 1.

- Check that electrodes are connected to this channel.
- Check whether the electrodes are old, worn, and/or the contact is poor: try using new electrodes.
- Try using the stimulation cable on a different channel. If the cable is still showing a fault, replace it.

Stimulation is not producing the usual sensation

- Check that all the settings are correct and ensure the electrodes are positioned properly.
- Change the positioning of the electrodes slightly.

The stimulation effect causes discomfort

- The electrodes are beginning to lose adhesion and no longer provide good contact on the skin.
- The electrodes are worn and need to be replaced.
- Change the positioning of the electrodes slightly.

The stimulator is not working

fig.15

- If an error screen appears while the equipment is in use, note the error number (in the example in **fig. 15**, error 1/0/0) and contact the consumer service stated and approved by Compex Médical SA.

VII PROGRAMMES

Programmes	Effects	Uses	Stimulation energies	mi Function if mi-SENSOR cable (sold separately) connected
Prevention category				
Ankle twist prevention	Increase in the contraction speed and strength of the lateral peroneal muscles	To avoid recurrence of ankle sprain, after the initial treatment allows a return to sports activities	Maximum tolerable energy (0-999)	mi-SCAN
Cramp prevention	Improvement of the circulatory system in order to prevent nocturnal cramps or stresses	Cramps involve fundamentally the muscles of the lower limbs, more particularly the calf muscles	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
Low back pain prevention	Improving the support to the lower back region by the abdominal and lower back muscles	To reduce the risk of pain occurring to the lower back	Maximum tolerable energy (0-999)	mi-SCAN

Programmes	Effects	Uses	Stimulation energies	mi Function if mi-SENSOR cable (sold separately) connected
Pain relief category				
Acute tendinitis	Pain relief by blocking the transmission of pain impulses	One or more times a day during the acute phase of the condition	Gradually increase the stimulation intensities until a very clear tingling sensation is felt	mi-TENS
Chronic tendinitis	Improvement of the tendon recovery phase by local increase of the blood circulation and antipain effect by the release of endorphin	After disappearance of the inflammatory signs during the chronic phase of the condition	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
Low back pain	Analgesic action through the release of endorphins Increases blood flow	Analgesic current specifically adapted to persistent low back pain (lumbar region)	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
Ankle twist	Pain relief by blocking the transmission of pain impulses	Several sessions per day interspersed with repeated applications of ice	Gradually increase the stimulation intensities until a very clear tingling sensation is felt	mi-TENS
Decontracture	Reduced muscular tension Relaxing and decontraction effect	To combat recent and localized muscular pain	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE

Programmes	Effects	Uses	Stimulation energies	mi Function if mi -SENSOR cable (sold separately) connected
Training category				
Core stabilization	Improvement of the stabilisation of trunk muscles	During training or maintenance throughout the sports season	Maximum tolerable energy (0-999)	mi -SCAN
Endurance	Increasing oxygen consumption by the muscles Improvement of performance to sustain an effort over a long period of time	For runners having a specific voluntary training volume less than 5 sessions per week	Maximum tolerable energy (0-999)	mi -SCAN
Strength	Increasing the maximum strength that the muscle is capable of developing Improving stride and better running economy	For runners having a specific voluntary training volume of more than 4 sessions per week	Maximum tolerable energy (0-999)	mi -SCAN
Strength trail specific	Increased strength of muscular contraction Improved running on hills and rough terrain	For runners preparing to test their running in harsh terrain	Maximum tolerable energy (0-999)	mi -SCAN
Long run optimisation	Marked local increase in blood flow Development of the blood capillary system in the muscle stimulated	After the standard long training session, to increase the benefits without imposing additional stress on the joints	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi -SCAN mi -RANGE
Overcompensation	Increase of the blood circulation in the muscle stimulated Development of the blood capillary system in the muscle stimulated	2 sessions per day on the last 4 days before the competition	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi -SCAN mi -RANGE

Programmes	Effects	Uses	Stimulation energies	mī Function if mī-SENSOR cable (sold separately) connected
Fartlek	Training and preparation of muscles for all kinds of muscular work (endurance, resistance, strength, explosive strength) through different working sequences	At the beginning of the season to "re-set" the muscles after a rest period and before more intensive and specific training During the season for persons who do not have any preferences for a single kind of performance and prefer to exercise their muscles with different workout methods	Maximum tolerable energy (0-999)	mī-SCAN

Programmes	Effects	Uses	Stimulation energies	mi Function if mi-SENSOR cable (sold separately) connected
Recovery category				
Active recovery	Marked increase in blood flow Accelerates elimination of waste from muscular contraction Endorphinic effect Relaxing and decontraction effect	To improve and accelerate muscle recovery after intensive exercise To be used during the first 3 hours after each session of intensive training or after a competition	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
Post-race recovery	Significant increase in the blood flow helping to drain the organic waste due to effort Analgesic and decontracting effects	To achieve the best possible and fastest muscular recovery after gruelling effort The particular rhythm of the session reduces the risk of muscular cramps appearing	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
Reduce muscle soreness	Significant increase in the blood flow helping to drain the organic waste due to effort Analgesic and decontracting effects	Reduced duration and intensity of muscular soreness To be used on the 2 or 3 days following the competition	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE

VIII EMC TABLE

(Electromagnetic Compatibility)

The Compex Runner needs special EMC precautions and must be installed and started according to the EMC information supplied in this manual.

All RF wireless transmission systems can affect the Compex Runner.

The use of accessories, sensors and cables other than those recommended by the manufacturer may result in stronger emissions or reduce the immunity of the Compex Runner.

The Compex Runner should not be used beside or stacked on top of any other equipment. If you must use it side by side or on top of another system, you should check that the Compex Runner works properly in the chosen configuration.

RECOMMENDATIONS AND DECLARATION BY THE MANUFACTURER CONCERNING ELECTROMAGNETIC EMISSIONS

The Compex Runner is intended for use in the electromagnetic environment specified below. The customer or user of the Compex Runner should ensure that it is used in this environment.

Emissions test	Compliance	Electromagnetic environment - guide
CISPR 11 RF emissions	Group 1	The Compex Runner uses RF energy only for its internal operation. Consequently, its RF emissions are very low and are unlikely to interfere with any adjacent electrical device.
CISPR 11 RF emissions	Class B	The Compex Runner is suitable for use in any establishment, including a private dwelling and a place connected directly to the low voltage mains supply which powers residential buildings.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ emission oscillations IEC 61000-3-3	Not applicable	

RECOMMENDED SPACING BETWEEN A PORTABLE AND MOBILE COMMUNICATION APPLIANCE AND THE COMPEX RUNNER

The Compex Runner is designed for use in an electromagnetic environment in which radiated RF waves are controlled. The buyer or user of the Compex Runner can contribute to preventing electromagnetic interference by maintaining a minimum distance between RF portable and mobile communication appliances (transmitters) and the Compex Runner according to the table of recommendations below and according to the maximum output power of the telecommunication appliance.

Maximum transmitter output power W	Spacing according to the frequency of the transmitter mCISPR 11		
	From 150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	From 80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	From 800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

In the case of transmitters whose maximum output power is not shown in the table above, the recommended spacing of d metres (m) can be calculated using the appropriate equation for the transmitter frequency, where P is the maximum output power of the transmitter in watts (W) as set by the transmitter manufacturer.

NOTE 1: At 80 MHz and at 800 MHz,

the spacing for high frequency amplitude is applied.

NOTE 2: These guidelines may not be appropriate for some situations. Electromagnetic wave propagation is modified by absorption and reflection due to buildings, objects and persons.

RECOMMENDATIONS AND DECLARATION BY THE MANUFACTURER -

The Compex Runner is designed for use in the electromagnetic environment stipulated below. The buyer or user of the Compex Runner must ensure it is used in this recommended environment.

Immunity test	Test level IEC 60601	Observance level
Electrostatic discharge (DES) IEC 61000-4-2	± 6 kV at the contact ± 8 kV in air	± 6 kV at the contact ± 8 kV in air
Rapid electrical transients in bursts of noise IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable System battery-powered
Surge currents IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Not applicable System battery-powered
Voltage trough, short power cuts and voltage variations on power supply input lines IEC 61000-4-11	< 5% VT (trough > 95% of VT) for 0.5 cycle < 40% VT (trough > 60% of VT) for 5 cycles < 70% VT (trough > 30% of VT) for 25 cycles < 5% VT (trough > 95% of VT) for 5 seconds	Not applicable System battery-powered
Magnetic field at mains frequency (50/60 Hz) IEC 61000-4-8	3 A/m	

NOTE: VT is the AC supply voltage before application of the test level.

RECOMMENDATIONS AND DECLARATION BY THE MANUFACTURER -

The Compex Runner is designed for use in the electromagnetic environment stipulated below. The buyer or user of the Compex Runner must ensure it is used in this recommended environment.

Immunity test	Test level IEC 60601	Observance level
RF immunity IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Signal lines Not applicable for power supplies of devices which are battery-powered
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz 10 V/m 26 MHz to 1 GHz	3 A/m 10 V/m

NOTE 1: From 80 MHz to 800 MHz, the spacing for high frequency amplitude is applied.

NOTE 2: These guidelines may not be appropriate for some situations.

a The field intensity from fixed transmitters, such as radio telephone base stations (cellular/wireless) and a mobile radio, amateur radios, AM and FM radio emissions and TV emissions cannot be predicted with any accuracy. It may therefore be necessary to consider an analysis of the electromagnetic environment of the site to calculate the electromagnetic environment coming from RF fixed transmitters. If the field intensity measured in the environment where the Compex Runner is located exceeds the appropriate RF observance level above, the Compex Runner should be monitored to ensure it is operating properly.

NOTE: VT is the AC supply voltage before application of the test level.

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© 12/2011 DJO, LLC