### Intended Use

**CardioPhase\* hsCRP** is an *in vitro* diagnostic reagent for the quantitative determination of C-reactive protein (CRP) in human serum, and heparin and EDTA plasma by means of particle enhanced immunonephelometry using BN\* Systems. In acute phase response, increased levels of a number of plasma proteins, including C-reactive protein, is observed. Measurement of CRP is useful for the detection and evaluation of infection, tissue injury, inflammatory disorders and associated diseases. High sensitivity CRP (hsCRP) measurements may be used as an independent risk marker for the identification of individuals at risk for future cardiovascular disease. Measurements of hsCRP, when used in conjunction with traditional clinical laboratory evaluation of acute coronary syndromes, may be useful as an independent marker of prognosis for recurrent events, in patients with stable coronary disease or acute coronary syndromes.

### Summary and Explanation

CRP is one of the 'acute-phase' proteins, the serum or plasma levels of which rise during general, nonspecific response to infectious and non-infectious inflammatory processes. CRP is synthesized in the liver and is normally present as a trace constituent of serum or plasma.

In various disease states resulting in tissue injury, infection or acute inflammation, CRP values may rise above normal to 20 to 500 mg/L<sup>1</sup>. As elevated CRP values are always associated with pathological changes the CRP assay provides useful information for the diagnosis, therapy and monitoring of inflammatory processes and associated diseases<sup>2,3,4</sup>. Increases in CRP values are non-specific and should not be interpreted without a complete clinical history. Studies have shown that measurement of CRP by high sensitivity assays is a strong independent predictor of risk for future cardiovascular and peripheral vascular diseases<sup>5</sup>. High sensitivity CRP measurements have also been shown to add to the predictive value of other markers used to assess the risk of cardiovascular and peripheral vascular disease<sup>6,7,8</sup>. Elevated CRP values determined by high sensitivity CRP assays may be indicative of the prognosis of individuals with acute coronary syndromes, and may be useful in the management of such individuals<sup>9</sup>.

### Principle of the Method

Polystyrene particles coated with monoclonal antibodies specific to human CRP are aggregated when mixed with samples containing CRP. These aggregates scatter a beam of light passed through the sample. The intensity of the scattered light is proportional to the concentration of the relevant protein in the sample. The result is evaluated by comparison with a standard of known concentration. The assigned value of CRP in N Rheumatology Standard SL is standardized against the international reference preparation BCR - CRM 470<sup>10,11</sup>.

### Reagents

#### Materials provided

**CardioPhase\* hsCRP Reagent, Code No. OQIY**

Five vials containing 5 mL each or three vials containing 2 mL each

#### Composition

**CardioPhase\* hsCRP Reagent** consists of a suspension of polystyrene particles coated with mouse monoclonal antibodies to CRP.

#### Preservatives:

**CardioPhase\* hsCRP Reagent:** Gentamicin 6.25 mg/L, Amphotericin 0.625 mg/L

#### Warnings and Precautions

- For *in vitro* diagnostic use.
- Each individual blood donation intended for use in the manufacture of the **CardioPhase\* hsCRP Reagent** was tested for HBSAg, anti-HCV, anti-HIV1 and anti-HIV2 by FDA required testing. Only donations with negative findings were used for manufacture. Nevertheless, since absence of infectious agents cannot be proven, all materials obtained from human tissue or body fluids should always be handled with due care, observing the precautions recommended for biohazardous material<sup>11</sup>.

#### Preparation of the Reagent

The **CardioPhase\* hsCRP Reagent** is ready-for-use as supplied, shake gently before first use.

#### Storage and Stability

**Shelf life at +2 to +8 °C:**

See expiry date on the label.

#### Stability once opened:

Four weeks if stored at +2 to +8 °C securely capped immediately after use; do not freeze.

#### On-board stability:

A minimum of five days, at eight hours per day, or comparable period of time.

**Note:** On-board stability may vary, depending on the BN\* System used and laboratory conditions. For further details, refer to the BN\* II and BN ProSpec® System Instruction Manual.

#### Materials Required but not Provided

BN\* System

N Rheumatology Standard SL, Code No. OQKZ

N/T Rheumatology Controls SL/1 and SL/2, Code Nos. OQDB and OQDC

Apolipoprotein Control Serum CHD, Code OUPH

N Supplementary Reagent/Precipitation, Code No. OUMU

N Diluent, Code No. OUMT

BN\* II Evaporation Stoppers (optional), Code No. OVLE

Additional materials and supplies as described in your BN\* System Instruction Manual.

#### Specimens

Suitable samples are human serum as well as heparin and EDTA plasma, either as fresh as possible (stored for no more than eight days at +2 to +8 °C) or stored frozen. Samples can be stored at below -20 °C for up to eight months, if they are frozen within 24 hours after collection and if repeated freeze-thaw cycles are avoided. Serum samples must be completely coagulated and, after centrifugation, must not contain any particles or traces of fibrin. Lipemic samples, or frozen samples, which become turbid after thawing, must be clarified by centrifugation (10 minutes at approximately 15,000 x g) prior to testing.

#### Procedure

Notes:

- Consult your BN\* System Instruction Manual for details regarding operation of the instrument.
- Allow reagents and samples to equilibrate to room temperature (+15 to +25 °C) before use on the BN\* 100 System. With a BN\* II or BN ProSpec® System reagents and samples stored at +2 to +8 °C can be used immediately.
- On the BN\* 100 System, samples should be run at approximately the same ambient temperature (maximum 3 °C deviation) as the measurements used for recording the reference curve.

#### Assay Protocols on the BN\* Systems

The assay protocols, for serum as well as plasma, are given in the BN\* System Instruction Manual and software of the instrument. All steps are performed automatically by the system.

#### Establishment of the Reference Curve

Reference curves are generated by multi-point calibration. Serial dilutions of N Rheumatology Standard SL are automatically prepared by the instrument using N Diluent. The standard dilutions are to be used within four hours. The reference curve is valid for four weeks and can be used beyond this period of time, as long as controls with corresponding method depending target values, e.g., N/T Rheumatology Controls SL/1 and SL/2 or Apolipoprotein Control Serum CHD, are reproduced within their respective confidence interval. If a different lot of reagent is used, a new reference curve must be generated. The exact measuring range depends upon the concentration of the protein in each lot of N Rheumatology Standard SL. For typical ranges refer to your BN\* System Instruction Manual.

#### Assay of Specimens

Samples are automatically diluted 1:400 (CRP1) or 1:20 in the **CardioPhase\* hsCRP** assay protocol (CRP2) with N Diluent. The diluted samples must be used within four hours. If the results obtained are outside the measuring range, the assay can be repeated using a higher or lower (only in the CRP1 assay protocol) dilution of the sample. Refer to the BN\* System Instruction Manual for information on repeat measurements using other dilutions.

#### Internal Quality Control

Assay N/T Rheumatology Controls SL/1 and SL/2 or Apolipoprotein Control Serum CHD after each establishment of a reference curve, the first use of a reagent vial as well as with each run of samples. The controls are to be assayed and evaluated as for patient samples. The assigned value and the confidence interval are listed in the Table of Assigned Values of the respective control. If a control value is outside the confidence interval, the determination must be repeated. If the repeated determination confirms the deviation, a new reference curve should be established. Do not release patient results until the cause of the deviation has been identified and corrected.

### Results

The results are evaluated automatically by the analyzer and are represented in mg/L or in a unit selected by the BN\* System user.

### Limitations and Interferences

Turbidity and particles in the samples may interfere with the determination. Therefore, samples containing particles must be centrifuged prior to testing. Lipemic or turbid samples, which cannot be clarified by centrifugation (10 minutes at approximately 15,000 x g), must not be used.

No interference was observed from bilirubin up to 600 mg/L, free hemoglobin up to 10 g/L and triglycerides up to 16 g/L.

Patient samples may contain heterophilic antibodies that could react in immunoassays to give a falsely elevated or depressed result. This assay has been designed to minimize interference from heterophilic antibodies. Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed. Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Due to matrix effects, inter-laboratory survey samples and control samples may yield results that differ from those obtained with other methods. It may therefore be necessary to assess these results in relation to method specific target values.

### AHA/CDC Expert Panel Recommendations<sup>5</sup>:

hsCRP levels should not be substituted for assessment of traditional cardiovascular risk factors. Application of management guidelines for acute coronary syndromes should not be dependent on hsCRP levels.

In patients with stable coronary disease or acute coronary syndromes, hsCRP measurement may be useful as an independent marker of prognosis.

When using the assay for risk assessment, patients with persistently unexplained, marked elevation of hsCRP (>10 mg/L) after repeated testing should be evaluated for non-cardiovascular etiologies.

The expert panel recommends against screening of the entire adult population for hsCRP as a public health measure.

Patients with evidence of active infection, systemic inflammatory processes or trauma should not be tested for cardiovascular disease risk assessment until these conditions have abated.

Application of secondary prevention measures should not depend on hsCRP determination, but rather an array of risk factors (global risk assessment).

Serial measurements of CRP should not be used to monitor effects of treatment.

Two separate CRP measurements (optimally two weeks apart) should be obtained before performing risk assessment, due to within-subject CRP variability.

Measurement of hsCRP is an independent marker of risk.

hsCRP levels may be useful in motivating patients to improve lifestyle behaviors.


### Interpretation of Results

Recent medical events resulting in tissue injury, infection or inflammation, which may cause elevated CRP levels should also be considered when interpreting results.

### Cardiovascular and Peripheral Vascular Disease

#### Risk Assessment:

Several studies have examined the utility of hsCRP measurements for cardiovascular disease risk prediction. In a subset consisting of 280 individuals participating in the Physician's Health Study, baseline samples were tested for a number of biomarkers, including hsCRP. The population was followed for nine years to assess the development of peripheral arterial disease (PAD). The following graph shows the relative risk (95% confidence interval) of developing future PAD as a function of quartile of hsCRP<sup>12</sup>.


The AHA/CDC Scientific Statement<sup>6</sup> provides the following risk assessment guidelines:

Risk	hsCRP (mg/L)
Low	< 1.0
Average	1.0 - 3.0
High	> 3.0

### Risk Stratification:

For patients with acute coronary syndromes, measurement of hsCRP may provide prognostic information. A value of CRP > 10 mg/L in the early period (6 - 24 hours after onset of symptoms), has been shown to be indicative of an increased risk for short term (30 days - 1 year) recurrent cardiac events<sup>5</sup>.

A substudy of 447 patients in the CAPTURE trial examined the clinical implications of elevated levels of CRP for risk stratification in patients with unstable angina. As shown in the following graph, patients with a CRP > 10 mg/L experienced a higher event rate (mortality or MI) than patients with a CRP < 10 mg/L<sup>13</sup>.

## Reference Interval

Expected values for healthy individuals as noted in the literature are typically  $\leq 3$  mg/L<sup>15</sup>. The normal range of CRP in the serum of 2147 apparently healthy individuals using the **CardioPhase\* hsCRP Assay** was found to be<sup>16</sup>:

90%	1.69 mg/L
95%	2.87 mg/L

A subset of the Stanislas Cohort was examined in this study<sup>17</sup>. The cohort subset used consisted of 1151 males and 996 females ranging in age from 5 to 71 years. All participants were of European ancestry, and free of previously diagnosed serious and/or chronic diseases (such as cancer or cardiovascular disease) and excluded individuals taking anti-inflammatory drugs or antibiotics.

As CRP is a nonspecific indicator for a wide range of disease processes, and as the reference intervals are affected by many factors that may differ for each population studied, each laboratory should determine its own reference interval.

## Specific Performance Characteristics

### Sensitivity

The analytical sensitivity of the assay is determined by the lower limit of the reference curve and therefore depends upon the concentration of the protein in the N Rheumatology Standard SL. A typical limit of detection for CRP is 0.175 mg/L for measurements performed using a sample dilution of 1:20.

A coefficient of variation (CV) of 7.6% was observed from ten replicates of a sample containing 0.41 mg/L of CRP.

### Specificity

No cross-reactivity is known.

### Precision

Measurements with **CardioPhase\* hsCRP Reagent** yielded the following coefficients of variation (CV):

Precision			
CRP Concentration (mg/L)	Intra-assay (n=20 each) CV (%)	CRP Concentration (mg/L)	Inter-assay (n=10 each) CV (%)
0.5	3.1	0.5	2.5
1.1	3.8	1.3	3.8
2.1	3.4	2.1	2.1
15	4.0	14	2.6
26	2.3	24	3.9
62	4.4	56	5.7

### Method Comparison

71 serum samples with CRP concentrations of up to approximately 500 mg/L were analyzed in parallel using the **CardioPhase\* hsCRP Reagent** (y) and the N Latex CRP kit (x). Regression analysis yielded the following regression equation:

$$y = 0.95 + 1.43 \text{ mg/L} \cdot x$$

**Note:** The values cited for specific performance characteristics of the assay represent typical values and are not to be regarded as specifications for the **CardioPhase\* hsCRP assay**.

## Bibliography

1. Kushner I, Rzewnicki DL. The acute phase response: general aspects. *Baillieres Clin Rheumatol* 1994; 8: 513-30.
2. Dixon JS, Bird HA, Sitton NG et al. C-reactive protein in the serial assessment of disease activity in rheumatoid arthritis. *Scand J Rheum* 1984; 13:39-44.
3. Hind CRH, Pepys MB. The role of serum C-reactive (CRP) measurement in clinical practice. *Int Med* 1984; 5:112-51.
4. van Leeuwen M, van Rijswijk MH. Acute phase proteins in monitoring of inflammatory disorders. *Baillieres's Clinical Rheumatology* 1994; 8:531-52.
5. Pearson TA, Mensah GA, Alexander RW, et al. Markers of Inflammation and Cardiovascular Disease: application to clinical and public health practice: A statement for healthcare professionals from the Centers for Disease Control and Prevention and the American Heart Association. *Circulation* 2003;107:499-511.
6. Ridker PM, Cushman M, Stampfer MJ, et al. Plasma concentration of C-Reactive Protein and risk of developing peripheral vascular disease. *Circulation* 1998; 97:425-8.
7. Ridker PM, Glynn RJ, Hennekens CH. C-reactive protein adds to the predictive value of total and HDL cholesterol in determining risk of first myocardial infarction. *Circulation* 1998; 97:2007-11.
8. Ridker PM, Rifai N, Rose L, et al. Comparison of C-reactive protein and low-density lipoprotein cholesterol levels in the prediction of first cardiovascular events. *N Engl J Med* 2002; 347: 1557-65.
9. Baudner S, Bienvu J, Blirup-Jensen S, et al. The certification of a matrix reference material for immunochemical measurement of 14 human serum proteins. CRM 470. Brussels: Community Bureau of Reference, Commission of the European Communities. BCR Information, Reference Materials. 1993 report EUR 15243 EN (ISSN 1018-5593): 1-172
10. Whicher JT, Ritchie RF, Johnson AM, et al. New international reference preparation for proteins in human serum (RPPHS). *Clin Chem* 1994; 40:934-8.
11. U.S. Department of Health and Human Services CDC. Biosafety in Microbiological and Biomedical Laboratories, HHS Publication (CDC) 93-8395; 1999; Section II; 8-16.
12. Graph adapted from Ridker PM, Stampfer MJ, Rifai N. Novel risk factors for systemic atherosclerosis: a comparison of C-reactive protein, fibrinogen, homocysteine, lipoprotein(a), and standard cholesterol screening as predictors of peripheral arterial disease. *JAMA* 2001; 285: 2481-5.
13. Graph adapted from Heeschen C, Hamm C, Brummner J, et al. Predictive value of C-reactive protein and troponin T in patients with unstable angina: A comparative analysis. *J Am Coll Cardiol* 2000; 35: 1535-42.
14. Tracy RP, Lemaitre RN, Psaty BM, et al. Relationship of C-reactive protein to risk of cardiovascular disease in the elderly. Results from the Cardiovascular Health Study and the Rural Health Promotion Project. *Arterioscler Thromb Vasc Biol* 1997; 17: 1121-7.
15. Macy EM, Hayes TE, Tracy RP. Variability in the measurement of C-reactive protein in healthy subjects: implications for reference intervals and epidemiological applications. *Clin Chem* 1997; 43: 52-8.
16. Cheniliot O, Henry J, Steinmetz J, et al. High sensitivity C-reactive protein: biological variations and reference limits. *Clin Chem Lab Med* 2000; 38:1003-11.
17. Siest C, Visvikis S, Herbeth B, et al. Objectives, design and recruitment of a familial and longitudinal cohort for studying gene-environment interactions in the field of cardiovascular risk: the Stanislas cohort. *Clin Chem Lab Med* 1998; 36: 35-42

BN ProSpec is a registered trademark of Dade Behring Marburg GmbH in the USA, in Germany and other countries.

\* BN and **CardioPhase** are trademarks of Dade Behring Marburg GmbH in the USA.

Dade Behring Marburg GmbH  
Emil-von-Behring-Str. 76  
D-35041 Marburg  
[www.dadebehring.com](http://www.dadebehring.com)

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\* BN e **CardioPhase** são marcas da Dade Behring Marburg GmbH nos EUA.
BN Pro Spec é uma marca registada da Dade Behring Marburg GmbH nos EUA, Alemanha e outros países.