HARTMANN

# Tensoval®

comfort

Instructions for use  
Instrucciones de uso  
Manual de utilização  
Οδηγίες χρήσεως  
Kullanım Kılavuzu  
تعليمات للاستعمال

 Please read this introduction carefully before use as correct measurement can only be carried out following these instructions. Should you have any further questions on our product or on measuring blood pressure, please contact our Customer Service (a list of addresses is given in the warranty document). On the Internet, you can visit our homepage on [www.hartmann.info](http://www.hartmann.info)

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Congratulations on your purchase of this HARTMANN quality product. Tensoval comfort is a fully automatic blood pressure monitor which uses intelligent measuring technology (Fuzzy Logic) to take rapid and reliable measurement of systolic and diastolic blood pressure as well as pulse rate using an oscillometric measuring method.

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
## 1. General information on blood pressure

The human heart beats approx. 60 to 80 times per minute. This causes blood to be pumped into the arterial vascular system and supplies the body with oxygen and necessary nutrients. A constant pressure known as blood pressure is required to enable the blood to flow into the smallest of blood vessels. In measuring blood pressure, two values are determined. The **systolic value** is the maximum pressure as the heart beats. **Diastolic value** is the

minimum pressure between two heart-beats. Blood pressure is expressed as, e.g. 120/80 in millimetres of a mercury column (mmHg).

Blood pressure changes constantly in each individual and thus creates the requirements necessary for the functional capacity of the body. Thus fluctuations in blood pressure are perfectly normal. However, if blood pressure values are constantly elevated at rest, **hypertension** or high blood pressure is present.

The World Health Organisation (WHO) has set out the following limit values:

Assessment	Systolic pressure	Diastolic pressure
Optimal	up to 120 mmHg	up to 80 mmHg
Normal	up to 130 mmHg	up to 85 mmHg
Normal limit values	130 – 139 mmHg	85 – 89 mmHg
Grade 1 hypertension	140 – 159 mmHg	90 – 99 mmHg
Grade 2 hypertension	160 – 179 mmHg	100 – 109 mmHg
Grade 3 hypertension	over 180 mmHg	over 110 mmHg

High blood pressure is one of the commonest causes of disability and death. And yet many people do not know they are suffering from high blood pressure. This is because the problem is only noticed when the conditions it leads to are at an advanced stage. These may be life-threatening, such as heart attack, stroke or kidney damage. Only regular blood pressure monitoring enables high blood pressure to be detected early. You have decided in favour of **blood pressure self-measurement** and this will contribute to your health care. This is the only way to ensure efficient protection against the dangerous results which may affect your heart and circulation.

## 2. Important notes for self-measurement

- Self measurement of blood pressure does not replace either treatment or **regular check-ups by your doctor**. The drug dose prescribed by your doctor may thus on no account be altered without consultation.
- Self-measurement of blood pressure must be conducted regularly and over a long period of time. With the assistance of the **blood pressure values recorded regularly by you in your blood pressure diary**, high blood pressure can be detected early and your doctor can initiate optimal treatment.

Individual readings are situation-related and thus are not useful.

- In order to compare values, it is important that you **always** take your blood pressure **at rest**. So you should relax for 5 minutes before taking a measurement.
- Even slight changes in internal and external factors (e.g. deep breathing, stimulants, talking, excitement, climatic factors) lead to **fluctuations in blood pressure**. This is why your doctor and pharmacist often obtain different readings.
- Measure your blood pressure regularly in the morning and evening **at about the same time**, as blood pressure fluctuates during the course of the day.
- You should **not move or talk** during the measurement process.
- Always measure blood pressure on the **naked upper arm**, and make sure that items of clothing rolled up do not impair circulation. Always measure blood pressure on the same arm and rest the arm, relaxed on a support.
- The normal cuff is for an **upper arm circumference between 22 and 32 cm**, the larger cuff is for an **upper arm circumference between 32 and 42 cm**. Outside these limits correct readings may no longer be guaranteed.
- When taking repeated measurements, **relax for at least one minute** between measurements as otherwise false

results are obtained.

- In severe **cardiac arrhythmias** (cases of disturbed heart rhythm), measurements should only be taken in consultation with a doctor. Due to the oscillometric measuring method, measurements taken with this device may in some instances be false or fail to produce a reading (Err).
- The monitoring of blood pressure is extremely important in **pregnancy**, as the pressure may be changed by the pregnancy. The results should, however, only be interpreted in consultation with the doctor.

### 3. Control displays and symbols

- ♥ flashes when the device is measuring and the pulse is being taken
- 🔋 change batteries
- Err measuring error, cf. Chap. 12
- ⬆ appears during inflation
- ⬇ appears during automatic checking
- M1 displays the saved reading for person 1

M2 displays the saved reading for person 2

protection from electric shock (type BF)


pay attention to operating instructions

### 4. Current supply

#### Battery operation

Open the battery cover on the underside of the device (see fig.) Insert batteries (see chap. 13, Technical data), ensuring correct polarity (“+” and “-”). Close the battery lid.

- We recommend the use of high-quality batteries as other batteries or accumulators could affect the precision of measurement. Never mix old and new batteries or batteries made by different

manufacturers.

- If you do not intend to use the device for a longer period of time, please remove the batteries.
- In the interests of environmental protection, used batteries should not be disposed of in household waste. Please observe valid regulations when disposing of the batteries or use municipal collecting bins.

#### Operation with mains adapter


There is a connecting socket on the rear side of the device for the mains adapter (output voltage: 6V DC/600mA). To ensure measuring precision, please use only a HARTMANN mains adapter which may be obtained from your pharmacist or specialist medical supplier.

### 5. Setting date and time

- After inserting the batteries the time function automatically appears. The flashing year key is then displayed. The year can be changed using the M1 (+) and M2 (-) buttons. Save the **year** by pressing the START/STOP button.
- Next the **month** is saved. The number on the right will flash. Use the buttons as above for setting the year. Repeat procedure for saving the **day, hour and minute**. Each time you change the batteries the date and time must be reset.

### 6. Applying the cuff

- Measurement should be taken on the **naked arm** which usually has the higher blood pressure value.
- When using the draw-clamp cuff thread the end of the cuff through the metal hoop if necessary, so as to form a loop. In this instance the Velcro fastener must be on the outside. Wrap the cuff round the upper arm, the tube lies in the centre of the elbow angle, pointing towards the hand. Take the free end of the cuff, wrap it firmly round the arm and close the Velcro fastener.
- Wrap the moulded cuff round the upper arm. The air tube of the cuff lies in the centre of the angle of the elbow, pointing towards the hand. Take the free end of the cuff, wrap it firmly round the arm and close the Velcro fastener.
- Check that the white strip is placed in the middle of the elbow, on the artery and that the lower edge of the cuff is approx. 2.5 cm from the bend of the elbow (see fig.).


- The cuff should be firm but not too tight. You should be able to push one or two fingers between the arm and the cuff. Note: Uneven wrapping of the cuff can lead to inaccurate results.
- Use the markings on the edge of the cuff to **check cuff size is correct**. The white arrow should lie within the red marking strip.
- Insert the cuff connector in the **cuff socket on the left side of the device**. NB: do not use the mains connection socket on the back of the device!

## 7. Measuring blood pressure

- It is recommended that blood pressure measurements are carried out whilst seated. Place your arm with the palm upwards on a support and make sure the cuff is at the level of the heart.
- Do not turn the device on until the cuff has been applied, otherwise the cuff

can become damaged through the resulting excess pressure.

- Press the **START/STOP button**. The appearance of all display segments followed by a flashing arrow pointing downwards, shows that the device is checking itself automatically and is ready for use. Subsequently, **automatic inflation** begins at approx. 190 mmHg. If the pump pressure is insufficient or if the measurement is interrupted, the device continues to pump at a rate of 30 mmHg until a high enough pressure is reached.
- If you require a higher inflation pressure, you can avoid having to repeat pumping by pressing the START/STOP button again shortly after inflation starts until the desired cuff pressure is reached. This should be approx. 30 mm Hg over the systolic value.
- **Important: do not move or talk throughout the entire measuring procedure.**
- As pressure in the cuff decreases, the heart symbol and the falling cuff pressure are displayed.
- A **beep** indicates the end of measurement. The **systolic and diastolic readings** appear simultaneously on the display, with the **pulse rate** beneath them (see fig.).


The time appears above the reading and M1 or M2 is displayed on the left. M1 represents the reading for the first person and M2 for the second. As long as the reading is displayed, you can apply the values to the respective person by pressing M1 or M2. Otherwise, if you do not apply them the reading is automatically saved for the person on the display.

- To switch off the device, simply press the START/STOP button. Otherwise the device will switch itself off automatically after 3 minutes.
- If you wish to stop measurement for any reason, simply press the START/STOP button. The inflation or measuring process is interrupted and an automatic fall in pressure occurs.

## 8. Memory function

- The memory recall is activated when the **"Memory Button"** is pressed, when the device is switched off. Press **M1** for the first person's recorded values and **M2** for the second person. The corre-

sponding symbol M1 or M2 will be displayed. First, the **average value** of all saved data will be shown for this person. An "A" appears on the display and the number of readings, from which the average value was calculated is displayed in the top right-hand corner (see fig.).

Press the memory button to access memory position No. 1

- Tensoval comfort can store up to **30 measurements**. The most recent measurement is always in memory position No.1 and all old measurement values shift down one position. When all memory positions are occupied the oldest value will be deleted each time.
- By **repeatedly pressing the memory-button** all stored values can be recalled one by one.
- The measurement and the relevant number of the memory position appear when you access stored values. The memory position, date and time appear in 2 – 3 second intervals.
- You can cancel the memories recall at

any time by pressing the START/STOP button. Otherwise the device will switch itself off after a few seconds.

- Should the power supply fail, e.g. when changing the batteries, the stored values will still be available.

clean the device. Do not use detergents or solvents.

- The cuff can be cleaned carefully with a moist cloth and mild, pH-neutral soap. Do not completely immerse the cuff in water.

### Deleting stored values

You can delete all stored data in M1 or in M2 for the respective person. To do this, press the corresponding memory button. The average value will then appear on the display. Press down on the memory button and after four seconds the display will flash. After eight seconds all data for the respective person will be deleted and only M1 or M2 will be displayed. If you release the button before the end of the eight-second time limit, no data will be deleted.

### 9. Maintenance

- Do not expose the device to extreme temperatures, humidity, dust or direct sunlight.
- This device consists of high quality precision components. Do not drop the device and protect it against hard knocks.
- The moulded cuff should not be folded or overstretched.
- Never open the device. Repairs may only be carried out by qualified individuals.
- Only use a soft, moistened cloth to

### 10. Checking of measurement function

We recommend that measurement function be checked at intervals of two years in the case of professional users, such as pharmacies, medical practices and clinics. You should also observe the national regulations, such as, in Germany, the Decree on Sales of Medical Devices. Measurement function checks can be carried out either by PAUL HARTMANN AG (see Service address), by competent authorities, or authorised maintenance providers against compensation.

#### Instructions for Calibration Mode:




Remove the batteries in order to switch to calibration mode. hold down on the START/STOP button and then insert the batteries. Next, release the button and, after a few moments, two zeros will appear one above each other.

Instructions on checking measurement function will be supplied on request to competent authorities, or authorised maintenance providers.

### 11. Warranty

This product comes with a **3-year warranty**. Details of warranty conditions and contact addresses can be found in the separate warranty document.

### 12. Error displays

Error which has occurred	Possible causes
Device will not switch on	No batteries, they have been inserted incorrectly or are dying. If necessary, check the mains adapter is connected correctly.
	Moving or speaking during the measurement.
	Cuff was incorrectly or loosely applied.
	Please repeat reading.
	Pressure in cuff exceeds 330 mmHg. Therefore, the pressure will be automatically released.
	Batteries are almost dead. Only a few more measurements are possible.
	Batteries are dead. Insert new batteries.

Error which has occurred	Possible causes
Cuff will not inflate	Cuff connection is incorrectly positioned in device.
Implausible reading	Incorrect cuff size. Cuff placed on top of clothing. Rolled up clothing impedes blood circulation. Moving, talking or excitement during measurement. Deep breathing during measurement process. Lack of relaxation before measurement. Stimulants taken just before measurement.

Switch the device off if an error symbol appears. Check the possible causes and note the instructions on self-measurement in Chapter 2. Relax for a minute and take the measurement again. Do not move or talk during the measurement.

### 13. Technical specifications

Measuring method:	oscillometric
Display range:	0 – 300 mmHg
Measuring range:	Systole (SYS): 50 – 250 mmHg Diastole (DIA): 40 – 180 mmHg Pulse: 40 – 160 beats / minute
Technical accuracy of measurement:	cuff pressure: +/- 3 mmHg Pulse: +/- 5 % of displayed pulse frequency
Clinical accuracy of measurement:	complies with the EN 1060 standard, part 3.
Power supply:	4 x 1.5 V alkaline Mangan Mignon (AA/LR06) batteries, or optional HARTMANN mains adapter
Battery capacity:	Tensoval comfort: > 1500 measurements Tensoval comfort large: > 1000 measurements
Inflation pressure:	ca. 190 mmHg
Automatic cut-out:	5 minutes after end of measuring process.
Cuff:	Normal cuff: 22 – 32 cm Large cuff: 32 – 42 cm
Pressure release valve:	electronically and pulse frequency controlled linear valve
Storage capacity:	2 x 30 measurements and average value
Operating conditions:	Ambient temperature: +10 °C to +40 °C Relative humidity: 15 – 90 %
Storage/transport conditions:	Ambient temperature: -20 °C to +50 °C Relative humidity: 15 – 90 %
Serial number:	see battery compartment

**14. Legal requirements and guidelines**

Tensoval comfort complies with the requirements of the EC directive, 93/42/EEC on Medical devices and bears the CE-marking.

The device complies with the EN 1060 standard relating to non-invasive blood pressure monitors Part 1/1995: General requirements and Part 3/1997: Additional requirements for electromechanical blood pressure measuring systems. Clinical testing of measurement accuracy was performed in Germany according to the testing protocol recommended by EN 1060-3 (DIN 58130: Central Hospital Reinkenheide, Prof. Anlauf, Germany, and according to ANSI/AAMI SP10-1992).

Manufacturer:  
PAUL HARTMANN AG, 89522 Heidenheim  
Germany

0 1 2 3

**15. Accessories and spare parts**

To ensure measurement accuracy, only use original HARTMANN accessories which may be obtained from your pharmacist or specialist medical supplier.

- HARTMANN mains adapter  
Code No.: **900 153**
- "Moulded cuff", normal  
Code No.: **900 166**  
Upper arm circumference 22 – 32 cm
- "Draw-clamp cuff" normal  
Code No.: **900 154**  
Upper arm circumference 22 – 32 cm
- "Draw-clamp cuff" large  
Code No.: **900 155**  
Upper arm circumference 32 – 42 cm

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DE – PAUL HARTMANN AG · 89522 Heidenheim
EG – PAUL HARTMANN Egypt S.A.E. · Giza, 12515
ES – Laboratorios UNITEX-HARTMANN S.A. · 08302 Mataró
GB – PAUL HARTMANN Ltd. · Heywood/Lancashire OL10 2TT
HK – PAUL HARTMANN Asia-Pacific Ltd. · Hong Kong
SG – PAUL HARTMANN Pte. Ltd. · Singapore 609 143