Tensoval[®] mobil

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Preliminary remarks

 Please read these instructions carefully before first use as correct blood pressure measurement depends on the appropriate use of the device.

These instructions for use are designed to instruct you, from the very start, in the individual steps of self-measurement of blood pressure using Tensoval mobil. You will thus receive important and helpful hints for producing reliable results for your personal blood pressure profile. Be sure to keep these instructions for use for future reference.

What is Comfort Air Technology?

Many automatic blood pressure devices use preset inflation pressures. Comfort Air Technology, on the other hand, automatically determines a person's individual systolic (upper) blood pressure during inflation and adjusts the inflation pressure by adding 30 mmHg to this individual systolic pressure.

For medical reasons and to ensure measuring accuracy, the inflation pressure used for blood pressure measurement should be about 30 mmHg greater than the expected systolic blood pressure.

What are the benefits of Comfort Air Technology?

A user whose systolic blood pressure is, for example, 130 mmHg and who is using a blood pressure device with a preset inflation pressure of, for example, 190 mmHg, may experience discomfort during blood pressure measurement because the difference between the systolic pressure and the device's inflation pressure is large.

When using Comfort Air Technology, the pressure difference would be only 30 mmHg. Automatically determining this person's individual systolic blood pressure during inflation rather than using a preset inflation pressure, this technology would only inflate to, for example, 160 mmHg.

This means that this blood pressure device only ever uses the inflation pressure that is really needed for accurate and gentle blood pressure measurements without causing user discomfort.

Comfort Air Technology gives you accurate readings and is simple and comfortable to use.

Table of Contents	Page
1. Introduction	160
2. General information on blood pressure	160
2.1 WHO and ISH classification of high blood pressure limits	160
2.2 Importance of self-measurement of blood pressure	161
2.3 Objectives of self-measurement of blood pressure	162
2.4 Regular blood pressure measurement	163
3. Getting ready for self-measurement	163
3.1 Inserting the batteries	163
3.2 Setting date and time	163
3.3 10 golden rules for blood pressure measurement	164
3.4 Correct posture and position	165
3.5 Applying the wrist monitor	166
4. Measuring blood pressure	167
5. Setting the memory function	168
5.1 Saving the measured values	168
5.2 Recalling the measured values	169
5.3 Deleting the measured values	170
5.4 Using the guest mode	170
6. Explanation of error displays	171
7. Significance of symbols on the device and cuff	175
7.1 Control displays	175
7.2 Symbols	175

	Page
8. Important notes	176
8.1 Drugs	176
8.2 Pregnancy	176
8.3 Diabetes, history of other medical conditions	177
8.4 Arrhythmias, heart rhythm disorders, cardiac pacemakers	177
8.5 Important notes for self-measurement	177
9. Maintenance of the device	178
10. Tensoval range of products	178
11. Warranty conditions	178
12. Technical data	180
13. Power supply, disposal notes and safety information	181
13.1 Batteries and disposal	181
13.2 Safety information	182
14. Legal requirements and guidelines	182
15. Instructions for the calibration check	183
16. Contact information for customer queries	184

1. Introduction

Congratulations on your purchase of this HARTMANN quality product. Tensoval mobil is a quality product for fully automatic self-measurement of blood pressure on the wrist. Requiring no preliminary settings, this device conveniently inflates automatically for easy, quick and reliable measurement of the systolic and diastolic blood pressures as well as the pulse rate. Using the HARTMANN Comfort Air Technology, this device automatically determines the maximum inflation pressure to be used for comfortable individualised blood pressure measurements. This device is designed to optimally support

you in taking blood pressure measurements. We wish you all the best for your health.

2. General information on blood pressure

2.1 WHO and ISH classification of high blood pressure limits

The World Health Organisation (WHO) and the International Society of Hypertension (ISH) have developed the following classification for blood pressure values:

World Health Organisation (WHO) blood pressure targets, 1999

Assessment	Systolic pressure	Diastolic pressure
Optimal	up to 120 mmHg	up to 80 mmHg
Normal	up to 130 mmHg	up to 85 mmHg
Normal limit values	130 – 139 mmHg	85 – 89 mmHg
Grade 1 hypertension	140 – 159 mmHg	90 – 99 mmHg
Grade 2 hypertension	160 – 179 mmHg	100 – 109 mmHg
Grade 3 hypertension	over 180 mmHg	over 110 mmHg

To determine your blood pressure you need to measure two values:

- **Systolic (upper) blood pressure:**
Is produced when the heart contracts and pumps blood into the blood vessels.
- **Diastolic (lower) blood pressure:**
This is the value measured when the heart muscle is dilated and again fills with blood.
- **Blood pressure readings are expressed in mmHg.**


Established hypertension (high blood pressure) is defined as repeated measurement of a systolic value greater than 140 mmHg and/or a diastolic value greater than 90 mmHg. Please note that this classification of blood pressure values is independent of age. Optimal blood pressure values have health benefits for all people.

There is no generally recognised definition of too-low blood pressure (hypotension). Readings of less than 100 mmHg systolic and less than 70 mmHg diastolic are considered too low. Please note that, unlike too-high blood pressure values, too-low blood pressure values are not usually expected to

be associated with health risks. However, if you are always feeling unwell, you should check with your doctor.

2.2 Importance of self-measurement of blood pressure

Constantly elevated blood pressure multiplies the risk for other health problems. Chief among these are the physical consequences of a heart attack or stroke including, for example, disability requiring nursing care, paralysis of one side of the body and organic heart or brain damage. In addition to other medical treatment measures, daily blood pressure monitoring will help to protect you from these adverse health outcomes.

2.3 Objectives of self-measurement of blood pressure

Your personal blood pressure profile is an important source of information. In case of drug treatment (e.g. in high blood pressure) your doctor can use your blood pressure profile to tailor your treatment to your particular needs. The better the drug treatment is tailored to your needs, the better you will feel, and the less you will suffer from side effects. Regular,

accurate blood pressure monitoring with Tensoval mobil will help you achieve this goal. Your blood pressure profile thus gives you excellent feedback on how successful your life-style changes are.

i Many people manage to lower their blood pressure through life-style changes (such as e.g. losing weight, dietary modification and getting more


exercise) to levels that do not require drug treatment.

2.4 Regular blood pressure measurement

Numerous factors including physical exertion, taking drugs or the time of day may have an impact on blood pressure. Blood pressure should therefore always be measured at the same time of day under similar conditions.

Our heart may beat up to 100,000 times a day, producing 100,000 different blood pressure values.


3. Getting ready for self-measurement

3.1 Inserting the batteries

Open the battery cover in the direction of the arrow on the upper

side of the device by applying gentle pressure on the housing notches.

Insert batteries (see chap. 12, "Technical data"), ensuring correct polarity ("+" and "-"). Close the battery lid.

3.2 Setting date and time

■ After inserting the batteries, the time function automatically appears. The flashing year is then displayed. The year can be changed using the M1 (+) and M2 (-) buttons. Store the year by pressing the START/STOP button. Next, the month is stored. The number on the right will flash. Use the buttons as above for setting the year. Repeat procedure for saving the day, hour and minute.

■ Each time you change the batteries the date and time must be reset.

3.3 10 golden rules for blood pressure measurement

1. Rest for approx. 5 minutes before measurement. Even deskwork increases blood-pressure by an average of approx. 6 mmHg systolic and 5 mmHg diastolic.

2. Do not consume any coffee or nicotine up to one hour before measurement.

3. Do not measure when you have a strong urge to urinate. A full bladder can lead to an increase in blood pressure of approx. 10 mmHg.

4. Take measurements from the naked wrist and while sitting upright.


5. In the case of using a wrist monitor, hold the cuff at the level of the heart during the measuring procedure. The cuff of an upper arm monitor automatically finds the correct level at the arm.

6. Do not talk or move during the measuring procedure. Talking increases the values by approx. 6 – 7 mmHg.

7. Wait at least one minute between two measurements, so that the vessels are relieved from pressure in preparation for a new measurement.


8. Enter values in the blood pressure diary: Note down the measured values,

together with any drugs taken, the date and time in your blood pressure diary.

9. Take measurements regularly. Even if your values have improved, you

should continue to check them for monitoring purposes.

10. Always take measurements at the same time of day. Because a person has

approx. 100,000 different blood pressure values every day, individual measurements have no significance. Only regular measurements at the same time each day over a long period of time allow a meaningful evaluation of blood pressure values.


Further notes:

- You should take your blood pressure in a quiet place, in a relaxed seated position. Measurement can be taken on the right or left arm. The arm giving higher readings should be used for long-term blood pressure monitoring.
- Do not take your blood pressure after taking a bath or exercising.

3.4 Correct posture and position

- Measurements can be taken on the right or left arm. The arm giving higher readings should be used for long-term blood pressure monitoring.
- The device must be placed at the level of the heart to give an accurate measurement result. To do this, place your arm on the table and bend it up so that the device is at the same level as the heart.


3.5 Applying the wrist monitor

Measurement should be taken on the naked wrist which usually has the higher blood pressure value. If you are unsure which arm to use, ask your doctor. The cuff should not be applied over a projecting bone, as it will otherwise not fit evenly around the wrist.

The device is permanently connected to the cuff, and the cuff should not be removed from the device. Now wrap the cuff round the wrist. The wrist monitor should be positioned axially on the inside of the wrist about 2 cm above the base of the wrist. The lettering on the upper side of the device points in your direction.


The cuff should be firm but not too tight.

Please note that inappropriate application of the cuff can lead to inaccurate readings. Also use the markings on the edge of the cuff to check if the cuff size is correct: The red dot on the silver-coloured strap should point on the red marking strip. If the red dot is outside the red marking strip, the cuff is too small for you. The Tensoval mobil cuff circumference ranges from 12.5 to 22.5 cm. People with larger wrist circumferences are recommended to use an upper arm blood pressure monitor instead because wrist monitors give inaccurate blood pressure readings in such individuals.


4. Measuring blood pressure

- We recommend carrying out blood pressure measurements whilst seated. Do not turn the device on until the cuff has been applied, otherwise the cuff can become damaged through the resulting excess pressure.
- Position the wrist with the strapped-on device and the palm of the hand pointing inwards at the level of the heart.
- Press the blue START/STOP button. The appearance of all display segments followed by a flashing arrow pointing downwards, shows that the device is checking itself automatically and is ready for use.

- Subsequently, automatic inflation begins. This is shown by an arrow pointing upwards and increasing pressure values. Tensoval mobil comes with Comfort Air Technology whereby the inflation pressure is only as high as required, thereby making comfortable measuring possible. If this inflation pressure is insufficient, the device continues to pump until a high enough pressure is reached.
- Once the required inflation pressure has been reached, three short beeps signal that the measurement starts. Then the measuring procedure begins with the release of the air.
- If you require a higher inflation pressure, you can avoid having to repeat pumping by pressing the blue START/STOP button again shortly after inflation starts and holding down on it until the desired cuff pressure is reached. This should be approx. 30 mmHg over the systolic (upper) value.
- **Important:** You should not move or talk throughout the entire measuring procedure!

- If you wish to stop measurement for any reason, simply press the blue START/STOP button. The inflation or measuring procedure is interrupted and an automatic fall in pressure occurs.
- As pressure in the cuff decreases, the heart symbol and the falling cuff pressure are displayed.
- A beep indicates the end of measurement. Then the systolic and diastolic blood pressure values appear simultaneously on the display, with the pulse rate beneath them (see fig.).

- Once the measurement is finished, the time appears above the measured values and M1 or M2 is displayed on the left. M1 represents the measured values

for a first person. Using M2, you can store the measured values for a second person (see 5.1 "Saving the measured values").

- In order to switch off the device, press the blue START/STOP button. Otherwise the device will switch itself off automatically after 3 minutes.

5. Setting the memory function

5.1 Saving the measured values

The device features two memory buttons, M1 and M2, to store readings for two different persons. M1 represents the measured values for a first person and M2 represents the measured values for a second person.

When a beep indicates the end of measurement, you may press M1 or M2 to apply the measured value to the respective person. This assignment can be made as long as the values are displayed. If you do not apply them, the measured value is automatically stored to the displayed memory.

5.2 Recalling the measured values

The device must be switched off to recall data from memory. Press M1 for the first person's stored values and M2 for the second person's values. The corresponding symbols M1 or M2 will be displayed. First, the mean value of all stored values for a person will be displayed. An "A" (for the English designation "Average") appears on the display and the number displayed in the top right-hand corner indicates the number of measurements, from which the average value was calculated (see fig.).

When pressing again, the most recently stored measured value will be displayed. By repeatedly pressing the memory button all stored values of the selected memory can be recalled one by one. The measured value and the number of the memory position are displayed when you access stored values. The number of the memory position, the date and time appear in 2 – 3 second intervals. Tensoval mobil can store up to 60 measured values for each person (M1 or M2). The most recent measured value is always in memory position No.1. When all memory positions are occupied the oldest value will be deleted each time.

! **Important:** The calculated mean value is based on all measured values stored for the respective person. If only two meas-

ured values are in memory, the mean value will be calculated from these two measured values. If only one measured value is in memory, no mean value will be calculated.

You can cancel the memory's recall of data at any time by pressing the START/STOP button. Otherwise the device will switch itself off automatically after approx. 30 seconds. Even if the power supply fails, e.g. when changing the batteries, the stored values will still be available.

5.3 Deleting the measured values

You can delete all stored data either in M1 or in M2. To do this, press the memory button of the corresponding memories M1 or M2. The average value will appear on the display. Hold down the memory button for a longer time. After four seconds, the display will flash and after further four seconds, all data will be deleted. On the display appears now M1 or M2.

Important: If you release the memory button before the end of the time limit, no data will be deleted. If there is only a single

measured value in memory, this stored value cannot be deleted. It is not possible to delete individual values.

5.4 Using the guest mode

When a third person uses Tensoval mobil, it is recommended that they use the guest mode. This one is used to store a measured value neither to the memory positions M1 nor to M2. This avoids bias to the average values and stored values of the two main users of the device. To take measurements in guest mode, start the measuring procedure by simultaneously pressing memory buttons M1 and M2. The START/STOP button need not be pressed.

At the end of measurement the display simultaneously shows both M1 and M2 above the measured values rather than either symbol alone. The measured value can thus not be applied to either person, and the measured values will not be stored.

To switch off the device, press the START/STOP button, also in guest mode. Otherwise the device will switch itself off automatically after 3 minutes.

6. Explanation of error displays

Error which has occurred	Possible causes	Remedy
Device will not turn on.	No batteries, they have been inserted incorrectly or are dead.	Check batteries and insert two identical, new batteries, if and when necessary.
Cuff will not inflate.	Cuff defective.	Return the device to the designated service centre for inspection.
	The measuring signals could not, or not correctly, be read. This can be caused by incorrect application of the cuff, moving, talking or by a very weak pulse.	Check the correct position of the cuff. Do not talk or move during the measuring procedure. Also observe the 10 golden rules.

Error which has occurred	Possible causes	Remedy
	<p>Cuff does not inflate or does not inflate fast enough. This can, among other things, be due to a too loosely applied cuff or to movement.</p>	<p>Apply the cuff so that it fits snugly around the wrist.</p>
	<p>Cuff loses air/has a leak.</p>	<p>Return the device to the designated service address.</p>
	<p>Air release during the measuring procedure is too fast or too slow. The cuff could have become undone or loosened itself. A movement during the measuring procedure is also a possibility.</p>	<p>Check the correct position of the cuff. Do not move during the measuring procedure.</p>
	<p>The pressure in cuff exceeds 300 mmHg. Therefore, an automatic fall in pressure occurs.</p>	<p>Please rest for at least a minute and take the measurement again.</p>
	<p>If the battery symbol flashes, the batteries are almost dead. Only a few more measurements are possible.</p>	<p>Keep new batteries of the same type handy (type AA LR03).</p>

Error which has occurred	Possible causes	Remedy
	If the battery symbol is permanently illuminated the batteries are dead and have to be replaced.	Insert new batteries of the same type (type AA LR03).
Implausible measured values.	Device not positioned at the level of the heart.	Position the wrist at the level of the heart, and repeat the measurement.
	Cuff too big or too small.	Use an upper arm blood pressure monitor.
	Cuff placed on top of clothing.	Apply cuff on the naked skin.
	Rolled-up clothing impedes blood circulation.	Wear loose clothing. Make sure that rolled-up sleeves do not impair circulation in the arm.
	Cuff wrongly applied.	Take note of the instructions and images showing how to apply the cuff correctly to the wrist.
	Cuff was not correctly inflated.	Check the correct position of the cuff.

Error which has occurred	Possible causes	Remedy
Implausible measured values.	Jewellery and/or a bracelet prevent(s) blood from circulating.	Take off your watch and/or any jewellery you may be wearing before measuring your blood pressure.
	Moving, talking or excitement during the measuring procedure.	Please take measurements in a relaxed position whilst seated. Do not talk or move during the measuring procedure.
	Lack of relaxation before taking a measurement.	Relax for 5 minutes before taking a measurement.
	Stimulants taken before measurement.	Please avoid alcohol/ nicotine and caffeine for one hour before taking a measurement.
Memory cannot be deleted.	Only one stored value is in memory. The memory cannot be deleted until two values have been stored.	Save a second value. Then delete the memory.
	The memory can only be deleted when the average value is displayed.	Press the memory button until the average value is displayed.

Switch the device off if an error symbol appears. Check the possible causes and note the 10 golden rules in chapter 3 and the instructions in chapter 8. Relax for a minute and take the measurement again. Do not move or talk during the measurement.

7. Significance of symbols on the device and cuff

7.1 Control displays

 Flashes when the device is measuring and the pulse is being taken.

 Replace batteries if the battery symbol is permanently illuminated.

 Measuring error, cf. chapter 6.

 Appears during inflation.

 Appears during automatic checking.

M1 Displays the stored measured values for person 1.

M2 Displays the stored measured values for person 2.

7.2 Symbols

 Pay attention to the operating instructions.

 Please note.

 Protection from electric shock (type BF).

 PAUL HARTMANN AG
D-89522 Heidenheim
Germany

 **CE**
0 1 2 3


 Electronic device disposal note.

 Disposal note.

8. Important notes

8.1 Drugs

Self-measurement of blood pressure does not replace treatment! So do not interpret your measured values on your own and do not use them for self-prescribed treatment. Take measurements as instructed by your doctor and have confidence in his diagnosis. Take drugs as prescribed by your doctor and never alter the dose on your own. Discuss the appropriate time for self-measurement of blood pressure with your doctor.

8.2 Pregnancy

Blood pressure may change during pregnancy. Regular blood pressure monitoring is particularly important if you have high blood pressure because the elevated blood pressure values may affect the development of the foetus. Check with your doctor whether and, if so, when you should carry out self-measurement of blood pressure.

8.3 Diabetes, history of other medical conditions

If you have diabetes, hepatic disorders or narrowed blood vessels (e.g. arteriosclerosis, peripheral arterial occlusive diseases (PAOD)), you should consult your doctor before carrying out self-measurement because altered measured values may occur in such cases. If you suffer from certain blood diseases (e.g. haemophilia), severely impaired blood flow, or if you take blood-thinning drugs, you should also ask your doctor before carrying out self-measurement.

8.4 Arrhythmias, heart rhythm disorders, cardiac pacemakers

- In the case of severe heart rhythm disorders (arrhythmias), measurements should only be taken in consultation with the doctor. Due to the oscillometric measuring method, in some cases, incorrect readings may be determined or no measurement results are obtained (Err).
- Altered measured values may occur in the case of self-measurement of blood pressure

performed by cardiac pacemaker wearers. The blood pressure monitor itself has no impact on the cardiac pacemaker. Please note that the displayed pulse rate is not suitable for checking the rate of cardiac pacemakers. Please check with your doctor whether self-measurement of blood pressure is advisable if you are wearing a cardiac pacemaker.

8.5 Important notes for self-measurement

- Individual readings are situation-related and thus are not useful.
- Even slight changes in internal and external factors (e.g. deep breathing, stimulants, talking, excitement, climatic factors) lead to fluctuations in blood pressure. This is why your doctor and pharmacist often obtain different readings.
- Always measure blood pressure on the same wrist and hold the blood pressure monitor at the level of the heart.
- Measurement can be taken on the right or left arm. The arm giving higher readings should be

used for long-term blood pressure monitoring.

- Patients who tend to bruise and/or are sensitive to pain on pressure, should only take blood pressure measurements after consulting a doctor.
- The cuff of the Tensoval mobil is for wrist circumferences between 12,5 and 22,5 cm. Outside these limits correct readings may no longer be guaranteed.

9. Maintenance of the device

- Do not expose the device neither to extreme temperatures nor to humidity, dust, or direct sunlight because this may lead to malfunction.
- This device consists of high-quality electronic precision components. Protect the device from knocks and do not immerse in water.
- Only use a soft, moistened cloth to clean the device. Please use neither a diluter nor alcohol, detergents or solvents. The cuff can be cleaned carefully with a lightly moistened cloth and mild, pH-neutral soap. Do not completely immerse the cuff in water

and do not remove it from the device.

- To protect it from external influences keep the blood pressure monitor in the storage box.

10. Tensoval range of products

Apart from this wrist monitor for blood pressure measurement, the HARTMANN range of products also includes upper arm blood pressure monitors. If you would like to learn more about other HARTMANN blood pressure devices, contact your specialist medical supplier (pharmacies or medical supply stores).

11. Warranty conditions

We give a **3-year warranty** on this high-quality device for measuring blood pressure from the day of purchase and in accordance with the following conditions. Claims must be made during the warranty period. The date of purchase may be documented by the appropriately completed and stamped warranty document or proof of purchase. Within the warranty period, HARTMANN shall replace or repair any faulty device components free

of charge which were caused by material or manufacturing errors. This does not extend the warranty period. This warranty is not applicable to damage caused by improper use or unauthorized interference. Parts that are subject to wear and tear (batteries, cuffs etc.) are excluded from the warranty. Claims for compensation are limited to the value of the goods; compensation for subsequent damages is expressly excluded. In warranty cases please send the device with cuff and together with the fully completed and stamped warranty certificate direct, or via your dealer to the Customer Services department for your country.

HK – PAUL HARTMANN
Asia-Pacific Ltd.
Hong Kong

ZA – HARTMANN South Africa
2194 Johannesburg

AE – PAUL HARTMANN
Middle East FZE
Dubai

AU – PAUL HARTMANN Pty. Ltd.
Level 6, 5 Ryder Boulevard
Rhodes, NSW 2138
Australia

GB – PAUL HARTMANN Ltd.
Heywood/Lancashire OL10 2TT

12. Technical data

Measuring method:	Oscillometric
Display range:	0 – 297 mmHg
Measuring range:	Systole (SYS): 50 – 250 mmHg Diastole (DIA): 40 – 180 mmHg Pulse: 40–160 beats/minute The displaying of values outside the measuring range cannot be guaranteed.
Technical measuring precision:	Cuff pressure: +/- 3 mmHg, Pulse: +/- 5 % of displayed pulse rate
Power supply:	2 x 1.5 V Mignon alkaline-manganese (AAA/LR03) batteries
Battery capacity:	> 1,000 measurements
Protection from electric shock:	Medical electrical equipment with an internal power supply (when using batteries)
Applied part:	Type BF
Comfort Air Technology:	Individually determined inflation pressure dependent on systolic blood pressure +30 mmHg.
Protection against harmful penetration of water or solid materials:	IPX0 (not protected)
Operating mode:	Continuous operation
Inflation pressure:	At least 150 mmHg

Automatic switch-off function:	3 minutes after end of measurement
Cuff:	12.5 – 22.5 cm
Pressure release valve:	Electronically controlled linear valve
Memory capacity:	2 x 60 measurements and mean value
Operating conditions:	Ambient temperature: +10 °C to +40 °C (+ 50 °F to + 104 °F) Relative humidity: 15 – 90 %
Storage/transport conditions:	Ambient temperature: -20 °C to +50 °C (- 4 °F to +122 °F) Relative humidity: 15 – 90 %
Serial number:	See battery compartment

13. Power supply, disposal notes and safety information

13.1 Batteries and disposal

- We recommend the use of high-quality batteries as other batteries or accumulators may result in a reduction in the measuring performance. Never mix old and new batteries or batteries made by different manufacturers.
- If the battery symbol is permanently illuminated you should

change the batteries. Please note that the symbol always looks “empty”.

- Remove the batteries from the device if it is not being used for a longer period.
- In the interests of environmental protection, exhausted batteries may not be disposed of in household waste.  Please observe the applicable waste disposal regulations or use public collecting bins.

13.2 Safety information

- The device is not waterproof!
- Do not leave the device unattended near toddlers or persons who cannot operate it themselves.
- Use the device for taking blood pressure measurements on the wrist only.
- Do not under any circumstances carry out blood pressure measurements on babies or toddlers.
- Do not expose the device to hard knocks or vibrations.
- Do not drop the device to the floor.
- The device must not be altered, dismantled, or repaired by the user.
- Do not excessively bend or fold the cuff.
- Never inflate the cuff when it is not properly applied to the wrist.
- Please do not apply the cuff over a wound, as this may result in further injuries.
- If you have had a mastectomy, do not carry out the measurement on the arm on the affected side of the body.
- Please note that the pressure built up by the cuff can lead to

temporary disruption to medical devices being simultaneously used on the same arm.

- If an intravenous treatment is being carried out or a venous catheter is present on the arm, blood pressure measurements can lead to injury. Never use the cuff on the arm on which these conditions apply.
- Please wait for one minute between two measurements.
- If you are carrying out the measurement on another person, please ensure that the use of the device does not result in persistent impairment of the blood circulation.

14. Legal requirements and guidelines

Tensoval mobil complies with the requirements of the EC directive 93/42/EEC on medical devices (Medical Device Directive MDD) and bears the CE mark.

The device complies, for example, with the European Standard EN 1060: Non-invasive blood pressure measuring devices, Part 1: General requirements and Part 3: Additional requirements for electro-mechanical

blood pressure measuring systems. Clinical testing of measuring accuracy was performed according to the European Standard EN 1060-4. Portable and mobile high-frequency and communication devices, such as telephone and mobile phone, can impair the functional capability of electronic medical devices. In compliance with European Standard EN 60601-1-2 further information can be supplied by HARTMANN on request.

15. Instructions for the calibration check

Each Tensoval mobil device has been carefully tested by HARTMANN for measuring precision, and been developed for a long useable service life. We recommend a calibration check at intervals of two years in the case of professionally used devices, for example in pharmacies, medical practices, or hospitals. You should also observe the national regulations determined by the legislator, such as, in Germany, the "Medizinprodukte-Betreiberverordnung" (Medical Device Operating Regulation). Calibration checks can

be carried out either by competent authorities or authorised maintenance providers against compensation.

Instructions for the calibration check:

A device function check can be carried out on people or using a suitable simulator. Calibration check involves testing for leak tightness of the pressure system and possible deviations of the pressure reading. Remove at least one battery in order to switch to calibration mode. Hold down on the START/STOP button and then insert the battery. Release the button after a few seconds and after a short moment, two zeros will appear one above each other on the display. Instructions on the calibration check will be supplied on request to competent authorities or authorised maintenance providers by HARTMANN.

16. Contact information for customer queries

AE – PAUL HARTMANN
Middle East FZE
Dubai

AU – PAUL HARTMANN Pty. Ltd.
Level 6, 5 Ryder Boulevard
Rhodes, NSW 2138
Australia

GB – PAUL HARTMANN Ltd.
Heywood/Lancashire OL10 2TT

HK – PAUL HARTMANN
Asia-Pacific Ltd.
Hong Kong

ZA – HARTMANN South Africa
2194 Johannesburg

Date of revision of the text:
2013-06