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ATTENTION

For greater safety:

Read and understand all the instructions contained in these instructions for use before installing or operating this equipment.

Note: These instructions for use must be read by all the operators of this equipment.

INDEX

PRESENTATION OF MANUAL	27
IDENTIFICATION OF THE EQUIPMENT	29
MODULES, ACCESSORIES, OPTIONALS AND CONSUMPTION MATERIALS.....	31
TECHNICAL SPECIFICATIONS	33
- Standards applied	34
- Content of accessible and non-accessible demarcations	34
- Electromagnetic emissions.....	34
- Packing symbols.....	39
- Product symbols.....	39
EQUIPMENT INSTALLATION	40
EQUIPMENT OPERATION	41
PRECAUTIONS, RESTRICTIONS AND WARNINGS	45
- Transport conditions, warehousing and operation	45
- Sensitiveness to environmental conditions foreseeable in normal situations of use	45
- Precautions and warnings "during the installation" of equipment.....	45
- Recommendations for the dental equipment maintenance.....	46
- Precautions and warnings "during the utilization" of the equipment	46
- Precautions and warnings "after" the use of equipment	47
- Precautions and warnings during the "cleaning and disinfection" of equipment.....	47
- Precautions in case of alteration in the functioning of equipment.....	47
- Recautions to be adopted against foreseeable or uncommon risks, related to the deactivation and abandoning of equipment.....	47
CORRECTIVE AND PREVENTIVE MAINTENANCE AND PRESERVATION.....	48
- Additional procedures for reuse	48
- Cleaning	48
- Disinfection	48
- Sterilization	49
- Preventive Maintenance	49
- Corrective Maintenance	49
UNFORESEEN EVENTS – SOLUTION OF PROBLEMS	50
WARRANTY OF EQUIPMENT	50
FINAL CONSIDERATIONS.....	50

IDENTIFICATION OF THE EQUIPMENT

Dear Customer

Congratulations. You have made a good choice when you decided to buy a GNATUS QUALITY product comparable to the best products available in the World. This manual is a general presentation of your product and it will give you important details to help you to solve possible problems.

Please, read it and keep this with you.

Equipment indication

The indications regarding the use of the Endus Apex locator are below:

- In routine situations of the dental treatment;
- For the detection of perforations, fractures and root reabsorptions;
- For the follow-up of the amount of work during the cleaning and modeling process of curve roots (dynamic odontometrics);
- Pregnant patients and patients that presented nausea during the radiographic measurements.

This equipment is for exclusive odontological use, and it should be handled and used by an expert professional (professional properly certified, according to local legislation of the country) observing the instructions contained in this manual.

It's an user obligation to use the equipment in perfect working conditions and to protect himself/herself, the patients and third parties against eventual dangers.

Physics Principle used by the equipment

It is well understood that the tooth works like a capacitor with accumulation of electric charges at the periodontal and on the inside of the root canal. The dentin works like an insulating for the propagation of electric current throughout all the extension of the root canal.

The apex locators work based on the principle of constancy of the electric current between the oral mucosa and the periodontal ligament. The electronic method takes into consideration the difference of electric conductivity of a metallic instrument on the inside of the root canal and the conductivity of the periapical tissue. The electric current existing on the root canal would complete the circuit at the moment that the electrode, file, touched on the tissue fluid, indicating the most apical portion of the root canal "the apical foramen".

IDENTIFICATION OF THE EQUIPMENT

Equipment Description


The Apex locator Endus is a high precision electronic equipment, modern and that should be used in endodontics to help with canal treatments.

It's used to measure the precise tooth length - (It determines the length of work) – up to where the instrumentation and filling of a tooth will take place, so that the endodontic treatment respects the biological limits of the dental structure.

Main aspects:

- Indicates precisely the cement dental junction.
 - It is not subject to interference from the anatomic structures adjacent to the tooth.
- Therefore, the fact of the canal being dry or full with electrolytes, blood and physiological brine does not affects the result of the measurement;
- Less time for obtaining the length of work;
 - It has a frontal LCD display with the following symbologies:
 - Tooth length meter;
 - Indicator when the file is placed on the canal of the dental root;
 - Sound indicator;
 - Energy consumption indicator;
 - It operates automatically, not needing manual adjustments;
 - Equipment functions with batteries AAA – 1,5V.


MODULES, ACCESSORIES, OPTIONALS AND CONSUMPTION MATERIALS

- 01 - On/Off Key
- 02 - Sound adjustment Key
- 03 - Key for the apex depth adjustment (deeper)
- 04 - Key for the apex depth adjustment (less deep)
- 05 - Indicator for the power level
- 06 - Sound level indicator: Mute, low, medium and high.
- 07 - Visualization of the file position
- 08 - Connection of the measurement cable
- 09 - Indicator of the level of advancement of the canal
- 10 - Meter of the depth of advancement of the canal
- 11 - APEX (ideal level of the root apex)
- 12 - OVER (apex exceeding area)

MODULES, ACCESSORIES, OPTIONALS AND CONSUMPTION MATERIALS

Accessories that accompany the product

- 01 - File support (02 units)
- 02 - Cable for measurement
- 03 - Stainless steel labial Clip (04 units)
- 04 - Batteries AAA 1,5V (05 units)

The use of any part, accessory or material neither specified nor foreseen in those use operation instructions and it is of the user's entire responsibility.

TECHNICAL SPECIFICATIONS

Technical characteristics of the equipment and its accessories:

Model
Endus
Classification of Equipment as per ANVISA
Class II
Classification of Equipment as per standard IEC 60601-1:
Protection against Electric Shock - Type B and Class II Equipment (IEC 60601-1)
Protection against harmful water penetration
IPX 0
Operation mode
Continuous operation
Batteries
AAA 1.5V (05 units)
Voltage
DC7,5V
Power consumption
<0.05W
Current
<10.μ A
Display mode
Display LCD 5"
Net weight
0,365 kg
Gross weight
0,787 kg
Time sounds
The device has 4 levels of sounds: 2,0mm to 1,0mm paused beeps / 1,0mm to 0,5 continuous / 0,5mm to 0,0mm short and frequent beeps and below 0.0mm (OVER) intermittent beeps.

TECHNICAL SPECIFICATIONS

The materials used to produce the equipment are Biocompatible.

Use of different cables, transducers and accessories from those specified may result in increased emissions or decreased immunity of the equipment.

Standards applied:

- IEC 60601-1:2005 - Equipamento Eletromédico - Parte 1: Prescrições gerais para segurança;
- EN ISO 980:2008 (Ed. 2) - Graphical symbols for use in the labelling of medical devices;
- EN ISO 14971:2012 - Medical devices - application of risk management medical devices;
- EN ISO 13485-2012 - Quality systems - medical devices;
- IEC 60601-1-2:2007 - Compatibilidade Eletromagnética;
- ISO 10993-1:2003 - Biogical evaluction of medical devices.

Content of accessible and non-accessible demarcations

FABRICANTE E RESPONSÁVEL PELA GARANTIA: MANUFACTURER AND RESPONSIBLE FOR THE WARRANTY: FABRICANTE Y RESPONSABLE POR LA GARANTIA:

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<0.05W

Electromagnetic emissions

Eletromagnetic emissions		
The equipment is made to be used in the electromagnetic environments specified below. The client or the user of the equipment must be sure that it is used in such environment.		
Emission test	Compliance	Eletromagnetic environment - Guide
RF emissions ABNT NBR IEC CISPR 11	Group 1	This equipment uses RF energy only for internal functions. However, its emissions are too low and it's unlikely to cause any interference in the equipments next to it.
RF emissions ABNT NBR IEC CISPR 11	Class B	This equipment is proper to be used in all establishments; including domestic settings and those directly connect to a public low voltage distribution which feeds domestic buildings.
Emissions of harmonics IEC 61000-3-2	Class A	
Fluctuation of Voltage / Emissions of flicker IEC 61000-3-3	As per	

34

TECHNICAL SPECIFICATIONS

Electromagnetic emissions

Recommended distances between portable and mobile RF communication equipments and the equipment			
The equipment is made to be used in an electromagnetic environment in which RF disturbances are controlled. The client or the user of the equipment may help preventing electromagnetic interference by keeping a minimal distance between mobile and portable RF communication equipment (transmitters) and the equipment, as recommended below, in accordance with the maximal voltage output of the communication equipment.			
Transmitter Maximum Output (W)	Separation distance according to transmitter frequency (M)		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 kHz to 800 MHz $d = 1,2\sqrt{P}$	800 kHz to 2,5 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters with a maximum nominal output power not listed above, the recommended d separation distance in meters (M) can be determined using an equation applicable to the frequency of the transmitter, where P is the transmitter maximum nominal output in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, is applied the separation distance for the higher frequency range.			
NOTE 2 These guidelines may not apply to all situations. The absorption and reflection from structures, objects and people affect the electromagnetic propagation.			

TECHNICAL SPECIFICATIONS

Electromagnetic emissions

Guidelines and manufacturer's declaration - electromagnetic immunity			
The equipment is made to be used in the electromagnetic environments specified below. The client or the user of the equipment must be sure that it is used in such environment.			
Immunity test	ABNT Test level NBR IEC 60601	Level of compliance	Electromagnetic environment Directives
Electrostatic discharge(ESD) IEC 6100-4-2	± 6 kV Contact ± 8 kV Air	± 6 kV Contact ± 8 kV Air	Floors should be wooden, concrete or ceramic. If the floor is covered with synthetic material, the relative humidity should be at least 30%
Quick electric transitory phases / train of pulses ("Burst") IEC 61000-4-4	± 2 kV in power supply lines ± 1 kV in input / output lines	± 2 kV in power supply lines ± 1 kV in input / output lines	It is advisable that the quality of the power supply should be that of hospital or typical commercial environment
Surges IEC 61000-4-5	± 1 kV lines (s) to lines (s) ± 2kV lines (s) to ground	± 1 kV lines (s) to lines (s) ± 2kV lines (s) to ground	It is advisable that the quality of the power supply should be that of hospital or typical commercial environment
Reduction, interruption and variance of voltage in power supply input lines IEC 61000-4-11	< 5% U_t (>95% drop in U_t) for 0,5 cycle 40% U_t (60% drop in U_t) for 5 cycles 70% U_t (30% drop in U_t) for 25 cycles < 5% U_t (>95% drop in U_t) for 5s	< 5% U_t (>95% drop in U_t) for 0,5 cycles 40% U_t (60% drop in U_t) for 5 cycles 70% U_t (30% drop in U_t) for 25 cycles < 5% U_t (>95% drop in U_t) for 5s	The recommended power supply quality is the same as used for commercial or hospital environment. If is required a continuous use during energy supply outages, it is recommended that the equipment be feed by an uninterruptible power supply or a battery.
Magnetic field in frequency of power supply (50/60Hz) IEC 61000-4-8	3 A/m	0,3 A/m	If an image distortion occurs, may be necessary place the equipment far from the supply frequency or to install a magnetic armour. The frequency magnetic field shall be measured at the installment place to assure that it is low enough.
NOTE U_t is the a.c. power supply voltage before the application of the test level			


TECHNICAL SPECIFICATIONS

Electromagnetic emissions

Guidelines and manufacturer's declaration - electromagnetic immunity			
The equipment is made to be used in the electromagnetic environments specified below. The client or the user of the equipment must be sure that it is used in such environment.			
Immunity test	ABNT test level NBR IEC 60601	Level of compliance	Electromagnetic Environment Directives
RF conducted IEC 61000-4-6	3 vrms 150 kHz up to 80 MHz	3 Vrms	It is advisable that portable and mobile RF communication equipment is not used near any part of the equipment, including cables, with a separation distance less than the one recommended, calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P} \text{ 80 MHz thru 800MHz}$ $d = 2,3\sqrt{P} \text{ 800 MHz thru 2,5MHz}$ Where P is the nominal maximum power of output of the transmitter in watts (W), as per the manufacturer of the transmitter, and d is the recommended separation distance in meters (m). It is advisable that the fiel intensity from the RF, transmitter as determined by means of electric inspection on-site, ^a is less than the level of compliance in each frequency range ^b. There may be interference near the equipment marked with the following symbol: 
RF radiated IEC 61000-4-3	3 V/m 88 MHz up to 2,5 GHz	3 V/m	
NOTE 1 At 80MHz and 800MHz, the highest frequency range applies. NOTE 2 These directives may not be applicable in every situation. The electromagnetic transmission is affected by the absorption and reflection of structures, objects and people. ^a The field intensities set by the fixed transmitters, such as radio base stations, telephones (mobile phone, wireless) land mobile radio, amateur radio, AM and FM radio transmissions and TV transmissions can not be predicted with accuracy. Due to the RF fixed transmitters is recommended to install an electromagnetic inspection at the local in order to evaluate the electromagnetic environment. If at the place where the equipment is be using the field intensity level exceeds the conformity level for the RF above, is recommended to observe if the operations are normal. Whether abnormal operations are observed, additional procedures shall be necessary such as reorientation or replace the equipment. ^b Whether above the frequency range of 150kHz to 80 MHz is recommended a field intensity below than 3 V/m.			

TECHNICAL SPECIFICATIONS

Dimensional (mm)

TECHNICAL SPECIFICATIONS

Packing symbols

It determines the maximum quantity of boxes which can be stacked during transportation and storage "as per packaging".

Packing to be transported and / or stored avoiding humidity, rains and wet floor.


Packing to be transported and / or stored with the harrows up.

The packing must be stored and transported away from direct sun light exposure.

Packing to be transported and / or stored with care (should not suffer drop and neither receive impact).

Temperature limit for the packing to be stored or transported.

Product symbols

On/Off Button

Apex depth adjustment key (deeper)

Apex depth adjustment key (less deep)

Sound level indicator

Battery level indicator

B type equipment

Careful : It indicates an important instruction for the operation of the product. Not following it can cause dangerous malfunctioning.

Note: It indicates useful information for operation of the product.

Important: It indicates an instruction of safety for operation of the product. Not following it, can lead to serious danger to the patient.


Warning - see the manual

Authorized representative in the European Community

EQUIPMENT INSTALLATION

Fig.01


To use the equipment you should proceed with the installation of the batteries. The batteries compartment is located on the back of the equipment (Fig.01).

Remove the back cover (A) and insert the 5 batteries AAA 1,5V (B).

Plug in the measurement cable (C) on the equipment, connect the other extremities of the cable to the file support (D) and to the labial clip (E), according to (Fig.02).

Warning: Check that the plug is firmly connected to the inlet. The wrong connection can interfere with the correct measurement.


Turn on the equipment (F) to certify that the installation was conducted correctly.

Fig.02

ALWAYS HOLD THE CONNECTORS TO CONNECT OR DISCONNECT THE CABLES.

EQUIPMENT OPERATION


After having installed the equipment correctly, place the labial clip on the lip of the patient and hold the file on its support. Introduce it in the canal, according to the picture below.

ENGLISH

Placement of the file on the support


Press the thumb in the direction of the arrow (A), hold the file (B) and release it (C).


Turn on the equipment by pressing the on/off button. To obtain precise measurements the instrument should be compatible with the dimension of the canal (diameter of the instrument and canal width). The equipment automatically detects the beginning of the measurements and tests the quality of the electric contact and the conductivity of the radicular canal.

EQUIPMENT OPERATION


Functioning of the LCD Display


1. Indicator of the level of canal advancement
2. Visualization of the file position.
3. Indicator of the sound level: Mute, low, medium and high.
4. Power level indicator.
5. Meter of the depth of canal advancement
6. APEX (ideal level of the root apex)
7. OVER (apex exceeding area)


0.C: Indicates that no measurement is being processed.


2.0mm: Beginning of the canal easurement where beeps will sound slowly from 2,0mm to 1,0mm.

1.0mm: Advancement of the canal measurement, where beeps will sound continuously from 1,0mm to 0,5 mm.

0.5mm: Indicates that the extremity of the file is near the apical region (APEX), where short and frequent beeps will sound 0,5mm to 0,0mm.

0.0mm: When the apical foramen is reached or outreached (APEX), intermittent beeps will sound and the word (OVER) will light and blink.

EQUIPMENT OPERATION

The equipment has 4 levels of sound adjustment : mute, low, medium and high, if you don't hear beeps, clean the labial clip and the instrument.

On the LCD Display the levels of the scale of the canal length are divided in the colors yellow, red and blue. As soon as the area is visible on the canal apex the yellow color will blink, alerting the beginning of the measurement.

The scale is divided in 10 levels, where each level represents 0,1 mm of the apex area. The spare yellows represent the narrow area of the root canal apex where they are defined in five levels and can be adjusted by the functions button, if needed.

It's possible to adjust the 0,1mm to 0,5mm scale through the control of the length of the narrow area. Generally the narrow area is 0,2-0,3mm distant to the canal apex.

When initiating the procedure the equipment will detect the beginning of the measurement and will light the drawing of the file on the LCD Display.

Push the file, rotating it slowly clockwise until the word "APEX" starts to blink. When it reaches the apex, rotate the file slowly counterclockwise until the meter reaches 0.5mm again. When the scale of the length of the narrow area on the LCD Display is on 0.5mm the file will have reached the desired position, i.e., at the root apex.



Position the rubber stop on the surface of the tooth as a reference point to determine the length of work of the root canal.

Reminding that as soon as the file reaches the 2,0 mm position from the canal apex, beeps will sound slowly from 2,0 to 1,0mm.


The duration of the beeps will raise according to the advancement, allowing the precise location of the tip of the instrument on the apical area. Automatically it will change to continuous range from 1,0mm to 0,5mm. After the 0,5mm scale the beeps will be short and frequent reaching the "APEX" ideal level, where this word will keep on blinking on the Display, when the tip of the instrument reaches the foramen.

CAUTION: If the instrument outreaches the apex, the word "OVER" (which represents the exceeding area of the apex) will keep on blinking followed by intermittent beeps.

The labial clip can be disconnected from the instrument and connected again during the measurement cycle, without affecting its functioning. It can be changed by one or another more adequate during the root canal treatment or when another canal is measured. In these cases the equipment automatically detects that a new measurement has been initiated.

ENGLISH

EQUIPMENT OPERATION

THE LENGTH OF WORK ON THE CANAL VARIES FROM CASE TO CASE DUE TO DIFFERENT FORMS OF THE TEETH AND OF THE CANALS OF THE DENTAL ROOTS.

AFTER THE WORKING LENGTH HAVING BEEN DEFINED TURN THE EQUIPMENT OFF AND DISCONNECT ALL THE CABLES.

YOU SHOULDN'T CONTINUE WITH THE MEASUREMENTS IF YOU DON'T HEAR THE BEEPS.

WHEN THE EQUIPMENT IS TURNED ON, THE SOUND LEVEL PREVIOUSLY SELECTED IS ACTIVATED. ACTIVATE THE SOUND ADJUSTMENT KEY TO CHANGE THE VOLUME OR TURN THE SOUND OFF.

AUTOMATIC TURNING OFF: THE EQUIPMENT WILL AUTOMATICALLY TURN OFF WHEN IT IS NOT BEING USED FOR MORE THAN 3 MINUTES.

PRECAUTIONS, RESTRICTIONS AND WARNINGS

Transport conditions, warehousing and operation

The equipment must be transported and stored with the following observations:

- Carefully, to not suffer falls and neither receive impacts.
- Protected from moist, do not expose to rain, water splashes or moist floor.
- Keep in a protected place from the rain and direct sun and its original package.
- When transporting, do not move it in irregular surfaces and protect the package from the direct rain and respect the maximum piling informed at the external part of the package.
- Room temperature range for transport or warehousing -12°C to $+50^{\circ}\text{C}$.
- Room temperature range recommended by Gnatus $+10^{\circ}\text{C}$ to $+35^{\circ}\text{C}$.

ATTENTION: The equipment must be only used with the power cable supplied by Gnatus, otherwise, it can cause the increasing of the emission or reducing of the equipment's immunity.

The equipment keeps its safety condition and efficiency, since kept (stored) as mentioned in this operation instruction. That way, the equipment will not lose or change its physical and dimensional characteristics.

Sensitivity to predictable environmental conditions in ordinary situations of use

- The equipment has been designed not to be sensible to interferences as magnetic fields, external electrical influences, electrostatic discharges, the pressure or the variation of pressure, since the equipment is installed and kept clean, conserved, transported and operated conforming this operation instruction.

Precautions and warnings “during the installation” of the equipment

- The equipment should only be installed by Gnatus authorized technical assistance or technicians.

- Position the equipment in a place where it will not get wet.
- Install the equipment in a place where it will not be damaged by pressure, temperature, humidity, direct sunlight, dust, salts, or sulfur compounds.
- The equipment should not be submitted to inclination, excessive vibrations, or blows (including during transportation and handling).
- This equipment was not planned for use in an environment where vapors, anesthetic mixtures inflammable with air, or oxygen and nitrous oxide can be detected.
- Before the first use and/or after long interruptions from work such as vacations, clean and disinfect the equipment; eliminate air and water deposited in the internal hoses.- Insert the batteries in the equipment before using it. Check that the power of the batteries is enough to assure that the result of the measurement is correct. When you are going to change the batteries, do not use new ones with old ones and don't mix alkaline batteries with manganese batteries. In case the equipment is not to be used for a long period, the batteries shall be removed, and this should apply in case there is a long time transportation of the equipment.

- Check that the measurement plug is firmly connected to the inlet. The inadequate connection can hinder the measurement.

PRECAUTIONS, RESTRICTIONS AND WARNINGS

These information also make part of the Manual of Installation and Maintenance of the equipment that can be found with the authorized Gnatus technician.

Recommendations for the dental equipment maintenance

Your Gnatus equipment has been designed and developed according to the standards of modern technology. Similarly to other kinds of equipment, it requires special care, which is many times neglected due to several reasons and circumstances. Therefore, here are some important reminders for your daily routine. Try to follow these simple rules, which will save you a lot of time and will avoid unnecessary expenses once they start making part of your working procedure.

Precautions and warnings “during the utilization” of the equipment

- The equipment should only be operated by duly enabled and trained technicians (Dental Surgeons, Capacitated Professionals).
- If any maintenance should be required, only use services of the Gnatus Authorized Technical Assistance.
 - The equipment has been manufactured to handle both continuous.
 - Although this equipment has been planned in accordance with the standards of electromagnetic compatibility, it can, in very extreme conditions, cause interference with other equipment. Do not use this equipment together with other devices very sensitive to interference or with devices which create high electromagnetic disturbance.
 - When the file accidentally touches the internal part of the canal of the dental root, the length meter of canal advancement will become a bit abnormal, but, after a few seconds, it will automatically go back to normal.
 - Do not use damaged file supports; a precise measurement cannot be made with a damaged file-holder.
 - Never use the equipment if the power indicator of the batteries is blinking. It'll not work properly if the batteries are worn out.
 - Only use files and file support with plastic cables. If the file has a metallic cable, there will be a current leakage when the fingers touch the cable, which will impede a precise measurement of the root canal. Even if the file cable is made of plastic, check to not touch the metallic part with the fingers.
 - Never use oxidized, deformed, discolored or with an abnormal aspect batteries.
 - Always hold the connectors to connect or disconnect the cables.
 - Replace the batteries as soon as the power level indicator begins to blink.

The equipment is not indicated for patients and dentist-surgeons bearers of a cardiac pacemaker.

Gnatus shall not be responsible for:

- Use of the equipment differing from that for which it is intended.
- Damages caused to the equipment, the professional and/or the patient by the incorrect installation and erroneous procedures of maintenance, differing from those described in these Instructions for use which come with the equipment or by the incorrect operation of it.

PRECAUTIONS, RESTRICTIONS AND WARNINGS

Precautions and warnings “after” the utilization of the equipment

- Turn off the main switch of the dental set when it is not in use for an extended period of time.
- Always maintain the equipment clean for the next operation.
- Do not modify any part of the equipment. Do not disconnect the cable or other connections without need.
- After using the equipment, clean and disinfect all the parts which may be in contact with the patient.
- The batteries shall not be kept inside the equipment in case it is not to be used for a long time. When you are not using the equipment, we advise you to remove the batteries, since they are corrosive and acid leakages can, with time, damage the equipment.

Precautions and warnings during the “cleaning and disinfection” of the equipment

- Avoid spilling water, even accidentally, or other liquids inside the equipment, which could cause short circuits.
- Do not use microabrasive material or steel wool when cleaning, or employ organic solvents or detergents which contain solvents such as ether, stain remover, etc.
- For the cleaning of the LCD Display use a soft and dry cloth.
- Do not sterilize other than in an autoclave (humid vapor), the parts appropriate to this process are: Labial clip and file support.

Precautions in case of alteration in the functioning of equipment

- If the equipment has any abnormality, check if the problem is related to any item listed in the topic of unforeseen events (failures, causes and solutions). If it is not possible to resolve the problem, turn off the equipment, remove the power supply cable from the socket and contact your representative (Gnatus).

Precautions to be adopted against foreseeable or uncommon risks, related to the deactivation and abandoning of equipment

In order to avoid environmental contamination or undue use of the Equipment after it has become useless, it should be discarded in the suitable place (as per the local legislation of the country).

- Pay attention to the local legislation of the country for the conditions of installation and disposal of residue.

CONSERVATION AND CORRECTIVE AND PREVENTIVE MAINTENANCE

Additional procedures for reuse

The equipment can be reused in undetermined, i.e. unlimited, quantities, only needing to be cleaned and disinfected.

Cleaning

Important: In order to execute cleaning or any type of maintenance, ensure that the equipment is disconnected from the electrical network.

The cleaning procedure should be executed at the start of the working day and after each patient. Always turn the main switch off before executing the procedures of daily maintenance.

• BODY OF THE EQUIPMENT, CABLE, TRANSDUCTOR AND HOSES

To clean the equipment, we recommend the use of **"BactSpray"** (Reg n° MS: 3.2079.0041.001-5) or any other similar product: Active component:

Benzalkonium chloride (tri-quaternary ammonium)
Solution 50%..... 0.329%

Chemical composition: Butyl Glycol, Decyl polyglucose, Sodium Benzoate, Sodium Nitrate, Essence, Deodorized Propane / Butane, demineralized Water.

For more information concerning cleaning procedures, see manufacturer's instructions.

NOTE: The registration at the Ministry of Health of the "BactSpray" is executed separately from the product described in this manual, as the "BactSpray" is not manufactured by Gnatus.

WARNING:

- In order to prevent risks and damages to equipment, make sure that the liquid does not enter into the equipment.
- The application of other solvent-based cleaning products or sodium hypochloride isn't recommended, because they may damage the equipment.

Disinfection

- Use clean and soft cloth dampened in alcohol 70% to disinfection of the equipment.
- Never use corrosive disinfectants or solvents.

Note: Use gloves and other systems of protection, during the disinfection.

CONSERVATION AND CORRECTIVE AND PREVENTIVE MAINTENANCE

Sterilization

The accessories that get in touch with the patient's mouth during the use of the equipment (labial clip and file support) should be sterilized in autoclave, between two patients in order to avoid cross contamination.

THIS EQUIPMENT IS NOT PROVIDED STERILIZED, IT SHALL BE CLEANED AND STERILIZED BEFORE ITS FIRST USE.

Preventive Maintenance

The equipment must suffer routinely measurements, following the current legislation of the country.

But, never with a period superior to 3 years.

For protecting your equipment, look for a Gnatus' technical assistance for periodic reviews as preventive maintenances.

Corrective Maintenance

Gnatus states that the supplying of the circuits' diagram, Part lists or any other information that permits the technical assistance by the user, can be requested, since previously agreed between the buyer and Gnatus.

In case of the equipment presents any abnormality; check if the problem is related to some of the listed items under the item Unpredictable (situation, cause and solution). If it's not possible to solve the problem, shutdown the equipment and call Gnatus' technical assistance.

UNFORESEEN EVENTS – SOLUTION OF PROBLEMS

⚠ Upon coming across any problem in operation, follow the instructions below to check and repair the problem, and/or get in touch with your representative.

Unexpectedness	Probable Cause	Solutions
- Equipment non-operating.	- Check if the batteries are perfectly positioned. - Check if the batteries are charged.	- Correctly position the batteries. - Replace the batteries.
- The alarm does not work.	- Check the sound volume. - The file might not have reached the point where the instrument activates the alarm.	- Adjust the volume deactivating the option "mute". - Reach the canal apex.
- The LCD is inoperative not presenting alterations to its part.	- First turn on the equipment and then manually connect the file.	- Correctly install the accessories.
- Irregular sound indicators.	- Involuntary contact of the file with the tooth walls.	- Restart the process.

EQUIPMENT'S WARRANTY

This equipment is covered by the warranty terms and norms contained in the Warranty Certificate that accompany the product.

FINAL CONSIDERATIONS

Among the care you have to take with your equipment, the most important is regarding of the spare parts replacement.

To ensure the lifetime of your device, only replace original spare parts from Gnatus. They have the assurance of the standards and technical specifications required by the Gnatus representative.

We call your attention to our authorized resellers' chain. Only this chain will keep your equipment constantly new, because it has trained technical assistant and specific tools for the correct maintenance of your device.

Whenever you need, demand the presence of a Gnatus' technician from the nearest resale, or ask through the Attendance Service GNATUS: + 55 (16) 2102-5000 / SAC: 0800-7015-054.

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Obelis S.A, Boulevard Général Wahis 53, 1030 Brussels, Belgium,