N Latex ASL

Intended Use

In vitro diagnostic reagents for the quantitative determination of Anti-Streptolysin O (ASL) in human serum by means of particle-enhanced immunonephelometry using the BN* Systems.

Summary and Explanation

Immunochemical assay of specific antibodies to streptococcal metabolites provides valuable information regarding past streptococcal infection (rheumatic fever, scarlet fever, tonsillitis, glomerulonephritis, etc.)(1 - 3). Rheumatic fever has, to be sure, become rarer in some regions; the increase in mild or sub-clinical cases requires a correspondingly more detailed evaluation by means of serological methods. Among the assays for antibodies to various streptococcal exoenzymes, the Anti-Streptolysin O assay should be given prominence as a sensitive parameter since it has been found to be elevated in about 80 to 85% of the cases. Moreover, up to now, international standardization has only been performed for ASL.

Principle of the Method

Polystyrene particles coated with streptolysin O are aggregated when mixed with samples containing human ASL. These aggregates scatter a beam of light passed through the sample. The intensity of the scattered light is proportional to the concentration of the relevant protein in the sample. The result is evaluated by comparison with a standard of known concentration.

Reagents

Materials provided

N Latex ASL, Code No. OUSS

N ASL Reagent, 3 vials for 2 mL each

N ASL Accelerator, 2 vials with 4 mL each or

N ASL Reagent, 4 vials for 3 mL each

N ASL Accelerator, 4 vials with 4 mL each

Composition

N ASL Reagent consists of a lyophilisate of polystyrene particles which are coated with streptolysin O.

N ASL Accelerator consists of a solution of polyethyleneglycol and detergents.

Preservative:

N ASL Reagent after reconstitution: Sodium azide < 0.6 g/L

N ASL Accelerator: Sodium azide < 1 g/L

Warnings and Precautions

- For in vitro diagnostic use.
- Reagents containing sodium azide must be handled with due caution.
Do not ingest or allow to contact skin or mucous membranes! If discarded into drain, flush with a large volume of water to prevent azide build-up. Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.
- Each individual blood donation for use in the manufacture of the N ASL Reagent was tested for HBsAg, anti-HCV, anti-HIV1 and anti-HIV2 by FDA-required testing. Only donations with negative findings are used for manufacture.
Nevertheless, since absence of infectious agents cannot be proven, all samples (e.g. patient serum) and products obtained from human blood should always be handled with due care, observing the precautions recommended for biohazardous material (4).

Preparation of Reagents

N ASL Reagent: The lyophilized content of a vial must be reconstituted in the volume of distilled water given on the label. The reagent may be used 15 minutes after reconstitution. Shake carefully before the first use.

N ASL Accelerator is ready-for-use.

Storage and Stability

Stability at +2 to +8 °C:

see expiry date on label;

Stability once opened:

4 weeks (reconstituted N ASL Reagent and N ASL Accelerator), if stored at +2 to +8 °C securely capped immediately after each use.

Do not freeze.

On-board Stability:

On-board stability may vary, depending on BN* System used and laboratory conditions. The reagents should not left opened on the BN* A, BN* 100 and BN* II System, but stored again at +2 to +8 °C securely capped immediately after use. Further details are inclosed in the BN* II and BN ProSpec® System Instruction Manual.

Materials required but not provided

BN* System

N Rheumatology Standard SL, Code No. OQKZ

N/T Rheumatology Controls SL/1 and 2, Code No. OQDB and OQDC

N Reaction Buffer, Code No. OUMS

N Diluent, Code No. OUMT

BN* II Evaporation Stoppers (optional), Code No. OVLE

Additional and supplies materials as described in your BN* System Instruction Manual.

Specimens

Suitable samples are human serum, either as fresh as possible (stored no more than 8 days at +2 to +8 °C) or stored frozen. Samples can be stored at below –25 °C for up to 3 months if they are frozen within 24 hours after collection and if repeated freeze-thaw cycles are avoided. Serum samples must be completely coagulated and, after centrifugation, must not contain any particles or traces of fibrin.

Lipemic samples or frozen samples which became turbid after thawing, must be clarified by centrifugation (10 minutes at approximately 15,000 x g) prior to testing

Procedure

Notes

- Consult your BN* System Instruction Manual for details regarding operation of the instrument.
- The reagents must not be used beyond the expiry date.
- The lyophilized reagent must not be used until properly reconstituted (at least 15 minutes after the addition of distilled water).
- Allow reagents and samples to equilibrate to room temperature (+15 to +25 °C) before use on the BN* A or BN* 100 Systems. With a BN* II or BN ProSpec® System reagents and samples stored at +2 to +8 °C can be used immediately.
- On the BN* A and BN* 100 Systems, samples should be run at approximately the same ambient temperature (maximum 3 °C deviation) as the measurements used for recording the reference curve.

Assay Protocol for the BN* Systems

The assay protocol is given in the BN* System Instruction Manual and software of the instrument. All steps are performed automatically by the system.

Establishment of the Reference Curve

Reference curves are constructed by multi-point calibration. Serial dilutions of the N Rheumatology Standard SL are automatically prepared by the instrument using N Diluent. The standard dilutions are to be used within 4 hours. The reference curve is valid for 2 weeks and can be used beyond this period, as long as the accuracy controls, e. g. N/T Rheumatology Control SL/1 and SL/2, are reproduced within their respective confidence interval. If a different lot of reagent is used, a new reference curve must be recorded.

The exact measuring range depends upon the concentration of the protein in each lot of N Rheumatology Standard SL. For typical figures refer to the respective BN* System Instruction Manual.

Assay of Specimens

Samples are automatically diluted 1:400 with N Diluent. The dilutions must be measured within 4 hours. If the results obtained are outside the measuring range, the assay can be repeated using a higher or lower dilution of the sample. Refer to BN* System Instruction Manual for information on repeat measurements using other dilutions.

Internal Quality Control

Assay N/T Rheumatology Controls SL/1 and SL/2 after each establishment of a reference curve, the first opening of a reagent vial as well as with each run of samples. The controls are assayed and evaluated as for patient samples. The assigned value and confidence interval are listed in the Table of Assigned Values of the respective control. If the result of the controls is outside the confidence interval, the determination must be repeated. If the repeated determination confirms the deviation, a new reference curve should be established. Do not release patient results until the cause of deviation has been identified and corrected.

Results

The results are evaluated automatically by means of a logit-log function.

Limitations of the Procedure

Turbidity and particles in the sample may interfere with the determination. Therefore, samples containing particles must be centrifuged prior to testing. Lipemic samples, which cannot be clarified by centrifugation (10 minutes at approximately 15,000 x g) and heat-inactivated samples must not be used. In the case of monoclonal gammopathies (e.g. Waldenström's macroglobulinaemia) the measurements yield false highly elevated ASL concentrations in some patients. Such samples should be assayed by another method.

Due to matrix effects, inter-laboratory survey samples and control samples may yield results that differ from those obtained with other methods. It may therefore be necessary to assess these results in relation to method-specific target values.

Reference Interval

As upper limit of the reference interval 200 IU/mL is accepted internationally. ASL concentrations depend upon the age of the patients (5), on the geographic location, and on the local frequency of streptococcal infections. In a study of 503 central European blood donors values up to 370 IU/mL (95th percentile) were yielded with N Latex ASL. In a further study performed in Japan on a population of 477 blood donors, values up to 240 IU/mL (95th percentile) were measured. Nevertheless, each facility should determine its own reference intervals since values may vary depending on the population studied.

Since an individual result provides no information about a change of the ASL concentration, an appropriate assessment is only possible if the assay is repeated after 1 or 2 weeks.

Specific Performance Characteristics

Sensitivity

The analytical sensitivity of the assay is determined by the lower limit of the reference curve and depends therefore on the concentration of the protein in the N Rheumatology Standard SL.

Specificity

The test is specific for anti-streptolysin O.

Precision

Two control samples with different ASL concentrations were analyzed in duplicate on 20 days each with N Latex ASL. Analysis of variance of the results (6) yielded the following coefficients of variation (CV):

Sample	Precision			
	Mean (IU/mL)	Intra Assay CV (%)	Inter Assay CV (%)	Total CV (%)
Control Serum 1	131	5.7	4.8	6.7
Control Serum 2	385	6.7	5.7	8.0

Method Comparison

64 serum samples with ASL concentrations ranging from 27 to 600 IU/mL were assayed on the BN* System with N Latex ASL (y) and in parallel with a hemolysis inhibition test (x) which was also calibrated against the WHO reference preparation. Comparison of the results by regression analysis yielded the following line: y = 0.80 x + 16.5 (IU/mL).

Note

The values cited for specific performance characteristics of the assay represent typical values and are not to be regarded as specification for N Latex ASL.

Bibliography

- Myrvik QN, Pearsall NM, Weiser RS. Fundamentals of Medical Bacteriology and Mycology. Philadelphia: Lea and Febiger, 1976.
- Hederstedt B, Holm S, Wadström T. Medical Laboratory, Behringwerke AG, Marburg, 1979; 7: 24-8.
- Dillon HC, Avery Reeves MS. Streptococcal immune responses in nephritis after skin infections. Am J Med 1974; 56: 333-46.
- U.S. Department of Health and Human Services CDC, Biosafety in Microbiological and Biomedical Laboratories, HHS Publication (CDC) 93-8395; 1999; Section II; 8-16
- Renneberg J, Söderström M, Prellner K, et al. Age-related variations in anti-streptococcal antibody levels. Eur J Clin Microbiol Infect Dis. 1989; 8: 792-5.
- National Committee for Clinical Laboratory Standards. Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline. NCCLS document EP5-A (ISBN 1-56238-368-X). NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898, USA, 1999.

BN ProSpec is a registered trademark of Dade Behring Marburg GmbH in the USA, Germany and other countries.

* BN is a trademark of Dade Behring Marburg GmbH in the USA.

Dade Behring Marburg GmbH
Emil-von-Behring-Str. 76
D-35041 Marburg
www.dadebehring.com

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BN ProSpec é uma marca registada da Dade Behring Marburg GmbH nos EUA, Alemanha e noutros países.
* BN é uma marca da Dade Behring Marburg GmbH nos EUA.