N Latex Myoglobin

Intended Use

In vitro diagnostic reagents for quantitative determination of myoglobin in human serum or heparinized and EDTA plasma by means of particle-enhanced immunonephelometry using the BN[®] Systems.

Summary and Explanation

Myoglobin is a heme-containing protein found primarily in the cells of striated muscle (skeletal muscle) as well as in cardiac muscle. When these muscles are damaged, myoglobin is released and appears in the blood. If cardiac muscle is damaged through, for example, inadequate perfusion following occlusion of a coronary artery (acute myocardial infarction), enzymes and other proteins are released into the vascular system due to the disruption of membrane function. Transported by the lymph, myoglobin appears in high concentrations in the peripheral blood 2 to 3 hours after symptom onset and thus reaches pathological levels 3 to 4 hours before the rise in myocardial markers such as CK-MB or troponin. Peak levels of myoglobin are reached after 6 to 9 hours while those of the aforementioned proteins do not reach peak values until after about 1 to 2 days. The steep increase in serum myoglobin concentration is followed by a rapid fall due to renal elimination so that, in cases without complications, the values measured about 24 hours after the infarct can have returned to within the normal range. Acute myocardial infarction can be ruled out with high probability if myoglobin levels in the patient blood are not found to be elevated 2 to 3 hours after symptom onset and if this is confirmed by a second determination approximately 5 hours later (1, 2). For suspected cases of acute myocardial infarction, the use of an early marker such as myoglobin is recommended as a diagnostic aid (1). Successful myocardial reperfusion by thrombolytic therapy will lead within approximately 90 minutes to a further release of cardiac proteins into the circulation (1). Elevations in blood myoglobin concentrations unrelated to myocardial infarction can occur, for example, after muscle trauma, resuscitation, myopathies, shock, strenuous exercise or in renal insufficiency involving decreased elimination. Extreme concentrations occur in rhabdomyolysis (3).

Principle of the Method

Polystyrene particles coated with specific antibodies to human myoglobin are aggregated when mixed with samples containing myoglobin; these aggregates scatter a beam of light passed through the mixture. The intensity of the scattered light is proportional to the concentration of the relevant protein in the sample. The result is evaluated by comparison with a standard of known concentration.

Reagents

Materials provided

N Latex Myoglobin, Code No. OWIA

N Myoglobin Reagent, 3 vials for 2.0 mL

N Myoglobin Standard (human), 3 vials for 1.0 mL

N Myoglobin Control (human), 3 vials for 0.5 mL

N Myoglobin Supplementary Reagent A, 3 vials of 0.5 mL

N Myoglobin Supplementary Reagent B, 1 vial of 1.6 mL

Composition and Standardization

N Myoglobin Reagent consists of lyophilized polystyrene particles coated with rabbit anti-human myoglobin antibodies.

N Myoglobin Standard consists of lyophilized, pooled human sera supplemented with human myoglobin. The myoglobin concentration for the reconstituted standard is given on the vial label.

N Myoglobin Control consists of lyophilized, pooled human sera supplemented with human myoglobin. The myoglobin concentration for the reconstituted control is given on the vial label.

The concentration of N Myoglobin Standard and N Myoglobin Control was calibrated with reference to highly purified human myoglobin (internal reference preparation of Dade Behring Marburg GmbH).

N Myoglobin Supplementary Reagent A consists of a buffered rabbit serum solution.

N Myoglobin Supplementary Reagent B consists of a detergent solution.

Preservatives: N Myoglobin Reagent after reconstitution: Gentamicin 6.25 mg/L, Amphotericin 0.625 mg/L

N Myoglobin Standard and

N Myoglobin Control after reconstitution: Gentamicin 100 mg/L, Amphotericin 1.25 mg/L

N Myoglobin Supplementary Reagents: Sodium azide <1g/L.

Warnings and Precautions

- For *in vitro* diagnostic use.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.
- Each individual blood donation for use in manufacture of N Myoglobin Reagent, N Myoglobin Standard or N Myoglobin Control was tested for HBsAg, anti-HCV, anti-HIV1 and anti-HIV2 by FDA-required testing. Only donations with negative findings are used for manufacture. Nevertheless, since absence of infectious agents cannot be proven, all samples (e.g. patient sera and plasma) and products (e.g. standard and control sera) obtained from human blood should always be handled with due care, observing the precautions recommended for biohazardous material (4).

Preparation of the Reagents

N Myoglobin Reagent: Dissolve the lyophilized contents of a vial in 2.0 mL of distilled water. The reagent may be used 15 minutes after reconstitution. The suspension should be carefully inverted to mix before its first use. Avoid vigorous shaking (foam formation).

N Myoglobin Standard: Dissolve the lyophilized contents of a vial in 1.0 mL of distilled water. Shake gently to mix. Let stand for 15 minutes before use.

N Myoglobin Control: Dissolve the lyophilized contents of a vial in 0.5 mL of distilled water. The solution should be carefully inverted to mix before its first use. Avoid vigorous shaking (foam formation). Let stand for 15 minutes before use.

N Myoglobin Supplementary Reagents: Pipette 0.5 mL of N Myoglobin Supplementary Reagent B into a vial of N Myoglobin Supplementary Reagent A and shake gently to mix.

Storage and Stability

Stability at +2 to +8 °C: The expiry date is given on the label;

Stability once opened:

1 week (reconstituted N Myoglobin Reagent, reconstituted N Myoglobin Standard, reconstituted N Myoglobin Control as well as the mixture of the N Myoglobin Supplementary Reagents),

6 weeks (N Myoglobin Supplementary Reagent B) if stored tightly closed at +2 to +8 °C directly after use. Do not freeze the reagents.

On-board Stability:

N Myoglobin Reagent and N Myoglobin Supplementary Reagent - at least 3 days (at 8 hours per day), or a comparable period of time.

Note: On-board stability may vary depending on the BN[®] System used as well as on laboratory conditions. For further details, refer to BN[®] II or BN ProSpec[®] System Instruction Manual.

The on-board stability of the N Myoglobin Control on the BN ProSpec[®] System is stated in the Instruction Manual of the system.

Materials required but not provided

BN[®] System

N Reaction Buffer, Code No. OUMS

N Diluent, Code No. OUMT

BN[®] II Evaporation Stoppers (optional), Code No. OVLE

Additional materials and supplies as described in the Instruction Manuals for the BN[®] Systems.

Specimens

Suitable assay specimens are human serum samples and heparinized and EDTA plasma samples, either as fresh as possible (stored for no more than 8 days at +2 to +8 °C) or stored in frozen condition. Samples frozen within 24 hours after collection can be used within 1 month if stored at below -20 °C and if repeated freeze-thaw cycles are avoided. Serum samples must have coagulated completely and, after centrifugation, must not contain any particles or traces of fibrin. Lipemic samples or turbid frozen samples which are turbid after thawing must be clarified by centrifugation (10 minutes at approximately 15,000 x g) prior to testing.

OWIA G11 E0541 (1653) H 1

Edition November 2003

Procedure

Notes

- For details regarding operation of the BN[®] Systems consult the relevant Instruction Manuals.
- The reagents must not be used beyond the stated expiry date.
- Only components (reagent, standard, control, supplementary reagents) from kits of the same lot number can be used together; do not combine kit components from kits of different lots.
- Reconstituted lyophilized reagents must not be used until after the stated reconstitution times (at least 15 minutes after the addition of distilled water).
- Reagents and samples should have reached room temperature (+15 to +25 °C) before assay on the BN[®] A and BN[®] 100 System. Reagents and samples stored at +2 to +8 °C can be used for immediate assay on the BN[®] II and BN ProSpec[®] System.
- On the BN[®] A and BN[®] 100 System, samples must be run at the same ambient temperature (maximum 3 °C deviation) as the measurements used for recording the reference curve.

Assay Protocol for the BN[®] Systems

The assay protocol for serum and plasma is given in the BN[®] System Instruction Manual and instrument software. All the steps necessary are performed automatically by the system.

Establishment of the Reference Curve

Reference curves are constructed by multi-point calibration. Serial dilutions of the N Myoglobin Standard are automatically prepared for this purpose using N Diluent. The standard dilutions must then be used within 4 hours. The assigned value is given on the respective vial label. On the BN ProSpec[®] System, it can be read into the system from the lot data diskette (Code No. OVLP). The reference curve is valid for one week.

A new reference curve must be recorded if a new reagent vial or a different lot of reagent is used, or if the lot of N Reaction Buffer is changed. If large numbers of samples are to be assayed, it is recommended that the volumes of the reagents needed for all assays in one week should be individually pooled and that the reference curve should then be established using the pooled reagents. The exact measuring range depends on the concentration of the protein in each lot of N Myoglobin Standard. For typical figures refer to the respective BN[®] System Instruction Manual.

Assay of Specimens

Samples are automatically diluted 1:20 with N Diluent. The dilutions must be assayed within 4 hours. If the results obtained are outside the measuring range, the measurement can be repeated using a higher or lower dilution of sample. Details on repeat measurements using further sample dilutions are described in the Instruction Manuals for the BN[®] Systems.

Internal Quality Control

N Myoglobin Control should be used after each establishment of a reference curve, after the first use of a reagent vial as well as in each series of samples. The control is assayed and evaluated in the same way as for the patient samples. The assigned value of the N Myoglobin Control is given on the vial label. The confidence interval is given by the assigned value ±20%. On the BN ProSpec[®] System, the assigned value can be read into the system from the lot data diskette (Code No. OVLP). The stated assigned value is intended for use as an assayed accuracy control. If used for precision control, the user should establish the control value and control limits in a preliminary period. If the result of the control measurements is outside of the confidence interval, the control determination must be repeated. If the repeat measurement confirms the deviation, a new reference curve should be established. Results for patients must not be released until the cause of the deviation has been identified and rectified.

Results

The results are evaluated automatically by means of a logit-log function.

Limitations of the Procedure

Interference by rheumatoid factors (up to 3 000 IU/mL) is suppressed by use of the Supplementary Reagent. However, in isolated cases, very high levels of rheumatoid factor can interfere with the assay. Such samples must be assayed with another method.

Turbidity and particles in the samples can interfere with the determination. Therefore samples containing particles must be centrifuged prior to testing. Heat-inactivated samples, as well as lipemic samples which cannot be clarified by centrifugation (10 minutes at approximately 15,000 x g), must be excluded from the assay.

After measuring samples with myoglobin levels of approximately > 30,000 µg/L, false-high values may be obtained in a myoglobin determination performed directly afterwards. If this is suspected, a repeat measurement should be performed.

Due to matrix effects, the results obtained for control samples and interlaboratory survey samples may differ depending on the assay method used. It may therefore be necessary to assess these results in relation to method-specific target values.

Reference Interval

The plasma and serum myoglobin levels in clinically healthy persons depend on age and sex and exhibit a wide range of distribution. The upper limit of the reference range (97.5th percentile) was established to be a myoglobin concentration of 70 µg/L (5).

Nevertheless, each laboratory should determine its own reference intervals as these are influenced by many factors which may vary depending on the population studied.

Specific Performance Characteristics

Sensitivity

The analytical sensitivity of the assay is determined by the lower limit of the reference curve and thus depends on the concentration of the protein in the N Myoglobin Standard.

Specificity

No cross-reactivity is known to occur with the antibody used.

Precision

The coefficients of variation (CV) obtained for determinations of 3 different myoglobin concentrations (about 85, 160, and 310 µg/L) with N Latex Myoglobin on the BN[®] System were 2.0%, 2.4%, and 2.9%, respectively, for within-series precision (n = 20) and 4.8%, 4.2%, and 5.0% for reproducibility from day-to-day (n = 15).

Method Comparison

In a study conducted in Europe, 117 serum samples (myoglobin concentrations up to about 8 400 µg/L) were comparatively assayed with N Latex Myoglobin on the BN[®] System (y) and with a commercial radioimmunoassay (x). An analysis of the correlation between the results yielded the following regression equation: y = 0.95 x – 4.33 µg/L.

Note: The values cited for specific performance characteristics of the assay represent typical values and are not to be regarded as specifications for N Latex Myoglobin.

Bibliography

- Wu AHB, Apple FS, Gibler WB, et al. National Academy of Clinical Biochemistry Standards of Laboratory Practice: recommendations for the use of cardiac markers in coronary artery diseases. Clin Chem 1999; 45: 1104-21.
- Plebani M, Zaninotto M. Diagnostic strategies in myocardial infarction using myoglobin measurement. Eur Heart J 1998; 19 Suppl N: N12-5.
- Beetham R. Biochemical investigation of suspected rhabdomyolysis. Ann Clin Biochem 2000; 37: 581-7.
- U.S. Department of Health and Human Services CDC, Biosafety in Microbiological and Biomedical Laboratories, HHS Publication (CDC) 93-8395; 1999; Section II; 8-16
- Chapelle JP. Akuter Myokardinfarkt: Diagnose und Überwachung durch Myoglobin-Bestimmung. Diagnose & Labor 1989; 39: 171–7.


BN ProSpec is a registered trademark of Dade Behring Marburg GmbH in the USA, Germany and other countries.

* BN is a trademark of Dade Behring Marburg GmbH in the USA.

 Dade Behring Marburg GmbH
Emil-von-Behring-Str. 76
D-35041 Marburg
www.dadebehring.com

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BN ProSpec é uma marca registada da Dade Behring Marburg GmbH nos EUA, Alemanha e noutros países.
* BN é uma marca da Dade Behring Marburg GmbH nos EUA.