N Latex Ferritin

Intended Use

In vitro diagnostic reagents for the quantitative determination of ferritin in human serum or heparinized plasma by means of particle-enhanced immunonephelometry using the BN* Systems.

Summary and Explanation

The serum ferritin concentration exhibits a direct correlation to the quantity of storage iron in the reticulo-endothelial system, so that its determination is used for diagnosis and monitoring of iron deficiency and iron overload (1 - 3). As iron deficiency can be detected prior to the appearance of anemia its diagnosis is particularly important to enable prevention of nutritionally-related anemia. Furthermore, ferritin determinations allow to monitor the efficacy of oral iron therapy through normalization of the storage iron levels (1).

Determination of serum ferritin also supports the diagnosis of iron-excess situations like hemochromatosis, and monitoring of therapeutic phlebotomy (2).

The serum level of ferritin is found to be elevated in patients with chronic inflammations, infections, neoplasias and liver disease. Specificity of ferritin determinations can be enhanced by combination with determination of soluble transferrin receptor (sTfR) (1,3).

Principle of the Method

Polystyrene particles coated with antibodies specific to human ferritin are aggregated when mixed with samples containing human ferritin. These aggregates scatter a beam of light passed through the sample. The intensity of the scattered light is proportional to the concentration of the relevant protein in the sample. The result is evaluated by comparison with a standard of known concentration.

Reagents

Materials provided

N Latex Ferritin, Code No. OQTH

N Latex Ferritin Reagent, 3 vials with 3 mL each

N Ferritin Supplementary Reagent A, 3 vials with 1 mL each

N Ferritin Supplementary Reagent B, 1 vial with 3 mL

Composition

N Latex Ferritin Reagent consists of a suspension of polystyrene particles coated with rabbit anti-human ferritin antibodies.

N Ferritin Supplementary Reagent A consists of a buffered solution of rabbit serum.

N Ferritin Supplementary Reagent B consists of a buffered solution of detergents.

Preservatives

N Latex Ferritin Reagent: Gentamicin 6.25 mg/L, Amphotericin 0.625 mg/L

N Ferritin Supplementary Reagents: Sodium azide < 1 g/L

Warnings and Precautions

- For *in vitro* diagnostic use
- Reagents containing sodium azide must be handled with due caution:
Do not ingest or allow to contact skin or mucous membranes. If discarded into drain, flush with a large volume of water to prevent azide build-up. Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.
- Each individual blood donation for use in the manufacture of the N Latex Ferritin Reagent was tested for HBsAg, anti-HCV, anti-HIV1 and anti-HIV2 by FDA-required testing. Only donations with negative findings are used for manufacture.
Nevertheless, since absence of infectious agents cannot be proven, all samples (e.g. patient serum and plasma) and products obtained from human blood should always be handled with due care, observing the precautions recommended for biohazardous material (4).

Preparation of Reagents

N Latex Ferritin Reagent is supplied ready-for-use. The suspension should be carefully shaken before the first use.

N Ferritin Supplementary Reagent: Pipette 1.0 mL N Ferritin Supplementary Reagent B into one vial of N Ferritin Supplementary Reagent A and mix gently.

Storage and Stability

Stability at +2 to +8 °C:

see expiry date on label;

Stability once opened:

4 weeks (N Latex Ferritin Reagent, mixture of the N Ferritin Supplementary Reagents) and

12 weeks (N Ferritin Supplementary Reagent B), if stored at +2 to +8 °C securely capped immediately after each use.

Do not freeze.

On-board Stability:

A minimum of five days, at eight hours per day, or comparable period of time.

Note: On-board stability may vary, depending on BN* System used and laboratory conditions. For further details, refer to BN* II and BN ProSpec® System Instruction Manual;

Materials required but not provided

BN* System

N Protein Standard SL (human), Code No. OQIM

N/T Protein Control SL/L (human), Code No. OQIN

N/T Protein Control SL/M (human), Code No. OQIO

N/T Protein Control SL/H (human), Code No. OQIP

N Diluent, Code No. OUMT

BN* II Evaporation Stoppers (optional), Code No. OVLE

Additional materials and supplies as described in your BN* System Instruction Manual.

Specimens

Suitable samples are human serum or heparinized plasma, either as fresh as possible (stored no more than 7 days at +2 to +8 °C) or stored frozen. Samples can be stored at below -20 °C for up to 3 months if they are frozen within 24 hours after collection and if repeated freeze-thaw cycles are avoided. Serum samples must be completely coagulated and, after centrifugation, must not contain any particles or traces of fibrin.

Lipemic samples, or frozen samples which became turbid after thawing, must be clarified by centrifugation (10 minutes at approximately 15,000 x g) prior to testing.

Procedure

Notes

- Consult your BN* System Instruction Manual for details regarding operation of the instrument.
- The reagents must not be used beyond the expiry date.
- Only components (reagent, supplementary reagents) from kits of the same lot may be used together; do not combine kit components from kits of different lots.
- Allow reagents and samples to equilibrate to room temperature (+15 to +25 °C) before use on the BN* A or BN* 100 Systems. With a BN* II or BN ProSpec® System reagents and samples stored at +2 to +8 °C can be used immediately.
- On the BN* A and BN* 100 Systems, samples should be run at approximately the same ambient temperature (maximum 2 °C deviation) as the measurements used for recording the reference curve.

OQTH G11 E0540 (131) H 1

Edition February 2002

Assay Protocol for BN* Systems

The assay protocol, for serum and plasma, is given in the BN* System Instruction Manual and software of the instrument. All steps are performed automatically by the system.

Establishment of the Reference Curve

Reference curves are constructed by multi-point calibration. Serial dilutions of the N Protein Standard SL are automatically prepared by the instrument using N Diluent. The standard dilutions are to be used within 4 hours. The reference curve is valid for 4 weeks and can be used beyond this period, as long as the accuracy controls, e. g. N/T Protein Controls SL/L, M and H, are reproduced within their respective confidence interval. If a different lot of reagent is used, a new reference curve must be recorded.

The exact measuring range depends upon the concentration of the protein in each lot of N Protein Standard SL. For typical figures refer to the respective BN* System Instruction Manual.

Assay of Specimens

Samples are automatically diluted 1:20 with N Diluent. The diluted samples must be measured within 4 hours. If the results obtained are outside the measuring range, the assay can be repeated using a higher or lower dilution of the sample. Refer to BN* System Instruction Manual for information on repeat measurements using other dilutions.

Internal Quality Control

Assay N/T Protein Controls SL/L, M and H after each establishment of a reference curve, the first opening of a reagent vial as well as with each run of samples. The controls are assayed and evaluated as for patient samples. The assigned value and confidence interval are listed in the Table of Assigned Values of the respective control.

If the result of the controls is outside the confidence interval, the determination must be repeated. If the repeated determination confirms the deviation, a new reference curve should be established. Do not release patient results until the cause of deviation has been identified and corrected.

Results

The results are evaluated automatically by means of a logit-log function.

Limitations of the Procedure

Interferences by rheumatoid factors (up to 2800 IU/mL) are suppressed by use of the N Ferritin Supplementary Reagent. With samples containing very high levels of ferritin, no antigen excess phenomenon has been observed up to the highest tested concentration of approximately 10,000 µg/L for measurements using the 1:20 dilution.

Turbidity and particles in the samples may interfere with the determination. Therefore, samples containing particles must be centrifuged prior to testing. Lipemic samples, which cannot be clarified by centrifugation (10 minutes at approximately 15,000 x g) must not be used.

Due to matrix effects, inter-laboratory survey samples and control samples may yield results that differ from those obtained with other methods. It may therefore be necessary to assess these results in relation to method-specific target values.

Reference Intervals

The ferritin concentration in the serum of healthy persons is age- and sex-dependent and exhibits a wide range of distribution. Samples from 456 healthy individuals from Central Europe were evaluated with N Latex Ferritin and the following data were obtained (5, plus additional calculations):

Men	(n=216):	20 - 290 µg/L
Women, premenopausal	(n=193):	4.5 - 170 µg/L
Women, postmenopausal	(n=47):	24 - 260 µg/L

Nevertheless, each facility should determine its own reference intervals since values may vary depending on the population studied.

Specific Performance Characteristics

Sensitivity

The analytical sensitivity of the assay is determined by the lower limit of the reference curve and depends therefore on the concentration of the protein in N Protein Standard SL.

Specificity

Investigation on crossreactivity of human spleen, liver, placenta and heart ferritin using N Latex Ferritin yielded the following results:

Ferritin source:	Crossreactivity:
Spleen	90.4 %
Liver	94.3 %
Placenta	103.4 %
Heart	7.8 %

Precision

Precision data were generated by evaluation of three levels of human serum pools ranging from ferritin concentrations of 36 to 620 µg/L in a manner consistent with NCCLS Guideline EP5-A (6). Coefficients of variation (CV) of intra-assay precision ranged from 1.0 to 4.6 %, while CV of inter-assay precision ranged from 1.2 to 3.1 % and total CV from 1.6 to 5.1 %.

Method Comparison

79 serum samples (ferritin concentrations from 12 to 440 µg/L) were analyzed comparatively using the N Latex Ferritin (y) versus a commercial ferritin nephelometric assay (x). Comparison of the results by regression analysis yielded the following regression equation: $y = 1.05 \times - 0.87$ (µg/L).

79 serum samples (ferritin concentrations from 12 to 440 µg/L) were analyzed comparatively using the N Latex Ferritin (y) versus a commercial microparticle enzyme immunoassay (x). Comparison of the results by regression analysis yielded the following regression equation: $y = 1.07 \times - 0.47$ (µg/L).

Note

The values cited for specific performance characteristics of the assay represent typical values and are not to be regarded as specification for N Latex Ferritin.

Bibliography

- Cook JD. Defining optimal body iron. Proc Nutr Soc 1999; 58: 489-95.
- Andrews NC. Disorders of iron metabolism. N Engl J Med 1999; 341: 1986-95.
- Sherwood RA, Pippard MJ, Peters TJ. Iron homeostasis and the assessment of iron status. Ann Clin Biochem 1998; 35: 693-708.
- U.S. Department of Health and Human Services CDC, Biosafety in Microbiological and Biomedical Laboratories, HHS Publication (CDC) 93-8395; 1999; Section II; 8-16
- Van den Bosch G, Van den Bossche J, Wagner C, et al. Determination of iron metabolism-related reference values in a healthy adult population. Clin Chem. 2001;47:1465-7.
- National Committee for Clinical Laboratory Standards. Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline. NCCLS document EP5 - A (ISBN 1-56238-368-X). NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898, USA, 1999.

BN ProSpec is a registered trademark of Dade Behring Marburg GmbH in the USA, Germany and other countries.

* BN is a trademark of Dade Behring Marburg GmbH in the USA.

 Dade Behring Marburg GmbH
Emil-von-Behring-Str. 76
D-35041 Marburg
www.dadebehring.com

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BN ProSpec é uma marca registada da Dade Behring Marburg GmbH nos EUA, Alemanha e noutros países.
* BN é uma marca comercial da Dade Behring Marburg GmbH nos EUA.