NycoCard[®] *HbA1c*

REF 1042184

GB	3 .. 6
DE	7 .. 10
FR	11 .. 14
IT	15 .. 18
ES	19 .. 22
PT	23 .. 26
NL	27 .. 30
GR	31 .. 34
PL	35 .. 38
CZ	39 .. 42
RU	43 .. 46
HU	47 .. 50
RO	51 .. 54
NO	55 .. 58
DK	59 .. 62
FI	63 .. 66

SYMBOLS/ SYMBOLE/ SYMBOLES/ SIMBOLI/ SÍMBOLOS/ SÍMBOLO/ SYMBOLEN/
 ΣΥΜΒΟΛΑ/ SYMBOLE/ SYMBOLY/ СИМВОЛЫ/ JELMAGYARÁZAT/ SIMBOLURI/
 SYMBOLER/ SYMBOLER/ SYMBOLIT

Conformity to the European directive 98/79/EC
 Konformität mit der europäischen Richtlinie 98/79/EC
 Conformité avec la directive européenne 98/79/CE
 Conformità con la direttiva europea 98/79/CE
 Conformidad con la directiva europea 98/79/CE
 Conformidade com a Directiva europeia 98/79/EC
 Voldoet aan de Europese richtlijn 98/79/EG
 Συμμόρφωση με την οδηγία 98/79/ΕΚ για ροιόντα IVD
 Zgodność z dyrektywą IVD 98/79/EC
 Ve shodě s předpisy evropské směrnice 98/79/ECC
 Соответствие с европейской директивой 98/79/EC
 Megfelel az európai 98/79/EC irányelvnek
 Conformitate cu directiva europeană 98/79/EC
 Overensstemmelse med EU-direktiv 98/79/EC
 Overensstemmelse med EU-direktiv 98/79/EC
 Yhdennukainen Euroopan direktiivin 98/79/EC kanssa

IVD

In Vitro Diagnostic Medical Device
 Medizinisches Gerät zur *In-vitro*-Diagnostik
 Dispositif médical de diagnostic *in vitro*
 Dispositivo medicale per la diagnostica *in vitro*
 Dispositivo médico de diagnóstico *in vitro*
 Dispositivo médico de diagnóstico *in vitro*
 Medisch hulpmiddel voor *in-vitro* diagnostiek
In Vitro διαγνωστικό ιατροτεχνολογικό προϊόν
 Urządzenie medyczne do diagnostyki *in vitro*
 Pro použití *In Vitro*
 Медицинское приспособление для диагностики
In Vitro
In vitro diagnosztikai orvostechnikai eszköz
In vitro Dispozitiv medical pentru diagnostic
 Medisinsk utstyr til *in vitro* diagnostikk
In vitro diagnostisk medicinsk udstyr
In vitro diagnostinen lääkinnäällinen laite

2°C
 8°C
 2...8°C
 2...25°C

Storage temperature
 Lagerungstemperatur
 Température de conservation
 Temperatura di conservazione
 Temperatura de almacenamiento
 Temperatura de conservação
 Beewaarttemperatuur
 Θερμοκρασία αποθήκευσης
 Temperatura przechowywania
 Teplota skladování
 Температура хранения
 Tárolási hőmérséklet
 Temperatura de depozitare
 Oppbevaringstemperatur
 Opbevaringstemperatur
 Säilytyslämpötila

Expiry date (year-month)
 Verfallsdatum (Jahr-Monat)
 Date de péremption (année-mois)
 Data di scadenza (anno-mese)
 Fecha de caducidad (año-mes)
 Data de validade (ano-mês)
 Uiterste gebruiksdatum (jaar-maand)
 Ημερομηνία λήξης (έτος-μήνας)
 Data ważności (rok-miesiąc)
 Datum expirace (rok-měsíc)
 Срок годности (год-месяц)
 Lejárata (év-hónap)
 Data expirării (an-lună)
 Utløpsdato (år-måned)
 Udløbsdato (år-måned)
 Viimeinen käyttöpäivämäärä (vuosi-kuukausi)

LOT

Lot number
 Chargen-Nr.
 Numéro de lot
 Numero di lotto
 Número de lote
 Número do lote
 Partijnummer
 Αριθμός Παρτίδας
 Numer partii
 Šarže
 номер партии
 Gyártási szám
 Numáru de lot
 Lotnummer
 Lot-number
 Eränumero

REF

Catalogue number
 Bestellnummer
 Numéro de référence du catalogue
 Numero di catalogo
 Número de catálogo
 Número do catálogo
 Catalogus nummer
 Αριθμός καταλόγου
 Numer katalogowy
 Katalogové číslo
 Номер по каталогу
 Katalógusszám
 Catalogul numărul
 Artikkelnummer
 Katalognummer
 Luettelonumero

Manufacturer
 Hersteller
 Fabricant
 Produttore
 Fabricante
 Fabricante
 Fabricant
 Κατασκευαστής
 Producent
 Výrobce
 Производител
 Gyártó
 Producatorul
 Producent
 Producent
 Valmistaja

PRODUCT DESCRIPTION**Intended use**

NycoCard[®] HbA1c is a rapid *in vitro* test for the measurement of glycated hemoglobin (% HbA1c) in human whole blood.

Test principle

NycoCard[®] HbA1c is a boronate affinity assay. The kit contains test devices with a porous membrane filter, test tubes pre-filled with reagent and a washing solution. The reagent contains agents that lyse erythrocytes and precipitate hemoglobin specifically, as well as a blue boronic acid conjugate that binds cis-diols of glycated hemoglobin. When blood is added to the reagent, the erythrocytes immediately lyse. All hemoglobin precipitates. The boronic acid conjugate binds to the cis-diol configuration of glycated hemoglobin. An aliquot of the reaction mixture is added to the test device, and all the precipitated hemoglobin, conjugate-bound and unbound, remains on top of the filter. Any excess of coloured conjugate is removed with the washing solution. The precipitate is evaluated by measuring the blue (glycated hemoglobin) and the red (total hemoglobin) colour intensity with the NycoCard[®] READER II, the ratio between them being proportional to the percentage of HbA1c in the sample.

Kit contents, 24 test kit

TD/Test Device 1 x 24 units
Plastic device containing a membrane filter.

R1/Reagent 1 x 24 x 0.2 mL
Glycinamide buffer containing dye-bound boronic acid and detergents.

R2/Washing Solution 1 x 2.0 mL
Morpholine buffered NaCl solution and detergents.

Material required (not supplied with the kit)

- Capillary tube or pipette (5 µL) for sample collection.
- Capillary tube holder
- Pipette (25 µL) and pipette tips
- NycoCard[®] READER II

Warnings and precautions

- IVD** For *in vitro* diagnostic use.
- Do not interchange components from different test kit lots.
 - Do not use test kits after the expiry date.
 - R1/Reagent and R2/Washing Solution contain sodium azide (0.05%), which is a toxic agent.
 - R2/Washing Solution contains morpholine (50 mmol/L), which is a toxic agent.

TEST CHARACTERISTICS**Analytical specificity**

NycoCard[®] HbA1c measures the total glycated hemoglobin (GHb), but reports a standardised HbA1c value.

Standardisation

NycoCard[®] HbA1c is traceable to the IFCC Reference Method for Measurement of HbA1c¹. HbA1c values are reported according to national recommendations. Please contact your local supplier if the national recommendation is unknown.

Measuring range

NycoCard[®] READER II, SN < 59681:
Measuring range: 3-18% HbA1c

NycoCard[®] READER II, SN ≥ 59681:
Measuring range: 4-15% HbA1c

Measuring interval: 0.1% HbA1c

Reference range²⁻⁴

For methods reporting DCCT traceable values, the upper limit of non-diabetic reference range is approximately 6% HbA1c.

Precision

In professional use, a coefficient of variation (CV) < 5% is usually obtained.

Limitations of the test

- Elevated amounts of glucose, bilirubin, lipids and fructosamine were added to blood samples with normal and elevated HbA1c values. No interference was obtained.
- Pre-glycated hemoglobin does not interfere with the test.
- Hemolysed samples with plasma Hb >3 g/100 mL will interfere with the test system.
- NycoCard® READER II corrects for Hb-concentrations in the range 6-18 g/100 mL.

STABILITY AND STORAGE

Unopened kits

The expiry date of the kit applies to storage at 2-8°C in the original container. Exposure to temperatures above 25°C or humidity above 70% should be avoided. Do not freeze.

Opened kits

R1/Reagent must be stored dark at 2-8°C. Equilibrate the R1/Reagent to room temperature (20-25°C) before use. Equilibration to room temperature can be achieved by holding the tube in a closed hand for 30 seconds. The reagent tube can be stored for maximum 6 hours at room temperature before use. Avoid direct sunlight.

TD/Test Device can be stored at room temperature (15-25°C). Store the test devices in the original bag and avoid humidity below 20% and above 70%. The TD/Test Device should have room temperature when used.

R2/Washing Solution can be stored at room temperature (15-25°C). The R2/Washing Solution should have room temperature when used.

Sample material

Blood samples (with anticoagulant) can be stored up to 10 days at 2-8°C before analysis. Avoid measuring hemolysed samples (see "Limitations of the test"). Do not freeze.

TEST PROCEDURE

Important procedural notes!

- Do not interchange components from different kits or kit lots.
- Bring the R1/Reagent to room temperature (20-25°C) before use.
- Do not touch the membrane with the pipette tip.
- Change the pipette tip between each pipetting step.

Sample material

Capillary blood and venous blood with or without anticoagulant (EDTA, heparin and NaF) can be used.

Internal quality control

A quality control material with NycoCard® HbA1c specific target values should be used to confirm the efficacy of the reagents and the correct performance of the test. Lyophilised control materials should not be used with the NycoCard® HbA1c assay.

Test procedure (illustrations on page 67)

1 Precipitate hemoglobin

Add 5 μ L whole blood to the test tube with R1/Reagent. Mix well. Incubate the tube for minimum 2 minutes, maximum 3 minutes. Use a timer.

Note! Equilibrate the R1/Reagent to room temperature (20-25°C) before use.

2 Apply sample

Remix to obtain a homogenous suspension. Apply 25 μ L of the reaction mixture to a TD/Test Device by holding the pipette approx. 0.5 cm above the test well. Empty the pipette quickly in the middle of the test well. Allow the reaction mixture to soak completely into the membrane (approx. 10 seconds).

Note! Avoid air bubbles.

3 Apply R2/Washing Solution

Apply 25 μ L R2/Washing Solution to the TD/Test Device. Allow the washing solution to soak completely into the membrane. Wait for minimum 10 seconds.

Note! Avoid air bubbles.

4 Read the test result

Read the test result within 5 minutes using the NycoCard® READER II.

Note! Follow the NycoCard® READER II instruction manual.

INTERPRETATION OF RESULT

Interpret the NycoCard® HbA1c test results with careful consideration to the patient's medical history, clinical examinations and other laboratory results. If the test result is questionable or if clinical signs and symptoms appear inconsistent with the

test result, re-test the sample or confirm the result using another method. Analyse control materials frequently to verify the performance of the NycoCard® READER II test system.

TROUBLE SHOOTING

Problem	Possible cause	Corrective action
The precipitate does not cover the whole surface of the test hole (white spots on the membrane surface).	The reaction mixture has been added too slowly or on the "wall" of the test well.	Retest the sample. Add the reaction mixture quickly in the middle of the test well.
	An air bubble has been formed when applying the reaction mixture to the test well.	Retest the sample and avoid air bubbles.
	Incorrect storage of TD/Test Device.	Retest the sample using a correctly stored TD/Test Device.
No precipitate on the TD/Test Device membrane.	Incomplete mixing before incubation, i.e., the capillary tube was not emptied.	Retest the sample. Mix the test tube more thoroughly.
The NycoCard® READER II displays "Hb conc. too low".	The hemoglobin concentration of the sample is below the measuring range.	Retest the sample and use 10 µL of whole blood: fill 2 capillary tubes simultaneously and drop both capillaries into one test tube with R1/Reagent. The NycoCard® READER II will display the correct test result.
The NycoCard® READER II displays "Reduce Hb conc.". .	The hemoglobin concentration of the sample is above the measuring range.	Mix the R1/Reagent from two test tubes (400 µL) and retest the sample (5 µL). The NycoCard® READER II will display the correct test result.
	The volume of R2/Washing Solution added to the test device is too small.	Retest the sample. Make sure that a correct volume of R2/Washing Solution is applied (25 µL).
Unexpected low results.	The R1/Reagent has not been equilibrated to room temperature (20-25°C) before use.	Retest the sample. Make sure that the R1/Reagent is equilibrated to room temperature before use.
	The R1/Reagent has been stored at room temperature/ not protected from sunlight.	Retest the sample with a test kit stored correctly.
	The sample material used was hemolysed.	Retest the sample using a fresh sample material.
The R2/Washing Solution does not soak into the TD/Test Device membrane.	The R1/Reagent has not been equilibrated to room temperature (20-25°C) before use.	Retest the sample. Make sure that the R1/Reagent is equilibrated to room temperature before use.
	The sample material used was hemolysed.	Retest the sample using a fresh sample material.

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1. Jeppsson, JO *et al.*, Approved IFCC Reference Method for the Measurement of HbA1c in Human Blood, *Clin Chem Lab Med* 2002; 40(1):78-89.
2. Little, RR *et al.*, Relationship of glycosylated hemoglobin to oral glucose tolerance. Implications for diabetes screening. *Diabetes* 1988; 37(1):60-64.
3. Little, RR *et al.*, Glycated haemoglobin predicts progression to diabetes mellitus in Pima Indians with impaired glucose tolerance. *Diabetologia* 1994; 37:252-256.
4. Rohlfing, CL *et al.*, Use of GHb (HbA1c) in Screening for Undiagnosed Diabetes in the U.S. Population. *Diabetes Care* 2000; 23(2):187-191.

DISTRIBUTORS

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Medinor ASA
Postboks 94 Bryn
Nils Hansensvei 4
0611 Oslo
Tlf: 22 07 65 00
Faks: 22 07 65 05
E-mail: medinor@medinor.no

Distributør i Danmark:

Medinor A/S
Park Allé 350 E
DK-2605 Brøndby
4000 Roskilde
Tlf.: +45 7015 1041
E-mail: info@medinor.dk

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Medinor Finland Oy Ab
Rajatorpantie 41 C
01640 VANTAA
puh. (09) 8520 2400
fax (09) 8520 2410
E-mail:
etunimi.sukunimi@medinorfinland.fi

Distributor in UK:

Axis-Shield UK
Ouse Road
Bicton Industrial Park
Kimbolton
Cambridgeshire PE 28 0LP
Tel: 01480 862100
Fax: 01480 862101
E-mail:
sue.youngusband@axis-shielduk.com

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Axis-Shield GmbH
Waldhofer Str. 101
69123 Heidelberg
Tel: +49 6221 825 750
Fax: +49 6221 825 751

A distributor list is available at www.axis-shield-poc.com

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