
3M ESPE

Sinfony™

Light-Curing Composite for the Fabrication of Crowns, Bridges, Inlays, and Onlays

- Ⓓ **Lichthärtendes Composite zur Herstellung von Kronen, Brücken, Inlays und Onlays**
- Ⓕ **Composite photopolymérisable pour la fabrication de couronnes, bridges et onlays**
- Ⓘ **Composito ibrido fotoindurente per la realizzazione di corone, ponti, inlay ed onlay**
- Ⓔ **Composite fotopolimerizable para la fabricación de coronas, puentes, inlays y onlays**
- ⒫ **Compósito fotopolimerizável para a fabricação de coroas, pontes, inlays e onlays**
- Ⓐ **Lichtuithardende composiet voor het vervaardigen van kronen, bruggen, inlays en onlays**
- ⒸⓇ **Φωτοπολυμεριζόμενη σύνθετη ρητίνη για την κατασκευή στεφανών, γεφυρών, ενθέτων και επενθέτων**
- Ⓒ **Ljushärdande komposit för framställning av kronor, broar, inlays och onlays**
- ⒻⒶ **Valokovettuva yhdistelmämuovi kruunujen, siltojen, inlayden ja onlayden valmistukseen**
- ⒹⓀ **Lyshærdende komposit til fremstilling af kroner, broer, inlays og onlays**
- Ⓐ **Lysherdende komposittmateriale til fremstilling av kroner, broer, inlays og onlays**

Instructions for Use
Gebrauchsinformation
Mode d'emploi
Informazioni per l'uso
Información de uso
Instruções de uso

Gebruiksinformatie
Πληροφορίες χρήσεως
Bruksanvisning
Käyttöohjeet
Bruksanvisning
Bruksveiledning

11/02
3M ESPE AG
Dental Products

ENGLISH

Sinfony™**Light-Curing Composite for the Fabrication of Crowns, Bridges, Inlays, and Onlays****Instructions for Use**

Table of Contents	Page
Product Description	1
Range of Application	2
Preparation	2
Prepare/Design Metal Framework	2
Bonding System	2
Preparation of the Model	2
Assembly of the Dispenser	2
Dosing from the Dispenser	2
Dosing Aid	3
Dispenser Holder	3
Shade Selection	3
Production of a Veneer	4
Opaker Processing	4
Intensive-Opaker Processing	5
Opaker-Dentin Processing	5
MAGIC-Intensive Color Processing	5
Building-up the Tooth Contours by Layering	5
Color Control	6
Polymerization	6
Shaping and Polishing	7
Correction/Repair	7
Correction after Shaping	7
Repair of Worn or Fractured Restorations	7
All-Composite Restorations	7
Guidelines for Preparatory Work	7
Fabrication of Full-Composite Inlays, Onlays, Crowns, Veneers/Facings	8
Preparatory Steps in the Laboratory for Trying-in the Restoration	8
Preparatory Steps for Cementation in the Dentist's Office	8
Adhesive Cementation of Full-Composite Restorations	8
Customization of Denture Teeth/	
Protemp 3 Garant Temporary Restorations	8
Denture Teeth	8
PROTEMP 3 GARANT Temporary Restorations	8
Customization of Dentures	9
Complete Dentures	9
Model Cast Dentures	9
Veneering of Glass-Fiber-Reinforced Frameworks	9
Vectris®	9
Stick Tech®	9
FibreKor®	10
DC-Tell®	10
Preventing Processing Errors	10
Framework Design	10
Opaker	10
Layering	10
Layering Scheme	10
Polymerization	10
Wall Thickness for Veneers	10
Veneer Cleaning	11
Restoration and Denture Care	11
Repair of Fixed Restorations	11
Notes	11
Incompatibilities	11
Storage and Shelf-life	11
Customer Information	11
Warranty	11
Limitation of Liability	11
Index	12

1. Product Description

Sinfony™, manufactured by 3M ESPE, is a light-curing hybrid composite that consists of very fine particles and is characterized by its high strength. The material is used for full coverage of frameworks made of metal or fiber-reinforced materials.

Sinfony is also suitable for fabrication of all-composite inlays, onlays, and individual crowns not requiring a framework. Other applications include the customization of plastic or ceramic stock teeth and temporary restorations made of Protemp®, manufactured by 3M ESPE.

Because of its special consistency, Sinfony can be used in a layering technique that is very similar to that used with ceramics. Because intermediate curing is possible, individual colors can be applied and fixed as desired.

The system consists of opaquer, dentin pastes, and incisal pastes for all 18 Vita-Classic shades, Intensive-Opaker, Opaque-Dentin, Transparent-Opal pastes, Magic colors, and accessories for processing. The combination of pastes and their designation are explained in the enclosed combination scheme.

All pastes are supplied in dispensers and are ready-for-use.

Sinfony Opaker is a **dual-curing** system consisting of a powder and a liquid component and can be cured with light in a rapid and reliable polymerization process. The ensuing process of self-curing ensures that the opaker hardens even in areas not exposed to light. A 0.1 mm layer of opaker is capable of forming a continuous coat.

Sinfony Intensive-Opakers are used for customizing veneers in cervical and interdental areas and extremely thin layers with an individual color. Intensive-Opaker IO1 pink is ideally suited for coating model cast retentions. The five Intensive-Opakers can be freely mixed with each other as well as with other opaques.

Sinfony Opaque-Dentin can be used for matching the color of extremely thin layers. The material can be freely mixed with all other Sinfony pastes. Five colors of Opaque-Dentin are available.

Sinfony Magic Intensive pastes are a set of 11 pastes for customized fabrication of restorations, denture teeth, and Protemp temporary restorations. These pastes can be used to produce a variety of effects, e.g. enamel cracks, calcareous spots, emphasis of cervical areas, and fissures. The 4 Transparent-Opal pastes are ideal for designing incisal edge translucence. Magic colors and Transparent-Opal pastes can be freely mixed with any Sinfony paste – with the exception of opaker.

Rocatec™ system, manufactured by 3M ESPE, is recommended as the adhesive system, as this system was specially designed for use in combination with Sinfony.

☞ The Instructions for Use of this product should not be discarded for the duration of product use.

2. Range of Applications

Sinfony:

Partial or complete veneering of

- crowns and bridges in repair applications
- telescope and cone crowns
- galvanic double crowns
- attachments
- implant supra-constructions
- adhesive bridges
- glass fiber-reinforced crowns and bridges

Fabrication of non-reinforced

- facings
- inlays and onlays
- all-composite single crowns

Customization of

- dentures
- denture teeth
- Protemp temporary restorations

Sinfony Opaquer is also used for:

- intraoral repair of ceramic or composite restorations with exposed metal parts (in combination with the Cojet™ system*, manufactured by 3M ESPE.

* Cojet is not available in some countries.

-
- ▶ Procedural step
 - 📌 Please note
 - 👉 Tip for processing
-

3. Preparation

Prepare/Design Metal Framework

Using Rocatec as the bonding system, the framework can be fabricated just like for veneering with ceramic. However, it is important to make sure to provide for a minimal wall thickness – see “Wall thickness for veneers”. The presence of an edge-to-edge bite, bruxism or other difficult bite situations exposing the restoration to special strains puts limitations on the fabrication with this technique. Under these circumstances, the framework must be shaped such that the veneer is sufficiently supported e.g. by a protective backing plate or some form of incisal edge protection.

However, if mechanical retentions are used instead of a bonding system, fabrication of the framework must comply with the guidelines of the composite veneering technique.

- ▶ Polish, clean, and de-grease the framework prior to veneering and treatment with Rocatec.

Bonding System

The recommended chemical adhesive system is Rocatec. This material bonds the veneering composite to the metal without leaving marginal gaps such that it is not required to use mechanical retentions.

The suitability of other silanization systems is beyond 3M ESPE's control. Therefore, 3M ESPE accepts no responsibility for the quality of bonding if materials other than Rocatec are used.

- 📌 For information on the Rocatec system, please refer to the respective Instructions for Use!

In the absence of a chemical bonding system:


- ▶ use common retention beads and shorten them to their equator;
- ▶ improve the mechanical adhesion by blasting surfaces to be veneered with aluminum oxide, grain size 100–250 µm;
- ▶ then clean the framework with alcohol in an ultrasonic bath.

Preparation of the Model

- ▶ Insulate all areas of the stone model exposed to Sinfony with Die Isolation, manufactured by 3M ESPE.
- ▶ Immediately re-seal the bottle after use.

- 📌 It is no problem if the liquid is turbid – this does not affect the insulation properties.

Assembly of the Dispenser

Dosing from the Dispenser


Option 1

- ▶ Dispense the paste from the dispenser by turning the front part of the dispenser clockwise, and apply directly.

Option 2


- ▶ Dose the paste onto a glass slab first, and then apply with a brush – size 1 – or a probe. Use the light protec-

tion lid to prevent premature curing of the paste on the glass slab.


Dosing Aid

The dosing aid can be used to consistently dispense the same quantity of the paste with each "click" of the aid. This mechanism facilitates the reproducible preparation of mixtures of identical composition.

Dispenser Holder

The Dispenser Holder can be used to keep up to 5 dispensers at hand with the color coding of the dispenser caps being well displayed.

Shade Selection

► Keep the suitable colors for the desired tooth shade at hand for the veneering procedure.

- ☞ Customized Sinfonycolor patterns are useful to select the tooth shade and all color shades directly in the patient. The Sinfony shade guide with single unfilled color trays is available from your distributor.
- ☞ Especially for fabrication of indirect inlays, it is recommended to generate pure dentin color patterns for shade selection directly in the patient.

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* Not available in all countries.

The assignment of the Opaque-Dentins can only be considered approximate, since the exact selection depends on the actual shade conditions:

- D01 + D02 for Vita-Classic shades A, B, and D, depending on brightness
- D03 for special orange effects, e.g. with shades from Ivoclar, Liechtenstein
- D02 + D04 for Vita-Classic shades, C and D, depending on brightness
- D05 for special dark gray-green effects, e.g. for older teeth
- ☞ Mixtures consisting of 1/3 Opaque-Dentin and 2/3 Dentin are useful for dental neck layering and coverage of small areas with small dabs of the material. Mixtures of this kind are less intensive and, thus, change the actual shade to a reduced degree.
- ☞ Customized shade trays simplify the selection of the Opaque-Dentins and assignment to the desired tooth shade. The Sinfony shade guide with unfilled color trays is available from your distributor.

4. Production of a Veneer

Silane the Rocatec-pretreated framework with ESPE™ Sil, manufactured by 3M ESPE. For further information, please refer to the Instructions for Use of the Rocatec system.

- Allow ESPE Sil to dry for **5 min**, and cover the framework with Sinfony Opaquer no more than **15 min** later.

Opaquer Processing

- Use the Combination Table to select the suitable opaquer powder for the desired color effect.
- Shake up the powder.
- Dose both the powder and the liquid into a groove in the ceramic plate at a weight ratio of 1:1.
 - ☞ Sinfony Opaquer can also be freely dosed. Overdosing the powder – producing a viscous-creamy consistency – accelerates, and underdosing the powder – producing a thin-creamy consistency – slows down the curing process.
 - Caution:** if the powder is strongly underdosed the opaquer will not harden!
 - ☞ For dosing the liquid, hold the bottle upside down in a vertical position and wait several seconds. This allows evenly-sized droplets to be dispensed from the bottle.
- ☞ Because the opaquer liquid is light-sensitive and volatile, the bottle must be immediately resealed after use!
- ☞ Do not use coated mixing pads, since any material abraded from the pad impairs the bonding to metal.
- ☞ Do not stir the opaquer on a flat surface to prevent the liquid from evaporating too rapidly, which would adversely affect the product properties.

- Use a plastic spatula and stir for **45 sec (!)**. The processing time is **5 min at 23°C/74°F** room temperature.
 - ☞ You must strictly comply with the specified mixing time to allow the catalyst components of the powder to completely dissolve in the liquid.
 - ☞ Do not use metal spatulas, since abraded material may cause changes in color.
 - ☞ Opaquer that has commenced to harden must not be diluted any further with opaquer liquid, since already formed polymer chains do not become re-dissolved. Doing so would adversely affect the bonding properties.
- Apply a thin layer of the opaquer with a clean disposable brush, made on behalf of 3M ESPE – similar to the procedure used for wash-baking ceramic veneers. Then use the Visio™ Alfa device, manufactured by 3M ESPE, for intermediate polymerization before subsequently applying a covering layer.
- For application and light-curing of the opaquer layer
 - choose option 1, if the framework is so small that the application time per each unit plus 10 sec of intermediate polymerization per each surface, taken together, does not exceed 5 min, which corresponds to the opaquer's processing time.
 - or
 - choose option 2, if the framework is so large that the application time per each unit plus 10 sec of intermediate polymerization per each surface, taken together, exceeds the 5 min of opaquer processing time.

Option 1

Apply the opaquer successively to each veneer surface and light-cure each surface separately for **2 x 5 sec** using the Visio Alfa device or for 10 sec using the Elipar™ unit manufactured by 3M ESPE. Repeat this sequential procedure until all veneer surfaces are completely covered by opaquer.

- ☞ Hold the opaquer surface as close as possible to the light exit, without directly contacting it, to ensure rapid and complete curing of the material.

Option 2

Apply the opaquer evenly to all veneer surfaces and polymerize by exposing to light for **7 min plus 10 sec of light + vacuum** using the Visio™ Beta Vario, manufactured by 3M ESPE.

- ☞ Polymerization with 1 min of light exposure plus 14 min light + vacuum would eliminate the inhibition layer and thus significantly reduce the strength of the chemical bond formed with the subsequent layer.
- Subsequently, apply another layer of opaquer ensuring that all metal is completely covered. Then repeat intermediate polymerization again according to option 1 or 2.
- Check that the entire veneer surface hardens properly. This can be done with a brush handle made of wood.
 - ☞ Make sure to touch the veneer surface only very gently, carefully avoiding to push the brush all the way through to the metal.

- ▶ Re-expose to light, if any opaquer adheres to the handle of the brush.
 - ❗ Incompletely cured opaquer areas will continue to harden with the subsequent Sinfony layer, but may cause the formation of streaks.
 - ❗ At this stage of the process, the opaquer surface must not be touched with the hands as this may adversely affect the chemical bond formed between the opaquer and the composite.
- ▶ Layering should be commenced once the veneer surface is completely covered and all opaquer completely polymerized.
 - ❗ After polymerization, the opaquer surface should not be left uncoated without a veneer for more than 1-2 hours.
 - ❗ Extended storage of the air-exposed opaquer surface without a veneer coat causes the reactive components to dry out and may lead to contamination and ultimately cause the quality of the bond between opaquer and Sinfony pastes to be reduced!

Intensive-Opaquer Processing

The translucence of this material is lower due to its large fraction of pigments. For this reason, mixtures containing Intensive-Opaquer must be given more time to polymerize.

- ❗ Intensive-Opaquers IO2-IO5 must not be processed!
- ▶ Add small amounts of Intensive-Opaquer (with the tip of the brush) to a small portion of ready-mixed opaquer of the respective tooth shade, until the individually desired color shade is attained.
- ▶ Cover the framework with opaquer using the common procedure, and apply the Intensive-Opaquer mixture especially in those places, in which customized staining is required.
- ▶ Mixtures with low Intensive-Opaquer content should be polymerized with the Visio Alfa device for **10 sec** = 2 exposure cycles.
- ▶ Mixtures with high Intensive-Opaquer content should be polymerized with the Visio Alfa device for **30 sec** = 6 exposure cycles, or with the Visio Beta Vario device using **7 min of light exposure and 10 sec of light plus vacuum**.

Intensive-Opaquer IO1 Pink is ideal for coating model cast retentions for which purpose it is mixed like normal opaquer, and then applied in the pure state.

- ▶ Polymerize Intensive-Opaquer IO1 Pink in the Visio Beta Vario device with **7 min of light exposure and 10 sec of light plus vacuum**; this prevents the formation of streaks in the denture plastics.


Opaque-Dentin Processing

Sinfony Opaque-Dentins can be combined with all Sinfony pastes.


- ❗ Because of their higher opacity, it is important to apply the pastes in layers of max. 1 mm to ensure complete curing of the material.

Magic Intensive Color Processing

- ❗ Mixtures with a high content of Sinfony Magic pastes are characterized by lower translucency! Therefore, mixtures of this type require 2-3 x longer light exposure = **10-15 sec** for polymerization with the Visio Alfa device.
- ▶ For customization of the veneer, dose Sinfony Magic pastes onto a glass slab.
- ▶ Use an instrument to mix the colors and the dentin pastes or apply very thin layers of max. 0.5 mm of the pure paste.
 - ☞ Intensive colors can be used both in internal layers and on the surface.
 - ☞ Standard mixtures allow the different color shades to be reproducibly dosed through the use of a dosing aid.
 - ☞ Magic colors, I9 Gingiva, I10 Red, and I11 Terra-cotta, are particularly well suited for generation of a gingiva-imitation effect.

Building-up the Tooth Contours by Layering

- ▶ There are two options available for building-up the tooth contours

Option 1:

Apply the pastes from the dispenser directly to the framework.

Option 2:

Dose the required quantities of paste onto a glass slab, and then use a size 1 brush or a probe for modeling.

- ❗ If a brush is used, first wet the brush with activator liquid and then streak out the brush to make sure that no air bubbles are included.

- ▶ Build up the entire restoration surface by surface in layers of max. 1 mm thickness, and light-cure each layer separately using the Visio Alfa device.
- ▶ Corrections can be made with a scalpel.
- ▶ Reseal the dispenser immediately after use to prevent the formation of clogs!

Reapplication of the Dispersion Layer

The dispersion layer must be reapplied if it was inadvertently touched by hand or became otherwise contaminated.

- ▶ Clean the contaminated surface, if any, with a lint-free cloth.
- ▶ Apply Symphony Activator with a clean short brush and lightly massage the substance into the surface. Do not cure with light!
- ▶ Continue layering as described under "Building-up the tooth contours".

Notes on Building-up the Tooth Contours

- ⓘ The tooth contours should be built up without any major delays to ensure optimal bonding and prevent contamination.
- ⓘ Symphony polymerized with the Visio Alfa device must not be subjected to grinding. The mechanical strain of grinding may cause cracks to form causing the material to chip off.

Building-up the Tooth Contours in 2 Steps


- ⓘ Building-up the tooth contours in 2 steps – i. e. 1., dentin layering followed by polymerization in the Visio Beta Vario device, grinding down; and 2., incisal layering – should only be done in "emergencies". Even exercising greatest care in this procedure, the bonding strength between the vacuum-cured dentin and the subsequently applied layers is only approx. 80% (on average) of the usual strength of the material. For this reason, continuous layering until the final tooth contour is attained followed by final polymerization is always preferable, since it affords superior bonding quality!
- ▶ However, if there is an emergency and the incisal layering or shape corrections must be made after polymerization in the Visio Beta Vario device, it is essential to comply with the procedure described under "Correction after shaping"!

Color Control


Sinfony polymerizes only if exposed to light. The yellow color of the polymerization photoinitiator is bleached only during the final polymerization in the Visio Beta Vario. Therefore, a definitive color control can be done only after final polymerization.

Polymerization

- ▶ Each layer must be polymerized separately (intermediate polymerization) for **5 sec** with the Visio Alfa device.
- ⓘ Applying and curing excessively thick layers without polymerizing the material in between causes tension to form, which may lead to chipping!
- ⓘ During intermediate polymerization, a dispersion layer remains on the composite's surface which is required for bonding to the subsequent layer. This dispersion layer must not be removed!
- ⓘ The quartz rod of the Visio Alfa device should be regularly cleaned. The cleaning procedure is described in the Operating Manual of the Visio Alfa device.

- ▶ Final polymerization is done in the Visio Beta Vario device without a stone model, using **1 min of light exposure plus 14 min of light plus vacuum**.
- ⓘ The photoinitiator loses its color during polymerization, and the Sinfony material attains the desired translucency.
- ⓘ The polymerization chamber of the Visio Beta Vario device must be kept clean. The reflector located behind the lamps should be regularly cleaned. The cleaning procedure is described in the Operating Manual of the Visio Beta Vario device under "Maintenance".
- ☞ If the polymerization chamber is strongly contaminated, soak respective areas in isopropanol or ethanol keeping the lid closed, and then remove the contamination with a cloth.
- ⓘ Acetone or acetic acid esters should not be used as these agents dissolve the rubber gasket.

Shaping and Polishing

Commence shaping and polishing with rotary instruments only after completing the final polymerization in the Visio Beta Vario device. Provide for aspiration during shaping and polishing!

- ⓘ The minimal thickness of the veneer should not be less than 0.8 mm.
- ⓘ Veneers on masticating surfaces should not be less than 1 mm thick (contra-indication).

For optimal shaping of the Sinfony material, the master set contains a shaping set of grinding and polishing tools manufactured by Hager & Meisinger – see package for ordering number – and Opal polishing paste from Renfert – ordering no. 520. These products can be procured directly from the manufacturers. The procedure described in the following can be used to produce a high gloss restoration.

Suitable Procedure using the Tools Provided in the Shaping Set

- ▶ Shape contours with cross-cut burs.
- ▶ Subsequently, smoothen the surface with a brown rubber polisher first, then with a green rubber polisher.
 - ⓘ The green rubber polishers are especially suitable for use at the transitions to the metal.
- ▶ Pre-polish all surfaces with the white sponge buff with no polishing agents.
 - ⓘ Do not use pumice powder!
- ▶ Use goat's hair brush and e.g. Acrypol (Bredent) for the interdental spaces.
- ▶ Finally polish the entire veneer to high gloss using the small cotton buffing disk and Opal high-gloss polishing paste.
- ▶ After completing polishing, clean the restoration with a medium hardness brush under tap water, or in an ultrasonic bath.
 - ⓘ **Do not use steam jet cleaning, since the different heat expansion coefficients of metal and composite may cause tension to form, which may lead to chipping!**

5. Correction/Repair

Correction after Shaping

- ⓘ In this context, please refer to the section, "Building-up the tooth contours in 2 steps"!
- ▶ The surface needing correction should be roughened, and oil-free air blown over the surface.
 - ⓘ Do not use steam, and do not clean with water!
- ▶ Apply ESPE Sil and allow to dry for 5 min.
- ▶ To reestablish a dispersion layer apply Sinfony Activator with a short clean brush and lightly massage the substance into the surface. Do not light-cure!
- ▶ Remove any excess material with a lint-free cloth, without wiping.
- ▶ Then layer as usual, polymerize individual layers with the Visio Alfa device, and conduct the final polymerization in the Visio Beta Vario. Then shape and polish.
 - ⓘ Subsequently applied material must not be ground and polished until after final polymerization in the Visio Beta Vario device.
- ▶ After polishing clean under tap water; do not use steam jet cleaning!

Repair of Worn or Fractured Restorations

- ▶ The surface needing correction should be roughened, and oil-free air blown over the surface.
 - ⓘ Do not use steam, and do not clean with water!
- ▶ Coat the area needing to be repaired with Rocatec Plus; for details please refer to the Instructions for Use of the Rocatec system.
- ▶ Apply ESPE Sil and allow to dry.
- ▶ Layer Sinfony as usual, polymerize individual layers with the Visio Alfa device, and conduct the final polymerization in the Visio Beta Vario. Then shape and polish.
 - ⓘ Subsequently applied material must not be ground and polished until after final polymerization in the Visio Beta Vario device.
- ▶ After polishing clean the work under tap water; do not use steam jet cleaning!

6. All-Composite Restorations

Guidelines for Preparatory Work

- Pronounced chamfer or circular shoulder
- Wall thickness must be at least 0.8 mm circular; 1.0-1.3 mm occlusal; 1.0 mm for veneers/facings
- Isthmus of inlays should be 1.5 mm (2.0 mm are preferable)
- Cavity walls 4° divergence versus occlusal

Fabrication of Full-Composite Inlays, Onlays, Crowns, Veneers/Facings

- ▶ To provide for precise control at the contacting points – especially in inlay and onlay fabrication – manufacture an un-sawn control model in addition to the master model.
- ▶ Block out undercuts of the master model with a light-colored or colorless wax.
 - ❏ Colored wax changes the color effect of the translucent composite.
 - ❏ Do not use light-curing materials, as these would bond to Sinfony.
- ▶ Mark the preparation limit.
- ▶ Wet the cavity/die twice with Ivoclar Separator* and then blow air over the surface.
 - * Ivoclar Separator is manufactured by Ivoclar, Liechtenstein
 - ❏ There must not be any pools remaining.
- ▶ Have at hand the required quantities of Sinfony pastes to be dispensed onto a glass slab.
 - ☞ During shaping, use the light protection lid to protect the pastes from exposure to light.
- ▶ Build-up the restoration with common instruments proceeding in layers of max. 1 mm thickness.
- ▶ Polymerize each layer separately with the Visio Alfa device for **5 sec.** Each surface of multi-surface restorations should be light-cured separately.
 - ❏ Layers in excess of 1 mm thickness may cause tension to form in the composite.
 - ☞ Cusps and ridges are easy to apply with thin probes or nerve broaches.
 - ☞ Coat amalgam-stained stumps with Opaque-Dentin; for standard layering use dentin, incisal pastes, and a sufficient quantity of Transpa pastes.
 - ☞ To provide for more stability especially for lifting-off MOD inlays, shape the approximal edges somewhat thicker without extending the length (!).
- ▶ Correct early contacts in the state of intermediate polymerization with a scalpel, do not subject the restoration to grinding!
- ▶ Lift the restoration off the model.
- ▶ Conduct the final polymerization in the Visio Beta Vario using **1 min of light exposure plus 14 min of light plus vacuum.**
- ▶ Shape according to the common procedure for veneers; please refer to “Shaping”.

Preparatory Steps in the Laboratory for Trying-in the Restoration

As a lift-off aid for the dentist, a sphere made from composite can be polymerized onto the restoration, which simplifies the trying-in and is easy to remove at a later time.

Preparatory Steps for Cementation in the Dentist's Office

Option 1

- ▶ Use a fine bur to roughen the basal surface of the restoration.
- ▶ Apply some composite bonding agent in accordance with the manufacturer's recommendations.

Option 2

- ▶ Coat the basal surface of the restoration with Cojet.
- ▶ Then apply ESPE Sil and allow to dry.
 - ❏ For information on the processing of 3M ESPE Cojet and ESPE Sil, please refer to the respective Instructions for Use!

Adhesive Cementation of Full-Composite Restorations

- ▶ Use an adhesive cement for cementation, e.g. Compolute™, manufactured by 3M ESPE. For information on the processing, please refer to the respective Instructions for Use!

7. Customization of Denture Teeth/ Protemp 3 Garant Temporary Restorations

Denture Teeth

- ▶ Roughen the surface of plastic or ceramic denture teeth with Rocatec Pre or rotary instruments.
- ▶ Pretreat with Rocatec Plus and ESPE Sil; please refer to the Instructions for Use of the Rocatec system.
- ▶ Layer Sinfony according to the common technique, then polymerize and shape.
 - ❏ For aesthetic reasons, it is recommended to coat the entire denture tooth with a thin layer, because transitions between different materials would otherwise become noticeable.

Protomp 3 Garant Temporary Restorations

Procedure for New Temporary Restorations

Sinfony can be applied equally well to the shaped surface or the unshaped clean inhibition layer of new restorations. Very good bonding is attained in either case.

- ▶ Remove any contamination – e.g. saliva or polishing dust – with alcohol, then dry with air. If the inhibition layer is clean or if the temporary restoration is shaped and cleaned, proceed as described in the following step.
- ▶ Used in the dentist's office, Sinfony should be dosed from the dispenser onto a mixing pad for hygienic reasons.
- ▶ Continue processing as described under “Building-up the tooth contours by layering”, “Polymerization”, and “Shaping and polishing”.

Procedure for Used Temporary Restorations

- ▶ Thoroughly roughen the surface of the temporary restoration with rotary instruments.
- ▶ Remove any contamination – e.g. saliva or polishing dust – with alcohol, then dry with air.
- ▶ Used in the dentist's office, Sinfony should be dosed from the dispenser onto a mixing pad for hygienic reasons.
- ▶ Continue processing as described under "Building-up the tooth contours by layering", "Polymerization", and "Shaping and polishing".

8. Customization of Dentures

Complete Dentures

Preparation of the Denture

- ▶ Fabricate the denture as usual.
- ▶ Areas to be customized should be thoroughly ground down or roughened.
- ▶ Pretreat with Rocatec Plus or Soft and ESPE Sil.

Customization of the Denture

Dentures can be customized as described for the following example:

- ☞ use I1 Ocean blue to imitate the bluish hue of interdental papillae,
- ☞ use I10 Red to reinforce the red color of gingiva,
- ☞ use I7 Vanilla for a hint of bony substance showing through,
- ☞ use I3 Chestnut or I11 Terracotta for the brown-red of gingiva,
- ☞ use I9 Gingiva to complete modeling.
- ▶ Select the desired colors for the individual patient.
 - ☞ All colors can be mixed with T1 neutral to generate delicate color effects.
 - 🔑 Use an instrument to mix the colors and make sure not to introduce any air bubbles.
- ▶ Apply the paste in layers of no more than 1 mm; apply progressively less material near the transition to the denture plastic.
- ▶ Light-cure each layer separately for **5 sec** with the Visio Alfa device.
- ▶ Once all layers of the denture are built-up, polymerize for **1 min of light exposure plus 14 min of light plus vacuum** using the Visio Beta Vario device.
- ▶ Shape with cross-cut burs and rubber polishers.
- ▶ Use sponge buff, goat's hair brush, and polishing paste to smoothen the surface.
- ▶ Use a cotton buffing disk and high-gloss polishing paste to polish to high gloss.

- ▶ After polishing, clean the denture with a brush and tap water or in an ultrasonic bath.
 - 🔑 Do not use steam jet cleaning!

Model Cast Dentures

- ▶ Apply a thin, but completely coating, layer of Opaquer pink to the retentions, and light-cure the opaquer for **7 min of light exposure plus 10 sec of light plus vacuum** using the Visio Beta Vario device.
 - ☞ If a Visio Beta device is not available for polymerization, bead polymer can be sprinkled onto the pink opaquer prior to completing the restoration with denture plastic. This allows the denture plastic to form a proper transition and prevents the formation of streaks.

9. Veneering of Glass-Fiber Reinforced Frameworks

Vectris®

Framework Preparation for max. 3-Unit Bridges with Sinfony Veneering:

- ▶ Fabricate the glass-fiber reinforced framework in accordance with the manufacturer's recommendations.

Preparation of the Framework:

Option 1:

- ▶ Blast the framework with 120 µm aluminum oxide at a max. pressure of 1 bar and then clean with steam.
- ▶ Apply Vectris wetting fluid and allow to act for 60 sec.
- ▶ Dry excess fluid with oil-free air.

Option 2:

- ▶ Blast the framework with Rocatec Soft in the Rocatecor™, manufactured by 3M ESPE.
- ▶ Then apply ESPE Sil and allow to dry for 5 min.

Veneering

- ▶ Use Sinfony Opaque-Dentin to fabricate both the crown margin and a thin coating over the entire framework. Individual layers must not exceed max. 1 mm in thickness, and need to be light-cured separately.
- ▶ Then build-up the veneer as usual, cure with light, and shape.

Stick Tech®

Framework Preparation for max. 3-Unit Bridges with Sinfony Veneering:

- ▶ Allow Sinfony Activator to permeate into Stick Tech® products prior to processing.
- ▶ Fabricate the glass fiber-reinforced framework in accordance with the manufacturer's Instructions for Use.
 - ☞ For adaptation of the glass fiber nets to the die, it is recommended to use transparent deep-drawing caps so that the net can be light-cured through the cap.

Preparation of the Framework

- The areas to be coated with veneer should be wetted and then wiped clean with a lint-free cloth. There must not be any pools of the substance.

Veneering

- Use Sinfony Opaque-Dentin to fabricate both the crown margin and a thin coating over the entire framework. Individual layers must not exceed max. 1 mm in thickness, and need to be light-cured separately.
- Then build-up the veneer as usual, cure with light, and shape.

FibreKor®

Framework Preparation for max. 3-Unit Bridges with Sinfony Veneering:

- Fabricate the glass fiber-reinforced framework in accordance with the manufacturer's Instructions for Use.
- Completely envelope the glass fibers with Sinfony making sure that the crown margin remains free of glass fibers.

Preparation of the Framework

- In framework fabrication with Sinfony there are no further preparatory steps required.

Veneering

- Build-up the Sinfony veneer as usual.

DC-Tell®

Framework Preparation for max. 3-Unit Bridges with Sinfony Veneering:

- Fabricate the glass fiber-reinforced framework in accordance with the manufacturer's Instructions for Use.

Preparation of the Framework

- Blast the framework with Rocatec Soft in the Rocatector.
- Then apply ESPE Sil and allow to dry for 5 min.

Veneering

- Use Sinfony Opaque-Dentin to fabricate both the crown margin and a thin coating over the entire framework. Individual layers must not exceed max. 1 mm in thickness, and need to be light-cured separately.
- Then build-up the veneer as usual, cure with light, and shape.

Vectris® is a registered trademark of Ivoclar, Liechtenstein
Stick Tech® is a registered trademark of Stick Tech Ltd. Oy, Finland
Fibrekor® is a registered trademark of Jeneric Pentron, USA
DC-Tell® is a registered trademark of DCS Dental AG, Allschwil, CH

10. Preventing Processing Errors

Framework Design

- ❗ The framework must provide sufficient support to the veneer.

- ❗ Difficult bite conditions can cause the veneer to shear off – partially or in full – if the framework design is faulty.

Opaquer

- ❗ Always apply the opaquer in very thin layers to ensure complete curing. If the opaquer layer is too thick, there may be a loss of bonding inside the opaquer layer.
- ❗ Do not use opaquer liquid that has become thickly viscous, because the material affords only reduced bonding strength at this consistency.
- ❗ The completely cured opaquer must not be touched by hand or wiped off, since this would adversely affect the chemical bond.
- ❗ After application of the opaquer do not allow the restoration to remain without a veneer for longer than 1 hour.

Layering

- ❗ Apply Sinfony in individual layers of max. 1 mm thickness, and polymerize layers separately to make sure that the material hardens completely.
- ❗ Multi-surface restorations and complete-veneer restorations should be built-up one surface at a time, and individual surfaces light-cured separately.

Layering Scheme

- ❗ Building-up the veneer make sure that the layer thickness is sufficient and a correct ratio of dentin paste to incisal paste is established in order to provide for an optimal color effect.

Polymerization

- ❗ Under no circumstance should the polymerization time with the Visio Alfa device be shorter than 5 sec per layer. Careful intermediate polymerization minimizes shrinking effects during polymerization and provides for a veneer with no internal tension.
- ❗ Hold the restoration as close as possible to the light exit to ensure optimal curing of the material.
- ❗ Keep the light exit of the Visio Alfa device and the polymerization chamber of the Visio Beta Vario device clean at all times, since this is a prerequisite of obtaining optimal light power.
- ❗ Under no circumstance should the polymerization time in the Visio Beta Vario device be shortened, since only the completely cured material possesses optimal material properties.

Wall Thickness for Veneers

- ❗ To prevent fracturing, the following wall thicknesses must be provided:
Veneers not involved in occlusion: 0.8 mm
Veneers of masticating surfaces 1.0 mm–1.3 mm
Facings 1.0 mm

Veneer Cleaning

- ⚠ Do not use steam jet cleaning, since the difference in heat expansion coefficients of metal and composite may cause tension in the bond to form which may ultimately lead to chipping.

11. Restoration and Denture Care

- ▶ Thoroughly brush restorations and dentures with a medium hardness toothbrush and a gel toothpaste with low abrasiveness.
 - ⚠ Strongly abrasive toothpaste should not be used as it may attack the polish.
 - ⚠ Do not place dentures in cleaning solutions, because the aggressive active ingredients of these solutions may attack the composite and destroy the polish. Although the bacteria present in the plaque are killed, cleaning solutions are not sufficient to remove all the existing plaque. Spots with remaining plaque are predilection sites for new plaque to form.
 - ⚠ Care instructions for patients can be obtained from 3M ESPE free of charge!

12. Repair of Fixed Restorations

- ▶ Repair defective fixed restorations with a filling composite and the suitable bonding agent.
 - ⚠ Repairs utilizing the Cojet system are particularly stable.

13. Notes

Sinfony veneer material must be polymerized with the Visio Alfa and the Visio Beta Vario light units. Since the suitability of other light units is beyond 3M ESPE's control, 3M ESPE accepts no liability for the quality of polymerization in cases, in which other light units are used.

14. Incompatibilities

In susceptible individuals, sensitization to the product cannot be excluded. Use of the product should be discontinued and the product completely removed, if allergic reactions are observed.

15. Storage and Shelf-life

Do not store Sinfony Pastes and Opaquer above 25°C/77°F. Do not use beyond the expiration date.

Large quantities of opaquer powder and liquid should be stored in a cool place, if they are to be stored for extended periods of time, e.g. in warehousing.

16. Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Index

Adhesive cement 8
 Blocking out 8
 Bonding quality 2
 Bonding strength 6
 Bonding system/Rocatec 2
 Bruxism (teeth grinding) 2
 Burs for shaping 7
 Chipping 6
 Coating of model cast retentions 5
 Color table 3, 6
 Control model 8
 Customized color trays 3, 4
 Denture cleaner 11
 Die isolation 2
 Difficult bite situation 2, 10
 Dispenser 2, 3
 Dispersion layer, reapplication 6
 Dispersion layer, smear layer 2, 6, 7, 10
 Dosing 2, 3
 Dosing aid 2
 Edge-to-edge bite 2
 Framework design 2
 Gel toothpaste 11
 Gingiva imitation effect 5
 Glass slab 2
 Heat expansion 11
 Indications 2
 Intensive pastes 1, 5
 Lift-off aid 8
 Light protection lid 2
 Magic color mixtures/pure usage 5
 Material strength 6
 Mechanical retentions 2
 Mixing – opaquer, pastes 1, 2, 4, 5
 Mixing spatula 4
 Mixtures containing Intensive-Opaquer 5
 Opaquer powder/liquid 4
 Other silanization systems 2
 Patient information 11
 Photoinitiator 6
 Polishing tools 7
 Polymerization chamber cleaning 6
 Pumice stone powder 7
 Reapplication of Sinfony 7
 Reflector cleaning 6
 Reproducible preparation of mixtures 3
 Shade guide, unfilled 3
 Shrinking during polymerization 10
 Silanization 4
 Stability for lifting-off 8
 Standard mixtures 5
 Steam jet 7
 Storage/Warehousing 11
 Streak formation 5
 Tension 6
 Thickness of layers 5
 Transitions between materials 8

Veneer support 2
 Veneer thickness 7, 10
 Visio Alfa cleaning 6
 Visio Beta Vario cleaning 6
 Vita-Classic shades 3
 Wall thickness of full-composites 7

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3M, ESPE, Cojet, Compolute, Elipar, Protemp, Rocatec, Rocatector, Sinfony and Visio are trademarks of 3M or 3M ESPE AG.
© 3M 2002. All rights reserved.