

Otto Bock®

QUALITY FOR LIFE

757L20

(D)	Otto Bock® Li-Ion Ladegerät	3
(GB)	Otto Bock® Lithium-Ion Battery Charger	6
(F)	Chargeur Li-Ion Otto Bock®	9
(I)	Caricabatteria al litio-ionio Otto Bock®	12
(E)	Cargador Li-Ion de Otto Bock®	15
(P)	Carregador de Acumuladores de Iões de Lítio Otto Bock®	18
(NL)	Otto Bock® Li-ion oplaadapparaat	21
(S)	Otto Bock® Li-Ion Laddningsaggregat	24
(DK)	Otto Bock® lithium-ion ladeapparat	27
(PL)	Ładowarka akumulatorów litowo-jonowych firmy Otto Bock®	30
(H)	Otto Bock® Li-ion akkutöltő	34
(CZ)	Nabíječka Li-Ion akumulátorů Otto Bock®	37
(GR)	Φορτιστής συσσωρευτών ιόντων λιθίου της Otto Bock®	40
(中国)	奥托博克锂电池充电器	44

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Before operating the Lithium-Ion Battery Charger, carefully read these operating instructions as well as the operating instructions for the 757B20/757B21 Otto Bock EnergyPack. Please pay particular attention to the Safety Precautions!

1 Application

The Otto Bock 757L20 Lithium-Ion Battery Charger is to be used **solely** for charging one or two 757B20/757B21 Otto Bock EnergyPacks.

2 Safety Precautions

Attention!

- Use the Charging Station only in combination with the 757L16-1 Power Supply Unit or with the 4X74 12V Charging Cable.
- Use the battery charger under dry conditions only.
- Do not use the battery charger if any damage is visible. If you discover any damage, send the battery charger immediately to an Otto Bock service facility.
- Do not allow objects or liquids to fall onto or spill into the battery charger.
- Do not expose the battery charger to excessive heat (avoid direct sunlight), excessive smoke, mechanical vibrations or shock.

3 Scope of Delivery

- 1 pce. Charging Station
- 1 pce. Power Supply Unit with EU plug
- 1 pce. US plug
- 1 pce. Operating Manual

4 Accessories

For using the battery charger with the 12V cigar lighter of a car, the 4X74 Charging Cable is optionally available.

5 Description and Function

Unlike NiCad batteries, the rechargeable lithium-ion battery EnergyPack has the benefits of a lower self-discharge, a higher cell potential and a higher energy density. Lithium-ion technology offers higher capacities, has no memory effect and requires considerably reduced charging times.

The 757L20 Otto Bock Lithium-Ion Battery Charger consists of a Charging Station (fig. 1, pos. 1), a 757L16-1 Power Supply Unit with EU plug (fig. 1, pos. 2), and an additional US plug (fig. 1, pos. 3). For exchanging the plug, please refer to figure 4.

The battery charger has an operating indicator (green LED, fig. 2, pos. 3) as well as four other LED indicators:

- Green LEDs (fig. 2, pos. 1 and 1a): Indicate the charging status of the EnergyPacks
- Yellow/red LEDs (fig. 2, pos. 2 and 2a): Indicate malfunctions of the Energy Packs or battery charger

Before charging, the battery charger conducts a self-test.

In a test procedure, short-circuited or damaged cells are identified and prepared for the CCCV process (constant current - constant voltage) during the conditioning phase of the charging process.

During the CCCV process, the EnergyPacks are charged with a continuously high current until a certain voltage level is reached. The current then automatically decreases to a lower level until it is switched off completely.

If the voltage of the EnergyPacks falls below a certain level, the battery charger automatically switches to trickle charge.

The charging process will also be interrupted when the temperature exceeds 50 °C/122 °F or falls below 0 °C/32 °F. This is indicated by the yellow LED. The charging process will resume as soon as the EnergyPack has reached room temperature.

6 Function

6.1 Self-Test

Connect the Charging Station to the 757L16-1 Power Supply Unit (see fig. 2, pos. 4) and plug into an electrical outlet. The operating indicator lights up (see fig. 2, pos. 3) showing the battery charger is ready for use.

Next the battery charger conducts a self-test, during which all other indicator lights will light up briefly.

6.2 Charging

Insert one or two EnergyPacks into the charging station (fig. 3).

The charging process begins when the indicator lights 1 and/or 1a light up.

The charging time for a depleted EnergyPack is approximately 3-1/2 hours. As soon as the EnergyPack is fully charged, the 1/1a indicator light(s) flash and the charging current is cut off.

If the EnergyPack voltage falls below a certain value due to self-discharge, the Otto Bock Lithium-ion Battery Charger automatically switches to trickle charging. The fully charged battery may be left in the battery charger without becoming damaged.

Note:

To avoid malfunctions, the electrical contacts of the EnergyPack and of the battery charger must be kept clean and free of grease. For cleaning please use a dry cotton pad stick.

6.3 Trouble-Shooting Tips

- The lithium-ion battery charger is plugged into an electrical outlet, but the operating indicator is not lit: ► Re-test charger function **at a different electrical outlet**.
- If the operating indicator still does not light up, the lithium-ion battery charger is defective and must be sent to an Otto Bock service facility.
- If the indicator lights 2 and 2a alternately flash red, the battery charger is defective and must be sent to an Otto Bock service facility.
- The lithium-ion battery charger is plugged into an electrical outlet, the operating indicator is lit, and two EnergyPacks are inserted, but only one of the indicators 1 or 1a is illuminated:

► Switch the locations of the EnergyPacks:

- a) If the same LED does not light up, the battery charger is defective and must be sent to an Otto Bock service facility for repair.

b) However, if the other LED now fails to light, the corresponding EnergyPack is defective and must be sent to an Otto Bock service facility for repair.

Other messages:

Battery Compartment	Indicator 1/1a green	Indicator 2/2a yellow	Indicator red	Function	Remarks
empty	off	off	off	okay	stand-by
empty	off	off	flashing	error	self-test – return charger
with EnergyPack	off	off	flashing	error	self-test – return charger
with EnergyPack	on	off	off	okay	charging
with EnergyPack	blinking	off	off	okay	EP is fully charged
with EnergyPack	off	on	off	error	temperature of EP > 50 °C: let cool down (see safety precautions)
with EnergyPack	off	off	on	error	EP is defective – return EP
empty	off	on	off	error	charger defective – return charger
empty	off	off	on	error	charger defective – return charger
empty	on	off	off	error	charger defective – return charger

7 Technical Data

Charging Time..... approx. 3-1/2 h
 Power supply..... 100 - 240 V
 Power supply frequency 47 - 63 Hz
 Operating temperature 0 - 50 °C (32-122°F)

8 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorized by the manufacturer (see instructions for use and catalogs). The manufacturer does not assume liability for damage caused by component combinations which it did not authorize.

The product may only be opened and repaired by authorized Otto Bock technicians.

9 Declaration of Conformity

Otto Bock as manufacturer with sole responsibility declares that the 757L20 Otto Bock Lithium-Ion Battery Charger conforms to the requirements of the 93/42 EEC Guidelines.

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Utility patent declared in ET, B, NL, L, D, F, GR, RI, I, MA, NA, P, E, TN, SCV.

Design Patent in China no. 99339427.2, in Danmark no. 2000 00795, in United Kingdom no. 2087 765, in Finland no. 2 1303, in Japan no. 1 094 598, in Canada no. 92 760 D , in Austria no. 38340, in Sweden no. 66 319, R.O.C. Design Patent no. 073025, in the USA no. D 438 170

Otto Bock Healthcare Products GmbH

Kaiserstraße 39 · 1070 Wien.Austria · Tel. +43 1 5269548 · Fax +43 1 5267985

e-mail: vertrieb.austria@ottobock.com · www.ottobock.at

Otto Bock has been certified by the German Society for the Certification of Quality Assurance Systems (DQS) in accordance with DIN EN ISO 13485, Reg. No. 070813 (Management System)