## Français - FR

Code technique : ST

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English - GB

### In vitro diagnostic reagent, for professional use only

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### CALIBRATION

ELITE Cholesterol LDL 2G Calibrator should be used for calibration. Its value is traceable to the reference method recommended by the CDC (Centers for Disease Control and Prevention).

The calibration frequency is specific for each analyser. (refer to § PERFORMANCE DATA).

### QUALITY CONTROL

Control sera ELITROL I and ELITROL II must be used to ensure adequate quality.

### PERFORMANCE DATA at 37 °C on analyser ELITech Clinical Systems Selectra XL

#### - Analytical range

Determined according to EP6-A protocol of CLSI.<sup>(6)</sup>

The reagent is linear from 15 to 700 mg/dL (0.39-18.10 mmol/L).

#### - Detection limit

Determined according to SFBC protocol<sup>(7)</sup>, the detection limit is equal to 0.3 mg/dL (0.01 mmol/L).

#### - Analytical sensitivity

The average variation of the analytical signal is 2.04 mA/mmol/L for mg/dL of LDL cholesterol (79 mA/mmol/L) for a light path of 1 cm.

#### - Precision

Determined according to EP5-A2 protocol of CLSI.<sup>(8)</sup>

N	Mean		Within-run CV (%)	Total
	mg/dL	mmol/L		
Level 1	80	103	2.66	1.6
Level 2	80	123	3.18	0.8
Level 3	80	164	4.24	2.0

#### - Correlation

A comparative study has been performed against a concurrent system (enzymatic colorimetric method - Selective detergent) on 53 human sera including 19 serum with a triglyceride concentration between 200 and 850 mg/dL. The sample concentrations ranged from 55 to 262 mg/dL (1.42 to 6.78 mmol/L).

The parameters of linear regressions are as follows :

Correlation coefficient: (r) = 0.9961

Linear regression:  $y = 0.9927x + 1.2 \text{ mg/dL}$  (0.03 mmol/L)

#### - Interferences<sup>(5)</sup>

According to SFBC recommendations, some studies have been performed to determine the level of interference from different compounds :

#### - Conjugated bilirubin:

No significant interference up to 22 mg/dL (376.3 µmol/L).

#### - Unconjugated bilirubin:

Negative bias from 27.5 mg/dL (470.4 µmol/L).

#### - Turbidity:

No significant interference up to 600 mg/dL (6.78 mmol/L) Triglyceride equivalent.

#### - Haemoglobin:

No significant interference up to 500 mg/dL (5 g/L).

In very rare cases, monoclonal gammopathies (multiple myeloma), in particular of IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results.<sup>(9)</sup>

Other compounds may interfere.<sup>(10,11)</sup>

#### - On board stability / calibration frequency

On-board stability: 21 days.

Calibration frequency: 21 days.

Make a new calibration when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.

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2. Naito, H.K., *Coronary artery disease and disorders of lipid metabolism*, Clinical Chemistry: Theory, Analysis, Correlation, 4<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C. (Mosby, Inc. eds. St Louis USA), (2003), 603.

3. Tietz, N.W., *Clinical guide to laboratory tests*, 3<sup>rd</sup> Ed., (W.B. Saunders eds. Philadelphia USA), (1995), 404.

4. Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), *Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP)*, JAMA, (2001), 285, 2486.

5. Bachorik, P.S., Ross, J.W., *Clin. Chem.*, (1995), 41, 1414.

6. *Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach*, Approved Guideline. CLSI (NCCLS) document EP6-A (2003), 23 (16).

7. Vassault A., et al., *Ann. Biol. Clin.*, (1986), 44, 666.

8. *Evaluation of Precision Performance of Quantitative Measurement Methods*; Approved Guideline—Second Edition. CLSI (NCCLS) document EP6-A (2002), 24 (25).

9. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*. Acta Clin Belg., (2004), 59, 263.

10. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AAC Press, (1997).

11. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AAC Press, (1995).

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