GreenLight MoXy™ Fiber Optic

- Product Insert
- Fiber Replacement Policy

CAUTION: Federal law (U.S.) restricts this device to sale by or on the order of a physician.

Manufactured by:

American Medical Systems
Innovation Center - Silicon Valley
3070 Orchard Drive
San Jose CA 95134-2011
USA
Customer Care:
U.S. Toll Free: 800 356 7600
Tel: +1 408 456 3585
Fax: +1 408 428 0512
www.AmericanMedicalSystems.com

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Made in USA

GreenLight MoXy™ is a trademark of American Medical Systems.
1008002-01 (AW, Rev 01, 10/12)

GreenLight MoXy™ Fiber Optic

REF 10-2400

These instructions detail the care and use of the GreenLight MoXy™ Fiber Optic. These instructions are not recommendations for the medical or surgical application of this device. Any physician using this device should be thoroughly familiar with the surgical procedure being performed prior to using the MoXy fiber. The MoXy fiber is designed for use with the AMS GreenLight XPS™ Laser System. Refer to the GreenLight XPS Operator's Manual, Safety and Professional Information Sections, for specific instructions concerning warnings, cautions, contraindications and clinical use of the laser.

Package Contents

- Individually packaged, sterile, disposable GreenLight MoXy Fiber Optic
- Fiber card
- Fiber Replacement Policy

Description

The GreenLight MoXy Fiber Optic features a side firing mechanism delivering up to 180W of 532nm light to tissue. The Fiber is designed/intended to be delivered to the desired treatment site via an endoscope/cystoscope. The fiber can be used for the surgical incision/excision, vaporization, ablation, hemostasis and coagulation of soft tissue. All soft tissue is included, such as skin, cutaneous tissue, subcutaneous tissue, striated and smooth tissue, muscle, cartilage, meniscus, mucous membrane, lymph vessels and nodes, organs and glands.

The fiber is supplied as a single-use, sterile product and is compatible with the GreenLight XPS Laser System. It can access the tissue in multiple planes. The fiber is a liquid cooled delivery device which enables the delivery of up to 180 watts of power and aids in maintaining a clear environment at the fiber cap.

Intended Use


The GreenLight MoXy Fiber Optic features a side firing mechanism delivering up to 180W of 532nm light to tissue. The fiber can be used for the surgical incision/excision, vaporization, ablation, hemostasis and coagulation of soft tissue. All soft tissue is included, such as skin, cutaneous tissue, subcutaneous tissue, striated and smooth tissue, muscle, cartilage, meniscus, mucous membrane, lymph vessels and nodes, organs and glands.

The GreenLight MoXy Fiber Optic will deliver 532nm laser energy from a compatible laser console to tissue during surgical procedures, including photoselective vaporization of the prostate for benign prostatic hyperplasia (BPH).

Specifications

Overall Length	3050 mm
Operating Length	360 mm
Core diameter	750 microns
Fiber tip O.D.	2.1 mm typical, 2.3 maximum
Compatible device	Most 22 to 24 Fr continuous flow cystoscopes

1. Metal cap	6. Beam alignment indicator
2. Laser firing point	7. Fiber handle
3. Blue triangle. Red octagon located on opposite side	8. Inlet flow tubing
4. Outer flow tubing	9. Luer lock connector
5. Control Knob	10. Optical connector

Contraindications

The fiber and GreenLight XPS Laser System should only be used by a qualified and trained surgeon.

The use of the Laser System is contraindicated for patients: whose general medical condition contraindicates surgical intervention, where appropriate anesthesia is contraindicated by patient history, where tissue (especially tumors) has calcified, for hemostasis of vessels over approximately two millimeters in diameter, where laser therapy is not considered the treatment of choice, uncontrolled bleeding disorders and coagulopathy, prostate cancer, acute urinary tract infection (UTI) or severe urethral stricture.

The fiber and GreenLight XPS Laser System is contraindicated in the presence of severe urethral strictures; however, the system can be used in the treatment of urethral strictures with proper cautions. A severe stricture is any stricture with visible narrowing via urethrography or ultrasonography, with near total obstruction that makes passage of instruments difficult or dangerous. Care should be taken to avoid injury to urethral tissue.

Potential Complications and Risks

The same Potential Complications and Risks that exist for conventional or traditional surgery also exist for laser surgery. These include, but are not limited to, the following:

Non-Thermal Complications and Risks: allergic reaction to medication, arrhythmia, aspiration, gas over-distension, hypertension, induced hemorrhage, infection, pain, perforation or pneumothorax

Thermal Complications and Risks: chills, edema, fever, induced hemorrhage, delayed hemorrhage. Leukocytosis, pain, perforation, ulceration, delayed healing, embolism, sepsis or stricture

As with conventional endoscopic treatment, the possibility of adverse reactions such as fever, chills, sepsis, edema and hemorrhage may occur following laser treatment. In extreme cases, death may occur due to procedural complications, concurrent illness or the application of the laser. Use caution when treating patients who have had difficulty with previous endoscopic procedures.

Other complications may include: abdominal bloating (intestinal gas), allergic reaction, bladder neck contracture, bleeding, chills, clot retention, constipation, deep venous thrombosis, delay in healing, diarrhea, dizziness, dysuria, edema, epididymitis, erectile dysfunction (ED), fatigue or weakness, fever, fluid overload/Hyponatremia, frequency, hematospermia, hematuria, hemorrhage, infection, impotence, loss of appetite, malfunction of laser fiber or console resulting in an injury or prolonged procedure, leukocytosis, nausea, nocturia, pain (abdominal pain unresponsive to NSAIDs, arm or leg pain, headache, back/low back pain, body aches, pelvic, penile), pelvic hematoma, penile urethral injury, perforation/injury of prostate or adjacent organ(s), profuse perspiration, not fever related, prostatitis, pulmonary embolus, retrograde ejaculation, sepsis, tissue sloughing, transfusion, ulceration, urethral stricture, urgency, urinary frequency or vomiting.

Consult the GreenLight XPS Operator Manual for additional complications and risk by indication.

Instructions for Use

Information regarding WARNINGS and PRECAUTIONS are contained in the shaded boxes through the instructions.

Laser System Set-up

1. Prior to use of the Laser System, consult the GreenLight XPS Operator Manual for information and instructions for use. Ensure that the user, the patient and all personnel in the room are wearing appropriate protective eye wear prior to using the laser.
2. Turn the laser system ON. Wait for self-test to be completed and the prompt "Insert Fiber Card" to appear.
3. Remove the sterile pouch from the box. Do not open the pouch.
4. Remove fiber card from the pocket on the outside of the pouch. Insert card into reader with image facing operator.
5. After the card is recognized by the laser system, there will be a prompt to "Attach a Fiber."

- The fiber card must remain inserted until the procedure is complete. Premature removal of the card will result in the discontinued use of the fiber during the procedure.

Prior to the Procedure

1. Prior to use, inspect the fiber package for damage.

- Do not use if the package is damaged. Damage to the package may result in damage to the fiber or compromised sterility.
- The fiber is a precision device. Damage to the fiber can result in the emission of uncontrolled laser energy. Do not excessively manipulate or bend the fiber.

To prevent drape fires or burns

- Do not wrap any portion of the fiber in drapes or cloth.
- Do not attach the fiber directly to surgical drapes with a crushing clamp such as a hemostat.
- Do not place or drop instrumentation onto the fiber.
- Do not step on the fiber. Do not bend the fiber where it connects to the laser system.

2. Using aseptic technique, open the pouch at the symbol.
3. Remove the fiber from the pouch and the protective tube.
Place the fiber onto the sterile field.
4. Inspect the fiber for damage.

- Do not use if the fiber appears damaged. Replace the fiber with a new one. Use of a damaged fiber may result in injury to the patient or injury to the user.
- Do not sharply bend the fiber.

5. Insert the Optical Connector of the fiber into the fiber port on the laser system and turn clockwise, until it locks (1/4 turn).
6. Place the laser in READY mode to activate the aim beam.

- Do not press the footswitch while checking the aim beam. Pressing the footswitch will activate the laser beam. Inappropriate activation of the laser beam may result in injury to the user, personnel in the room and/or the patient.

7. Position the distal end of the fiber on a non-reflective sterile surface and turn slightly until the aim beam is visible.
8. Return the laser to the STANDBY mode.

- Do not use if the red colored aim beam is not visible. Replace the fiber with a new one. Use of a damaged fiber may result in injury to the patient and/or injury to the user.

Performing the Procedure

1. Ensure that ALL persons (user, patient and personnel) in the room are wearing appropriate protective eye wear prior to using the laser.
 - Do not press the footswitch without having the appropriate protective eye wear. Pressing the footswitch may activate the laser beam. Activation of the laser beam when not wearing protective eye wear may result in eye injury to the user, personnel in the room and/or the patient.
2. Connect the tubing from the sterile saline for irrigation solution to the fiber luer lock connector. Position the saline solution at least 106 cm higher than the level of the endoscope/cystoscope.
3. Open the flow valve. Fluid should be observed flowing out of the metal cap window.
 - Failure to position the saline solution at least 106 cm or not opening the flow valve will decrease the amount of fluid across the fiber end. Decreased fluid flow may limit the lasing time of the fiber.
4. Insert and advance the endoscope/cystoscope to the treatment area using standard medical practice.
5. Insert the fiber into the endoscope/cystoscope until the distal end of the fiber is visualized. A blue triangle and a red octagon are located on fiber to provide an extension reference.
6. Laser light is emitted opposite the blue triangle. **The blue triangle must be visible when firing the laser.**
 - Do not fire the laser when the red octagon is visible as this may result in damage to the cystoscope or unintended tissue.
 - The laser beam is emitted at an angle. Avoid forward deflection of the beam at the bladder neck as damage to the ureteral orifices may occur.
7. Place the laser in READY mode. This will enable the footswitch control.
8. Determine the direction of the laser beam by using the red colored aiming beam. The laser beam is in-line with the raised indicator on the circular knob.
 - Do not depress the footswitch to activate the working beam unless the aim beam is visible on the intended tissue.
 - Do not depress the footswitch to activate the working if the aim beam exits straight from the end of the fiber. Replace the fiber.
9. With the aim beam on the targeted tissue, depress the footswitch and begin lasing.

The treatment times will vary based on distance to tissue, the power settings, duration of beam application and other factors. Use settings that are as low as possible for the desired treatment effect.

 - Do not activate the laser if the blue triangle is not visible. Damage may occur to the fiber and/or endoscope/cystoscope and require replacement of the damaged devices.
 - Do not bury the tip of the fiber into the tissue. Reduced irrigation fluid flow may result in damage to the fiber.
 - Do not retract, dissect or probe tissue with the tip of the fiber. Damage may occur to the fiber tip.
 - Do not press the fiber end into the tissue. This creates stress between the fiber and the opening (edge) of the cystoscope. Damage to the fiber may result.
 - Do not bend the fiber at sharp angles. Damage may occur to the fiber.
10. The fiber can be used up to an energy or time limit. The user is alerted by an audible tone and a screen message when less than 50,000 Joules or 5 minutes remain. Once the energy or time limit is reached, to continue the procedure, a new fiber card and fiber must be inserted.

Troubleshooting and Procedural Considerations

Condition	Steps
A sudden increase in bubbles at the fiber tip	1) Increase the amount of irrigation through the cystoscope. 2) If increased irrigation does not eliminate the bubbling FOLLOW the steps described above. 3) If the problem persists, discontinue use and replace with a new fiber.
Lack of bubbles Inefficient vaporization due to the sweep speed and distance from tissue	1) Adjust sweep rate of the laser beam for even tissue removal. 2) Ensure that the distance from the fiber to the tissue is approximately 2 mm (1 to 3 mm). 3) Adjust the output of the laser.
Accumulation of debris in the metal cap window or on the fiber	NOTE: The liquid cooling design of the MoXy fiber eliminates the need to clean the fiber during recommended use. ▪ Clean only when tissue or other foreign material is forced into the window. Unnecessary cleaning may damage the fiber. 1) Turn the laser to the STANDBY mode 2) Remove the fiber from the endo/cystoscope 3) Gently clean the metal window opening with wet sterile gauze or pointed swab 4) Insert the fiber into the endo/cystoscope and resume the procedure. ▪ Failure to place the laser in the standby mode may result in user, personnel or patient exposure to laser energy and potential injury. ▪ Excessive or rough cleaning of the fiber tip may result in damage to the fiber. Accumulation of debris may result in over heating of the fiber cap leading to premature fiber degradation and/or fiber failure.

After the Procedure

1. With the endoscope/cystoscope still in place, pull the fiber into the scope approximately 5 cm. Confirm the amount of desired tissue is removed. Observe for any bleeding and restart irrigation and coagulate as necessary.
2. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

How Supplied and Storage

Store at room temperature (4°– 40° C, 39.2° to 104° F) and keep dry.

WARNING: Contents supplied STERILE.

Do not use if sterile barrier is damaged. If damage is found, call your AMS representative.

For single patient use only. Do not reuse, reprocess or resterilize.

Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Warranty

The Warranty for this product covers only defects in the materials and workmanship. The GreenLight MoXy Fiber Optic is a disposable device and is intended to be used in a single procedure and then discarded. Should you need to return this product, call the AMS Customer Care Center or the appropriate international distributor for "Return Authorization" (contact information is listed on the cover page). Instructions for returning product to AMS are detailed in the Fiber Replacement Policy form included with each device product package.

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