Web:www.esacrom.com E-Mail:esacrom@iol.it www.esacrom.it

Esacrom srl

Manuale operativo Micromotore

Versione 1/1 del 05/05/00

Challenger 40

Web:www.esacrom.com E-Mail:esacrom@iol.it www.esacrom.it

Indice

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Web:www.esacrom.com E-Mail:esacrom@iol.it www.esacrom.it

EQUIPAMENT

- Pedal-operated power unit complete with power
 1.5 mm setscrew L-wrench -
- Micromotor handpiece complete with lead -
- Soft rubber micromotor unit holder -
- Collet adapter for 1.6 mm dia. Burrs
- Collet wrench -
- Use and maintenance manual -
- Certificate of guarantee

CAUTION

IMPORTANT!!!

ESACROM has compiled this booklet to help you handle and operate the equipment better.

ESACROM cannot accept any responsibility for damage resulting from incorrect use or modifications made to the equipment.

A copy of this booklet accompanies each machine. To ensure correct machine operation, please make sure you follow carefully the instructions provided. The machine must never be used for purposes other than those indicated in the operating instructions and leaflet.Non-compliance with the instructions contained in this booklet will result in the immediate termination of the guarantee. On all requests, always indicate the date of purchase, the machine model and the serial number shown on the plate located on the back of the machine. Removal of this plate will result in the immediate termination of the guarantee.

INTRODUCTION

This manual refers to the micro-motor with pedal control model 'MI 400'.

ESACROM congratulates and thanks you for your choice.

The following instructions will help you get to know and use the machine better. Always make sure the certificate of guarantee is attached to the machine and that this has been completed and stamped by the stockists. A copy of this guarantee should remain with the stockists and with ESACROM. Before unpacking the machine, use the following check list to ensure no component parts are missing. After carefully removing all parts from their packing, place these on a table to ensure they have not been damaged during transport. If any parts are found to be damaged, contact your stockists immediately. Also make sure that all packing residue is removed from the machine. Now follow the installation instructions in this booklet.

PRECAUTIONS

IMPORTANT -The micromotor is a dental laboratory instrument.

IMPORTANT -No other possible uses exist for this instrument outside the dental sphere.

The instrument must be used by qualified personnel who are familiar with the uses for which the instrument is intended. This instrument provides high-speed circular movement to tools called burrs. If it is used incorrectly or for purposes other than those for which it was intended in dental laboratories, it can become dangerous for the user and for those in the immediate vicinity. Never attempt to stop micromotor rotation with a foreign body or, even worse, with your hands. This would be hazardous for the user and could damage the instrument. It should always be remembered too that this instrument has a power rating of 220V so that all necessary precautions must be taken during handling. Always remove the burr from the micromotor after use to prevent this representing a hazard in the event of accidental pressing of the foot-pedal.

DESCRIPTION

Challenger 40 is the result of close cooperation with modern dental laboratories and of a careful study of current technological innovations in the electro-mechanical and electronic fields. The machine comprises a handpiece operating in low-voltage and connected by means of a power cable to a foot-pedal which, in turn, is connected to the electric network. The handpiece is made of two different materials; blu galvanized alluminium for the rear part and black derling for the front part. Challenger is a dental laboratory machine which operates at 40,000 rpm, making it suitable for modern machining techniques and the use of high-speed cutting burrs. The machine is distinguished by the type of main electric motor component which operates according to the latest 'brushless' technique used in modern numeric control industrial processes. Unlike traditional motors, the above technology allows the motor to operate with direct current without those two typical components, the graphite brush and the commutator, which are subject to wear with consequent maintenance problems. Special features of machine mechanics are the high-speed low-noise bearings, an adjustable and easily replaceable tool chucking system supplied in the standard size of 2.35 mm and, on request, 3 mm. The chucking system or collet is adjusted to within a maximum concentricity discrepancy, with respect to the rotation axis, of

ESACROM S.R.L. P.IVA 02007591205 R.E.A Bo n° 64274 CAP SOC.30.000.000 I.V.

Web:www.esacrom.com E-Mail:esacrom@iol.it www.esacrom.it

0.02 mm. The foot-pedal is covered in shockproof plastic and rests on rubber supports for better stability. It features a large, easy-to-use control lever and a sophisticated electronic control system which allows choosing between two speed ranges by either pressing the pedal to the right or left. It also features torque compensation adjustment which recovers, at low speeds, 80% of maximum torque. The machine is equipped with dual current limiter which protects the electric motor in the event of overloads.

The power lead connecting the electric motor to the foot-pedal can be disconnected at both ends for easy replacement.

INSTALLATION

Though during packing and preparation for shipment, all precautions have been taken to prevent damage during transport, it is important to check the condition of the machine after removing it from the packing. In the event of any damage, immediately notify the supplier. It is best not to throw away the packaging as this may be useful if the machine has to be sent back for checking or overhauling. The foot-pedal should be located under the work bench within easy reach of the operator to ensure a comfortable working position. Care should be taken that there are no objects (leads, tubes, platforms, etc) near the pedal which might interfere with operation and prevent the rubber supports from resting correctly on the floor with consequent pedal instability. Place the micromotor on the rubber holder on the work bench, in an area where there is plenty of room for manoeuvre. Check the power lead between the pedal and the handpiece. This must be as straight as possible, without any bends. It is best not to use pulley lead winding systems. Connect the plug of the micromotor lead to the round socket at the back of the pedal, being careful to match up the position references on the plug and the socket. Fasten the plug in position with the locking ring to prevent bad contacts.

<u>IMPORTANT</u> - Connect the pedal power lead to a power socket. Always ensure the power ratings are the same as those shown on the plate at the back of the machine.

Make sure the power socket is correctly earthed. At this point, power can be provided to the foot-pedal by activating the switch on the rear of the pedal. This will light up. After fitting a burr on the micromotor, the machine is ready for use. To operate, press the pedal lever.

USE AND MAINTENANCE

The micromotor can be used for all kinds of dental laboratory machining operations where suspended handpieces, fixed motors and micromotors are called for.Before starting the micromotor, fit a burr in the collet by turning the ring nut clockwise. Always remember to tighten the ring nut (turn anti-clockwise).


Fig. 1 - Opening and closing the collet

Always ensure that the burrs used are in good condition, especially the high-speed rotation burrs. If these are badly balanced or have a crooked shank, the considerable centrigugal force could cause breakages with

Web:www.esacrom.com E-Mail:esacrom@iol.it www.esacrom.it

consequent danger for the user or other people in the vicinity. The rotation speed of the burrs must be calculated and adjusted in relation to their mass and construction features. For example, burrs in carborundum grain or similar, of large size and with inserted shanks cannot rotate at maximum speed. The speed limit for such burrs is 18,000 rpm. Never exceed the speeds recommended by the manufacturers. To facilitate speed control, the pedal has been built so that, by turning the pedal lever to the right, rotation speed can be adjusted between minimum and maximum. By turning the lever to the left, rotation speed can be adjusted from minimum to a maximum of 30,000 rpm.

The micromotor has a forced-air cooling system which sucks air in from the central part of the handpiece and blows it out at the back near the lead connection; never completely obstruct the suction area when holding the handpiece. To achieve best performance in terms of useable power, it is important to 'feel' the effort being required from the motor. During use, it is normal to pass from situations of minimum effort to situations of maximum effort depending on the pressure applied by the burr to the workpiece. If the effort required from the motor exceeds maximum levels, power will cut out automatically. It is always best to lighten pressure on the burr periodically while working so the micromotor can return to maximum torque.

MAINTENANCE PROGRAMME

And now a few words about maintenance. The handpiece is equipped with high-speed pre-lubricated bearings. **Do not lubricate these bearings**, as this would only prevent correct operation.

- To clean the micromotor, never use compressed air. A soft cloth or brush are quite sufficient, except in the case of the specific maintenance operations indicated.
- Periodically check the air ducts. Air should always exit from the rear of the motor. Remove any material which might obstruct these ducts.
- Remove the dust protection cap at the end of the handpiece and clean it with compressed air. Subsequently brush with a soft brush the area around the shaft, holding the handpiece upright with the motor on top.
- Check the screw of the collet ring nut is tight.
- Check the collet is working properly and that it is clean.

If the collet does not hold the burr correctly during operation, this must be adjusted as follows:

- 1. Take off the cap pleced on the front part of the handpiece.
- 2. Turn the ring nut clockwis so the collet opens.
- 3. Insert the triangular collet wrench in the collet.
- 4. Insert the 1.5 mm collet wrench in the hole of the locking microscrew and, holding the triangular screw steady, turn anti-clockwise.
- 5. At this point the collet is unfastened and, by turning it anticlockwise with the triangular wrench, can be removed for cleaning. After removing with the 1.5 mm wrench, further loosen the microscrew for greater field of adjustment.
- 6. Replace the collet in its housing after cleaning it with a cloth.
- 7. For better adjustment, use a small size burr with shank in good condition.
- 8. Insert the burr inside the collet and tighten the latter as much as possible. This is the maximum adjustment position.
- 9. From this position, turn the collet 180° anticlockwise to create a margin of play.
- 10. At this point, turn the ring nut anticlockwise so the collet closes.
- 11. Check the rotation concentricity of the burr by operating the micromotor. As an instrumental check cannot be carried out, assess the truing by resting the burr on a hard surface to see if there are any vibrations.
- 12. If vibrations do occur, alter the position of the collet with respect to the shaft by turning it a few degrees anticlockwise and repeat operation 10. To find the best position, it may ne necessary to repeat the operation a few times. Always be careful not to loosen the collet to much.
- 13. Once the best position has been found, lock the collet in position and then, with the triangular wrench, hold it steady, (being careful not to move the collet with respect to the shaft as this would jeopardize previous operations) and with the 1.5 mm wrench fasten the internal microscrew (reverese operation to number 3).

Web:www.esacrom.com E-Mail:esacrom@iol.it www.esacrom.it


Fig. 3 - Collet adjustment

TECHNICAL FEATURES

Voltage 220 V. +/- 10% (110 V on request)

Frequency 50/60 Hz **Absorbed power** 150 Watt

Fuse 1.25 Amp (delayed in glass 5 x 20 mm)

Minimum speed 3500 rpm **Maximum speed** 40000 rpm

SPARES

List of spare parts subject to wear and replaceable by user -

Plasticized cap complete with dust protection rings	Code	
2.35 mm dia. Collet	Code	
3.00 mm dia. Collet	Code	
Adapter for 1.6 mm collet	Code	
Ring nut screw	Code	
Ring nut screw ring	Code	
Complete lead for micromotor	Code	
Fuse	Code	
Foot-pedal power lead	Code	