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Torino - Turin, 22/04/2011

Amm. Delegato - *Managing Director*
Alberto CATELLA

C.I.E. DENTALFARM s.r.l.
 Via Susa, 9/a - 10138 TORINO - ITALY

 **SERVIZIO TECNICO COMMERCIALE - TECHNICAL COMMERCIAL SERVICE**
 (+39) 011/4346588
 **SERVIZIO ASSISTENZA TECNICA - AFTER-SALE TECHNICAL SERVICE**
 (+39) 011/4346632
 **FAX (+39) 011/ 4346366**
 E-mail: info@dentalfarm.it
 Website: www.dentalfarm.it

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SOLAR LED 3x+12x is a precision instrument which has been manufactured in compliance with the highest optical and mechanical features in order to ideally satisfy the daily requirements of the Dental Laboratory.

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Plafoniera - Roof lamp :

| lampade - lamps | tensione - voltage |
|--|----------------------------|
| 60 LED su circuito stampato (5W ca.) 60 LED sources on printed circuit (5W approx.) | 220 V – 240 V 50 Hz |

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EXPLODED DRAWING AND SPARE PART LIST SOLAR 3x+12x LED

| POS. | CODE | DESCRIPTION |
|------|---------|------------------------------------|
| 1 | 7003013 | LENS PROTECTION CAP |
| 2 | 7003019 | MAGNIFYING LENS (3 + 12 DIOPTERS) |
| 3 | 7002023 | SWITCH |
| 4 | --- | HEX KEY FOR FRICTION ADJUSTING |
| 5 | 7002022 | ELECTRONIC BALLAST |
| 6 | 7003024 | MIDDLE ARM JOINT WITH FRICTION |
| 7 | 7003023 | ROOF LAMP BALL JOINT |
| 8 | 7003021 | COMPLETE LED RING |
| 9 | 7003015 | TRANSPARENT LAMP PROTECTION SHIELD |

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Gently remove the lamp from its packaging. The appliance is equipped with a balanced base to be positioned on the bench, near the working area. The arm features two joints allowing to move the roof lamp and to adjust its height. The joint in the middle of the arm is equipped with a friction adjuster (use the hex key incorporated in the lower part of the lamp base). As ideal place for the installation, choose a well ventilated, clean and dry location, with temperature ranging between 0°C and 40°C and relative humidity not exceeding 90% (indoor use only).

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Before connecting the electrical cable, be sure that the rated voltage stated on the label corresponds to the mains voltage and that the socket is duly supplied with earthed wire. The section of the feeding cable supplied with the unit suits to the highest absorbed current (specified in the technical features) and it is one meter long. Should you use any extension, be sure that the section of the cable matches to its length and that it is at least of the same size as the original cable.

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The roof lamp is fitted with a protective cover ensuring that the lens is kept clean; the cover must be closed in idle mode. Switch on and tune the lamp on the object to be enlarged or lightened, then adjust the focusing by gently lifting or lowering the lamp. Never force the mechanisms and do not try to force the movements of the middle joint; if needed, adjust its friction as explained above. Do not forget to switch off the lamp at the end of the work.

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Use a soft cloth to clean the lens, eventually use a liquid, antistatic and non-abrasive cleansing agent. The other parts of the lamp can be cleaned with any non-aggressive detergent or disinfectant agent specific for the plastic and the enameled metallic surfaces.

Operations to be avoided: to clean the lens surface with your hands, to misuse and apply unnecessary loads to the arm, to clean the appliance with volatile solvents or abrasive cleansing agents, to try to recondition the lamp by yourselves.

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Replacement of the lens – unscrew the 5 screws fixing the transparent plastic screen and lift the LED ring; the lens is now free and can be removed from its seat.

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Replacement of the complete LED ring – unscrew the 5 screws fixing the transparent plastic screen, lift the LED ring and unweld the two electric contacts from the printed circuit.

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Replacement of the electronic ballast:- unscrew the 4 screws hidden under the supporting feet, unscrew the two screws fixing the plastic case and open it; after replacement of the part, restore all electric connections one by one.

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

Never use the unit for any use differing from the operating purposes it has been designed for and which are clearly illustrated in the present instruction manual.

In case of demolition, comply scrupulously with the regulations in force in the Country where the operation takes place, always considering that you are handling an electrical device.

C.I.E. Dentalfarm Srl reserves the right to modify at any time and without prior notice any features of its products.

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