DAKO 5410

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SARATOGA S.P.A.

Manufacturing of dental surgeries and laboratories equipments

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INDEX

Aim of manual.....	8
Identification of the machine manufacturer	8
Technical information.....	8
Safety devices	8
Warranty.....	9
Technical data.....	10
Packaging and unpackaging	10
Loading and unloading.....	10
Installation instruction	10
Control devices descriptions	11
Function and use instructions	12
General instructions	12
Ordinary maintenance.....	12
Extraordinary maintenance.....	12
Precaution.....	12
Spare parts (explosion diagrams).....	13
Declaration of conformity	15

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AIM OF THE MANUAL

This manual has been written by the manufacturer and it is an integral part of the machine.


Its information are addressed to the user and contain all Safety Recommendations.

Before using this mixer, it is recommended that you read this manual carefully, especially if it is the first time, to ensure you know the control devices, their function and position. It is also advisable to make performance tests.

This manual must be kept for future references.

IDENTIFICATION OF THE MACHINE MANUFACTURER

On the machine you can find the label as shown in pic. 1.

pic. 1

TECHNICAL INFORMATION

The DAKO 5410 vacuum mixer has been planned and built to mix coatings, casts and plasters in dental laboratories in accordance with all Safety Recommendations (pic.2).

SAFETY DEVICES

On the machine are some Fixed Guard: to avoid access to the electronic and mechanical parts.

WARRANTY

Approved by UNIDI

1. Through the present document the manufacturer certifies a correct construction of the product and assures that first materials have been used; all necessary approvals and trials have been made so that the product is conform to the regulations in force.

The guarantee on this product is valid for a period of 12 months starting on the day of delivery to the purchaser. The purchaser approves the guarantee regulations by filling in and signing his part and returning it to the manufacturer. The guarantee only covers those parts to be substituted or replaced because of construction faults and doesn't cover labour expenses, travelling indemnity of technicians, transport expenses, packing etc.

The guarantee excludes damages or failures caused through a bad maintenance, a wrong feeding supply, negligences, unskillfulness or other causes not imputable to the manufacturer.

The guarantee excludes as well failures caused through omitted maintenance due to the negligences of the user.

The guarantee doesn't include any kind of refund versus persons or objects for damages, direct or indirect, due to eventual inefficiency of the equipment.

2. The guarantee forfeits automatically in case that the equipment has been repaired, modified or tampered by the purchaser himself or by unauthorized third parties.
3. For interventions under guarantee, the purchaser must address himself to the selling-agent, or the by him indicated technical-service or direct to the manufacturer.
The guarantee gives the right to substitution of the defective part, but excludes the substitution of the entire equipment.
4. In case of contestation concerning the application of the guarantee regulations about the quality or conditions of the delivered equipments, the purchaser will in no way suspend or delay the payment of the price or the rates of price.
5. No refund can be requested by the purchases for the temporary stop of the equipment.
6. The guarantee forfeits if:
 - a) the equipment presents damages of dropping; exposure to flames; overturning of liquid; natural disasters; or any other causes not imputable to fabrication defaults;
 - b) the installation wasn't correct;
 - c) the connection to the main supply was wrong rated feeding voltage;
 - d) the serial no has been taken of, cancelled or altered.
7. The components, to be substituted under guarantee must be sent back to the manufacturer for reimbursement following carefully his shipment instructions. For no reason we will accept the returning of equipments or spare parts of them without having given written authorization before. In case of non observation of these conditions the manufacturer reserves himself the right to send the goods back to the sender.
8. **Due to tax reasons the spare parts will be granted under guarantee only if we have received back the guarantee coupon properly filled in with all customer's data.**

Warranty regulations

1. The non restitution of the guarantee certificate involves the immediate withdrawal of the certificate itself.
The guarantee doesn't cover labour expenses and travelling indemnity which will always be on charge of the purchaser.
2. The payment of the labour expenses, travelling indemnity, and right of appeal bill will be have to fulfilled immediate at the presentation of the bill itself. Not fulfilling the payment of the bill involves immediate withdrawl of the guarantee itself.
3. Neither the manufacturer or the dental deposit are obliged to provide substitutive equipment during the repair period.
4. For any other cases not foreseen in the present guarantee certificate and regulations reference is made to the italian civil code.

TECHNICAL DATA

Overall dimensions (lxdxh)	23x28x44cm
Weight	16kg
Voltage	230V
Power	300W
Paddling speed	100-500rpm/min
Standard bowl capacity	650cc
Little (R1) bowl capacity	250cc
Maximum vacuum	720mmHg
Fuse rating	5A

PACKAGING AND UNPACKAGING

The packaging is as follows:

- polystyrene housing K10;
- polystyrene plate;
- hard cardboard box.

LOADING AND UNLOADING

While moving the machine avoid absolutely any kind of bumping, dropping or tilting: they could seriously damage the machine. Obviously, the manufacturer is not responsible for damages caused by droppings, improper use and maintenance which are not in strict accordance with the manufacturer's instructions contained in this manual.

Packaging waste disposal must be done with respect to the environment and the law in force.


INSTALLATION INSTRUCTIONS

Connect the power supply plug into a 220V AC outlet.

It is duty of the user to ensure that the electric plant is in accordance with the safety laws in force. It is particular important to make sure that the grounding is well efficient.

Furthermore, it is important to verify the network voltage: in case the voltage is too low (lower than 210V), the engine will start too slow and it might be necessary to install a voltage stabilizer.

CONTROL DEVICES DESCRIPTION (PIC. 2)

pic. 2

1. Fusesholder.
2. 220 V AC power supply plug.
3. General power switch: the machine turns on and the display 7 lights up.
4. Fast clutch: it connects the bowls to the vacuum pump.
5. Pump (ON/OFF button): turns on and off the vacuum pump.
6. Vacuum exhaust (VAC/OFF button): by pressing and keeping it pressed the air will go in the bowls and therefore eliminates the vacuum.
7. Pump warning light: it lights up when the pump is on.
8. Timer display: it indicates the set time for the paddling cycle. During the paddling the time will decrease up to 0. Time is indicated in seconds.
9. On by rotating that knob we can vary the motor speed : 100 - 500 rpm.
10. Paddling time: by rotating this knob you will set the paddling time and the display 8 will show it. By rotating the knob to the right the time will increase up to a maximum of 99 seconds whereas by rotating it to the left the time will decrease up to 1 seconds.
11. START pushbutton: by pushing it, you start the paddling cycle.
12. Jointed coupling: it transmits the rotation from the engine to the spattle spindle.
13. Bowls lid.
14. Spattle.
15. Bowl.
16. Vacuum gauge.

FUNCTION AND USE INSTRUCTIONS

- Turn the machine on with the general power switch 3.
- Remove the bowls 15 from the machine and place in it the needed amount of coating and water.
- Close the upper bowl with the lid 13.
- Lift the coupling 12 and leave it at its high position.
- Place the bowl on the machine centering it in the special reference.
- Let the jointed coupling 12 go down and make it fit on the bolt at the top of the bowl.
- Hook up the clutch 4 and close the nearby slow exhaust manual valve by rotating it to the right.
- With the knob 10 set the paddling time: the display 7 will show the set time.
- With the knob 9 set the motor speed (100-500 rpm)
- Start the paddling cycle with the pushbutton 11 (START).
- Start the vacuum pump with the blue button 5 (PUMP): the vacuum gauge pointer 16 will reach in few seconds values close to 70 mmHg. If this doesn't happen it means that the bowl is not properly closed .
- As the set paddling cycle time elapses the engine will stop.
- Turn off the pump (PUMP button 5) eliminate the vacuum in the bowl by pressing and keeping pressed button 6 (VAC/OFF) or by opening the slow exhaust valve located on the bowls lid.
- Remove the bowl from the machine by lifting the coupling 12 and disconnecting the clutch 4.
- Clean carefully the bowl, the paddle and the bowl lid: if necessary lubricate the OR with silicone grease.

GENERAL INSTRUCTIONS

It is advisable to prepare the mixture of the coating, with the right dosage, in a separate container and perform a manual quick pre-paddling: this is to avoid introducing in the bowl powder and water fully separated with the risk that the powder could be sucked out by the vacuum pump and the water could deposit at the bottom.

Avoid placing too much coating in the bowl so not to block the vacuum pipes: stay 2 centimetres minimum below the bowl upper edge.

ORDINARY MAINTENANCE

All the packings of the machine must be kept well cleaned and lubricated with silicone grease and checked they are always in perfect conditions.

To clean the packings never use any solvents that could damage them.

EXTRAORDINARY MAINTENANCE

To repair the machine or replace some of its parts, the user should consult directly the manufacturer or only specialized technicians.

Absolutely avoid to open any of the fixed guards without proper precautions.

Disconnect the power supply plug from the outlet before any extraordinary maintenance operation.

PRECAUTIONS

Inside the machine there are components subjects to 220V network voltage.

Before performing any operation inside the machine make sure you turned off the general switch 3 and unplugged the power supply 2 (same precautions must be taken when substituting the fuses).

If the user doesn't strictly follow the rules contained in this manual the manufacturer is lifted from any guarantee and responsibility.

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Application of the Directives E.E.C. 89/392 - E.E.C. 89/336 - E.E.C. 73/23

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Conforms with the following Directives E.E.C. 89/392 - E.E.C. 89/336 - E.E.C. 73/23.

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