Zanotti Medical Research

APPARECCHIATURA PER APPLICAZIONE DEI TRATTAMENTI TERMICI
UNIT FOR COLD AND HOT THERAPY
EQUIPEMENT POUR L'APPLICATION DE TRAITEMENTS THERMIQUES
WÄRMEAUSTAUSCHGERÄT FÜR KÄLTE UND WÄRMEBEHANDLUNGEN
EQUIPO PARA LA APLICACIÓN DE TRATAMIENTOS TÉRMICOS

MS

MANUALE DI USO E MANUTENZIONE
USE AND MAINTENANCE MANUAL
NOTICE DU MODE D'EMPLOI ET D'ENTRETIEN
GEBRAUCHS-UND WARTUNGSHANDBUCH
MANUAL DE USO Y MANTENIMIENTO

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CONTENTS

1. Safety recommendations	Pag 18
2. Plates	Pag 29
3. Description of the unit	Pag 20
4. Electric controls - Operation and starting	Pag 21
4.1 Control panel of unit AA (Cold-Heat)	Pag 21
4.2 Operation and starting of unit AA (Cold-Heat)	Pag 22
4.3 Control panel of unit BA (Only-Cold)	Pag 24
4.4 Operation and starting of unit BA (Only-Cold)	Pag 25
5. Start up	Pag 27
5.1 Application (Installation)	Pag 27
5.2 Warnings	Pag 27
5.3 Dimensions	Pag 27
5.4 Safety measures and recommendations	Pag 27
5.5 Cleaning	Pag 28
6. Connecting the unit to external power sources	Pag 28
7. Connection to the mains	Pag 28
8. At start up	Pag 28
8.1 Checks	Pag 28
9. Maintenance and repairs	Pag 29
9.1 Glycol level	Pag 29
10. Routine maintenance	Pag 29
10.1 Special maintenance	Pag 29
10.2 Operations to be carried out by qualified technicians	Pag 29
10.3 Troubleshooting	Pag 29
11. How to order spare parts	Pag 30
12. How to dispose of the packing	Pag 30
13. How to dispose of the unit	Pag 30

Thank you for choosing this Zamar product.

Please read these instructions carefully. They provide details and advice on the correct method of installing, using and maintaining this unit, in order to use it in the best way and obtain its maximum performance.

1) SAFETY RECOMMENDATIONS

When installing and using the unit please follow the recommendations listed here below.

- Installation shall be carried out in strict compliance with the diagrams and instructions supplied by the manufacturer.
- Damages due to improper connections are excluded from any warranty.
- The mains electricity supply where the unit is installed must meet with the relevant standards in force.
- Maintenance should be effected by trained personnel and by the manufacturer, in compliance with all regulations and provisions in force.

Attention!

The unit is not suitable for working in explosive environments.

Therefore the use of the unit in an explosive-dangerous atmosphere is absolutely forbidden.

Attention!

The unit is not suitable for working in corrosive environments. In this event it should be protected with appropriate means.

When maintenance involves operations on the refrigerating circuit, empty the system and let it reach the atmospheric pressure.

Warning!

Do not discharge the refrigerant in the atmosphere. Specialised technicians using suitable equipment must recover it.

- The quantity and type of the refrigerant to be charged is indicated on the data plate.
- Do not use refrigerants of a different kind, especially inflammable fluids (for example hydrocarbons) or air.
- Do not modify or alter the refrigerating circuit or its components (for example: welding on compressor body).
- The final user shall protect the system from external fire hazards.

Warning: Use of the unit as well as heat/cold applications, regulation of treatment

temperature and length shall be carried out by qualified personnel or on the basis of instructions supplied by the same.

2) PLATES

GB

ZAMAR		CE
ZAMAR srl Str. Naz. Cisa Km 162,34/A 46029 Suzzara (MN) ITALY		
Anno di fabbricazione Year of manufacture	2010	
Modello Model	MS385BA00F	
Matricola Serial number	0809102A	
Tensione Voltage	V/ph/Hz 230-1-50	
Assorbimento Absorpt	A 2.5	W 440
Potenza Compressore Nominal horse power	Kw 0.39	
Refrigerante Refrigerant	R404A	Kg 0.29
Massa Mass	Kg 23.3	
Isolamento classe 1- parte applicata tipo B F 		
Class insulation 1- Application part type B F		

Technical plate


- 1) Year of manufacture
- 2) Zamar unit code
- 3) Serial number
- 4) Voltage
- 5) Absorption
- 6) Rated power of compressor
- 7) Refrigerant: Type, Quantity
- 8) Mass of the unit

Refrigerant plate

Warning plates
Attention: cold or hot parts

Attention: moving parts

3) DESCRIPTION OF THE UNIT

This unit is used for the local application of heat and cold therapy and is suitable for treating injuries on different parts of the body.

The unit includes three main parts:

A) Main body consisting of compressor, expansion element, heat pump, heat exchanger and control panel. These elements are contained in a portable case for an easy transport (picture 1). Refrigerant used: R404A.

B) Heat exchanger of different sizes, made of polyurethane (PU), consisting of a coil-shaped PU hose placed in a form-fitting treatment pad. Inner part in hypoallergenic material. Equipped with elastic strips in Velcro™ for correct positioning (picture 2).

C) PVC insulated hoses to supply specially treated liquid from the unit to the exchanger (picture 3).

The unit can produce “cold” only (unit MS---B-00F) or “cold” and “heat” (unit MS---A-00F) and convey them wherever required through the hoses connecting the case and the heat-exchanging bags. Tight quick couplings are used to the purpose.

The control panel of the unit is used to select the desired temperature treatment required.

Picture 1 =Case

MS365/385_ _00F

MS599_ _00F

Picture 2 =Heat-exchanging bag


Picture 3 = Connecting hoses

4) ELECTRIC CONTROLS - OPERATION AND STARTING

4.1) Control panel of unit

MS_ _ _ A _ 00F (Cold-Hot)

- **A** “display”: shows liquid temperature. During programming it shows work set points.
- **B** “display”: shows the pause/work time left (minutes). During programming it shows the length of each work/pause time (minutes).
- **C** “display”: it shows the number of cycles, during programming it shows the number of set cycles.
- **D** “display”: dedicated to the massage function, it shows pump times (adjustable from 2 to 9 seconds).
- **E** “button”: to increase temperature set points and times at programming stage.
- **F** “button”: to decrease temperature set points and times at programming stage.
- **G** “button”: to increase/decrease pump work times in massage function.
- **H** “button”: to increase/decrease pump pause times in massage function.
- **I** “button”: chose the degrees in Celsius or Fahrenheit
- **L** “button”: to enable/disable the massage function. When led 1 is ON, massages is enabled.
- **M** “button”: it allows to set point value of Cold-Hot.
- **N** “button”: to select the operating mode desired. On starting it is in Cold position, Press it to reach the following operating modes: “Only Cold” again.
- **O** “button”: it allows to change the working time in Cold-Hot.
- **P** “button”: it allows to enter programming, number of cycles in repetition of the set program.
- **Q** “button”: it allows to change the pause time in Cold-Hot.
- **R** “button”: ON-OFF programs.

- led 1: massage function on.
- led 2: program selection mode.
- led 3: number of cycles mode.
- led 4: unit connected to the mains supply; flashing light (stand-by), fix light (working).
- led 5: set point for Cold-Hot mode being reached.
- led 6: Cold-Hot mode, work time in course
- led 7: Cold-Hot mode, pause time in course
- led 8: Cold mode, on.
- led 9: Hot mode, on.
- led 10: continuous cold mode.
- led 11: pump alarm.
- led 12: no-liquid alarm
- led 13: high pressure alarm.
- led 14: temperature alarm.
- led 15: temperature in Celsius.
- led 16: temperature in Fahrenheit.
- led 17: setup **ON** of the pump.
- led 18: pump operating.
- led 19: setup **OFF** of the pump.

- Switch **0/1** on the side of the unit: to switch the unit on/off. To start operation press button **"R"** on the front.

4.2) OPERATION AND STARTING

The unit is controlled by a thermostat (A) by which the required temperature can be set, according to the different treatments.

Two temperatures (COLD/HOT) can be set; work times can be alternated with pause periods. A cycle can be repeated as many times as desired.

A) Setting temperatures

Make sure that the connecting cable has been correctly placed in the socket. Press starting switch (C) on the left side of the unit, near the socket.

After a short pause the unit is switched on. If necessary set the temperature as required

Above set points can be changed by pressing buttons as follows:

- 1 – press button **M**: led **5** flashes and the factory-set temperature for "cold" mode appears on display **A**;
- 2 - change the set point as required using buttons **E** (up) and **F** (down): the new set point is stored in memory automatically;
- 3 – press button **M**: led **5** flashes and the factory-set temperature for "hot" mode appears on display **A**;
- 4 – repeat the operations described in point 2.

B) Start up

After setting the desired temperature (see **A**) the unit is ready to start the working cycle.

- Press button **N** to select the working mode:

Led **8** on = Cold

Led **9** on = Hot

Leds **8** and **9** on = Cold-Hot

Led **10** on = Cold Continuous

- Press button **L** (led 1 lights up) to start the "massage" function.

- When the massage function is active (led 1 on), work (**ON**) and pause (**OFF**) times of the pump can be set using **H** and **G** buttons: in this way you can obtain the massage pressure desired. Display **D** show the set times (in seconds).

After above operations you can set/change the work and pause times as well as the number of cycles.

- Press button **O** "work time – cold mode". When led **6** flashes, the set time (in minutes) appears on display **B**.

- The set times can be changed by pressing buttons **E** (up) and **F** (down).

- Repeat above operations for "pause time "cold-hot" (button **Q**), "work time – hot mode".

- After setting the work / pause times, press button **P** "cycles": led **3** flashes and the number of pre-set cycles appears on display "**C**". Use buttons **E** and **F** to set the number of times by which the cycle will be repeated. A complete cycle corresponds to the total work and pause times.

- Set "--" to have the cycle repeated continuously.

To start the unit press button **R**: now the unit is operating and it will carry out the set cycle by the set times automatically, as follows:

Start • work in "cold" mode • pause in "cold" mode • work in "hot" mode • pause in "hot" mode

C) Operating modes

- When pressing "**R**" the unit always starts in the "cold" mode first.

- The green led **5** near button **M** (cold mode set point) on the control panel lights up and it remains on until the set point for "cold" mode is reached. At this stage the unit starts counting the work time in "cold" mode and the green led 6 lights up. The same occurs for the "hot" mode.

- Display **B** shows the number of cycles to be repeated; during the work and pause times the display shows the work/pause minutes left.

After programming the unit is operating.

- Display "**A**" shows the liquid temperature.

- Display "**C**" shows the number of missing cycles.


- Green led "**18**" on: the pump is working in continuous mode. When the Massage function is active, the green led "**18**" flashes; during the pause times it is off.

- Green led corresponding to buttons: when on, the corresponding mode is active.

4.3) Control panel of unit

MS_ _ _ B_ _ 00F (Cold)

GB

- **A** “display”: shows liquid temperature. During programming it shows work set points.
- **B** “display”: shows the pause/work time left (minutes). During programming it shows the length of each work/pause time (minutes).
- **C** “display”: it shows the number of cycles, during programming it shows the number of set cycles.
- **D** “display”: dedicated to the massage function, it shows pump times (adjustable from 2 to 9 seconds).
- **E** “button”: to increase temperature set points and times at programming stage.
- **F** “button”: to decrease temperature set points and times at programming stage.
- **G** “button”: to increase/decrease pump work times in massage function.
- **H** “button”: to increase/decrease pump pause times in massage function.
- **I** “button”: chose the degrees in Celsius or Fahtenheit
- **L** “button”: to enable/disable the massage function. When led 1 is ON, massages is enabled.
- **M** “button”: it allows to set point value of Cold.
- **N** “button”: to select the desired operating mode: when pressing ON the unit starts in Cold mode ; pressing the button again it changes to Continuous Cold.
- **O** “button”: it allows to change the working time in Cold.
- **P** “button”: it allows to enter programming, number of cycles in repetition of the set program.
- **Q** “button”: it allows to change the pause time in Cold.
- **R** “button”: ON-OFF programs.

- led 1: massage function on.
- led 2: program selection mode.
- led 3: number of cycles mode.
- led 4: unit connected to the mains supply; flashing light (stand-by), fix light (working).
- led 5: set point for "cold" mode being reached.
- led 6: "cold" mode, work time in course
- led 7: "cold" mode, pause time in course
- led 8: "cold" mode, on.
- led 9: N.A.
- led 10: continuous cold mode.
- led 11: pump alarm.
- led 12: no-liquid alarm
- led 13: high pressure alarm.
- led 14: temperature alarm.
- led 15: temperature in Celsius.
- led 16: temperature in Fahrenheit.
- led 17: setup **ON** of the pump.
- led 18: pump operating.
- led 19: setup **OFF** of the pump.

- Switch **0/1** on the side of the unit: to switch the unit on/off. To start operation press button **"R"** on the front.

4.4) OPERATION AND STARTING

The unit is controlled by a thermostat (A) by which the required temperature can be set, according to the different treatments.

Cold temperatures can be set; work times can be alternated with pause periods.

A cycle can be repeated as many times as desired.

A) Setting temperatures

Make sure that the connecting cable has been correctly placed in the socket. Press starting switch (C) on the left side of the unit, near the socket.

After a short pause the unit is switched on. If necessary set the temperature as required

Above set points can be changed by pressing buttons as follows:

- 1 – press button **M**: led **5** flashes and the factory-set temperature for "cold" mode appears on display **A**;
- 2 - change the set point as required using buttons **E** (up) and **F** (down): the new set point is stored in memory automatically;

B) Start up

After setting the desired temperature (see **A**) the unit is ready to start the working cycle.

- Press button **N** to select the working mode:

Led **8** on = Cold

Led **10** on = Cold Continuous

- Press button **L** (led 1 lights up) to start the “massage” function.

- When the massage function is active (led 1 on), work (**ON**) and pause (**OFF**) times of the pump can be set using **H** and **G** buttons: in this way you can obtain the massage pressure desired. Display **D** show the set times (in seconds).

After above operations you can set/change the work and pause times as well as the number of cycles.

- Press button **O** “work time – cold mode”. When led **6** flashes, the set time (in minutes) appears on display **B**.

- The set times can be changed by pressing buttons **E** (up) and **F** (down).

- Repeat above operations for “pause time” (button **Q**), “pause time”.

- After setting the work / pause times, press button **P** “cycles”: led **3** flashes and the number of pre-set cycles appears on display “**C**”. Use buttons **E** and **F** to set the number of times by which the cycle will be repeated. A complete cycle corresponds to the total work and pause times.

- Set “--” to have the cycle repeated continuously.

To start the unit press button **R**: now the unit is operating and it will carry out the set cycle by the set times automatically, as follows:

Start • work in “cold” mode • pause in “cold” mode • work in “hot” mode

C) Operating modes

- When pressing “**R**” the unit always starts in the “cold” mode first.

- The green led **5** near button **M** (cold mode set point) on the control panel lights up and it remains on until the set point for “cold” mode is reached. At this stage the unit starts counting the work time in “cold” mode and the green led **6** lights up.

- Display **B** shows the number of cycles to be repeated; during the work and pause times the display shows the work/pause minutes left.

After programming the unit is operating.

- Display “**A**” shows the liquid temperature.

- Display “**C**” shows the number of missing cycles.

- Green led “**18**” on: the pump is working in continuous mode. When the Massage function is active, the green led “**18**” flashes; during the pause times it is off.

- Green led corresponding to buttons: when on, the corresponding mode is active.

5) START UP

After above operations the unit is ready to start working.

The following operating modes are available:

MS---A-00F "Cold-Heat": Cold, Continuous Cold, Heat, Cold-Heat.

MS---B-00F "Cold": Cold, Continuous Cold.



5.1) Application (Installation)

- 1) The unit should be positioned correctly on a flat surface to prevent the liquid from leaking out. We recommend that the unit should be elevated and that a purpose built shelf or trolley should be used. The unit should not be used and operated on the floor. (containing the coils).
- 2) Place the bags where the treatment is required.
- 3) Switch on and start the unit.
- 4) Connect the hoses with the unit and the bags.
Follow this sequence of operations so that a correct pressure is maintained in the whole circuit.
- 5) When the treatment has been completed carry out operations contrary to the ones described above.
- 6) When the unit is off and the pipes are connected do not hold the treatment bags above the level/high of the unit or it may overflow.
- 7) If you want to switch off the unit without disconnecting the hoses and the bags, make sure that the latter are at a lower level than the unit.


5.2) Warnings

There are no warnings or warning labels on the unit.

5.3) Dimensions

	MS3	MS5
A	585	600
B	305	1088
C	266	353

5.4) Safety measures and recommendations

The unit has been assembled taking all measures necessary to ensure user's safety.

- Glycol-liquid leakage: should any glycol-liquid leak out of the unit during the treatment (because of a breakage of the coil applied on the part to be treated) act as follows:

- a) disconnect the coil immediately and stop the unit;
- b) take the actions described below, keeping in mind that the liquid can be harmful if swallowed.

First aid measures – general information

- **Inhalation** – in case of difficulty in breathing after inhaling vapours/aerosol, have the person transported in the open air. Seek medical attention.
- **Contact with skin** – rinse with plenty of soap and water.
- **Contact with eyes** – rinse with plenty of water for 15 minutes, keeping eyes wide open.
- **Ingestion** - rinse mouth with plenty of water. Administer 50 ml of pure ethanol in a drinkable concentration.
- **Note for the doctor** - treat symptomatically.

Accidental release measures (spillage of stocked liquid)

- Do not discharge into sewers and/or streams.
- Cleaning method: in case of spillage or leakage dam with earth, sand, sawdust, cloths, etc. and dispose of according to the provisions in force. Should large quantities of liquid be spilled, put the product into containers and dispose of according to the provisions in force.
- Additional information: high slip hazard, due to product spillage.

Disposal considerations

- Disposal procedures: discharge or incinerate according to local provisions. Non-washable containers should be disposed of following the same procedures as for the material contained.


Transport information

- Not classified as dangerous in any transport regulation (ADR, RID, ADNR, IMDG/GGVSee, ICAO/IATA)

Regulatory information

Directive 1999/45/EC:

- Symbol of danger
Xn - Harmful
- Risk phrases:
R 22: Harmful if swallowed
- Safety phrases:
S2: Keep out of reach of children
S24/25: Avoid contact with skin and eyes
S.46: If swallowed, seek medical advice immediately and show container or label.

5.5) Cleaning

The unit does not require any particular cleaning operations.

Clean the unit periodically following the recommendations supplied here below:

- frequency: at least once a month

- method:

1. disconnect the electrical socket from the power supply;
2. unscrew the securing screws on the case;
3. with a bristle brush dust battery and pipes;
4. be careful not to touch any electrical part;
5. check that the glycol-liquid tank does not show any leak (visual control).

6) CONNECTING THE UNIT TO EXTERNAL POWER SOURCES

Attention!

Before connecting the unit, make sure that voltage and frequency of mains correspond to those on the plate of the unit and that voltage is $\pm 10\%$ of rated value.

7) CONNECTION TO THE MAINS

After checking that the unit is in good condition, connect it to the mains: insert the supply cable in the socket at the back of the unit and then introduce the plug into the mains socket.

8) AT START UP

8.1) Checks

Before starting the unit, check that it is connected correctly.

Should the unit have been opened, make sure that:

- no tools have been left inside;
- assembly has been carried out correctly;
- there is no gas leakage;
- the control panel has been properly fixed.

9) MAINTENANCE AND REPAIRS

Routine maintenance is crucial for obtaining longer life, perfect working conditions and high efficiency of the unit as well as ensuring the safety features provided by the manufacturer.

9.1) Glycole

In order to guarantee a good performance of the thermic wraps, it is necessary to check periodically the glycol level within the relevant tank; in case of lack of this liquid it is simply necessary to do a filling up until the max level; we remind you that when glycol lacks, an automatic no-liquid alarm announces this droop.

Warning: it is extremely important to use as refrigerant the product suggested by ZAMAR. It is absolutely avoided to use any other kind of liquid.

10) ROUTINE MAINTENANCE

Good operation of the unit requires the condenser to be cleaned periodically (frequency of cleaning depends on the environment where the unit is installed).

Switch off the unit and clean it by blowing air from the inside outwards. Should no air jet be available, use a long haired brush and work on the outside of the condenser.

Warning - Disconnect the unit before working on it.

10.1) Special maintenance

Special maintenance is only required in case of troubles and only qualified technicians should carry it out.

10.2) Operations to be carried out by qualified technicians or by the manufacturer

The following operations require special technical skills and knowledge and should be carried out by qualified technicians or by the manufacturer exclusively.

Under no circumstances is the user allowed to:

- replace electrical components,
- work on the electric equipment,
- repair mechanical parts,
- work on the refrigerating system,
- work on the control panel, ON/OFF and emergency switches,
- work on safety and protection devices.

10.3) Troubleshooting

During operation following troubles may occur:

- the unit stops because there is no voltage:
check voltage in the mains and fuse conditions. if they are correct, have the unit inspected by a qualified technician.
- glycol-liquid leaking out.
identify the leaking point. if it is inside, have the unit inspected by a qualified technician.

Warning

Removal of protections during machine operation is absolutely forbidden. They have been developed by the manufacturer to safeguard the user's safety.

The following electrical safety devices can also interrupt operation of the unit:

- high pressure alarm
- no-liquid alarm
- pump alarm
- temperature alarm

In these cases:

HAVE THE UNIT INSPECTED BY A QUALIFIED TECHNICIAN.

11) HOW TO ORDER SPARE PARTS

When ordering spare parts refer to the serial number written on the unit data plate.

Warning

Worn parts should be replaced only by qualified personnel or by the manufacturer.

12) HOW TO DISPOSE OF THE PACKING

Wooden, plastic, polystyrene packing shall be disposed of according to the regulations in force in the country where the unit is used.

13) HOW TO DISPOSE OF THE UNIT

Do not discharge scrapped components in the environment. They should be disposed of by companies dealing with special waste collection and recovery, according to the regulations in force in the country where the unit is used.

Warning

Do not discharge the refrigerant in the atmosphere. Companies dealing with special waste collection and recovery should dispose of it.

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ZAMAR S.r.l.