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English

CONTENTS

1	-	General information	20
1.1	-	Using and storing the manual	20
1.1.1	-	Symbols used in this manual	20
1.2	-	Description of staff	21
2	-	Description of appliance	22
2.1	-	Description of appliance and foreseen use	22
2.2	-	Environmental conditions	22
3	-	Safety	23
3.1	-	General warnings	23
3.1.1	-	Precautions for users	23
3.2	-	Contraindications	24
4	-	Use	25
4.1	-	Instructions and precautions	25
4.1.1	-	Application precautions	25
4.1.2	-	Specifications	26
4.1.3	-	Reference to standards	26
4.2	-	Tanning procedures	27
4.2.1	-	Face and Décolleté treatment (figs. 3 and 4)	27
5	-	Maintenance	28
5.1	-	Gun maintenance	28
5.1.1	-	Cleaning the spray gun	28
5.2	-	Compressor maintenance	30
5.3	-	Maintenance of the work area	30
5.4	-	Troubleshooting	31

1 - General information

Congratulations on purchasing the new Spray Tan Face Compressor tanning system!

Spray Tan Face Compressor appliance uses an innovative application system comprising a concentrate of cutting-edge technology packed into a compact machine.

You are about to enjoy the outstanding benefits of a spray tan. Whether you have come across spray tanning for the first time or whether you have already used it before, some fundamental spraying guidelines will help you achieve excellent and successful results. Carefully read this information and follow the simple instructions to get the most out of your new Spray Tan Face Compressor tanning system.

What are the advantages of this appliance?

Spray Tan Face Compressor exploits MicroCore technology to guarantee a precise, uniform and intense tan.

This elegantly designed, practical and silent-running system is suitable for use in all operating environments.

The HP (high-pressure) technology allows tanning lotion to be applied with unmatched precision.

Spray Tan Face Compressor transfers more than 80% of the tanning lotion onto face and décolleté compared with the approximately 30% achieved by turbine systems.

Additionally, application time is much shorter: 15/20 seconds for face and décolleté application. Approximately half the amount of lotion is used compared with a traditional turbine system.

There's no need to wait before getting dressed after application as the mixture of air and lotion vaporises to dry the skin instantaneously.

Lotion application is totally personalised.

1.1 - Using and storing the manual

This manual provides all the information required to use the appliance correctly, independently and safely.

This manual has been written for the owner of the appliance and his/her operators.

The manual is divided into chapters, paragraphs and subparagraphs: the index page is an easy way of finding all subjects of interest.

The material contained in this document is provided exclusively for information purposes and may be modified without notice.

Keep this manual and all attached documents in good condition, legible and complete; store it near the appliance in an accessible place the whereabouts of which is known to all operators.

1.1.1 - Symbols used in this manual

This symbol marks information and warnings which, if not observed, may damage the appliance or compromise operator and/or user safety.

1.2 - Description of staff

This manual has been written for operators who are qualified to use the appliance.

Operators must never perform operations reserved for specialised technicians.

The manufacturer declines all liability for damage deriving from failure to observe the above requirements.

- **Operator qualified to start and operate the appliance:**
specialised person capable of working with the appliance under normal conditions by using the relative controls. He/she must also be able to operate with the appliance under “maintenance” in order to perform simple routine maintenance operations (cleaning, product loading), and start or reset the appliance following an unscheduled stop.

Operators must carefully read the instructions contained in this manual before working on the appliance.

Before starting to use the appliance, the operator must be perfectly familiar with all its controls and characteristics. All the instructions and warnings contained in this manual and the use and maintenance manual must be observed.

- **Users:**
persons subject to tanning sessions, in normal psycho-physical conditions, capable of correctly carrying out the instructions given by the operator responsible for starting and using the appliance.

Women should not be treated (without medical authorisation) during the first three months of pregnancy.

2-Description of appliance

2- Description of appliance

2.1 - Description of appliance and foreseen use

The new Spray Tan Face Compressor (fig.1) has been exclusively designed for tanning using special products supplied by Quadra Medical or by one of its authorised distributors. Other uses or the use of different products invalidate the warranty and relieve the manufacturer and Quadra Medical of all liability.

The appliance uses an innovative application system comprising a concentrate of cutting-edge technology packed into a compact machine. Spray Tan Face Compressor exploits MicroCore technology to guarantee a precise, uniform and intense tan.

The HP (high-pressure) technology allows tanning lotion to be applied with unmatched precision.

This elegantly designed, practical and silent-running system is suitable for use in all operating environments.

Fig. 1

2.2 - Environmental conditions

Do not use the appliance in the presence of:

- humidity;
- dust or dirt.

Use the appliance in well-ventilated areas to allow perfect cooling and to prevent the sprayed product from depositing in the work area or stagnating in the atmosphere.

3 - Safety

3.1 - General warnings

IMPORTANT: before using the appliance, carefully read this operating manual and the compressor operating manual and scrupulously follow the technical operating instructions and advice contained in them. The operator must be perfectly familiar with the position and functions of all control devices and the characteristics of the appliance.

This appliance complies with current safety regulations. Improper use, however, may injure people or damage property.

The director of the centre must instruct staff as regards accident risks, operator and user safety devices, the general accident prevention regulations established in international directives and in the legislation in force in the country where the appliance is used.

The appliance may only be started and operated by trained personnel.

Operators must always scrupulously observe the accident prevention regulations in force in the country where the appliance is used.

Make sure that the ducts and air intakes of the appliance are not obstructed.

Spray Tan Face Compressor may only be used with non-flammable spray solutions approved by Quadra Medical.

Always check the indications and warnings on the spray pack.

Always spray the product in a well-ventilated area.

WARNING: TO PREVENT THE RISK OF ELECTROCUTION, DO NOT EXPOSE THE APPLIANCE TO WATER OR OTHER LIQUIDS.

Never leave the appliance unguarded. Keep it out of reach of children or people who are unfamiliar with the appliance.

3.1.1 - Precautions for users

The operator must make sure that users are in sufficiently good psycho-physical condition to use the appliance and undergo the requested tanning cycle.

Do not use the appliance on women during their first 3 months of pregnancy (without medical authorisation), without first checking the warnings on the spray pack.

Users must observe the following recommendations given by the operator:

- if you have had waxing or depilation treatment, wait at least 12 hours before tanning;
- keep your eyes closed during treatment;
- apply a balsam to your lips;
- do not inhale the product (use nasal filters or hold your breath);
- keep your mouth closed so as not to swallow the product;
- if breast-feeding, protect the nipple area with barrier cream to prevent your child from swallowing product residues;
- cover any open cuts or injuries;
- if your skin is very sensitive, perform a patch-test to check for allergy reactions.

If your skin becomes irritated or swollen, consult a doctor.

If in doubt, consult a doctor.

- on women during their first 3 months of pregnancy (see paragraph 3.1.1);
- in insufficiently ventilated areas.

3.2 - Contraindications

The appliance may not be used:

- for purposes other than those indicated in paragraph 2.1 "Description of appliance and foreseen use";
- for purposes other than those shown on the spray pack;
- after incorrect installation;
- by untrained staff;
- exposed to the weather;
- with liquids other than tanning lotions especially designed for the appliance;

4. - Use

4.1 - Instructions and precautions

- 1 - Once the machine is unpacked, it should be checked over, visually for any signs of damaged during transportation, and defected must be reported to the point of purchase immediately.
- 2 - Place your air compressor on a flat level surface, in a well ventilated area, not in contact with anything other surface it sits on, and the switch position off. Ensure compressor is in a safe working position before connecting power lead to power point. Plug the compressor into a mains supply making sure the voltage is correct.
- 3 - Connect Airbrush hose to pressure regulator outlet with rapid connector
- 4 - Auto- stop function. The maximum pressure preset at the factory is 4 bar/60 PSI. The compressor will automatically cut-off each time the air line reaches the preset range of pressure and will re-start when pressure reach 3 bar/40 PSI
Without proper air tight seals the Auto-stop will not operate
- 5 - Regulate the working pressure (with the air tank full) according to your working necessities, proceeding as follows: Maximum recommended pressure for tanning 2 bar / 30 PSI.

To **INCREASE PRESSURE** turn regulator filter knob (on the filter) clockwise  until the desired pressure is reached, indicated on the pressure gauge.

To **DESCREASE PRESSURE** turn regulator filter knob (on the filter) anti-clockwise  until the desired pressure is reached, indicated on the pressure gauge.

Trial and errors will determinate the best operating pressure for various makes of Airbrush with various type of colour.

- 6 - If the machine will become very HOT during and immediately after use, stop it till is cold (about one our)
- 7 - The motor is thermally protected, automatically it shuts-off when it becomes overheated or lack of voltage. Turn power point "OFF" and turn on again after cooling down.

**If no pressure gauge/regulator is used, do not use this machine to inflate tires of bicycle, car or truck since it may explode.

**For your own safety read instruction before operating this compressor.

4.1.1 - Application precautions

- 1 - If the air compressor does not work when restarted during operation, release the air pressure, and turn-on again. Otherwise may result overloading and damaging the motor.
- 2 - This air compressor has a thermal switch that automatically shuts-off the compressor when it becomes overheated or lack of voltage. Continuous working time max 20 minutes and 20 minutes stop. Prolonged use could cause the motor and compressor cylinder irreparable damage.
- 3 - Do not touch the air compressor during or after use as motor may be hot.
- 4 - Do not use any lubricant, which may damage the air compressor.
- 5 - Use the air compressor in a dry, clean well-ventilated area.
- 6 - Never spray flammable or toxic materials.
- 7 - If the supply cord is damaged, it must be replaced by the manufacturer or its service agent or a competent person in order to avoid hazard.

- 8 - The compressor is provided with an automatic filter for discharging condensate.
A few drops of water may appear under the machine.

4.1.2 - Specifications (see tab. 1)

Specifications	
Voltage	220/240V
Horse power	70 W
Maximum pressure	4 Bar/60 PSI
Flow	16 LPM/0,57 CFM
Tank	LT 1,2-Gal.0,32
Weight	9 Kg
Size	28x18,5x23 cm
Automatic operation	YES
Pressure reducer with filter	YES
Automatic motor save	YES
Noise	50 dB
Exit	1x1/8
Tab. 1	

4.1.3 - Reference to standards

The Face Tan Spray compressor is built in compliance with the current European directives.

		Quadra Medical Srl 30061 - Bassano del Grappa (VI) - ITALY via Strada Cartigliana 143/B tel. +39/0424/33966 - fax +39/0424/33979	
Modello/Model: SPRAYTAN FACE COMPRESSOR		Numero di serie/Serial Number: 9999999	
Alimentazione/Power Supply: 230V 1N~+PE 50Hz 70W 0,3A		Compressore/Compressor 1,2 l - 16 l/min - 4 bar 50 dB(A)	
Mese-Anno/Month-Year: ott-09		 	
Fusibili quadro/Panel Fuses: 1x2A aM			

Fig. 2

4.2 - Tanning procedures

- Remember that the ideal distance of the spray gun from the skin is approximately 12-15 cm.

At the end of the day, dismount the spray gun, wash the parts in hot water, dry them and remount the gun (see paragraph 5.1.1 "Cleaning the spray gun"). Pour any product remaining in the cup into a different container from the original one. Never pour product residue into the original container.

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- Adjust the atomiser vertically so that the air flows horizontally.
- Treat the face starting from above the head and apply three short up and down strokes, making clean and fluid movements. Never point the liquid directly on the face at the beginning of application (to avoid staining).
- Treat the décolleté with a "V" movement, without forgetting the area behind the neck. With the chin raised, apply to the neck. Depending on the phototype of the skin, apply in several stages to obtain a more intense colour and a darker skin.

Remember that any errors can be corrected immediately by removing the lotion with a damp towel and treating the part again.

Fig. 3

Fig. 4

5- Maintenance

5.1 - Gun maintenance

5.1.1 - Cleaning the spray gun

Make sure to eliminate the tanning solution from the spray gun at the end of every day.

Everyday maintenance prolongs the lifetime of the spray gun.

WARNING: TO PREVENT THE RISK OF ELECTROCUTION, DO NOT EXPOSE THE APPLIANCE TO WATER OR OTHER LIQUIDS.

To maintain the spray gun correctly, carry out the following operations in the indicated order:

- 1 - Disconnect the spray gun from the air hose using the relative bayonet coupler (ref. 1 fig. 5).
- 2 - Unscrew the cap (ref. 1 fig. 6) of the pressure cup and remove any residue tanning liquid (ref. 2 fig. 6).
- 3 - Rinse the pressure cup under running water (fig. 7). **Take care not to get other components of the appliance wet in order to avoid the risk of electrocution.**
- 4 - Pour water (possibly hot) into the pressure cup and screw up the cap (fig. 8).
- 5 - Connect the spray gun to the air hose using the relative bayonet connector (ref. 2 fig. 5).
- 6 - Perform a no-load spray cycle to remove all product residues remaining inside the gun (fig. 9); spray the water into a towel or sink until it comes out clean.
- 7 - Unscrew the cap (ref. 1 fig. 9) of the pressure cup and remove any excess water.

Screw the cap (ref. 1 fig. 9) of the pressure cup back on and spray until all the water has left the spray gun.

- 8 - Make sure to keep the threads of the pressure cup and gun body clean.
- 9 - Clean the outside of the spray gun and the pressure cup. Store the spray gun facing upwards, possibly in the work area.


Fig. 5



Fig. 6

Fig. 7

Fig. 8



Fig. 9

5.2 - Compressor maintenance

It is advisable to eliminate regularly (every day) any condensate that gathers inside the air tank, proceeding as follows:

- 1 - Unscrew the cap of the air tank (ref. 1 fig. 10) and let the compressed air flow out slowly.
- 2 - Tilt the compressor slightly so as to assist the drainage of any water residue.
- 3 - Screw the cap back onto the air tank (ref. 1 fig. 10).

5.3 - Maintenance of the work area

STATUS OF APPLIANCE

When maintaining the work area, the appliance must be switched off by turning the general switch to “O” (OFF) and disconnecting the power cable from the mains socket.

To keep the spraying area looking professional, clean the floor and the walls at the end of each working day to remove all spray residues.

Clean with a cloth soaked in water and a mild non-flammable disinfectant.

WARNING: TO PREVENT THE RISK OF ELECTROCUTION, DO NOT EXPOSE THE APPLIANCE TO WATER OR OTHER LIQUIDS.

5.4 - Troubleshooting (see tab. 2)

PROBLEM	PROBABLE CAUSE	SOLUTION
The spray gun drips.	Impurities between the needle and the nozzle interrupt the hermetic seal.	Perform a spray gun cleaning cycle (see paragraph 5.1.1 "Cleaning the spray gun"). Contact the manufacturer if the problem persists.
Drop- or oval-shaped jet. 	Dirt on the tip of the nozzle or in the air circuit.	Carefully clean the outside of the spray gun with a damp cloth. If the problem persists, perform a spray gun cleaning cycle (see paragraph 5.1.1 "Cleaning the spray gun").
The material comes out in spurts, or bubbles in the pressure cup.	<ul style="list-style-type: none"> - The nebulised air enters the liquid tank through the liquid channel. - The air nozzle is not completely tight, the air circuit is obstructed. The seat is faulty or the nozzle set is damaged. 	<ul style="list-style-type: none"> - Tighten the various parts. - Clean the various parts. - Contact the manufacturer for detailed instructions.
<hr style="width: 10%; margin-left: auto; margin-right: 0;"/> Tab. 2		

That'so
cosmetics

QUADRA MEDICAL srl

Quadra Medical Srl

Strada Cartigliana, 143/b - Bassano del Grappa (VI) Italy - Tel. +39 0424 33988 - Fax +39 0424 33979
info@qmed.it - www.qmed.it - info@thatso.it - www.thatso.it

Cod. 0452700