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
Phone number: [Empty]

Established date: 01/01/2008 0:16:27 AM

Measurement unit: Metric system Imperial system


Weight: 130 Kilogram

Height: 170 Centimeter

Default date: 01/01/2008 0:16:27 AM

Save & Analysis

Appointment Preliminary Search BT-505PC Introduction Instruction Backup EXIT

Body Composition Analyzer BT-505PC

Preliminary

Name: John

ID number: 0123456789

Sex: Female

Birth date: Jan 1, 1970

Medical record number: [Empty]

Phone number: [Empty]

Established date: 01/01/2008 0:16:27 AM

Measurement unit: Metric system Imperial system


Weight: 130 Kilogram

Height: 170 Centimeter

Default date: 01/01/2008 0:16:27 AM

Height: 170 cm
 Weight: 130 kg
 METRIC UNITS
 Fat: 15.5%
 Muscle: 15.1%
 Bone: 11.1%
 Water: 68.3%
 Fat-free mass: 111.5 kg
 Fat-free mass index: 37.4
 Fat-free mass percentage: 84.5%
 Fat-free mass ratio: 0.845
 Fat-free mass ratio squared: 0.714
 Fat-free mass ratio cubed: 0.603
 Fat-free mass ratio to the power of 4: 0.509
 Fat-free mass ratio to the power of 5: 0.428
 Fat-free mass ratio to the power of 6: 0.359
 Fat-free mass ratio to the power of 7: 0.299
 Fat-free mass ratio to the power of 8: 0.251
 Fat-free mass ratio to the power of 9: 0.211
 Fat-free mass ratio to the power of 10: 0.177
 Fat-free mass ratio to the power of 11: 0.148
 Fat-free mass ratio to the power of 12: 0.123
 Fat-free mass ratio to the power of 13: 0.101
 Fat-free mass ratio to the power of 14: 0.082
 Fat-free mass ratio to the power of 15: 0.066
 Fat-free mass ratio to the power of 16: 0.053
 Fat-free mass ratio to the power of 17: 0.042
 Fat-free mass ratio to the power of 18: 0.034
 Fat-free mass ratio to the power of 19: 0.028
 Fat-free mass ratio to the power of 20: 0.023

Appointment Preliminary Search BT-505PC Introduction Instruction Backup EXIT

Body Composition Analyzer BT-505PC

Preliminary

Name: John

ID number: 0123456789

Sex: Female

Birth date: Jan 1, 1970

Medical record number: [Empty]

Phone number: [Empty]

Established date: 01/01/2008 0:16:27 AM

Measurement unit: Metric system Imperial system


Weight: 130 Kilogram

Height: 170 Centimeter

Default date: 01/01/2008 0:16:27 AM

Fat: 15.5%
 Muscle: 15.1%
 Bone: 11.1%
 Water: 68.3%
 Fat-free mass: 111.5 kg
 Fat-free mass index: 37.4
 Fat-free mass percentage: 84.5%
 Fat-free mass ratio: 0.845
 Fat-free mass ratio squared: 0.714
 Fat-free mass ratio cubed: 0.603
 Fat-free mass ratio to the power of 4: 0.509
 Fat-free mass ratio to the power of 5: 0.428
 Fat-free mass ratio to the power of 6: 0.359
 Fat-free mass ratio to the power of 7: 0.299
 Fat-free mass ratio to the power of 8: 0.251
 Fat-free mass ratio to the power of 9: 0.211
 Fat-free mass ratio to the power of 10: 0.177
 Fat-free mass ratio to the power of 11: 0.148
 Fat-free mass ratio to the power of 12: 0.123
 Fat-free mass ratio to the power of 13: 0.101
 Fat-free mass ratio to the power of 14: 0.082
 Fat-free mass ratio to the power of 15: 0.066
 Fat-free mass ratio to the power of 16: 0.053
 Fat-free mass ratio to the power of 17: 0.042
 Fat-free mass ratio to the power of 18: 0.034
 Fat-free mass ratio to the power of 19: 0.028
 Fat-free mass ratio to the power of 20: 0.023

Appointment Preliminary Search BT-505PC Introduction Instruction Backup EXIT

Body Composition Analyzer BT-505PC

Analysis Chart

Who's body fat content: weight: 15.5%

Who's body fat content: series: stage: 2%

Who's body fat content: lean: 84.5%

Who's body fat content: issue: content: 36.8%

Who's body fat content: moisture: content: 831 (33.8%)

Who's body fat content: metabolic: rate: 146 cal

Suggested idea: fat: content: 15.5%

Suggested idea: body: weight: 68 kg

Actual body: weight: 48 kg

Print

Name: John
 ID number: 0123456789
 Height: 170 Centimeter
 Sex: Female
 Age: 33
 Date of birth: 01/01/2008 0:16:27 AM
 Fat-free mass: 111.5 kg
 Fat-free mass index: 37.4
 Fat-free mass percentage: 84.5%
 Fat-free mass ratio: 0.845
 Fat-free mass ratio squared: 0.714
 Fat-free mass ratio cubed: 0.603
 Fat-free mass ratio to the power of 4: 0.509
 Fat-free mass ratio to the power of 5: 0.428
 Fat-free mass ratio to the power of 6: 0.359
 Fat-free mass ratio to the power of 7: 0.299
 Fat-free mass ratio to the power of 8: 0.251
 Fat-free mass ratio to the power of 9: 0.211
 Fat-free mass ratio to the power of 10: 0.177
 Fat-free mass ratio to the power of 11: 0.148
 Fat-free mass ratio to the power of 12: 0.123
 Fat-free mass ratio to the power of 13: 0.101
 Fat-free mass ratio to the power of 14: 0.082
 Fat-free mass ratio to the power of 15: 0.066
 Fat-free mass ratio to the power of 16: 0.053
 Fat-free mass ratio to the power of 17: 0.042
 Fat-free mass ratio to the power of 18: 0.034
 Fat-free mass ratio to the power of 19: 0.028
 Fat-free mass ratio to the power of 20: 0.023

Appointment Preliminary Search BT-505PC Introduction Instruction Backup EXIT

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INTRODUCTIONS

The accepted method of measuring total body water is based on the dilution of isotopes of water Deuterium or Tritium as it is metabolized in the human subject over time. And the traditionally estimate body composition is to use electrical impedance plethysmography.

However, there is a need for a total body composition measurement technique that is rapid, safe, noninvasive and sufficiently accurate and convenient to permit its use in body composition. Recently, we used a new approach to human body composition estimation that appears to fulfill many of these requirements.

This new technique entails measuring total body electrical resistance impedance. The amount of water, fat free mass in man can be accurately estimated by measuring total body electrical resistance. Specifically, four electrodes are placed over metacarpals and metatarsals where a 50 KHz current is introduced. The detected signal level when corrected for subject height will be an index of total body water and fat free mass.

BODY COMPOSITION

The body is made up of both lean and fat tissue, as well as water.

Together they equal the weight you see on your weight scales.

Lean tissue includes muscle, bone, vital organs, etc; both fat and lean tissues contain water. To have a healthy, well-functioning body, it is extremely important that the lean and fat tissue weights are properly balanced.

Now the BT-905PC will help you to manage your nutrition, fitness and weight. It not only can tell if you are at the right weight but also tell you whether or not you are nutritionally healthy.

HOW CAN THE BODY COMPOSITION ANALYZER HELP ME GET AND STAY FIT?


Listed below are some of the many ways BT-905PC will help maximize your personal level of fitness-whether you want to lose weight, build muscle, tone and trim, or simply be sure you are maintaining your current level of fitness

- 1) It allows you to check the effectiveness of your weight loss program by telling you exactly where your weight loss is coming from-fat, lean, or water.
- 2) It establishes a baseline for your body composition so that nutritional and fitness improvements can be monitored on a regular basis.
- 3) It provides possible answers to weight "plateaus" -a common plague of dieters.
- 4) It helps you monitor improvements in muscle mass as a result of exercise and/or diet, and it warns you of muscle loss caused by inactivity, the aging process, or illness.
- 5) It help explain "mysterious" weight gains or losses associated

FEATURES

- 1) Results - FAT%, LEAN %, Water, BMR, Target Weight and Target FAT%.
- 2) Accurate - The results accurate same as the instruments cost thousands of US dollars.
- 3) Time saving - Testing and results take less than 2 minutes.
- 4) Private - You can perform the test almost anywhere.
- 5) Portable - Small size and light weight.
- 6) Warranty - One year limited warranty.
- 7) Durable.

INDICATORS AND CONTROLS

Power Switch - Power ON/OFF control.

Lead Connector - There are two jacks for two lead connectors to plug in. The jacks are on one side of the machine.

Sensor Cable - Each lead wire has two color red and black clips, each clip will connect to one sensor pad.

Key Pad - Input data age, height, weight.

Sensor Pads - Connect the machine and body to get signal.

DC-IN Jack - For DC 9V adaptor.

GETTING YOURSELF READY FOR TESTING

To be sure your test results are accurate...

You must follow the guidelines listed below AS EXACTLY AS POSSIBLE to be sure that you will consistently obtain the most valid and reliable results possible you're your BT-905PC Instrument.


- 1) **Meals.** Wait at least 5 hours after a meal before performing the test.
- 2) **Hydration.** Do not consume large amounts of water just prior to testing. Some medications - particularly diuretics - may also affect test results.
- 3) **Coffee.** Do not drink coffee for 4 hours before testing. Coffee is a diuretic and may affect the hydration level of your body.
- 4) **Alcohol.** Do not drink alcoholic beverages for 24 hours prior to testing. Alcohol is a diuretic. In your bloodstream, it also conducts electricity in such a way as to interfere with an accurate test.
- 5) **Exercise.** Do not test your body composition for at least 12 hours after exercise. This is because exercise has a significant effect on your body's water content.
- 6) **Weight.** Check your weight on a reliable pair of scales immediately before testing. This is particularly important for those who tend to be quite lean. The scientific formulas used to accurately calculate your percents of lean and fat are dependent on an accurately calculate your percents of lean and fat are dependent on an accurately scale weight. Even using yesterday's weight is not close enough to assure a valid test.

Just before the test ... Cleanse your skin.

Using a washcloth with soap and water or a cotton pad with rubbing alcohol, cleanses the areas where you will be placing the sensor pads. (See drawing on next page). Scrub the areas until the skin becomes pink but not sore or irritated. This will remove all dead skin, perspiration, and any excess oils that could prevent good contact with the sensor pads. Dry the skin with a clean cloth.

SENSOR PAD PLACEMENT

(Right foot & right hand or usually use side)


Stick on the sensor pads.

Sensor pads should be placed on one hand and one foot, the side use usually, as shown in the accompanying diagrams.

On the back of the hand, place a sensor pad on the third knuckle of the middle finger.

Place a second pad on the arm just behind the wrist at the crease observed when the hand is bent back. The edge of the sensor pad should just touch the crease.

The pad on the foot should be applied just behind the place where the big toe and the first toe join. The position for the sensor pad at the ankle is found by placing your thumb and forefinger on the two protruding bones on either side of the ankle and then bringing the fingers forward to the midway point between the two bones.

Place the sensor pad at the bend of the ankle, in line with the shinbone.

Refer to diagrams for proper placement, and be sure to be consistent each time you test.

Connect the sensor cable clips to the sensor pads.

There are two long sensor cable wires.

At the ends of each of these wires are two clips - a red clip and a white clip (a total of four clips).

On the foot - Using either of the main wires - the white clip should be attached to the sensor pad near the big toe; the red clip should be attached to the sensor pad on the ankle.

On the hand, the white clip on the remaining wire should be attached to the sensor near the knuckle; the red clip should be attached to the sensor at the wrist.

- It may help to remember that the red clips are always placed closer to your heart than the white clips.

GET YOUR BODY IN POSITION FOR THE TEST

For testing, lie flat on your back on a hard, flat surface (Preferably uncarpeted floor or table). Do not touch your arms to your body while testing. Arms should be at your side, palms flat on the floor, a comfortable distance from your trunk. Legs should be unbent and slightly apart so that your body is perfectly symmetrical when the display counts down for your test.


BODY COMPOSITION ANALYZER SYSTEM INSTALLATION

Installation procedure

- 1) Turn on the computer.
- 2) Put the CD of the Body Composition Analyzer software into the CD-ROM drive.
- 3) The screen will display main windows file. In the main windows, it will display the Body Composition Analyzer content. Click the "**setup.exe**", the main unit can automatically carry on the installation.


- 4) At the same time, the screen will appear "Welcome to the BT-905PC Setup Wizard" chart. This will install BT-905PC on your computer. It is recommended that you close all of other applications before continuing. Click "**Next**" to continue, or "**Cancel**" to exit Setup.


- 5) Setup Information. If not first time install, before install the program, please backup the data file, because when install the new program will cover the old data file. After finish the installs, please use the icon "**Export/Import**" to restore the old data to Body Composition Analyzer system. If the first time install, you can install directly.


The preset pathway of Body Composition Analyzer system is “**C:\Program files\ Bodily composition analyzer\BT 905- Data\BT-905.mdb**”. When you are ready to continue with Setup, click “**Next**”.

6) **Destination Location.** Setup will install BT-905PC into the following folder “**C:\Program files\BT-905PC**”. To continue, click “**Next**”. If you would like to select a different folder, click “**Browse**”.


7) **Select Star Menu Folder.** Setup will create the program’s shortcuts in the following Start Menu folder “**BT-905PC**”. To continue, click “**Next**”. If you would like a different folder, click “**Browse**”.


8) **Select Additional Tasks.** Select the Additional tasks you would like Setup to perform while installing BT-905PC “**Create a desktop icon**” or/and “**Create a Quick Launch icon**”. Then click “**Next**”.


Set the BT-905PC as a shortcut to the desktop icon or Quick Launch icon, user enter the WINDOWS desktop after the computer is on. Click the icon and enter the process of Body Composition Analyzer system.


9) Ready to Install. Setup is now ready to begin installing BT-905PC on your computer. Click “Install” to continue with the installation, or click “Back” if you want to review or change any settings.

10) Installing and Extracting files. Please wait while setup installs BT-905PC on your computer. If the screen display “Error”, it is recommended that you keep the existing file. Click “Yes” to keep the existing file and to continue install, or “No” to copy file and to continue install.

11) Completing the BT-905PC Setup Wizard. Setup has finished installing BT-905OC on your computer. The application may be launched by selecting the installed icon "Launch BT-905PC.EN", then click "Finish" to exit setup.

Note: There is no need to restart the computer, it can be operate right the way.

CONNECTION DIAGRAM OF BODY COMPOSITION ANALYZER INSTRUMENT

- 1) Connect the RS-232 wire to RS-232 TRANSMIT jack or USB wire to USB TRANSMIT jack of BT-905PC main unit.
- 2) Connect the 9PIN (RS-232) end or USB port to the computer.
- 3) Connect one end of adaptor into the DC-IN jack and plug the other end to 110VAC or 220VAC socket.
- 4) Connect the sensor cable clips to the sensor pads.

OPERATION PROCEDURE

Connection PC and BT-905PC

- 1) Connect PC and BT-905PC.
- 2) Turn on the computer and waiting for the processing into the WINDOW desktop.
- 3) Connect the RS-232 wire to RS-232 TRANSMIT jack or USB wire to USB TRANSMIT jack of BT-905PC main unit.
- 4) Connect the 9PIN (RS-232) end or USB port to the computer.
- 5) Connect one end of adaptor into the DC-IN jack and plug the other end to 110VAC or 220VAC socket.
- 6) Connect the sensor cable clips to the sensor pads.

Stick on the sensor pads

There are two long sensor cable wires. At the ends of each of these wires are tow clips - a red clip and a white clip (a total of four clips).

Sensor pads should be placed on one hand and one foot, the side use usually, as shown in the accompanying diagrams.



Operation the Body Composition Analyzer system

- 1) Switch ON the BT-905PC and warm up the unit for at least 10 minutes before testing.
- 2) Move the mouse to the desktop “BT-905PC” icon and continue click left key of the mouse twice into the Body Composition Analyzer system go to main menu.
- 3) If any [error] message appears , please check the electrodes in proper positions, make sure lead wires connectors are firmly seated and restart the BT-905PC main unit.
- 4) Once again for the Body Composition Analyzer system.

Appointment Steps

- 1) After executing the  “Body Composition Analyzer” go to main menu, click the





Appointment



- 2) Or after you enter the sub-main page, you may also click the **Appointment** down at the function bar there.



3) Next, enter “Appointment” page, please click **Add** or choose the examinee and measurement unit, and then input height and weight.




4) Next, Click **Save & Analysis** and wait for approximately one and a half min for processing.


5) If any “error” message appears, please click **Save & Analysis** once again for the analysis chart.



6) Click **Print** to print out the analysis chart.

Preliminary Steps



1) After executing the "Body Composition Analyzer" go to main menu, click the **Preliminary**.





2) Or after you enter the sub-main page, you may also click the **Preliminary** down at the function bar there.




3) Next, enter “**Preliminary**” page, please input basic personal data, height and weight, and select measurement unit.


4) Next, Click **Save & Analysis** and wait for approximately one and a half min for processing.

5) If any “error” message appears , please click **OK** , then click again **Save & Analysis** for the analysis chart.


6) Click **Print** to print out the analysis chart.


Read Data Steps / Search Steps



1) After executing the  “Body Composition Analyzer” go to main menu, click the

2) Or after you enter the sub-main page, you may also click  **Search** on the at the function bar there.

3) Next, enter “Search” page. You could search by Name, ID and medical history number accordingly.

Body Composition Analyzer -- Search --

BT-800PC **Body Composition Analyzer**

Search

Name: Medical record number: Measure unit: Imperial system:

ID number: Phone number: Weight: Kilograms:

Sex: Establishment: Height: Centimeter:

Birth date: Detect date: Detect case:

Name	ID number	Medical record	Sex	Phone number	Birth date	Date
Isoko	0123456789		Female	0301807	03/18/73	01/20/08
Isoko	0123456789		Female	0301807	03/18/73	01/20/08
Ry	9		Male	01/18/55	01/18/55	01/20/08
Ry	04		Female	01/18/54	01/18/54	01/20/08
Isoko	1		Female	01/18/55	01/18/55	01/20/08
Isoko	099		Male	01/18/57	01/18/57	01/20/08
Isoko	345674		Male	01/18/75	01/18/75	01/20/08
Isoko	100		Female	01/18/81	01/18/81	01/20/08
Isoko	112		Female	01/18/81	01/18/81	01/20/08

Save & Analysis Record

First Record Previous Next Last Record

Appointment Preliminary Search BT-800PC Introduction Instruction Backup EXIT

4) After search, click

Record

Body Composition Analyzer -- Search --

BT-800PC **Body Composition Analyzer**

Search

Name: Medical record number: Measure unit: Imperial system:

ID number: Phone number: Weight: Kilograms:

Sex: Establishment: Height: Centimeter:

Birth date: Detect date: Detect case:

Name	ID number	Medical record	Sex	Phone number	Birth date	Date
Isoko	111		Male	01/18/81	01/18/81	01/20/08
Isoko	312		Male	01/18/81	01/18/81	01/20/08
Isoko	300		Male	01/18/81	01/18/81	01/20/08
Isoko	Shoek		Male	05/18/75	05/18/75	01/20/08

Save & Analysis **Record**

First Record Previous Next Last Record

Appointment Preliminary Search BT-800PC Introduction Instruction Backup EXIT

5) All history records will be listed of the specified.

Body Composition Analyzer -- Search --

BT-800PC **Body Composition Analyzer**

Search

Name: Medical record number: Measure unit: Imperial system:

ID number: Phone number: Weight: Kilograms:

Sex: Establishment: Height: Centimeter:

Birth date: Detect date: Detect case:


Detect date	Age	Height	Weight	fat FAT	muscle FAT	fat LEAN	muscle LEAN	Protein Water
01/20/2008 0:00:01	13	178	80	11.3	52.7	15.7	37.0	81

Analysis Chart Multiple Comparison



First Record Previous Next Last Record

Appointment Preliminary Search BT-800PC Introduction Instruction Backup EXIT

6) Choose the record and click **Analysis Chart**.

7) The animated analysis chart will appear according to your appointed record.

Add Data Steps

After “Search Steps” 1) to 4), click **Add**. Please follow the “Preliminary Steps” 3) to 6).

Modify Data Steps



1) After “Search Steps” 1) to 4), please click **Modify** and choose the examinee and measurement unit, and then input height and weight.

Delete Data Steps





After “Search Steps” 1) to 5), please click **Delete** to delete data.

Multiple Comparison Steps


1) Multiple history: after “Search Steps” 1) to 4), click **Multiple Comparison**

2) Choose 2-3 data from the records and click **Multiple Comparison**.



3) Click **Print** to print out the analysis chart.


Backup Steps


1) After executing the  “Body Composition Analyzer” go to main menu: click the .



2) Or after you enter the sub-main page, you may also click the **Backup** at the function bar there.


3) Next, enter **Backup** page. In the path of source data, the system will automatically bring into the path of database located.

4) When file exporting, please choose source database and purpose route in advance, then click **Export**. The dialogue window with message of file export successfully will appear when completing file exporting.


5) When file importing, please choose your source database and target route in advance. The target route is at C:\Program Files\Body Composition Analyzer\BT-905PC-Data. Then click **Import**.

6) While importing the file, the importing file will automatically replace the current database and file name will be changed to a system use file name.

7) When the file import is done, a dialogue window with message file import successfully will appear.

Note: If the BT-905PC machine idled for more than 5 minutes, the LCD will show “+++++++” and enter the power saving mode. You need to turn OFF the machine first and then turn ON again to begin a new analysis. It doesn't need to warm up the unit for extra 10 minutes when you re-start the machine and continue to operating the test.

NOTE

- 1) If you use English units, when you input your height the units is in inches, for example: If you are 5 feet 1/2 inch tall, you would key in 60.5 inches tall. If you are an even 6 feet tall, you should key in 7.2 inches. In the event that you are an even number of inches tall, be sure to include a zero as your final entry.
- 2) To be sure your results are accurate, do not guess at your weight or rely on your memory. Please indicate your weight just before your test, and use the same scale each time you test, if possible.
- 3) If you accidentally press a wrong number or letter, do not hit the **Save & Analysis**, then the correct response, and then click the. If you press a wrong number of letters and click **Save & Analysis**, you must restart the entire test.

INTERPRETING THE RESULTS

lbs FAT

Pounds of FAT, the actual weight of fat on your body. There are approximately 3500 calories stored in each pound of body fat.

lbs LEAN

The pounds of tissue in your body mostly muscle and vital organs represent the amount of your weight that is not fat.

Pints of water

A measure of the pints of water in your body can be very helpful in explaining rapid weight loss or gain, A Pints of water weights about one pound. Most of the body's fluid is contained in the lean body mass, which normally consists of 70~75% water.

% FAT

The percent of fat weight in total body weight.

%LEAN

The percent of tissue weight in total body weight.

BMR

Basal metabolic rate refers to the number of calories your body uses each day just to “**stay alive**” - breathing, pumping blood, thinking, and so on. It doesn’t include calories your burn by exercising.

Your basal metabolic rate is directly related to the amount of lean tissue in your body. The more lean tissue you have, the higher your metabolic rate - in other words, the more calories you will “**use up**” just by being alive.

Target Weight

Target weight is simply a weight goal to aim for under optimal conditions. It is based on the assumption that all the weight you lose will be fat and that lean tissue will be neither lost nor gained. Under an ideal condition, all weight loss would be fat and a small amount of lean would be increased. This is why your target weight may change slightly as you alter your diet and exercise plan. With a healthy lifestyle, eventually your actual weight and target weight should match.

Target Percent Fat

Target percent fat is a value which represents a “normal” percent for your sex and age.

It does not necessarily represent the “best” percent fat for you. For example, very fit individuals often fall well below these target values. The importance of the target percent fat is that it is a worthy goal for those who tend to be over fat and is generally considered healthy. This value can be seen in the display only when Program, Target Weight is entered.

TAKING CARE OF YOUR BT-905PC

Do

- 1) Cleaning - Periodically you may wipe BT-905PC with a damp cloth without injury to the internal components of the instrument.
- 2) When removing lead connectors, grasping them by the plug, not the cord.

Don't

- 1) Don't expose to excessive amounts of moisture and definitely don't immerse in water.
- 2) Don't leave in direct sunlight for extended periods of time.
- 3) Don't drop or hit against other objects. Treat it as you would other sophisticated electronic devices such as a calculator or recorder.

TROUBLE SHOOTING

No Display

If you turn on you BT-905PC and there is no display in the window, you should

- 1) **Check power:** If the above attempts fail, return the machine to the manufacturer.
- 2) **“error” message:** This error message pop up results from one of several possibilities. The most likely cause is failure to follow protocol closely enough. Carefully re-read the sections in your manual entitled “**Getting Yourself Ready for Testing**” and “**Operate Procedure**” Then retest.
- 3) Recheck your exact scale weight and your height; then reperform the test. This is absolutely necessary for thin individuals in particular. Weight fluctuates from day to day enough to produce invalid data if not entered correctly.
- 4) A Very small percent of the extremely thin and tall population may not be possible to measure using bioresistance methods. This would be a very rare situation, and the Invalid Data message is probably due to improper testing protocol rather than to being immeasurable lean. Be sure to follow protocol.
- 5) Dehydration. Body water levels that are unusually low can cause unusual readings. If you are an endurance athlete or if you take diuretics, this may invalidate your data. To receive a valid test, be sure you are stable in your water consumption.

Check Probes

This message indicates that the electrical signals from your body are not strong enough for your BT-905PC to read and interpret.

- 1) Check the sensor pads. Make sure they fit snugly against your skin. Did you remember to cleanse your skin thoroughly before application?
- 2) Check the cable connector to be sure it is properly inserted in the jack.
- 3) Check the sensor cable clips to be sure each one is clipped firmly onto the sensor pads.

When these checks have been completed, retake the test.

SPECIFICATION

Power source: operated by 9VDC adaptor

Size: 66x176x208 mm

Disposal: *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.*

For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product.

This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty.

The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use.

GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc.

The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.