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MAINTENANCE
&
USERS MANUAL

STUFE ED
INCUBATORI
52-80-120-240-400 Lt

OVENS AND
INCUBATORS
52-80-120-240-400 Lt

FALC INSTRUMENTS s.r.l.

Treviglio (BG) Italy

INDICE

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INDEX

1.0 SAFETY.....	page 1
1.1 GENERAL.....	page 1
1.2 OPERATOR QUALIFICATION.....	page 1
1.3 INSTRUMENT SAFETY SYSTEMS.....	page 1
2.0 PRELIMINARY OPERATIONS.....	page 1
1.1 VERIFICATION OF TRANSPORTATION DAMAGES.....	page 1
2.2 IN CASE OF DAMAGE.....	page 2
2.3 CLEANING THE INSTRUMENT.....	page 2
3.0 INSTRUMENT CHARACTERISTICS.....	page 2
4.0 CONTROLS.....	page 4
5.0 INSTALLATION AND USE PROCEDURES.....	page 5
5.1 PHYSICAL CHARACTERISTICS.....	page 5
5.2 STARTING.....	page 5
6.0 MANUAL REGULATION FOR AIR VENTILATION AND ITS USE.....	page 6
7.0 VENTILATION.....	page 6
8.0 FORESEEN AND UNFORESEEN USES.....	page 7
9.0 ACCESSORIES.....	page 7
10.0 MAINTENANCE.....	page 8
10.1 INSTRUMENT MAINTENANCE.....	page 8
10.1.1 ORDINARY MAINTENANCE.....	page 8
10.1.2 PROGRAMMED MAINTENANCE.....	page 8
10.2 CLEANING.....	page 8
10.2.1 CLEANING TOOLS AND PRODUCTS.....	page 8
10.3 TRANSPORT.....	page 8
10.4 STORAGE.....	page 8
11.0 PROBLEM SOLVING	page 9
12.0 MANUFACTURER SPARE PARTS	page 10
13.0 DECLARATION OF CALIBRATION	page 10
14.0 ELECTRICAL SCHEME.....	page 11
15.0 DEMOLITION.....	page 14
15.1 OPERATOR QUALIFICATION.....	page 14
15.2 DEACTIVATION AND DEMOLITION	page 14
15.3 DEACTIVATION AND DEMOLITION PROCEDURE	page 14
15.4 POST-DEMOLITION RESIDUAL RISKS	page 14

This manual must be kept by user for the whole technical life of the instrument in order to be immediately available in case of need. This manual must always accompany the instrument, even if it is sold to another user.

Keep this manual in a well protected and easy accessible place in order to permit a quick consultation in case of need. We also suggest to keep this manual in a clean place in order to avoid the formation of stains that could make it difficult to read the text.

1.0 SAFETY

1.1 GENERAL

Falc Instruments declares that the engineering and construction of this instrument comply with the criteria and notions that satisfy the basic safety requirements under the Low Voltage Directive 2006/95/CE.

A meticulous risk analysis, carried out by the manufacturer, made it possible to eliminate most of the risks that were connected to both foreseen and reasonably predictable instrument use conditions.

Full documentation of all the adopted safety measures can be found in the instrument technical file, held by the manufacturer.

The manufacturer recommends that users strictly observe all instructions, procedures and recommendations contained in this manual, and stick to the current workplace safety regulations, even when using the protection devices.

FALC INSTRUMENTS declines any liability for any damages to people, animals or objects, resulting from an improper use of the instrument or from use by unqualified operators.

1.2 OPERATOR QUALIFICATION

All the instrument installation operations must be carried out exclusively by trained, qualified and authorized personnel, after studying the complete instruction manual.

1.3 INSTRUMENT SAFETY SYSTEMS

The instrument has 2 different types of safety systems: a digital one which is set automatically 2°C or 5°C over the set point and a mechanical one which is activated when the instrument temperature exceeds the maximum limit of 120°C or 350°C.

- FIRST SAFETY SYSTEM (digital safety system)

When temperature is 2°C (ICT models) or 5°C (STF models) higher than the set point, the led  turns on automatically.

This safety system disconnects automatically as soon as the temperature returns to normal. This safety system is also used to keep the internal temperature of the chamber constant and to prevent thermal imbalances when opening the instrument door.

- SECOND SAFETY SYSTEM (mechanical safety system)

When temperature exceeds 120°C (ICT) or 350°C (STF-F e STF-N) the warning light  turns on and the heating system stops working.

Wait until the internal temperature of the chamber decreases and press the SAFETY button to restart the instrument.

The instrument will start working again.

ATTENTION → Please contact the technical assistance if the warning signals appear again.

For further information read chapter “PROBLEM SOLVING” on page 9.

FALC INSTRUMENTS declines any liability for any damages to people, animals or objects, resulting from failure to comply with precautions contained in this manual or from failure to use the prescribed PPD (personal protection devices).

2.0 PRELIMINARY OPERATIONS

2.1 VERIFICATION OF TRANSPORT DAMAGES

Carefully unpack the instrument from the shipping container and perform a visual inspection in order to check its conditions. Any deformation of visible parts of the instrument indicates that it has suffered damage during transportation that could compromise normal working. Please pay particular attention to the conditions of the following parts of the instrument:

- Instrument structure.
- Controls and warning lights.
- Make sure all screws and bolts are tightened properly.

2.2 IN CASE OF DAMAGE

Damages due to transport must be attributed to the carrier and immediately communicated to the manufacturer or his representative. Send a claim letter to the carrier.

2.3 CLEANING THE INSTRUMENT

Get rid of any external dust or dirt that may have accumulated during transport.

3.0 INSTRUMENT CHARACTERISTICS

COMMON CHARACTERISTICS

Capacity		52Lt	80Lt	120Lt
Internal dimensions LxWxH	mm	390x350x390	450x400x450	500x470x520
External dimensions LxWxH	mm	705x570x610	775x610x675	825x680x740
Supplied shelves	n°	2	2	2
Allowed number of shelves	n°	6	7	8
Maximum load per shelf*	Kg	20		
Maximum load	Kg	50	75	100
Weight	Kg	55	60	68

Capacity		240Lt	400Lt
Internal dimensions LxWxH	mm	600x500x800	700x580x1000
External dimensions LxWxH	mm	925x710x1030	1025x790x1230
Supplied shelves	n°	2	2
Allowed number of shelves	n°	12	15
Maximum load per shelf*	Kg	20	
Maximum load	Kg	120	
Weight	Kg	80	100

*Weight was calculated distributing it uniformly over the shelf length.

Electric supply 220/240V - 50/60Hz	Safety class 3.1	Protection class IP 42
---------------------------------------	---------------------	---------------------------

- Steel structure painted in oven with antacid epoxy dust.
- Internal chamber made of INOX steel.
- Door with double locking system.
- Ø45mm backside hole for smoke discharge or cables passage.
- Ø10mm upper side hole for thermometer insertion.
- PID microprocessor-based temperature control.
- Timer can be set from 0 to 12 hours with possibility to set it to infinite time.
- Incubators are provided with a second glass door which allows users to inspect the internal chamber without causing any thermal imbalances.

FORCED VENTILATION OVENS

Model		STF-F 52Lt	STF-F 80Lt	STF-F 120Lt	STF-F 240Lt	STF-F 400Lt
Working temperature	°C	from ambient+5 to +300				
Accuracy	70°C	±1,0	±1,0	±1,0	±1,5	±1,5
	150°C	±1,5	±1,5	±1,5	±2,0	±2,0
	300°C	±2,0	±2,0	±2,5	±2,5	±3,0
Heating power	W	1200	1400	1600	2600	3400
Code		720.2090.05	720.2090.08	720.2090.12	720.2090.24	720.2090.40

NATURAL VENTILATION OVENS

Model		STF-N 52Lt	STF-N 80Lt	STF-N 120Lt	STF-N 240Lt	STF-N 400Lt
Working temperature	°C	from ambient+5 to +300				
Accuracy	70°C	±1,5	±1,5	±2,0	±2,5	±3,0
	150°C	±2,5	±2,5	±3,0	±3,5	±3,5
	300°C	±3,0	±3,0	±3,5	±3,5	±4,0
Heating power	W	1200	1400	1600	2600	3400
Code		720.2080.05	720.2080.08	720.2080.12	720.2080.24	720.2080.40



INCUBATORS


Model		ICT 52Lt	ICT F 80Lt	ICT 120Lt	ICT 240Lt
Working temperature	°C	from ambient+5 to +80			
Accuracy	37°C	±0,4	±0,4	±0,8	±0,8
Heating power	W	300	350	400	650
Code		720.2070.16	720.2070.18	720.2070.20	720.2070.22

SUPPLEMENTARY GRATING SHELVES/SPARE PARTS

Capacity	52Lt	80Lt	120Lt	240Lt	400Lt
Code	720.2091.05	720.2091.08	720.2091.12	720.2091.24	720.2091.40

4.0 CONTROLS

①	Regulates the quantity of air passing through the back flue.
Display	Shows current or set temperature.
	Turns on when resistances are working or during the setting phase.
	Turns on when temperature is lower than the set temperature.
	Turns on when temperature is equal to the set temperature.
	Turns on when temperature is higher than the set temperature.
	Turns on during the setting phase or when temperature is 2°C or 5°C higher than the set point temperature.
	Confirms that changes are made.
	Increases the set temperature.
	Decreases the set temperature.
	Turns on when the instrument is working.
	Turns on when the fan is working.
	Turns on in case of alarm.
②	Sets the time.
	Serves to manually restart the second safety system.
③	Regulates the fan speed.
IG	Serves to turn the instrument on and off.
④	Fuse door.

For 240Lt and 400Lt models.

①

5.0 INSTALLATION AND USE PROCEDURES

5.1 PHYSICAL CHARACTERISTICS

All operations described in this manual must be carried out by trained, qualified and authorized personnel (see chapter 1 SAFETY).

1. Place the instrument in a suitable location. (Risk evaluation in accordance with the decree law 626/94 of the Italian legislation). Please consider the instrument volume and make sure the following conditions are satisfied:

- The instrument is placed in the immediate vicinity of a power source.
- There is sufficient space around the instrument for normal use, maintenance and any peripheral equipment.

The instrument must be placed in a covered location in order to prevent any contact with atmospheric agents. For a safe working and maintenance of the instrument, a good lighting is required. A normal ambient light will prevent any risk on work connected to shade areas. For a normal functioning of the instrument, the allowed temperature goes from a minimum of +5°C to a maximum of +40°C with a maximum relative humidity of 80% for temperatures up to +31°C decreasing linearly to 50% at +40°C.

2. Before connecting the instrument to the power supply, please check the voltage data.

5.2 STARTING

1) ELECTRICAL CONNECTION

The instrument can be powered by using only one supply point, placed on the back side of the device.

2) STARTING TEST

Before proceeding to use this instrument in a regular and continuative way, verify its perfect functioning by carrying out a complete test cycle.

3) INSTRUMENT SWITCHING ON

- Insert the plug of the power supply cable in the network socket 220/240V-50/60Hz.
- Move the **IG** general switch into the “**I**” position.
- Rotate the timer knob to the desired time.
- The display turns on and the instrument starts heating up.



4) SET PHASE

- The SET POINT value appears by pressing button  or .

If no button is pressed within 4 seconds, the display starts showing the current temperature again.

- Press button  or  to set the desired temperature value.
- To save the set temperature value press button .
- The instrument starts showing the room current temperature again.
- Rotate the knob **③** to regulate the fan speed (only model STF-F).

In the central position, under the display, there is a three LED deviation indicator: one central green LED and two lateral red LEDs. There are 5 possible signals:

	→ The temperature is 1,5°C lower than the SET POINT.
	→ The temperature is within -1,0°C of the SET POINT.
	→ The temperature is ± 0,5°C from the SET POINT.
	→ The temperature is within +1,0°C of the SET POINT.
	→ The temperature is 1,5°C higher than the SET POINT.

5) INSTRUMENT SHUTDOWN PROCEDURE

- Move the **IG** general switch into the “**O**” position.
- Remove the containers from the instrument.
- **ATTENTION** → handling of hot objects.
- Disconnect the instrument from the power supply.

NOTE

If, for any reason, a treatment has to be interrupted before the end of a cycle, do **NOT** force the knob to move the **TIMER** into the “**O**” position.

ATTENTION

FALC INSTRUMENTS declines any liability for any damages to people, animals or objects, resulting from failure to comply with precautions contained in this manual or from failure to use the prescribed PPD (personal protection devices). (See chapter SAFETY on page 1).

6.0 MANUAL REGULATION FOR AIR VENTILATION AND ITS USE

52Lt, 80Lt and 120Lt models

This air vent-hole has a diameter of 45mm and is placed on the back side of the instrument. It is possible to regulate this opening using the knob placed on the control panel. Turning the knob clockwise the air vent-hole closes while turning the knob anticlockwise the air vent-hole opens. This opening can be connected to an intake system or used as a cable entrance.


240Lt and 400Lt models

This air vent-hole has a diameter of 45mm and is placed on the back side of the instrument. It is possible to regulate this opening using the lever placed on the control panel. Pulling the lever upward the air vent-hole closes while pulling the lever down the air vent-hole opens. This opening can be connected to an intake system or used as a cable entrance.

7.0 VENTILATION

To ensure a good working of the instrument and a good air circulation inside the chamber, arrange the materials leaving some space between one test piece and another and make sure the vent for the air circulation of the engine is not obstructed.

8.0 FORESEEN AND UNFORESEEN USES

☺ FORESEEN USES

- 👍 These instruments are engineered and constructed for incubation (ICT) or desiccation (STF) treatments.
- 👍 These instruments are engineered and constructed to treat substances that are normally used in laboratories.
- 👍 These instruments are engineered and constructed to be used in laboratories for professional, industrial and didactic purposes.
- 👍 These instruments can be used in the field of biochemistry, microbiology, haematology and diagnostics.


☹ UNFORESEEN USES

- 👎 All uses are unforeseen unless specifically indicated in the FORESEEN USES.
- 👎 Use of the instruments in environmental conditions different from those indicated.
- 👎 Use of the instruments in explosive gas atmospheres.
- 👎 Electromedical use.
- 👎 Water jet cleaning.
- 👎 Acids and solvents must not come into direct contact with the instruments.

ATTENTION

FALC INSTRUMENTS declines any liability for any damages to people, animals or objects, resulting from failure to comply with precautions contained in this manual.

9.0 ACCESSORIES

INOX steel perforated plate

Capacity	52Lt	80Lt	120Lt	240Lt	400Lt
Code	720.2091.55	720.2091.58	720.2091.62	720.2091.74	720.2091.80

Collapsible oven table

In steel plate, painted in oven with antacid epoxy dust

Capacity	52Lt - 80Lt	80Lt - 120Lt
LxPxH	800x650x850	900x750x800
Code	720.2091.92	720.2091.93

Oven trolley

In steel plate, painted in oven with antacid epoxy dust, with swivel wheels (Ø150) and brake

Capacity	240Lt – 400Lt
LxPxH	1000x800x500
Code	720.2091.98

Note: all measures are expressed in mm

Oven table

Aluminium profile table with shelf in anodized aluminium

Capacity	52Lt - 80Lt	80Lt - 120Lt
LxPxH	800x650x850	900x750x800
Code	720.2091.95	720.2091.96

10.0 MAINTENANCE

10.1 INSTRUMENT MAINTENANCE

Before embarking on any maintenance work, check that the instrument is switched off and is not connected to the power supply and make sure all its parts are cold.

Maintenance operations must be carried out exclusively by trained, qualified and authorized maintenance personnel, after studying and understanding all the instructions contained in this section.

10.1.1 ORDINARY MAINTENANCE

Ordinary maintenance operations can be carried out by the user. These are cleaning operations and periodic or precautionary inspection operations that ensure a safe use of the instrument.

10.1 PROGRAMMED MAINTENANCE

Programmed maintenance operations can be carried out exclusively by authorized personnel. These are periodic or precautionary inspections and interventions.

If any problem persists after the intervention of qualified and trained personnel, please contact the technical assistance.

10.2 CLEANING

ATTENTION

Switch off the instrument and disconnect it from the power supply before proceeding with any cleaning operation.

At the end of every day, proceed with an accurate general cleaning of the instrument. Please pay particular attention to the following areas:

- Instrument chamber.
- Instrument surrounding surface.
- Work area flooring.

10.2.1 CLEANING TOOLS AND PRODUCTS

Cleaning operations can be carried out with a normal detergent and water. Use isopropanol only if dirt is difficult to remove.

ATTENTION → Do not use any solvents as these may damage the varnish and the synthetic materials. In particular avoid the use of petrol, nitro-perchlorate and trichlorethene.

10.3 TRASPORT

If transportation is carried out on road vehicles with an open container, make sure the instrument, its outfit and the packing are completely covered with a plastic wrapping to protect them from atmospheric agents.

10.4 STORAGE

All the conditions contained in this section must be respected during the storage period of the instrument:

- Altitude up to 2000m.
- Allowed temperatures: from +5°C to +40°C.
- Maximum relative humidity of 80% for temperatures up to 31°C decreasing linearly to 50% at +40°C.
- Suitable protection from atmospheric agents.
- There is sufficient space to safely execute all the lifting and transport operations.
- The horizontal support base capacity is higher than the instrument weight.
- DO NOT PUT ANY OBJECTS ON THE INSTRUMENT.

11.0 PROBLEM SOLVING

In case of need, the instrument can only be repaired by the manufacturer (Falc Instruments) or by trained and authorized personnel.

Spare parts codes are indicated in the following chapter “MANUFACTURER SPARE PARTS”.

Type of failure	Problem solving
The instrument does not turn on	<p>Fuse control:</p> <ol style="list-style-type: none"> 1) Disconnect the instrument from the power supply. 2) Open the small drawers located below the tag. 3) Remove the fuses and, in case they are blown, replace them with new fuses with the same amperage and voltage rating as the previous ones. 4) Close the small drawers.
The instrument does not heat	<ol style="list-style-type: none"> 1) Verify that the thermoregulator parameters are set correctly (contact the technical assistance). <i>If the problem persists, replace the static relay:</i> 2) Disconnect the instrument from the power supply. 3) Disassemble the lateral panel using a screwdriver or another similar tool. 4) Unscrew the two self-tapping screws that hold the static relay. 5) Disconnect the relay from the rest of the incubator. 6) Connect the new static relay paying attention to the polarity. 7) Use the two self-tapping screws to connect the new static relay. 8) Assemble the lateral panel. <p style="text-align: center;">OTHERWISE</p> <p>The failure is caused by the malfunction of the instrument resistances. For further information, please contact the technical assistance.</p>
The instrument keeps heating	<ol style="list-style-type: none"> 1) Verify that the thermoregulator parameters are set correctly (contact the technical assistance). <i>If the problem persists, replace the static relay:</i> 2) Disconnect the instrument from the power supply. 3) Disassemble the lateral panel using a screwdriver or another similar tool. 4) Unscrew the two self-tapping screws that hold the static relay. 5) Disconnect the relay from the rest of the incubator. 6) Connect the new static relay paying attention to the polarity. 7) Use the two self-tapping screws to connect the new static relay. 8) Assemble the lateral panel.
The display shows ---- uuuu oooo	Verify that the probe is properly connected to the instrument and is working correctly.
The display shows ErEP	Press the “set” button.

All the failures described above can be attributed to an incorrect use of the instrument, to the technical life and to the daily wear and tear of the components.

If the authorized and trained personnel can not solve a particular problem, please contact the technical assistance.

12.0 MANUFACTURER SPARE PARTS

Spare part	Manufacturer code
Thermoregulator	TR0013
Static relay for 52Lt - 80Lt - 120Lt	RT0004
Static relay for 240Lt - 400Lt	RT0005
Probe	SN0006
Mechanical timer 0/12h + ∞	TR0020

STF-F	52Lt	80Lt	120Lt	240Lt	400Lt
FUSE	8 A	10 A	12 A	16 A	20 A
CODE	PF0024	PF0025	PF0026	PF0027	PF0028

STF-N	52Lt	80Lt	120Lt	240Lt	400Lt
FUSE	8 A	10 A	12 A	16 A	20 A
CODE	PF0024	PF0025	PF0026	PF0027	PF0028

ICT	52Lt	80Lt	120Lt	240Lt
FUSE	3.15 A	4 A	6.3 A	8 A
CODE	PF0020	PF0021	PF0023	PF0024

Note: each instrument is provided with two fuses.

	52Lt	80Lt	120Lt	240Lt	400Lt
RESISTANCE	600W	700W	800W	1300W	1700W
CODE	ER0035	ER0036	ER0037	ER0038	ER0039

Note: each instrument is provided with two resistances.

13.0 DECLARATION OF CALIBRATION

Calibration was performed using a thermometric chain calibrated in accordance with the SIT centre certified thermometric chain. (It is possible to ask the FALC to release a copy of the SIT certification).

Time periods to reach different temperatures are:

Model		°C		
STF-F		70	150	300
52Lt	min	10	22	50
80Lt	min	10	22	50
120Lt	min	12	25	60
240Lt	min	15	28	65
400Lt	min	20	30	70

Model		°C		
STF-N		70	150	300
52Lt	min	15	25	60
80Lt	min	15	25	65
120Lt	min	20	35	70
240Lt	min	40	50	80
400Lt	min	50	65	100


Model		°C	
ICT		70	150
52Lt	min	40	60
80Lt	min	45	90
120Lt	min	60	100
240Lt	min	70	130

Time periods to reach different temperatures are approximate and when they were calculated, the chamber was empty.

14.0 ELECTRICAL SCHEME


STF-F MODEL

NR	COMPONENTS	NR	COMPONENTS
A	FUSES	I	PT100 PROBE
B	FUSE DOOR	L	R.P.M. REGULATOR
C	GREEN SWITCH	M	GREEN WARNING LIGHT
D	0/12h E ∞ TIMER	N	RELAY
E	YELLOW WARNING LIGHT	O	ARMOURED RESISTANCE
F	RED WARNING LIGHT	P	ENGINE
G	SAFETY THERMOSTAT	Q	STATIC RELAY
H	TM-FALC THERMOREGULATOR		


STF-N MODEL

NR	COMPONENTI	NR	COMPONENTI
A	FUSES	I	PT100 PROBE
B	FUSE DOOR	L	R.P.M. REGULATOR
C	GREEN SWITCH	M	\
D	0/12h E ∞ TIMER	N	RELAY
E	YELLOW WARNING LIGHT	O	ARMOURED RESISTANCE
F	RED WARNING LIGHT	P	\
G	SAFETY THERMOSTAT	Q	STATIC RELAY
H	TM-FALC THERMOREGULATOR		

ICT MODEL

NR	COMPONENTI	NR	COMPONENTI
A	FUSES	I	PT100 PROBE
B	FUSE DOOR	L	R.P.M. REGULATOR
C	GREEN SWITCH	M	\
D	0/12h E ∞ TIMER	N	RELAY
E	YELLOW WARNING LIGHT	O	ARMOURED RESISTANCE
F	RED WARNING LIGHT	P	\
G	SAFETY THERMOSTAT	Q	STATIC RELAY
H	TM-FALC THERMOREGULATOR		

15.0 DEMOLITION

15.1 OPERATOR QUALIFICATION

Qualified operators are all those who have read and understood all contents included in chapter SAFETY on page 2 and in chapter DEMOLITION on page 14.

15.2 DEACTIVATION AND DEMOLITION

Once the instrument has reached the end of its technical and operative life, it must be deactivated. Deactivation and demolition must be carried out in a manner that allows raw materials that compose the instrument to be recycled for other purposes.

Deactivation and demolition of the instrument must be carried out in a safe way.

ATTENTION

FALC INSTRUMENTS declines any liability for any damages to people, animals or objects, resulting from recycling single parts of the instrument for functions or assembly situations which are different from the original ones.

15.3 DEACTIVATION AND DEMOLITION PROCEDURE

Deactivation and demolition operations must be performed exclusively by properly trained and equipped personnel.

-  Turn off the instrument
-  Disconnect the power supply
-  Disassemble the following parts:
 - Electrical parts.
 - Plastic components.
 - Resin support feet.
-  In case the instrument has to be moved, please refer to chapters INSTALLATION AND USE PROCEDURES on page 5 and TRANSPORT AND STORAGE on page 8.

ATTENTION

The instrument is made with non-biodegradable materials. The instrument is considered to be a RAEE (disused electrical or electronic device), according to the following directives:

2002/95/CE

2002/96/CE

2003/108/CE

Deposit the instrument at an authorized disposal site.

15.4 POST-DEMOLITION RESIDUAL RISKS

If all the instructions indicated in chapter DEACTIVATION AND DEMOLITION on page 14 are followed meticulously, there are no post-demolition residual risks.

Via G. Monzio Compagnoni, 2
24047 Treviglio (BG) – ITALY
Tel. +39-0363.304660
Fax +39-0363.304395
e-mail: falc@falcinstruments.it
web site: www.falcinstruments.it

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FALC INSTRUMENTS s.r.l.

SUBJECT: GUARANTY

For a period of 24 months from the date of delivery to the carrier. With the word "guaranty" we engage ourselves to substitute any and every faulty instrument, product components or part, due to faults or defect of our manufacturing process, not due to transport or misuse, as soon as we are informed and the item sent back to us for verification. The item should be sent back to us at your charge with a report explaining and describing the fault. When the problem has been verified and confirmed, the item shall be substituted or changed repaired and sent back to the customer free of charge. In a case the item has been tampered with, misused, or used not following the instructions, or in case wrong repair trials have been done on it, there shall be no guaranty whatsoever.

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The following unit complies with the safety requirements outlined by the EC DIRECTIVE concerning the guidelines for ELECTROMAGNETIC COMPATIBILITY (2004/108/EC) and LOW VOLTAGE REGULATIONS (2006/95/EC).

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The units are manufactured in compliance with the regulations

CEI EN 61000-6-2

CEI EN 61000-6-4

CEI EN 61010-1

CEI EN 60355-1

CEI EN 61326-1

Il responsabile legale / legal responsible

Scotti Francesco