ISO 9001 and ISO 13485 Certified  
Reg. No. 207267 MP21

# MD-150T

## Dual Halogen Illuminator Operator Manual

6018 Bowdendale Avenue  
Jacksonville, FL 32216  
Customer Service: 904 731-5869 FAX 904 733 4832  
Toll Free 800-684-6404

# TABLE OF CONTENTS

1. INTRODUCTION
2. WARNINGS
3. SPECIFICATIONS
4. OPERATING ELEMENTS, SYMBOLS AND FUNCTIONS
  - 4.1 FRONT PANEL
  - 4.2 REAR PANEL
5. INSTALLATION
  - 5.1 SETTING UP LIGHT SOURCE
  - 5.2 CONNECTING LIGHT CABLE
6. OPERATION
  - 6.1 POWERING UP LIGHT SOURCE
  - 6.2 LIGHT BRIGHTNESS CONTROL
7. CLEANING
8. DISINFECTION AND STERILIZATION
  - 8.1 DISINFECTING CONTROL UNIT
  - 8.2 STERILIZING LIGHT CABLE
  - 8.3 UL
9. MAINTENANCE, SERVICING & REPAIR
  - 9.1 LAMP REPLACEMENT
  - 9.2 LIMITED WARRANTY
10. END OF PRODUCT LIFE
11. TROUBLESHOOTING
12. AUTHORIZED REPRESENTATIVE

# 1. INTRODUCTION

**Congratulations on the purchase of your new MD-150T Dual Halogen illuminator!**

This user-friendly Halogen Illuminator is a high efficiency light source utilizing state-of-the-art superior illumination technology. It offers a variety of features such as:

- 3400°K color temperature
- Quiet operation
- Turret which adapts to various types of light guides
- Easy lamp replacement

In short, you have chosen the best and we would like to make sure you receive the optimal results with your new Halogen Illuminator by using it correctly.

This Operator Manual will help you to install the device and optimally integrate it with other components of your system. It will also instruct you how to operate the Halogen Illuminator, how to keep it clean, and sterilize it. It will give you maintenance and service guidelines as well as recommendations for best performance results.

## 2. WARNINGS

To prevent fire or electric shock, do not open or expose the light source unit to rain or moisture. Refer all servicing to qualified personnel only.

Not suitable for use in presence of flammable anesthetic mixture with air or with oxygen or nitrate oxide.

This product should be used only with type BF endoscopic instruments which have been certified according to IEC 60601-1 for medical equipment and IEC 60601-2-18 for endoscopic equipment.

This symbol indicates type BF equipment

### Caution

This product is not provided as sterile.

All devices connecting to the Illuminator must be classified as medical equipment. Additional information processing equipment connected to the Illuminator, a Medical System and the operator must determine that all equipment complies with the appropriate end-product standards (such as IEC 60950 or IEC 60065 and the Standard for Medical System, IEC 60601-1-1).

### Caution

Light source can cause permanent eye damage if viewed directly with unprotected eye. To reduce the chance of eye damage, set the intensity control **always** to the minimum level and insert the fiber optic cable into the unit before turning on the power.

### Caution EQUIPMENT CONNECTION

The fiberoptic cable must be a **NON-CONDUCTIVE CABLE**. It should not have conductive shielding or any other conductive connection between the patient and equipment. Such connection will impair safety of the equipment.

It must be rinsed free of soaking/disinfectant solution and dried before plugging into the light source receptacle. Ensure the optical surface is clean before engaging into the light source.

### 3. SPECIFICATIONS

Item	Specification
Lamp type	EJA Halogen
Power	150 Watt
Color temperature	3400°K
Lamp life	40 hours (typical)
Lamp replacement	Open lamp compartment door on bottom of unit
Light guide adapter	Turret type to fit your choice of four: Storz, ACMI, Wolf, Olympus
Brightness control	Mechanical dimming
Input voltage	120 AC / 230AC 50/60 Hz
Power consumption	345 watt max
Regulatory Approvals	UL2601, IEC 601.1 and CSA 601.1
Equipment Class	Class 1, BF-type
Mode of Operation	Continuous operation
Water Resistant	Not Protected Equipment, IPX0
Operating Environment Temperature Relative Humidity Air Pressure	+10° to +40° C (50° to 104° F) 30 to 85% 700 to 1060 hPa
Storage Environment Temperature Relative Humidity Air Pressure	-20° to +60° C (-4° to 140° F) 0 to 95% 700 to 1060 hPa
Dimensions	355(14)W x 130 (5)H x 250 (10)D mm
Weight	6.8kg (15Lbs)

## 4. OPERATING ELEMENTS, SYMBOLS AND FUNCTIONS

### 4.1 FRONT PANEL


Figure 1. Light Source Front Panel

No.	Name	Function
1	Turret	Accepts light source end fitting of light guide cable
2	Intensity Control	Manually controls the light output
3	Lamp switch	Turns unit and lamp ON and OFF

## 4.2 REAR PANEL


Figure 2. Light Source Rear Panel

No.	Name	Function
1	Circuit breaker	Over current protection. Type, resettable
2	Voltage selector switch	Set the voltage from 115 VAC to 230 VAC
3	AC main inlet	Accepts AC power cords
4	Serial number label	Serial number and power requirements
5	Venting grids	Allows adequate ventilation and cooling of light source
6	Fan	Fan 110VAC, airflow for cooling of unit
7	Connector	For potential equalization

## 5. INSTALLATION

### 5.1 SETTING UP THE LIGHT SOURCE

Place the light source on a stable surface (cart, counter, stand, etc.).

**NOTE**

*Avoid places where the light source may be splashed with liquid.*

*Absolutely DO NOT use in any environment with explosive or flammable gases.*

*DO NOT block the venting grids of the light source.*

Make sure the power switch is in OFF position.

Connect AC power cord to the power inlet located on the rear panel of the light source. Make sure power cord is completely seated into the power inlet. Secure power cord by tightening the screw with a Phillips Head screwdriver.

**CAUTION**

*Use only cords provided with the light source.*

**CAUTION**

*Connecting and disconnecting the power cord from the light source is the only means of connecting and disconnecting the unit from the mains supply.*

Plug the AC power cord into a wall outlet using the three (3) prong plug supplied with the unit.

**CAUTION**

*To prevent electric shock, connect power cords of peripheral equipment through medical isolation transformers.*

**NOTE**

*When using medical isolation transformer, be sure to check the transformer power ratings.  
Make sure that the power cord is connected to the main power with three prong plug (USA use UL2601-1 rated isolation transformers and/or power strips only).*

### 5.2 CONNECTING THE LIGHT CABLE

Connect the light cable to endoscope and then plug the light cable endfitting into the appropriate port of the turret on the front panel.

## 6. OPERATION

### 6.1 POWERING UP THE LIGHT SOURCE

To operate the light source: Attach the light source to the instrument or headlamp.

Turn on the power switch. The fan and lamp will start.

Each power switch controls the individual lamp as marked, i.e., Lamp 1, Lamp 2

### 6.2 LIGHT BRIGHTNESS CONTROL

The light brightness adjustment enables the user to obtain a brighter or darker illumination of the object of observation.

Adjust the light intensity by turning the intensity control knob.

Each lamp has its own control as marked

## 7. CLEANING

**NOTE** *Always disconnect the power cord before cleaning the system.*

The control unit can be cleaned with any cleaning agents used for external cleaning of electric equipment, according to instructions given by the manufacturer of the cleaning solution.

Do not allow excessive moisture or liquids to come in direct contact with the control unit.

Do not use cleaning agents that are not permitted for use with plastics, i.e., ammonia, acetone, salty acids (HCl), etc.

Do not allow cleaning agents or liquids to enter the control unit outlets.

## 8. DISINFECTION AND STERILIZATION

### 8.1 DISINFECTING THE UNIT

**NOTE** *Always disconnect the power cord before cleaning the system.*

Use any disinfectant agents which are commonly applied while disinfecting surfaces of electric medical equipment. Such disinfectant agents usually arrive in the form of sprays or damp cloths.

Follow the instructions given by the manufacturer of the disinfectant solution.

### 8.2 STERILIZING THE LIGHT CABLE

The light cable should be sterilized according to the sterilizer manufacturer's specification.

WITH RESPECT TO ELECTRIC SHOCK,  
FIRE MECHANICAL HAZARD ONLY  
IN ACCORDANCE WITH UL-2601-1  
98EF

WITH RESPECT TO ELECTRIC SHOCK, FIRE MECHANICAL  
AND OTHER SPECIFIED HAZARDS ONLY  
IN ACCORDANCE WITH CAN/CSA C22.2 NO. 601.1  
(CAN/CSA 601.2XX, IF APPLICABLE)  
MEDICAL EQUIPMENT CERTIFIED FOR CANADA  
98EF

## 9. MAINTENANCE, SERVICING AND REPAIR

Performance of preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious; thus, enhancing the instrument's reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer.

Defective items of equipment are to be serviced and repaired exclusively by persons authorized by the manufacturer. All repair work shall employ original manufacturer's parts only.

### 9.1 LAMP REPLACEMENT

**CAUTION** *Always disconnect power cord and turn main switch off before lamp replacement.*

#### Use lamp replacement EJA type, rated 150 WATT, 21 Volts

Turn light source off. Allow light source to cool completely. Disconnect power cord from outlet. Use a flathead screwdriver or coin, turn the large screw located on lamp compartment door counter-clockwise one-quarter turn. Manually pull the screw upward. This will open the door of the light source. You should now see the inside portion of the lamp compartment. To remove the lamp, push the lever toward the back of the unit; this will raise the lamp. Remove the lamp and replace with the new lamp. Close the lamp door and secure with the large screw. Plug the unit into the outlet and turn the power on.

## 9.2 LIMITED WARRANTY

Your fiberoptic product(s) carry a 3-year warranty from the date of purchase on workmanship and all defects of material, excluding replaceable sweatbands, headband, fiberoptic cable and lamps. Should your product prove to have such defects within three years of the purchase, SUNOPTICS TECHNOLOGIES® will repair or replace the product or component part without charge. Should your fiberoptic product(s) need servicing under this warranty, please contact SUNOPTICS TECHNOLOGIES® for return authorization documentation. You should carefully pack unit in a sturdy carton and ship it to the factory. Please include a note describing the defects, your name, telephone number and a return address. Warranty does not cover equipment subject to misuse, accidental damage, normal wear and tear or if transferred to a new owner without authorization from SUNOPTICS TECHNOLOGIES®. This warranty gives you specific legal rights and you may also have other rights that vary from state to state.

**POST WARRANTY REPAIRS:** You may return your product(s) for repair, shipping prepaid to the factory, with \$10.00 to cover handling. Your product will be inspected and an estimate of repair charges will be submitted to you for approval. Payment must be received before repairs are completed.

- In the U.S 800-684-6404 (toll free)
- FAX 904 733-4832

## 10. END OF PRODUCT LIFE

In accordance with the European Waste from Electrical and Electronic Equipment (WEEE) directive, we encourage our customers to recycle this product whenever possible. Disposal of this unit must be performed in accordance with the applicable local environmental regulations.

In the US a list of recyclers in your area can be found at: <http://www.eiae.org/>.

Please contact customer service to issue a return authorization to return product to manufacturer at the end of product life.

## 11. TROUBLESHOOTING

<u>Problem</u>	<u>Solution</u>
<b>The power is on, but lamp will not ignite.</b>	A. Check that the AC power cord is properly connected. B. Check the unit circuit breaker. If necessary, reset. C. Check the lamp connection. D. Replace the lamp (refer to 9.1)

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