Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

# COBRA Fusion Ablation System

Catalog Numbers 700-001, 700-002, 700-003, 700-004

Patent No. 5,651,780; 5,769,847; 6,106,522; 6,129,724; 6,387,092; 6,447,506; 6,471,699; 6,500,172; 7,115,122; 7,226,448; 7,335,196; 7,413,568; 7,542,807; 7,837,684; 7,957,820; Additional Patents Pending

## Instructions for Use

### TABLE OF CONTENTS

<u>English</u> .....	2
<u>Dutch</u> .....	10
<u>French</u> .....	20
<u>German</u> .....	30
<u>Italian</u> .....	40
<u>Spanish</u> .....	50

---

**Estech, Inc.**  
2603 Camino Ramon Suite 100  
San Ramon, CA 94583 USA  
Phone: +1 925-866-7111  
Fax: +1 925-866-7117  
E-mail: [info@estech.com](mailto:info@estech.com)

European Authorized Representative:  
**Medical Partners BV**  
Groenendijk 39  
4926 RE Lage Zwaluwe  
The Netherlands  
Phone: +31.168.48.37.32  
Fax: +31.168.48.47.90

## ENGLISH

# Instructions for Use

Carefully read all instructions prior to use. Observe all contraindications, warnings and precautions noted in these directions. Failure to do so may result in patient complications.

ESTECH relies on the physician to determine, assess, and communicate to each patient all foreseeable risks of the procedure.


### System Description:

The Cobra Fusion™ System is comprised of four main components: the COBRA Fusion 50 or 150 Ablation System, the Magnetic Retriever System, and the Pacing/Recording Adapter Cable.


**Figure 1. COBRA Fusion System: Left to Right: (1) COBRA Fusion 150 ablation catheter, (2) Magnetic Introducer, (3) Accessory Magnetic Retriever, (4) Alternate Accessory Ball-tipped Magnetic Retriever, (5) Accessory Pacing/Recording Adapter Cable**

The COBRA Fusion ablation catheter includes an integrated suction stabilizer designed to engage tissue under vacuum such that constant contact between catheter electrodes and tissue to be ablated is maintained throughout the procedure. The Introducer and Retriever System are designed to facilitate introduction and advancement of the catheter to the desired anatomical position. The accessory cable allows the catheter to be connected to a temporary pacing/recording device.

**Figure 2: Features of the COBRA Fusion Ablation Catheter**

The Cobra Fusion™ comprises a flexible distal ablation catheter, designed to conform to the specific anatomy of the tissue to be ablated. The distal section of the catheter allows for one to six 25 mm electrodes to be activated corresponding to the numeric indicators on the device. Any combination of active electrodes may be used. The user can optionally direct RF energy from the active RF electrodes toward the integrated indifferent electrode (Bipolar mode) or toward the indifferent electrode pads placed on the patient's back (Monopolar mode) by selecting the corresponding mode on the handle as shown below.

**Monopolar Mode:** Active electrode delivers RF energy to indifferent electrode pads on patient's back when this mode is selected by depressing this button.


**Bipolar Mode:** Active electrode delivers RF energy to the integrated indifferent electrode when this mode is selected by depressing this button.

The Cobra Fusion™ is connected to the Electrosurgical Unit (ESU) by a cable permanently attached to the handle of the Catheter. Directions for use of the ESU may be found in the ESU Operator's Manual.

The distal ablation catheter includes an integrated vacuum stabilizer, which is connected to a vacuum source using the included accessory tubing. Accessory tubing consists of a 3-way stopcock and tubing with a male luer connector on one end and vacuum source connector on opposite end. The stopcock/tubing attaches to the tubing that exits the catheter handle. A three foot segment of tubing for connecting a fluid canister to a vacuum source is also provided.

The Cobra Fusion™ comprises a magnet at its distal end which couples to the magnetic connector at the proximal end of the Introducer. This allows the user to temporarily attach the Introducer, guide the Catheter to the desired anatomical position, and then decouple the Introducer. The Cobra Fusion™ Introducer can be used as needed to advance the Catheter to the desired anatomical position. The Introducer has a curved shape that is straightened out with an inserted stainless steel stylet for initial use. After the straightened distal end of the Introducer is advanced to the desired location of the heart, the stylet is removed while the Introducer is advanced. The proximal magnetic connector of the Introducer is then attached to the distal magnetic connector of the Cobra Fusion™ Catheter. The distal end of the introducer can be retrieved by using either the ball-tipped retriever or the flexible retriever provided in the accessory retriever kit. Once coupled, pulling on the Introducer/Retriever will bring the catheter into alignment with the target tissue. Once in place, the Introducer is detached from the ablation catheter.

When the Cobra Fusion™ is connected to an auxiliary temporary external pacemaker it can be used to provide transient cardiac pacing, sensing, recording, and stimulation for the assessment of electrical isolation / conduction block of ablation lesions in the surgical treatment of arrhythmias. The COBRA Fusion must be disconnected from the ESU and connected to a temporary external pacemaker using the accessory cable provided.

## Indications

*The ESTECH Cobra Fusion Ablation System is intended to ablate cardiac tissue during cardiac surgery using radiofrequency (RF) energy when connected directly to the Estech Electrosurgical unit (ESU).*

*The ESTECH Cobra Fusion Pacing/Recording Adapter Cable may be used for temporary cardiac pacing, sensing, recording, and stimulation during the evaluation of cardiac arrhythmias during surgery when connected to a temporary external cardiac pacemaker or recording device.*

**Sterilization Method:** The sterilization method is EtO and provides sterility assurance level  $10^{-6}$  in compliance with the obligatory requirements of ISO 11135.

## Contraindications

- Local or systemic infection
- Do not apply suction over an artery, large vein or aneurismal tissue

## Warnings

- The device(s) should be used by physicians thoroughly trained in the techniques of invasive surgical procedures and in the specific approach to be used
- The operator should keep the set temperature and power limit as low as possible to achieve the desired end effect. This minimizes excessive thermal damage to tissue, potential for collateral damage to adjacent tissue not intended for ablation, and potential for char, thrombus, or steam pop formation or occurrence. **Recommended temperatures should not be exceeded.**
- Care should be taken to assure that the Catheter is not in contact with tissue other than that to be ablated to avoid inadvertent tissue damage
- Care should be taken when using the Catheter in the proximity of vascular or nerve tissue to avoid inadvertent tissue damage.
- Care should be taken to thermally isolate the tissue to be ablated when anatomically possible to avoid damage to unintended tissues or structures
- Following RF ablation, visual inspection of underlying tissues should be routinely performed to rule out the presence of inadvertent tissue damage
- Care should be taken to ensure that the Catheter is not in contact with other surgical instruments, staples or other objects while coagulating. Inadvertent contact with objects while coagulating could lead to conduction of RF energy or heat and unintentional ablation of tissues in contact with that object
- Care should be taken when positioning the Catheter to prevent perforation or other damage to adjacent tissue. Do not force the Introducer or Catheter during advancement if snagged
- When using the push button magnetic decoupling feature, ensure that the distal end of the accessory shafted device is directly facing the magnet to avoid applying the ejector rod to tissue.
- Turn vacuum control stopcock open to atmosphere and off to vacuum source prior to removal of suction attachment from tissue.
- Care should be taken when positioning Suction Stabilizer to prevent perforation or other damage to adjacent tissue during the application of a vacuum.
- Take care not to occlude vacuum lumen or puncture the Suction Stabilizer as this may result in loss of vacuum.
- Do not exceed recommended vacuum limit. Excessive vacuum may cause bruising and/or hematoma.

- Inadvertent application of vacuum or ablation over an artery may constrict or occlude the artery resulting in infarction.
- Read auxiliary device manual and observe instructions and warnings. This device can be used with a temporary pacemaker.
- Do not connect the pacing/recording adapter cable to supply mains (line voltage) operated equipment without verifying isolation of the connected equipment to EN60601-1-1. Supply mains operated equipment may introduce dangerous leakage currents into the heart.
- Avoid the possibility of unintentional contact between the device tip electrodes or pacing lead connectors and any conductive surface contact. The pacing lead connectors should only be connected to a temporary pacemaker or to a recorder designed for safely recording biopotentials.

## Precautions


- Contents supplied **STERILE** using an ethylene oxide (EtO) process. Carefully inspect the package prior to use for any breach of the sterile barrier or damage to the contents. If the sterile barrier integrity is compromised or the contents damaged, **DO NOT USE** and contact your ESTECH representative. Use of product with a compromised sterile barrier may lead to patient injury.
- For single use only. Do not reuse, reprocess or re-sterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another, patient injury, illness or death.
- Do not attempt to operate the system before thoroughly reading ESU Operator's Manual and Instructions for Use.
- Care should be taken when manipulating the device. While the distal portion of the catheter is designed to be flexible to conform to the anatomy of the area to be ablated, excessive or rough shaping of the suction stabilizer, including bending and axial torque, may damage internal components of the device. When using surgical instruments, only grasp at the distal or proximal ends, away from the electrodes. Do not grasp the electrodes or spine of the suction stabilizer with instruments.
- If using a TEE probe, care should be taken to withdraw the TEE probe prior to ablation to avoid compressing the esophagus against the left atrium during ablation
- Many variables, including patient anatomy, pathology, and surgical techniques, may influence procedural outcomes. Patient and procedure selection is a responsibility of the medical professional.
- Avoid attaching the device or ablating over extremely thin, fragile, or aneurysmal tissue.
- Pacemakers and implantable cardioverter/defibrillators can be adversely affected by RF signals and magnetic fields. Refer to the manufacturer's Directions for Use
- Dispersive Indifferent Patch (DIP) electrodes used with the ESTECH system should be applied carefully according to the manufacturer's directions. Poor or incomplete contact of the DIP electrodes may result in skin burns. The use of DIP electrodes which meet or exceed ANSI/AAMI HF-18 requirements is recommended
- The risk of igniting inflammable gases or other materials is inherent in the application of RF power. Precautions must be taken to restrict flammable materials from the area where tissue ablation is performed
- Electromagnetic interference (EMI) produced by the ESU during the delivery of RF power may adversely affect the performance of other equipment
- Bare electrode tips and lead wire connectors constitute a direct path to the heart. Small currents (approximately 10  $\mu$ amps) may be sufficient to cause fibrillation. Only trained personnel should handle the COBRA Fusion and COBRA Fusion Pacing/Recording Adapter Cable and connecting pins.

## Storage

- Store in a cool, dry place

**Directions for Use (Ablation)**

1. Inspect the System and all packaging materials carefully. Open the package using aseptic technique
2. Ensure 2 DIP electrodes are in good contact, according to the manufacturer’s directions
3. Connect the ESU cable into the black receptacle on the face of the ESU.
4. If required, attach 3-foot line from the vacuum canister to the vacuum regulator connection
5. Connect the blue connector on one end of the 12-foot vacuum line to the vacuum canister. Connect the male fitting of the three-way stopcock vacuum tube of the catheter handle. Turn the stopcock to the off position with the “Off” indicator pointed in the direction of the vacuum source.
6. If using the Introducer, insert the stylet into the grey curved tubing to straighten it and provide pushability for advancement.
7. Advance the Introducer to the desired anatomical location. By withdrawing the stylet the Introducer tubing will advance in a preformed curve. The tubing will curve opposite the axial dark stripe.
8. Remove the stylet and couple the proximal magnetic fitting of the Introducer tubing to the distal magnetic fitting of the Catheter. Ensure that the dark stripe of the Introducer tubing is oriented opposite the open side of the suction stabilizer.
9. If using the Magnetic Retriever, advance toward the Magnetic Introducer. Magnets at the distal end of the Introducer and Retriever will couple automatically when in close proximity.
10. By pulling on the Retriever, the Introducer and connected Catheter are brought into alignment with the target tissue. Once in place, detach the Introducer from the Catheter. The accessory Positioning devices can then be coupled to the distal end of the catheter to aid in catheter manipulation and placement.
11. Turn vacuum regulator on to -500 mmHg to affix device to tissue. Turn vacuum control stopcock on to vacuum source. Push stabilizer against tissue to complete seal. Allow vacuum to build up prior to activating RF energy. Maintain manual pressure as required to ensure catheter/tissue contact and seal integrity.
12. Select the mode of ablation to enable or disable the integrated indifferent electrode. When ablating in bipolar mode, indifferent electrode pads should be disconnected from the ESU receptacles.

13. Select up to 3 electrodes to be activated simultaneously on the ESU.
14. Activate selected electrodes by depressing the RF on/off switch on the ESU.
15. Radiofrequency energy may be discontinued by depressing the RF on/off switch on the ESU.
16. When the cycles are complete, turn vacuum control stopcock to off position prior to removing the Stabilizer. Remove the device in the opposite direction it was advanced.
17. Upon completion of procedure, disconnect device from ESU and discard after use.

**Ablation Time (seconds) COBRA Fusion at 50W/Electrode, -500 mmHg of vacuum**

Tissue Thickness	60°C		70°C	
				
5mm	60	30	60	30
7mm	120	120	120	90

*Times based on mode of energy delivery selected by surgeon*

### **Directions for Use (Pacing/Recording Mode)**

1. Inspect the System and all packaging materials carefully. Open the package using aseptic technique
2. Connect the COBRA Fusion to the temporary external pacemaker using the non-sterile accessory COBRA Fusion Pacing/Recording Adapter Cable. Connect the Cobra Fusion cable into the female receptacle of the adaptor cable. Connect the two pins to the temporary external pacemaker or recording device.
3. If required, attach 3-foot line from the vacuum canister to the vacuum regulator connection.
4. Connect the blue connector on one end of the 12-foot vacuum line to the vacuum canister. Connect the male fitting of the 12-foot vacuum line to the three-way stopcock. Turn the stopcock to the off position with the "Off" indicator pointed in the direction of the vacuum source.
5. If using the Introducer, insert the stylet into the curved tubing to straighten it and provide pushability for advancement.
6. Advance the Introducer to the desired anatomical location. By withdrawing the stylet the Introducer tubing will advance in a preformed curve. The tubing will curve opposite the axial dark stripe.
7. Remove the stylet and couple the proximal magnetic fitting of the Introducer tubing to the distal magnetic fitting of the Catheter. Ensure that the dark stripe of the Introducer tubing is oriented opposite the open side of the suction stabilizer.
8. If using the Magnetic Retriever, advance toward the Introducer. Magnets at the distal end of the first and second introducers will couple when in close proximity.
9. By pulling on the Retriever, the Introducer and connected Catheter are brought into alignment with the identified anatomical sites for pacing/recording. Once in place, the Introducer is detached from the Catheter. The accessory positioning devices can then be coupled to the distal end of the catheter to aid in catheter placement.
10. The pacing/recording leads are located within the distal end of the suction chamber of the Catheter. Turn vacuum regulator on to -500 mmHg to affix device to tissue. Turn vacuum control stopcock on to vacuum source. Push stabilizer against tissue to complete seal. Allow vacuum to build up prior to activating device. Maintain manual pressure as required to ensure catheter/tissue contact
11. Turn on the temporary external pacemaker. Refer to its IFU for appropriate settings and procedure.
12. The device will pace, sense, record, or stimulate when the auxiliary device is on, according to its intended use.
13. Upon completion of procedure, disconnect catheter from pacing/recording adapter cable. Discard catheter after use. The pacing adapter is reusable and should be cleaned and stored in accordance with hospital procedures.

### **Complications**

The following potential risks or discomforts may be associated with electrosurgical procedures. The frequency and severity of these events can vary, and may necessitate additional medical intervention, including surgery. Strict adherence to the forgoing instructions before use will help reduce the incidence of complications.

Allergic reaction, Arrhythmias, Cardiac or respiratory arrest, Cardiac valve damage, Chest pain, Damage to vessel intima or cardiac ultrastructures, Death, Embolus, Hematoma / ecchymosis, Hemorrhage, Infarction, Infection, Perforation, Pericardial effusion, Pericarditis, / pleuritis, Pseudoaneurysm, Pulmonary edema, Sinus or AV node injury, Stroke, Tamponade, Thrombosis, Vasovagal reaction.

### **Disposal**

After use, dispose of product and packaging in accordance with hospital, administrative and/or local, state, federal and international laws and regulations.

## Warranty and Limitations

ESTECH warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond ESTECH's control directly affect the instrument and the result obtained from its use. ESTECH's obligation under this warranty is limited to the repair or replacement of this instrument and ESTECH shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. ESTECH neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. ESTECH assumes no liability with respect to instruments reused, reprocessed or re-sterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such instrument.

### For Technical Information Contact:

**Endoscopic Technologies, Inc. dba Estech**  
 2603 Camino Ramon Suite 100  
 San Ramon, CA 94583 USA  
 Phone: +1 925-866-7111  
 Fax: +1 925-866-7117  
 E-mail: [info@estech.com](mailto:info@estech.com)

European Authorized Representative:  
**Medical Partners BV**  
 Groenendijk 39  
 4926 RE Lage Zwaluwe  
 The Netherlands  
 Phone: +31.168.48.37.32  
 Fax: +31.168.48.47.90

## Graphic Symbols for this Device Labeling

Catalogue Number 	Caution-See instructions for Use 	Manufacturer/Company Address 	Sterilized using Ethylene Oxide 
Batch Code 	Single Use Only 	Expiration Date 	Authorized Representative in the European Union 
Quantity 	Date of Manufacture 	Non-Pyrogenic 	Does not contain Natural Rubber Latex 
Contains no di (2-ethylhexyl) phthalate (DEHP) released from polyvinyl chloride (PVC) 	"Consult Instructions For Use" 	CE Mark with Notified Body - class II and III medical devices comply with EC Directive 93/42/EEC 	"Caution: U.S. Federal law restricts this device to sale by or on the order of a physician or other licensed practitioner." 

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