Medtronic

TRANSVENE™ 6937

TRANSVENE-SVC Unipolar, Endocardial, Superior Vena Cava Lead

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Medtronic, Transvene

Contents

| | | |
|----|--|---|
| 1 | Description | 3 |
| 2 | Indications | 3 |
| 3 | Contraindications | 3 |
| 4 | Warnings and precautions | 3 |
| 5 | Potential adverse events | 4 |
| 6 | Implant procedure | 4 |
| 7 | Specifications (nominal) | 7 |
| 8 | General warning | 8 |
| 9 | Medtronic disclaimer of warranty | 8 |
| 10 | Service | 8 |
| 11 | Explanation of symbols on package labeling | 9 |

1 Description

The Medtronic Transvene-SVC Model 6937 unipolar, endocardial, superior vena cava lead is designed for the delivery of cardioversion and defibrillation therapies. It is intended for single long term use in the superior vena cava or for single acute use in the central venous circulatory system. It is designed to be used with a compatible Medtronic implantable arrhythmia management device or with a compatible non-implantable Medtronic defibrillator implant support device along with at least one other Medtronic defibrillation electrode.

Silicone insulation surrounds the low resistance conductor coil and the coil electrode is a platinum alloy. The lead also has a DF-1¹ connector.

1.1 Package contents

Leads and accessories are supplied sterile. Each package contains the following items:

- 1 lead with 1 radiopaque anchoring sleeve², stylet, and stylet guide
- 1 vein lifter
- extra stylets
- product documentation

1.2 Accessory descriptions

Anchoring sleeve – An anchoring sleeve secures the lead to prevent it from moving and protects the lead insulation and conductors from damage caused by tight sutures.

Stylet – A stylet provides additional stiffness and controlled flexibility for maneuvering the lead into position. Each stylet knob is labeled with the stylet diameter and length.

Stylet guide – A stylet guide facilitates stylet insertion into the lead.

Vein lifter – A vein lifter facilitates lead insertion into a vein.

2 Indications

The lead has application for patients in which arrhythmia management devices are indicated or for acute internal cardioversion of atrial arrhythmias.

Prior to lead implant, it is strongly recommended that patients undergo a complete cardiac evaluation, which should include extensive electrophysiologic testing. Also, extensive electrophysiologic evaluation and testing of the safety and efficacy of the proposed pacing, cardioversion, or defibrillation therapies are recommended during and after the implantation of the system.

3 Contraindications

Please refer to an appropriate device manual.

4 Warnings and precautions

Line-powered and battery-powered equipment – An implanted lead forms a direct current path to the myocardium. During lead implant and testing, use only battery-powered equipment or line-powered equipment specifically designed for this purpose to protect against fibrillation that may be caused by alternating currents. Line-powered equipment used in the vicinity of the patient must be properly grounded. Lead connector pins must be insulated from any leakage currents that may arise from line-powered equipment.

Diathermy – People with metal implants such as pacemakers, implantable cardioverter defibrillators (ICDs), and accompanying leads should not receive diathermy treatment. The interaction between the implant and diathermy can cause tissue damage, fibrillation, or damage to the device components, which could result in serious injury, loss of therapy, and/or the need to reprogram or replace the device.

Electrophysiologic testing – Prior to lead implant, it is strongly recommended that patients undergo a complete cardiac evaluation, which should include electrophysiologic testing. Also, electrophysiologic evaluation and testing of the safety and efficacy of the proposed pacing, cardioversion, or defibrillation therapies are recommended during and after the implant of the system.

Single use – The lead is for single use only.

Inspecting the sterile package – Inspect the sterile package with care before opening it.

- Contact a Medtronic representative if the seal or package is damaged.
- Do not store this product above 40 °C.
- Do not use the product after its expiration date.

¹ DF-1 refers to the International Connector Standard (ISO 11318) whereby pulse generators and leads so designed are assured of a basic mechanical fit.

² Two radiopaque anchoring sleeves are provided with leads 85 cm or longer.

Sterilization – Medtronic has sterilized the package contents with ethylene oxide before shipment. This lead is for single use only and is not intended to be resterilized.

Handling the lead – Handle the lead with care at all times.

- Protect the lead from materials that shed particles such as lint and dust. Lead insulators attract these particles.
- Handle the lead with sterile surgical gloves that have been rinsed in sterile water or a comparable substance.
- Do not severely bend, kink, or stretch the lead.
- Do not use surgical instruments to grasp the lead or connector pins.
- Do not immerse the lead in mineral oil, silicone oil, or any other liquid, except blood, at the time of implant.

Handling the stylet – Handle the stylet with care at all times.

- Do not use a sharp object to impart a curve to the distal end of the stylet.
- Do not use excessive force or surgical instruments when inserting the stylet into the lead.
- Avoid overbending or kinking the stylet.
- Use a new stylet when blood or other fluids accumulate on the stylet. Accumulated blood or other fluids may damage the lead or cause difficulty in passing the stylet into the lead.

Necessary hospital equipment – Keep external defibrillation equipment nearby for immediate use during acute lead system testing, the implant procedure, or whenever arrhythmias are possible or intentionally induced during post-implant testing.

Magnetic resonance imaging (MRI) – Do not use magnetic resonance imaging (MRI) on patients who have this device implanted. MRI can induce currents on implanted leads, potentially causing tissue damage and the induction of tachyarrhythmias.

Concurrent devices – Output pulses, especially from unipolar devices, may adversely affect device sensing capabilities. If a patient requires a separate stimulation device, either permanent or temporary, allow enough space between the leads of the separate systems to avoid interference in the sensing capabilities of the devices. Previously implanted pulse generators and implantable cardioverter defibrillators should generally be explanted.

Chronic repositioning or removal of a lead – Proceed with extreme caution if a lead must be removed or repositioned. Chronic repositioning or removal of a lead may not be possible because of the fibrotic tissue. If a lead is abandoned, it should be capped to avoid transmitting electrical signals. Return all removed leads, unused leads, or lead sections to Medtronic for analysis.

See additional cautions in the section titled Section 6, “Implant procedure”, page 4.

Second anchoring sleeve – Leads that are 85 cm or longer include 2 anchoring sleeves. Use both anchoring sleeves to assure adequate fixation.

Connector compatibility – Although the lead conforms to the International Connector Standard DF-1, do not attempt to use the lead with any device other than a commercially available

implantable defibrillator system with which it has been tested and demonstrated to be safe and effective. The potential adverse consequences of using such a combination may include, but are not limited to, undersensing cardiac activity and failure to deliver necessary therapy.

5 Potential adverse events

The potential adverse events related to the use of transvenous leads include, but are not limited to, the following patient-related conditions that can occur when the lead is being inserted or repositioned:

- air embolism
- cardiac or venous perforation
- cardiac tamponade
- constrictive pericarditis
- endocarditis
- fibrillation or other arrhythmias
- heart wall or vein wall rupture
- hemothorax
- infection
- pneumothorax
- tissue necrosis

Other potential adverse events related to the lead include, but are not limited to, the following conditions:

- insulation failure
- lead conductor or electrode fracture
- lead dislodgement
- poor connection to the device, which may lead to oversensing, undersensing, or a loss of therapy

6 Implant procedure

Proper surgical procedures and sterile techniques are the responsibility of the medical professional. The following procedures are provided for information only. Some implant techniques vary according to physician preference and the patient’s anatomy or physical condition. Each physician must apply the information in these instructions according to professional medical training and experience.

6.1 Opening the package

Use the following steps to open the sterile package and inspect the lead:

1. Within the sterile field, open the sterile package and remove the lead and accessories.
2. Inspect the lead. Leads that are shorter than 85 cm should have 1 anchoring sleeve on the lead body. Leads that are 85 cm or longer should have 2 anchoring sleeves on the lead body.


6.2 Using the stylet guide and stylets

The lead is packaged with the stylet guide attached to the connector pin and a stylet already inserted.

A stylet provides additional stiffness and controlled flexibility for maneuvering the lead into position. Stylets vary in stiffness to accommodate a physician's preference for lead and stylet flexibility.

If the stylet guide has been removed, replace it by gently pushing it as far as possible onto the connector pin (Figure 1). Insert a stylet through the stylet guide and into the lead body. If a slight curve is needed for the stylet, refer to the procedure in Section 6.4.

Figure 1. Stylet guide attachment.

Caution: To avoid damage to the lead or body tissue, do not use excessive force or surgical instruments to insert a stylet into the lead. To avoid lead tip distortion, the stylet should always remain fully inserted into the lead during lead introduction and while advancing the lead, especially through tortuous veins, that may cause the stylet to “back out” of the lead. When handling a stylet, avoid overbending, kinking, or blood contact. If blood is allowed to accumulate on a stylet, passage of the stylet into the lead may be difficult.

6.3 Selecting an insertion site

Caution: Use care when handling the lead during insertion.

- Do not severely bend, kink, or stretch the lead.
- Do not use surgical instruments to grasp the lead or connector pins.


Insert the lead using the following techniques:

1. Select a site for lead insertion. The lead may be inserted by venotomy through several different venous routes, including the right or left cephalic vein or the external or internal jugular vein. Use the cephalic vein whenever possible to avoid lead damage in the first rib or clavicular (thoracic inlet) space.

Cautions:

- Certain anatomical abnormalities, such as thoracic outlet syndrome, may also precipitate pinching and subsequent fracture of the lead.
- When using a subclavian approach, avoid techniques that may damage the lead.
- Place the insertion site as far lateral as possible to avoid clamping the lead body between the clavicle and the first rib (Figure 2).

Figure 2.


1 Suggested entry site

- Do not force the lead if significant resistance is encountered during lead passage.
 - Do not use techniques such as adjusting the patient's posture to facilitate lead passage. If resistance is encountered, it is recommended that an alternate venous entry site be used.
2. Insert the tapered end of a vein lifter into the incised vein and gently push the lead tip underneath and into the vein (Figure 3).

Note: A percutaneous lead introducer (PLI) kit may be used to facilitate insertion. Refer to the technical manual packaged with an appropriate percutaneous lead introducer for further instructions.

Figure 3.

³ If additional stylets are needed, Medtronic recommends using the same type of Medtronic stylet that is packaged with the lead.

6.4 Positioning the lead

Caution: Use care when handling the lead during positioning.


- Do not severely bend, kink, or stretch the lead.
- Do not use surgical instruments to grasp the lead or connector pins.

To place the electrode in the superior vena cava, advance the tip such that the distal portion of the electrode is approximately 3 to 5 cm proximal to the superior vena cava right atrial junction.

For added control in maneuvering the lead tip, the distal end of the lead may be slightly curved by inserting a gently curved stylet as shown in Figure 4

Caution: Do not use a sharp object to impart a curve to the distal end of the stylet. Imparting a curve to the stylet can be accomplished with a smooth-surface, sterile instrument (Figure 4).

Figure 4.

The electrode for the Transvene-SVC must be as far as possible from other defibrillation electrodes. Also, the cardiac mass between the electrodes should be as large as possible. Insufficient separation between the electrodes may result in reduced defibrillation efficacy.

Accurate positioning of the electrode is essential for defibrillation efficacy. Use fluoroscopy (P-A and lateral position) to ensure that the distal end of the electrode is properly positioned in the superior vena cava.

6.5 Taking electrical measurements and defibrillation efficacy measurements

Caution: Prior to taking electrical or defibrillation efficacy measurements, move objects made from conductive materials, such as guide wires, away from all electrodes. Metal objects, such as guide wires, can short a lead and an active implantable device, causing electrical current to bypass the heart and possibly damage the implantable device and lead.

In order to demonstrate reliable defibrillation efficacy, it is recommended that fibrillation be initiated and successfully terminated. Please refer to the appropriate device manual for proper instructions.


6.6 Anchoring the lead

Use the triple groove anchoring sleeve to secure the lead from moving and to protect the lead insulation and conductor coil from damage caused by tight ligatures (Figure 5 and Figure 6).

Anchor the lead with nonabsorbable sutures.

Caution: Tabs on anchoring sleeves are provided to minimize the possibility of the sleeve entering the vein. Do not remove the tabs (Figure 5). If using a large diameter percutaneous lead introducer (PLI) sheath, extreme care should be taken to prevent passage of the anchoring sleeve into the PLI lumen and/or the venous system.

Figure 5. Triple groove anchoring sleeve with tabs.

1 Anchoring sleeve tab

With a triple groove anchoring sleeve, all three of the grooves should be used with the following procedure (Figure 6).

The triple groove anchoring sleeve is situated at the connector end of the lead. Partially insert the anchoring sleeve into the vein.


Use the most distal suture groove to secure the anchoring sleeve to the vein.

Use the middle groove to secure the anchoring sleeve to the fascia and lead. First, create a base by looping a suture through the fascia underneath the middle groove and tying a knot. Continue by firmly wrapping the suture around the middle groove and tying a second knot.

Use the third and most proximal groove to secure the anchoring sleeve to the lead body.

Caution: Do not use the anchoring sleeve tabs for suturing.

Figure 6. Triple groove anchoring sleeve secured to the lead and fascia using three grooves.

Tie the sutures securely but gently to prevent damage to the triple groove anchoring sleeve.

Caution: Do not secure the ligatures so tightly that they damage the vein or lead. Do not tie a ligature directly to the lead body (Figure 7). During anchoring, take care to avoid dislodging the lead tip.

Figure 7. Do not secure the sutures too tightly and do not tie a suture to the lead body.

6.7 Connecting the lead

Use the following steps to connect the lead to an implantable device:


1. Carefully remove the stylet and stylet guide. When removing the stylet and stylet guide, grip the lead firmly just below the connector pin to prevent dislodgement.
2. Insert the lead connectors into the lead connector block. Consult the product documentation packaged with the implantable device for instructions on proper lead connections.

6.8 Placing the device and lead into the pocket

Caution: Use care when placing the device and leads into the pocket.

- Ensure that the leads do not leave the device at an acute angle.
- Do not grip the lead or device with surgical instruments.
- Do not coil the lead. Coiling the lead can twist the lead body and may result in lead dislodgement (Figure 8).


Figure 8.

Use the following steps to place the device and leads into the pocket:

1. To prevent undesirable twisting of the lead body, rotate the device to loosely wrap the excess lead length (Figure 9).

Figure 9.

2. Insert the device and leads into the pocket.
3. Before closing the pocket, verify sensing, pacing, cardioversion, and defibrillation efficacy.

6.9 Post-implant evaluation

After implant, monitor the patient's electrocardiogram until the patient is discharged. If a lead dislodges, it usually occurs during the immediate postoperative period.

Recommendations for verifying proper lead positioning include x-rays and pacing/sensing thresholds taken at pre-hospital discharge, 3 months after implant, and every 6 months thereafter.

In the event of a patient death, explant all implanted leads and devices and return them to Medtronic with a completed Product Information Report form. Call the appropriate phone number on the back cover if there are any questions on product handling procedures.

7 Specifications (nominal)

7.1 Detailed device description


Table 1. Specifications (nominal)

| Parameter | | Model 6937 |
|--------------------------------|-------------|---|
| Type | | Unipolar |
| Position | | Superior vena cava |
| Fixation | | None. Tip is placed in the superior vena cava |
| Length | | 20–110 cm |
| Connector | Unipolar: | DF-1 |
| Materials | Conductors: | Multifilar MP35N Composite |
| | Insulation: | Silicone |
| | Electrode: | Platinum alloy |
| Electrical shadow area | | 568 mm ² |
| Exposed electrode surface area | | 160 mm ² |
| Unipolar resistance | | 1.5 Ω (58 cm) |
| Diameters | Lead body: | 1.8 mm |

Table 1. Specifications (nominal) (continued)

| Parameter | Model 6937 |
|------------------------------------|---------------------|
| | Tip: 2.3 mm |
| Lead introducer (recommended size) | 2.3 mm (7.0 French) |

Figure 10. Model 6937 lead components

- 1 Lead tip diameter: 2.3 mm²
- 2 Electrode length: 8 cm
- 3 Electrode surface area: 160 mm²; electrical shadow area: 568 mm²
- 4 Insulator: Silicone
- 5 Lead length: 20–110 cm
- 6 Lead body diameter: 1.8 mm
- 7 Lead body diameter: 3.2 mm

8 General warning

Medtronic implantable leads are implanted in the extremely hostile environment of the human body. Leads are necessarily very small in diameter and must still be very flexible, which

unavoidably reduces their potential performance or longevity. Leads may fail to function for a variety of causes, including, but not limited to: medical complications, body rejection phenomena, allergic reaction, fibrotic tissue, or failure of leads by breakage or by breach of their insulation covering. In addition, despite the exercise of all due care in design, component selection, manufacture, and testing prior to sale, leads may be easily damaged before, during, or after insertion by improper handling or other intervening acts. Consequently, no representation or warranty is made that failure or cessation of function of leads will not occur or that the body will not react adversely to the implantation of leads or that medical complications (including perforation of the heart) will not follow the implantation of leads or that the lead will, in all cases, restore adequate cardiac function.

9 Medtronic disclaimer of warranty

Although the Medtronic Transvene-SVC Model 6937 Implantable Defibrillation Lead, hereafter referred to as "Product" has been carefully designed, manufactured and tested prior to sale, the Product may fail to perform its intended function satisfactorily for a variety of reasons. The warnings contained in the Product labeling provide more detailed information and are considered an integral part of this DISCLAIMER OF WARRANTY. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the Product. Medtronic shall not be liable for any incidental or consequential damages caused by any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort or otherwise.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

10 Service

Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed on the back cover.

11 Explanation of symbols on package labeling

Table 2. Explanation of symbols on package labeling

Refer to the package labels to see which symbols apply to this product.






| | |
|--|--|
|  | Conformité Européenne (European Conformity). This symbol means that the device fully complies with European Directive AIMD 90/385/EEC. |
|  | Date of manufacture |
|  | Manufacturer |
|  | Authorized representative in the European community |
|  | Use by |
|  | Serial number |
|  | Reorder number |
|  | Dangerous voltage |
|  | Open here |
|  | Attention: see accompanying documents |
|  | Lot number |
|  | Ethylene oxide sterilized |
|  | Do not reuse |
|  | Lead length |
|  | A-V distance |
|  | Storage temperature |
|  | Maximum storage temperature |
|  | Implantation date |

Table 2. Explanation of symbols on package labeling (continued)

| | |
|---|--|
|  | For your patient file/for your registration card |
|  | With guide wire |
|  | Without guide wire |
|  | Lead introducer |
|  | Insertion diameter |
|  | J-shaped |

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© Medtronic, Inc. 2010
The following are trademarks or registered trademarks of Medtronic in the United States and possibly in other countries:
Medtronic B.V.
Medtronic Australasia Pty Ltd