01

# EDI

**LAMPADA DENTALE**  
DENTAL LIGHT

MANUALE D'USO  
INSTRUCTION MANUAL  
MODE D'EMPLOI  
GEBRAUCHSANLEITUNG  
MANUAL DE USO

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Electrical connection diagrams (\*)

---

(\*) Illustrations and connection diagrams at end of booklet.

---

---

## ***SYMBOLS USED IN THIS MANUAL***

### ***DANGER***

The paragraphs marked with this symbol contain instructions which must be observed carefully to avoid damage to the device, the operator and the patient.

### ***WARNINGS***

These instructions advise that attention should be paid to avoid situations which could damage the device.

### ***PROHIBITION***

This icon shows what to avoid doing to avoid damage to the device.

### ***SUGGESTIONS***

This icon signals information permitting use of the device in the most effective manner.

Disposal container. Dispose the device observing the standards about the differentiated collection of electrical material.

Device classe II.

### ***ATTENTION***

Consult the enclosed documentation.

## ***SAFETY RULES***

The dental operating lamp, used to light the operating site, should be used exclusively by medical staff in a dental practice.

- The device should be cleaned before use (see the section “Cleaning the device”).
- The package of the lamp is suitable to protect the same against external agents.
- Do not use the lamp if any of its parts are damaged.
- Do not use the lamp without the front screen: risk of projection of parts or irritation.
- The device should be installed exclusively by qualified staff.
- The dental operating lamp should be installed on a specific control and power device such as dental units, or equipped with electrical system compliant with standard IEC 60364-1 and the “national installation standards for electrical installations in medical locations”.
- The device should be installed with an omnipolar circuit breaker compliant with standard IEC 61058-1.
- Should the lamp be installed on a dental unit, the installer should ensure its compliance with standard IEC 60601-1.
- Make sure the supply voltage specified on the data plate matches the mains voltage.
- Do not clean the plastic parts of Edi lamp with disinfectants containing:  
AMMONIUM HYDROXIDE - SODIUM HYDROXIDE - METHYLENE CHLORIDE - METHYL ALCOHOL.

Failure to comply with this provision may lead to:

- risk of damage and subsequent rupture of plastic parts.

Should you have any doubt, contact The Customer Care Faro.

- Do not carry out any maintenance intervention on the lamp when the power is on; disconnect the power cord from the mains before proceeding.
- Do not insert in the slots the head of the lamp, objects or equipment that could come into contact with parts under voltage.
- Do not fix nor point the light beam directly to patient’s eyes especially in the case of patients at increased risk of eye injury (e.g. children with eye conditions). In this case, always use suitable protective equipment and adopt due caution: Faro suggests you use safety goggles BLUE-BAN or BABY BLUE-BAN (for children).

- Do not spray detergents-disinfectants directly on the head.
- Do not obstruct the openings next to the fan.
- Do not use the product with accessories that are not specified in this manual.
- Data plate affixed to the rear arm of the device.
- The device is not suitable for installation in explosive atmosphere or rich in oxygen.
- The articulated arm and head joints enable correct positioning of the light beam. Do not overload the arms and the joints with shocks on the limit switch.
- The dentist must use disposable safety equipment on lamp handles or guarantee their sterilisation.
- Unsuitable performance of the device does not compromise patient's safety.
- Divide the materials according to the type (ferrous, rubber, plastic etc.). Dismantling and disposal of materials should be carried out in accordance with the standards in force in your country, contacting special waste disposal centres, if required, that are acknowledged and authorized.
- Do not leave small parts of the equipment unattended or at reach of exposed persons (children) as they can be hazardous.
- Nominal arm load 5 Kg.
- Do not carry out maintenance interventions or replacements other than those specified in the manual. Any intervention not specified in this manual may compromise the safety level provided by the device.
- Carry out exclusively the operations given in this manual; centre for any other issue, contact the technical support.
- FARO does not allow any change to the product unless specifically authorized in writing; failure to comply with this rule will nullify the compliance with safety standards and the warranty.
- The device of the original package can be transported or stored for 15 weeks if the following conditions are met:
  - ⇒ **Room temperature from  $-20^{\circ}\text{C}$  to  $+70^{\circ}\text{C}$**
  - ⇒ **Relative humidity from 10% to 90%**
  - ⇒ **Atmospheric pressure from 500 to 1060 mBar**

# SAFETY RULES

## REQUIREMENTS FOR ELECTROMAGNETIC COMPATIBILITY

This section contains specific information concerning conformity of the product with the standard IEC 60601-1-2: 2007. The MAIA dental lamp is an electrical medical device which requires special precautions as regards: electromagnetic compatibility, and which must be installed and put into service in accordance with the electromagnetic compatibility information provided. Mobile and portable RF communication equipment (mobile phones, radio transceivers, etc.) may influence the medical system. The use of accessories, transducers and cables sold by the manufacturer of the equipment and the system as replacement parts may result in an increase in emissions or a decrease in the immunity of the equipment or systems.

Manufacturer's guidelines and statement – Electromagnetic emissions		
The lamp EDI is designed to function in the electromagnetic environment specified below. The client or user must ensure its use in the said environment.		
Emission tests	Compliance	Electromagnetic environment - Guidelines
RF Emission CISPR15	Compliant	The lamp EDI uses RF energy only for its internal function. Therefore its RF emissions are very low and most likely do not cause any interference in neighbouring electronic devices.
RF Emission CISPR15	Compliant	The lamp EDI is fit for use in all buildings, including domestic ones and those directly connected to the public low voltage supply network that feeds buildings for domestic use.
Harmonic emission	Class C	
Voltage fluctuations/flicker emission	Compliant	

Recommended distances between portable and mobile radiocommunication devices and the dental unit			
The lamp EDI is designed to function in an electromagnetic environment in which irradiating RF disturbances are under control. The client or operator of the unit can contribute toward preventing electromagnetic interferences by ensuring a minimum distance between mobile and portable RF communication devices (transmitters) and the dental unit, as recommended below, depending on the maximum output power of the radiocommunication devices.			
Maximum nominal output power of the transmitter W	Distance for transmitter frequencies (m)		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters whose maximum nominal power is not listed above, the recommended distance d in metres (m) can be calculated by using the applicable equation for the transmitter frequency, with P as maximum nominal output of the transmitter in Watts (W), depending on the manufacturer. Notes: The highest frequency interval is applied at 80 MHz and 800 Mhz. These guidelines might not apply to all situations. Electromagnetic propagation is influenced by absorption and reflection of structures, objects and persons.			

## ELECTROMAGNETIC IMMUNITY

<b>Manufacturer's guidelines and statement – Electromagnetic immunity</b>		
The lamp EDI is designed to function in the electromagnetic environment specified below. The client or user must ensure its use in the said environment.		
<b>Immunity test</b>	<b>Compliance</b>	<b>Electromagnetic environment - Guidelines</b>
Electrostatic discharge (ESD) IEC/EN61000-4-2	± 6kV contact ± 8kV air	The floor must be in wood, concrete or ceramic. If the floor is covered with synthetic material, relative humidity should be at least 30%.
Electrical fast transient/burst IEC/EN61000-4-4	± 2kV power supply ± 1kV for input/output lines	The quality of supply network voltage should be typical of commercial or hospital environments.
Surge IEC/EN61000-4-5	± 1kV differential mode ± 2kV common mode	The quality of supply network voltage should be typical of commercial or hospital environments.
Voltage dips, short interruption and voltage variation IEC/EN61000-4-11	< 5% Ut for 0,5 cycle 40% Ut for 05 cycle 70% Ut for 25 cycle <5% Ut for 5 sec.	The quality of supply network voltage should be typical of commercial or hospital environments. If the user of the lamp EDI requires continuous use even without a supply network, use an uninterruptible power supply.
Power frequency magnetic field IEC/EN61000-4-8	3A/m	Level of magnetic field at the network frequency typical of commercial or hospital environments.
Conducted immunity IEC/EN61000-4-6	3Vrms 150kHz to 80MHz (for non life-supporting equipment)	Portable and mobile RF communication devices should not be used near any part of the dental unit, including cables, unless they comply with recommended distances calculated with the applicable equation for transmitter frequency. Recommended distances: $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ from 80 Mhz to 800 MHz $d = 2.3\sqrt{P}$ from 800 Mhz to 2.5 GHz P is the maximum nominal power issued by the transmitter in Watts (W) depending on the manufacturer of the transmitter, and d is the recommended distance in metres (m). The intensity of the fixed RF transmitter field, as established in an electromagnetic investigation of site a, could be less than the compliance level of each frequency interval. There can be interference near devices marked with the following symbol: 
Conducted immunity IEC/EN61000-4-6	3Vrms 80MHz to 2.5GHz (for non life-supporting equipment)	
<p>Note: Ut is the power-line voltage</p> <p>Note 1: The highest frequency interval is applied at 80 MHz and 800MHz.</p> <p>Note 2: These guidelines might not apply to all situations. Electromagnetic propagation is influenced by absorption and reflection of structures, objects and persons.</p> <p>a) ISN bands (industrial, scientific and medical) between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz and 40.66 MHz to 40.70 MHz.</p> <p>b) Compliance levels in ISN bands between 150 kHz and 80 MHz and 80 MHz to 2.5 GHz present a decreasing probability of portable transmission devices causing interference if inadvertently taken to the patient area. Therefore, an additional 10/3 factor has been incorporated into the formula used to calculate the distance between transmitters.</p> <p>c) Field intensities for fixed transmitters such as base stations for radiotelephones (mobiles and cordless) and cellular mobile radios on land, CB user equipment, AM and FM transmitters and TV transmitters cannot be theoretically estimated with precision. To establish an electromagnetic environment caused by fixed RF transmitters, an electromagnetic investigation of the site should be considered. If field intensity measured at the site of use of the dental unit exceeds the aforementioned applicable compliance level, normal function of the lamp should be monitored. If any abnormal performance is noticed, additional provisions such as a different orientation or position of the lamp might be necessary.</p> <p>d) The field intensity in an interval of frequencies from 150 kHz to 80 MHz should be less than 3 V/m.</p>		

## ***FEATURES***

The “edi” Lamp is supplied in the following models.

### ***STANDARD***

- *S/TS Lamp (with or without Transformer) with fixed handle and mechanical switch.*
- *S/TS Lamp (with or without Transformer) ceiling model*
- *Head and accessories in the FARO standard colors.*

### ***VARIANTS***

- *With detachable handle for sterilization*
- *With “infrared sensor” switch*

### ***APPLICATIONS***

- *Can be mounted on dental unit, wall, floor or ceiling*
- **Fan** incorporated in head to hold temperature within range ensuring compliance with safety requirements of Std 60601-1.
- **Light source** consisting of a halogen bulb whose light is reflected by a parabolic reflector.
- **The parabolic reflecting** surface is glass with dichroic treatment and has a geometry assuring a regular, uniform light spot at all levels of intensity. The reflector is designed with approximately 600 facets to distribute light uniformly in the operating area without shadows or obscuration by the operator.
- **Adjusting light intensity** by means of an adjustable screen. This system keeps color temperature of 5000 K constant.
- **Front Screen** protects operator and patient from any ultraviolet radiation.
- **Detachable Handles** allow sterilization.
- **Infrared sensor** permits to switch On/Off the lamp without direct contact to avoid the risk of cross contamination
- **Maintenance** is facilitated by application of new technology allowing for safety, human engineering and hygiene requirements.

## INSTALLATION

### Mod. S/TS “edi” Dental Lamp

Make sure all the components listed below are included in the package received.

- S/TS Dental Lamp (in the requested version)
- Wrench for removing Front Screen
- Instruction Manual

**Clean the device before use** (see section Cleaning the device).

**This device should be installed by specialized technicians.**

**Power should always be cut off during installation.**

### LAMP ASSEMBLY, DENTAL CHAIR VERSION

- Install the lamp by inserting the terminal lamp pin into the specific hole on the dental chair.

	SAFE WORKING LOAD	MINIMUM BREAKING LOAD
Long arm 855 mm	2.92 Kg	23.5 Kg.
Long arm 550 mm	2.56 Kg	20.5 Kg

### LAMP WITH TRANSFORMER


- It should be powered directly to the mains corresponding to the characteristics written on the data plate or the technical specifications given in the manual.
- Make sure the master switch comply with standard IEC/EN 61058-1.

### LAMP WITHOUT TRANSFORMER

- It should be powered by alternate current at low voltage (17 V AC) using a safety transformer compliant with EN 60601-1.
- The final medical installation should be declared compliant with IEC/EN 60601-1 by the installer.
- Make sure that upstream the lamp is installed a fuse of T 8A L 250 V. Warning: the fuse is not supplied with the lamp.

### MOUNTING THE LAMP ON CEILING-WALL-FLOOR

- The applications are not supplied with the lamp.

- To mount every single application please refer to the relative instruction manual

## ***OPERATING INSTRUCTIONS***

The Edi lamp is designed to illuminate the oral cavity during intervention procedures in the dentist's surgery. The optical parameters are obtained by positioning the lamp at a distance of 700 mm from the subject. At different distances the light spot may lose its optimum parameters in terms of size, homogeneity and illuminance. It is therefore advisable to work as close as possible to this distance.

### **Mod. S/TS "edi" Dental Lamp**

Use of the lamp is limited to the following operations:

- **Turning On**

⇒ Mechanical. Use Switch "A"

(fig.1)

⇒ By Sensor. Place hand within not more than 3 cm distance of Sensor "B"(fig 2)

- **ON / OFF**  
**Adjusting Light Intensity**

⇒ Use Wheel "C" to adjust light intensity (fig 3).

### **Assembling forked lower handle (on request)**

To insert forked lower handle proceed as follows:

- Remove front screen (fig.4 pag.20)
- Remove plugs by acting on points "M" (fig.10)
- Insert forked lower handle by pressing on points "L" (fig.9).  
Release it when fitted in the holes.

**PLUGS OR FORKED LOWER HANDLE SHOULD BE ALWAYS  
FITTED IN THE HEAD WHILE LAMP IS WORKING**

**FOR FIXED HANDLES USE DISPOSAL COVER**

## ***MAINTENANCE***

### **Replacing Front Screen (fig 4)**

To remove Front Screen from Lamp proceed as follows:

- Turn screw on side of Screen counterclockwise using Tool supplied.
- Pull out Screen.

To reinstall Front Screen perform above operations in reverse.

### **Replacing Bulb (fig 5/6)**

**(mod.Halogen Lamp GY 6.35 17 V 95 W)**

To remove Bulb proceed as follows:

- Remove Front Screen (fig 4).
- Pull out Screen while pressing sides (fig 5).
- Withdraw Bulb “D” by pulling it upward (fig 6).
- Fit new Bulb. It must be as specified in the Data Sheet or on Data Plate.

**The new Bulb should not be touched with the hands; handle it with the protective shell. If it is accidentally touched with the fingers wipe away the finger prints with cotton soaked in alcohol.**

### **Replacing Fuse in Transformer model (fig 7).**

The Lamp with Transformer has two Fuses of the same power. To replace them proceed as follows.

- Make sure Power Cord is disconnected.
- Unscrew Plugs “E” located on Transformer case.
- Withdraw Fuses, check for burn out, replace if necessary.

**New Fuses must have same rating as shown on data plate and in Data Sheet.**

## **REGULAR CHECKS**

### **• *Yearly***

- Check mechanical soundness of joints and plastic elements.
- Legibility of label.

### **• *Every two years***

- Insulation test in accordance with EN 62353.
- Test for leakage in the casing.
- Check maximum and minimum illuminance.

## **CLEANING**

The handles are not supplied sterile. Sterilize them before use. The handles can be sterilized with saturated steam under pressure according to the following sterilisation cycles with saturated steam under pressure for up to maximum 200 sterilizations:

Sterilization temperature C°	Minimum holding time min	Press kPa
121	15	207
134	3	308

To remove sterilizable Handles proceed as follows.

### **Side Handle open (fig 8)**

Turn knob “G” counterclockwise and pull until completely unhooked. To fit, reverse the operation.

### **Forked lower Handle (fig 9)**

Press points “L”, then pull out. To fit, reverse the operation.

### **Cleaning Front Screen and Parabolic Reflector**

To remove Screen proceed as explained in **Replacing Front Screen** on page 20 fig.4. Cleaning should be done with a soft cloth moistened with a glass-cleaning liquid,

 **Abrasive or detergent materials such as alcohol, trichloroethylene, gasoline (petrol), turpentine or the like should not be used on any part of the Lamp.**

 **For cleaning plastic parts in the Lamp EDI DO NOT USE detergents-disinfectants containing: AMMONIUM HYDROXIDE - SODIUM HYDROXIDE - METHYLENE CLORIDE - METHYL ALCOHOL. Failure to comply with this indication might cause: • deterioration risk of plastic parts with consequent breakage. In case of doubt, please contact Customer Care Faro.**

---

## ***TROUBLESHOOTING***

---

**The Lamp does not light.**

- Make sure Electric Power Cord is correctly plugged in.
- Check condition of Fuses.
- Check condition of Bulb.
- If the problem is not solved, call our Service Department.

**Light intensity is considerably dimmed**

- Clean Parabolic Reflector or Front Screen with soft cloth moistened in a glass-cleaning detergent.
- If light remains dim, call our Service Department.

## DATA SHEETS

	Version with transformer	Version without transformer
Supply voltage	<b>230 V</b>	<b>17 V</b>
Frequency	<b>50 / 60 Hz</b>	<b>50 / 60 Hz</b>
Input power	<b>105 W</b>	<b>105 W</b>
Input current	<b>0,5 A</b>	<b>6,2 A</b>
Halogen lamp	<b>17 V 95 W</b> (mod Halogen Lamp GY 6.35)	<b>17 V 95 W</b> (mod Halogen Lamp GY 6.35)
Transformer	<b>230 V-17 V</b>	//
Spotlight dimensions (1)	<b>60 mm x 180 mm</b>	<b>60 mm x 180 mm</b>
Maximum illuminance (1)	<b>25000 lux</b>	<b>25000 lux</b>
Colour temperature (1)	<b>5000 K</b>	<b>5000 K</b>
Maximum temperature accessible surfaces	<b>&lt; 71°C under normal operating conditions</b> (plastic materials)	<b>&lt; 71°C under normal operating conditions</b> (plastic materials)
Fuses	<b>T 1,6A L 250V</b>	<b>T 8A L 250V (*)</b>
Fan noise level	<b>&lt; 40 dB (A)</b>	<b>&lt; 40 dB (A)</b>
Protection against electrical hazard	<b>Equipment class II</b>	<b>Equipment class II</b>
Type of operation	<b>Continuous</b>	<b>Continuous</b>
Environmental conditions for use	<ul style="list-style-type: none"> <li>● Temperature from 10°C to 40°C</li> <li>● Relative humidity from 30 to 75%</li> <li>● Max altitude: 2000 m</li> <li>● Atmospheric pressure from 800 to 1060 mba</li> </ul>	<ul style="list-style-type: none"> <li>● Temperature from 10°C to 40°C</li> <li>● Relative humidity from 30 to 75%</li> <li>● Max altitude: 2000 m</li> <li>● Atmospheric pressure from 800 to 1060 mbar</li> </ul>

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**GB****GUARANTEE CERTIFICATE**

Faro offers the final customer a **12 month** guarantee starting from the date of purchase.

Repairs under guarantee must be performed at FARO; expenses and transport risks are at the risk of the purchaser.

**Repair under guarantee is considered valid only when:**

- **all sections of the certificate have been filled in and sent in advance to FARO by Fax (039.6010540).**

The guarantee covers faults due to the bad quality of the material or manufacturing defects; in the case of valid claims, the guarantee covers free repair or replacement. **Claims for damages and/or interest are excluded.** The guarantee is not considered valid, at the sole discretion of FARO, if the fault is due to tampering, damage, incorrect use, improper maintenance and normal wear and tear.

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The logo for FARO, featuring the word "FARO" in a stylized, red, serif font. The letters are bold and have a slightly irregular, hand-drawn appearance. The "F" and "A" are particularly prominent.

DAL 1948: ESPERIENZA  
E RINNOVAMENTO

DAL 1948: ESPERIENZA  
E RINNOVAMENTO

 **FARO S.p.A.**

via Faro, 15 - 20876 Ornago (MB) - Italy  
Tel. +39 039.68781 - Fax +39 039.6010540  
**www.faro.it - comm.italia@faro.it - export@faro.it**

**FARO FRANCE**

Za Tgv Coriolis - 71210 Monchanin - France  
Tel. +33 385.779680 - Fax +33 385.779688  
**www.farofrance.com - farofrance@farofrance.com**

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Gewerbepark Heideckhof Heideckstr. 179  
D-47805 Krefeld - Germany  
Tel. +49 2151.936921 - Fax +49 2151.936933  
**www.faro.it - info@faro-dental.de**

Azienda  
Certificata

CERT. 9124.FAR2

CERT. 9120.FAR1

FARO SpA si riserva il diritto di modificare, senza preavviso, le caratteristiche indicate nel presente manuale.  
FARO SpA reserves the right to change the specifications of this equipment without notice.  
FARO SpA se reserve le droit de modifier, sans préavis, les caractéristiques dans ce manuel.  
FARO SpA behält sich recht vor, jederzeit stillschweigend technische oder bauliche Änderung vorzunehmen.  
FARO SpA se reserva el derecho de modificar sin aviso previo la características incluidas en el presente manual de uso.