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## **REC**

### ***Piastra riscaldante in ceramica***

**F20700080 (230V - 50÷60Hz)**

**F20710080 (115V - 60Hz)**

## **Prefazione:**

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**F20710080 (115V - 60Hz)**

## Preface:

Velp Scientifica thanks you for choosing the Ceramic plate **REC**.  
The **REC** model is designed for general use and for all applications which require a set heating plate temperature.

## Safety instructions:

Before using the unit, please read the operating manual supplied with the apparatus carefully.

Warning ! Hot surface

Do not dispose of this equipment as urban waste

In order to prevent possible risk of electric shocks, fire and personal injury when the unit is being used, basic safety measures must always be taken, including:

- The heating plate can reach a temperature of up to 540°C. Care must be taken during the heating phase and also when the appliance is switched off as it takes time to cool down.
  - The Heating Ceramic Plate must not be used:
    - in hazard areas
    - in bain-marie
    - to heat combustible liquids that have a low combustion temperature
  - The vessel must be smaller in diameter than the ceramic plate to prevent any risk of tipping over.
  - If the heating function is used the vessel must be made of a suitable material to withstand the foreseen temperature – the use of Pyrex vessels is recommended.
  - The following precautions must be taken before using the instrument:
    - position the instrument on a level surface.
    - position the instrument at least at 30 cm from the wall
    - an operator must be present all the times the instrument is in use
  - The liquid must not come into contact with the electric power cable or with the electrical parts inside the instrument.
  - Check that the power supply corresponds to the rating plate on the rear of the unit and distant from the heating plate.
- 
- Only correctly earthed three-terminal power cables must be used.
  - Do not use the unit if it is not working correctly. In case of malfunctioning, contact your nearest assistance centre.
  - Personal protective equipment must be compatible with the possible risks posed by the material being processed and the glass parts.
  - Follow the cleaning instructions described in this manual

**This unit must only be used for laboratory applications**

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

This unit has been designed and produced in compliance with the following standards.

Safety requirements for electrical apparatus for:

Measurement and control and for laboratory use	CEI EN 61010-1
Electrical equipment for laboratory use	UL 3101-1
General requirement – Canadian electrical code	CAN/CSA-C22.2

***N.B.:*** *The manufacturer is committed to constantly improving the quality of the products and reserves the right to modify the characteristics without prior notice.*

# CONTENTS

<b>1</b>	<b>INTRODUCTION</b>	<b>1</b>
1.1	Parts included	2
1.2	Instrument description	2
1.3	Knobs and functions	3
1.4	Information regarding construction materials	3
<b>2</b>	<b>ASSEMBLY AND INSTALLATION</b>	<b>4</b>
2.1	Electrical connection	4
2.2	Start-up	4
<b>3</b>	<b>OPERATING CONTROLS</b>	<b>4</b>
<b>4</b>	<b>END OF WORK CYCLE OPERATIONS</b>	<b>5</b>
<b>5</b>	<b>MAINTENANCE</b>	<b>5</b>
5.1	Cleaning	5
5.2	Disposal	5
5.3	Error messages	6
<b>6</b>	<b>ACCESSORIES AVAILABLE ON REQUEST</b>	<b>7</b>
<b>7</b>	<b>SPARE PARTS</b>	<b>7</b>
<b>8</b>	<b>TECHNICAL FEATURES</b>	<b>8</b>
<b>9</b>	<b>WIRING DIAGRAM</b>	<b>9</b>
<b>10</b>	<b>WARRANTY</b>	<b>10</b>

---

# 1. Introduction

---


The **REC** heating ceramic plate, is used in the laboratory for general use and for all the applications that require a precise regulation of the heating plate temperature.

The ceramic heating plate **(1)** makes it possible to obtain higher temperatures than the traditional aluminium heating plate and a faster sample heating. Ceramic is an inert and very hard material and it is resistant to almost any type of chemical or mechanical aggression; its surface characteristics remain unaltered over time.

The temperature of heating plate **(2)**, indicated on the digital display **(3)**, are set using the relative knobs on the frontal panel.

The “Hot Plate” light **(4)** turns on when the temperature of the heating plate is higher than 50°C.

**NOTE:** If the heating function is used the vessel must be made of a suitable material to withstand the foreseen temperature.

## 1.1 Parts included

Check the unit is complete after unpacking.  
The table below shows the parts included.

Parts included		Qty
1	Heating ceramic plate	1
2	Power cable	1
3	Operating manual	1

## 1.2 Instrument description

The structure of the unit, made out of die-cast aluminium, combined with its innovative design and low profile, allows the operator to work in extreme comfort and safety under all conditions.

The new design combines ergonomics with a high Electrical Protection rating (IP42) according to the rules CEI-EN 60529.

The control panel that contains the electronic parts, is completely separate from the heating plate in order to increase the safety rating during the use. The inclination of the knobs has been carefully studied in order to facilitate the use by the operator.


The temperature of the heating plate can be set within a range of 5 – 540°C and with a 5°C resolution, using the knob on the front panel.

The temperature set and the control of the function run are shown on the display.

A temperature probe in direct contact with the heating plate takes continuous temperature readings whilst an electronic closed-loop feedback system ensures that the set temperature is maintained.

## 1.3 Knobs and Functions

The knob is positioned on the front panel, it is easily accessible and enough far from the heating plate in order to ensure maximum operator safety as well as to safeguard the electronic components inside the unit.

<p><b>TEMPERATURE</b></p> 	<p>Turn the heating knob/switch on the front panel in order to start the heating function. It is possible to select heating plate temperatures from 5 to 540°C with a 5°C resolution. The set temperature is indicated on the display.</p>
<p><b>“HOT PLATE”LIGHT</b></p> 	<p>The HOT PLATE light is always on when the temperature of the heating plate exceeds 50°C.</p> <p>When the heating function is turned Off (Heating knob on OFF) and the heating plate temperature is over 50°C, the HOT PLATE light flashes and OFF is indicated on the display until the temperature of the heating plate falls below 50°C.</p> <p><b><i>N.B.</i></b> This warning is not active if the instrument is not connected to the electric power supply.</p>

## 1.4 Information regarding construction materials

Housing	Die-cast aluminium
Paint	PET
Front control panel	PET
Feet	EPDM
Heating plate	Ceramic

If any liquid being processed falls onto the above materials, clean them immediately to avoid chemical corrosion.

---

## 2. Assembly and installation

---

### 2.1 Electrical connection

After unpacking the instrument, place the unit on a laboratory bench. Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the cable supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and easy to reach.

### 2.2 Start-up

After the connecting the instrument to the power supply, place a flask containing the sample.

---

## 3. Operating controls

---

To start the heating function turn the heating knob “Heating °C” on the front panel. Temperatures from 5 to 540°C can be selected at 5° intervals. The temperature setting is shown on the display.

The HOT PLATE light remains lit when the temperature of the heating plate exceeds 50°C.

When heating is turned Off (Heating knob OFF) and the heating plate temperature is over 50°C, the HOT PLATE indicator flashes and OFF is shown on the display until the temperature of the heating plate falls to below 50°C.

***N.B.*** *This warning is not active if the instrument is not connected to the power supply.*

---

---

## **4. End of work cycle operations**

---

---

At the end of the work cycle turn off the heating knobs.

---

---

## **5. Maintenance**

---

---

No routine or extraordinary maintenance is necessary apart from periodically cleaning of the unit as described in this manual (chap.5.1).

In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors.

Should the fuse need changing, disconnect the unit from the power supply and open the cover on the back of the unit to access the two fuses 5A (230V) or 8A (115V).

### **5.1 Cleaning**

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable and non-aggressive detergent to clean the unit.

### **5.2 Disposal**

The final disposal of the unit or of its components must be carried out in compliance with the directives and laws in force in the country concerned.

### 5.3 Error messages

If the display shows an error message please contact Velp Scientifica's technical service department.

When the display shows an error message, the heating function stops automatically.

**NOTE:** To remove the error message disconnect the instrument from the power supply.

Message	Problem
AL1	Thermocouple is out of range
AL2	Excessive heating time

## 6. Accessories available on request

Description:	Code No
Hemispheric bowl for 250ml flasks	A00001071
Hemispheric bowl for 500ml flasks	A00001072
Hemispheric bowl for 1000ml flasks	A00001073
Support rod	A00001069

A00001056	A00001060	A00001061
		
A00001071	A00001072	A00001073
		
A00001069		


## 7. Spare parts

Description:	Code No.
Knob	10002097
Foot	10000239

## 8. Technical features

Power supply	V/Hz	230/50-60 115/60
Total electric power	W	770
Dimensions	mm (wxhxd)	205x96x335
	Inch (wxhxd)	8x3.7x13.2
Weight	Kg	3.15
	Lbs	7
Temperature range	°C	5-540
	°F	41-1004
Temperature control		Electronic probe
Temperature selection		Display reading 5°C intervals
Max load-bearing capacity of plate	Kg	25
	lbs	55
Ceramic plate dimensions	mm	200x200
Construction material		Die cast aluminium
Environmental temp. range	°C	+5...+40
	°F	+41...+104
Storage temp. range	°C	-10...+60
	°F	+14...+140
Max. humidity	%	Max 80
Operating mode		Continuous
Pollution degree CEI EN61010-1		2
Installation category CEI EN61010-1		2
Fuse	A	2x5A (2x8A – 115V/60Hz)

## 9. Wiring diagram

- 1) Main board
- 2) Thermocouple
- 3) Temperature potentiometer
- 4) Display board
- 5) Heating element
- 6) Fuses rit. 5x20mm n°2x5A (n°2x8A-115V)

---

## 10 Warranty

---

The unit is guaranteed against production defects for 25 months from the date on the delivery note and the serial number on the unit.

In accordance with this guarantee VELP SCIENTIFICA undertakes to repair any units resulting as faulty due to the quality of the materials used or poor workmanship.

Costs for returning the unit will be met by the Customer.

Units rendered faulty due to inexpert handling/use or carelessness by the Customer's personnel will not be replaced or repaired under warranty.

Repairs not covered by the warranty will be charged to the Customer.

### **Exclusions:**

The guarantee will be considered null and void for faults resulting from:

- inexperience and carelessness of the Customer
- repairs, maintenance or replacement of parts carried out by personnel or Companies not approved by the manufacturer

## ***Declaration of conformity***

*We the manufacturer*            *VELP SCIENTIFICA s.r.l.*  
*Address*                            *Via Stazione, 16*  
    *USMATE (Milan)*  
    *Italy*

*under our responsibility declare that the product is manufactured in conformity with the following standards:*

*EN 61010-1 (2001)*  
*EN 61326-1 (1997) + A1 (1998) + A2 (2001) + A3 (2003)*  
*2002/95/CE (RoHS)*

*and satisfies the essential requirements of the following directives:*

*Machines directive 98/37/CEE*  
*Low voltage directive 2006/95/EC*  
*Electromagnetic compatibility directive 2004/108/EC*  
*plus modifications and that the documents listed in annex V are available at Velp's offices as foreseen by the machine directive.*

## **Thank you for having chosen a VELP product !**

Since 1983 Velp has offered to professionals in the sector a range of sophisticated and reliable equipment with excellent operating capacity and backed by high levels of know-how at competitive prices.

Velp works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification. Instruments are built according to the International norms IEC 1010-1 and to the rules of the CE mark.

Let us present our product Lines:

### **Environment**

C.O.D.

B.O.D. – B.O.D. Sensor

REFRIGERATED THERMOSTATS

INCUBATORS

JAR TEST

METALS IN SLUDGES AND SEDIMENTS

NEPHELOMETRIC TURBIDITY

NITROGEN IN DIFFERENT FORMS

PERFORMANCE OF ANAEROBIC DIGESTERS

PHENOLS

LEACHABILITY TESTS

### **Stirring**

MAGNETIC STIRRERS

HEATING MAGNETIC STIRRERS

MULTIPLACE MAGNETIC STIRRERS

LIGHTED MAGNETIC STIRRERS

HEATING PLATES

VORTEX

STIRRERS

### **Food and feed**

DIGESTERS

STEAM DISTILLING UNITS

SCRUBBER

CATALYST TABLETS

NITROGEN-FREE WEIGHING BOATS

EXTRACTORS FOR RAW FIBER

EXTRACTORS FOR DIETARY FIBER

EXTRACTORS FOR SELECTIVE EXTRACTION OF SOLUBLE PRODUCTS BY SOLVENTS (FOR FAT)

### **Pumps**

PERISTALTIC PUMPS

RECIRCULATING WATER VACUUM PUMP

**Fisher Scientific**

Part of Thermo Fisher Scientific

**Belgique / België**

tél 056 260 260

fax 056 260 270

[bioblock.belgium@thermofisher.com](mailto:bioblock.belgium@thermofisher.com)

[www.bioblock.be](http://www.bioblock.be)

**España**

tfno 91 515 92 34

fax 91 515 92 35

[bioblock.ventas@thermofisher.com](mailto:bioblock.ventas@thermofisher.com)

[www.es.fishersci.com](http://www.es.fishersci.com)

**France**

tél 03 88 67 53 23

fax 03 88 67 85 11

[bioblock.vente@thermofisher.com](mailto:bioblock.vente@thermofisher.com)

[www.bioblock.com](http://www.bioblock.com)

**Italia**

tel 02 953 28 258

fax 02 953 27 374

[fisher.vendite@thermofisher.com](mailto:fisher.vendite@thermofisher.com)

[www.it.fishersci.com](http://www.it.fishersci.com)

**Netherlands**

tel 020 487 70 00

fax 020 487 70 70

[nl.info@thermofisher.com](mailto:nl.info@thermofisher.com)

[www.fishersci.nl](http://www.fishersci.nl)

**Suisse / Switzerland**

tel 056 618 41 11

fax 056 618 41 41

[info.ch@thermofisher.com](mailto:info.ch@thermofisher.com)

[www.ch.fishersci.com](http://www.ch.fishersci.com)