KARTELL spa

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Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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This unit must be used for laboratory applications only.

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use **IEC/EN 61010-1**

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Electrical equipment for laboratory use **UL 61010-1**

General requirement - Canadian electrical code **CAN/CSA-C22.2 No.61010-1**

The manufacturer reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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Do not use with explosive and dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it.

It is also responsibility of the user to use safety substances for cleaning or decontaminating, which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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Contents / Indice

1.	INTRODUCTION.....	3
2.	ASSEMBLY AND INSTALLATION.....	3
	2.1 ELECTRICAL CONNECTIONS.....	3
	2.2 START-UP.....	3
	2.3 INFORMATION ABOUT CONSTRUCTION MATERIALS.....	3
3.	OPERATING CONTROLS.....	3
4.	MAINTENANCE.....	4
	4.1 CLEANING.....	4
5.	TECHNICAL DATA.....	4
6.	ACCESSORIES.....	4
7.	SPARE PARTS.....	4
8.	WIRING DIAGRAM.....	5
9.	WARRANTY.....	5
10.	DECLARATION OF CONFORMITY CE	5

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Mixing by spinning can be carried out quickly and safely by placing the tube containing the sample on the mixing cup head. The TK3S can be used in automatic (TOUCH) or continuous mode (CONT.): in both cases the mixing speed can be adjusted using the speed selector knob (1).

When the instrument is switched on, it is always set in TOUCH mode. Mixing starts when the test-tube tube is placed on the mixing cup head (4). For continuous mode press the mode select button (2), mixing starts automatically at the selected speed, no pressure is required on the cup head.

The led (5) indicates that the instrument is running in TOUCH mode.

The press-on fitting means that the mixing cup head can easily be replaced by any of the various accessories available for test-tubes of different sizes.


Fig. 1

2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- TK3S Shaker complete with mixing cup head
- Adapter plug
- Power supply 100-240V/12V no plug
- Instruction manual

2.1 Electrical connections

After having unpacked the instrument, place the unit on the laboratory bench.

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the transformer supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

2.2 Start-up

Check that the mixing cup head is correctly positioned. Also check that the speed regulation knob is turned to "0 rpm". Use the on-off switch on the front of the unit to turn the unit on (position "ON"). The start-up and the control of the speed are effected with the knob. When switched on the unit starts up in "SENSOR" mode with the indicator lights on.

2.3 Information about construction materials

Housing	Zinc alloy/Polymer	Feet	MQ/MVQ
Paint	PET	Platform	PP
Frontal control panel	PET	Foamed part for platforms	PE
Cup head	SEBS		

3. Operating controls

Connect the unit to mains and turn it on using the on-off switch. Select the operating mode using the mode select button on the front of the unit.

Fig. 2

When the **Touch** mode is selected, mixing starts automatically when the tube is lightly pressed down on the mixing cup head as highlighted in Figure 2.

In **Continuous** mode the mixing cup head runs continuously at the set speed without pressure on the mixing cup head. However, it is important to grip the tube firmly during mixing.

CONTINUOUS-TOUCH

Pushing the mode select button, the respective modes are activated:

- **Touch** (indicator light on) allows automatic operation mode even with different accessories.
- **Continuous** (indicator light off). The instrument runs in continuous mode with various accessories.

SPEED REGULATION KNOB

The knob on the front of the unit can be used for quick precision regulation of mixing speed from 0 to 3000 rpm.

ON-OFF SWITCH

The on-off switch turns the unit on and off. If the switch is in the "OFF" position the unit is off; if the switch is in the "ON" position the unit is on.

Always turn the unit off after use.

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors.

To replace the cup head, pull it upwards and off. The other push-fit accessories can then be inserted.

4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

5. Technical data

Power supply input	AC 100-240V ($\pm 10\%$); 50-60 Hz ; 0.5A
Power supply output	DC 12V ; 1.25A
Dimensions (WxHxD)	150x130x165 mm (5.9x5.1x6.5 in)
Weight	2.7 Kg (6.2 lb)
Speed range	0-3000 rpm
Spin diameter	4.5 mm (0.177 in)
Max. weight on the mixing cup head/ Accessories	0.5 kg (1.1 lb)
Permitted operation	Continuous
Operation mode	Continuous or Touch
Level of electrical protection IEC/EN60529+A1	IP 42
Noise level	<< 80 dBa
Admitted temperature	Storage: -10...+60 °C (+14...+140 F°) Working: 5°...+40 °C (+41...+104 °F)
Admitted humidity	Max. 80%
Power	15 W
Overvoltage category	II
Pollution degree IEC/EN61010-1	2
Max altitude	2000 m

6. Accessories

- A00000012 Foam stand for 19 microvials 1.5 ml - Eppendorf®
- A00000013 Customizable soft foam top
- A00000014 Foam stand 5 test tubes Ø16mm
- A00000015 Foam stand for microtiter
- A00000016 Small rubber supporting plate Ø 50mm
- A00000017 Big foam supporting plate Ø 94 mm
- A00000019 Foam stand 4 test tubes Ø29mm


Do not exceed a speed of 800 rpm when any one of the above accessories is in use. High speeds may cause the sample to spill.

	A00000012	A00000013	A00000014
A00000015	A00000016	A00000017	A00000019

7. Spare Parts

10005195	Mixing cup head TX, ZX, RX	10003082	UE plug for power supply 40001099
10002097	Knob 24D black	10003083	USA plug for power supply 40001099
10002564	Foot 28Dx5H antivibration	10003084	UK plug for power supply 40001099
40001099	Switching 100-240V/12V without plug	10003085	AU plug for power supply 40001099

8. Wiring Diagram

- | | |
|-----------------------------|--|
| 1. Transformer | 5. Electric motor |
| 2. Socket of the instrument | 6. Potentiometer for speed regulation |
| 3. ON-OFF switch | 7. Micro-switch |
| 4. Electronic board | 8. Front control panel with button and indicator light |

9. Warranty

The warranty is effective from the date shown on the delivery note and refers to the registration number of the individual article.

10. Declaration of conformity

We address KARTELL spa
Via Delle Industrie 1, 20082 NOVIGLIO (MI), Italy

under our responsibility declare that the product is manufactured in conformity with the following standards:

EN 61010-1 (2001)
EN61010-2-051 (2003)
EN 61326-1 (2006)
2011/65/EU (RoHS)
2002/96/EC (WEEE)

and satisfies the essential requirements of the following directives:

Machinery directive 2006/42/EC
Low voltage directive 2006/95/EC
Electromagnetic compatibility directive 2004/108/EC
plus modifications and that the documents listed in annex I are available at the manufacturer's offices as foreseen by the machinery directive.

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