KARTELL spa

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This unit must be used for laboratory applications only.

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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This unit has been designed and manufactured in compliance with the following standards:

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UL 61010-1

CAN/CSA-C22.2 No.61010-1

The manufacturer reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The heating plate can reach a temperature of up to 370 °C.

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The magnetic fields are active on magnetic or metallic pieces (magnetic supports, CDs, floppy discs, cardiac stimulators, magnetic cards). Keep these pieces far away from the magnetic stirrer plate and the magnetic stir bars.

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Any container used must be smaller in diameter than that of the working surface to prevent any risk of tipping over and must be compatible with the temperature used. Normally, Pyrex glass containers are used.

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Observe the following precautions before using the instrument:

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Contents / Indice

EN

1.	INTRODUCTION.....	3
2.	ASSEMBLY AND INSTALLATION.....	3
3.	START-UP	3
3.1	USING THE THERMOREGULATOR VERTEX 2239.....	3
4.	END-OF-WORK OPERATIONS	4
5.	MAINTENANCE	4
5.1	CLEANING.....	4
6.	TECHNICAL DATA	4
7.	ACCESSORIES	4
8.	SPARE PARTS	4

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1. Introduction

The application of new technology has led to the creation of these modern magnetic stirrers for laboratory use whose basic features are safety, reproducibility of results, high performance and energy saving.

The structure of the unit is made out of pressure die-cast aluminium treated with epoxy resins offering a high resistance to the many chemical aggressions typical of the laboratory environment.

The meticulously studied design of the pressure die-cast casing combines aesthetic fluidity with ergonomics.

A wide central recess on both sides offers a good grip and facilitates bench-top handling whilst the handle moulded directly into the rear of the die-cast structure can be used to transport the instrument once it has been disconnected from the power supply.

The control panel is easily accessible, is positioned at a safe distance from sources of heat and is protected against possible damage caused by liquid spills thanks to a special run-off groove.

The pressure die-cast structure is designed so that accidental liquid spills cannot reach the internal parts of the unit.

The electrical connections are recessed into the rear of the instrument and offer optimum electrical protection of the external electrics in the case of liquid spills.

The engineering of the unit facilitates internal inspection if technical assistance is required.

The heating plate is made out of **cast aluminum** is treated with a particular process which ensures:

- Optimum heat distribution and a high specific power thanks to the circular configuration
- Optimum temperature homogeneity across the plate
- High resistance to thermal stress, thermal shock, chemical aggressions, scratches and surface abrasions.

NOTE: using the heating plate at high temperatures may cause discoloring. This does not alter the thermal, mechanical and chemical resistance of the plate in any way.

Magnetic stirring is generated by the Patented Composite Magnet driven by an asynchronous mono-phase brushless motor in alternating current which offers a virtually unlimited duration.

The stirring speed is adjusted by an impulse microprocessor which makes it possible to obtain optimum speed progression. Speeds of from 50 to 1.200 rpm can be selected using the speed-selector knob.

The heating plate temperature is controlled by an analog bulb thermoregulator fitted directly onto the plate itself. The heating plate temperature and all the set temperatures are thus controlled constantly.

2. Assembly and installation

The unit is delivered in special shock-resistant packaging and comes complete with: Power cable

Check that the values on the rating plate correspond to those of the power supply.

Connect the unit to the power supply using the power cable supplied.

Turn the unit on using the mains switch, the switch will light-up indicating that the unit is powered.

3. Start-up

Turn the unit on using the mains switch and place a beaker with magnetic stir bar on the heating plate.

Select the stirring speed using the speed-selector knob.

Stirring speeds of from 50 to 1200 rpm can be selected.

Select the heating plate temperature using the temperature control knob.

Temperatures of between 50 and 370 °C can be selected.

3.1 Using the thermoregulator Vertex 2239

Screw the threaded support rod (A00001069) into its seat on the back of the instrument and fasten the Vertex 2239 thermoregulator onto the support rod. Place the temperature probe in the receptacle making sure that it is completely immersed in the sample.

Connect the two instruments (TK 23 and Vertex 2239) by plugging the Vertex 2239 into the dedicated socket on the back of the TK 23.

Select the operating temperature required on the VERTEX 2239 thermoregulator.

Turn the temperature control knob on the front panel of the TK 23 to maximum.

The heating magnetic stirrer always has primary control of the heating plate temperature.

When using the Vertex 2239 thermoregulator always select the max temperature on the TK 23.

The temperature control function of the heating plate can also be used as a safety thermostat.

In this case the maximum temperature of the heating plate will not exceed the temperature setting on the TK 23 meaning that a longer heating time is required in order to reach the Vertex 2239 thermoregulator temperature setting.

4. End-of-work operations

When the unit is not in use it is advisable to turn off the mains switch, the light will go out.

5. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. It may be necessary to replace the fuses. The unit is fitted with two fuses located in the socket on the back of the unit. To replace one or more of the fuses, disconnect the unit from the power supply and use a screwdriver to lift the fuse box cover.

5.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a mild, non-flammable detergent when the plate is cold.

6. Technical data

Dimensions WxHxD	mm	165x115x280
Weight	Kg	2,6
Overall power	W	630
Power of the heating plate	W	600
Diameter of the heating plate	mm	155
Programmable temperature range	°C	50 – 370
Type of temperature control		Analog
Overtemperature protection		Yes
Stirring power	W	30
Stirring capacity	H ₂ O liter	20
Programmable speed range	rpm	50 - 1200
Type of motor control		electronic
Temperature range	°C	+5...+40
Storage temperature range	°C	-10...+60
Max humidity	%	80
Level of electrical protection CEI EN60529	IP	42
Pollution degree CEI EN61010-1		2

7. Accessories

F208B0063KA	Electronic Vertex 2239 -10...+300 °C
40000781*	Spiral electric cable complete of plug
A00001069	Support rod
A00001056	Magnetic stir bar diam. 6x35mm
A00001061	Magnetic stir bar diam. 9,5x60 mm
A00001060	Magnetic stir bar diam. 10x40mm
A00001071	Hemispheric bowl for 250 ml flasks
A00001072	Hemispheric bowl for 500 ml flasks
A00001073	Hemispheric bowl for 1000 ml flasks

* For connection to other contact thermometers other than Vertex 2239

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