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## 1. INTRODUCTION

**MEDEL ELITE** blood-pressure monitor is a automatic, digital blood pressure measuring device. The device offers a very high measurement accuracy and has been designed to provide the maximum user-friendliness.

### ATTENTION!

Before using, please read through this instruction manual carefully. For further questions on the subject of blood-pressure and its measurement, please contact your doctor.

### ATTENTION!

#### 1.1 Important information about self-measurement

- Do not forget: self-measurement means control, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.
- The pulse display is not suitable for checking the frequency of heart pacemakers! In cases of cardiac irregularity (**I.H.B.**), measurements made with this instrument should only be evaluated after consultation with the doctor.

#### 1.2 Electromagnetic interference

The device contains sensitive electronic components (Microcomputer). Therefore, avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These can lead to temporary impairment of the accuracy.

## 2. THE VARIOUS COMPONENTS OF THE BLOOD-PRESSURE MONITOR( fig.1)

- |   |                                  |
|---|----------------------------------|
| 1. LCD  | 10. Average value symbol         |
| 2. O/I button (Start / Stop)                            | 11. Heartbeat symbol             |
| 3. Memory button (MEM)                                  | 12. Battery symbol               |
| 4. Average button (T3)                                  | 13. Systolic blood pressure      |
| 5. Air socket   | 14. Diastolic blood pressure     |
| 6. AC/DC adapter socket                                 | 15. Pulse display                |
| 7. Battery compartement                                 | 16. Inflation / Deflation symbol |
| 8. Standard Cuff ( type B applied )<br>for 22~30 cm arm | 17. Date / time display          |
| 9. Memory symbol  | 18. Stored data                  |

#### Accessories ( not included) - Cuffs:

- **XL:** for arm circumference 30 - 42 cm (optional accessory)
- **XXL:** for arm circumference 42 - 48 cm (optional accessory)

**NOTE:**

*If your arm Cuff broken or not functional, Please use a new Cuff. If the new arm Cuff does not include an air plug, Please continue to the old air plug on the new arm Cuff*

**Power:**

**4 AA** Alkaline Batteries

batteries are not suitable for this monitor

Adapter (not included)

### **3 PUTTING THE BLOOD-PRESSURE MONITOR INTO OPERATION**

---

#### **3.1. Inserting the batteries**

Insert the batteries. (pic.2)

1. Remove the battery cover from the battery compartment.
2. Insert four **AA** powerful batteries into the compartment and ensure each battery is in the proper direction.
3. Replace the battery cover.

#### **ATTENTION !**

- After the battery warning symbol "" appears, the device is blocked until the batteries have been replaced.
- Please use "**AA**" Long-Life **1.5 V** batteries.
- Never leave any low battery in the battery compartment since they may leak and cause damage to the unit.
- If the blood-pressure monitor is left unused for long periods, please remove the batteries from the device.

#### **3.2 Clock and date adjustment**

1. Once you install the battery or turn off the monitor, it will enter Clock mode, and LCD will display time and date.
2. While the monitor is in Clock Mode, press the **O/I** and **MEM** button simultaneously for **2** seconds, the month will blink at first ( **pic. 3**). Press **MEM** button to set the month. Keep on pressing the button **MEM**, the number will increase fast.
3. Press the button **O/I** to confirm and to pass setting the day.
4. Repeat the instruction above to set hour and minutes .
5. The monitor will turn off automatically after **1** minute of no operation with the time and date unchanged.
6. Once you change the batteries, you should readjust the time and date.

### 3.3 Connecting Tube

Insert the connector into the socket of the cuff for cuff on the left side of the meter (**Pic.1**)

## 4. HOW TO TAKE PROPER MEASUREMENTS

---

### 4.1 Before the measurement:

- Relax for about **10** minutes. avoid eating, smoking as well as all forms of exertion directly before the measurement. All these factors influence the measurement result.
- Remove any garment that fits closely to your upper arm.
- Measure always on the same arm (normally left).
- Attempt to carry out the measurements regularly at the same time of day, since the blood-pressure changes during the course of the day.

### NOTE:

*Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is of extraordinary importance. The cuff size is dependent upon the circumference of the arm (measured in the centre). The permissible range is printed on the cuff. If this is not suitable for your use, please contact your dealer.*

### ATTENTION!

- Only use **MEDEL** original cuffs!
- A loose cuff or a sideways protruding air-pocket causes false measurement values.
- With repeated measurements, blood accumulates in the respective arm, which can lead to false results.
- Correctly executed blood-pressure measurements should therefore first be repeated after a **2-3** minute pause or after the arm has been held up in order to allow the accumulated blood to flow away .

### 4.2. Fitting the cuff

1. Pass the end of the cuff (with sewn-in rubber stopper) through the metal stirrup so that a loop is formed. The Velcro closer must be facing outwards. (Ignore this step if the cuff has already been prepared.) (**Pic 4**)
2. Push the cuff over the left upper arm so that the tube points in the direction of the lower arm. (**Pic 5**)
1. Position the cuff on the arm as illustrated in the picture (**Pic. 6**). Make certain that the lower edge of the cuff lies approximately **2 to 3** cm above the elbow and that the rubber tube leaves the cuff on the inner side of the arm )
2. Tighten the free end of the cuff and close the cuff with the closer. (**Pic 7**)

 **ATTENTION!**

The mark (ca. 3 cm long bar) must lie exactly over the artery which runs down the inner side of the arm.

- There must be no free space between the arm and the cuff as this would influence the result. Clothing must not restrict the arm. Any piece of clothing which does (e.g. a pullover) must be taken off.
- Secure the cuff with the Velcro closer in such a way that it lies comfortably and is not too tight. Lay the arm on the table (palm upwards) so that the cuff is at the same height as the heart. Make sure that the tube is not kinked. **(Pic 8-9).**
- Remain seated quietly for two minutes before you begin the measurement.

**NOTE:**

*If it is not possible to fit the cuff to the left arm, it can also be placed on the right one. However all measurements should be made using the same arm.*

**4.3 Measuring procedure**

After the cuff has been appropriately positioned, the measurement can begin

- Press **O/I** button. All display symbols appear for **1** second. If the monitor has stored results, the LCD will momentarily display the most recent one. If no result has been stored, all the digits show **0**. Then **0** will appear on LCD on the DIA column. **(pic. 10)**. The monitor is ready for measurement.
- Then the cuff inflates automatically, "▼" symbol is shown. **(pic. 10)**
- When appropriate pressure is reached, inflation stops and the cuff starts to deflate, "▼" symbol is shown. **(pic. 11)**
- Once a pulse is detected, the "♥" symbol blinks. **(pic. 11)**
- When the measurement completes, the systolic / diastolic pressure and pulse are shown simultaneously. **(pic. 12)**
- In case the measurement shows irregular heartbeat (**I.H.B.**), the "♥" symbol is displayed. In this case the results of the measurement have deviated from a normal situation. Then the measurement needs to be repeated. In the majority of cases, this is not cause for concern. If the symbol "♥" appears regularly (for example several times a week during different daily measurements) it is advisable to discuss it with your medical advisor. **(pic. 12).**

#### 4.4 How to use the memory function:

- In Clock mode ( Off mode) , press **MEM** button. The number of stored data is showed (**pic14**), sudden the latest result will be displayed with date and time (**pic15**).
- Press **MEM** button repeatedly to review the results measured .
- Press **T3** button, the average value of the last **3** measurements (Done within a period of **15** minutes) will be displayed, with **AVG** symbol blinking. (**pic.16**) .
- To display stored data press again **MEM** button.

#### NOTE:

When displaying the stored results, the monitor will turn off automatically after **10** second of no operation. You can also press the button **O/I** to turn off the monitor manually..

#### 4.5 To delete all values stored in the memory:

In Clock mode ( Off mode) press **MEM** button , after the latest data stored will be displayed press again and hold **MEM** button for **3** sec. All results will be deleted display. Press the button **MEM** or **O/I**, the display will turn off. (**pic.17**)

## 5. ERROR MESSAGE/MALFUNCTION

Follow the below recommended actions and press **O/I** button again to initiate another measurement.

Error code	Possible cause(s)	Recommended Action
Er 0	Pressure system is unstable before measurement	Don't move and try again.
Er 1	Fail to detect systolic pressure	
Er 2	Fail to detect diastolic pressure	
Er 3	Pneumatic system blocked or cuff is too tight during inflation	Apply the cuff correctly and try again

<b>Er 4</b>	Pneumatic system leakage or cuff is too loose during inflati	Apply the cuff correctly and try again
<b>Er 5</b>	Cuff pressure above 300mmHg	Measure again after five minutes. If the monitor is still abnormal, please contact the local distributor or the factory.
<b>Er 6</b>	More than 3 minutes with cuff pressure above 15 mmHg	
<b>Er 7</b>	Problem in internal memory device	
<b>Er 8</b>	Device parameter checking error	
<b>Er A</b>	Pressure sensor parameter error	

### 5.1 Other possible malfunctions:

<b>Problem</b>	<b>Possible cause(s)</b>	<b>Recommended Action</b>
Only the symbol "  " flashes	Batteries are drained.	Replace all batteries with new ones.
The blood pressure values displayed are abnormally high or low.	The cuff is not wrapped correctly or not at the level of the heart.	Wrap the cuff correctly and rise your hand so that the cuff is at the same level of your heart.
	Too much stress is applied on your shoulder or arm.	Relax yourself and make measurement.
	You move your arm or arm muscles during measurement.	Remain still and do not move /contract the muscles during measurement.
The <b>(I.H.B.)</b> symbol "  " is displayed, but my heart beat rate should be normal.	You move your arm or arm muscles during measurement.	Remain still and do not move /contract the muscles during measurement.

## 6. MAINTENANCE AND REPAIRS

---

In case of failure, address to qualified personnel authorized by **MEDEL Group S.p.A.** Never open the device. There are no user-serviceable parts within. The device requires no lubrication or maintenance.

## 7. CARE, MAINTENANCE AND STORAGE

---

1. For cleaning only use a soft, lightly moistened cloth.
2. Do not use petroleum spirits, thinners or similar solvents.
3. Do not wash the cuff.
4. Take care to ensure the cuff is not damaged by sharp or pointed objects such as scissors, knives, etc.
5. Take the batteries out if the device is not used for a long period.
6. In case the device becomes malfunction due to the influence of strong electromagnetic fields, power off the device and repeat the procedure again.
7. Disposal of the device (and used batteries) should be carried out according to the national regulations for the disposal of electronic products.

## 8. REFERENCE TO STANDARDS

---

Device standard: Device corresponds to the requirements of the European standard for non-invasive blood-pressure monitor.

- EN1060-1 Non-invasive sphygmomanometers - General requirements.
- EN1060-3 Non-invasive sphygmomanometers - Supplementary requirements for electro-mechanical blood pressure measuring systems.
- EN1060-4 Non-invasive sphygmomanometers - Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers.
- EN 60601-1
- ANSI / AAMI SP10, NIBP - requirements.
- EN 60601-1-2 Electromagnetic compatibility:

**CE 0123** In conformity with European Directive on Medical Devices 93/42/EEC

**Separate collection for waste electrical and electronic equipment.** (Applicable to European Union countries and countries with separate collection systems) This symbol on the product or its literature indicates that the product complies with the Electrical and electronic equipment and should not be disposed of as household waste. The user is responsible for the transfer of the equipment at end of life with the appropriate collection facilities, subject to the penalties provided for by current legislation on waste. For more detailed information collection systems available contact your local waste disposal service.

## 9. PRODUCT SPECIFICATION

---

Display :	Large digital backlight display
Technology:	Oscillometric method
Classification :	Internally powered equipment, with type B applied part
Measurement range:	Pressure : 0 to 295 mmHg Heart rate : 40 to 180 times/minute
Accuracy :	pressure : +/- 3 mmHg Heart rate : +/- 5%
Memory :	Maximum 60 sets of data (Systolic / Diastolic pressure, Heart beat)
Power supply:	4 AA 1.5V (batteries)
Accessories:	Standard cuff: for arm circum. 22 - 30 cm
Optional Accessories:	a) AC/DC Adapter 6Vdc; 600 mA b) XL cuff: for arm circum. 30 - 42 cm c) XXL cuff: for arm circum. 42 - 48 cm
Operation conditions:	Temperature : +5°C to +40° C; Relative Humidity: < 90 % pressure: 80 to 105 kPa
Transport/Storage conditions:	Temperature : -20°C to +55°C ; Relative Humidity : < 95% pressure: 80 to 105 kPa
Weight of main unit :	426 g (without batteries)
External dimensions:	171 X 131 X 62 mm

**NOTE:** Specifications are subject to change without prior notice. Meanings of various symbols.

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© Medel Group S.p.A. 2010