...for your health

DISPLAY

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Instruction manual

EN

1. INTRODUCTION	25
1.1. MEDEL DISPLAY features	25
1.2. Important information about self-measurement	25
2. IMPORTANT INFORMATION ON THE SUBJECT OF BLOOD-PRESSURE AND ITS MEASUREMENT	26
2.1. How does high/low blood-pressure arise?	26
2.2. Which values are normal?	26
2.3. What can be done, if regular high/low values are obtained?	27
3. THE VARIOUS COMPONENTS OF THE BLOOD-PRESSURE MONITOR	29
4. BATTERY INSTALLATION / REPLACEMENT	31
4.1. How to Install Battery	31
5. PROPER USE OF THE ARM CUFF	32
6. QUICK REFERENCE GUIDE	32
7. TAKE A MEASUREMENT	34
8. USE THE MEMORY FUNCTION	35
8.1. How to Use the Memory Function	35
8.2. The display method of memory values before turning off the unit when a measurement completed	35
8.3. The display method of memory values when the unit is powered off	35
9. EXCEPTIONAL SITUATIONS	37
10. MAINTENANCE AND REPAIRS	39
11. CARE AND MAINTENANCE, RECALIBRATION	40
12. REFERENCE TO STANDARDS	41
13. TECHNICAL SPECIFICATIONS	42
14. WARRANTY	43

1. INTRODUCTION

1.1. MEDEL DISPLAY features

MEDEL DISPLAY blood-pressure monitor is a fully automatic, digital blood-pressure measuring device for use on the upper arm, which enables very fast and reliable measurement of the systolic and diastolic blood-pressure as well as the pulse frequency by way of the oscillometric method of measuring.

The device offers a very high and clinical tested measurement accuracy and has been designed to provide the maximum user-friendliness.

Thanks to the large MEDEL DISPLAY LCD display, both operational status and cuff pressure are visible during measurement.

Before using, please read through this instruction manual carefully and then keep it in a safe place. For further questions on the subject of blood-pressure and its measurement, please contact your doctor.

Attention!

1.2. Important information about fuzzy-measurement

- Do not forget: self-measurement means control, not diagnosis or treatment. Unusual values must always be discussed with your doctor.
Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.
- The pulse display is not suitable for checking the frequency of heart pacemakers!
- In cases of cardiac irregularity (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with the doctor.

Electromagnetic interference:

The device contains sensitive electronic components (Microcomputer). Therefore, avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These can lead to temporary impairment of the measuring accuracy.

2. IMPORTANT INFORMATION ON THE SUBJECT OF BLOOD-PRESSURE AND ITS MEASUREMENT

2.1. How does high/low blood-pressure arise?

The level of blood-pressure is determined in a part of the brain, the so-called circulatory centre, and adapted to the respective situation by way of feedback via the nervous system.

To adjust the blood-pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. The latter is effected by way of fine muscles in the blood-vessel walls.

The level of arterial blood-pressure changes periodically during the heart activity: during the 'blood ejection' (Systole) the value is maximal (systolic blood-pressure value), at the end of the heart's 'rest period' (Diastole) it is minimal (diastolic blood-pressure value).

The blood-pressure values must lie within certain normal ranges in order to prevent particular diseases.

2.2 Which values are normal?

Blood pressure is too high if at rest, the diastolic pressure is above 90mmHg and/or the systolic blood-pressure is over 160mmHg. In this case, please consult your doctor immediately. Long-term values at this level endanger your health due to the associated advancing damage to the blood vessels in your body.

Should the systolic blood-pressure values lie between 140mmHg and 160mmHg and/or the diastolic blood-pressure values lie between 90mmHg and 100mmHg, likewise, please consult your doctor. Furthermore, regular self-checks will be necessary.

With blood-pressure values that are too low, i.e. systolic values under 105mmHg and/or diastolic values under 60mmHg, likewise, please consult your doctor.

Even with normal blood-pressure values, a regular self-check with your blood-pressure monitor is recommended. In this way you can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, please

keep a record of the level of your blood pressure by carrying out regular self-measurements at specific times of the day. Show these values to your doctor. **Never use the results of your measurements to alter independently the drug doses prescribed by your doctor.**

Table for classifying blood-pressure values (units mmHg):

Range	Systolic Blood-pressure	Diastolic Blood-pressure	Measures
Hypotension	lower than 100	lower than 60	Check with your doctor
Normal range	between 100 and 140	between 60 and 90	Self-check
Mild hypertension	between 140 and 160	between 90 and 100	Consult your doctor
Moderately serious hypertension	between 160 and 180	between 100 and 110	Consult your doctor
Serious hypertension	higher than 140	higher than 110	Urgent! See your doctor
Specific systolic hypertension	higher than 140	lower than 90	Consult your doctor

Further information

- If your values are mostly standard under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called "labile hypertension". Please consult your doctor if you suspect that this might be the case.
- Correctly measured diastolic blood-pressure values above 120mmHg require immediate medical treatment.

2.3. What can be done, if regular increased/low values are obtained?

- Please consult your doctor.
- Increased blood-pressure values (various forms of hypertension) are associated long- and medium term with considerable risks to health. This concerns the arterial blood vessels of your body, which are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can be the result.

Furthermore, with long-term continuously increased blood-pressure values, the heart will become structurally damaged.

- c) There are many different causes of the appearance of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organic malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are measures which you can take, not only for reducing a medically established high blood pressure, but also for prevention. These measures are part of your general way of life:

A) Eating habits

- Strive for a normal weight corresponding to your age. Reduce overweight!
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

B) Previous illnesses

- Follow consistently any medical instructions for treating previous illness such as:
 - Diabetes (Diabetes mellitus)
 - Fat metabolism disorder
 - Gout

C) Habits


- Give up smoking completely.
- Drink only moderate amounts of alcohol.
- Restrict your caffeine consumption (Coffee).

D) Physical constitution:


- After a preliminary medical examination, do regular sport.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your sporting activities. He will advise you regarding the type and extent of types of sport that are possible for you.

3. THE VARIOUS COMPONENTS OF THE BLOOD-PRESSURE MONITOR

Main Unit


Display

Accessorial

Standard Cuff:
Fit for 22cm~32mm range
of upper arm perimeter
(type BF applied part)


⚠ If your arm cuff is broken or not functional, Please use a new cuff.

⚠ A new arm cuff does not include an air plug. Please continue to use the old air plug on the new arm cuff.

Power:

4 AA Alkaline Batteries (not included).

Adapter (not included).

**Accessories (not included)**


Small: for arm circumference 15-22 cm (available as optional accessory)
X-Large: for arm circumference 32-48 cm (available as optional accessory)

Instruction Manual:


Please keep the instruction manual well after using.

4. BATTERY INSTALLATION / REPLACEMENT**4.1. How to Install Battery**

1 Remove the battery cover from the battery compartment.


2 Insert four AA powerful batteries into the compartment and ensure each battery is in the proper direction.

3 Replace the battery cover.


⚠ The battery cover will be closed tight when you hear a clicking sound.

**⚠ Caution!**

- Use four identical 1.5V alkaline batteries type AA LR6.
- Replace all of the batteries if the symbol of low battery appears on the display.
- Never leave any low battery in the battery compartment since they may leak and cause damage to the unit.

5. PROPER USE OF THE ARM CUFF

- 1 Make sure the air plug is inserted properly in the main unit.

6. QUICK REFERENCE GUIDE


- 1 Please take off thick clothing and keep yourself quiet for 5~10 minutes before taking a measurement.


- 2 Remove all clothes from your upper arm allowing the cuff to fit directly on the skin.


- 3 Pass the end of the cuff through the metal, and keep the hose outward.

- 4 The hose runs down the inside of your arm. The bottom of the cuff should be approximately 1~2 cm above your elbow.

- 5 Wrap the cuff properly around your upper arm. Not more than 2-finger-room allowed between the cuff and your arm.

- 6 Have your upper arm relax, keep your palm upward and your fingers curving naturally, then turn on the unit and start to measure.

7. TAKE A MEASUREMENT

1 Apply the cuff as described.

2 Press the **START/STOP** button.

All display symbols appear on the screen. The cuff starts to inflate automatically.

3 Inflation stops automatically and the measurement is started.


The heartbeat symbol flashes once a pulse is detected.

4 The blood pressure and pulse rate are displayed when the measurement is completed.

The cuff is deflated automatically, and the results are stored in the memory with date and time.

The **IHB** symbol will be displayed if arrhythmia is detected.

5 Press the **START/STOP** button to turn off the unit off.


The unit will automatically turn off after 2 minutes if you have forgotten to turn it off.

8. USE THE MEMORY FUNCTION


8.1. How to Use the Memory Function

- The unit stores the blood pressure and pulse rate in the memory after completing a measurement every time. It can automatically store 60 sets of measurement values. The oldest record is deleted to save the latest measurement value when more than 60 set.
- The unit also calculates an average reading based on the values of the latest 3 times measurement taken. You can get the average measurement value when you press the memory button for the first time. Press the memory button for the second time, you can get the latest measurement value, so on and so forth.
- Press the left "MEM" button. The average value symbol will display on the screen with the average reading.

8.2. The display method of memory values before turning off the unit when a measurement completed.

8.3. The display method of memory values when the unit is powered off.

Unit Setting (kPa or mmHg)

Long press the memory button for 3~4 seconds until the units displayed when the unit powered off.

Memory Values Cleaning out

When reading the measurement values, press and hold the "MEM" button, then press START/STOP button for about 3-4 seconds.

9. EXCEPTIONAL SITUATIONS

Error Indicators

The following symbol will appear on the display when measuring abnormal.

Symbol	Problem	Correction
Err	The course of inflating appears error.	Wrap the cuff properly.
		Make sure that the air plug is properly inserted in the unit. Inflate over again after ensuring.
		Do not move and talk while measuring. Measure over again according to correct way.
Lo	When the batteries power are too low.	Replace all of the worn batteries with new ones.

Eliminate Fault

Problem	Causes and Solutions
No power.	Replace all the worn batteries with 4 new ones.
No display appears on the display screen.	Check the battery installation for proper placement of the batteries polarities.
Measurement values appear too high or too low.	Blood pressure varies constantly. Many factors including stress, time of day, how you wrap the cuff, may affect your blood pressure. Review the sections "Proper Use of the Arm Cuff" and "Take a Measurement".

Other possible malfunctions and their elimination

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

Malfunction	Remedy
The display remains empty when the instrument is switched on although the batteries are in place.	1) Check whether the batteries are installed with the correct polarity and correct if necessary. 2) If the display is unusual, remove the batteries and then them, or exchange them for new.
The pressure does not rise although the pumping ball is pressed.	1) Check the connection of the cuff tube and connect properly if necessary. 2) The cuff could be damaged, change it.
The instrument frequently fails to measure the blood pressure values, or the values measured are too low (too high).	1) Fit the cuff correctly on the arm. 2) Before starting the measurement make sure that the cuff is not fitted too tightly and that there is no rolled-up sleeve exerting pressure on the arm above the measuring position. Take off articles of clothing if necessary. 3) Measure the blood pressure again in complete peace and quiet.
The blood pressure values measured differ from those measured by the doctor.	1) Record the daily development of the values and consult your doctor about them.

Further information

The level of blood-pressure is subject to fluctuations even with healthy people. Important thereby is, that **comparable measurements always require the same conditions (Quiet conditions)!**

If, in spite of observing all these factors, the fluctuations are larger than 15mmHg, and/or you hear irregular pulse tones on several occasions, please consult your doctor.

Should the device not start working properly, address immediately to an authorized technical service centre or to MEDEL CUSTOMER SERVICE.

10. MAINTENANCE AND REPAIRS

In case of failure, address to qualified personnel authorized by Medel.

Never open the device. There are no user-serviceable parts within. The device requires no lubrication or maintenance.

11. CARE AND MAINTENANCE

Care

- Keep the unit in the storage case when not in use.
- Clean the unit with a soft dry cloth. Do not use any abrasive or volatile cleaners.
- Never immerse the unit or any components in water.

Maintenance


- Do not clean the body and cuff with naphtha, thinner or gasoline etc.


- Do not wet the cuff or attempt to clean the cuff with water.

- Store the unit in a clean and dry location. Do not subject the unit to extreme hot or cold temperature, humidity and direct sunlight.

- Remove the batteries if the unit will not be used for 3 months or longer. Always replace all of the batteries with new ones at the same time.

Use the unit consistent with the instructions provided in this manual.
Use only authorized parts and accessories.

12. REFERENCE TO STANDARDS

Device standard: Device corresponds to requirements of the European standard for non-invasive blood-pressure monitor.

-EN 1060-1 Non-invasive sphygmomanometers - Part 1: General requirements

-EN 1060-3 Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro-mechanical blood-pressure measuring systems.

-EN 1060-4 Non-invasive sphygmomanometers - Part 4: Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers.

-EN 60601-1

-ANSI/AAMI SP10

Electromagnetic compatibility: Device fulfils the stipulations of the European standard IEC60601-1-2

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CE 0123

The stipulations of the EU-Guidelines 93/42/EEG for Medical Products Class IIa have been fulfilled.

Correct disposal of the product (electric and electronic waste)

(Applicable in the countries of the European Union and in the countries using systems for separate collection of waste)

The symbol on the product or its manuals indicates that the product complies with the directive WEEE (waste electrical and electronic equipment) and must not be disposed of among domestic waste. At the end of its lifecycle, the user shall dispose of the equipment in specific collection areas. Otherwise, the sanctions foreseen by the current laws about waste are applicable. For detailed information concerning available collection systems, refer to the local waste disposal services.

ATTENTION: Read the instructions

13. TECHNICAL SPECIFICATIONS

Classification:	Internally powered equipment, with type BF applied part
Display:	LCD Digital Display
Measuring Principle:	Oscillometric Method
Measurement Localization:	Upper Arm
Measurement Range:	Pressure: 0 to 299 mmHg (0 ~ 39.9 kPa) Pulse: 40 to 180 Pulses/min Pressure: ± 3 mmHg (0.4 kPa) Pulse: $\pm 5\%$ of reading
Accuracy:	Pressure: ± 3 mmHg (0.4 kPa) Pulse: $\pm 5\%$ of reading
LCD Indication:	Pressure: 3 Digits Display of mmHg/kPa Pulse: 3 Digits Display Symbol: Memory/ Average/ IHB/ Heartbeat/ Low Battery, etc.
Inflation:	Automatic by Internal Pump
Deflation:	Automatic by Air Valve
Rapid Air Release:	Automatic by Air Valve
Memory Function:	Memory 60 Sets of Measurement Values
Power Source:	4 AA Alkaline Batteries
Automatic Power-Off:	In 2 Minutes
Main Unit Weight:	Approx. 330g (Batteries not included)
Accessories:	Standard Cuff
Optional Accessories:	a) AC/DC adaptor 6Vdc ; 600mA b) Small cuff: for arm circum. 15-22 cm c) X-Large cuff: for arm circum. 32-48 cm
Operating Environment:	Temperature: $+10^{\circ}\text{C} \sim +40^{\circ}\text{C}$ Humidity: 15 ~ 90%RH (Noncondensing) Barometric Pressure: 105 ~ 80 kPa
Storage/transport Environment:	Temperature: $-20^{\circ}\text{C} \sim +50^{\circ}\text{C}$ Humidity: 10 ~ 90%RH (Noncondensing) Barometric Pressure: 105 ~ 80 kPa

Technical modifications reserved!

14. WARRANTY CERTIFICATE

Warranty conditions

- The device is covered by a 2 years warranty from date of purchase for any defects in the design or in the materials used.
- The warranty consists in the replacement and/or repair free of charge of originally defective components.
- The warranty does not cover the accessories supplied with the device and those parts subject to normal wear and tear.
- The device must be repaired by authorized technical service centres only.
- The device must be sent to the authorized service centre for repair within 8 days after the defect is noticed.
- Transport costs for the device are at user's charge.
- Any repair out of warranty conditions are at user's charge.
- The warranty does not cover failures resulting from misuse, abuse of the product or in case the damage is not due to manufacturer's liability (accidental fall, rough transportation, etc.).
- The warranty does not involve any compensation for any kind of damages to persons or things, either direct or indirect, occurred when the device is not operating.
- The warranty is valid from the date of purchase certified by the ticket or the invoice which have to be peremptorily attached to the warranty coupon.
- The lack of warranty coupon properly filled in and validated by a certification of purchase voids the warranty.

**COUPON TO BE RETURNED IN CASE OF REPAIRS****THE WARRANTY IS VALID ONLY IF
ACCOMPANIED BY THE TICKET/INVOICE**

Product type: _____

Model: _____

Serial No.: _____

Date of purchase: _____

BUYER DATA

Name: _____

Address: _____

Phone number: _____

Defect description: _____

Signature for acceptance of the above mentioned warranty conditions

☐ I hereby authorise use of the above information in accordance with Privacy Law 675/96.

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© Medel group S.p.A. 2009