# Snowden-Pencer™

- en 1.0 Switch-Blade™ Reusable Handle and Disposable Scissor Tips
- fr Poignée réutilisable et pointes de ciseaux jetables 1.0 Switch-Blade™
- de Produkte – 1.0 Switch-Blade™ Wiederverwendbarer Griff und Einweg-Scherenspitzen
- it Impugnatura riutilizzabile e punte di forbice monouso 1.0 Switch-Blade™
- es Mango reutilizable 1.0 Switch-Blade™ y puntas de tijera desechables
- fi tuotteet – 1.0 Switch-Blade™ uudelleenkäytettävä kahva ja kertakäyttöiset leikkauskärjet

**USA** Rx Only

**STERILE**

For Single Use Only

Snowden-Pencer and Switch Blade are trademarks or registered trademarks of CareFusion Corporation or one of its subsidiaries.

© 2012 CareFusion Corporation, or one of its subsidiaries.  
All rights reserved.

**CareFusion**  
1500 Waukegan Road  
McGaw Park, IL 60085 U.S.A.

CF36-1803 / 26-2910-A • 2012-04

**CareFusion**

## Catalog Numbers

**Table 1 (Reusable handles - non-sterile)** SP90-1050, SP90-1150, SP90-1250

**Table 2 (Single use scissor tips - sterile)** E89-5100B, E89-5300B, 89-5105B, 89-5200B, 89-5300B

**Table 3 (Accessories)** 88-9860, 88-9861

## Indications For Use

Endoscopic devices are designed to transmit cutting and dissecting force through delicate mechanisms to delicate working tips through a cannula.

Electrosurgical devices are intended to be connected to Electrosurgical Units via cables to allow monopolar cutting and coagulation.

## How Supplied

The switch blade handles are supplied as non sterile. Cleaning and sterilization must occur prior to use. (See table 1)

The switch blade scissor tips are supplied as sterile. (See table 2)

Do not use if single use sterile package is damaged.

## Limitations on Reprocessing

Repeated processing on the reusable devices has minimal effect. End of life is normally determined by wear and damage due to use.

The switch blade scissor tips are for single use only. Reuse may result in a non-functional product or contribute to cross contamination.

## Warnings

Devices shall be used in accordance with these instructions for use. Read all sections of this insert prior to use. Improper use of this device may cause serious injury. In addition, improper care and maintenance of the device may render the device non-sterile prior to patient use and cause a serious injury to the patient or health care provider.

When not in use, store in a location that is isolated from the patient.

Prior to use of the device in patients who have cardiac pacemakers or other active implants, a possible hazard exists because interference with the action of the device or damage to the device may occur. Consult device manufacturer's instructions for use. In case of doubt, qualified advice should be obtained.

Gas embolism may result from over-insufflation of air, inert gas prior to HF surgery or laser assist gas.

The Snowden Pencer instrument is a monopolar device and should only be connected to a generator that is compatible with monopolar devices.

Do not exceed the 0.5kVp maximum rating specified for this device. Adjust the electrosurgical generator to a setting appropriate for the maximum peak output voltage.

A dispersive electrode, such as a grounding pad, should be used with the device and generator to prevent burns/injury to the patient when using the monopolar device. Ensure the pad is properly sized and properly placed and always recheck the pad when patient is moved. It is recommended that the device is used with an HF generator that contains contact quality monitoring (return electrode monitoring) with a signal to indicate there is contact to the patient.

Conductive fluids, (e.g., blood or saline) in direct contact with an active electrode may carry electrical current or heat, which may cause unintended burns to the patient.

Do not use Snowden Pencer devices with energy in the presence of combustible/explosive gases or other flammable anesthetics, such as nitrous oxide (N<sub>2</sub>O) and oxygen. Non-flammable agents used for cleaning and disinfecting, or as solvents of adhesive should be allowed to evaporate before the application of HF surgery. Attention should be called to the danger of ignition of endogenous gases. Some materials, for example cotton, wool, and gauze, when saturated with oxygen may be ignited by sparks produced in normal use of the HF surgical equipment.

There is a risk of pooling of flammable solutions under the patient or in body depressions such as the umbilicus, and in body cavities such as the vagina. Any fluid pooled in these areas should be mopped up before HF surgical equipment is used.

To reduce capacitive coupling, the device should only be activated when in position to deliver energy to the target tissue.

Activating the electrosurgical unit while simultaneously utilizing the suction/irrigation may alter the path of the electrical energy away from target tissue.

Improper connections of accessories may result in inadvertent accessory activation or other potentially hazardous conditions.

Connect adaptors and accessories to the electrosurgical unit only when the unit is off. Failure to do so may result in an injury or electrical shock to the patient or operating room personnel.

Start with the lowest possible power setting on the electrosurgical generator. To achieve the desired cutting and coagulation, check the patient circuit and then gradually increase the power setting. **Note:** the output power selected should be as low as possible for the intended purpose.

Do not activate the electrosurgical unit until the device has made contact with the patient. Interference produced by the operation of the HF surgical equipment may adversely influence the operation of other electronic equipment.

Use of low frequency power could result in neuromuscular stimulation.

## Cautions

Always use the device in field of view of the camera and use caution when using other devices so as not to damage device insulation.

Smoke-plume extraction may be necessary during electrosurgical cutting or coagulation.

Device (including insulation) and device accessories must be inspected prior to use to ensure integrity. In particular, electrode cables and endoscopically used accessories should be checked for possible damage to the insulation.

To prevent the possibility of electrical shocks or burns, do not use devices with breaks in the insulation.

When used in conjunction with laser equipment, always use caution (i.e. eye protection).

If there are any variations between these instructions and either your facility's policies and/or your cleaning/sterilization equipment manufacturer's instructions, those variations should be brought to the attention of the appropriate responsible hospital personnel for resolution before proceeding with cleaning and sterilizing your devices.

Use of device for a task other than that for which it is intended will usually result in a damaged or broken device.

## Examples

1. Use of a delicate dissector as a grasper.
2. Use of a delicate scissor to cut suture.
3. Use of a dissector to remove clips.
4. Use of a 5mm grasper or dissector instead of a 10mm claw extractor forceps to remove excised tissue through cannula.

Prior to use, inspect device to ensure proper function and condition. Do not use devices if they do not satisfactorily perform their intended function or have physical damage.

Inspect insulation. Any interruptions in the coating may compromise the safety of the device. To prevent the possibility of electrical shocks or burns, do not use devices with breaks in the insulation.

Avoid mechanical shock or overstressing the devices. Close distal ends prior to insertion or removal through cannulas.

Devices are designed to be held with one finger and thumb in ring handles. If force is applied in any other fashion (such as holding handles in a whole hand pistol grip), the mechanism can be easily overstressed resulting in damage or breakage.

Only the cleaning and sterilization processes which are defined within these instructions for use have been validated.

Always use caution when inserting or removing devices through cannula. Lateral pressure on the device during removal can damage the working tip, shaft of the device and/or insulation. Be sure the tips are closed and the device is pulled straight out until completely clear of cannula to avoid catching the valve assemblies in cannulas or dislodging the cannula.

Use only neutral pH (6-8) detergent solutions.

## Rated Voltage

The maximum recurring peak voltage for devices applicable to this IFU is 0.5kVp.

**Note:** Failure of the electrosurgical generator may result in an unintended increase in output power.

## Electrode Cable Compatibility

All of Snowden-Pencer monopolar electrosurgical devices are compatible with the Unipolar Bovie Cord with universal plug (88-9199) or equivalent.

Pre-processing instructions

Initiate cleaning of device within 2 hours of use.

Transport devices via the institutions established transport procedure.

Remove excess gross soil as soon as possible after use by rinsing or wiping the device.

All devices must be processed in the completely open position (i.e. jaws, flush-ports etc).

Devices should not be disassembled.

## Manual Cleaning

1. Ensure all pre-processing instructions are followed prior to cleaning.
2. Prepare the enzymatic / neutral pH detergent solution, utilizing tap water with a temperature range of 27°C to 44°C (81°F to 111°F), per manufacturer's instructions.
3. Place device in the open/relaxed position, with flush port open, and completely immerse in the detergent solution and allow device to soak for a minimum of 5 minutes. Actuate all movable parts during the initiation of the soak time.
4. Using a soft bristled brush, remove all visible soil from the device. Actuate device while brushing, paying particular attention to hinges, crevices and other difficult to clean areas. **Note:** It is recommended that the detergent solution is changed when it becomes grossly contaminated (bloody and/or turbid).
5. For lumen devices, use a soft bristled brush with a brush diameter and length that is equivalent to lumen diameter and length. Scrub the lumen (i.e. angulated/nonangulated positions) until no visible soil is detected in the lumen rinsing step below.
6. For lumen devices, place the device into the open/relaxed position with the distal tip pointed down. Flush the device with a minimum of 50 ml of detergent solution utilizing a temperature range of 27°C to 44°C (81°F to 111°F), by using the flushing port located on the handle/shaft. Repeat the flush process a minimum of 2 times (i.e. total of 3 times), ensuring all fluid exiting the lumen is clear of soil.

7. For lumen devices, if visible soil is detected during the final lumen flush, re-perform brushing and flushing of the lumen.
8. Rinse the device by completely immersing in tap water with a temperature range of 27°C to 44°C (81°F to 111°F), for a minimum of 30 seconds to remove any residual detergent or debris.
9. For lumen devices, following the rinsing step, place the device into the open/relaxed position with the distal tip pointed down. Flush the device with a minimum of 50 ml of tap water utilizing a temperature range of 27°C to 44°C (81°F to 111°F), by using the flushing port located on the handle/shaft. Repeat the flush process a minimum of 2 times (i.e. total of 3 times).
10. Dry the device with a clean, lint-free towel.
11. For lumen devices, manipulate the device to allow rinse water to drain from the lumen.
12. Visually examine each device for cleanliness.
13. If visible soil remains, repeat cleaning procedure.

#### Automatic Cleaning

1. Ensure all pre-processing instructions are followed prior to cleaning.
2. Clean the devices via the Automatic Cleaning parameters below.

PHASE	MINIMUM RECIRCULATION TIME	WATER TEMPERATURE	DETERGENT TYPE AND CONCENTRATION (IF APPLICABLE)
Pre-Wash <sup>1</sup>	15 Seconds	Cold Tap Water 1°C - 16°C (33°F - 60°F)	N/A
Enzyme Wash	1 Minute	Hot Tap Water 43°C - 82°C (110°F - 179°F)	<ul style="list-style-type: none"> <li>• Detergent: pH-neutral/enzymatic detergent</li> <li>• Concentration: Per the detergent manufacturer's recommendations</li> </ul>
Wash 1	2 Minute	Tap Water 43°C - 82°C (110°F - 179°F)	<ul style="list-style-type: none"> <li>• Detergent: pH-neutral cleanser</li> <li>• Concentration: Per the detergent manufacturer's recommendations</li> </ul>
Rinse 1	15 Seconds	Tap Water 43°C - 82°C (110°F - 179°F)	N/A
Pure Rinse	10 Seconds	Purified Water 43°C - 82°C (110°F - 179°F)	N/A
Drying	N/A	N/A	N/A

3. Manipulate the device to allow rinse water to drain from the lumen.
4. If visible moisture is present dry the instrument with a clean, lint-free towel.
5. Visually examine each instrument for cleanliness.
6. If visible soil remains, repeat cleaning procedure.

#### Inspection/Maintenance

Proper care and handling is essential for satisfactory performance of any surgical device. The previous cautions should be taken to ensure long and trouble-free service from all your surgical devices. Inspect devices before each use for broken, cracked, tarnished surfaces, movement of hinges, and chipped or worn parts. If any of these conditions appear, do not use the device. Return devices to an authorized repair service center for repair or replacement.

Before sterilizing, lubricate the device with instrument milk or a steam permeable/ water soluble lubricant, following the lubricant manufacturer's instructions. Check for smooth action of all moving parts.

Immediately proceed to recommended steam sterilization. Do not store instrument prior to sterilization.

Let devices drip dry for three (3) minutes before packaging for sterilization.

#### Packaging

Devices can be loaded into dedicated packaging systems. Sterilization wrap material must be cleared for the applicable sterilization modality by your country's regulatory body. Use in accordance with packaging manufacturer's sterilization instructions being sure to protect jaws and cutting edges from damage.

#### Sterilization

All devices must be processed in the completely open and disassembled configuration. Note that applicable device disassembly should not require any mechanical tooling (i.e. screwdriver, pliers etc.) unless otherwise indicated.

All devices shall be positioned to allow sterilant contact of all surfaces.

All devices with concave surfaces shall be configured so that water pooling does not occur.

#### Prevacuum Steam Sterilization Parameters

Minimum Preconditioning Pulses: 3

Minimum Temperature: 132°C (270°F)

Minimum Exposure Time: 3 minutes

Minimum Dry Time: 30 minutes

Sterilization Configuration: Wrapped (2 layer 1-ply or 1 layer 2-ply)

## Gravity Steam Sterilization Parameters

Minimum Temperature: 132°C (270°F)

Minimum Exposure Time: 15 minutes

Minimum Dry Time (lumens): 45 minutes

Sterilization Configuration: Wrapped (2 layer 1-ply or 1 layer 2-ply)

## Storage


After sterilization, devices must remain in sterilization packaging and be stored in a clean, dry environment.

## Assembly Instructions

Assembly (Applicable to catalog numbers SP90-1050, SP90-1150, SP90-1250, E89-5100B, E89-5300B, 89-5100B, 89-5105B, 89-5200B, 89-5300B)


### Assembly

1. Open sterile packaging and remove scissor tips.
2. Ensure scissor tips are closed and protective guard is on.
3. Grasp device shaft behind tip of device.
4. Open handle to fully extend center rod. (This is best accomplished with shaft vertical and handle end down).
5. Slide rod into scissors slot.
6. Push scissor tips down (toward device) and hand tighten scissors by turning clockwise.
7. Remove protective guard. Device is ready for use.

### Disassembly

1. With blades closed, replace protective guard on scissor tips.
2. Grasp device shaft below scissors with shaft vertical and handle end down.
3. Grasp black area of scissor tips and loosen by turning counterclockwise.
4. Slide scissors off of device rod.
5. Discard used scissor tips in accordance with hospital policy.

## Warranty

CareFusion guarantees every surgical device bearing the Snowden-Pencer brand name to be free of functional defects in workmanship and materials when used normally for its intended surgical purpose. Any Snowden-Pencer device proving to be defective will be replaced or repaired at no charge. This device carries a lifetime warranty against manufacturer defects, 3 year warranty against wear.

**IEC Classification:** BF Applied Part

## Repair Service

Regardless of age, if any Snowden Pencer device needs service, return it to an authorized repair service center. For repairs outside the U.S., please contact your local distributor.

**Note:** All devices being returned for maintenance, repair, etc. must be cleaned and sterilized per these instructions prior to shipment.

## Contact Information

### CareFusion

1500 Waukegan Road  
McGaw Park, IL 60085 U.S.A.

800-323-9088

[www.carefusion.com](http://www.carefusion.com)

For Domestic inquiries email: [GMB-VMueller-Cust-Support@carefusion.com](mailto:GMB-VMueller-Cust-Support@carefusion.com)

For international inquiries email: [GMB-SIT-International-Team@carefusion.com](mailto:GMB-SIT-International-Team@carefusion.com)

## Other Resources

To learn more about sterilization practices and what is required of manufacturers and end users, visit [www.aami.org](http://www.aami.org), [www.aorn.org](http://www.aorn.org), or [www.iso.org](http://www.iso.org)

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