# Ventricular Drainage System

Sterile For Single Use Only

STERILE EO

CE 0123

NL850-500V

## Description

The Ventricular Drainage System, Model NL850-500V (see Figure 1), provides a sterile fluid path for the drainage of cerebrospinal fluid (CSF) from an intraventricular catheter to a collection bag. System components facilitate ventricular drainage, fluid injection, CSF sampling and intracranial pressure (ICP) monitoring when used in conjunction with Camino Ventricular Pressure Monitoring Kits.

The 4-way stopcock provides control of the CSF drainage system. By turning the stopcock to various positions, you will be controlling the direction of the fluid flow.

The collection bag is calibrated to provide approximate volumetric measurement of the fluid collected. Integral with the collection bag are a 50 ml rigid drip chamber, 30 cm H<sub>2</sub>O and 20 mmHg height scales used for positioning the collection bag. The proper placement is critical, as the height of the bag affects the ICP. A removable 1-way check valve is located on the bag side of the 4-way stopcock. Exiting the bottom of the collection bag is a tube with clamp and self-sealing injection/sampling port.

## Indications

The major indication for external drainage is the management of hydrocephalic shunt infections. If an internal shunt is not indicated, treatment of other cerebral conditions such as pre-operative or post-operative pressure

monitoring may also require external drainage of CSF to control increased intracranial pressure.

## Contraindications

This device is not designed, sold, or intended for use except as indicated.

## Instructions for Use

### Alignment

- The ICP is controlled by the height of the collection bag above the tip of the catheter; therefore, proper placement of the bag is critical.

- When using the cord for alignment, secure the cord of the drainage bag to an IV pole or other appropriate support. Utilize the height scale on the bag to maintain proper ICP. For example, to maintain an ICP of approximately 15 mmHg, the 15 mmHg mark on the height scale should be located at the approximate level of the catheter tip in the patient's ventricles.

- When using the pole clamp (NL850 PC), attach the pole clamp to the IV pole. Adjust the pole clamp height so the extended pointer is located approximately at the patient's ventricles (the top of the patient's ear) and secure in place. Then attach the bag to the pole clamp by sliding the burette down through the burette holder until it stops (see Figure 1). Utilize the height scale on the bag to maintain the proper ICP. For example, to maintain an ICP of approximately 15 mmHg, the 15 mmHg mark on the bag should be aligned with

the pointer, which is approximately at the same level as the patient's ventricles. After alignment, ensure the burette holder is locked in place.

**It is critical that neither the patient nor the collection bag be raised or lowered accidentally. Height changes should only be made by qualified personnel on the orders of the physician.**

### Set-Up

- Aseptic technique should be used during set-up and use of the ventricular drainage system.

- If desired, the 1-way check valve can be removed from the system by disconnecting the luer locks on both sides of the valve. Then, re-connect the tubing directly to the 4-way stopcock.

Note: The 1-way valve helps to control the direction of the fluid flow and reduce the chance of fluid backflow.

### Priming

- The external portion of the Ventricular Drainage System must be primed before connection to the catheter. Using aseptic technique, fill a syringe with sterile saline. Before puncturing, the port must be aseptically cleaned to prevent possible contamination.

- Close the 4-way stopcock to the extension tube and inject saline through the 4-way stopcock injection/sampling port. Prime the drainage tube until fluid

102575-702-02

enters the drip chamber. Then, close the 4-way stopcock to the drainage tube and prime the extension tube. Ensure that all air bubbles are eliminated and all connections are free from leaks.

#### Connection to Patient

- Introduce and position the ventricular catheter into the ventricle. Placement of these catheters may be accomplished through a variety of surgical techniques; therefore, the surgeon is best advised to use the method which his/her own practice and training dictate to be best for the patient.

- Remove the dust cap from the end of the extension tube and verify that the cap is removed from the ventricular catheter.

- Connect the extension tube to the ventricular drainage port. Open the 4-way stopcock and verify that the tubing clamps are open.

- The system is now ready for measurement of ICP and drainage of CSF, when indicated.

#### Draining CSF

- To drain CSF, turn the 4-way stopcock so the "off" arrow points to the sampling port. To measure small volumes, ensure that the one-way stopcock at the bottom of the 50 ml drip chamber is closed. When the drip chamber becomes full, empty it by opening the one-way stopcock. To calculate flow rate, assume an approximate conversion of 18 drops/ml in the drip chamber.

#### Injection and Sampling

- To inject fluid and to sample CSF, close the 4-way stopcock to the drainage tube. This allows the injection/sampling port to communicate with the catheter only. Before puncturing, the port must be aseptically cleaned to prevent possible contamination. After sampling or injection, re-open the 4-way stopcock to the drainage tube. Sampling may also be taken from the bottom of the collection bag.

#### Transport of Patient

- When patient requires transport, while connected to the Ventricular Drainage System, **temporarily** close the tube clamp on the drip chamber vent. This

will prevent wetting of the vent's filter element. Open the clamp after the patient is transported. The tube clamp must be open to allow proper operation of the vent.

#### Emptying of collection bag

- Before puncturing, the port must be aseptically cleaned to prevent possible contamination.

- When the collection bag becomes full, empty it by inserting a sterile 20G or smaller needle into the injection/sampling port at the bottom of the collection bag and withdraw fluid with a syringe.

#### Removal/Replacement of collection bag

The collection bag can be replaced in one of two locations:

##### Location 1

- Close the 4-way stopcock to the drainage tube and close the four tube clamps.

- Prepare the new collection bag by closing its tube clamps and removing the drainage tube, 4-way stopcock and extension tube.

- Disconnect and remove the full collection bag from the drainage tube.

- Secure the cord of the replacement collection bag to an IV pole or other appropriate support and connect the bag to the existing drainage tube.

- Align and prime system as described above.

- Resume drainage by opening the 4-way stopcock to the drainage tube and opening the tube clamps on the collection bag.

##### Location 2

- Close the 4-way stopcock to the drainage tube and close the four tube clamps.

- Prepare the new collection bag by closing its tube clamps and removing the extension tube and 4-way stopcock.

- Disconnect and remove the full collection bag and drainage tube.

- Secure the cord of the replacement collection bag to an IV pole or other appropriate support and connect the bag to the existing 4-way stopcock.

- Align and prime system as described above.

- Resume drainage by opening the 4-way stopcock to the drainage tube and opening all tube clamps.

#### Do Not Resterilize

The External CSF Drainage System components are disposable devices and cannot be resterilized.

#### Warnings

Hydrocephalic patients with cerebrospinal fluid drainage systems must be kept under close observation for signs and symptoms of changing intracranial pressure due to shunt failure. These signs and symptoms may vary from patient to patient. Increasing intracranial pressure is characterized by headache, vomiting, irritability, listlessness, drowsiness, other signs of deterioration of consciousness and nuchal rigidity. In the infant, increased scalp tensions at the anterior fontanelle and congestion of scalp veins will be noted.

Proper placement of the collection bag is critical, as its height controls the ICP. It is essential that neither the patient nor the collection bag be raised or lowered accidentally. Height changes should only be made by qualified personnel on the orders of the physician. CSF will not drain from the ventricles if the clear plastic tubing is occluded. If shunting is the primary function, reopen the tube clamps promptly after bag replacement and resume CSF drainage to maintain the intracranial pressure.

Physicians who infuse drugs into this device must be assured of the compatibility of the drugs, with polycarbonate, polyvinyl chloride, polyethylene, ABS and acrylic. Some change in the drug potency may occur when in contact with these materials.

#### Precautions

Prior to surgery, all prospective patients or their representatives should be informed of the possible complications associated with this product.

When the ICP monitoring port of the ventricular catheter is not in use, the protective cap must be in place to help prevent contamination. Aseptically clean all injection/sampling ports before using.

This device is intended for **SINGLE USE ONLY; DO NOT RESTERILIZE.**

**Integra NeuroSciences makes no claim for or representation as to the**

**performance characteristics of this product if it is used in conjunction with components of other manufacturers.**

**Complications**

Complications which may result from the use of this product include the risks associated with the medications and methods utilized in the surgical procedures, as well as the patient's degree of intolerance to any foreign

object implanted in the body. The principal complications associated with cerebrospinal fluid drainage are: catheter obstruction, functional failure of the system, infection or intracranial hypotension/hypertension.

Ventricular catheters may be obstructed by particulate matter such as blood clots, fibrin, or brain fragments. If not properly located in the lateral ventricle, the tube may become embedded in the


Figure 1. System Description

ventricular wall or choroid plexus. Less commonly, the catheter may be obstructed by the excessive reduction of ventricular size to slit-like proportions.

Functional failure of the drainage system due to separation of its component parts can result in serious complications. Ventricular catheters may migrate into the lateral ventricles.

Infection is a common and serious complication with drainage systems and is most frequently caused by skin contaminants. Septicemia, which occurs most frequently in debilitated infants, can result from infections anywhere in the body and may develop with few or no symptoms. Septicemia may also occur as a result of a wound infection. The presence of a foreign body (i.e., the catheter system) may trigger ventriculitis or a dormant meningitis.

Lesions developing from the breakdown of skin or tissue around the drainage system may also serve as foci of serious infections. In the event of an infection, replacement of the drainage system is indicated in addition to the appropriate antibiotic therapy.

Excessive lowering of intracranial pressure may result in complications, particularly in the infant. These include: subdural hematomas, markedly sunken fontanelles, overriding of cranial bones and conversion of a communicating to a noncommunicating hydrocephalus due to obstruction of the Aqueduct of Sylvius.

Failure of the drainage system may be evidenced by continuing symptoms of increased ICP.

## Product Information Disclosure

**Integra NeuroSciences has exercised reasonable care in the choice of materials and manufacture of this product. Integra NeuroSciences excludes all warranties, whether expressed or implied by operation of law or otherwise, including, but not limited to any implied warranties of merchantability or fitness for a particular purpose. Integra NeuroSciences shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from use of this product. Integra NeuroSciences neither assumes nor authorizes any other person to assume for it any other or additional liability or responsibility in connection with this device.**

## Product Order Information

All products can be ordered through your Integra NeuroSciences Neuro Specialist or customer service representative or by contacting :

Integra NeuroSciences  
311 Enterprise Drive  
Plainsboro, NJ 08536 USA  
Telephone: 1-800-654-2873  
Outside the US: 1-609-275-0500  
Fax: 609-275-5363

or

Integra NeuroSciences  
Newbury Road, Andover  
Hampshire SP10 4DR England  
Tel: +44(0) 1264-345-700  
Fax: +44 (0) 1264-332-113

**Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. Do not use if the package has been opened or damaged.**

### Symbols Used On Labeling

See instructions for use

Expiration date

Do not reuse after opening

Lot number

STERILE EO

Sterile unless package is opened or damaged. Method of sterilization-ethylene oxide.

0123

Product complies with requirements of directive 93/42/EEC for medical devices.

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© Copyright 2002 Integra LifeSciences. All rights reserved.