






Otto Bock®

QUALITY FOR LIFE

Ⓓ	Freistehbarren	<i>Bedienungsanleitung</i> 3
ⒼⒷ	Standing Frame	<i>Instructions for Use</i> 12
Ⓘ	Standing	<i>Istruzioni d'uso dello standing</i> 21
Ⓔ	Bipedestador	<i>Instrucciones de uso</i> 29
Ⓢ	Ståstativ	<i>Bruksanvisning</i> 38

Inhalt

Seite

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Contents

	Page
Preface	13
Intended Use	13
Range of Application	13
Liability	14
CE Conformity	14
Delivery / Preparing the Standing Frame for Use	14
Safety Tips	15
 Accessories	 16
Manual Strap Re-winder	16
Electric Strap Re-winder	16
Electric Strap Re-winder with Battery	16
Swivelling Pelvic Pads	17
Support Bar with Pelvic Pads	17
Upper Body Support.....	17
 Adjustability / Assembly Instructions	 17
Applying the User Belt	17
Approaching the Standing Device.....	17
Connecting the Straps with the User Belt	17
Getting Up / Standing / Sitting down	18
 Servicing / Maintenance / Cleaning	 18
 Technical Data	 19
Standing Frame for Adults.....	19
Standing Frame for Children	19
Swivelling Pelvic Pads	19
Upper Body Support.....	19
Electric Strap Re-winder	19
Electric Strap Re-winder with Battery	20
Table Battery Charger	20
Accessories.....	19

Preface

The Standing Frame is a quality product offering versatile use in everyday life.

Before using the Standing Frame **please read these Instructions for Use completely and thoroughly**. They contain important information for a safe and reliable use of this device. Please keep these Instructions for Use safe for later reference.

Delivery / Preparing the Standing Frame for Use

See page 14

Safety Tips

See page 15

Accessories

See page 16. Features additional components for this standing device, which may broaden the application possibilities and increase the user comfort.

Adjustability / Assembly Instructions

See page 17. Describes the possibilities of adjusting the standing device to your individual needs.

The design, as described in these Instructions for Use, is subject to technical alterations without notice.

If your Standing Frame is equipped with an electric strap re-winder make sure that the voltage of your mains supply corresponds with the indications on the type plate.

When adjusting the table height of the Standing Frame make sure that the emergency/safety cut-out device has sufficient distance to the handrail's cross tube to ensure proper function of the automatic emergency/safety cut-out mechanism (it might be necessary to adjust the electric strap re-winder higher).

Intended Use

The Standing Frame is designed solely for people who are unable to walk or who have a walking impediment. The Standing Frame can be used by the patient without assistance or with the assistance of an attendant.

The Standing Frame is intended to be used as part of a standing therapy programme. A level of supervision appropriate to the user's abilities should be applied at all times. **DO NOT EXCEED HEIGHT AND WEIGHT LIMITS SET OUT IN THESE INSTRUCTIONS**. The Standing Frame is not intended to position a user at any angle other than vertical and should not be positioned at any reclined angle. Do not use the Standing Frame on any inclined surface.

Range of Application

The Standing Frame is an important part of many therapy regimes. The act of standing helps to stretch muscles and tendons to prevent contractures, improves circulation, extends the abdominal cavity for respiration and digestive processes to function easily and helps in the evacuation of the upper urethral tract. It allows safe, stable, and independent standing.

Although there are several configurations of the Standing Frame the basic principle behind their use is the same. By blocking the knees against the knee pads on the Standing Frame, positioning the feet on the footplate and lifting the pelvis, the user can be safely raised from a sitting position into the standing position in stages depending on the user ability.

The Standing Frame is a free-standing device and can be used on any flat, solid ground without additional attachment to a wall.

The transport wheels allow easier movability (change of place living room, terrace, etc.).

Standing Frame for Adults

The device can be used by persons having a body height between approx. 140 cm and 185 cm (for persons with a body height of up to approx. 200 cm an extension of the table and handrail bars is required).

- Adjustable handrail bars within a range of approx. 80 cm and 105 cm.
- Gas cylinder assisted height adjustment of the table board within a range of approx. 90 cm and 130 cm.
- The automatic strap device can be shifted by approx. 20 cm independently from the table board.
- Vertically adaptable knee pads by approx. 20 cm.
- Individual fixation of the feet in four possible positions using hook & loop straps and detachable heel restraints.

To provide individual fitting the following items should be considered:

- Height and weight of the body (max. user weight 120 kg).

- Physical and psychological constitutions.
- Home environment.

Standing Frame for Children

The device can be used by persons having a body height between approx. 100 cm and 140 cm.

- Adjustable handrail bars within a range of approx. 65 cm and 85 cm.
- Height adjustable table board within a range of approx. 85 cm and 105 cm.
- The automatic strap device can be shifted by approx. 20 cm independently from the table board.
- Vertically adaptable knee pads by approx. 12 cm.
- Individual fixation of the feet in four possible positions using hook & loop straps and detachable heel restraints.

To provide individual fitting the following items should be considered:

- Height and weight of the body (max. user weight 70 kg).
- Physical and psychological constitutions.
- Home environment.

Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

The manufacturer is not responsible for damages caused by components and spare parts not approved by the manufacturer. Repairs must be carried out exclusively by authorised dealers or by the manufacturer.

CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

Delivery / Preparing the Standing Frame for Use

To remove the transport packaging turn the carton box by 90° (fig. 1) and put it on its longitudinal side so that the transport wheels are at the bottom/floor. Pull the Standing Frame out of the carton box (fig. 2). Pull up to the standing position stepping on the tube between the transport wheels. For moving the Standing Frame (fig. 3), also step on the tube between the transport wheels and tilt back the Standing Frame then pull/wheel to required area for use.

Components of the Standing Frame (fig. 4):

- 1 Switch lever (for manipulation while standing)
- 2 Table board
- 3 Switch lever (for manipulation while sitting)
- 4 Clamping levers
- 5 Knee pads
- 6 Transport wheels
- 7 Toe hook & loop straps
- 8 Rubber bumpers
- 9 Foot restraints
- 10 Footplate
- 11 Cross tube
- 12 Manual strap re-winder
- 13 Handrail bars

In general, your Standing Frame will be delivered completely assembled. If you have ordered your standing device with accessory components, a few simple steps will suffice to prepare it for use.

- ☞ If your Standing Frame is equipped with a manual strap re-winder or swivelling pelvic pads, you need to insert the long switch lever (fig. 4, item 1) in the bush provided and tighten the screw. (The black ball of the switch lever must be at the top).

- ☞ If your Standing Frame is equipped with an electric strap re-winder, place your standing device near to a wall socket and connect the mains plug (fig. 5, item 1) to the wall socket. Once the strap re-winder has been connected to the mains supply the indicator lamp lights up and the device is ready for use.
- ☞ Before using the Standing Frame, it must be adapted to the individual user.
Especially for Standing Frames with electric strap re-winder (mains supply or battery), correct adjustment of the table height and the switch-off button of the emergency/safety cut-out device is an important safety function.
To adjust the table height, release the clamping levers (fig. 4, item 4) and adjust the table so that the table board pad is located (in accordance with the individual requirements) between the sternum and the upper edge of the pelvis of the user when in the standing position. The table board pad must never be set below the upper edge of the pelvis! After adjustment retighten the clamping levers; do not overtighten clamping levers (see section servicing/maintenance/cleaning).
To adjust the switch-off button (fig. 5, item 2) of the emergency/safety cut-out mechanism, loosen the lock nut (fig. 5, item 3) with a 10 mm spanner and rotate the switch-off button on the thread of the switch-off bar. The switch-off button is set correctly if the motor is not yet switched off when the desired standing position is achieved but automatic emergency/safety cut-out is about to occur immediately. After any adjustment has been carried out make sure to re-tighten the lock nut (fig. 5, item 3).

Caution must be taken if the Standing Frame is to be used by several users!

Important

After an emergency/safety cut-out has occurred, both the re-wind and unwind functions are shut off. The switch-off bar must be brought into its original position (with no pressure being applied to the button).

- ☞ If your Standing Frame is equipped with an electric strap re-winder with battery: The integrated battery has already been pre-charged. However, it is recommended to charge the battery once again before using the device. Please note that the batteries will develop their full capacity only after several charging processes have been carried out.
The regulated battery charger for the electric strap re-winder with battery works fully automatically. Connect the secondary line to the strap's charging receptacle (observing correct pole connection). Plug the charger into the mains electricity supply; once connected, the charging process will start and the red LED lights up. When the red LED goes out (with the mains electric supply still connected), the battery is fully charged, disconnect the electric mains supply.
- ☞ If your Standing Frame is equipped with a thoracic pad: For transport reasons, the thoracic pad is **not** mounted in the correct way. Release the two clamping levers of the support bar with pads and slide the support bar backwards until the right side slides off the right handrail.
Turn the support bar in counter-clockwise direction approx. 90°. Close the left clamping lever of the support bar with pads. Pull the indexing plungers of the thoracic pad height adjustment and remove the thoracic pad. Now insert the pad from above into the support bar's holding bracket and fasten the clamping levers, **do not overtighten** (see section servicing/maintenance/cleaning).
- ☞ Please **ensure** when mounting the pad that the indexing plungers of the adjustment device have popped into position.

Safety Tips

- ➔ Check the belt is positioned correctly on the user before standing. The lip along the bottom edge is to fit under the user bottom. Then attached to the user belt is a hook & loop closure which fastens around the user and two car seat belt type buckles which join the user belt to the lifting mechanism (these straps should run slightly upwards when in standing position).
- ➔ To avoid increased pressure on the knees when standing up, make sure to have a sitting position in which the knee angle is at least 90°.
- ➔ If your Standing Frame is equipped with 4 wheels (standard with the CHILDREN'S Standing Frame), apply all wheel locks before use!
Do not move the Standing Frame while a person is using it!
- ➔ Speak with your physician or therapist regarding your standing training programme, especially before using the Standing Frame for the first time.
- ➔ **Pay attention** to safety! Use the Standing Frame only in the presence of an attendant.
- ➔ Standing Frame with electric strap re-winder: We recommend you inspect the emergency/safety cut-out mechanism for correct function on a weekly basis! To do this, press the switch-off button during operation of the strap re-

winder; this will cause the electric re-winder to stop. If this fails to happen there is a malfunction, and an authorised specialist must check the device before it is used again.

→ The air hose of the suck and blow control must always be laid and placed in a way that any unwanted functions are not possible (by the user, attendant or a third person). Malfunctions can occur if the air hose is kinked or pinched. If the suck and blow control is not required, the air hose can be pulled off the connection piece.

→ In case of a power failure (with electric strap re-winder)

If there is a power failure while using the electric strap re-winder, at least one attendant is required to release the user from the Standing Frame in order to help the user back into the wheelchair. The following description starts from the assumption that the attendant is right-handed. For a left-handed attendant the instructions must be carried out from the opposite side.

Swing the support bar out of the way (see section support bar with pelvic pads). Position the wheelchair to the left of the Standing Frame as close as possible (for use later) but allowing enough space for the attendant to stand behind the user.

The attendant stands behind the user.

Open the hook & loop closure.

With the right arm, grasp under the user's axilla and around his thorax.

With the left hand open the strap closure of the user belt - it may be necessary to relieve the pressure on the user belt by pushing the user slightly forward.

With the left hand grasp under the user's buttocks to support him when gliding back into the wheelchair.

Accessories

The Standing Frame is designed as a modular assembly system, which means that certain authorised accessories may be fixed to your standing device.

The following options and accessories are available for your standing device.

Manual Strap Re-winder

Important: The manual strap re-winder has been designed for a user who has good shoulder and arm strength enabling them to stand up and sit down on their own.

The strap tensioning system takes up the slack in the sling as the user lifts themselves to standing. For standing up, bring the switch lever (3) into the "Release" position, thus allowing to unwind the straps. Connect the straps with the user belt. Bring the switch lever into "Standing" position. Now, the mechanical system of the manual strap re-winder allows only re-winding of the straps.

This system does not actually assist the lift but makes self-lifting safer, as the sling will lock in any position if it is put under tension, allowing the user to rest or re-position their arms during the lift.

For sitting down again, bring the switch lever (1) into the "Release" position. The locked straps will be released only after the tension is relieved by a small forward movement of the hips.

Caution: When disconnecting the straps from the belt ensure to hold on to the straps as they re-wind. Do not let go as this may cause damage to the mechanism and cause a malfunction!

The manually operated Standing Frame is suitable for lower lesion spinal injuries and paraplegics with good shoulder and arm strength who are able to push on the support bars.

Electric Strap Re-winder, Type KR213167 (fig. 5)

The electric strap re-winder has been designed as an accessory to the Standing Frame for users with very little arm strength.

This accessory can be operated with a hand control switch (item 4) or with a suck and blow control (item 5). By pressing the membrane key button on the hand switch or slight blowing or sucking on the air hose, the electric motor will re-wind or unwind the straps.

Due to the modular design it is possible to carry out a retrofit to an existing Standing Frame with an electric strap re-winder later.

Electric Strap Re-winder with Battery, Type K213167

The electric strap re-winder with battery has been designed as an accessory to the Standing Frame for users with very little arm strength and for use of the Standing Frame independently from a mains supply.

This accessory can be operated with a hand control switch (fig. 5, item 4) or with a suck and blow control (fig. 5, item 5). By pressing the membrane key button on the hand switch or slight blowing or sucking on the air hose, the electric motor will re-wind or unwind the straps.

Due to the modular design there is no problem to retrofit a Standing Frame with an electric strap re-winder with battery or vice versa with a manual strap re-winder later.

☞ For changing from manual to electric strap re-winder, please proceed as follows:

Remove the two securing screws (fig. 6, item 1) on the two guide bars. Stand behind the device and support the strap re-winder with your forearms to avoid uncontrolled gliding off when pulling the indexing plungers (fig. 7, item 1). Now let the strap re-winder glide downwards on the guide bars and then remove from the frame.

☞ **Important:** The Standing Frame can only be used with one type of strap re-winder. A combination of manual and electric strap re-winder is not allowed! The strap re-winder may only be exchanged when the device is not in use by any person and when it is disconnected from the mains supply.

Swivelling Pelvic Pads

The Standing Frames KR 020034 and KR 020033 are equipped with swivelling pelvic pads instead of the electric, battery, manual strap re-winder and user belt. These swivelling pads have been designed primarily for users who are able to get into the standing position in one single movement.

By swivelling in the left and right pelvic pads towards the body (swivelling one pad will move both pads simultaneously) the pelvic pads lock automatically into position, the user will now be stabilised in the pelvic area and will be safely supported in the standing position.

☞ **Attention:** The swivelling pelvic pads cannot be combined with the strap re-winders!

Support Bar with Pelvic Pads (fig. 8, item 1)

The support bar with pads provides more positive pelvic control for the user and must be used in conjunction with the lifting strap.

Before using the Standing Frame, release the two clamping levers (fig. 8, item 3) of the support bar and slide the support bar backwards on the handrails. Please note the left-hand rail is slightly longer, this is to allow the support bar to slide off the right rail and swing downwards out of the way for transferring. A small securing screw at the end of the left handrail bar prevents the support bar from pulling off completely.

After the user has stood up, an attendant slides the support bar back on the right handrail bar and fixes it into position with the clamping levers; **do not overtighten** (see section servicing/maintenance/cleaning).

The pelvic pads can be individually adjusted laterally to centralise the pelvis. The user position is held with indexing plungers. Pull the plungers out and slide the pelvic pads along the bar to the desired position. Release the plungers into the nearest set of holes and ensure they have popped into position. The aim is to position the pads symmetrically.

Thoracic Pad/Thoracic Lateral Pads and Head Support (fig. 8 item 2)

A range of posteriorly (behind user) positioned supports are available to provide additional postural control. The thoracic pad prevents spinal extension, the thoracic lateral pads provide lateral (sideways) support of the chest and the head support provides limited head control.

To assemble these parts, the support bar (fig. 8 item 1) with two lateral pads (art. no. RR02.010.00.000) is required. After the user has stood up, an attendant slides the support bar with the two lateral pads from behind over the handrail bars and fixes it with the clamping levers; **do not overtighten** (see section servicing/maintenance/cleaning). Now individual adaptation of the thoracic pad is possible.

Adjustability / Assembly Instructions

Applying the User Belt (fig. 9)

Before approaching the standing device with your wheelchair apply the user belt. If you have a standing device with castors, first apply the locks! The user belt should be applied so that its lower edge is clamped between the user's buttocks and the wheelchair seat. Close the abdominal strap.

Approaching the Standing Device (fig. 10)

Laterally swing away the wheelchair's legrests or footplate and place your feet into the hook & loop straps for the toes. The foot restraints can be removed for this. Approach the standing device with the wheelchair until the front wheels contact the footplate of the Standing Frame in the cutouts provided for this purpose. Apply the wheelchair brakes.

Connecting the Straps with the User Belt (fig. 11)

Bring the switch lever (item 1) into the position "Release". It is now possible to unwind the straps against a light resistance. Connect the strap closure tongues on the left and right with the corresponding strap closures of the user belt. Bring the switch lever (item 1) into the position "Standing". Now it is no longer possible to unwind the straps.

Getting Up / Standing / Sitting Down (fig. 12-14)

After grasping the handrails, tablesides or table board, the user can pull up or push themselves up. The manual strap re-winder allows the user to stand up in stages, since the manual strap re-winder locks and holds the user in any position when in the functional position "Standing". If the user has not the necessary strength in their arms, an attendant should support the user.

- ☞ Place the right hand under the user's right axilla, left hand under their buttocks.
- ☞ The strap re-winder allows the user to correct his grasp already during standing up, since their body weight is held in any position. Caution: Proper positioning of the belt is required, the lip along the bottom edge is to fit under the user bottom. Then attached to the user belt is a hook & loop closure which fastens around the user and two car seat belt type buckles which join the user belt to the lifting mechanism (these straps should run slightly upwards when in standing position).
- ☞ To get back from the standing position into the wheelchair, first turn the switch lever to "Release" (fig. 11, item 2). In spite of the "Release" position, the user belt is still locked. To unlock the belt, the pressure on both straps must be relieved. Caution: The user is no longer held by the strap system. With the user holding on to the table or handrails guide them to a sitting position within the waiting wheelchair.

Servicing / Maintenance / Cleaning

Your Standing Frame has been provided with the CE sign.

The manufacturer herewith guarantees that this medical product as a whole conforms to the requirements of 93/42/EEC Guideline. Should you have problems with your Standing Frame, please contact your healthcare provider.

- ☞ **Keep castors free from dirt, hair and other loose fibres.**
- ☞ Do not attempt to move the Standing Frame without releasing the brakes on all four castors.
- ☞ Clean pads promptly after soiling. Avoid heavy brushing of pads with stiff bristles. Clean upholstery with a damp cloth and mild detergent.
- ☞ Do not clean any part of the Standing Frame with bleach or solvent cleaners.
- ☞ Clean metal work with a slightly damp cloth.
- ☞ Clean woodwork with standard domestic furniture polish.
- ☞ Always store the Standing Frame away from excessive moisture as this may lead to rust and corrosion.
- ☞ Caution must be taken through use, transport and storage. Do not to drop or knock the Standing Frame.
- ☞ Do not store other equipment on top of the Standing Frame or its accessories as this may cause damage.
- ☞ Do not store the Standing Frame where it may block access to doors, walkways or fire exits.
- ☞ **If the webbing straps are heavily frayed or the lifting sling shows signs of tearing, fraying or other damage,** take the Standing Frame out of service immediately and contact the supplier.
- ☞ **Regularly check the hand levers** which control tray and handrail height for signs of wear, such as the thread looking rounded and highly polished, and replace them as necessary to ensure that they lock up securely.
- ☞ **Do not overtighten clamping levers** as this may cause damage to the threads and metal work.

For reasons of technical safety we recommend that you have the product inspected at least once a year by an authorised specialist. The device should be checked for wear and tear, damage and functional safety.

Technical Data of the Standing Frame for Adults / Children as well as of Accessories

1. Standing Frame for Adults

Width	77.5 cm
Length	107 cm
Height	90 cm - 130 cm
Weight	45 kg

2. Standing Frame for Children

Width	70 cm
Length	85 cm
Height	85 cm - 105 cm
Weight	43 kg

3. Swivelling Pelvic Pads

Width	64 cm - 102 cm
Length	95 cm
Height	14 cm
Weight	14 kg

4. Upper Body Support

Width	65 cm
Height	100 cm
Depth	35 cm
Adjustment range	32 cm - 58 cm
Weight with support bar	12 kg

5. Electric Strap Re-winder

Rated voltage	230 V
Rated power	160 W
Type	B
Protective system	II
Operating time 50%	SD=10 min
Maximum user weight	120 kg
Length	33 cm
Width	43.5 cm
Height	9 cm
Weight	11.5 kg

Recommended Accessories for the Electric Strap Re-winder

Article description	Article number
User belt small with lip	RR02.039.00.000
User belt medium with lip	RR02.040.00.000
User belt large with lip	RR02.041.00.000

6. Electric Strap Re-winder with Battery

Rated voltage	24 v
Rated power	160 W
Type	B
Protective system	–
Operating time 50%	SD = 10 min
Maximum user weight	120 kg
Length	33 cm
Width	43.5 cm
Height	9 cm
Weight	8 kg (without battery)

7. Table Battery Charger for the Electric Strap Re-winder with Battery

Input voltage	230 V + 6% cable with Euro-plug
Output voltage	27.6 V
Output current	0.8 A
Secondary connection	approx. 2 m stranded wire, bared and tinned, positive pole marked
Indicator	red LED, which goes out after transition to conservation charging
Temperature range	-20° C to +40° C
Degree of interference suppression	"N" according to DIN VDE 88 75
Protective system	IP 20
Length	1.43 cm
Width	7.8 cm
Height	5.8 cm

Regulated table battery charger for maintenance-free lead batteries (model DL 24-0,8)

The linearly regulated battery charger serves the charging of the maintenance-free lead batteries in the electric strap re-winder.

The output voltage is regulated to 2.3 V pile $\pm 1\%$ and is provided with a current limitation. In case of overcharge or short circuit the device is electronically regulated to a lower value or is thermally shut off. The charging process is indicated by an LED which goes out after transition to conservation charging.

The battery charger is located in a plastic table housing in accordance with the regulations being valid on international level.

Note: The values indicated are measurements which have been theoretically determined. We reserve the right to technical alterations and deviations by 5%.

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Otto Bock HealthCare LP
Two Carlson Parkway North, Suite 100
U.S.A.-Minneapolis, Minnesota 55447
Phone +1 800 328 4058 · Fax +1 800 962 2549
e-mail: usa.customerservice@ottobock.com
www.ottobockus.com

Asia/Pacific

Otto Bock Australia Pty. Ltd.
Suite 1.01, Century Corporate Centre
62 Norwest Boulevard · Norwest Business Park
AUS-Baulkham Hills NSW 2153
Tel. +61 2 88182800 · Fax +61 2 88182898
e-mail: healthcare@ottobock.com.au
www.ottobock.com.au

Beijing Otto Bock Orthopaedic Industries Co. Ltd.
HengXiang Tower · No.15 Tuanjihu South Road
Chaoyang District · Beijing 100026 · P.R.China
Tel. +86 10 85986880 · Fax +86 10 85980040
e-mail: china@ottobock.com.cn · www.ottobock.com.cn

Otto Bock Asia Pacific Ltd.
Suite 3218, 32/F., Sun Hung Kai Centre
30 Harbour Road, Wanchai · Hong Kong
Tel No. +852 2598 9772 · Fax No. +852 2598 7886
e-mail: info@ottobock.com.hk

Otto Bock HealthCare India Pvt. Ltd.
Behind Fairlawn Housing Society
Sion Trombay Road
Chembur · IND-Mumbai 400 071
Tel. +91 22 2520 1268 · Fax +91 22 2520 1267
e-mail: information@indiaottobock.com
www.ottobockindia.com

Otto Bock Japan K. K.
Yokogawa Building 8F · 4-4-44 Shibaura, Minato-ku
J-Tokyo 108-0023
Tel. +81 3 3798-2111 · Fax +81 3 3798-2112
e-mail: ottobock@ottobock.co.jp
www.ottobock.co.jp

Otto Bock Korea HealthCare Inc.
Beakyoung B/D 2FL · 37-22, Samsung-dong
Gangnam-gu · ROK-Seoul 135-090
Tel. +82 2 577-3831 · Fax +82 2 577-3828
e-mail: info@ottobockkorea.com
www.ottobockkorea.com

Otto Bock South East Asia Co. Ltd.
1741 Phaholyothin Road,
Kwaeng Chatuchark, Khet Chatuchark,
T-Bangkok 10900
Tel. +66 2 930 3030 · Fax +66 2 930 3311
e-mail: obsea@ottobock.co.th
www.ottobock.co.th

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Otto Bock Mobility Solutions GmbH has been certified by the German Society for the Certification of Quality Assurance Systems (DQS) in accordance with DIN EN ISO 9001 standard, reg. no. 779 (management system)

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Otto Bock HealthCare GmbH
Otto Bock Manufacturing Königsee GmbH