HAEMATOKRIT 210

(DE) Bedienungsanleitung	5
(EN) Operating Instructions	17
(FR) Mode d'emploi	29
(IT) Istruzioni per l'uso	42


Fig. 1


Fig. 2

**EG-Konformitätserklärung
EC Conformity Declaration
Déclaration de conformité CE
Dichiarazione di conformità alle norme CEE**

Andreas Hettich GmbH & Co. KG • Föhrenstraße 12 • D-78532 Tuttlingen • Germany

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Contents

1	Use according to specification	18
2	Residual risks	18
3	Technical specifications	18
4	Notes on safety	19
5	Symbol meanings	20
6	Delivery checklist	20
7	Initial operation	20
8	Opening and closing the lid.....	21
8.1	Opening the lid	21
8.2	Closing the lid	21
9	Control and display elements.....	21
9.1	Symbols on the control panel	21
9.2	Keys and setting options.....	21
10	Setting the brake step	22
11	Centrifugation.....	22
11.1	Centrifugation with preselected time	22
11.2	Continuous operation.....	22
11.3	Short-time centrifugation	23
12	Relative centrifugal force (RCF)	23
13	Centrifugation of materials with higher density.....	23
14	Haematocrit measurement	23
14.1	Haematocrit capillaries - Overview.....	23
14.2	Filling of the haematocrit capillaries	24
14.3	Sealing and centrifugation of the haematocrit capillaries	24
14.3.1	Standard capillary tubes, Order No.: 2074, and mylar-reinforced capillary tubes, Order No.: 1072....	24
14.3.2	Self-sealing capillary tubes (Order No. 1071)	24
14.4	Loading the rotor	24
14.5	Determination of the centrifugation time.....	25
14.6	Reading the haematocrit value	25
15	Emergency release	25
16	Maintenance and servicing.....	25
16.1	Cleaning the centrifuge	25
16.2	Cleaning and disinfecting the haematocrit rotor	26
16.3	Procedure in case of broken haematocrit capillaries.....	26
16.4	Rotor removal / installation.....	26
17	Faults	27
18	Change mains input fuse.....	28
19	Acceptance of the centrifuges for repair.....	28
20	Disposal	28
21	Anhang / Appendix	55
21.1	Rotoren und Zubehör / Rotors and accessories.....	55

1 Use according to specification

The machine presented here is a medical product (laboratory centrifuge) according to the IVD guideline 98/79/EG. The centrifuge is used to separate substances or substance mixtures with a density of max. 1.2 kg/dm³. This also includes substances and substance mixtures of human origin. The centrifuge is only intended to be used for this purpose. A different use or application over and above this is deemed not in accordance with the specifications. The company Andreas Hettich GmbH & Co. KG undertakes no liability for damages resulting therefrom.

Belonging to the application according to specification is also the observance of all references contained in the Instruction Manual and compliance with the inspection and maintenance works.

2 Residual risks

The machine is constructed according to the state of the art and the recognized technical safety regulations. Improper use and handling can result in dangers to life and limb of the user or third parties and impairments to the machine or to other material assets. The machine is only to be used for the specified applications and only in an impeccable technical safety condition.

Disturbances that can interfere with the safety are to be immediately rectified.

3 Technical specifications

Manufacturer	Andreas Hettich GmbH & Co. KG D-78532 Tuttlingen	
Model	HAEMATOKRIT 210	
Type	2104	2104-01
Mains voltage ($\pm 10\%$)	208 - 240 V 1~	100 - 127 V 1~
Mains frequency	50 - 60 Hz	50 - 60 Hz
Connected load	250 VA	220 VA
Current consumption	1.35 A	2.2 A
Max. capacity	24 x 2.2 ml	
Allowed density	1.2 kg/dm ³	
Speed (RPM)	13000	
Force (RCF)	16060	
Kinetic energy	850 Nm	
Obligatory inspection (BGR 261)	no	
Ambient conditions (EN 61010-1)	Indoors only Up to 2000 m above sea level 2°C to 40°C Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C. II 2	
Device protection class	I	
Not suitable for use in explosion-endangered areas.		
EMC	EN 55011, Group 1, Class B EN 61000-3-2 EN 61000-3-3 EN 61000-6-1	
		FCC Class B ----
Noise level (dependent on rotor)	≤ 58 dB(A)	
Dimensions	275 mm 330 mm 247 mm	
- Width	275 mm	
- Depth	330 mm	
- Height	247 mm	
Weight	approx. 10 kg	

4 Notes on safety

No claim under guarantee will be considered by the manufacturer unless the above instructions have been adhered to.

- Before the initial operation of your centrifuge you should read and pay attention to the operating instructions.
Only personnel that has read and understood the operating instructions are allowed to operate the device.
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner.
These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate.
 - However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge should be installed on a good, stable base.
- Before using the centrifuge absolutely check the rotor for firm placement.
- When the centrifuge is running, according to IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either heat up for at least 3 hours in the warm room before being connected to the mains, or run hot for 30 minutes in the cold room.
- Only the rotors and accessories approved by the manufacturer for this device may be used (see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").
- The centrifuge rotor may only be loaded in accordance with the chapter "Loading the rotor".
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge may not be operated in explosion-endangered areas.
- The centrifuge must not be used with:
 - inflammable or explosive materials
 - materials that react with one another producing a lot of energy.
- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures.
For hazardous substances centrifuge containers with special screw caps must strictly be used. In addition to the screw cap centrifuge containers, for materials in hazard category 3 and 4 a biosafety system must be used (see the World Health Organisation's "Laboratory Biosafety Manual").
In a biosafety system, droplets and aerosols are prevented from escaping by a bioseal (packing ring).
If the hanger of a biosafety system is used without the lid, the packing ring must be removed from the hanger in order to prevent the packing ring from being damaged during the centrifugation run. Damaged packing rings must not be used to seal the biosafety system.
Without the use of a biosafety system the centrifuge is not microbiologically sealed in the sense of the EN 61010-2-020 standard.
For further details of available biosafety systems see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". If in doubt, you should obtain relevant information from the manufacturer.
- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.
- Rotors, suspensions and accessories that possess traces of corrosion or mechanical damage or if their term of use has expired may not be used any longer.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Andreas Hettich GmbH & Co. KG company are allowed to be utilised.
- The following safety regulations apply:
IEC 61010-1 and IEC 61010-2-020 as well as their national deviations.
- The safe operation and reliability of the centrifuge can only be guaranteed if:
 - the centrifuge is operated in accordance with the operating instructions,
 - the electrical installation on the site where the centrifuge is installed conforms to the demands of IEC stipulations,
 - prescribed tests to BGV A1, BGR 261 are carried out by an expert.

5 Symbol meanings

Symbol on the machine:

Attention, general hazard area.

Before using the centrifuge implicitly read the operating instructions and pay attention to the safety relevant references!

Symbol in this document:

Attention, general hazard area.

This symbol refers to safety relevant warnings and indicates possibly dangerous situations.

The non-adherence to these warnings can lead to material damage and injury to personal.

Symbol in this document:

This symbol refers to important circumstances.

Symbol on the machine and in this document:

Symbol for the separate collection of electric and electronic devices according to the guideline 2002/96/EG (WEEE). The device belongs to Group 8 (medical devices).

Applies in the countries of the European Union, as well as in Norway and Switzerland.

6 Delivery checklist

The following items and accessories are delivered with the centrifuge:

- 1 Connecting cable
- 2 Fuses
- 1 Notes on moving the equipment safely
- 1 Operating instructions
- 1 Hexagonal pin-type spanner

The rotor(s) and associated accessories are included in the delivery in the quantity.

7 Initial operation

- Remove the transportation safety device from the bottom of the housing, see sheet "Transportation safety device"
- Position the centrifuge in a stable and level manner in a suitable place. During set-up, the required safety margin of 300 mm around the centrifuge is to be kept according to IEC 61010-2-020.

When the centrifuge is running, according to IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

Substances can be discharged through the ventilation opening on the rear side of the centrifuge. The device is to be placed in such a way that the air stream is not directed towards people.

- Do not place any object in front of the ventiduct.
Keep a ventilation area of 300 mm around the ventiduct.
- Check whether the mains voltage tallies with the statement on the type plate.
- Connect the centrifuge with the connection cable to a standard mains socket. For connection ratings refer to Chapter "Technical specifications".
- Turn on the mains switch. Switch position "I".
The last used centrifuge data will be displayed.
- Open the lid.
- Remove the transportation safety device in the centrifuge, see sheet "Transportation safety device".

8 Opening and closing the lid

8.1 Opening the lid

The lid can only be opened when the centrifuge is switched on and the rotor is at rest. If it cannot be opened under these circumstances, see the section on "Emergency release".

- Rotate the turning handle on the front panel to the left. The symbol "L" (lid open) illuminates in the rotation indicator ↗.
- Open the lid.

8.2 Closing the lid

Do not bang the lid shut.

- Place the lid and lightly press down the front edge of the lid.
- Turn the handle on the front panel to the right . The symbol "—" (lid closed) illuminates in the rotation indicator ↗.

9 Control and display elements

See figure on page 2.

Fig. 2: Display and control panel

9.1 Symbols on the control panel

Rotation indicator. The rotation indicator lights up and rotates anticlockwise while the rotor is turning.

When the rotor is stationary, the status of the lid is displayed by symbols in the rotation indicator:

Symbol L : Lid open
Symbol — : Lid closed

Operator errors and occurring faults are indicated on the display (see Chapter "Faults").

9.2 Keys and setting options

RPM x 100

- Speed
A numeric value of 500 RPM up to the maximum rotor speed can be set. For maximum rotor speed, see chapter "Anhang/Appendix, Rotor und Zubehör/Rotor and accessories". Preset in steps of 100 (RPM = displayed value x 100).
If the key ▲ or ▼ is kept pressed, the value changes with increasing speed.
- Display the brake step.

t

- Running time
- Preset from 1 - 99 minutes, in 1 minute steps
- Continuous operation "--"
- Braking steps 0 or 1. Step 1 = short run-down time, Step 0 = long run-down time.
If the key ▲ or ▼ is kept pressed, the value changes with increasing speed.

START

- Start centrifugation run.


STOP

- End centrifugation run.
The rotor runs down with the preselected brake step.
- Save the brake step.

- Short-time centrifugation.
The centrifugation run occurs while the key
- Display the brake step.

10 Setting the brake step

- Switch off the mains switch.
 - Keep the key ▲ beneath the speed indicator and the key ▾ pressed simultaneously.
 - Switch on the mains switch and release the keys again.
- The speed indicator shows the machine version and the time indicator shows the set brake step: e.g.:

If the machine version and brake step are not displayed, press the ▲ key under the speed indicator until they are displayed.

The machine version is set by the manufacturer and cannot be changed.

- Set the desired brake step with the keys ▲ ▾ beneath the time indicator.
Step 1 = short run-down time, Step 0 = long run-down time.
For run-down times, see chapter "Anhang/Appendix, Rotor und Zubehör/Rotor and accessories".
- Press the key **STOP** to save the setting.

11 Centrifugation

When the centrifuge is running, according to IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

If the permissible weight difference is exceeded within the rotor loading, the drive switches off during the run-up time, and error -3- is displayed (see chapter "Faults").

The centrifugation run can be interrupted at any time by pressing the key **STOP**.

The time and speed can be changed during the centrifugation run, with the keys ▲ ▾.

If the key ▲ or ▾ is kept pressed, the value changes with increasing speed.

After a centrifugation run, the display flashes until the cover is opened or a key is pressed.

If the symbol "—" (lid closed) and "L" (lid open) flashes alternately in the rotation indicator ↗, operation of the centrifuge can only be continued after opening the lid.

- Switch on the mains switch (switch position "I").
- Load the rotor and close the centrifuge cover.

11.1 Centrifugation with preselected time

- Set the desired speed with the keys ▲ ▾ beneath the speed indicator.
- Set the desired time with the keys ▲ ▾ beneath the time indicator.
- Press the key **START**. The rotation indicator ↗ appears while the rotor is turning.

The time is displayed in minutes. The last minute is counted down in seconds.

When the time is displayed in minutes, a point flashes next to the number.

- After expiry of the time or if the centrifugation run is interrupted by pressing the key **STOP**, the rotor runs down with the set brake step.

During the centrifugation run the rotor speed and the remaining time are displayed.

11.2 Continuous operation

- Set the desired speed with the keys ▲ ▾ beneath the speed indicator.
- Set the time to zero with the key ▾ beneath the time indicator. "--" is displayed.
- Press the key **START**. The rotation indicator ↗ appears while the rotor is turning. The time count starts from 0.

The first minute is counted up in seconds, and then the time is displayed in minutes.

When the time is displayed in minutes, a point flashes next to the number.

- Press the key **STOP** to end the centrifugation run. The rotor runs down with the set brake step.

During the centrifugation run the rotor speed and the expired time are displayed.

11.3 Short-time centrifugation

- Set the desired speed with the keys ▲▼ beneath the speed indicator.
- Keep the key ⌂ pressed. The rotation indicator ↗ appears while the rotor is turning. The time count starts from 0.

The first minute is counted up in seconds, and then the time is displayed in minutes.
When the time is displayed in minutes, a point flashes next to the number.

- Release the key ⌂ again to end the centrifugation run. The rotor runs down with the set brake step.

During the centrifugation run the rotor speed and the expired time are displayed.

12 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$\text{RCF} = \left(\frac{\text{RPM}}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad \text{RPM} = \sqrt{\frac{\text{RCF}}{r \times 1,118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge.

For more on the centrifugal radius see the chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

13 Centrifugation of materials with higher density

The rotors are designed to centrifuge substances up to a maximum mean homogenous density of 1.2 kg/dm³ when rotating at the stated speed.

Denser substances must be centrifuged at lower speed.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}\text{)} = \sqrt{\frac{1.2}{\text{Greater density}}} \times \text{Rated speed}$$

e.g.: RPM 4000, density 1.6 kg/dm³

$$\text{n}_{\text{red}} = \sqrt{\frac{1.2}{1.6}} \times 4000 = 3464 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

14 Haematocrit measurement

14.1 Haematocrit capillaries - Overview

Standard capillary tubes, Order No.: 2074

The classical heparinized haematocrit capillary tube made of glass. Tried and tested for years.

Mylar-reinforced capillary tubes, Order No.: 1072

The heparinized glass capillary is coated with plastic. Should such a capillary break the glass splinters are retained in this coating. This means that the risk of injury for the user is reduced.

Self-sealing capillary tubes, Order No. 1071

These capillary tubes, which are also heparinized, no longer need to be sealed. They are already equipped with a seal. At first, the seal material is air-permeable so that the capillary can be filled. The seal will only become impermeable when it comes into contact with blood. Therefore, there will be no leaks during centrifugation.

14.2 Filling of the haematocrit capillaries

The capillaries have different volumes according to type.

The following are the exact data:

Capillary:	Standard	Mylar-reinforced	Self-sealing
Order No.:	2074	1072	1071
Total volume without seal (μl):	49.6	74.9	74.9
Internal diameter (mm):	0.92	1.12	1.12
Overall length (mm)	75	75	75

Therefore, more blood is needed for the capillaries No. 1072 and 1071 to be filled. The dimensions, however, are still in compliance with the requirements of DIN 12846 (standard for haematocrit capillaries) and the accuracy of the haematocrit values is not affected.

For the haematocrit measurement the capillaries are always filled up to $\approx 75\%$:

Capillary:	Standard	Mylar-reinforced	Self-sealing
Order No.:	2074	1072	1071
Volume when filled up to 75%(μl)	37.2	56.2	56.2

14.3 Sealing and centrifugation of the haematocrit capillaries

14.3.1 Standard capillary tubes, Order No.: 2074, and mylar-reinforced capillary tubes, Order No.: 1072

- When filling the capillary ensure that the end of the capillary opposite the filler opening remains dry.
- Seal at the dry end of the capillary:
Insert the capillary tube vertically into the seal up to the point where the end of the capillary touches the bottom of the seal plate.
Slightly tilt the capillary tube to the side and remove it from the seal.
- Place the capillary tubes with the lute side facing outwards in the haematocrit rotor and apply the rotor lid and click into place.
- Centrifuge.
- Determine the haematocrit value.

14.3.2 Self-sealing capillary tubes (Order No. 1071)

- Fill the capillary tube.
- Let the column of blood settle in the direction of the seal.

The blood must touch the seal!

- Tap the seal end of the capillary three times on the laboratory table!
This serves to intensify the contact with the seal material and to enhance the sealing process. Only in this manner the capillaries will be sealed off!
- Place the capillary tubes with the lute side facing outwards in the haematocrit rotor and apply the rotor lid and click into place.
- Centrifuge.
- Determine the haematocrit value.

14.4 Loading the rotor

- Load opposite positions equally.
- Put lid onto rotor and click into place.
- Place the lid and lightly press down the front edge of the lid.
- Turn the handle on the front panel to the right.

If the symbol is illuminated in the rotation display, the lid has been closed properly.

14.5 Determination of the centrifugation time

The following approximation formula is used for the determination of the centrifugation time:

$$\frac{100000}{RCF} = \text{time [min]}$$

Example:

Given an RCF of 16060 $\frac{100000}{16060} \approx 7 \text{ min}$ of centrifugation time.

For the calculation of the RCF see chapter "Relative Centrifugal Force (RCF)".

14.6 Reading the haematocrit value

Perform the set-up as follows:

1. hold the haematocrit rotor.
2. move the evaluation disk (lid of the haematocrit rotor) until the 100%-line of the evaluation disk is in line with the upper margin of the liquid column.
3. Read the haematocrit value at the upper margin of the erythrocyte column.

Repeat steps 1 to 3 until all capillaries have been analysed.

15 Emergency release

The lid cannot be opened in the event of a power failure. An emergency release has to be executed by hand.

For emergency release disconnect the centrifuge from the mains.
Open the lid only during rotor standstill.

See figure on page 2.

Switch off the mains switch (switch position "0").

Look through the window in the lid to be sure that the rotor has come to a standstill.

Position the centrifuge on the table sufficiently far forward so that the hole (Fig. 1, A) for the emergency release on the base of the casing can be reached.

- Insert the hexagonal pin driver into this hole, push it upwards and at the same time rotate the lid-lock handle to the left.
- Open the lid.

16 Maintenance and servicing

Pull the mains plug before cleaning.

Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Cleaning agents and disinfectants which lie in the pH range 5 – 8 are to be utilised. Alkaline cleaning agents with a pH value > 8 are to be avoided.
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

Haematocrit rotor, rotor lid, padded shells, and sealing rings are not autoclavable.

16.1 Cleaning the centrifuge

- Clean the centrifuge housing and the centrifuging chamber regularly, using soap or a mild detergent and a damp cloth if required. For one thing, this serves purposes of hygiene, and it also prevents corrosion through adhering impurities.
- In the event of condensation water formation, dry the centrifugal chamber by wiping out with an absorbent cloth.
- If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately. For surface disinfection we recommend Bacillol® manufactured by Bode Chemie in Hamburg or Biocidal ZF™ from the company WAK-Chemie Medical GmbH in Steinbach.
- Lightly rub the rubber seal of the centrifuge chamber with talcum powder or a rubber care product after each cleaning.

16.2 Cleaning and disinfecting the haematocrit rotor

If despite all precautions a leakage of haematocrit capillaries during centrifugation should occur, rotor and lid must be cleaned according to the following instructions:

1. Carefully remove padded shells or sealing ring from the haematocrit rotor.
2. Soak haematocrit rotor, lid and padded shells or sealing ring in cold water, until the blood residues are completely dissolved.
3. Perform a suitable disinfection. For disinfection we recommend Helipur® H plus N from the company B. Braun Melsungen. For the removal of radioactive material we recommend decon neutracon® from the company Decon Laboratories Limited. The application of disinfection substances on the lid is to be avoided as the inscription can dissolve.
4. Rinse the components in cold water and dry them.
5. They can now be replaced into the centrifuge.

16.3 Procedure in case of broken haematocrit capillaries

Risk of injury!


Best work with two pairs of gloves worn on top of each other. Wear safety goggles and mask, as glass splinters could be blown off!

1. Carefully remove lid of haematocrit rotor.
2. Remove larger pieces of broken capillaries using pincers.
3. Remove rotor. Carefully and slowly remove the affected padded shells or the sealing ring with pincers and exchange them.
4. Clean and disinfect as described above.

16.4 Rotor removal / installation

Removal of haematocrit rotor:

1. Open the lid and disconnect the centrifuge from the power supply.
2. Remove the rotor lid by pressing the press lock.
3. Unscrew the fixing nut with the supplied spanner.
4. Pull the haematocrit rotor off the motor shaft.

Installation of the haematocrit rotor:

1. Insert the haematocrit rotor vertically onto the motor shaft. The motor shaft dog has to fit in the rotor slot.
2. Tighten the fixing nut with the supplied spanner.

The fixing nut must always be firmly tightened to the motor shaft.

17 Faults

If the fault cannot be eliminated with the help of the fault table, please inform Customer Service.

Please state the type of centrifuge and the factory serial number. Both values are visible on the centrifuge type plate.

Perform a MAINS RESET:


- Switch off the mains switch (switch position "0").
- Wait at least 10 seconds and then switch on the mains switch again (switch position "I").

Fault	Display	Cause of fault	Remedy
No display	---	No voltage Mains input fuses defective.	<ul style="list-style-type: none"> - Check distribution voltage. - Check mains power input fuse, refer to Chapter "Change mains input fuse". - Mains switch ON.
System reset	- 2 -	Power failure during the centrifugation run. (The centrifugation run was not finished.)	<ul style="list-style-type: none"> - When stationary, open lid and press START key. - Repeat the centrifugation run if necessary.
Balance error	- 3 -	The rotor is unevenly loaded.	<ul style="list-style-type: none"> - Open lid when rotor is stationary. - Check the loading of the rotor, see chapter "Loading the rotor". - Repeat the centrifugation run.
Communication	- 4 -	Fault in control unit or power unit.	<ul style="list-style-type: none"> - Perform a MAINS RESET when the rotor has been stationary.
Overload	- 5 -	Motor or motor control defective.	<ul style="list-style-type: none"> - Perform a MAINS RESET when the rotor has been stationary.
Overvoltage Undervoltage	- 6 - - 8 -	Supply voltage outside tolerance (see Technical Data).	<ul style="list-style-type: none"> - Perform a MAINS RESET when the rotor has been stationary. - Check supply voltage.
Excess temperature	- 9 -	Excess temperature switch in motor has triggered.	<ul style="list-style-type: none"> - When rotor is stationary, open lid using emergency unlocking (see Emergency Unlocking chapter). - Allow motor to cool down.
Controller watchdog	- C -	Fault in control unit.	<ul style="list-style-type: none"> - Perform a MAINS RESET when the rotor has been stationary.
Lid error	- d -	Error in lid locking or lid closure.	<ul style="list-style-type: none"> - Perform a MAINS RESET when the rotor has been stationary.
Short circuit	- E -	Short circuit in control unit / power unit.	<ul style="list-style-type: none"> - Perform a MAINS RESET when the rotor has been stationary.
Configuration	- F -	Incorrect machine version	<ul style="list-style-type: none"> - Inform Customer Service.

18 Change mains input fuse

Switch off the mains switch and separate the centrifuge from the mains!

The fuse holder (A) with the mains input fuses is located next to the mains switch.

- Remove the connecting cable from the machine plug socket.
- Press the snap-fit (B) against the fuse holder (A) and remove.
- Exchange defective mains input fuses.

Only use fuses with the rating defined for the type. See the following table.

- Reinsert the fuse holder until the snap-fit clicks shut.
- Reconnect the centrifuge to the mains supply.

Model	Type	Fuse	Order no.
HAEMATOKRIT 210	2104	T 3,15 AH/250V	E997
HAEMATOKRIT 210	2104-01	T 5 AH/250V	E914

19 Acceptance of the centrifuges for repair

If the centrifuge is returned to the manufacturer for repair, it must be decontaminated and cleaned to protect persons, environment and material.

We reserve the right to accept contaminated centrifuges.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

20 Disposal

When you are disposing of the device, the respective statutory rules must be observed.

Pursuant to guideline 2002/96/EC (WEEE), all devices supplied after August 13, 2005 may not be disposed as part of domestic waste. The device belongs to group 8 (medical devices) and is categorized in the business-to-business field.

The icon of the crossed-out trash can shows that the device may not be disposed as part of domestic waste.

The waste disposal guidelines of the individual EC countries might vary. If necessary, contact your supplier.

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All rights reserved. No part of this publication may be reproduced without the written prior permission of the copyright owner.